



# METABOLIC METHODS

*Clinical Procedures in the Study of  
Metabolic Functions*

By

C FRANK CONSOLAZIO

Chief of Biochemistry United States Army  
Medical Nutrition Laboratory  
Chicago Illinois

ROBERT E JOHNSON, M D, D Phil (Oxford)

Professor and Head of the Department of Physiology  
University of Illinois Urbana Illinois

and

EVELYN MAREK, M A

Biochemist, United States Army  
Medical Nutrition Laboratory  
Chicago Illinois

ILLUSTRATED

ST LOUIS  
THE C V MOSBY COMPANY  
1951

COPYRIGHT 1951 BY THE C V MOSBY COMPANY

*Press of*  
*The C V Mosby Company*  
*St Louis*

THIS BOOK IS DEDICATED  
TO OUR COLLEAGUES PAST AND  
PRESENT AT THE HARVARD FATIGUE  
LABORATORY AND THE MEDICAL  
NUTRITION LABORATORY





## PREFACE

---

This book presents in detail methods which have proved to be useful in our experience during fifteen years of research on various aspects of mammalian metabolism particularly that of human beings in health and in disease. Three somewhat different types of work have posed methodological problems: first metabolic balance studies, second routine laboratory analysis of tissues, food and body fluids, and third studies out in the field. A common characteristic of all three is the accumulation of large numbers of specimens on which a wide variety of analyses must be performed. The most satisfactory methods in all three cases are those which use the least amount of samples and which can be run routinely in large numbers. Sometimes it is desirable to sacrifice some specificity in favor of simplicity, and sometimes it is necessary to sacrifice some accuracy for the sake of speed and reproducibility. We have described in detail only those methods which have been used routinely by us long enough to clarify the precautions essential for success.

We have omitted intentionally any description of isotopic, ultracentrifuge and electrophoretic techniques for the reason that as yet they are not widely used routinely in metabolic balance studies. This omission will doubtless have to be rectified in future editions.

The manner of presentation adopted in this book has evolved over many years of training technicians and graduate students. We have found that a uniform detailed presentation following a single basic outline is a great saving of time for the instructor, especially when examples of calculations and a listing of important precautions are also included.

Some of the material in this book was taken from two previous compilations:

"Syllabus of Laboratory Methods from the Harvard Fatigue Laboratory" by W. V. Consolazio, S. M. Horvath and D. B. Dill (1941) and

"Laboratory Manual of Field Methods for Biochemical Assessment of Metabolic and Nutritional Condition" by R. E. Johnson, F. Sargent, C. F. Consolazio and P. F. Robinson (1945).

We are of course indebted to many people for help and advice during the preparation of this book. Our former colleagues D. B. Dill, H. S. Belding, R. C. Darling and F. Sargent have given valuable assistance. We are particularly indebted to Mr. Carl Gordon for producing or redrawing all of our charts. Almost all members of the Medical Nutrition Laboratory helped at one time or another, especially the whole secretarial staff. In particular Mrs. I. DeCordial worked hard in preparing many drafts and Mrs. L. Vrhel

## 10 Preface

Mrs I DeCordial and Miss E Bielat were solely responsible for the final typography of the book. Members of the Statistical Branch assisted in many ways and Captain J Rielly MSC rendered great service in the assembling of the manual. The members of the Biochemistry Branch (Dr G Weigand Mr J Bell Corporal T Hutton Corporal R Powell and Private R Golden) helped out on many occasions and Private H Rosen drafted the sections on the Beckman spectrophotometer and flame photometer. The section on microbiological procedures was prepared by Captain L Crowley MSC Mrs B Enwall and Miss L Blanksma. Miss M Maloney assisted in the preparation of the sections on clinical laboratory procedures.

In a book such as this most of the work is not original. Permission has been obtained from authors and publishers for the reprinting of certain charts and graphs. For the use of illustrations of apparatus we wish to thank Dr Peter Scholander Swarthmore University Pennsylvania Dr Ancel Keys Laboratory of Physiological Hygiene University of Minnesota Mr W V Consolazio Office of Naval Research Washington D C Dr D B Dill Medical Division Edgewood Arsenal Maryland and Dr H S Belding of the Climatic Research Laboratory Lawrence Mass.

We have received willing permission from numerous manufacturers to describe their products and have been supplied by them with photographs of their instruments and with circuit diagrams. In the appropriate places in the text their contributions are acknowledged.

C F CONSOLAZIO  
R E JOHNSON  
EVELYN MAREK

Urbana Illinois and  
Chicago Illinois

# CONTENTS

|  | SECTION I   | PAGE |
|--|-------------|------|
| INTRODUCTION   | -           | 17   |
|  | SECTION II  |      |
| THE COLLECTION AND STORAGE OF SPECIMENS  | -           | 19   |
| General Considerations   | 19          |      |
| Food (Whole Diets in the Field)  | 19          |      |
| Whole Blood and Serum (In the Field)   | 20          |      |
| Urine (In the Field)   | 2           |      |
| Tissue   | 25          |      |
| Feces  | 26          |      |
|  | SECTION III |      |
| INSTRUMENTATION  | -           | 27   |
| General Principles of Spectrophotometry  | 7           |      |
| Spectrophotometer Coleman Junior Model 6   | 33          |      |
| Beckman Photoelectric Quartz Spectrophotometer (Model DU)  | 38          |      |
| Fluorometry—General Principles   | 50          |      |
| Photofluorometry—The Coleman Photo fluorometer Model 10B   | 51          |      |
| Fluorometric Methods for Qualitative Use and for Quantitative Measurements in the Field—Visual Fluorometry | 57          |      |
| Flame Photometry—The Beckman Flame Photometer Attachment   | 59          |      |
| Balances   | 65          |      |
| Volumetric Analysis  | 69          |      |
| Instrumentation Miscellaneous References   | 79          |      |
|  | SECTION IV  |      |
| BIOCHEMICAL PROCEDURES   | -           | 80   |
| A MINERALS   | -           | 82   |
| Total Water and Total Solids in Food Tissue Cells Serum and Feces  | 8           |      |
| Total Base in Serum Urine Cells and Tissue   | 83          |      |
| Sodium in Serum Urine Cells Tissues Diets and Stools (Gravimetric Method)                                  | 8           |      |
| Potassium in Serum Urine Cells Tissue Diets and Stools   | 91          |      |
| Calcium in Serum Urine Food Tissue and Feces   | 93          |      |
| Magnesium in Serum   | 97          |      |
| Inorganic Iron in Serum  | 99          |      |
| Copper in Serum  | 100         |      |
| Chloride in Serum and Urine  | 105         |      |
| Inorganic Phosphorus in Serum Urine Food and Feces   | 109         |      |
| Sulfur in Urine  | 112         |      |
| Minerals Miscellaneous References  | 115         |      |
| B PROTEIN AND NITROGEN COMPOUNDS   | -           | 120  |
| Total Nitrogen (Micro Kjeldahl) in Serum Urine Food and Feces  | 125         |      |
| Serum Proteins (Biuret Method)   | 1           |      |
| Albumin and Globulin in Serum  | 1           |      |
| Whole Blood Hemoglobin and Serum Protein (Copper Sulfate Specific Gravity)                                 | 127         |      |
| Protein in Urine   | 133         |      |
| Fibrinogen in Plasma   | 134         |      |
| Nonprotein Nitrogen in Serum and Whole Blood   | 135         |      |
| Creatinine in Urine  | 136         |      |
| Creatine in Urine  | 139         |      |
| Urea Nitrogen  | 140         |      |
| Amino Acid Nitrogen in Urine   | 14          |      |
| Uric Acid  | 143         |      |
| Uric Acid (Uricase Method)   | 145         |      |
| Paper Chromatography   | 147         |      |
| Proteins and Nitrogenous Compounds Miscellaneous References  | 151         |      |
| C CARBOHYDRATES AND DERIVATIVES  | -           | 154  |
| Blood Glucose (Colorimetric)   | 154         |      |
| Glucose (Volumetric)   | 156         |      |
| Lactic Acid  | 157         |      |
| Glycogen in Tissue   | 161         |      |
| Alcohol in Blood and Urine   | 163         |      |
| Total Bisulfite Binding Capacity   | 165         |      |
| Carbohydrates Miscellaneous References   | 166         |      |

|  | PAGE |
|--|------|
| D FATS AND DERIVATIVES   | 163  |
| Total Lipids in Serum 168 Total Cholesterol 169 Free and Combined Cholesterol in Blood 171 Lipid Phosphorus in Blood and Serum 173 Acetone Bodies 175 Fats and Derivatives Miscellaneous References 178  |      |
| E VITAMINS   | 179  |
| Vitamin A and Carotenoids (Carr Price) 179 Vitamin A and Carotenoids (Dichlorohydrin) 181 Thiamine in Urine Food Blood and Feces 184 N-Methyl nicotinamide in Urine 189 Riboflavin in Urine Blood Food and Feces 190 Ascorbic Acid in Whole Blood Serum and Urine 195 Ascorbic Acid in Tissue 199 Rapid Field Methods for Ascorbic Acid, Riboflavin Thiamine and N-Methyl nicotinamide 201 Vitamins and Nutrition Miscellaneous References 210 |      |
| F HORMONES   | 215  |
| Ketosteroids in Urine 215 Hormones Miscellaneous References 218  |      |
| G ENZYMES  | 221  |
| Lipase in Serum 221 Phosphatase Alkaline and Acid 2-2, Cholinesterase in Serum and Cells 224 Amylase in Serum 226 Enzymes Miscellaneous References 227   |      |
| H PIGMENTS   | 229  |
| Hemoglobin (Sanford Sheard) 229 Total Hemoglobin (Evelyn Malloy) 230 Oxyhemoglobin Methemoglobin Sulfhemoglobin 231 Pigments Miscellaneous References 233  |      |
| I HYDROGEN ION CONCENTRATION   | 234  |
| Indicators 234 Buffer Solutions 234 Glass Electrodes 235 Measurement of pH in Solutions 243 Hydrogen Ion Concentration Miscellaneous References 243  |      |

## SECTION V

|  |     |
|--|-----|
| MICROBIOLOGICAL PROCEDURES   | 244 |
| Pteroylglutamic Acid (Folic Acid) 244 Nicotinic Acid 247 Biotin 249 Pantothenic Acid 250 Pyridoxine (Vitamin B) 251 Riboflavin 253 Tryptophane 255 Agar Plate Assay for Biotin 256 Microbiology Miscellaneous References 259 |     |

## SECTION VI

|  |     |
|--|-----|
| PHYSIOLOGICAL MEASUREMENTS   | 261 |
| A BLOOD GASES  | 261 |
| Calculation of Gas Volumes 261 Techniques for Obtaining Arterial Blood and for Handling Blood Anaerobically 261 The Van Slyke Neil Manometric Apparatus 265 Oxygen in Whole Blood (Van Slyke) 267 Carbon Dioxide in Whole Blood and Serum (Van Slyke) 270 Oxygen and Carbon Dioxide in Single Samples of Blood (Van Slyke) 272 Oxygen Capacity (Van Slyke) 273 Carbon Monoxide in Blood (Total Heme Pigments) 275 Oxygen in Blood (Roughton Scholander) 277 Dissolved Nitrogen in Blood (Edwards Scholander Roughton) 283 Carbon Monoxide in Blood (Scholander Roughton) 285 Carbon Dioxide in Blood (Scholander Roughton) 291 Combined Carbon Dioxide and Oxygen (Scholander) 293 Manometric Methods Miscellaneous References 299 |     |

|   | PAGE |
|---|------|
| <b>B PHYSICAL CHEMISTRY OF BLOOD GASES</b>  | 301  |
| Technique for Equilibrating Arterial and Venous Blood Against Gas Mixtures of Known Composition 301 Dissociation Curves of Blood 305 Miscellaneous Calculations of the Physicochemical Properties of Blood 310  |      |
| <b>C ESTIMATION OF RESPIRATORY GASES</b>  | 310  |
| The Haldane Apparatus for Measuring Carbon Dioxide and Oxygen 31 The Scholander Micrometer Gas Analyzer 30  |      |
| <b>D RESPIRATORY METABOLISM</b>   | 317  |
| Collection of Expired Air—The Tissot Gasometer 37 Collection of Expired Air—The Douglas Bag Technique 331 Respiratory Metabolic Calculations 333 Basal Metabolic Rate 337 Metabolism During Work—The Treadmill 341  |      |
| <b>E PULMONARY MEASUREMENTS</b>   | 346  |
| Alveolar Air in Rest and Work 346 Functional Residual Air 351 Cardiac Output (Grollman's Method) 353 Pulmonary and Circulatory Measurements Miscellaneous References 356  |      |
| <b>F PLASMA VOLUME BLOOD VOLUME EXTRACELLULAR VOLUME</b>  | 358  |
| Plasma and Blood Volumes (Dye Method) 358 Total Circulating Red Cell Blood and Plasma Volumes (Carbon Monoxide) 361 True Blood Volume 363 Extracellular Fluid (Thiocyanate Space) 364 Plasma Volume Blood Volume and Extracellular Fluid Miscellaneous References 366   |      |
| <b>SECTION VII</b>  |      |
| <b>STUDIES IN THE FIELD</b>   | 367  |
| General Considerations 367 Laboratory and Clinical Methods 367 Treadmill Test of Physical Fitness for Hard Muscular Work 369 Field Tests of Physical Fitness 371 Surveys of Food Consumption 375 Heat Balance—Clothing 38 Meteorological Measurements 381 Preservation of Animal Tissues for Histological Studies 389 |      |
| <b>SECTION VIII</b>   |      |
| <b>TECHNIQUES FOR METABOLIC WARDS</b>   | 391  |
| General Considerations 39 Derived Data—Theoretical Balances 393 Derived Data—Theoretical Weight 394   |      |
| <b>SECTION IX</b>   |      |
| <b>CLINICAL LABORATORY PROCEDURES</b>   | 396  |
| <b>A MISCELLANEOUS REFERENCES TO TEXTBOOKS</b>  | 396  |
| <b>B HEMATOLOGY</b>   | 397  |
| Calculation of Mean Corpuscular Volume Mean Corpuscular Hemoglobin Content and Mean Corpuscular Hemoglobin Concentration 397 Calculation of Sedimentation Rate 399 Direct Eosinophil Count 399 Hematology Miscellaneous References 401  |      |

## 14 Contents

|   | PAGE |
|---|------|
| C TOLERANCE TESTS - - - - -   | 403  |
| Glucose Tolerance Tests 403 Vitamin Tolerance Tests 405   |      |
| D KIDNEY FUNCTION - - - - -   | 409  |
| Titratable Acidity 409 Phenolsulfonphthalein Test 410 Renal Clearance Tests<br>Miscellaneous References 411   |      |
| E LIVER FUNCTION - - - - -  | 419  |
| Bromsulfalein Test 412 Serum Bilirubin 413 Urobilinogen in Urine and Feces<br>415 Liver Function Miscellaneous References 417   |      |
| F FUNCTIONAL TESTS OF ENDOCRINE GASTRIC AND OTHER SYSTEMS - - - -   | 419  |
| Pituitary Adrenocorticotrophic Hormone (ACTH) Test 419 References to<br>Endocrine Gastric and Other Functional Tests 420  |      |
| G SUBSTANCES OF OCCASIONAL CLINICAL IMPORTANCE MISCELLANEOUS REFERENCES - -   | 421  |
| SECTION X   |      |
| STATISTICAL METHODS - - - - -   | 42   |
| References to Standard Works 422  |      |
| SECTION XI  |      |
| MISCELLANEOUS DATA - - - - -  | 424  |
| Atomic Weights 424 Conversion Between Different Systems of Weights and<br>Measures 425 Elements Listed in Order of Atomic Number, 428 Composition of<br>Body Fluids 429 |      |

## **METABOLIC METHODS**

*Clinical Procedures in the Study of Metabolic Functions*





## SECTION I

### INTRODUCTION

---

The purpose of this manual is to assist in the carrying out of metabolic studies whether in the laboratory or in the field. Perhaps a brief description of how the book is organized will facilitate its use.

A standard form of presentation has been adopted for all detailed descriptions of methods. The major headings presented in each case are (a) title of method (b) a key reference (c) a brief outline of the principle of the method (d) a detailed list of apparatus that will be required (e) a detailed list of reagents and how to prepare them (f) a step by step description of the method (g) the general equation for calculating the results (h) an example of a typical set of results calculated according to the general equation and (i) precautions which we have found to be important in using the method. Those sections that do not deal with specific methods have no fixed pattern the presentation being appropriate to the subject matter.

The first major section of the book deals with the collection and storage of specimens this being a primary concern in all metabolic studies. The second major section deals with instrumentation and includes detailed descriptions of those instruments which will be used in the various procedures presented later in the book. The third main section presents biochemical procedures covering the categories of minerals protein and nitrogen compounds carbohydrates and derivatives fats and derivatives vitamins hormones enzymes pigments and hydrogen ion concentration. The fourth section deals with microbiological procedures which are applicable to routine analysis of vitamins and amino acids. The fifth section is concerned with physiological measurements blood gases physical chemistry of blood gases estimation of respiratory gases measurement of respiratory metabolism pulmonary measurements and the estimation of plasma volume blood volume and extracellular fluid.

Starting with the sixth major section the presentation changes somewhat. The sixth section itself deals with studies in the field in which the authors have had a wide and varied experience. (Field physiology poses many interesting methodological problems which have to be met in a manner not always as rigorous as is customary in the laboratory.) The seventh section discusses briefly some of the techniques which we have found useful in metabolic wards. The last major section deals with clinical laboratory procedures emphasis being placed on functional tests. This section does not attempt to be comprehensive in fields of clinical pathology that have been covered well in numerous other textbooks.

Two minor sections finish the text of this manual. A list of references is given to standard statistical treatises. Finally a series of miscellaneous tables of useful data is appended.

Throughout the book extensive use is made of D Ocagne alignment charts. We have found that these charts have two great virtues in metabolic studies: first they speed up the computation of the voluminous data which are always collected, second they minimize errors in computation. All of the alignment charts presented herein are used in the same way: the final answer is read from the middle line when a string or ruler is stretched straight between the values listed on the right and left hand lines.

Since no two workers will always agree on the merits of individual methods we have included at the end of many sections an eclectic bibliography of other useful methods. These bibliographies are by no means complete, but they do give some suggestions on where to go for further information.

The index contains references to substances, methods, authors, tables and charts. It is felt that a single index is easier to use than several separate indices.

## SECTION II

# THE COLLECTION AND STORAGE OF SPECIMENS

### A GENERAL CONSIDERATIONS

The handling of tissue food body fluids and excreta will depend upon the analyses required. Below are listed techniques that have proved useful in metabolic and field studies when ordinarily we have required hematology and analyses for nitrogen minerals vitamins of the B complex and ascorbic acid.

### B FOOD (WHOLE DIETS IN THE FIELD)

#### Apparatus

- 1 Diet scales
- Food chopper
- 3 Waring Blendor one quart capacity
- 4 Suitable pots and pans
- 5 Two mixing cylinders of 1 liter capacity
- 6 Deep freeze unit or equivalent
- 7 Brown bottles with screw top of 100 ml capacity

#### Reagents

- 1 If ascorbic acid is to be estimated a 10% solution of oxalic acid in water is needed for stabilizing the vitamin

#### Procedure

- 1 Weigh the total day's diet including fluid
- Sampling technique is critical in all kinds of food analysis. Obtain a truly representative collection of the diet. Mix this well by hand or in a mechanical grinder.
- 3 Weigh accurately on diet scales a sample of about 500 grams
- 4 Mix this with approximately 1 liter of water (If vitamin C is to be estimated use 100 ml of 10% oxalic acid and 800 ml of water)
- 5 Homogenize the mixture in the Waring Blendor
- 6 Dilute the homogenate to exactly 1000 ml
- 7 Mix well by pouring from one cylinder to another
- 8 Fill a 50 ml brown bottle about 75 full and store in deep freeze

#### Calculation

$$\text{gm nutrient/day} = (\text{mg nutrient/ml final mixture}) \times \frac{1000}{1000} \times \frac{\text{gm total daily diet}}{\text{gm sample used}}$$

## 20 Collection and Storage of Specimens

### Example

Total day's diet (including fluid) weighed 2750 gm

Sample taken for mixing was 500 gm

Final volume of mixture was 2000 ml

The final mixture contained 20 mg N/ml

Therefore

$$\text{gm N/day} = (20) \times \frac{2000}{1000} \times \frac{2750}{500} = 220$$

### Precautions

- 1 Sampling of the diet is critical especially when whole chunks of meat are used.
- 2 Complete homogeneity of the mixture is essential. For the highest accuracy it must be possible to pipette the final mixture.
- 3 In deep freezing bottles will break if they are too full. Two thirds is about right.
- 4 Computations are made easier if the total day's diet can be mixed for analysis not just a sample of it.

## C WHOLE BLOOD AND SERUM (IN THE FIELD)

A properly trained crew using convenient apparatus and facilities can easily draw 70 specimens of blood in one hour

Apparatus (assuming 100 subjects)

*For drawing blood*

- 1 Convenient table at least six feet long and comfortable chair for subject
- 2 Kit containing 20 ml syringes 15 ml centrifuge tubes sharp No 20 needles sponges alcohol 80% by volume and bandages
- 3 Racks containing  $\frac{3}{4}$  oz vials prepared with heparin as described below
- 4 Dishpan containing clean water
- 5 Small pan containing isotonic saline solution made as described below
- 6 Autoclave or other field sterilizer

*For preparing and storing serum*

- 1 Centrifuge. Clay Adams Senior Clinical Angle Head Eight Place AC DC Centrifuges are thoroughly satisfactory
- 2 Wooden applicators
- 3 Two stainless steel 6 inch No 20 needles for transferring samples
- 4 Racks of  $\frac{3}{4}$  oz vials for storing samples
- 5 At least a dozen 10 ml syringes

### Reagents

- 1 Saline solution. About 2 liters for every 100 subjects. It is most conveniently prepared from commercial salt tablets designed for the preparation of physiological saline solution. The usual tablet weighs 10 gm and should be dissolved in 118 ml of water.
- 2 Heparin. Potent preparations of the Connaught Laboratory are the most satisfactory. There should be solutions of two different strengths. The first is used when it is desired to keep 10 or more ml of blood from clotting. The heparin is dissolved in hot physiological saline solution in such a strength that 1 drop will hold 10 ml of blood for 24 hours. The second or weak strength solution is designed to be dried in  $\frac{3}{4}$  oz vials to hold 2 ml of whole blood for 24 hours. It is prepared in plain water in such a strength that 1 drop will hold 2 ml of blood. When 100 subjects are being examined  $\frac{3}{4}$  oz vials are placed in their racks one or two days before the test and one drop of weak heparin is placed in each vial and is allowed to dry.

The tops are screwed on and adhesive tape labels are attached. The solutions of heparin once made are stable in cool places for at least 2 months. When dried in vials the heparin keeps indefinitely. As an example of making up heparin solutions a recent batch of Connaught Laboratory heparin was shipped in vials each containing 100 mg. For strong heparin 100 mg were dissolved in 14 ml of physiological saline and for weak heparin 100 mg were dissolved in 70 ml of distilled water. It is obvious that these amounts vary from batch to batch.

### 3 Alcohol 80% for sterilizing

## Procedure

### 1 Drawing blood

- a. On the night before needles are sharpened and sterilized syringes are checked for chips and vials are labelled.
- b. The six foot table is arranged so that the subject sits on a chair at one end with his right arm straightened over the table conveniently for one man to draw blood from an antecubital vein. Racks containing centrifuge tubes and heparinized vials are placed near the subject's arm. Dishes containing sponges in alcohol are placed conveniently within reach of the man drawing blood. A roll of 50 needles that have been carefully wrapped in gauze and sterilized is opened with just the first needle showing. As needles are needed the roll is gradually opened so as not to contaminate all of them. At the other end of the table are placed a dishpan containing clean water and the pan containing physiological saline solution. The space in the middle of the table is reserved for syringes.
- c. A crew of 4 men engages in these activities. The subjects are commanded to line up single file. Each man comes up in turn sits in the chair and straightens out his arm over the table. One man behind the chair applies venous stasis by hand and another man commences to draw blood. A third man on the side opposite the chair labels a centrifuge tube and a vial with the subject's identification. When about 15 ml of blood have been drawn the needle is removed from the arm, a sponge moistened with alcohol is placed in the antecubital space and the subject is released with instructions to maintain firm pressure by hand for 5 minutes. In our experience the mere flexing of the elbow in some instances is not enough to prevent hematomas.
- d. The man who has labeled the vials and tubes holds the vial in such a position that the operator can easily introduce about 2 ml of blood into the heparinized vial and then screws on the top and mixes the blood and heparin thoroughly by inversion. DO NOT SHAKE. At the same time the operator introduces the rest of the blood into the 15 ml centrifuge tube down the side gently in order to prevent hemolysis.
- e. The fourth member of the crew washes dirty syringes and needles by receiving them from the operator. He is stationed at the far end of the table. A sterile needle is fitted on a clean syringe and the operator begins to draw blood from a new subject who is being prepared by the second member of the crew.
- f. After the syringe and needle are washed in clean water the needle is placed in the basin of alcohol and the syringe is prepared for use by 2 washings of about 5 ml each with physiological saline solution. The residual saline solution in the syringe is in large part removed by vigorous pumping of the barrel. The amount of saline now remaining has no significant effect on 15 ml of blood. The clean syringe is placed on the table within convenient reach of the operator.

### 2 Preparation of serum

- a. The centrifuge tube containing blood is allowed to stand for at least 15 minutes so that clotting is well advanced. It may stand up to 2 hours without harm.

## 22 *Collection and Storage of Specimens*

- b A clean wooden applicator is inserted carefully down the side of the tube to the bottom. The applicator is moved around the whole circumference of the tube touching the side for the full length of the tube at all times. This step breaks strands of fibrin and is essential for obtaining a good yield of serum.
- c The centrifuge tube is placed in the centrifuge and is centrifuged at high-speed for about 10 minutes.
- d Just before step c vials are labeled and numbered to correspond with the samples in the centrifuge. It has proved unsatisfactory to label many vials at once owing to mix ups in transferring samples of serum when this is done.
- e With the centrifuge tube held in the left hand and a 10 ml syringe fitted with a long needle in the right, serum is removed and is placed in a vial.
- f Samples of whole blood and serum in carrying cases should be stored in a convenient cool place.

### Precautions

#### 1 *In drawing blood*

- a The man washing syringes must be sure that needles do not become plugged that chipped syringes are removed from circulation and especially that there is no excess saline solution left in syringes that he has cleaned.
- b Before beginning to draw blood it is advisable to brief the whole group of subjects on the desirability of maintaining firm pressure after blood has been drawn.
- c Be sure that heparinized tubes are completely dry before use.

#### 2 *In preparing serum*

- a Hemolysis is best prevented by making sure that centrifuge tubes are dry, that the needle is removed from the syringe before blood is introduced into the centrifuge tube and that blood is introduced gently down the side.
- b Make quite sure that there is adequate padding in the bottoms of the brass centrifuge cups and all centrifuge tubes to be used fit loosely in the cups without their lips touching the brass. If this is not done there will be annoying breakage.

## D URINE (IN THE FIELD)

When a hundred or more subjects are to be studied simultaneously, the work is facilitated by the use of disposable vessels, convenient storage bottles and oxalic acid as a stabilizing agent.

### Apparatus

- 1 About 120 unwaxed cardboard drinking cups. (These must be tested for substances interfering with subsequent estimations. Dixie cups No 2108 8 oz are good.)
- 2 Measuring cylinders of 250 ml capacity.
- 3 Carrying racks containing 1 oz round short wide mouth amber bottles with plastic screw tops. Each new batch of tops must be tested for substances interfering with estimations.
- 4 Porcelain spatula for delivering approximately 100 mg of oxalic acid into storage bottles.

### Reagents

- 1 Oxalic acid CP powdered

### Procedure

- 1 On the preceding night add to each storage bottle approximately 100 mg of oxalic acid by means of a calibrated small spatula. The functions of the oxalic acid are (a) to stabilize ascorbic acid and (b) to bring the acidity of the

specimen to pH 3 to 5 which is the range of maximum stability of thiamine riboside and 1-methylnicotinamide

The night before labels of adhesive tape are placed on the tops of the storage bottles and the necks are lined up

- 3 It is usual in our field surveys to collect fasting specimens in the morning. The subjects are aroused at 4:30 A.M. and are lined up at a latrine just about 5:00 A.M. At command they empty their bladders into the latrine. They are then ordered to drink at least one half pint of water in order to insure diuresis. They are also ordered to urinate in the next hour and a half only at command. Light activity is allowed and we customarily draw blood during this time as described in Calove.
- 4 Just before 6:30 A.M. each subject is issued a paper cup with his identification on it. At the command all subjects urinate into their cups and are instructed to empty their bladders completely using more than one cup if necessary. After urination the subjects place their cups in order in a convenient place for subsequent measuring.
- 5 The time at the beginning and end of collection is recorded. Subjects unable to urinate at command must be detained until they can urinate and the times for such men have to be noted particularly carefully.
- 6 It is convenient to have a crew of 2 men for measuring and storing the urine. The first man measures the volume of each specimen, reports it to the recorder and pours an aliquot of about 25 ml into the brown urine bottle, pouring excess urine into a waste jar and throwing paper cups into a garbage can. The second man records the volumes in a notebook. The third member of the crew is responsible for putting identification numbers on the 1 oz bottles to correspond with the sample from the cup serving on tops and putting the specimens in order in a rack.
- 7 The racks of specimens are placed in order in their carrying boxes and are stored in any convenient cool place.

#### Calculations

- 1 It is desirable to record urinary excretions of vitamins in terms of milligrams or micrograms excreted per hour. One column in the notebook should therefore show ml of urine excreted per hour calculated as follows:

$$\text{ml urine per hour} = \frac{\text{ml urine collected in total period}}{\text{time in hrs. and decimal fraction of hr.}}$$

#### Example

Bladder emptied 5:05 and urine collected until 6:06

ml collected 110

Therefore

$$\text{ml urine/hr} = \frac{110}{1\frac{11}{60}} = \frac{110}{1.18} \approx 91$$

The calculations are facilitated by the line chart (Fig. 1)

#### Precautions

- 1 Always have spare cups on hand in case a few men have a large diuresis.
- 2 Never put oxalic acid in the paper cups because it quickly makes the bottoms leak.
- 3 Pouring from the cups is facilitated by pinching the cup and making it into a spout.
- 4 In our experience trying to put tops on the cups is more bother than it is worth.
- 5 If calcium is to be measured the urine samples must be preserved with a layer of toluene. Oxalic acid will precipitate out calcium in the urine.



# VOLUME OF URINE PER HOUR FROM TOTAL VOLUME AND TIME OF COLLECTION




Fig 2

## E TISSUE

### Apparatus

- 1 Torsion balance (Roller Smith 25 gm capacity)
- 2 Mortar and pestle
- 3 Deep freeze unit
- 4 A variety of pipettes
- 5 A number of 15 ml round bottom centrifuge tubes
- 6 A number of one ounce brown bottles
- 7 A number of tissue grinders (Potter & R and Elvehjem C A. A Modified Method for the Study of Tissue Oxidation J Biol Chem 114 495 504 June 1936)
- (The apparatus consists of a Pyrex tube of 30 ml capacity graduated to 25 ml and fitted with a revolving glass pestle that has a 15 cm glass shaft on its upper end. The shaft fits by means of a rubber stopper to the shaft of a 1/50 hp electric motor that has a built in reducing gear of 8 1/2 ratio and a suitable rheostat.)
- 8 Centrifuge

### Reagents

- 1 If ascorbic acid is to be estimated a 5% solution of metaphosphoric acid in water or a 10% stannous chloride solution made up in 5% metaphosphoric acid is needed for stabilizing the vitamin.
- 2 Carbon dioxide solid (when unstable compounds are involved)
- 3 Trichloroacetic acid a 4% solution in water is the most generally applicable deproteinizing agent

### Procedure

- 1 For specific analyses refer to the full details of the method as presented in the appropriate sections
- 2 As a general tool for preparing solutions for routine analysis of water soluble compounds such as lactic acid the tissue homogenizer has been of considerable use
  - a Weigh tissue on torsion balance and drop into the bottom of the homogenizer tube (If an unstable system such as the glycolytic system of muscle is involved the tube should be in dry ice for quick freezing)
  - b Add approximately 10 ml of the appropriate deproteinizing solution. (For unstable systems this should be near the freezing point)
  - c Insert the rotor of the homogenizer and grind the tissue until no discrete particles can be seen in the tube
  - d Remove the rotor washing it down into the tube with a stream of deproteinizing solution delivered from a pipette
  - e Dilute the contents of the tube exactly to the desired volume with deproteinizing fluid.
  - f Centrifuge 10 minutes at moderate speed
  - g Use aliquots of the supernatant fluid for analysis
- 3 Under many circumstances samples of tissue can be weighed placed in a closed container which will prevent evaporation and stored in the deep freeze until required for analysis.
- 4 For some estimations the tissue homogenizer is unsuitable and a mortar and pestle has to be used. Skin and bones cannot be ground in the tissue homogenizer

### Calculation

See specific methods

### Precautions

See specific methods.

**F FECES****Apparatus**

- 1 A deep freeze unit
- 2 A refrigerator
- 3 A Waring Blender one quart capacity
- 4 A 2000 ml graduated cylinder
- 5 Large open mouth jars of 3 to 4 liter capacity
- 6 A number of 2 ml and 5 ml calibrated 1/2 inch BD Lok Syringes (These are recommended because of their wide bore)
- 7 A number of No 14 needles to fit the above syringes (These needles should be sawed off at one inch)

**Reagents**

- 1 Carmine in gelatin capsules containing 600 mg each
- 2 Caprylic alcohol C P

**Procedure—Preparation of Subject**

- 1 At the time to start the collection period usually of 3 days in this laboratory, give the subject one carmine pill by mouth. Seventy two hours later give him another carmine pill by mouth.
- 2 He should collect his feces as described below in a large jar until carmine first appears in a stool. Then he should collect all stools thereafter in a separate wide mouth jar until the second dose of carmine first appears. This marks the end of the 3 day collection.

**Procedure—Preparation of Specimen**

- 1 The sample is collected into a large jar with a wide mouth. Up to three days the sample is always kept in a refrigerator.
- 2 At the end of that time approximately 1500 ml of water is added and the lid of the jar stoppered tightly.
- 3 Add 1 ml of caprylic alcohol to prevent foam and shake vigorously.
- 4 Pour portions gradually into a Waring Blender and homogenize completely.
- 5 Add all the homogenized portions to a 2000 ml graduated cylinder and dilute exactly to 2000 ml.
- 6 Pour all into a large bottle or flask mix by shaking and pour aliquots into two 250 ml brown bottles. Stopper well.
- 7 Deep freeze until used for analyses.
- 8 For individual analyses refer to the appropriate section.
- 9 For almost all kinds of analyses aliquots may be measured accurately from a 1/2 inch BD Lok Syringe which delivers accurately to 0.1 ml if properly handled.

**Calculations**

$$\text{If the collection period was three days } \left. \begin{array}{l} \text{gm substance per day} = \\ (\text{gm per ml final mixture}) \times \frac{2000}{3} \end{array} \right\}$$

**Precautions**

- 1 Fill bottles only 2/3 full as in freezing the expansion will cause the bottles to crack if full.
- 2 Step 5 is critical. Unless the specimen is homogenized well enough to be pipetted accurately, results will be unreliable. No discrete particles must remain in the final mixture.
- 3 Some subjects produce feces the color of which masks carmine. In such cases charcoal or malachite green should be substituted for carmine.

## SECTION III

# INSTRUMENTATION

---

### A GENERAL PRINCIPLES OF SPECTROPHOTOMETRY

#### 1 Transmission

The ratio of the intensity of the exit beam to that of the incident beam is called the transmission of the tube

$$\text{Transmission} = \frac{B}{A} \quad (1)$$

Transmission is an expression of the total effect of the intercepting transparent body on the light beam including the reflection losses as well as the actual absorption of the body itself. Thus transmission is a proper factor for expressing the performance of a light filter since it is a measure of its total performance. It is not a term ordinarily used in spectrochemistry because the companion expression transmittance is much more convenient. Both transmission and transmittance are mentioned here in order that they may be recognized as different measurements and not regarded as different terms for the same measurement.

#### 2 Transmittance

Suppose now that a test tube of water is further modified by the addition of a portion  $C$  of a soluble colorant. The intensity of the exit beam will of course fall from the original value to some lower value. The ratio of the intensity of this exit beam with  $C$  added to that with  $C'$  absent is called the transmittance ( $T_c$ ) of  $C$ .

$$T_c = \frac{B}{B'} \quad ( )$$

It is very important to notice that  $T_c$  is an expression dependent only on the nature and amount of the added constituent independent of reflection losses and of absorption losses from all other constituents of the solution.

In spectrochemical procedures this independence of measurements is achieved by the use of optically identical cuvettes. One of these is filled with a solution that is either free from or else contains a known proportion of the considered constituent. This cuvette of solution is called the reference. Another of the identical cuvettes is filled with the solution under test containing an unknown proportion of constituent but otherwise optically identical with the reference. This latter is termed the 'sample'. The spectrophotometer is used to determine the difference in concentration  $C'$  between the

reference and sample by measuring the relative intensities of the exit beam first with the reference and then with the sample in position

It should now be clearly understood that we are always dealing with the transmittance of the difference in concentration of a constituent between reference and sample and not with absolute concentration. It will be evident also that transmittance measurements are based on the premise of optical identity and this fact must be given due consideration in the selection and handling of cuvettes and in deciding on the proper reference solution

### 3 Selecting the Proper Reference Solution

In clinical chemistry the reference solution will usually consist of a blank identical in composition with the sample except for the concentration of the considered constituent. In some instances distilled water may prove an acceptable blank but its use should be the exception rather than the rule because the possibility of introducing light absorbing constituents with reagents is usually very real and the slight gain in convenience seems hardly worth the cost. It is certain that a reference consisting of all ingredients possibly including a known proportion of the considered constituents is by far the safest selection and as explained above it is perfectly permissible to regard the considered constituent  $C$  as the difference in concentration between reference and sample. The practice of employing an air reference is to be condemned.

In the development of new procedures as well as in the routine applications of existing methods too much stress cannot be placed on the very obvious but frequently overlooked fact that no spectrochemical analysis can be more reliable than its reference.

### 4 Color and Wave Length

Referring again to a cuvette of colorant the transmittance ( $T$ ) of this solution will be less than that of pure water because it allows less light to pass. But let us assume further that the colorant is hemoglobin. Why does the solution appear red when viewed by white light?

White is actually a mixture of light rays of many different colors and only when these colors are combined in the proper proportions does the light appear white. Thus when this cuvette of solution is viewed against white light it appears red because the hemoglobin has the peculiar ability to stop the blue and green light rays much more effectively than the red. The exit beam thus contains a much greater proportion of red light and so appears colored red. This property of a substance to interfere unequally with the passage of different colored light is termed selective absorption.

It will now be apparent that because of selective absorption the  $T$  of a colored solution will depend on the color of the light used for the measurement and that the presence of a constituent will be most evident when the  $T$  measurement is made with light of that color most strongly absorbed by the constituent. This then is the outstanding advantage of the spectrophotometer as compared to a colorimeter that color most strongly absorbed by the con

stituent being measured can be readily determined and the sensitivity and reliability of the analysis greatly improved by then operating at this wave length the relative effect of interfering constituents decreasing accordingly

With a filter photometer or colorimeter it is customary to express the color of the light used for the  $T$  measurement by a name such as 'blue green' or by an arbitrary number designating the type of filter used to produce the colored light. But such terms are only approximate color expressions and for absorption spectrochemistry it is both more convenient and more exact to express the color of the light in terms of wave length ( $\lambda$ ). This numerical expression of color is particularly convenient because it makes possible a graphical expression of the relationship between  $\lambda S$  and  $T$  and this spectral transmittance ( $S-T$ ) curve strikingly indicates that particular wave length where the considered constituent absorbs most strongly and can most readily be measured.

The wave length scale of most spectrophotometers is calibrated in millimicrons (1 millimicron = 10 Angstrom units =  $10^{-6}$  millimeter) and the wave length of the measuring light is controlled by adjusting the  $\lambda$  knob of the instrument until the required wave length appears at the index of the  $\lambda$  dial. The beginner should acquaint himself with the relationship between color and wave length by placing a white paper in the instrument's cuvette well and then noting the color of the light projected on this paper with various settings of the  $\lambda$  dial.

## 5 Preparation of Spectral Transmittance Curves

The development of the spectral transmittance ( $S-T$ ) curve of the considered ingredient should be one of the first steps in the rational development of any new spectrochemical method. Obviously all such measurements should be made in an environment identical with that contemplated for the analysis and it need hardly be added that the final method must incorporate control measures for pH or other environmental factors that can affect the color characteristics of the considered constituent or of any other constituents of the solution on which the spectrophotometric measurements will depend. Usually the initial study will contemplate also the  $S-T$  curves of possible interfering constituents and with such a group of curves before him the analyst can then intelligently select that wave length where the proposed method will be most sensitive to the concentration of the constituent being measured and least sensitive to interference.

The  $S-T$  curve is also of great importance in deciding the optimum wave length because if the photometric measurements are made with monochromatic light having the same wave length as the characteristic absorption band of the considered constituent the data will almost invariably obey the Lambert Beer Law with precision. Indeed it is only when measurements are made with such monochromatic light that this law may be expected to hold exactly and while exact observance is not imperative to exact analysis it certainly is both convenient and reassuring. With properly developed and applied  $S-T$  curves it will be found possible in almost every instance to select

a proper wave length and so regulate the chemistry of the analytical method that the Lambert Beer Law applies well within clinical precision.

There do occur instances when it may prove desirable to work at some wave length other than that at which the considered constituent has maximum absorption and here again the S T curve is a chart from which one may reach logical conclusions concerning the optimum wave length. An example of this sort is the spectrochemical measurements of pH using the dye methyl orange. If the beginner is not already familiar with the S T curves of the acid and basic forms of this indicator it is suggested that he prepare them using the same concentration in each instance. Then note that the wave length of maximum pH sensitivity does not coincide with either point of maximum absorption. Because the preparation of S T curves is so simple and their significance so fundamental and applicable every operator should familiarize himself with their use by examining several typical colored solutions preparing their S T curves and then deciding the wave length best suited for the spectrochemical measurement of the colored constituent in each instance.

The S-T curve is prepared by measuring the transmittance of a single sample at a series of wave lengths within the range of the spectrophotometer, using a suitable reference. The data are then plotted on cross section paper conveniently with transmittance (T) as ordinate and expressed as percent and with wave length ( $\lambda$ ) as abscissa and expressed as millimicrons ( $m\mu$ ). Usually the points first plotted will disclose regions of low T values and additional points will then be located in these regions until the minima are established to the desired exactitude. Most spectrophotometers are capable in most instances of indicating a minimum with an accuracy of 1 or 2  $m\mu$ , simply by taking a series of closely spaced T readings in the region of the critical value.

The principal purpose of the S T curve is to indicate that wave length at which the considered constituent has minimum transmittance and the exactness of this indication is highly dependent on the amount of the constituent that the measuring beam is caused to penetrate. If the sample solution is too dilute or the cuvette too small the amount of constituent intercepting the light beam may be insufficient to produce an exactly measurable response. Conversely if the sample is too concentrated or the cuvette too wide the S-T curve may be so flattened at the region of minimum transmittance that the optimum wave length cannot be exactly ascertained. It is important to remember that the curve will be most indicative when the sample concentration and depth are so adjusted that the minimum T value will fall between 5 to 90 per cent. Thus good practice demands that the first step in the preparation of the S T curve be the selection of a suitable cuvette size and solution concentration.

This brief discussion of S T curves is presented primarily because their significance is too little recognized or at least too little utilized. Actually one of the outstanding advantages of the spectrophotometer is that it allows the development of these curves and thereby a simple and intelligent approach to the development of new spectrochemical methods and a more exact application of old ones.

## 6 The Lambert-Beer Law

The absorption of light by a substance is due to the transmission of the individual atoms or molecules from one energy state to another. The spectrophotometer is an instrument capable of measuring the intensity of light absorbed by a given substance at a given wave length making possible identification and quantitative analysis of substances in solution.

The most satisfactory operation of spectrophotometers is with systems which obey Beer's Law dealing with light absorption. The law can be stated thus: (a) all layers of a solution each of a given thickness absorb the same fraction of light incident upon them, and (b) all layers of solution each containing a given number of molecules absorb the same fraction of light incident upon them. Thus the intensity of light absorbed depends both on concentration of solute and length of light path through the solution.

If  $I_0$  is intensity of incident light  
 $dI$  is intensity of light absorbed  
 $c$  is concentration in moles per liter and  
 $dl$  is the thickness of a layer

$$\frac{dI}{I} \propto dl \quad \text{and} \quad \frac{dI}{I} \propto dc$$

$$\int_I^{I_0} \frac{dI}{I} \propto \int_0^l dl \quad \text{and} \quad \int_I^{I_0} \frac{dI}{I} \propto \int_0^c dc$$

$$\ln \frac{I_0}{I} \propto l \quad \text{and} \quad \ln \frac{I_0}{I} \propto c$$

combining

$$\ln \frac{I_0}{I} \propto lc \quad \text{or} \quad \ln \frac{I_0}{I} = klc$$

$$\text{or } \log \frac{I_0}{I} = \epsilon l \quad \text{or } D = \epsilon l \quad (\text{Equation 1})$$

where  $D$  = optical density or  $\log \frac{I_0}{I}$   
 $\epsilon$  = molar extinction coefficient and  
 $c$  = concentration in moles per liter

In this way using the last equation if  $D$  is determined experimentally for a solution of known concentration and lightpath at a given wave length the  $\epsilon$  can be calculated. Since  $\epsilon$  is independent of  $c$  and  $l$  when Beer's Law holds  $\epsilon$  for any other solution can be calculated from density readings. To test the applicability of Beer's Law optical densities of solutions of the same thickness ( $l$ ) are plotted against concentration. Since  $\epsilon l = \frac{D}{c}$ ,  $\epsilon l$  should be the slope of the graph which should be a straight line.

If it is desired to use transmittance instead of optical density Equation 1 transforms into the following expression:

$$\epsilon = K \log I/T \quad (\text{Equation 2})$$

or more simply

$$\epsilon = -K \log T \quad (\text{Equation 3})$$



In words Equation 3 states that the concentration of solute is proportional by a negative constant to the log of the transmittance

## 7 The Concentration Transmittance Calibration Graph

The importance of the Lambert Beer Law is stressed here first because it is the most simple and second because it is the only rational means for translating photometer readings to expressions of the corresponding concentration of the sample. The law also frequently allows the simultaneous determination of two or even three different constituents in the same sample without mutual interference.

The relationship between  $C$  (Concentration) and  $T$  (Transmittance) is exactly expressed by the Lambert Beer Law in practically every instance where the requirements of this law are fully met. These four requirements are

- a The  $T$  measurement must be made with monochromatic light
- b The  $T$  measurement must be made at a wave length corresponding to a region of the constituent's  $S T$  curves where  $T$  is substantially constant i.e. at the middle of a region of maximum absorption or on a 'flat' portion of the curve. The measurement should not be made on the 'side' of an absorption band except under unusual circumstances
- c The reference must be so selected that  $C = 0$  when  $T = 100\%$
- d The nature of the sample solution must be such that its  $T$  responds only to changes in  $C$ .

Equation 3 significantly states that any concentration  $C$  of a given solution is proportional to the log of its corresponding transmittance  $T$ . This means that the concentration transmittance ( $C T$ ) graph of this relationship will be a straight line if plotted on semi log coordinates and such a  $C T$  graph constitutes a calibration by which the concentration corresponding to any particular  $T$  value may be conveniently determined.

The equation further states that this specified straight line will always intersect the point ( $C = 0, T = 100\%$ ). Accordingly it is possible to prepare the  $C T$  calibration graph of an analytical method without calculations. To do this determine the  $T$  of one solution of known concentration  $C$  plot this known point on semi log coordinates and then draw a straight line intersecting this point and the point ( $C = 0, T = 100\%$ ). This line is the  $C T$  calibration graph of the method.

The precision of this method of plotting will depend on how faithfully the analytical method follows the Lambert Beer Law and it is always advisable to prepare and measure the  $T$  of several solutions of known concentration within the contemplated range of the method. When plotted these several points serve both to confirm the validity of the law and to discover erroneous or careless technique. The points should fall on a straight line.

Let it be repeated the technique of almost any determination can be so developed that the data follow the Lambert Beer Law with precision provided

that measurements are made with a true spectrophotometer. When these conditions are met the construction of the corresponding C T calibration graph is fully as simple as it sounds.

## 8 Optical Density

It is sometimes convenient in spectrochemistry to express the light stopping ability of a sample in terms of its optical density D when

$$D = 2 - \log \text{ percentage transmittance}$$

The advantage of this method of expression is that D values are directly proportional to corresponding C values. For example

The optical density of a solution (let us say copper carbamate that has been extracted from serum) is designated as D ml the serum blank as Ds and the reagent blank as Dr. Copper concentration is obtained by the following calculation

- a K is determined from the calibration curve
- b  $D \text{ ml} - (D_s + D_r) = D_{cu}$
- c  $D_{cu} \times K \times 100 = \text{Copper micrograms per 100 ml}$

# B SPECTROPHOTOMETER COLEMAN JUNIOR MODEL 6

## 1 General Principles of the Spectrophotometer

### a. Spectrophotometer

A spectrophotometer as its name implies is really two instruments in a single case a spectrometer and a photometer. A spectrometer is a device for producing colored light of any selected color and when used as a part of a spectrophotometer such a unit is usually called a monochromator and is usually calibrated to express the color of the monochromatic light it produces in terms of wave length.

A photometer is a device for measuring the intensity of light and when combined in a spectrophotometer is used to measure the intensity of the monochromatic beam produced by the associated monochromator. Usually the photometric measurement is made both with a reference and a control sample interposed in the light beam the difference between the two intensity measurements then being a measure of the transmittance of the sample at the wave length of the light.

### b The Monochromator

The monochromator of the Junior Clinical Spectrophotometer is illuminated by a 6500 exciter lamp. This lamp is powered by a 6 volt storage battery (or suitable transformer) and is turned on and off by a switch mounted integral with the fine control knob of the instrument. The filament of this lamp is a thin straight coil about 0.3 inches long and exactly positioned with respect to the prefocus base of the lamp. This base registers exactly in a special lamp socket mounted on the instrument's lamp arm so that the filament thus occupies a very exact position in the monochromator optical system. It is very important to note that exciter lamps must accordingly be very carefully handled to avoid possible bending of the prefocus base and the calibration of the instrument verified each time the exciter lamp is disturbed. Light from the exciter lamp filament enters the double lens system of the monochromator and were this an ordinary optical system would be caused to form a sharp white image at the plane of the monochromator exit slit. However there is interposed between the two condensing lenses a diffraction grating and this grating has the

peculiar ability to bend the transmitted light, the amount of the bending being proportional to the wave length. As a result of this grating the anticipated white image is actually "smeared" and there appears instead at the plane of the exit slit a true and brilliant spectrum. Most of this spectral image falls on one or the other of the slit sidewalls and is stopped; that narrow portion of the spectrum falling directly on the slit is not intercepted and so emerges alone from the opposite side of the slit as a monochromatic beam.




Fig. 7.—The Coleman Junior Spectrophotometer Model 6A (Reproduced by permission of the manufacturers.)

As indicated the position of the lamp filament may be altered by moving the lamp arm under compulsion of the wave length cam and this cam is rigidly fixed to the wave length dial. Thus, as the calibrated  $\lambda$  dial is rotated the lamp filament is swung across the optical axis of the monochromator resulting in a corresponding shift of the spectral image. The design of the instrument is such that that portion of spectral light actually passing through the slit always has exactly that wave length value currently indicated on the  $\lambda$  dial; i.e., light of the desired wave length is delivered from the exit slit of the monochromator simply by adjusting the instrument's  $\lambda$  dial to read that wave length. The position and curvature of the arc through which the filament moves is such that the waveband passing the slit is always in sharp focus.

#### c The Photometer

The photometer of this instrument consists of a barrier layer photocell connected by a suitable control circuit to a projection galvanometer illuminated by a 6500° gal



**SCHEMATIC DIAGRAM OF JUNIOR CLINICAL SPECTROPHOTOMETER  
OPTICAL AND ELECTRICAL SYSTEMS**

vanometer lamp and projecting its index image on a galvanometer scale mounted in the top of the instrument case. Mounted integral with the photocell is a correction filter selected largely to compensate for the very great differences in the photocell response to light of different wave lengths. It is to be understood that this compensation is by no means complete and that when the  $\lambda$  dial is turned the position of the index may change markedly. The method of operation of the galvanometer controls is evident from the circuit diagram (Figs 2 and 3).

It will be noted that the galvanometer of the Junior Clinical spectrophotometer is double calibrated with a linear transmittance (T) scale printed in black and with a logarithmic density (D) scale appearing in red.

#### d. The Optimum Transmittance Range

Errors in position of the galvanometer index due to instrument limitations and errors of the operator in reading the galvanometer scale together constitute the absolute instrument error, a value entirely independent of the nature of the analytical method and of the manipulative skill of the operator. The effect of this absolute error on the precision of an analysis is called the effective instrument error.

The magnitude of the effective error corresponding to a given absolute value is surprisingly dependent on the transmittance of the sample being tested and as a consequence the efficacy of an analytical procedure may well depend quite as much on the dilution ratio selected as on the verity of the chemical reactions involved. This is a fact the analyst too frequently disregards in the development of his new procedures although its observance is always well repaid by more exact and consistent data.

## 2 Setting Up and Calibration of the Instrument

### a. Mounting and Connecting

- (1) Place the instrument on a reasonably level surface free from extreme vibration mounting one round sponge rubber pad under each instrument foot to minimize any slight vibrational effects that might otherwise disturb exact galvanometer readings.
- (2) Connect the twisted black and red pair of wires to either a 6 volt storage battery or to a constant voltage transformer delivering not less than .5 watts at 6 volts. The polarity of this connection is immaterial.
- (3) If a storage battery is used always connect its positive (+ red) terminal to the positive (+ red) lead from the battery charger and connect its negative (- black or unmarked) terminal to the negative (- colorless) charger lead. Connect the line cord plug of the battery charger to the power source indicated on the charger instrument plate. Do not attempt to charge the battery during a test because variations of charger potential with line voltage changes will cause unstable readings.
- (4) Preferably locate the instrument in a position where strong light cannot fall on its panel or into the cuvette well. If this condition cannot be met then always place the 6312 light shield over the cuvette well during every reading.

### b. Adjusting the Mechanical Galvanometer Zero

The Junior is equipped with a very sensitive but very rugged taut suspension galvanometer. The zero adjustment for this element is located under the raised housing just to the left of the cuvette well. If the spectrophotometer is not disturbed or its leveling altered then this galvanometer adjustment remains very stable. The purpose of the raised housing is to protect the galvanometer adjusting lever from any accidental dislocation.

With the instrument set up and connected as specified in the preceding section first adjust the galvanometer zero as described below. Always verify the galvanometer zero adjustment before undertaking an analysis and if necessary rectify. It may con

veniently be checked at any time during an analysis and this procedure is recommended during prolonged routine examinations where extreme precision is demanded.

- (1) Turn the instrument on by rotating the fine knob clockwise until its switch clicks and the  $\lambda$  dial is illuminated.
- (2) Desensitize the galvanometer circuit by rotating the coarse knob to its extreme counterclockwise position and the fine knob as far counterclockwise as possible without turning off the switch.
- (3) Darken the photocell. If a cuvette adapter is available the photocell is conveniently darkened by placing the adapter in the well with its key 90° clockwise to the keyway of the well. In this position the adapter is not fully seated and also its window is not in the path of the light beam. An alternate method is to place a slip of opaque paper in the well so that the monochromator beam is completely stopped, covering the well with the light shield. With the photocell darkened there is no electrical input to the galvanometer which may now be caused to read exactly zero on the black scale by moving the galvanometer adjusting lever. To verify the galvanometer zero perform steps 1° and 3.
- (4) With a pencil point or other convenient tool reach under the housing and move the galvanometer adjusting lever until the galvanometer index reads exactly zero on the black  $\%T$  scale i.e. at the left hand end. Do not alter this adjustment during an analysis except to correct.

#### Verifying the Wave Length Calibration

Obviously this operation is not necessary yet confirmation of the  $\lambda$  calibration is such a quick and simple task that it may well be included in the initial setup. It must be undertaken if the exciter lamp is disturbed. A 6400 d dym um standard filter is required for this test.

The  $\lambda$  calibration is verified as follows:

- (1) Turn on the instrument and allow to warm up for 5 minutes before proceeding.
- (2) Set the wave length scale to read 610 exactly.
- (3) Insert the 6400 d dym um calibration filter in the cuvette well but with its key at right angles to the cuvette well keyway. With this angular positioning the photocell is completely darkened.
- (4) Set the galvanometer to mechanical zero as described in section d above.
- (5) Remove the 6400 calibration filter, cover the well opening and then adjust the coarse and fine knobs until the galvanometer index line coincides exactly with 100% T on the black galvanometer scale. (Tap the instrument if the index tends to stick off scale.)
- (6) Now place the previously cleaned 6400 calibration filter in cuvette well with its key entering the well keyway. Make sure that the filter is now fully seated in well.
- (7) Read the new position of the galvanometer index line on the  $\%T$  galvanometer scale. If this measured value of T checks within  $\pm 3$  scale divisions of the value of T engraved on the 6400 calibration filter the  $\lambda$  calibration is correct within  $\pm 2$  millimicrons and satisfactory.

If by any chance the  $\lambda$  calibration is found to have been shifted, in transit refer to section (d) below for full calibration details.

When the  $\lambda$  calibration is found correct the galvanometer zero properly adjusted and the instrument mounted and connected as specified it is ready for operation.

#### d. Calibration of Instrument When the $\lambda$ Calibration Is Not Correct

- (1) Set the instrument on its side, unscrew the four rubber feet and remove the sheet metal bottom to expose the lamp arm F and adjustments F' and G (see Fig. 3).
- (2) Loosen set screw F' on arm F thereby freeing the calibration screw G.

vanometer lamp and projecting its index image on a galvanometer scale mounted in the top of the instrument case. Mounted integral with the photocell is a correction filter selected largely to compensate for the very great differences in the photocell response to light of different wave lengths. It is to be understood that this compensation is by no means complete and that when the  $\lambda$  dial is turned the position of the index may change markedly. The method of operation of the galvanometer controls is evident from the circuit diagram (Figs 2 and 3).

It will be noted that the galvanometer of the Junior Clinical spectrophotometer is double calibrated with a linear transmittance (T) scale printed in black and with a logarithmic density (D) scale appearing in red.

#### **d. The Optimum Transmittance Range**

Errors in position of the galvanometer index due to instrument limitations and errors of the operator in reading the galvanometer scale together constitute the absolute instrument error, a value entirely independent of the nature of the analytical method and of the manipulative skill of the operator. The effect of this absolute error on the precision of an analysis is called the effective instrument error.

The magnitude of the effective error corresponding to a given absolute value is surprisingly dependent on the transmittance of the sample being tested and as a consequence the efficacy of an analytical procedure may well depend quite as much on the dilution ratio selected as on the verity of the chemical reactions involved. This is a fact the analyst too frequently disregards in the development of his new procedures although its observance is always well repaid by more exact and consistent data.

## **2 Setting Up and Calibration of the Instrument**

### **a. Mounting and Connecting**

- (1) Place the instrument on a reasonably level surface free from extreme vibration mounting one round sponge rubber pad under each instrument foot to minimize any slight vibrational effects that might otherwise disturb exact galvanometer readings.
- (2) Connect the twisted black and red pair of wires to either a 6 volt storage battery or to a constant voltage transformer delivering not less than 25 watts at 6 volts. The polarity of this connection is immaterial.
- (3) If a storage battery is used always connect its positive (+ red) terminal to the positive (+ red) lead from the battery charger and connect its negative (- black or unmarked) terminal to the negative (- colorless) charger lead. Connect the line cord plug of the battery charger to the power source indicated on the charger instrument plate. Do not attempt to charge the battery during a test because variations of charger potential with line voltage changes will cause unstable readings.
- (4) Preferably locate the instrument in a position where strong light cannot fall on its panel or into the cuvette well. If this condition cannot be met then always place the 6312 light shield over the cuvette well during every reading.

### **b. Adjusting the Mechanical Galvanometer Zero**

The Junior is equipped with a very sensitive but very rugged taut suspension galvanometer. The zero adjustment for this element is located under the raised housing just to the left of the cuvette well. If the spectrophotometer is not disturbed or its leveling altered then this galvanometer adjustment remains very stable. The purpose of the raised housing is to protect the galvanometer adjusting lever from any accidental dislocation.

With the instrument set up and connected as specified in the preceding section first adjust the galvanometer zero as described below. Always verify the galvanometer zero adjustment before undertaking an analysis and if necessary rectify. It may con

veniently be checked at any time during an analysis and this procedure is recommended during prolonged routine examinations where extreme precision is demanded.

- (1) Turn the instrument on by rotating the fine knob clockwise until its switch clicks and the  $\lambda$  dial is illuminated.
- (2) Desensitize the galvanometer circuit by rotating the coarse knob to its extreme counterclockwise position and the fine knob as far counterclockwise as is possible without turning off the switch.
- (3) Darken the photocell. If a cuvette adapter is available the photocell is conveniently darkened by placing the adapter in the well with its key 90° clockwise to the keyway of the well. In this position the adapter is not fully seated and also its window is not in the path of the light beam. An alternative method is to place a slip of opaque paper in the well so that the monochromator beam is completely stopped covering the well with the light shield. With the photocell darkened there is no electrical input to the galvanometer which may now be caused to read exactly zero on the black scale by moving the galvanometer adjusting lever. To verify the galvanometer zero perform steps 1 and 3.
- (4) With a pencil point or other convenient tool reach under the housing and move the galvanometer adjusting lever until the galvanometer index reads exactly zero on the black % T scale 18 at the left hand end. Do not alter this adjustment during an analysis except to correct.

#### c. Verifying the Wave Length Calibration

Ordinarily this operation is not necessary yet confirmation of the  $\lambda$  calibration is such a quick and simple task that it may well be included in the initial setup. It must be undertaken if the exciter lamp is disturbed. A 6400 didymium standard filter is required for this test.

The  $\lambda$  calibration is verified as follows:

- (1) Turn on the instrument and allow to warm up for 5 minutes before proceeding.
- (2) Set the wave length scale to read 610 exactly.
- (3) Insert the 6400 didymium calibrating filter in the cuvette well but with its key at right angles to the cuvette well keyway. With this angular positioning the photocell is completely darkened.
- (4) Set the galvanometer to mechanical zero as described in section b above.
- (5) Remove the 6400 calibrating filter, cover the well opening and then adjust the coarse and fine knobs until the galvanometer index line coincides exactly with 100% T on the black galvanometer scale (tap the instrument if the index tends to stick off scale).
- (6) Now place the previously cleaned 6400 calibrating filter in cuvette well with its key entering the well keyway. Make sure that the filter is now fully seated in well.
- (7) Read the new position of the galvanometer index line on the % T galvanometer scale. If this measured value of T checks within  $\pm 3$  scale divisions of the value of T engraved on the 6400 calibrating filter the  $\lambda$  calibration is correct within 2 mill microns and satisfactory.

If by any chance the  $\lambda$  calibration is found to have been shifted in transit refer to section (d) below for full calibrating details.

When the  $\lambda$  calibration is found correct the galvanometer zero properly adjusted and the instrument mounted and connected as specified it is ready for operation.

#### d. Calibration of Instrument When the $\lambda$ Calibration Is Not Correct

- (1) Set the instrument on its side, unscrew the four rubber feet and remove the sheet metal bottom to expose the lamp arm E and adjustments F and G (see Fig. 3).
- (2) Loosen set screw F on arm E thereby freeing the calibrating screw G.



- (3) Now if the value of  $T$  measured at step (7) is less than that engraved on the calibrating filter turn the calibrating screw "G" clockwise about  $\frac{1}{8}$  turn for each 5 percent difference in  $T$  (or vice versa) Then lock in the new position by slightly tightening on screw F
  - (4) Repeat steps c (1) to c (7) If the measured and specified values of  $T$  still fail to agree within 3% repeat steps d (1-4)
  - (5) When adjustment is satisfactory tighten screw F" firmly and then again repeat steps c (1) to c (7) as a final check
- Note When replacing the sheet metal bottom make sure that its drain hole is below the cuvette well

### 3 Instructions for Routine Use of the Instrument

- a Insert in the cuvette well a cuvette adapter of the proper size to accept the type of cuvette specified in the contemplated analytical method Standard calibrated cuvettes may be procured from the company
- b Turn on the instrument by rotating the fine knob clockwise until the switch clicks and the  $\lambda$  dial is illuminated.
- c Verify the galvanometer zero setting and readjust if necessary
- d Adjust the  $\lambda$  dial to that wave length specified in the contemplated analytical method
- e Wipe clean and then insert at the proper angle in the cuvette well a cuvette containing a sufficient volume of the reference solution (This volume will depend on the size of the cuvette)
- f Adjust the coarse and fine knobs until the galvanometer index reads that value specified in the contemplated analytical method. (Usually this will be 100%  $T$  if the black transmittance scale is used or 0 if the red density scale is utilized.)
- g Remove the reference cuvette Wipe clean and then insert at the proper angle a similar cuvette containing a portion of the sample solution
- h Read the position of the galvanometer index on the same (red or black) galvanometer scale as was used for the initial adjustment (f)
- i Refer to tables curves or scale reading for the actual concentration of the sample solution as directed in the contemplated analytical method

## C BECKMAN PHOTOELECTRIC QUARTZ SPECTROPHOTOMETER (MODEL DU)

### 1 General Construction


#### a Electrical System

The Beckman Spectrophotometer is essentially a potentiometer Light passing through the quartz prism monochromator and the solution falls on a light sensitive photo tube inducing in it a voltage proportional to the intensity of the light incident upon it This voltage disturbs the balanced potentiometer which must be rebalanced. This rebalancing is done with a variable resistance calibrated in %  $T$  and  $D$  In operation the circuit is first balanced with a blank or pure solvent in the light path then with the solution in the light path If the potentiometer reading of the solvent is arbitrarily set at zero  $D$  (100%  $T$ ) then the balanced potentiometer with the solution in the light path will read directly in %  $T$  or  $D$


#### b Optical System

The optical system of the spectrophotometer is shown schematically in Figs 4 and 5 An image of the light source A is focused by the condensing mirror B and diagonal mirror C on the entrance slit at D The entrance slit is the lower of two slits





SCHEMATIC DIAGRAM OF OPTICAL SYSTEM



SCHEMATIC WIRING DIAGRAM  
BECKMAN MODEL DU QUARTZ SPECTROPHOTOMETER

## 42 Instrumentation

windows are available with 10 mm 25 mm and 50 mm light paths and may be used for measurements on liquids over the entire wave length range of the instrument. With the exception of the demountable cells and certain gas cells all of the above cells are of fused construction and are not cemented

### g Phototubes

Two phototubes are furnished with each instrument mounted in a sealed desiccated compartment adjoining the cell compartment. A sliding rod brings either tube into position and simultaneously switches their electrical connections. A switch coupled to the phototube aperture shutter permits the dark current of the phototube to be checked at any time without removing cells or changing any control setting.

A cesium oxide phototube with high sensitivity from 1000 millimicrons to 600 millimicrons and a blue ultraviolet sensitive phototube in ultraviolet transmitting envelope with high sensitivity from 6.5 millimicrons to 400 millimicrons are supplied with the Model DU Spectrophotometer.

### h Spectral Range

The wave length scale is calibrated from 200 millimicrons in the ultraviolet region through the entire visible range and in the infrared to 2000 millimicrons. As no single phototube exhibits maximum sensitivity over the entire spectral range each Beckman instrument is equipped with two phototubes, one having maximum sensitivity in the red infrared region and the other having maximum sensitivity in the blue ultraviolet region. Either tube can be instantly switched into use.

At present the short wave length limit of operation is governed by several practical considerations such as the limiting sensitivity of phototubes, the intensity of radiation from the light source, radiation losses at reflection surfaces and the opacity of fused silica or special glasses used in the construction of absorption cells and phototubes. The upper wave length limit of operation is governed primarily by the limiting sensitivity of available photoelectric tubes. Full sensitivity performance is guaranteed for all instruments to 1000 millimicrons and satisfactory operation to 1100 millimicrons may reasonably be expected.

### i Sensitivity

The phototube and amplifier circuits of the Beckman Photoelectric Spectrophotometer are unusually stable and sensitive. A useful current amplification of  $5 \times 10^7$  permits full scale readings on a rugged milliammeter with as little as  $2.5 \times 10^{-8}$  watts of radiant energy. Because of this high sensitivity slit openings can be narrowed to secure maximum resolution. Measurements can be made over most of the spectral range with nominal band widths of less than one millimicron.

### j Accuracy

Precision optical and mechanical parts permit slit openings and wave length settings to be reproduced with great fidelity. Wave length settings are reproducible to 0.05 millimicrons in the ultraviolet and 0.5 millimicrons in the infrared. Optical parts are securely mounted in a rigid casting and are not affected by vibration, dirt or fumes.

Stray light, often an unsuspected source of error in spectrophotometric measurements, has received special attention. Instead of a replica grating the crystal quartz prism is used. Carefully designed baffles and optical surfaces of the highest quality keep diffuse reflection at a minimum. The total stray light effect over all but the extreme ends of the spectral range is less than 0.1%.

For measurements between zero and 10% transmittance the transmittance scale is expanded by merely throwing a switch so that the entire transmittance scale reads 0 to 10% instead of 0 to 100%.

## 2 Setting Up and Calibrating the Instrument

### a. Power Requirements

The power required for operation of the Beckman DU Quartz Spectrophotometer, other than that supplied by the dry batteries installed in the instrument is supplied from a 6 volt storage battery. To insure the tube stability a battery of approximately 100 ampere hours rating (i.e. a heavy duty truck battery) is recommended.

For operation securely attach the battery clamps marked - + and +6 to the corresponding voltage terminals of the storage battery. If two batteries are used in parallel it is not necessary to parallel the 2 volt terminals. This apparatus requires 6 Burgess dry cell batteries #5510 (Beckman #911) each  $7\frac{1}{2}$  volts.

### b. Selector Switch

The main instrument battery switch has four positions: OFF, CHECK, 10, 01. In the OFF position all battery circuits within the monochromator and those to the phototube compartment are disconnected. (Separate switches are provided for the tungsten lamp and hydrogen discharge lamp.) The CHECK position provides a very convenient means of adjusting the instrument for 100% transmittance on the solvent without having to set the transmittance scale to 100% because the transmission slide wire is replaced with a fixed resistor circuit calibrated to be the equivalent of 100% transmittance with the selector switch in 10 position. Measurements on the unknown solutions are usually made with the switch on 10. The transmittance scale is then read 0% to 110% transmittance and the density scale infinity to zero. The useful part of the density scale is from 0 to 1.0. Position 01 gives a tenfold increase for sensitivity in range from 0 to 11% and is used for samples having less than 11% transmittance. In this case the transmittance scale is read 0 to 11% transmittance and the density scale 10 to infinity with the readable calibration being 10 to 30 density units (read by adding 10 to the density scale reading). The 01 position also permits measurements with extremely narrow spectral band widths.

### c. Phototube Selection

The Model DU is equipped with ultraviolet sensitive and red sensitive phototubes for the range 200 to 1000 millimicrons. The proper phototubes to use is determined by the spectral range to be investigated. The red sensitive phototube required for use above 650 millimicrons (but usable to 400 millimicrons) is in position when the knob on the phototube housing is pushed IN. The ultraviolet sensitive phototube for use below 650 millimicrons is in position when the knob is pulled OUT as far as possible. The two electronic tubes of the amplifier circuit are also mounted in the phototube compartment.

### d. Light Source

The tungsten lamp is used from 330 to 1000 millimicrons and the hydrogen discharge lamp is used from 200 to 50 millimicrons. The lamp housings may be interchanged by pulling up the locking slide on the left side of the lamp housing. The power for operation of the tungsten lamp is supplied from the storage battery through the spectrophotometer circuit. A special electronic power supply unit Part 20-0000 is necessary to provide current for operation of the hydrogen discharge lamp.

### e. Filter Slide

The filter slide is located in the mounting block between the slit plate and cell compartment and is positioned by means of the knob on the front of the instrument. The front position knob pushed IN is blank and is used for measurements in the range 400 to 1000 millimicrons—also with the hydrogen tube. The second position contains a Corning No. 9463 Red Purple Correx A filter and is used in the range of 400 to 700 millimicrons with tungsten lamp. The third position is blank and may be fitted with special filter lens or screen for various applications.

## 44 Instrumentation

### f Sensitivity Control

Maximum accuracy on the percent transmittance measuring circuit is attained when the sensitivity knob is used 1 to 3 turns from its clockwise limit and requires comparatively wide slit openings. Even these slit openings give a spectral band width of less than 1 millimicron over most of the spectral range of the instrument. With the sensitivity knob in this position transmittance readings may be easily made to 0.1%. At this high sensitivity needle fluctuations may be encountered. The fluctuations may appear objectionable to the user but do not effect the accuracy to which measurements can be made. They may be reduced at some sacrifice of accuracy by operating with the sensitivity knob near its counterclockwise limit or at reduced speed of operation but no sacrifice in accuracy by installing IR /0348B Condenser Assembly. The maximum resolution of the monochromator which may be realized with the energy measuring system is obtained when the sensitivity knob is near its counter clockwise limit, or when the switch position 0.1 is used because under these conditions less light is required for balance permitting narrow slits to be used.

### g Characteristics and Use of Absorption Cells and Holder

The *Corex* cells regularly supplied with the DU Spectrophotometer are used in the wave length region from 300 to 1000 millimicrons. Cells made of silica glass must be used in the lower ultraviolet region to 220 millimicrons. Silica cells are not limited in use to the ultraviolet region but may be used through the entire spectral range of the instrument.

All absorption cells of both *Corex* and silica are filled with distilled water and tested for percent transmittance in the spectrophotometer. This transmission of *Corex* cells must be at least 75% of that of an air light path and this transmission of silica cells must be at least 70% of that of an air light path.

All standard absorption cells are matched for light transmission at various wave lengths. Matching is done in a spectrophotometer with the cells filled with distilled water. *Corex* cells are matched to within 1.0% transmittance at 800, 500, 400 and 350 millimicrons and to within 2% at 3.0 millimicrons. Silica cells are matched to within 0.5% transmittance at 350, 300, 2.0 and 240 millimicrons and to within 3% at 2.0 millimicrons. All cells approach perfect matching at longer wave lengths.

These cells are also matched with respect to the mean light path length to within 0.5% in that portion of each cell through which the beam passes. Perfect matching of the light path length is not usually necessary providing it is known. Each 10 mm cell is marked giving its light path length to the nearest thousandth of a centimeter.

The absorption cells should be thoroughly cleansed with distilled water or other suitable solvent before use. Do not habitually use hot concentrated acids which might etch the polished surfaces. Remove cells from the holder for cleaning and filling and before replacing carefully wipe the outside surfaces. For the most precise work however the outside surfaces should not be wetted or touched during a series of measurements.

It is recommended that the cell holder be used with the corner springs toward the phototube and that it be removed from the compartment for loading and unloading. It is also recommended that one position in the holder be consistently reserved for the cell containing the solvent. One of the ground glass sides of each cell should be marked and it should always be oriented the same way in its holder position.

When the holder containing filled cells is placed on the slide in the cell compartment make certain it is seated properly with knob pushed IN. Each cell is moved into the light beam by pulling the slide knob to the correct detent stop. The cell compartment cover must be in place when measurements are being made.




Fig. 1



**h Nominal Band Width**

Fig 6 gives the band width versus wave length' relationship from which the required slit openings for a given spectral band width can be determined for any wave length of interest. If it is desired to make measurements with the very minimum spectral band width, the instrument switch may be used in the 01 position and the slide wire set to 100% transmittance for balancing with the solvent or "blank" in the light beam. The sensitivity knob should be near its counterclockwise limit. In this procedure the transmission scale is read 0 to 110%—not 0 to 11%.

For measurements at a predetermined spectral band width the necessary slit width must be calculated with the consideration that optical aberrations tend to increase that slit width by approximately 0.02 mm. If  $X$  is equal to the actual slit width to be used at one wave length corresponding to the nominal band width it may be calculated as follows:

$$\begin{aligned} W_D &= \text{Spectral band width per 1 mm slit opening} \\ W_S &= \text{Spectral band width to be used for measurement} \\ X &= \frac{W_S - 0.02 W_D}{W_D} \end{aligned}$$

For measurements at wide spectral band widths the light intensity to the phototube is decreased by inserting a fine screen of proper mesh size (determined experimentally) in position 3 of the filter slide.

**i Effects of Band Widths on Transmission or Reflectance Measurements**

The apparent percent transmission at any wave length will usually depend on the resolution or effective band width and will include such influences as the slope of the transmission curve within the wave band being transmitted and the variation of the photoreceiver. It is usually easy to determine if these effects are significant by increasing or decreasing the slit width by a factor of 2 or more. It is suggested that whenever reliable values of percent transmission must be obtained for instance when intercomparing the results between instruments or when changing from the red sensitive to the ultraviolet sensitive phototube that the narrowest possible slit widths be used and the above check for the effective resolving power be made.

**j Calibration of Cells**

Because each set of cells is matched at the factory with respect to light transmission it is seldom necessary to apply individual cell transmission corrections. However when using shorter wave lengths or where extreme precision is necessary the cells should be calibrated. The calibration of each cell in a group of cells to be used should be made with respect to one cell selected as the reference cell. This calibration should be checked at several points in the wave length region to be used. The cells must be calibrated at each wave length of use when filled with solvent. Do not calibrate unfilled cells because variation in reflection from unwetted internal cell surfaces will alter standardization values.

Clean and fill each cell with distilled water or other solvent. Place the cells in the cell compartment and replace compartment cover. Set wave length dial at the wave length desired (see Sections c d e g). Turn the lamp ON, set selector switch to 10 and the transmission dial to zero density. Turn shutter switch OFF, and balance meter needle with dark current knob. Turn shutter switch ON, rotate sensitivity knob to within one turn of its clockwise limit and approximately balance the needle with the slit dial. Balance the needle accurately with the sensitivity knob and record the density reading. Slide each of the cells into the light beam and balance the needle with the transmission dial, recording the reading as density.

Select one cell as the reference cell to contain the solvent throughout all measurements. Subtract its density reading algebraically from the values for all cells. Record these resulting values as  $D$ .

Subsequent experimental density measurements  $D$  are corrected to actual sample density  $D$  by algebraically adding the correction  $D$

Example

|                         |                                     |
|-------------------------|-------------------------------------|
| $D = -0.03$             | Calibrated cell correction.         |
| $D = 1.5$               | Measured density of sample and cell |
| $D = 1.5 - 0.03 = 1.47$ | Actual sample density               |

#### k. Calculation of Solution Concentration

For quantitative absorption measurements where the unknown component is the only absorbing material at the wave length of measurement it is necessary to consider the concentration of the solution and light path length of the cell. Calculations of this sort may be carried out as follows:

Calculate  $k$  (extinction coefficient) at each wave length of measurement from measurements made on a solution of known concentration as follows:

$$D - k \cdot c \cdot d = \left( -\log \frac{T}{100} \right) \quad \begin{array}{l} c = \text{concentration of sample in solution} \\ d = \text{light path of the cell} \\ D = \text{corrected optical density} \\ T = \text{corrected percent transmittance} \end{array}$$

Using the value of  $k$  concentrations of the unknown in other solutions may be calculated from this same formula by substituting the proper density value and light path length of the cell.

#### l. Focusing of the Tungsten Lamp

When a new tungsten lamp is received in its housing from the factory only the final focusing adjustment if any should be needed. Extensive refocusing is seldom necessary but may be carried out as follows:

- (1) Remove cells and holder from compartment, unscrew the four knurled thumb screws on the end of the phototube housing and remove the housing and absorption cell compartment. Unscrew the four screws that fasten the mounting block (containing entrance mirror and filter slide) to the monochromator and remove the block, tilting the top away from the end plate in order not to drop the filter slide. Note mirror mounting and adjustment but DO NOT TOUCH ANY MIRROR SURFACES.
- (2) Turn on the current and adjust the lamp and mirror screws so that the light beam comes through the hole in the front of the lamp housing. Hold a piece of white paper in the light beam to locate the crescent shaped focal point which should appear opposite the monochromator lower (entrance) window with the center of the image about 5/16" out from the window. To raise or lower the light beam turn the upper lamp adjusting screw. If this adjustment is insufficient tilt the condensing mirror slightly.
- (3) The condensing mirror in the lamp housing may be extensively adjusted by means of the three small screws at the left side. The screw in the bushing rotates the mirror from side to side. The two clamping screws above and below bushing are loosened to tilt the mirror or to slide it toward or away from the lamp and must be tightened after adjustments are made. CAUTION: Do not release screws more than two or three revolutions as there is danger of releasing the mirror mounting.
- (4) To direct the light beam toward or away from the monochromator end plate turn the mirror adjusting screw in the bushing. To change the focal point of the image turn lower lamp adjusting screw. If this adjustment is insufficient, shift the whole mirror assembly toward the lamp to lengthen the focal distance or vice versa. After the focal point adjustments are made, replace filter slide and mounting block securely.

- (5) Set slit opening at 2 mm and wave length scale at about 800 millimicrons. Push filter slide knob in. Through the exit slit the light beam should appear as a magnified (red) image of the lamp filament. If it does not a slight adjustment of either the condensing or entrance mirror should bring out the horizontal filament lines. To adjust entrance mirror insert a screw driver into the hole in the front of the mounting block and turn the screw at the back of mirror. After a good visual adjustment has been made replace the absorption cell compartment and phototube housing.
- (6) Pull phototube knob out. Push filter slide knob in. Turn the instrument switch to 'CHECK,' the shutter switch to 'OFF' and then adjust the dark current knob to zero the needle. Rotate sensitivity knob fully clockwise. Set slits to 2 mm. Turn shutter switch ON. Rotate wave length knob toward shorter wave length (around 300 millimicrons) to zero the needle. Turn the upper lamp adjusting screw *slightly* until movement of the needle to the left indicates that more light is entering the monochromator. Rebalance the instrument by rotating wave length knob to shorter wave lengths to zero the needle. Then again adjust the lamp screw until needle is deflected to the left. Continue this procedure of alternate adjustment of screw and wave length knob until further adjustment of the screw does not cause deflection of the needle to the left. Next follow this adjusting procedure with the lower lamp adjusting screw then with the condensing mirror screw and finally with the entrance mirror screw. Then repeat the entire procedure several times for each of the four screws until the wave length scale cannot be set to any lower value. The instrument should then be in optimum focus.

#### m Focusing of Hydrogen Lamp

- (1) Attach the hydrogen lamp housing to the monochromator mounting plate.
- (2) Unscrew the four knurled thumbscrews on the end of the phototube compartment and detach the phototube and absorption cell compartments. Remove the four machine screws in the corners of the filter slide block and detach it holding the filter slide in one hand to prevent it from being damaged by falling out.
- (3) Turn on the hydrogen lamp. Hold a piece of white paper directly in front of the entrance window (lower) of the monochromator and at right angles to the light beam coming from the condensing mirror in the lamp house. Using a small screw driver turn the adjusting screw in the bushing on the left side of the lamp house until the crescent shaped image on the paper falls about 2 mm from the surface of the monochromator mounting plate. Using the upper screw on the back of the lamp house raise or lower the image until it is approximately on a level with the entrance window. Focus the beam with the lower screw on the back of the lamp house so that the crescent is narrowest opposite the entrance window as determined by moving the paper back and forth along the light beam. If the range of the adjusting screws is insufficient to make the above adjustments it may be necessary to readjust the position of the hydrogen lamp in its holder or to loosen the two mirror clamping screws on the left side of the lamp house and readjust the mirror position. CAUTION DO NOT TOUCH THE MIRROR SURFACE.
- (4) Replace the filter slide block being careful to assemble the filter slide so that its detent functions properly.
- (5) Push the filter slide in open the slits to a width of 20 mm and set the wave length scale for 450 millimicrons. Placing the eye at the hole in the filter slide block look through the exit slit into the interior of the monochromator. If the focusing adjustment is approximately correct an image of the quartz prism will be seen in the collimating mirror illuminated with an intense blue light. Again adjust the lamp house screws, also the entrance mirror adjustment screw (which

is accessible through the hole in the front of the filter slide block) until the illumination is most intense and almost completely fills the aperture of the prism.

- (6) I assemble the phototube and absorption cell compartments. I place the selector switch in the 'Check' position turn the sensitivity knob fully clockwise push in the filter slide and pull out the phototube slide. With the shutter 'Off' adjust the dark current control in the usual manner and then place the shutter in the 'On' position. With the slits set at  $\infty$  mm rotate the wave length scale toward the shorter wave lengths until the galvanometer is brought to balance. Readjust the screw on the left side of the lamp house to make the galvanometer needle go off the scale toward the left. Rebalance the galvanometer by turning the wave length scale to a shorter wave length. Repeat this process so long as the galvanometer needle can be moved to the left by adjusting the screw. Then using the other two lamp house adjustment screws continue the process until optimum focus is reached i.e. when turning any screw in either direction causes the galvanometer needle to go to the right.
- (7) When the optimum focus is reached by the procedure of Paragraph 6 the wave length scale should read less than  $\infty$ 30 millimicrons. A high wave length reading indicates (1) one or more of the aluminized mirrors in the light path has deteriorated (2) the radiation output of the hydrogen lamp is below standard or (3) the phototube sensitivity is low. The mirrors most subject to deterioration are the concave mirror in the lamp house and the small flat entrance mirror. These may be replaced or returned for resurfacing. (Inspection under visible light does not necessarily indicate reflecting quality in the ultraviolet.)

#### n. Wave Length Calibration

Because of its numerous intense spectral lines a mercury arc in an ultraviolet transmitting envelope is a useful radiation source for checking the wave length scale of the monochromator.

The hydrogen discharge lamp may be used to check the wave length scale at the hydrogen alpha line (6563 millimicrons) which is emitted with moderate intensity. Using a source of spectral lines the wave length scale may be checked as follows:

- (1) Place selector switch in '10' position and the sensitivity control in the middle of its range. Push the filter slide in and select the appropriate phototube. Adjust for dark current in the usual manner.
- (2) Set the transmission scale at about 3% and place shutter switch in the 'On' position.
- (3) Adjust the wave length scale to give maximum galvanometer deflection to the left. If necessary reduce slit width to keep the galvanometer needle on the scale.
- (4) If the wave length scale reading at maximum deflection differs from the known wave length of the spectral line by more than 0.3 millimicron, adjustment may be made as follows:

Set the wave length scale to read the correct wave length of the spectral line. Remove the snap button on the left hand end of the monochromator and turn the adjusting screw to give maximum galvanometer deflection. A 5/16" socket wrench is most convenient for turning the screw although a screw driver may be used. This adjustment has the effect of changing all readings on the wave length scale by the same angular amount. Consequently in terms of millimicrons it produces a larger effect at the long wave lengths where the scale is compressed. After making the adjustment repeat the wave length check to make certain that the best adjustment has been obtained.

#### o. Photometric Calibration

The photometric scale is checked by using a series of four filters the transmittances of which have been accurately determined by the National Bureau of Standards.

## 50 Instrumentation

These four filters are designated as carbon yellow," "cobalt blue" copper green and selenium orange. The reported transmittances of these filters are in the range from 390 m $\mu$  to 50 m $\mu$ . The temperature of the filters should be held as close to 25 °C as possible.

## 3 Operating the Instrument

- a Securely attach the battery clamps -1 +2 and +6 to the corresponding voltage terminals of the storage battery
- b Turn the lamp switch ON "
- c Set selector switch to CHECK
- d Rotate wave length knob until desired scale value is shown under hairline
- e Select proper phototube and filter if required.
- f Turn shutter switch OFF
- g Rotate dark current knob to zero the meter needle. This adjustment must be repeated occasionally (for highest accuracy between each transmission reading) Remove cell compartment cover and put filled cells in place
- h Make sure that the cells and holder are seated properly then replace the compartment cover and place the solvent or blank ' in the light path
- i Turn shutter switch ON. This switch must be turned gently otherwise the phototube housing may be jarred resulting in needle shift or unsteadiness
- j Rotate slit knob to approximately zero the needle or to give the desired slit width
- k Rotate the sensitivity knob to zero the needle accurately. For highest accuracy this adjustment should be checked between each transmission reading
- l Operate the slider knob to place an unknown sample in the light beam
- m Set selector switch to 1.0. If the transmission is less than 10% set selector switch to 0.1
- n Rotate transmission knob to zero the needle
- o Record the transmittance or density reading
- p Place the next unknown sample in the light path and zero the needle by rotating transmission knob

## D FLUOROMETRY—GENERAL PRINCIPLES

### References

- Radley J. A. and Grant J. Fluorescence Analysis in Ultraviolet Light ed. 3 London 1939 Chapman and Hall Ltd
- Dake H. C. and De Ment J. Fluorescent Light and Its Applications New York 1941 Chemical Publishing Company
- De Ment J. Fluorescent Chemicals and Their Application New York 1940 Chemical Publishing Company

### 1 Definition

As the term is used in biochemical analysis 'fluorescence' means the emission of visible light by a substance when it is irradiated with ultraviolet light. If emission of visible light occurs after the removal of the ultraviolet light the afterglow is termed 'phosphorescence'.

### 2 Laws of Fluorescence

Fluorometry is based on the general laws of photochemistry and there are three specific 'laws of fluorescence'.

- a. Energy must be absorbed by a luminescent system before emission can occur (the absorption law of De Ment)

- b In a fluorescent system the emitted wave lengths are longer than the absorbed wave lengths that is the energy emitted is less than the energy absorbed (the emission law of Stokes)
- c The absorption of radiation by a luminescent system is a quantum process involving one quantum per absorbing unit (e.g. atom or center) the yield being unity only in the ideal case (one aspect of the photochemical equivalence law of Einstein)

### 3 General Principles of Fluorometry

All quantitative fluorometry is based on a fundamentally simple system. A source of ultraviolet light casts a beam on the solution to be examined. Visible light is emitted and is measured either visually (against a comparison standard) or photoelectrically.

In general it is necessary to conduct fluorometry with dilute solutions which are measured against calibration curves prepared under conditions identical with those of the analysis. The reason for this is that the absolute intensity of fluorescence is highly variable depending upon many important factors including temperature, time of exposure, nature of instrument, traces of impurities, pH and quenching. 'Quenching' is a complex phenomenon of which there are at least three kinds: (1) interference of a molecule's fluorescence by other organic molecules, (2) interference by inorganic molecules, (3) and reabsorption of visible light by other molecules of the same compound as those being measured. Fluorometry by its nature is not as capable of extreme accuracy and reproducibility as is spectrophotometry. Every kind of analysis is strongly influenced by the nature, purity, history and condition of the material under analysis.

## E PHOTOFLUOROMETRY—THE COLEMAN PHOTOFLUOROMETER MODEL 12B

### 1 General Construction of the Photofluorometer

This instrument consists basically of: (a) a mercury vapor lamp source of ultraviolet light, (b) a series of filters to transmit a beam of ultraviolet light of the appropriate wave length, (c) a cuvette to hold the solution to be analyzed in the beam of ultraviolet light, (d) at right angles to the beam of ultraviolet light where it passes through the solution, a filter to transmit the beam of visible light in the appropriate range of wave length, (e) a photo cell in the path of the visible light, (f) an amplifying system and a direct reading galvanometer to measure the potential from the photocell. (See Fig 7 for general construction.)

### 2 Setting Up and Adjusting the Instrument

#### a. Location

Place the instrument in any convenient location and plug into 110-120 volt 60 cycle supply line, first insuring that the 'STD' knob is in the extreme counterclockwise position. Do not attempt to operate the instrument except at the frequency and voltage specified. It is wise to place rubber pads under the four feet to minimize vibration conducted through the desk.

These four filters are designated as carbon yellow, cobalt blue, copper green and selenium orange. The reported transmittances of these filters are in the range from 390 m $\mu$  to 750 m $\mu$ . The temperature of the filters should be held as close to 25°C as possible.

### 3 Operating the Instrument

- a. Securely attach the battery clamps -1, +2 and +6 to the corresponding voltage terminals of the storage battery.
- b. Turn the lamp switch ON.
- c. Set selector switch to CHECK.
- d. Rotate wave length knob until desired scale value is shown under hairline.
- e. Select proper phototube and filter if required.
- f. Turn shutter switch OFF.
- g. Rotate dark current knob to zero the meter needle. This adjustment must be repeated occasionally (for highest accuracy between each transmission reading). Remove cell compartment cover and put filled cells in place.
- h. Make sure that the cells and holder are seated properly, then replace the compartment cover and place the solvent or blank in the light path.
- i. Turn shutter switch ON. This switch must be turned gently, otherwise the phototube housing may be jarred, resulting in needle shift or unsteadiness.
- j. Rotate slit knob to approximately zero the needle, or to give the desired slit width.
- k. Rotate the sensitivity knob to zero the needle accurately. For highest accuracy this adjustment should be checked between each transmission reading.
- l. Operate the slider knob to place an unknown sample in the light beam.
- m. Set selector switch to 1.0. If the transmission is less than 10%, set selector switch to 0.1.
- n. Rotate transmission knob to zero the needle.
- o. Record the transmittance or density reading.
- p. Place the next unknown sample in the light path and zero the needle by rotating transmission knob.

## D FLUOROMETRY—GENERAL PRINCIPLES

### References

- Radley J. A. and Grant J. *Fluorescence Analysis in Ultraviolet Light* ed. 3. London 1939 Chapman and Hall Ltd.
- Dake H. C. and De Ment J. *Fluorescent Light and Its Applications* New York 1941 Chemical Publishing Company.
- De Ment J. *Fluorescent Chemicals and Their Application* New York 1941 Chemical Publishing Company.

### 1 Definition

As the term is used in biochemical analysis, "fluorescence" means the emission of visible light by a substance when it is irradiated with ultraviolet light. If emission of visible light occurs after the removal of the ultraviolet light, the afterglow is termed "phosphorescence".

### 2 Laws of Fluorescence

Fluorometry is based on the general laws of photochemistry and there are three specific 'laws of fluorescence'.

- a. Energy must be absorbed by a luminescent system before emission can occur (the absorption law of De Ment).

- b In a fluorescent system the emitted wave lengths are longer than the absorbed wave lengths that is the energy emitted is less than the energy absorbed (the emission law of Stokes)
- c The absorption of radiation by a luminescent system is a quantum process involving one quantum per absorbing unit (e.g. atom or center) the yield being unity only in the ideal case (one aspect of the photochemical equivalence law of Einstein)

### 3 General Principles of Fluorometry

All quantitative fluorometry is based on a fundamentally simple system. A source of ultraviolet light casts a beam on the solution to be examined. Visible light is emitted and is measured either visually (against a comparison standard) or photoelectrically.

In general, it is necessary to conduct fluorometry with dilute solutions which are measured against calibration curves prepared under conditions identical with those of the analysis. The reason for this is that the absolute intensity of fluorescence is highly variable depending upon many important factors including temperature, time of exposure, nature of instrument, traces of impurities, pH and 'quenching'. Quenching is a complex phenomenon of which there are at least three kinds: (1) interference of a molecule's fluorescence by other organic molecules, (2) interference by inorganic molecules, (3) and reabsorption of visible light by other molecules of the same compound as those being measured. Fluorometry by its nature is not as capable of extreme accuracy and reproducibility as is spectrophotometry. Every kind of analysis is strongly influenced by the nature, purity, history and condition of the material under analysis.

## E PHOTOFLUOROMETRY—THE COLEMAN PHOTOFLUOROMETER MODEL 12B

### 1 General Construction of the Photofluorometer

This instrument consists basically of: (a) a mercury vapor lamp source of ultraviolet light, (b) a series of filters to transmit a beam of ultraviolet light of the appropriate wave length, (c) a cuvette to hold the solution to be analyzed in the beam of ultraviolet light, (d) at right angles to the beam of ultraviolet light where it passes through the solution, a filter to transmit the beam of visible light in the appropriate range of wave length, (e) a photocell in the path of the visible light, (f) an amplifying system and a direct reading galvanometer to measure the potential from the photocell. (See Fig 7 for general construction.)

### 2 Setting Up and Adjusting the Instrument

#### a. Location

Place the instrument in any convenient location and plug into 110-120 volt 60 cycle supply line, first insuring that the ETD knob is in the extreme counterclockwise position. Do not attempt to operate the instrument except at the frequency and voltage specified. It is wise to place rubber pads under the four feet to minimize vibration conducted through the deck.



## 52 Instrumentation

### b Preliminary Adjustments

Before making tests with this instrument it is essential that both its mercury vapor lamp and electronic amplifier be turned on and allowed to reach equilibrium. Proceed as follows:

- (1) Plug extension cord into a 110 120 volt AC supply line




Fig 7—The Coleman Photofluorometer Model 1-B (Reproduced by permission of the manufacturers)

- (2) Insert an oblong panel opening #1 the proper filter for the measurement involved e.g. use filter B for thiamine or use filter B<sub>2</sub> when considering riboflavin. These two particular filters pass the 365 m  $\mu$ m crown and 4.6 m  $\mu$ m crown mercury bands respectively. The filtered beams are then projected to focus within the  $\frac{3}{4}$  round cuvette. Filter B is also fitted with a "90°" transmission screen to reduce the light intensity to a level at which significant bleaching of theochrome does not occur.
- (3) In oblong panel opening #2 insert the corresponding filter using PC-1 when measuring thiamine and filter PC-2 when examining riboflavin. The purpose of

these second filters is to screen the phototube from the waveband passed by the corresponding first filter while allowing the fluorescent light of the sample to enter the phototube and be measured.

- (4) Energize the instrument by rotating the STD knob clockwise until it just clicks. This turns on the amplifier and also starts the cooling fan and excites the mercury vapor lamp.
- (5) After 9 seconds continue rotation of the STD knob to extreme clockwise position. This adjusts the amplifier to maximum sensitivity. (The STD knob thus serves both to turn the instrument on and off as well as to adjust its sensitivity.)
- (6) Depress button B and while it is in the depressed position adjust knob BAL until the instrument meter reads approximately zero. The purpose of this adjustment is to condition the electronic amplifier to the condition of the A batteries. The adjustment is not at all critical but should be checked occasionally during prolonged tests, always again checking BLK as in (7).
- (7) Adjust the BLK knob until the instrument meter reads zero. This knob serves initially to balance the amplifier to zero as indicated above and may subsequently be used automatically to correct for fluorescence of blanks during actual analysis of fluorescent samples.
- (8) Now allow at least ten minutes for the instrument to stabilize and then repeat steps (6) and (7). (Note that the mercury vapor lamp once turned off has to cool down before it will light again.)

### 3 Making Routine Measurements

For every run of unknowns certain steps are necessary. These are (a) setting sensitivity to a fixed value against a stable reference standard (b) reading a reagent blank and (c) reading the unknowns. The following example is for thiamine by the thiochrome method. Other substances may require other stable reference standards with an emission similar to that of the unknown.

- a. Complete adjustments (1) through (8) in section 2 b above.
- b. Be sure that enough standard cuvettes are available to carry out all readings expeditiously. (This fluorometer accepts ordinary 34 mm round test tubes and is so designed optically that the tube diameter is not critical. Use not less than 8 ml of sample and be sure that tubes are covered to exclude light during all tests.)
- c. Place a cuvette of stable reference solution in cuvette well, cover with cap and then open the lamp shutter by fully depressing button P. The beam of filtered UV light projected into the cuvette excites fluorescence and the resulting light entering the phototube causes the instrument meter to register. With P fully depressed adjust the STD knob until the meter reads a value consistent with the concentration of the sample to be examined. This operation constitutes the sensitivity adjustment. Note the value of the reference standard and immediately before testing each sample first insert the reference solution and then adjust the meter to read the value of the reference standard when button P is depressed. This stable reference solution is necessary because the instrument's sensitivity is dependent in part on the intensity of the mercury vapor lamp which in turn depends on the voltage and frequency of the power supply. Since the latter factors may vary with time it is advisable to check the sensitivity of the photofluorometer immediately before each test by inserting in the cuvette well a stable reference solution, e.g. quinine sulfate when examining thiamine. A solution containing 0.05 mg of quinine sulfate dissolved in one liter of N/10 sulfuric acid is set at 60 on the scale. It is recommended that this stable reference solution be utilized solely as a means of checking and maintaining the instrument sensitivity.

at the selected value. The practice of assigning a definite value to the ratio of instrument responses to quinine and to vitamin is unnecessary and of questionable merit. A much more reliable calibration consists of running a calibration curve with standard vitamin solutions.

- d. Insert a cuvette containing a reagent blank, cap the cuvette, fully depress button P, and read the meter as soon as the pointer has come to rest. (The glassware, solvent, and the other chemicals used for preparing the sample may themselves exhibit some measure of fluorescent response, and the activity of these 'blank' constituents must therefore be considered in all measurements. In practice the blank solution is prepared by combining in the same manner and amounts all of the ingredients of the corresponding sample solution except the considered constituent. In the determination of thiamine, the fluorescent compound is thiochrome, formed by oxidizing thiamine with potassium ferricyanide. In this instance the 'blank' solution is prepared exactly the same as the sample, using water in place of the sample to be analyzed. Thus the blank is slightly fluorescent. The corresponding meter deflection caused by the thiochrome is equal to the difference between the meter deflection of the blank and the deflection occasioned by the corresponding sample.)
- e. Set the sensitivity of the instrument with the stable reference standard as in step (c) above.
- f. Now read an unknown. Insert a cuvette containing a solution of the considered constituent (e.g. thiochrome for thiamine) of approximately the concentration to be considered in subsequent tests. Cap the tube, fully depress button P, and then read the meter as soon as its pointer comes to rest. (The best accuracy is obtained when the range of concentration to be examined gives a convenient deflection, and it is good practice not to attempt measurements at very high concentrations, since the fluorescence-concentration curves depart from linearity at high concentrations. For example, the recommended range of concentration for thiamine is not more than 2 gamma oxidized to thiochrome in 10 ml of isobutanol. Higher concentrations should be diluted to this range.)
- g. Before all subsequent readings, set the sensitivity of the instrument by means of the stable reference standard.

## 4 Computation of Results

### a. Calibration Curves

We recommend the use of a calibration curve derived from samples of standards run through the analytical procedure. (The reason for this recommendation is that fluorescence may be affected markedly by slight alterations in analytical conditions. Hence computation of a standard  $K$ , i.e. units of solution concentration per meter scale unit, is inherently liable to error. In addition, the fluorescence galvanometer reading relation is not always linear.)

### b. Calculated Constant

We do not recommend this calculation for the reasons stated above. However, it is given here for the sake of completeness.

- (1) Let  $K$  = units of solution concentration per scale division. Let  $A$  and  $B$  be scale readings for two known concentrations of the substance under consideration,  $S_1$  and  $S_2$  respectively. Then  $K = \frac{B - A}{S_2 - S_1}$ .
- (2) Let  $K$  = constant as defined above.  
 $U$  = scale reading for unknown  
 $B$  = scale reading for blank

Then

$$\text{Concentration of unknown} = K \times (U - B)$$

## 5 Adjustment and Replacement of Parts

### a. Adjustment of Mercury Vapor Lamps Autotransformer to Match Exciting Voltage

The H 4 mercury vapor used in the Coleman Model 1. photofluorometer is adjusted to operate at 120 volts when shipped. If with this adjustment it is connected to a higher supply voltage the intensity of the lamp will be very markedly increased and there will be an appreciable reduction of the lamp life. In instances when it is necessary to use this lamp in conjunction with a line or voltage regulator delivering other than 120 volts the connections on the autotransformer must be changed to fit the alternate voltage. This adaption can be made by the following procedure:

- (1) Remove the two screws which secure the cover on the small autotransformer.
- (2) On the terminal board inside this transformer will be found terminals marked 110, 115 and 120. The instrument is shipped with a green lead connected to the 120 terminal.
- (3) If the supply voltage is between 108 and 113 volts move this lead to the 110 terminal. (This same connection will also assure maximum possible brilliance when it is necessary to operate at less than 108 volts.) If the supply voltage is between 113-118 volts connect the green lead to the terminal marked 115. If the supply voltage is above 118 volts leave the green lead connected to the terminal marked 120.
- (4) After making the proper connection be sure that all terminals are screwed tight and replace the cover. The instrument is now ready for operation.
- (5) Transformers are individually adjusted to insure the maximum possible voltage control and as a result the output voltage of a given transformer may exceed 115 volts. Therefore when using the transformer be sure to ascertain its actual output voltage which is given on a plate attached to the front of the unit and then adjust the lamp autotransformer to correspond.

### b. Replacement of Mercury Vapor Lamp

The lamps used in the Model 1 is an AH 4 medium high pressure quartz mercury arc. The life expectancy of this lamp is in the neighborhood of 1000 hours. When failure occurs as evidenced by unsteadiness of illumination, obscuring of the walls of the central quartz capillary and finally refusal to start when lamp has cooled replacement should be made with another lamp of exactly the same type. Specify serial number of instrument when ordering.

- (1) Remove the five screws around the edges of the panel.
- (2) Place the instrument at the edge of the table with the left hand side of the box far enough off the table so that the hole in the bottom from which the rubber covered cable to the transformer issues is just clear of the edge of the table. This will free the cable and facilitate removal of the panel.
- (3) Pull the panel straight upward until checked by the battery cable. Bring the back edge of the panel up and forward and lay it on its face on the top of the box which should be moved back onto the table.
- (4) The mercury lamp is secured by two spring retaining straps which hold it to an aperture plate. Connections are made by twisted leads which connect to a terminal strip at the side of the chassis adjacent to the lamp.
- (5) Disconnect the leads from the terminal strip.
- (6) Loosen the nuts which hold down the lamp retaining straps. It is not generally necessary to remove the straps. With the straps loose slide the lamp endwise till clear of the straps.
- (7) Insert a new lamp into the straps with the small conical tip on the inner quartz capillary pointed away from the aperture plate and the two support rods inside the lamp approximately equidistant from the plate. Center the lamp by sliding

endwise so that only the transparent part of the capillary is in line with the three-quarter inch hole in the aperture plate the metalized end caps being clear of the hole

- (8) Push down the retaining straps and pull the nuts down gently not tightly A gentle pressure here will hold the lamp firmly, excessive pressure may shatter the outer tube of the lamp
- (9) Slide the spade lugs on the connecting leads under the terminal screws and tighten Be sure that spade lugs do not touch each other or touch the rivets holding the terminal strip to the chassis
- (10) Replace the panel by repeating instructions 1, 2 and 3 in reverse (Note Be sure battery cable does not crowd over so as to interfere with fan or mercury lamp)

#### c Replacement of Vacuum Tube

The vacuum tube used in the Model 12 has been carefully selected and rigorously tested before incorporation in the instrument Although a standard radio tube of the correct type may be used in this instrument it is advisable to secure replacement tubes through a Coleman dealer When ordering a replacement tube always specify the serial number of the instrument Symptoms of tube failure are (1) inability to reach the proper adjustments even though the batteries check well, and (2) when filament burns out continuous off scale reading to the right not responsive to position of the BAL or BLK knobs When it becomes necessary to replace the vacuum tube

- (1) Remove panel as described in previous sections
- (2) The vacuum tube is surrounded by a perforated metal shield with a removable end cap Remove this cap exposing the grid connection to the tube Pull the grid connection off the tube Remove the perforated shield by pulling it up off the tube The tube may now be removed from its socket
- (3) After inserting the new tube be sure to remove any conducting film existing on outside of bulb between grid terminal and base by wiping gently with tissue or a clean cloth before replacing the shield Replace the grid connection before placing the end cap on the shield being careful to avoid again touching glass
- (4) Replace panel as described in previous section

#### d Replacement of Batteries

The Model 12B requires the following batteries for its operation

- One 45 v B Battery type 5308
- One 3 v A Battery type 2F2H
- One 90 v C Battery type Z60BP
- One 6 V Bucker Battery type F4X

The life expectancy of the 2F2H should be above 400 hours of service that of the F4X about 1000 hours The other two batteries should be good for over a year Failure of the 2F2H or the #5308 will be shown by inability of the BAL knob to bring the meter up to zero while failure of the F4X will be shown by the inability of the BAL knob to bring the meter down to zero The Z60BP should be replaced when its voltage drops below 50 volts When it becomes necessary to change batteries

- (1) Remove the five screws around the edge of the panel and disassemble as described in previous sections
- (2) Disconnect battery to be replaced. It is wise also to check the voltage of all the batteries and replace any which show a voltage less than three fourths of their rated voltage
- (3) Remove screws securing hold-down strap of battery to be replaced. Remove strap thus freeing battery (The Z60BP is held down by a wooden block instead of a strap)

- (4) Replace battery with identical battery. Replace strap and reconnect battery. To facilitate this the following table gives the color coding of the battery cable

| COLOUR |                       | CONNECTION |
|--------|-----------------------|------------|
| Gray   | -A                    | ( F-H)     |
| Black  | +A                    | ( F H)     |
| Blue   | -B                    | ( 308)     |
| Yellow | +2- $\frac{1}{2}$ V B | (5308)     |
| Red    | +45 V B               | (5308)     |
| Brown  | -C                    | (260BP)    |
| Green  | +C                    | (260BP)    |

The leads to the F4X battery are separate comprising a twisted red and black pair the red being plus

- (5) Replace panel as described in previous section

## F FLUOROMETRIC METHODS FOR QUALITATIVE USE AND FOR QUANTITATIVE MEASUREMENTS IN THE FIELD—VISUAL FLUOROMETRY

### Principle

Qualitative fluorometry is best done by eye because it is necessary to distinguish fine shades of color. No photoelectric fluorometer of which we are aware is suitable for qualitative use.

For many quantitative purposes in the field visual fluorometry is equally as satisfactory as photoelectric fluorometry. The chief reason is that under field conditions line voltage is often extremely erratic. This impairs the accuracy of photoelectric fluorometers but has no effect on the accuracy of visual fluorometry.

### Apparatus

- 1 A suitable darkened room for reading unknown solutions. This need not be completely light tight. We have used successfully x-ray rooms, photographic dark rooms, bath rooms, abandoned kitchen cupboards and U.S. Army command post tents.
- 2 A source of ultraviolet light consisting of a General Electric Type R4 mercury vapor bulb enclosed in a box with louvers for air cooling and a Wood filter for absorbing visible light. A General Electric Autotransformer Type 59G18 is installed between the bulb and any suitable source of 110 volt AC. In the field we have used simple motor generators without line stabilization.
- 3 The fluorometric comparator consisting of a thin wooden rack containing 10 Corning glass stoppered thin walled 1 ml test tubes outside diameter 10 mm and length 75 mm. Holes are bored so that these tubes which contain standard solutions are flanked by empty holes to contain 10 x 75 mm thin walled test tubes containing the unknown solution. The rack is constructed so as to slide conveniently in a direction perpendicular to the ultraviolet ray and 8 cm from the Wood filter.
- 4 Suitable beakers, delivery pipettes and other usual appurtenances including a large number of 10 x 75 mm thin walled test tubes. One special pipette is of great assistance. This is an automatic air syringe pipette the use of which is described in the section on volumetric apparatus.
- 5 A small spotlight in the darkened room in such a position that a recorder can easily take down notes without disturbing the operator who is doing the fluorometry.

**Reagents**

- 1 These are different for different estimations but in all cases thoroughgoing reagent blanks must be run on each new batch and if there is a perceptible blank the reagents must be either replaced or purified until the blank is negligible
- 2 Standard solutions are prepared by running solutions of known concentration through the whole method. The final solutions are placed in the glass stoppered 2 ml test tubes up to a convenient level usually 1 ml, and the wooden racks are clearly labeled in fluorescent paint or ink with a number equivalent to the concentration of unknown substance in 100 ml of original solution. This procedure automatically takes care of calculations and also includes any correction necessary for systematic losses in the methods. The preparation of standard solutions will be described under individual methods

**Procedure**

- 1 Each method is carried out as described in detail in subsequent sections
- 2 As the final step in each method the compound to be measured is extracted into an appropriate solvent and the solvent layer is separated from the aqueous layer, either by simple standing or by centrifuging
- 3 An automatic syringe pipette fitted with a long stainless steel needle is stopped at approximately 1 ml. The tubes containing samples in which the solvent layer has been separated from the aqueous layer are placed in convenient racks in front of the operator. Other racks containing clean 75 x 10 mm test tubes in positions corresponding to the unknowns are similarly placed
- 4 Approximately 1 ml of solvent is transferred to a 75 x 10 mm test tube. Excess solvent in the automatic syringe pipette is removed by several up and down motions of the plunger and the next sample is transferred. The transfer pipette need be rinsed with blank solvent only when it is expected that there will be very wide variations in concentration from unknown to unknown e.g. if fasting and loaded samples are being run at the same time
- 5 The rack containing 75 x 10 mm test tubes is now transferred to the darkroom after the notebook is appropriately ruled and after the ultraviolet bulb has been turned on for at least 10 minutes to allow thorough warming up
- 6 After entering the darkroom allow about 5 minutes for the eyes to become adapted to the dark
- 7 Comparison is made by inserting the unknown solution from empty hole in the comparator until a match is obtained. It is usually wise to commence reading with the unknown in position between two standard solutions both of which are unquestionably brighter than the unknown and then to move the unknown down the rack for final matching
- 8 The recorder should record the value for the unknown in the appropriate place in the notebook. It will be recalled that these readings represent without further calculation the concentration of the original unknown solution

**Calculations and Examples**

- 1 These will be discussed individually in individual methods

**Precautions**

- 1 Some of the standard solutions are themselves unstable and should be checked daily or even oftener against freshly prepared standards. This is particularly true of riboflavin. Quinine is apparently stable indefinitely under usual conditions but we have seen quinine standards deteriorate in very hot weather
- 2 The standard solutions and unknown solutions must be at approximately the same temperature during reading because fluorescence is markedly affected by changing the temperature of the solution. The temperature will be the same if reagents, samples and standards are kept on the same bench

## G FLAME PHOTOMETRY—THE BECKMAN FLAME PHOTOMETER ATTACHMENT

### 1 General Principles of Flame Photometry

The purpose of flame photometers now commercially available is to introduce inorganic substances into a flame at a uniform rate so that light produced thereby will be constant. The flame used is very hot so that the spectral lines are excited much more strongly than with a cool flame and more lines of more elements can be excited to a usable intensity than with an ordinary cool flame. In this way it is possible to make instantaneous spectrophotometric measurements whereas the usual method requires a spectrogram followed by a densitometer measurement and the results have to be interpreted with the aid of sensitometry measurements made on each plate at all wave lengths used. A flame source has certain advantages over the customary arc and spark sources for chemical analysis. Unnecessary lines are not excited, the temperature is lower so the continuous background is lower, the temperature and rate of feed of specimen can be held constant over a long period allowing numerous individual line measurements to be made in succession, the flame and rate of feed can be exactly reproduced at any time, being extended in area the flame does not require precise optical alignment and slight sidewise motions do not alter the luminosity and especially the specimen to be analyzed is introduced into the region of excitation at a uniform rate with the more volatile and less volatile constituents in perfect admixture. (In other sources the more volatile constituents vaporize first so that photographic or integrating methods must be used and the sample must all be consumed for accurate determinations.)

Of the several good instruments now on the market we shall describe in detail only one, the flame photometer attachment of the Beckman spectrophotometer.

### 2 General Construction of the Beckman Flame Photometer

#### a. Photometer Attachment

With the flame photometer attachment the Beckman Type DU Spectrophotometer becomes an emission spectrophotometer and percent transmittance becomes light intensity. Basically there are four separate parts (see Fig. 8).

- (1) The Beckman Type DU Spectrophotometer for measuring the intensity of emitted light. Some modifications are required as described below.
- (2) A burner to produce emission from the sample. The burner of this instrument burns illuminating gas with air or oxygen, the mixture emerging from 30 ports. Through the burner seven ducts carry fine droplets of mist containing the specimen under analysis. The burner is of stainless steel and is water cooled. A water cooled chimney which may be extended out of doors surrounds the flame and protects against heat and fire hazard. The burner and chimney are arranged as an integral unit with gas, oxygen, and water cooling lines going to the control unit and two vapor ducts going to the two sprayers. One of the ducts is for unknown specimens, the other for standard solutions.



**Reagents**

- 1 These are different for different estimations but in all cases thoroughgoing reagent blanks must be run on each new batch and if there is a perceptible blank the reagents must be either replaced or purified until the blank is negligible
- 2 Standard solutions are prepared by running solutions of known concentration through the whole method. The final solutions are placed in the glass stoppered 2 ml test tubes up to a convenient level usually 1 ml and the wooden racks are clearly labeled in fluorescent paint or ink with a number equivalent to the concentration of unknown substance in 100 ml of original solution. This procedure automatically takes care of calculations and also includes any correction necessary for systematic losses in the methods. The preparation of standard solutions will be described under individual methods.

**Procedure**

- 1 Each method is carried out as described in detail in subsequent sections
- 2 As the final step in each method the compound to be measured is extracted into an appropriate solvent and the solvent layer is separated from the aqueous layer either by simple standing or by centrifuging
- 3 An automatic syringe pipette fitted with a long stainless steel needle is stopped at approximately 1 ml. The tubes containing samples in which the solvent layer has been separated from the aqueous layer are placed in convenient racks in front of the operator. Other racks containing clean 75 x 10 mm test tubes in positions corresponding to the unknowns are similarly placed
- 4 Approximately 1 ml of solvent is transferred to a 15 x 10 mm test tube. Excess solvent in the automatic syringe pipette is removed by several up and down motions of the plunger and the next sample is transferred. The transfer pipette need be rinsed with blank solvent only when it is expected that there will be very wide variations in concentration from unknown to unknown, e.g. if fasting and loaded samples are being run at the same time
- 5 The rack containing 75 x 10 mm test tubes is now transferred to the darkroom after the notebook is appropriately ruled and after the ultraviolet bulb has been turned on for at least 10 minutes to allow thorough warming up
- 6 After entering the darkroom allow about 5 minutes for the eyes to become adapted to the dark
- 7 Comparison is made by inserting the unknown solution from empty hole in the comparator until a match is obtained. It is usually wise to commence reading with the unknown in position between two standard solutions both of which are unquestionably brighter than the unknown and then to move the unknown down the rack for final matching
- 8 The recorder should record the value for the unknown in the appropriate place in the notebook. It will be recalled that these readings represent without further calculation the concentration of the original unknown solution

**Calculations and Examples**

- 1 These will be discussed individually in individual methods

**Precautions**

- 1 Some of the standard solutions are themselves unstable and should be checked daily or even oftener against freshly prepared standards. This is particularly true of riboflavin. Quinine is apparently stable indefinitely under usual conditions but we have seen quinine standards deteriorate in very hot weather
- 2 The standard solutions and unknown solutions must be at approximately the same temperature during reading because fluorescence is markedly affected by changing the temperature of the solution. The temperature will be the same if reagents, samples and standards are kept on the same bench

## G FLAME PHOTOMETRY—THE BECKMAN FLAME PHOTOMETER ATTACHMENT

### 1 General Principles of Flame Photometry

The purpose of flame photometers now commercially available is to introduce inorganic substances into a flame at a uniform rate so that light produced thereby will be constant. The flame used is very hot so that the spectral lines are excited much more strongly than with a cool flame and more lines of more elements can be excited to a usable intensity than with an ordinary cool flame. In this way it is possible to make instantaneous spectrophotometric measurements whereas the usual method requires a spectrogram followed by a densitometer measurement and the results have to be interpreted with the aid of sensitometry measurements made on each plate at all wave lengths used. A flame source has certain advantages over the customary arc and spark sources for chemical analysis. Unnecessary lines are not excited, the temperature is lower so the continuous background is lower, the temperature and rate of feed of specimen can be held constant over a long period allowing numerous individual line measurements to be made in succession, the flame and rate of feed can be exactly reproduced at any time, being extended in area the flame does not require precise optical alignment and slight sidewise motions do not alter the luminosity and especially the specimen to be analyzed is introduced into the region of excitation at a uniform rate with the more volatile and less volatile constituents in perfect admixture. (In other sources the more volatile constituents vaporize first so that photographic or integrating methods must be used and the sample must all be consumed for accurate determinations.)

Of the several good instruments now on the market we shall describe in detail only one, the flame photometer attachment of the Beckman spectrophotometer.

### 2 General Construction of the Beckman Flame Photometer

#### a. Photometer Attachment

With the flame photometer attachment the Beckman Type DU Spectrophotometer becomes an emission spectrophotometer and percent transmittance becomes light intensity. Basically there are four separate parts (see Fig. 8):

- (1) The Beckman Type DU Spectrophotometer for measuring the intensity of emitted light. Some modifications are required as described below.
- (2) A burner to produce emission from the sample. The burner of this instrument burns illuminating gas with air or oxygen; the mixture emerging from 30 ports. Through the burner seven ducts carry fine droplets of mist containing the specimen under analysis. The burner is of stainless steel and is water cooled. A water-cooled chimney which may be extended out of doors surrounds the flame and protects against heat and fire hazard. The burner and chimney are arranged as an integral unit with gas, oxygen, and water cooling lines going to the control unit and two vapor ducts going to the two sprayers. One of these ducts is for unknown specimens, the other for standard solutions.

- (3) Sprayers are made of Pyrex glass and are conveniently located to facilitate filling and washing and observation of the specimens
- (4) The control unit has the following on its panel a gas manometer to read the gas pressure at the burner nozzle an oxygen gauge to read the oxygen pressure at the burner nozzle an air pressure gauge to read the pressure at the sprayer nozzles The air valve is a cut off valve only it is used either open or closed The air pressure regulator controls the air pressure reducing it from 30 to 50 lbs/sq in at input to 15 at output The gas valve controls and cuts off the gas The oxygen valve controls and cuts off the oxygen




Fig 8—The flame photometer attachment for the Beckman Type DU Spectrophotometer  
(Reproduced with the permission of the manufacturers)

### 3 Setting Up the Instrument

#### a. Alterations in the Beckman DU Spectrophotometer

The conversion of the spectrophotometer into an emission measuring instrument by addition of the flame attachment involves four small changes in the spectrophotometer

- (1) The flame used in emission is not as steady in intensity as the hydrogen or tungsten lamp used in absorption measurements. The galvanometer fluctuations therefore increase considerably and must be minimized or smoothed out. This is done by changing the 2000 megohm resistor in the phototube compartment to a 10 000 megohm resistor. The change of resistors gives the galvanometer five times the previous time lag enabling it to be read with a high degree of reproducibility.

- (2) An aluminum shield is placed over the phototube compartment to protect it from the heat of the flame
- (3) The spectrophotometer has to be raised onto the platform of the control panel.
- (4) When the flame attachment is working the entrance mirror of the spectrophotometer has to be adjusted to an angle optimal for measuring the emitted light

#### b Preliminary Adjustments and General Precautions

- (1) Place the combined unit of control panel flame tower and sprayers on a rigid table convenient to sources of gas air water and electricity To avoid erratic results the room in which the instrument is to work must be as free of dust and drafts as possible
- (2) Place the spectrophotometer on the base of the control panel and its left end flush with the right margin of the gas manometer
- (3) The hoses are attached securely Scotch tape or small clamps serve well The air line must be securely clamped both on the input side and the exits to sprayers The air pressure must be turned off (gauge reading zero) before removing a hose from the sprayer Otherwise a large air flow will occur and water from the high pressure bubbling towers may be carried through. A drip column has been added to catch a reasonable amount of such water However it can only be drained by removing the back of the cabinet
- (4) The water is turned on at about one liter per minute Too high a flow is objectionable Water to drain should be hot enough but not boiling with burner running If too much is run through the bubbling towers may become over filled.
- (5) Turn on gas low just enough to light and light it When it is lit turn gas up to the fiducial marks on the manometer Always remember that water must be turned on before the flame is lit and the flame must be turned out before the water is turned off
- (6) Turn oxygen on to about 2 inches of water and check gas pressure on the manometer In case the gas pressure rises too high the flame may back fire as all oxygen flames do This is harmless in the present burner the small mixing chamber in the base of the burner has low pneumatic resistance to the outside air and no damage will result

## 4 Procedure for Routine Analyses

This laboratory prefers to use direct comparison between unknown samples and aqueous standard solutions rather than an internal standard of lithium added to each unknown specimen Some of the reasons for this choice are to be found in the discussion section below To make a run of analyses

- a. Prepare a group of samples the approximate concentrations of which are known properly diluted to bring them within the linear C-T range Standards are prepared from a stock standard so that they also be in the correct range The special conditions of slit opening flame temperature etc as determined by trial and error are produced in the instrument and the samples and standards are put into the 5 ml beakers provided
- b Turn switch to check
- c Turn on water
- d Plug in atomizer heater
- e Light burner low at first then high
- f Turn on oxygen
- g Allow 5 minutes for warming up

## 62 Instrumentation

- h Turn on atomizer, read an unknown sample on the "C %T" scale and immediately after two standards whose concentrations i.e., "C %T" readings lie on either side of the sample Repeat a number of times until all readings are steady and reproducible and good checks are obtained.
- i Repeat this process with the rest of the unknown samples

## 5 Computations

If the graph of "C %T" is linear in the range used then for each sample the concentration can be calculated from the simple ratio (standard galvanometer reading)/(sample galvanometer reading) In this case the ratios of sample reading with either of the standards bracketing it and the reading from the "C %T" standard curve should check. If a determination is done in which it is not possible to obtain linear readings the simplest method of computation is to use a calibration curve prepared from standards When graphing is resorted to, a semilogarithm or logarithmic plot can be used, the criterion being which gives the greater degree of linearity

## 6 Discussion of Sources of Error in Flame Photometry

### a Metals Which Can Be Measured

In clinical chemistry about all the flame photometer is really useful for is estimating sodium, potassium and calcium as shown in Table 1 Other alkali metals and alkaline earths can be measured readily, but for other metals such as copper and iron chemical methods are far more reliable than the flame photometer Since sodium, potassium and calcium are by far the most common metals measured with the flame photometer the rest of this discussion is restricted to them

### b Sources of Instrumental Error

Readings for any given sample in the flame photometer may be affected by many variables among which are quality of gas rate of water flow through chimney and burner sensitivity setting of the instrument width of slit wave length setting and pressures of gas oxygen, and air The adjustment of these various components depends not only upon the element to be analyzed but also upon other interfering substances which may be present Therefore no hard and fast rules can be formulated for the operation of the instrument but a few general precautions can be stated to help in obtaining consistent results

- (1) In our experience ordinary city gas is too variable in pressure and in quality to be completely satisfactory We have found it necessary to purchase tanks of gas of constant source and composition By means of a suitable reducing valve this gas is led into the instrument With some kinds of such gas it may be necessary to adjust the burner
- (2) Water flow should be as constant as possible This can be accomplished by using a reservoir with an overflow thus keeping a constant pressure of water at all times
- (3) Even with all variables held constant the readings of a given sample will rise continuously for about  $\frac{1}{4}$  hour However, if the standards and samples are read immediately one after the other the determination can be made with good accuracy without waiting for the readings to become exactly constant
- (4) The lines radiated by each element vary in intensity those of the alkali metals in general being most intense As intensity increases the slit width can be decreased to obviate interference by other elements When this is done, however

TABLE 1  
SENSITIVITY OF FLAME PHOTOMETER WITH DIFFERENT METALS IN HOT FLAME  
(GAS OXYGEN)

| ELEMENT<br>AND<br>COMPOUND        | SLIT WIDTH<br>MILLIMETERS | LINES USED<br>ANGSTROM UNITS | CONCENTRATION               |                             |
|-----------------------------------|---------------------------|------------------------------|-----------------------------|-----------------------------|
|                                   |                           |                              | MEASURABLE<br>PARTS/MILLION | DETECTABLE<br>PARTS/MILLION |
| Sodium                            | 0.102                     | 5890                         | 0.5                         | 0.1                         |
| NaCl                              |                           | 5896                         |                             |                             |
| Potassium                         | 0.102                     | 7665                         | 0.5                         | 0.1                         |
| KCl                               |                           |                              |                             |                             |
| Lithium                           | 0.10                      | 6708                         | 0.5                         | 0.1                         |
| Li <sub>2</sub> SO <sub>4</sub>   |                           |                              |                             |                             |
| Cesium                            | 0.10"                     | 851                          | 0.5                         | 0.1                         |
| Cs <sub>2</sub> CrO <sub>4</sub>  |                           |                              |                             |                             |
| Calcium                           | 0.102                     | 618                          | 3                           | 1                           |
| CaCl                              |                           | 603                          |                             |                             |
| Manganese                         | 0.10                      | 4031                         | 2.3                         | 1                           |
| MnCl                              |                           | 4034                         |                             |                             |
| Copper                            | 0.102                     | 348                          | 30                          | 10                          |
| CuSO <sub>4</sub>                 |                           | 3274                         |                             |                             |
| Magnesium                         | 0.10"                     | 2852                         | 100                         | 25                          |
| MgCl                              |                           | 3838                         |                             |                             |
| Cobalt                            | 0.10"                     | 3520                         | 100                         | 25                          |
| CoCl <sub>2</sub>                 |                           | 3502                         |                             |                             |
| Nickel                            | 0.10"                     | 355                          | 150                         | 40                          |
| NiCl                              |                           | 3493                         |                             |                             |
| Boron                             | 0.105                     | 5481                         | 150                         | 50                          |
| H <sub>2</sub> BO <sub>3</sub>    |                           | 5440                         |                             |                             |
| Iron                              | 0.10"                     | 3735                         | 00                          | 50                          |
| Fe(NO <sub>3</sub> ) <sub>3</sub> |                           | 3737                         |                             |                             |

The results were obtained with aqueous solutions of pure single compounds, with no other ions interfering.

the wave length must be set very carefully. This can be done by setting the wave length approximately with the instrument in operation and then turning it slightly to the right or left until the galvanometer shows maximum deflection.

- (5) The capillary leading into the atomizer is very fine. It has a tendency to become clogged gradually without the operator noticing it. This results in gradually falling readings. Therefore it is necessary to blow out the capillary after every few determinations. This can be done simply remove the stopper holding the capillary and with the air turned on hold a finger over the inside end for a few moments and replace. Recreational sprayers of this type were adopted for several reasons: they will spray at continuous rates for a long time; they require only small samples (5 ml in the present model); they are economical of sample with a turnover of about 5 ml per minute; the residue of the sample is left unchanged in concentration; and they are easy to clean and are chemically inert. If the operator desires the present sprayers can be used for a continuous flow of sample; attach drain tubes to the drain plugs and trickle the sample directly to the sprayer.
- (6) The concentration of the element analyzed should be such that a straight line relationship will exist between  $C$  and  $C^2T$ . For the alkali metals this range is about 0.5-150 mEq/liter.

#### c. Theoretical Considerations Important for Accuracy

- (1) *Wave Length* The wave length selected for each element depends on interfering substances. Therefore although 556 mμ is the most intense calcium line it cannot be used due to the intense sodium line at 589. Potassium is best read at 767 mμ.

- (2) *Spectral Line Interference* While the lines at which the three elements are read are quite widely separated and do not interfere with each other directly there are several indirect types of interference. For example while the alkali metal radiation is of an atomic nature (and therefore quantized into discrete values) the calcium radiation is in large part of the molecular type. This leads to a band spectrum of which the line at  $422\text{ m}\mu$  is merely the 'band head'. Therefore calcium is an interfering substance always to be taken into account. Also while potassium and sodium separately have no radiation at each other's wave lengths each has a large effect in the presence of the other. For example in a urine determination when the ratio of sodium to potassium is about 4:1 the reading of potassium at  $767\text{ m}\mu$  is increased considerably by the sodium present when read in a hot flame.
- (3) *Characteristics of the Flame* Since the calcium spectrum is of a molecular nature i.e. due to vibration and rotation it will increase in intensity with increasing flame temperature. The hottest possible flame must be used in its determination. However potassium radiates more strongly at a relatively low temperature due probably to a favorable tendency toward primary ionization at low temperature. Moreover at low temperature sodium does not affect potassium to any measurable extent. Sodium loses some of its intensity at low temperature but this is offset by the decreased interference due to any band spectra present. The spraying of water reduces the luminosity of the flame as compared with a no-spray flame doubtless due to lowering of the temperature. Since water is always sprayed with a sample the background (no salt) must be read with a distilled water spray not just the bare flame. With a spray of distilled water there is observed a continuous spectrum of rather constant intensity from  $2500\text{ \AA}$  to the deep red except for a hump near  $3100\text{ \AA}$ . This hump is presumed to be a band spectrum perhaps of water since it rises with the addition of water to bare flame. Resolved at narrow slits it shows two peaks at approximately  $3067\text{ \AA}$  and  $3091\text{ \AA}$ . Its height appears to serve as a fiducial point for checking the adjustment of the flame from time to time. Its intensity has not shown alteration with the spraying of a salt that gives no flame line in that region.
- Although we recommend a hot (gas oxygen) flame it is possible to use a cool (gas air) flame. Under these conditions it is necessary to use wider slits and to sacrifice selective and sensitivity except in the case of potassium. To increase luminosity a mirror of thick stainless steel has been mounted back of the flame. This mirror helps when using a cool flame or a feeble radiating element.
- (4) *Standards* While any salts of sodium or potassium can be used in preparing standards this is not true of calcium for which the spectrum will vary with the anion used. In urinalysis the standard must be made up with the calcium in very nearly the same form as it occurs in urine and it is necessary to prepare a standard containing calcium phosphate, chloride and sulfate.
- (5) *Concentrations* Samples should be of such a concentration of the ion under analysis that they fall into a range in which the  $C\%T$  relation is linear. For sodium and potassium this range is from  $0.5\text{ mEq/liter}$  to about  $15\text{ mEq/liter}$ . For calcium however the linear range is at much higher concentrations so high that to obtain it and to get a measurable radiation the sample must be read undiluted. This leads to difficulties due to differences in viscosity between the standards and sample and must be allowed for.
- (6) *Slit Width* Generally speaking a wide slit will increase interference possibilities. Therefore the slit is made as small as is practicable. For calcium a weakly radiating substance a slit of about  $0.2\text{ mm}$  must be used to get a measurable response. For potassium in a cool flame a slit of about  $0.8\text{ mm}$  is used and for sodium in a moderately hot flame a slit of about  $0.4\text{ mm}$  is used. With the Beckman DU Spectrophotometer the intensity of light received by the photo-

cells from a line is directly proportional to the slit width. However, received from a continuum is proportional to the square of the slit width. Therefore, packing out lines against a continuous background it is as narrow a slit as is feasible.

## H BALANCES

Standard textbooks cover very well the calibration of analytical balances. The torsion balance and damped balances in metabolic studies that a brief discussion of them is included here.

### 1 The Roller Smith Precision Torsion Balance Model B

#### a General Construction of the Instrument

In this balance the inside end of a spiral spring is attached to the axis of rotation. The outside end is attached to an index pointer (see Fig. 9). Balances of this type for maximum capacity and sensitivity may be purchased. We have found that for rapid accurate weighing of small amounts of material especially when sensitivity is 0.2 per cent but the final accuracy of any weight measurement depends on so many factors that it is impossible to say that the instrument has an inherent percentage accuracy. With the vernier the scale may be read to 0.1 scale divisions.

#### b Setting Up and Calibrating the Balance

- (1) Set the instrument on a table high enough so that the operator is level with the balance pointer. It is a good idea to lay the balance on a leveling case.
- (2) Be certain that the table is firm.
- (3) Adjust the leveling screws until the bubble of the spirit level is centered.
- (4) Release the beam by rotating the beam clamp lever until it points forward or back.
- (5) Adjust the index lever until the zero on the vernier coincides with zero on the dial.
- (6) The balance pointer should now coincide with its zero line. If not, turn the zero adjuster knob on the back of the case until the beam is horizontal.
- (7) The balance is now ready to operate. The index lever can be shifted to positions by loosening the knurled nut.
- (8) The pan and compensating hook of the same weight are selected.
- (9) The balance is calibrated with the pan in place on the beam. It is adjusted to zero without both the pan and compensatory hook. A standard weight should be purchased with each balance.

#### c Operating the Instrument

- (1) Before any weight is placed on the hook the index lever is rotated to the vernier nearly to the highest value on the dial so that the upper limit of travel is reached.
- (2) Hang the material to be weighed on the hook and rotate the index lever until the beam drops down and the balance pointer again indicates zero.
- (3) The index vernier then indicates the weight directly on the dial.

#### d Disassembling the Roller Smith Balance

- (1) Remove the front spindle knob taking off the index lever.
- (2) Remove the small plate on the left of the case which will be removed.






Fig 3—The Roller Smith Precision Balance Model B (Reproduced with the permission of the manufacturers)

- (3) The back plate screws should now be taken out and the mechanism removed from the rear taking care that the beam hooks are not injured when the mechanism is withdrawn from the case
- (4) Jewels and pivots in the balance should never be oiled as oil at these points will in time collect dust and become gummy thus introducing friction into the bearings. The central knob which rotates the spindle should always be perfectly free to rotate and when an instrument is being overhauled it is a good plan to lubricate this spindle which has a bearing in the cross arm. The lubrication should be done with a little tallow.

## 2 Damped Balances

### a. General Principles

For the sake of speed without sacrifice of accuracy it is convenient to use damped balances either magnetic or air damped. We prefer the air damped balance in which an inner and outer damping cylinder are never in contact with each other and are affected only by air friction. The action is so adjusted that the balance remains stationary after only one or two oscillations.

### b. General Construction of the Sartorius Air Damped Balance

This balance has a multiple rider device for weights from one gram down to 10 mg. For simplicity in reading an illuminated projection reader is switched on automatically when the beam is released making values from 0.1 to 10 mg easy to read. All parts of the reader are so arranged that opening of the doors has no effect on its adjustments. Small errors in the zero point may be corrected by means of a special knob which is operated from the outside. This balance has a capacity of 100 gm with a sensitivity of 0.1 mg.

### c. Adjusting the Balance

- (1) Level the case with the pillar on the bench by means of the leveling screws on the legs.
- (2) Place the luted lock knob entirely into the hollow body bar in the middle of the plane and the two lock bolts into the holes made for this purpose. These are on the left and right hand side of the pillar.
- (3) Switch off the lock knob and if the control bar is inside the pillar the lock mechanism rises and drops regularly.
- (4) Fit the beam on the center agate knife edge with the balance in the locked state and then check it.
- (5) Lock the balance again and fit the suspensions so that the holes are in front and side. Then hang the pans in place.
- (6) On the right and left hand side of the pillar top are two very light aluminum bells that have two gaps at the bottom.
- (7) Remove the bottom and put them into the outside cylinder in order to hang them on the respective suspension hooks.
- (8) After hanging the inside cylinders shut the bells with the previously removed bottoms but before fixing the pans it is advisable to make sure that the marks on each bell are in the right position.
- (9) Make sure when releasing the balance that the inside cylinders do not touch the outside ones. Then be certain that the balance rests at zero. If not use the correction screws at either the left or right hand side of the beams.
- (10) The balance is then put in order by turning the knob outside the case.
- (11) The annular riders for external mechanical application of the weights are hung from the heaviest to the lightest that is from 100 to 10 mg on the double hook.

### d. Calibrating the Balance

- (1) Before starting to use the balance set at zero by turning the sensitivity screws on the left or right hand side of the beam.



- ( ) Place a 10 mg weight on the right pan release the balance and check the weight on the glass screen indicator It should read exactly 10 mg If not adjust the calibration screws until the reading is correct
- (3) Then repeat the process by placing the weight on the left pan
- (4) If the oscillations are too small raise the sensitivity screw in the center of the beam and if too violent lower the screw
- (5) Be sure the pivot is not twisted otherwise the zero point will be off

#### e Operation of the Instrument

The operation is the same as for any analytical balance

## I VOLUMETRIC ANALYSIS

### Reference

Peffer E L and Mulligan G C U S Bureau of Standards Testing of Glass Volumetric Apparatus Circular C434 Washington 1941 U S Government Printing Office

### 1 Definition of Volumetric Analysis

By volumetric analysis is meant quantitative analysis by the use of a measured amount of standard solutions The two aspects discussed below are first volumetric apparatus and second the preparation and standardization of common solutions

### 2 Volumetric Apparatus

#### a Pipettes

- (1) *Common Pipettes* Many kinds of pipettes are used of which some of the more common kinds are illustrated in Fig 10 Calibrated pipettes differ in their accuracy the common commercial kind (e.g. Exax ) having the widest tolerance and those calibrated to Bureau of Standards specifications the closest tolerance (e.g. Normax ) Table 1 lists permitted tolerance for pipettes burettes volumetric flasks and graduated cylinders All pipettes are calibrated at 20°C and are marked either TD (to deliver) or TC (to contain) TD (to deliver) does not specify how delivery is to be made Hence there are two kinds of TD pipettes (a) pipettes with frosted tips must be blown out and (b) pipettes with no frosted band on the stem must be allowed to drain against the side of the vessel for the time specified on the pipette Bureau of Standards pipettes have been calibrated at the Bureau of Standards and carry individual calibration certificates Pipettes calibrated by the company (e.g. Normax ) are guaranteed to fall within Bureau of Standards tolerances When pipettes are marked TC (to contain) it is necessary either to use a nonwetting liquid such as mercury or to wash the pipette out pooling the wash fluid with the original sample
- ( ) *Syringe Pipettes* It is very convenient to use syringe pipettes for sampling (See Fig 11 Procurable from Macalaster Dicknell Co Cambridge Mass) These syringe pipettes were developed by Krogh and are calibrated in the same manner as ordinary pipettes They can be adjusted to any volume from 0.1 ml to 100 ml For the analysis of food and feces we have found it convenient to calibrate a Luer Lok syringe because of the large bore at the tip of the syringe The syringe makes the sampling much easier

#### b Burettes

- (1) *Common Burettes* Burettes are procurable with different tolerances of accuracy (Table 2) but there is no real economy in purchasing other than those with the greatest accuracy (e.g. Normax ) since few laboratories use many burettes.

- (2) *Microburettes* When titrating small amounts of solution (as in serum potassium analysis when 0.05 ml of 0.1 N sodium thiosulfate is needed) we have found it convenient to use a Rehberg microburette (See Figs 1° 13. Procurable from Macalaster Bicknell Co., Cambridge, Mass.) This is a capillary burette with a capacity of 0.1 or 0.2 ml and is graduated into fractions of 0.001 ml. At the lower end of the burette the capillary tubing bore becomes larger and is filled with mercury so that a turn of the screw causes the mercury to deliver aqueous

## VARIETIES OF PIPETTES





Fig 10—Common kind of volumetric pipettes.




## SYRINGE PIPETTES

C E GORDON


Fig 11

SPRING TYPE  
JOINT



TO AIR  
PRESS

SCALE  
cm.




SCREW TYPE  
MERCURY PLUNGER

REHBERG MICRO BURETTE

G. E. BORDON

# KROGH BURETTE SYRINGE TYPE



G. E. GORDON


Fig. 12.



solutions with the greatest of accuracy. The tip of the burette is on a greased swivel joint that makes titrating convenient.

- (3) *Syringe Screw Burette* This was popularized by Krogh. (For an illustration see Fig. 13.) This burette has a capacity of 2 or 3 ml and is graduated to 0.01 ml. In titrations a twist of the screw delivers any fraction of a drop with great precision.

## MICROMETER BURETTE



—G.E. GORDON—

Fig. 14—Scholander's micrometer burette. (Redrawn with permission of the author.)

### 3 Calibration of Volumetric Apparatus

For most purposes volumetric apparatus is calibrated by weighing mercury or water delivered under the conditions of use of the apparatus. At the same time the temperature of the fluid is measured and the specific gravity of the mercury or water is found in standard tables. The volume delivered then equals the weight of the fluid divided by the density at the given tem

perature For quick and extremely accurate calibration of micro apparatus especially pipettes and burettes the Scholander micrometer burette is excellent (Scholander P F 'Microburette' Science 95 177 178 Feb 1942 See Fig 14)

### a. Principle of the Scholander Micrometer Burette

Mercury is displaced by the spindle of a machinist's micrometer The great precision of this kind of micrometer enables volumetric readings to be made with an accuracy of 0.0003 ml

### b. Construction of a Micrometer Burette (Fig 14)

- (1) Procure a machinist's micrometer with jaws of capacity 1 to 2 inches The best procurable micrometer should be purchased since only good ones have complete linearity of scale
- (2) Replace the anvil with a set screw (#1 Fig 14)
- (3) The lower part of the burette is a cylindrical reservoir with a ground plane thick bottom and plane ground opening and an internal diameter just enough to clear the spindle
- (4) A fiber ring (#2) is placed on the closed end of the reservoir and next to it is placed a steel disk (#3) with a center punch to receive the set screw
- (5) A fiber washer (#4) greased with a heavy stop-cock grease is placed at the open end
- (6) Fasten the reservoir in the micrometer frame Mercury will not leak through the spindle bearing when it is properly greased
- (7) If a small spindle is used the fiber washer (#4) must be reamed to a perfectly tight fit around the spindle which will then be proof against leakage of mercury
- (8) The micrometer frame is clamped to a convenient stand and by movement of the spindle mercury is sucked in to fill the reservoir
- (9) When the spindle is screwed out far enough to be level with the fiber washer air bubbles can be readily brought up into the funnel-shaped opening to the capillary For proper functioning of the burette it is essential that the whole mercury containing system be free of air The reservoir is air free when there is no visible bouncing when the mercury column is moved quickly from one place to another in the dry burette The last traces of air are removed by applying vacuum suction at the burette tip while the micrometer spindle is moved in and out

### c. Calibration of the Microburette

- (1) Note that the precise size and shape of the glass part are immaterial It is only the mercury expelled from the micrometer that is actually measured
- (2) Bring the micrometer close to 10 making sure that mercury fills the burette exactly to the tip Read
- (3) Into a weighed vessel force mercury until the micrometer reads close to 0 Read again
- (4) Measure the temperature of the burette by laying a thermometer alongside
- (5) Reweigh the vessel
- (6) From standard tables find the density of mercury at the observed micrometer temperature
- (7) Compute the calibration constant of the burette from the equation

$$\text{Milliliters per scale division} = \frac{(\text{mg Hg})}{(\text{density Hg})} \div (\text{Reading 2 minus Reading 1})$$

- (8) For greatest accuracy recalibrate by successive tenths of an inch with weighings each time

**d Calibration of Apparatus With the Scholander Micrometer Burette**

- (1) Note that for calibration of other apparatus the temperature of the micrometer and of the apparatus calibrated must be the same but that the absolute temperature is immaterial. This is so because the initial calibration of the micrometer burette gives calibration constant equating scale divisions of the micrometer with volume.
- (2) As an example let us assume that a serological pipette of 0.1 ml capacity is to be calibrated.
  - (a) Place tightly on the tip of the burette a piece of strong capillary rubber tubing extending about 2 cm beyond the burette tip.
  - (b) Make sure that the micrometer is filled almost to 10.
  - (c) Insert firmly the mouth end of the pipette and clamp it almost horizontal, with the tip end slightly higher than the mouth end. (If the pipette is clamped upright subsequent pressure of the mercury column might give erroneous readings or even spillage.)
  - (d) Manipulate micrometer until mercury column is exactly at top mark of pipette.
  - (e) Read micrometer.
  - (f) Force mercury column exactly to tip of pipette.
  - (g) Read micrometer again.
  - (h) Calculate volume of pipette as follows

$$\text{ml in pipette} = (\text{Reading 1 minus Reading 2}) \times (\text{Micrometer calibration})$$

**4 Standard Solutions****a. Definitions**

- (1) The *atomic weight* of an element is the relative weight of its atom compared to the common isotope of oxygen as 16.
- (2) The *molecular weight* of a compound is the sum of the atomic weights of all the atoms composing the molecule.
- (3) The *equivalent weight* or *combining weight* of an element or compound is that weight which will react with or release 8 units of oxygen or any equivalent thereof. Equivalent weight depends upon the specific reaction under consideration.
  - (a) The *gram equivalent* is the equivalent weight in grams.
  - (b) The *milliequivalent (mEq)* is the equivalent weight in milligrams.
- (4) A *molar solution* contains one gram molecular weight in 1 liter e.g. a molar solution of sulfuric acid contains 98.08 gm of sulfuric acid diluted to a liter of solution.
- (5) A *normal solution* contains one gram equivalent weight of a substance in a liter of solution. The gram equivalent weight depends upon the reaction for which the substance is used. E.g. a normal solution of sulfuric acid will contain  $\frac{98.08}{2} = 49.04$  gm of the acid diluted to one liter with water.
- (6) A *normal reducing solution* contains in one liter of solution one gram of oxidizable hydrogen or its equivalent in other reducing substances. Oxalic acid ( $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ ) contains two hydrogen atoms oxidizable by potassium permanganate ( $\text{KMnO}_4$ ) so that its normal reducing solution contains  $\frac{1}{2} \times 63 = 31.5$  gm per liter.
- (7) One liter of a *normal oxidizing solution* will oxidize one gram atom of hydrogen or equivalent reducing substance. A normal solution of  $\text{KMnO}_4$  for use in the reaction with oxalic acid contains  $\frac{158.03}{5} = 31.606$  gm per liter.

## 5 Preparations of Standard Solutions

### a. Hydrochloric Acid $\text{HCl}$ , mol. wt. 36.47

- (1) Preparation of 1 000 N Solution Mix 500 ml of water and 50.0 ml of concentrated hydrochloric acid (sp. gr. 1.2) At 25 °C adjust density to 1.096 with more acid or water as needed. Distill at the rate of 3-4 ml/minute and collect the last quarter of the distillate for use. Record the barometric pressure during distillation. Prepare exactly 1 000 N  $\text{HCl}$  by use of the following table.

| Barometric Pressure<br>mm Hg | $\text{HCl}$<br>Concentration<br>% by weight | Required for One<br>Liter 1 000 N $\text{HCl}$<br>gm |
|------------------------------|--|--|
| 730                          | 0.93   | 179.56   |
| 740                          | 0.69   | 170.77   |
| 750                          | 0.45   | 179.98   |
| 760                          | 0.1  | 180.19   |
| 770                          | 0.197  | 180.41   |

Prepare standard hydrochloric acid by diluting the proper weight of  $\text{HCl}$  to volume. Eg. Barometric pressure is 760. To prepare exactly 1 000 N  $\text{HCl}$  weigh out exactly 180.19 gm of distillate and dilute to 1 liter with water.

- (2) Standardization This is primary standard and needs no standardization.

### b. Potassium Hydrogen Phthalate $\text{C}_8\text{H}_5(\text{COOH})(\text{COOK})$ , mol. wt. 204.14

- (1) Preparation of 0.1000 N Solution Obtain pure dry salt from the U. S. Bureau of Standards. Dilute 0.414 gm to 1 liter with water. When titrating with this solution use phenolphthalein for indicator. Before the end point is reached boil the aliquot being titrated so as to remove dissolved  $\text{CO}_2$ .

- (2) Standardization This is a primary standard and needs no standardization.

### c. Sodium Hydroxide $\text{NaOH}$ , mol. wt. 40.01

- (1) Stock Concentrated Solution For many purposes it is desirable to keep on hand a concentrated solution. In a 500 ml Pyrex flask or bottle mix approximately 90 gm  $\text{NaOH}$  and 60 ml water. Cool and stopper with a rubber not glass stopper. This keeps indefinitely but should not be exposed to air except as necessary. Carbonates are practically insoluble at this concentration and will settle out. The supernatant fluid is approximately 1.5 N in  $\text{NaOH}$ .

- (2) Preparation of 1 000 N Sodium Hydroxide Boil approximately 1400 ml of distilled water to remove  $\text{CO}_2$  and cool to room temperature. Decant or filter through asbestos 60-80 ml of concentrated stock solution. Dilute 50 ml of this to exactly one liter with boiled distilled water.

- (3) Standardization Using phenolphthalein as indicator and standard potassium acid phthalate in the burette pipette 20 ml portions of sodium hydroxide into 50 ml Erlenmeyer flasks. Add approximately 15 ml of the standard acid bring to boil and titrate to end point. If alkali is over 1 000 N add water according to the following equation.

To 1000 ml of  $\text{NaOH}$  add  $(1000) \times (\text{actual normality} - 1.000)$

Repeat standardization until exactly 1 000 N is reached. Store in Pyrex with rubber not glass stopper. Avoid exposure to air.

### d. Oxalic Acid $\text{C}_2\text{H}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ mol. wt. 126.06

- (1) Preparation of 0.1000 N Solution Dissolve exactly 6.303 gm of pure dry oxalic acid in water and dilute to 1 liter. This keeps indefinitely.

- (2) Standardization Titrate against hot sodium hydroxide with phenolphthalein as indicator.

## 76 Instrumentation

### d. Calibration of Apparatus With the Scholander Micrometer Burette

- (1) Note that for calibration of other apparatus the temperature of the micrometer and of the apparatus calibrated must be the same, but that the absolute temperature is immaterial. This is so because the initial calibration of the micrometer burette gives calibration constant equating scale divisions of the micrometer with volume.
- (2) As an example let us assume that a serological pipette of 0.1 ml capacity is to be calibrated.
  - (a) Place tightly on the tip of the burette a piece of strong capillary rubber tubing extending about 2 cm beyond the burette tip
  - (b) Make sure that the micrometer is filled almost to 10
  - (c) Insert firmly the mouth end of the pipette and clamp it almost horizontal, with the tip end slightly higher than the mouth end. (If the pipette is clamped upright subsequent pressure of the mercury column might give erroneous readings or even spillage.)
  - (d) Manipulate micrometer until mercury column is exactly at top mark of pipette
  - (e) Read micrometer
  - (f) Force mercury column exactly to tip of pipette
  - (g) Read micrometer again
  - (h) Calculate volume of pipette as follows

$$\text{ml in pipette} = (\text{Reading 1 minus Reading 2}) \times (\text{Micrometer calibration})$$

## 4 Standard Solutions

### a. Definitions

- (1) The *atomic weight* of an element is the relative weight of its atom compared to the common isotope of oxygen as 16.
- (2) The *molecular weight* of a compound is the sum of the atomic weights of all the atoms composing the molecule.
- (3) The *equivalent weight* or *combining weight* of an element or compound is that weight which will react with or release 8 units of oxygen or any equivalent thereof. Equivalent weight depends upon the specific reaction under consideration.
  - (a) The *gram equivalent* is the equivalent weight in grams
  - (b) The *millicquivalent (mEq)* is the equivalent weight in milligrams
- (4) A *molar solution* contains one gram molecular weight in 1 liter. e.g. a molar solution of sulfuric acid contains 98.08 gm of sulfuric acid diluted to a liter of solution.
- (5) A *normal solution* contains one gram equivalent weight of a substance in a liter of solution. The gram equivalent weight depends upon the reaction for which the substance is used. E.g. a normal solution of sulfuric acid will contain  $\frac{98.08}{2} = 49.04$  gm of the acid diluted to one liter with water.
- (6) A *normal reducing solution* contains in one liter of solution one gram of oxidizable hydrogen or its equivalent in other reducing substances. Oxalic acid ( $\text{H}_2\text{C}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$ ) contains two hydrogen atoms oxidizable by potassium permanganate ( $\text{KMnO}_4$ ) so that its normal reducing solution contains  $\frac{126}{2} = 63$  gm per liter.
- (7) One liter of a *normal oxidizing solution* will oxidize one gram atom of hydrogen or equivalent reducing substance. A normal solution of  $\text{KMnO}_4$  for use in the reaction with oxalic acid contains  $\frac{158.03}{5} = 31.606$  gm per liter.

## 5 Preparations of Standard Solutions

a. Hydrochloric Acid  $\text{HCl}$ , mol. wt. 36.47

- (1) Preparation of 1 000 N Solution Mix 500 ml of water and 500 ml of concentrated hydrochloric acid (sp gr 1.2). At  $20^{\circ}\text{C}$  adjust density to 1.096 with more acid or water as needed. Distill at the rate of 3-4 ml/minute and collect the last quarter of the distillate for use. Record the barometric pressure during distillation. Prepare exactly 1 000 N  $\text{HCl}$  by use of the following table.

| Barometric Pressure<br>mm Hg | $\text{HCl}$<br>Concentration<br>% by weight | Required for One<br>Liter 1 000 N $\text{HCl}$<br>gm |
|------------------------------|--|--|
| 730                          | 0.93   | 179.56   |
| 740                          | 0.92   | 179.77   |
| 750                          | 0.91   | 179.98   |
| 760                          | 0.90   | 180.19   |
| 770                          | 0.89   | 180.41   |

Prepare standard hydrochloric acid by diluting the proper weight of  $\text{HCl}$  to volume. Eg. Barometric pressure is 760. To prepare exactly 1 000 N  $\text{HCl}$  weigh out exactly 180.19 gm of distillate and dilute to 1 liter with water.

- (2) Standardization This is primary standard and needs no standardization.

b. Potassium Hydrogen Phthalate  $\text{C}_8\text{H}_5\text{O}_4\text{K}$ , mol. wt. 204.14

- (1) Preparation of 0.1000 N Solution Obtain pure dry salt from the U. S. Bureau of Standards. Dilute 0.414 gm to 1 liter with water. When titrating with this solution use phenolphthalein for indicator. Before the end point is reached boil the aliquot being titrated so as to remove dissolved  $\text{CO}_2$ .
- (2) Standardization This is a primary standard and needs no standardization.

c. Sodium Hydroxide  $\text{NaOH}$ , mol. wt. 40.01

- (1) Stock Concentrated Solution For many purposes it is desirable to keep on hand a concentrated solution. In a 500 ml Pyrex flask or bottle mix approximately 200 gm  $\text{NaOH}$  and 300 ml water. Cool and stopper with a rubber not glass stopper. This keeps indefinitely but should not be exposed to air except as necessary. Carbonates are practically insoluble at this concentration and will settle out. The supernatant fluid is approximately 17 N in  $\text{NaOH}$ .
- (2) Preparation of 1 000 N Sodium Hydroxide Boil approximately 1400 ml of distilled water to remove  $\text{CO}_2$  and cool to room temperature. Decant or filter through asbestos 60-70 ml of concentrated stock solution. Dilute 59 ml of this to exactly one liter with boiled & distilled water.
- (3) Standardization Using phenolphthalein as indicator and standard potassium acid phthalate in the burette pipette 0.1 ml portions of sodium hydroxide into 50 ml Erlenmeyer flasks. Add approximately 15 ml of the standard acid bring to boil and titrate to end point. If alkali is over 1 000 N add water according to the following equation:
- $$\text{To 1000 ml of NaOH add } (1000) \times (\text{actual normality} - 1.000)$$
- Stop at standardization until exactly 1 000 N is reached. Store in Pyrex with rubber not glass stopper. Avoid exposure to air.

d. Oxalic Acid,  $\text{C}_2\text{H}_2\text{O}_4 \cdot 2\text{H}_2\text{O}$  mol. wt. 126.06

- (1) Preparation of 0.1000 N Solution Dissolve exactly 6.303 gm of pure dry oxalic acid in water and dilute to 1 liter. This keeps indefinitely.
- (2) Standardization Titrate against hot sodium hydroxide with phenolphthalein as indicator.

## 78 Instrumentation

### e Sulfuric Acid, $\text{H}_2\text{SO}_4$ mol. wt. 98.08

- (1) Preparation of 1 N Solution Add about 500 ml water to a 1 liter volumetric flask Add approximately 29 ml of concentration sulfuric acid (*sp gr 1.84*) Dilute to 1 liter with water
- (2) Standardization Titrate against standard sodium hydroxide The concentration should be near 1 N

### f Benzoic Acid $\text{C}_6\text{H}_5\text{COOH}$ mol wt 122.05

- (1) Preparation of 0.1000 N Solution Weigh out exactly 12.205 gm dry acid into a one liter volumetric flask Add about 500 ml warm 95% ethyl alcohol, dissolve and dilute to 1 liter
- (2) Standardization Titrate against hot sodium hydroxide using phenolphthalein as indicator

### g Silver Nitrate $\text{AgNO}_3$ mol wt. 169.89

- (1) Preparation of 0.1000 N Solution Dissolve 16.989 gm pure fused silver nitrate in water add about 100 ml of concentrated nitric acid and dilute to 1 liter Keep in the dark in a brown bottle
- (2) Standardization This may be used as a primary standard without standardization However if standardization is desired a gravimetric procedure is recommended as follows
  - (a) Make an approximately 1% solution of HCl by diluting 100 ml of concentrated HCl with 900 ml of water
  - (b) Running triplicates prepare 250 ml beakers with exactly 50 ml of silver nitrate solution
  - (c) With stirring add 75 ml of HCl small portions at a time
  - (d) Bring to boil and allow to cool and settle in the dark
  - (e) Saving all washings filter through weighed sintered glass filter Wash with 5 successive 50 ml portions of water containing a few drops of nitric acid and test the last washing with silver nitrate for presence of chloride If negative wash with two 50 ml portions of plain water
  - (f) Record volume and temperature of total washings
  - (g) Dry precipitate to constant weight at 110 C
  - (h) Calculation

$$\text{Normality of AgNO}_3 = \frac{(\text{gm ppt}) + \left( \frac{\text{total ml washings}}{1000} \times \text{loss/liter} \right)}{\left( \frac{\text{ml AgNO}_3 \text{ used}}{1000} \right) \times (143.34)}$$

- (i) Solubility of  $\text{AgCl}$  in 1% HCl is
  - 0.000 gm  $\text{AgCl}/1$  at 20 C and
  - 0.0027 gm  $\text{AgCl}/1$  at 30 C

### h Potassium Thiocyanate $\text{KCNS}$ mol wt 97.16

- (1) Preparation of 0.10 N Solution Dissolve 9.8 gm KCNS in water and dilute to 1 liter
- (2) Standardization To 20 ml of 0.1000 N  $\text{AgNO}_3$  add 5 ml concentrated  $\text{HNO}_3$  5 ml of water and 1 ml of saturated aqueous ferric alum solution Titrate with KCNS to the first permanent salmon pink color

### i Sodium Oxalate $\text{Na}_2\text{C}_2\text{O}_4$ mol wt. 133.99

- (1) Preparation of 0.1000 N Solution Weigh out exactly 6.700 gm of pure dry salt (Bureau of Standards) Dissolve and dilute to 1 liter with water
- (2) Standardization This is a primary standard and requires no standardization.

**j. Potassium Permanganate  $\text{KMnO}_4$ , mol. wt. 158.03**

- (1) Preparation of 0.1 N Solution Dissolve approximately 3.3 gm pure  $\text{KMnO}_4$  in 1 liter of water. This must age. Let stand in dark several months and filter through asbestos.
- (2) Standardization To 20 ml of 0.1000 N sodium oxalate in a 500 ml Erlenmeyer flask add 60 ml of water and 3 ml of concentrated  $\text{H}_2\text{SO}_4$ . Heat to  $60^\circ\text{C}$  and titrate to first pink persisting 30 seconds. Be sure to run adequate numbers of blanks.

**k. Potassium Iodate  $\text{KIO}_3$ , mol. wt. 214.03**

- (1) Preparation of 0.1000 N Solution Weigh out exactly 3.06 gm pure dry  $\text{KIO}_3$  dissolve in water and dilute to 1 liter.
- (2) Standardization This is a primary standard and need no standardization.

**l. Potassium Dichromate  $\text{K}_2\text{Cr}_2\text{O}_7$ , mol. wt. 294.18**

- (1) Preparation of 0.1000 N Solution Dissolve exactly 4.904 gm of pure salt in water and dilute to exactly 1 liter.
- (2) Standardization This may be used as a primary standard.

**m. Sodium Thiosulfate  $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ , mol. wt. 248.19**

- (1) Preparation of 0.1000 N Solution Dissolve 4.819 gm of pure salt in approximately 100 ml water and add 1 ml of approximately 1 N  $\text{NaOH}$ . Dilute to 1 liter with water.
- (2) Standardization against Iodate To 5 ml of standard 0.1000 N iodate in an Erlenmeyer flask add 10 ml of 10%  $\text{KI}$  and 0 ml of 1 N  $\text{H}_2\text{SO}_4$  or  $\text{HCl}$ . Titrate with thiosulfate until the solution is pale yellow. Add a few drops of a 1% aqueous starch solution (1%) and titrate to disappearance of blue color.
- (3) Standardization against Permanganate Dissolve approximately 1 gm of  $\text{KI}$  in 10 ml of water. Add approximately 3 ml of concentrated  $\text{HCl}$  and pipette exactly 5 ml of 0.1000 N  $\text{KMnO}_4$ . Titrate with thiosulfate until bright yellow, then add  $\frac{1}{2}$  to 1 ml of 1% starch and titrate to disappearance of the blue color.
- (4) Standardization with Dichromate In a 60 ml Erlenmeyer flask mix approximately 300 ml of water, 5 gm of  $\text{KI}$ , 4 gm of  $\text{NaHCO}_3$  and then excess of concentrated  $\text{HCl}$ . Add exactly 10 ml of 0.1000 potassium dichromate and let stand in the dark 10 minutes. Titrate with thiosulfate using  $\frac{1}{2}$  to 1 ml of 1% starch as indicator.

**n. Iodine  $\text{I}_2$ , mol. wt. 253.81**

- (1) Preparation of 0.1 N Solution In approximately 60 ml of water dissolve 30 gm of pure potassium iodide and 13 gm of iodine (resublimed). Dilute with water to 1 liter.
- (2) Standardization In a 100 ml Erlenmeyer flask mix exactly 10 ml of iodine and approximately 100 ml of water. Titrate with 0.1000 N thiosulfate to faint yellow. Add approximately 1 ml of 1% starch solution and titrate to disappearance of blue color.

**J INSTRUMENTATION MISCELLANEOUS REFERENCES****Balance**

Lowry O H. A Simple Quartz Torsion Balance. *J Biol Chem* 152: 229, 24 (Feb) 1944

**Burette**

Longwell R and Hill R M. A Modified Redox Burette for Titrating Solutions Which React With Mercury. *J Biol Chem* 112: 313-314 (Dec) 1936

Scholander I F, Edwards G A and Irving L. Improved Microtiter. *Bull J Biol Chem* 148: 493-500 (June) 1943

Scholander I F. Microburette. *Science* 95: 171-18 (Feb) 1944

Pehberg I B. A Method of Microtitration. *Biochem J* 19: 1-6 (Feb) 1925



## 80 Instrumentation

### Electron Microscope

Wyckoff, R W G *Electron Microscopy Technique and Applications* New York 1949  
Interscience Publishers, p 248

### Electrophoresis

Abramson H A Moyer L S and Gorin M H *Electrophoresis of Proteins* New York  
1949 Reinhold Publishing Corp

Alberty R A *An Introduction to Electrophoresis I. Methods and Calculations* J Chem.  
Ed 25 496 433 (Aug) 1948

Blux G and Pedersen K O *Electrophoresis and Ultracentrifugation of Lipide Free Hu  
man Serum* Acta Chem Scandinav 1 (5) 511 590 1947

Jager B V Schwartz, T B Smith E L, Nickerson M and Brown D M *Comparative  
Electrophoretic and Chemical Estimation of Human Serum Albumin an Evaluation  
of Six Methods* J Lab & Clin Med. 35 76 86 (Jan) 1950

Labhart H and Staub H *Mikro Elektrophorese*, Helv Chim Acta 30 ( ) 1954 1964  
1947

Lang E R Van Winkle Q and France W G *Electrophoresis in an Ultracentrifugal  
Field* J Colloid Sci 2 315 332 (June) 1947

Longworth L G *A Modification of the Schlieren Method for Use in Electrophoretic  
Analysis* J Am Chem Soc 61 549 530 (Feb) 1939

Longworth L G *The Quantitative Interpretation of the Electrophoretic Pattern of  
Proteins* J Physical & Colloid Chem. 51 171 183 (Jan) 1947

Moore D H and White J U *A New Compact Tiselius Electrophoresis Apparatus* Rev  
Sci Instr 19 700 706 (Oct) 1948

Popyak G and McCarthy, E I *Osmotic Pressure of Experimental and Human Lipaemia  
Sera Evaluation of Albumin Globulin Ratios With the Aid of Electrophoresis*  
Biochem J 40 789 803 1946

Stern K G *Amnco Stern Electrophoresis Apparatus* American Instrument Company  
Silver Springs Md Bull No 4175 (June) 1949

### Flame Photometer

Domingo W R, and Klyne W *Photoelectric Flame Photometer* Biochem J 45 400  
408 1949

### Gas Analysis

Scholander P F and Evans H J *Microanalysis of Fractions of a Cubic Millimeter of  
Gas* J Biol Chem 169 (3) 551 560 (Aug) 1947

Scholander P F and Irving L *Micro Blood Gas Analysis in Fractions of a Cubic  
Millimeter of Blood* J Biol Chem 169 (3) 561 569 (Aug) 1947

### Infrared Spectrometer

Gore R C *Infrared Spectrometry of Small Samples With the Reflecting Microscope*  
Science 110 710 711 (Dec) 1949

### Manometric Techniques

Dixon M *Manometric Methods as Applied to the Measurement of Cell Respiration and  
Other Processes* ed 2 London 1943 Cambridge University Press

Holmes E *The Metabolism of Living Tissues* New York 1937, The Macmillan Company

Laser H and Rothschild L *A New Manometric Method for Determination of Respira  
tory Quotients* Biochem J 45 598 612 1949

Umbreit W W Burris R T and Stauffer J F *Manometric Techniques and Related  
Methods for the Study of Tissue Metabolism* 1st ed., 4th Printing Minneapolis  
(Feb) 1948 Burgess Publishing Company

### Mass Spectrometer

Hipple J A and Shephard M *Mass Spectrometry* Anal Chem 21 3 36 (Jan.) 1949

Kerwin L *New Type Mass Spectrometer* Rev Scient Instruments 21 96 97 (Jan) 1950

Near A O *Mass Spectrometer* Rev Scient Instruments 18 398 1947

Washburn H W *Mass Spectrometer* U S Nav Bull Suppl (March April) 1948.

Wiley H F *A Mass Spectrometer for Isotope Ratio Determinations* Science 110 51  
355 (Oct) 1949

### Microrespirometer

Scholander P F *Volumetric Microrespirometers* Rev Sc Instr 13 3 33 (Jan) 194

### Oximeter

Milikan G A *The Oximeter an Instrument for Measuring Continuously the Oxygen  
Saturation of Arterial Blood in Man* Rev Scient Instruments 13 (10) 434 444 (Oct)  
1942

- Hemingway A and Taylor C B Laboratory Tests of the Oximeter With Automatic Compensation for Vasomotor Changes *J Lab & Clin Med.* 29 987-991 (Sept.) 1944  
 Comroe J H Jr., and Walker P Normal Human Arterial Oxygen Saturation Determined by Equilibration With 100% O<sub>2</sub> in Vivo and by the Oximeter *Am J Physiol* 152 365-371 (Feb) 1948

#### Pipettes

- Shohl A. T A Pipet for Microanalyses *J Am Chem. Soc* 60 417 (Feb) 1938

#### Polarography

- Kolthoff I M and Lingane J J Polarography (Revised Reprint), New York 1948 Interscience Publishers

#### Respirometer

- Scholarer P F Volumetric Plastic Micro Respirometer, *Rev Scient Instruments* 21 378-380 (Apr) 1950

#### Ultracentrifuge

- Pedersen K O Ultracentrifugal Studies on Serum and Serum Fractions *Uppsala Sweden* 1945, Almquist and Wiksells  
 Pickels E G The Ultracentrifuge Practical Aspects of the Ultracentrifugal Analysis of Proteins *Chem Rev* 30 341-355 (June) 1947  
 Svedberg T., and Pedersen K O The Ultracentrifuge, London 1940 Oxford University Press

#### Ultrafiltration

- Lavielles P H Anaerobic Ultrafiltration *J Biol Chem* 100 26-45 (Aug) 1937

#### Ultramicroanalysis

- Kirk P L Quantitative Ultramicroanalysis, New York 1950 John Wiley & Sons Inc.

## 80 Instrumentation

### Electron Microscope

Wyckoff R W G Electron Microscopy Technique and Applications, New York, 1949  
Interscience Publishers p 248

### Electrophoresis

Abramson H A, Mover L S and Gorin M H Electrophoresis of Proteins New York  
1942 Reinhold Publishing Corp

Alberty R A An Introduction to Electrophoresis I Methods and Calculations J Chem  
Ed 25 466-439 (Aug) 1948

Blax G and Peder sen K O Electrophoresis and Ultracentrifugation of Lipids Free F.  
man Serum Acta Chem. Scandinav 1 (5), 511-520 1947

Jager B V, Schwartz T B, Smith E L, Nickerson M and Brown D M Comparative  
Electrophoretic and Chemical Estimation of Human Serum Albumin: an Evaluation  
of Six Methods J Lab & Clin Med 35 78-86 (Jan) 1950

Labhart H, and Staub H Mikro-Elektrophorese, Helv Chim Acta 30 (7) 1941-1951  
1947

Lang E P, Van Winkle Q and France W G Electrophoresis in an Ultracentrifugal  
Field J Colloid Sci 2 310-337 (June) 1947

Longworth J G A Modification of the Schlieren Method for Use in Electrophoretic  
Analysis J Am Chem Soc 61 59-30 (Feb) 1939

Longworth J G The Quantitative Interpretation of the Electrophoretic Pattern of  
Proteins J Physica & Colloid Chem 51 171-183 (Jan) 1946

Moore D H, and White J U A New Compact Tiselius Electrophoresis Apparatus Rev  
Sci Instr 19 1007-06 (Oct) 1948

Popjak G and McCarthy E F Osmotic Pressure of Experimental and Human Lipoprotein  
Sera: Evaluation of Albumin-Globulin Ratios With the Aid of Electrophoresis,  
Biochem J 40 789-803 1946

Stern K G Aminco Stern Electrophoresis Apparatus American Instrument Company  
Silver Springs Md, Bull No 2175 (June) 1949

### Flame Photometer

Domingo, W R and Klyne W Photoelectric Flame Photometer Biochem J 45 400  
408 1949

### Gas Analysis

Scholander P F and Evans H I Microanalysis of Fractions of a Cubic Millimeter of  
Gas J Biol Chem 169 (3) 501-560 (Aug) 1947

Scholander, P F and Irving L Micro Blood Gas Analysis in Fractions of a Cubic  
Millimeter of Blood, J Biol Chem 169 (3) 561-569 (Aug) 1947

### Infrared Spectrometer

Gore, R C Infrared Spectrometry of Small Samples With the Reflecting Microscope  
Science 110 710-711 (Dec) 1949

### Manometric Techniques

Dixon M Manometric Methods as Applied to the Measurement of Cell Respiration and  
Other Processes ed. London 1943, Cambridge University Press

Holmes F The Metabolism of Living Tissues New York 1937 The Macmillan Company

Laser, H and Rothchild L A New Manometric Method for Determination of Respiratory  
Quotients Biochem J 45 593-612 1949

Umbreit, W W, Burris R T, and Stauffer J F Manometric Techniques and Related  
Methods for the Study of Tissue Metabolism 1st ed., 4th Printing, Minneapolis,  
(Feb) 1949 Burgess Publishing Company

### Mass Spectrometer

Hipple J A and Shephard M Mass Spectrometry Anal Chem 21 37-36 (Jan) 1949

Kerwin J A New Type Mass Spectrometer Rev Scient Instruments 21 96-97 (Jan) 1950

Near A O Mass Spectrometer Rev Scient Instruments 18 308 1947

Washburn H W Mass Spectrometer, U S Nav Bull Suppl (March-April) 1948

Wile, H F A Mass Spectrometer for Isotope Ratio Determination Science 110 334  
335 (Oct) 1949

### Microrespirometer

Scholander P F Volumetric Microrespirometers Rev Sci Instr 13 32-33 (Jan) 1944

### Oximeter

Milikan G A The Oximeter an Instrument for Measuring Continuously the Oxygen  
Saturation of Arterial Blood in Man Rev Scient Instruments 13 (10) 404-444 (Oct)  
1942

- Hemingway A., and Taylor G B Laboratory Tests of the Oximeter With Automatic Compensation for Vasomotor Changes *J Lab & Clin. Med.* 29 987 991 (Sept) 1944  
 Cumro, J H Jr, and Walker P Normal Human Arterial Oxygen Saturation Determined by Equilibration With 100% O<sub>2</sub> in Vivo and by the Oximeter *Am J Physiol* 152 365 371 (Feb) 1948

#### Pipettes

- Shohl A. T A Pipet for Microanalyses *J Am. Chem. Soc* 60 417 (Feb) 1938

#### Polarography

- Kolthoff I M and Lingane J J Polarography (Revised Edition), New York 1946, Interscience Publishers

#### Respirometer

- Scholander P F Volumetric Plastic Micro Respirometer *Rev Scient Instruments* 21 783 0 (Apr) 1950

#### Ultracentrifuge

- Pedersen K O Ultracentrifugal Studies on Serum and Serum Fractions Upsala Sweden, 1945 Almqvist and Wiksells  
 Pickels E O The Ultracentrifuge Practical Aspects of the Ultracentrifugal Analysis of Proteins *Chem Rev* 36 341 355 (June) 1946  
 Svedberg T and Pedersen K O The Ultracentrifuge London 1940 Oxford University Press

#### Ultrafiltration

- Lavietes P H Anaerobic Ultrafiltration *J Biol Chem* 120 267 275 (Aug) 1937

#### Ultramicroanalysis

- Kirk P L Quantitative Ultramicroanalysis New York 1940 John Wiley & Sons Inc.

# SECTION IV

## BIOCHEMICAL PROCEDURES

---

### A MINERALS

#### 1 Total Water and Total Solids in Food, Tissue, Cells, Serum and Feces

##### Reference

Association of Official Agricultural Chemists ' Official and Tentative Methods of Analysis Total Solids Method I Official ' ed. 5, Washington D C 1940 p 71

##### Principle

The difference in weight between a sample before and after drying at 110 C overnight is the amount of water in the sample The residue represents the total solids of the sample

##### Apparatus

- 1 A number of weighing bottles with lids
- 2 An analytical balance
- 3 An oven maintained at 110 C
- 4 A desiccator

##### Reagents

None

##### Procedure

- 1 Weigh accurately a dry weighing bottle that has been cooled in a desiccator for at least one half hour
- 2 Add the sample and weigh again.
- 3 Open the dish and leave it in a 110 C oven overnight
- 4 Cool in the desiccator close and weigh again.

##### Calculation

$$1 \text{ gm water } \% = \frac{(\text{Weight of bottle before drying} - \text{weight of bottle after})}{\text{gm sample taken}} \times 100$$

$$2 \text{ gm solids } \% = \frac{(\text{Weight of bottle after drying} - \text{initial weight of bottle})}{\text{gm sample taken}} \times 100$$

##### Example

|                              |            |
|------------------------------|------------|
| Weight of weighing bottle    | 10.9586 gm |
| Weight with sample of muscle | 12.5542 gm |
| Weight after drying          | 11.3507 gm |
| Therefore, weight of sample  | 1.5956 gm  |
| weight of residue            | 0.3921 gm  |
| and weight of water          | 1.2035 gm  |

$$\text{gm water } \% = \frac{1.2035}{1.5956} \times 100 = 75.4$$

$$\text{gm solids } \% = \frac{0.3921}{1.5956} \times 100 = 24.6$$

## Precautions

- 1 In the presence of unsaturated fat in large concentration the weight of the residue of a specimen sometimes reaches a minimum and then rises again as oxygen is absorbed. It is desirable to make weighings as often as every six hours to determine the minimum.
- 2 In the presence of a high concentration of water preliminary drying on a steam bath is necessary to prevent frothing which often occurs when a sample is placed directly into an oven.


## 2 Total Base in Serum, Urine Cells and Tissue

### Reference

Consolazio W V and Talbott J H The Determination of Total Base in Biological Material by Electrodialysis J Biol Chem. 132 753 67 (Feb) 1940

### Principle

The basic ions are electro dialysed into mercury with which they amalgamate. The amalgam is decomposed with standard acid and the excess acid back titrated with standard alkali.



### Apparatus

- 1 Total base machine with accessories (Macalaster Bicknell Company Cambridge Massachusetts). An individual unit is illustrated in Fig 15. This includes two millimeters a pilot lamp a rheostat a shunt and a switch. It is operated at 110 or 20 volts DC the possibility of using lower voltages has not been tested fully. It draws 65 milliamperes at full load. The anode is a platinum electrode fused in a Pyrex cup and connected by a metal socket directly to the source of power. The anode acid conductor of Adair and Keys (Adair G S and Keys A B J Physiol 81 16 1934) has been eliminated. Anode cups should not be cleansed with cleaning solution. Precipitation of chromium on platinum may introduce serious errors. It is recommended that the cups be brushed with water filled with concentrated nitric acid and allowed to stand for several hours. Several rinses with distilled water remove all adherent base.
- 2 #300 plain transparent cellophane (Dupont)
- 3 Rehberg microburette 0 ml
- 4 Slow speed rotating motor for applying collodion
- 5 Membrane tubes. These are prepared as follows. Cut cellophane into 3 inch squares and soak for 5 minutes in distilled water. Apply the wet squares to the open end and opposite the constriction. Now with the thumb and the forefinger pull the wet membrane taut by applying pressure along the tube wall and away from the former opening. The membrane should always be stretched across the grain since this axis has greatest elasticity. A good procedure to follow is to stretch the membrane slightly cross grain before applying it to the membrane tube. Allow the membrane to dry 1<sup>1</sup>/<sub>2</sub> hours. The excess cellophane is now trimmed with a sharp razor or scalpel about 1/4 inch down the tube wall. The partially made membrane tubes are now rotated parallel to the floor about 20-30 rpm and while they are rotating a thin coat of collodion is applied. Starting at the edge of the tube pour a continuous stream about 1/4 inch from the tube end. The collodion sets in a minute or two. Care should be taken that collodion does not flow over the membrane face. After drying (3-5 hours) the finished membrane tubes are rinsed thoroughly with distilled water and stored in distilled water in the refrigerator. Membrane tubes prepared in this fashion have

## ELECTRODIALYSIS UNIT FOR TOTAL BASE




Fig 16

been used successfully two months later. To remove the membranes from dried membrane tubes the collodion seal should be soaked at least 15 minutes in water. It then merely becomes necessary to run the fingernail down the collodion seal and remove the old membrane in one piece.

6 Krogh syringe pipettes 0.2 ml capacity and 17 ml capacity

### Reagents

- 1 Water glass distilled over alkaline permanganate. It must meet the specifications for conductivity water.
- 2 Mercury base free. Prepare this by cleaning with 5% nitric acid with subsequent thorough washing with distilled water.
- 3 Stock solution of benzoic acid ( $\text{C}_6\text{H}_5\text{COOH}$ ) 1.000 N. Dry benzoic acid by placing it in a vacuum desiccator at least 48 hours. Accurately weigh 122.047 gm. dissolve in 100 ml absolute EtOH and make up to liter with absolute EtOH. If stored in the cold it can be kept indefinitely.
- 4 Benzoic acid 0.0 N made from stock solution. Pipette exactly 5 ml of 1.000 N benzoic acid into a 50 ml volumetric flask. Dilute to about 40 ml with glass distilled water with persistent agitation reprecipitated benzoic acid will dissolve. Dilute to mark.
- 5 Sodium hydroxide 0.200 N (Kahlbaum CO free).
- 6 Phenolphthalein 0.10% (Lamotte) dilute 1:1 with EtOH.
- 7 Collodion (Mallinckrodt) CP.
- 8 Standard solution of mixed ions. Dissolve and dilute to exactly one liter with glass distilled water: 8.190 gm NaCl, 0.73 gm KCl, 0.018 gm  $\text{CaCl}_2 \cdot \text{H}_2\text{O}$  and 0.095 gm  $\text{MgCl}_2$ . The resulting solution contains per liter 140 mEq Na, 5 mEq K, 5 mEq  $\text{Ca}^{++}$  and 5 mEq  $\text{Mg}^{++}$  making a total base of 150 mEq per liter.
- 9 Concentrated nitric acid.

### Procedure

- 1 Wash anode cups with soap. Dreft. Rinse and leave filled with concentrated  $\text{HNO}_3$ . Rinse anode cups with distilled water several times then once with glass distilled water.
- 2 From a burette add to each cup approximately 0.5 ml mercury and with a broken tip pipette approximately 10 ml of glass distilled water.
- 3 Add 0.2 ml serum, cells or urine or weighed tissue approximating 0.1 gm (the quantity used depends upon the base content). Rinse several times with water in anode cup keeping the tip of the Krogh (automatic) pipette below the surface of the water to prevent splashing. Then flush out pipette with 2 ml of water. With each run have 0.2 ml of standard total base solution which has to be used for calculating the blank.
- 4 Always rinse membrane tubes several times with distilled water and once with glass distilled water before use. To each membrane tube add approximately 1 ml of mercury.
- 5 Add exactly a measured 17 ml of 0.00 N benzoic acid. The amount added depends on base unknown. (17 ml is equivalent to 170 mEq of base/liter of unknown when 0.2 ml of unknown is used.)
- 6 Set membrane tubes in solution in anode cups submerged about 0.5 cm.
- 7 Close circuit by attaching all electrode contacts and turn on current. After letting it run low for 1 hour advance the rheostat and let it run overnight on high. Dialysis is complete when milliammeter indicates no current or not more than 3 milliamperes for all 9 units.
- 8 Remove contacts and lift membrane tube above solution in cups.
- 9 Rinse off electrodes in membrane tubes with 0.2 ml glass distilled water.
- 10 Add 1 drop of phenolphthalein.
- 11 Titrate to first sign of pink color with 0.000 N NaOH using Rehberg burette.



### Calculation

1 Blank

In this method the reagent blank is best computed from a standard solution not from a run of distilled water.

$$\text{mEq blank} = (\text{mEq acid} - \text{mEq alkali}) - (\text{ml standard} \times \text{mEq total base per ml standard})$$

2 Unknown

mEq total base/liter unknown =

$$(\text{mEq acid} - \text{mEq alkali} - \text{mEq blank}) \times \frac{1000}{\text{ml unknown}}$$

### Example

17 ml acid used titration of 0.2 ml of unknown is 0.01 ml NaOH titration of 0.2 ml of standard is 0.0034 ml NaOH

$$1 \text{ Blank} = (17 \times 02) - (0.004 \times 0^{\circ}) - (0.152 \times 0^{\circ}) = 0.003$$

$$2 \text{ Unknown} = (17 \times 02) - (0.021 \times 03) - (0.003) \times \frac{1000}{0.2} = 135 \text{ mEq/liter}$$

3 This computation is made more expeditious by the use of Table 3

TABLE 3 TOTAL BASE

Millicivalents Total Ba s per Liter of Serum Whole Blood Cells or Urine When 0.6 ml  
Samples Are Analyzed When 17 ml of 0.02 N Benzoic Acid Is Used in the Cathode  
Cups During Electrolysis and When A Means  
Millicivalents (NaOH Plus Blank) in Back Titration

[illegible]

### Precautions

- 1 In order to avoid erratic results only glass distilled water should be used.
- 2 Water gives irregular blanks therefore use one standard during each group of analyses to obtain the blank correction.

## 3 Sodium in Serum, Urine Cells Tissues Diets and Stools (Gravimetric Method)


### A. IN SERUM

#### References

- Butler A. M. and Tothill E. An Application of the Zinc Uranyl Acetate Method for Determination of Sodium in Biological Material *J Biol Chem* 93 171 180 (Sept) 1931 193
- Consolazione W. V. and Dill D. B. The Determination of Sodium *J Biol Chem* 197 587 59 (Feb) 1941

#### Principle

The ash from biological material is treated with uranyl zinc acetate and the triple salt of sodium uranyl zinc acetate is precipitated quantitatively. The precipitate is estimated gravimetrically.



#### Apparatus

- 1 Jena or Pyrex sintered glass filters (fine)
- 2 Platinum or quartz crucibles 5 ml capacity
- 3 Small electric stirring motor with platinum wire
- 4 Muffle oven set at 600 C
- 5 Electric centrifuge or good analytical filter paper (Whatman #40) and funnels
- 6 Syringe pipettes 10 ml and 5 ml
- 7 Analytical balance
- 8 An oven set at 110 C
- 9 Rubber stoppers #6 solid
- 10 Watch glasses
- 11 Source of suction

#### Reagents

- 1 Uranyl zinc acetate (UZA)
  - a. 80 gm of sodium free uranyl acetate  $\text{UO}_2(\text{C}_2\text{H}_3\text{O}_2)_2 \cdot 2\text{H}_2\text{O}$   
45 ml of acetic acid ( $\text{CH}_3\text{COOH}$ ) 30% by volume  
Water to make up to 500 gm
  - b. 50 gm of zinc acetate  $\text{Zn}(\text{CH}_3\text{COO})_2 \cdot 2\text{H}_2\text{O}$   
23 ml of acetic acid 30% by volume  
Water to make up to 500 gm

Warm both solutions on steam bath. Mix while hot let stand ~4 hours and seed with sodium uranyl zinc acetate. Stir well and filter before using.

Ethyl alcohol (EtOH) 95% saturated with sodium uranyl zinc acetate. Filter before use.

- 3 Sulfuric acid concentrated
- 4 Ethyl ether anhydrous ( $\text{C}_2\text{H}_6\text{O}$ )
- 5 Calcium hydroxide powdered (sodium free)
- 6 Hydrochloric acid one N (10 ml concentrated HCl in 100 ml of water)
- 7 Sodium uranyl zinc acetate (triple sodium salt)

**Procedure**

- 1 Pipette exactly 1 ml of serum into a platinum or quartz crucible
- 2 Add one drop of concentrated sulfuric acid and char at 110 C for at least 15 hours
- 3 Ash in muffle oven (600 C) overnight
- 4 To the bottom of an accurately weighed sintered glass filter fit a #6 rubber stopper
- 5 Add approximately 15 ml of freshly filtered UZA.
- 6 To the ashed serum add 5 to 7 drops of one N hydrochloric acid.
- 7 By means of a medicine dropper transfer the sample quantitatively to the weighed crucible while mechanically stirring the UZA
- 8 Wash remaining serum ash with three 0.2 ml portions of approximately 0.1 N HCl transferring it to the crucible
- 9 Now wash the crucible with 2 ml of UZA.
- 10 Stir mechanically for a total of 5 minutes. Then raise the stirring rod with 2 ml of UZA cover the crucible with a watch glass and let stand for 30 minutes at room temperature
- 11 Filter the ingredients in crucible by suction
- 12 Using capillary pipettes, wash 5 times with 2 ml portions of freshly filtered 95% alcohol saturated with the triple salt of UZA.
- 13 Then wash 5 times with 2 ml portions of absolute ether
- 14 Dry by vacuum. Clean outside of crucible with a clean cloth and place in a desiccator for 1 hour
- 15 Weigh

**Calculation**

- 1 If 1.0 ml of serum was used  $\text{mEq Na/liter} = (\text{weight of precipitate in grams} - \text{blank in grams}) \times 650$
- 2 To convert mEq Na to mg Na, multiply by 23

**Example**

- 1 Weight of precipitate was 0.2220 gm and weight of blank was 0.0010 gm,  $\text{mEq Na/liter serum} = (0.2220 - 0.0010) \times 650 = 143.0 \text{ mEq/l}$
- 2  $\text{mg Na} = 143.0 \times 23$
- 3 This calculation is assisted by Table 4.

**Precautions**

- 1 In washing out the precipitate be sure that no water is added as the sodium zinc uranyl acetate precipitate dissolves very readily

TABLE 4 SERUM SODIUM  
When 1 ml of serum is analyzed,  $\text{mEq Na/liter} =$   
 $(\text{gm ppt} - \text{gm blank}) \times 650$  or  
 $(\text{mg ppt} - \text{mg blank}) \div 0.65$

| MG PPT<br>MINUS<br>MG<br>BLANK          | 0   | 1   | 2   | 3   | 4   | 5   | 6   | 7   | 8   | 9   |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| MILLIEQUIVALENTS SODIUM PER LITER SERUM |     |     |     |     |     |     |     |     |     |     |
| 190                                     | 1.4 | 124 | 1.5 | 126 | 1.6 | 128 | 1.7 | 130 | 1.8 | 132 |
| 200                                     | 1.5 | 131 | 131 | 132 | 133 | 134 | 135 | 136 | 137 | 138 |
| 210                                     | 137 | 138 | 139 | 140 | 141 | 142 | 143 | 144 | 145 | 146 |
| 220                                     | 143 | 144 | 145 | 146 | 147 | 148 | 149 | 150 | 151 | 152 |
| 230                                     | 150 | 151 | 152 | 153 | 154 | 155 | 156 | 157 | 158 | 159 |
| 240                                     | 156 | 157 | 158 | 159 | 160 | 161 | 162 | 163 | 164 | 165 |
| 250                                     | 163 | 164 | 165 | 166 | 167 | 168 | 169 | 170 | 171 | 172 |

- 2 Make sure that the filters are dry before weighing
- 3 The accuracy of the method must be checked against a standard solution (see total base)

## B SODIUM IN URINE

Apparatus See serum

Reagents See serum

### Procedure

- 1 To approximately 10 ml of urine add approximately 1 gm of powdered Ca(OH)
- 2 Shake manually at intervals for 1 hour
- 3 Centrifuge or filter through hardened filter paper

TABLE 5 SODIUM IN URINE

When 2 ml of urine are used for analysis mEq Na/liter urine =

(gm ppt - gm blank)  $\times$  3<sup>95</sup> or

(mg ppt - mg blank)  $\times$  0.35

| MG PPT<br>MINUS<br>MG BLANK | 0                                       | 1     | 2     | 3     | 4     | 5     | 6     | 7     | 8     | 9     |
|-----------------------------|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|
|                             | MILLIEQUIVALENTS SODIUM PER LITER URINE |       |       |       |       |       |       |       |       |       |
| 0                           | 0.0                                     | 0.3   | 0.7   | 1.0   | 1.3   | 1.6   | 0     | 3     | 2.6   | 9     |
| 10                          | 3.3                                     | 3.6   | 4.0   | 4.2   | 4.6   | 4.9   | 5.3   | 5.6   | 5.9   | 6     |
| 20                          | 6.5                                     | 6.8   | 7.2   | 7.5   | 7.8   | 8.1   | 8.5   | 8.8   | 9.1   | 9.4   |
| 30                          | 9.8                                     | 10.1  | 10.4  | 10.7  | 11.1  | 11.4  | 11.7  | 1.0   | 1.4   | 12.7  |
| 40                          | 13.0                                    | 13.3  | 13.7  | 14.0  | 14.3  | 14.6  | 1.0   | 1.3   | 1.6   | 15.9  |
| 50                          | 16.3                                    | 16.6  | 16.9  | 17    | 17.6  | 17.9  | 18    | 18.5  | 18.9  | 19.2  |
| 60                          | 19.5                                    | 19.8  | 0.2   | 0.5   | 0.8   | 21.1  | 15    | 21.8  | 2.1   | 4     |
| 70                          | 2.8                                     | 23.1  | 3.4   | 3.7   | 4.1   | 4.4   | 24.7  | 2.0   | 3.4   | 25.7  |
| 80                          | 26.0                                    | 26.3  | 26.7  | 27.0  | 27.3  | 7.6   | 28.0  | 8.3   | 28.6  | 8.9   |
| 90                          | 9.3                                     | 9.6   | 9.9   | 30.2  | 30.6  | 0.9   | 31.2  | 31.5  | 31.9  | 32    |
| 100                         | 3.5                                     | 30.8  | 3.3   | 33.5  | 33.8  | 34.1  | 34.5  | 34.8  | 3.1   | 35.4  |
| 110                         | 35.8                                    | 36.1  | 36.4  | 36.7  | 37.1  | 37.4  | 37.7  | 38.0  | 38.4  | 38.7  |
| 120                         | 39.0                                    | 9.3   | 9.7   | 40.0  | 40.3  | 40.6  | 41.0  | 41.3  | 41.6  | 41.9  |
| 130                         | 4.3                                     | 42.6  | 4.9   | 43.2  | 43.6  | 43.9  | 44.2  | 44.5  | 44.9  | 45    |
| 140                         | 45.5                                    | 45.8  | 46    | 46.5  | 46.8  | 47.1  | 47.5  | 47.8  | 48.1  | 48.4  |
| 150                         | 48.8                                    | 49.1  | 49.4  | 49.7  | 50.1  | 50.4  | 50.7  | 51.0  | 51.4  | 51.7  |
| 160                         | 5.0                                     | 2.3   | 5.7   | 53.0  | 53.3  | 53.6  | 54.0  | 4.3   | 54.6  | 54.9  |
| 170                         | 55.3                                    | 55.6  | 55.9  | 56    | 56.6  | 56.9  | 57.2  | 57.5  | 57.9  | 58.2  |
| 180                         | 58.5                                    | 58.8  | 59    | 9.5   | 59.8  | 60.1  | 60.5  | 60.8  | 61.1  | 61.4  |
| 190                         | 61.8                                    | 6.1   | 6.4   | 6.7   | 63.1  | 63.4  | 63.7  | 64.0  | 64.4  | 64.7  |
| 200                         | 6.0                                     | 65.3  | 65.7  | 66.0  | 66.3  | 66.6  | 67.0  | 67.3  | 67.6  | 67.9  |
| 10                          | 68.3                                    | 68.6  | 68.9  | 69.6  | 69.6  | 69.6  | 70    | 70.5  | 70.9  | 71    |
| 0                           | 71.0                                    | 71.8  | 7.2   | 7     | 7.8   | 73.1  | 73.5  | 73.8  | 74.1  | 74.4  |
| 230                         | 74.5                                    | 75.1  | 5.4   | 75.7  | 76.1  | 76.4  | 6     | 77.0  | 7.4   | 77.7  |
| 40                          | 78.0                                    | 78.3  | 78.7  | 9.0   | 79.3  | 79.6  | 80.0  | 80.3  | 80.6  | 80.9  |
| 250                         | 81.3                                    | 81.6  | 81.9  | 8.5   | 8.6   | 8.9   | 83.2  | 83.5  | 83.9  | 84.2  |
| 260                         | 84.0                                    | 84.3  | 85    | 85.5  | 85.8  | 86.1  | 86.5  | 86.8  | 87.1  | 8.4   |
| 70                          | 87.8                                    | 88.1  | 88.4  | 88.7  | 89.1  | 89.4  | 89.7  | 90.0  | 90.4  | 90.7  |
| 280                         | 91.1                                    | 91.3  | 91.7  | 9.0   | 9.3   | 9.6   | 93.0  | 9.3   | 93.6  | 93.9  |
| 290                         | 94.3                                    | 94.6  | 94.9  | 95.2  | 95.6  | 95.9  | 96    | 96.3  | 96.9  | 97    |
| 300                         | 97.5                                    | 97.8  | 98    | 98    | 98.8  | 99.1  | 99.5  | 99.8  | 100.1 | 100.4 |
| 310                         | 100.8                                   | 101.1 | 101.4 | 101.7 | 102.1 | 10.4  | 10    | 103.0 | 103.4 | 103.7 |
| 320                         | 104.0                                   | 104.3 | 104.7 | 105.0 | 105.3 | 105.6 | 106.0 | 106.3 | 106.6 | 106.9 |
| 330                         | 107.3                                   | 107.6 | 107.9 | 108   | 108.6 | 108.9 | 109.2 | 109.5 | 109.9 | 110   |
| 340                         | 110.5                                   | 110.8 | 111   | 111.5 | 111.8 | 112.1 | 11.5  | 11.8  | 113.1 | 113.4 |
| 350                         | 113.8                                   | 114.1 | 114.4 | 114.7 | 115.1 | 115.4 | 115.7 | 116.0 | 116.4 | 116.7 |
| 360                         | 117.0                                   | 117.3 | 117.7 | 118.0 | 118.3 | 118.6 | 119.0 | 119.3 | 119.6 | 119.9 |
| 370                         | 1.03                                    | 120.6 | 120.9 | 1.1   | 1.16  | 1.19  | 1     | 1.5   | 1.9   | 1.3   |
| 380                         | 1.35                                    | 1.38  | 1.4   | 124.5 | 1.48  | 1.51  | 1.55  | 1.58  | 126.1 | 1.64  |
| 390                         | 126.8                                   | 1.71  | 1.74  | 127.7 | 128.1 | 1.84  | 128.7 | 1.90  | 129.4 | 129.7 |

## 92 *Biochemical Procedures*

- 5 To the ash add 0.7 ml of 20%  $\text{H PtCl}$
- 6 Agitate the tube for  $\frac{1}{2}$  to one minute or until complete solution occurs and  $\text{K PtCl}$  begins to precipitate
- 7 Add approximately 5 ml absolute  $\text{EtOH}$  and mix contents thoroughly with stirring rod. Rinse the rod with 1.2 ml absolute  $\text{EtOH}$ .
- 8 Let the tubes stand for 30 minutes
- 9 Centrifuge gently for 10 minutes at about 1500 r p m
- 10 Aspirate or decant supernatant liquid. (Save it for platinum recovery if desired)
- 11 Wash contents with 5 ml of absolute  $\text{EtOH}$  mix with stirring rod and rinse rod as before
- 12 Centrifuge again for 2.3 minutes at 1500 r p m, aspirate the supernatant liquid and save residues
- 13 Place tubes in 70 C water bath
- 14 Bring temperature to boiling and slowly maintain for 5.10 minutes until  $\text{EtOH}$  is completely volatilized. If  $\text{EtOH}$  is not completely evaporated it will reduce the iodoplatinate to free platinum black
- 15 Leave the tubes in hot water and add approximately 4.5 ml of water to each.
- 16 Dissolve salt ( $\text{K PtCl}$ ) completely. If the salt is not dissolved free platinum will be formed upon addition of  $\text{KI}$ .
- 17 Add approximately 1 ml 2 N  $\text{KI}$  mix and let stand 5 minutes in hot water
- 18 Then cool and titrate with 0.2 ml Rehberg burette using 0.1000 N  $\text{Na}_2\text{S}_2\text{O}_3$  until complete disappearance of red color. The remaining solution will be lemon yellow. Save the residue for recovery of the platinum present

### Calculation

If 1 ml serum used  $\text{mEq K/liter serum} = \text{ml of } \text{Na}_2\text{S}_2\text{O}_3 \times \text{normality} \times 1000$

To convert  $\text{mEq potassium per liter}$  to  $\text{mg potassium per 100 ml}$  multiply  $\text{mEq/liter}$  by 3.9 (i.e. equivalent weight of potassium  $\times 1/10$ )

### Example

Titration figure was 0.0530 ml of 0.1000 N sodium thiosulfate and the blank was 0.0030 ml

$$\text{mEq K/liter of serum} = (0.0530 - 0.0030) \times 0.1000 \times 1000 = 5.0$$

### Precautions

- 1 Make certain that all the alcohol is evaporated before adding the water
- 2 Then make certain that the precipitate is dissolved in the water before adding the potassium iodide. (If not free platinum will be formed upon addition of the  $\text{KI}$ )
- 3 Be sure that the ash is white otherwise muffle again
- 4 It is wise to check the accuracy of the method with a standard solution (see total base)

## B POTASSIUM IN RED BLOOD CELLS

### Reference

Consolazione W V (Unpublished)

- 1 Pipette exactly 0.1 ml of cells with automatic syringe pipette into quartz or Vycor tube as for serum
- 2 Add 1 drop of saturated sodium sulfate ( $\text{Na}_2\text{SO}_4$ )
- 3 Add 1 drop of  $\text{H}_2\text{SO}_4$ . Dry at 110 C
- 4 In order to ash add 1 drop of concentrated  $\text{HCl}$  and place in 100 C oven until dry. This step prevents occlusion of K salts on the iron, since all ferric oxide is converted to ferric chloride
- 5 Now proceed as for Serum, steps 4.18

## C POTASSIUM IN URINE

- 1 Into a quartz or Vycor tube pipette 0.2 ml of urine  
Add 1 drop of sodium sulfate solution
- 3 Add 1 drop of sulfuric acid solution
- 4 Proceed as for Serum steps 3-18

## D POTASSIUM IN DIET STOOLS AND TISSUE

- 1 Use approximately 0.5 to 1 gm for food feces or tissue. Ash in muffle.  
To the proper amount of HCl extract of ash add 1 drop of  $\text{Na}_2\text{SO}_4$ . Do not add  $\text{H}_2\text{SO}_4$ .  
Dry at 110°C. Do not ash
- 4 Proceed as for Serum steps 5-18
- 5 Pyrex tubes may be used here in place of quartz tubes since no ashing is employed after step

## 5 Calcium in Serum Urine Food Tissue and Feces


### A. IN SERUM

#### References

- Kramer B. and Tisdall F. F. A Simple Technique for the Determination of Calcium and Magnesium in Small Amounts of Serum. *J Biol Chem* 47: 475-481 (Aug.) 1921
- Clark E. P. and Collip J. B. A Study of the Tisdall Method for the Determination of Blood Serum Calcium With a Suggested Modification. *J Biol Chem* 63: 461-464 (March) 1925

#### Principle

Calcium is precipitated as the oxalate which is titrated in hot acid solution against standard potassium permanganate



#### Apparatus

- 1 A quantity of 15 ml conical centrifuge tubes
- A water bath with boiling water
- 3 A 3 or 5 ml microburette
- 4 An electric centrifuge
- 5 Accurately calibrated 1 and 2 ml syringe pipettes

#### Reagents

- 1 Phenol red approximately 0.04% in water
- 2 Ammonium oxalate approximately 45 gm per liter
- 3 Ammonium hydroxide approximately 0.35 N i.e. approximately 2 ml of concentrated  $\text{NH}_4\text{OH}$  diluted to 100 ml with water
- 4 Sulfuric acid approximately N prepared by diluting 3 ml concentrated  $\text{H}_2\text{SO}_4$  to 100 ml with water
- 5 Standard potassium permanganate 0.01 N. This is best prepared from a well aged 0.1 N standard (3.161 gm  $\text{KMnO}_4$  per l) and should be standardized weekly against oxalate as follows:
  - a A standard 0.1 N solution of sodium oxalate is prepared by diluting exactly 6.700 gm of pure dry  $\text{Na}_2\text{C}_2\text{O}_4$  to exactly one liter
  - b A standard 0.01 N solution is prepared by diluting exactly 10 ml of the 0.1 N solution to exactly 100 ml.

## 94 *Biochemical Procedures*

c In conical tubes place exactly 2 ml of the 0.01 N standard and approximately 2 ml of 1 N H<sub>2</sub>SO<sub>4</sub>. Heat in the water bath at 60-70 C and titrate while hot against the dilute KMnO<sub>4</sub>.

d Calculation    Normality of dilute KMnO<sub>4</sub> =  $\frac{2}{\text{Titration figure in ml}} \times 0.01$

6 Ammonia    Approximately a 2% solution

7 Calcium Standard (Primary Standard)

0.250 gm of dry calcium carbonate (pure Iceland spar) in a small amount of water in a 1000 ml flask. Add some dilute HCl (1:10) to dissolve the crystals and then add a few ml of acid in excess. Dilute to 1000 ml with water and mix. This solution is stable indefinitely and contains 0.1 mg of Ca per ml. It corresponds to serum with a Ca content of 10 mg/100 ml.

### Procedure

- 1 Pipette exactly 2 ml of serum into a conical tube
- 2 Add approximately 1 ml of ammonium oxalate and 1 drop of phenol red. Mix by rotation of the tube
- 3 Add 2% ammonia drop by drop until the indicator turns red
- 4 Stand at room temperature for at least 10 minutes for precipitation and then centrifuge at high speed for 10 minutes. Do not use an angle head centrifuge
- 5 Using a fine capillary and slow suction, suck off and discard the supernatant fluid taking great care not to lose any of the precipitate. Do not try to decant the supernatant
- 6 Pipette approximately 3 ml of 2% NH<sub>4</sub>OH down the sides of the tube and stir the precipitate by twirling the tube by hand.
- 7 Centrifuge at high speed for 3 minutes, suck off the supernatant fluid and repeat steps 6 and 7 until no color of phenol red remains in the precipitate. Wash 4 times in this way.
- 8 Tap the tube sharply to loosen the precipitate and add 2 ml of 1 N H<sub>2</sub>SO<sub>4</sub>
- 9 Place tubes in boiling water until the precipitate dissolves
- 10 Titrate with 0.01 N KMnO<sub>4</sub> to the first permanent pink color lasting 30 seconds. Tube and contents must never become cooler than 70 C during the titration.
- 11 With each run set a reagent blank consisting of 2 ml of 1 N H<sub>2</sub>SO<sub>4</sub> and titrate as in Steps 9 and 10. The blank should not be above 0.05 ml if reagents are satisfactory.

### Calculation

$$\frac{\text{mg Ca/100 ml serum} = (\text{ml titration} - \text{ml blank})}{2} \times \text{Normality KMnO}_4 \times 10.04 \times 100 = (\text{ml titration} - \text{ml blank}) \times \text{Normality KMnO}_4 \times 1002$$

### Example

Normality of KMnO<sub>4</sub> was 0.0095. Titration was 1.10 ml. Blank was 0.03 ml.  
Therefore mg Ca per 100 ml serum = (1.10 - 0.03) (0.0095) (1002) = 10.2

### Precautions

- 1 The phenol red turns reddish purple in Step 3, never true purple
- 2 Take great care not to lose precipitate in Steps 5, 6 and 7
- 3 The titration must be carried out above 70 C. When the first drop of KMnO<sub>4</sub> is added it persists for some time. As soon as some manganous salt is formed, the titration proceeds more rapidly to the end point.

- 4 Wash at least 4 times with ammonia. The amount and type of washing should be determined for each laboratory. It depends upon shape of tubes, efficiency of centrifuge and probably the quality of reagents especially distilled water. In our laboratory we get very erratic results with only two washings.
- 5 Be sure that the distilled water is of constant quality free of calcium.

## B. CALCIUM IN URINE

### Reference

(Same as for serum)

### Principle

Same as in serum calcium. The presence in urine of magnesium, iron and phosphate in significant amounts makes attention to the acidity during precipitation very important.

### Apparatus

(Same as for serum)

### Reagents

- 1 Phenol red, sulfuric acid and standard potassium permanganate as in serum calcium.
- 2 Ammonium hydroxide 4N. Dilute 9 ml of concentrated  $\text{NH}_4\text{OH}$  to 100 ml with water.
- 3 Acetic acid. Dilute 10 ml of glacial acetic acid to 100 ml with water.
- 4 Sodium acetate saturated. Approximately 100 gm are added to 100 ml of water.
- 5 Oxalic acid 1N. Dissolve 63 gm of acid in one l of water.

### Procedure

- 1 Use freshly collected urine without added stabilizers particularly without added oxalate. Pipette exactly 1 ml of urine into a conical tube.
- 2 Add 1 drop of phenol red and then 4 N  $\text{NH}_4\text{OH}$  drop by drop until the indicator is pink.
- 3 Quickly add 10<sup>cc</sup> acetic acid drop by drop until indicator turns yellow.
- 4 Add 1 ml of oxalic acid and 1 ml of sodium acetate.
- 5 Mix well and stand 1 1/2 hours at room temperature.
- 6 Proceed exactly as in Steps 4 through 10 in serum calcium.

### Calculation

$$\begin{aligned} &\text{mg Ca per 100 ml urine} = \\ &\frac{(\text{ml titration} - \text{ml blank})}{\text{ml urine}} \times (\text{normality KMnO}_4) \times 20.04 \times 100 = \\ &(\text{ml titration} - \text{ml blank}) \times \text{normality KMnO}_4 \times .004 \end{aligned}$$

### Example

Normality of  $\text{KMnO}_4$  was 0.0096. Titration was 1.55 ml. Blank was 0.05 ml.  
Therefore mg Ca per 100 ml urine =  $(1.55 - 0.05) (0.0096) (.004) = .88$

### Precautions

- 1 Use fresh urine before precipitation has taken place. If the urine must be allowed to stand, add 4 ml of glacial acetic acid per 100 ml of urine and correct appropriately in the calculations.
- 2 Other precautions as in serum calcium.
- 3 Wash at least 4 times with ammonia. The exact number of washings has to be determined for each laboratory.



- 4 This method is not very specific Under some conditions as during therapy with sulfaguanidine, urine contains substances which interfere with the estimation Under such conditions the urine has to be ashed with sulfuric acid 3 drops of concentrated sulfuric acid to every 2 ml of urine

### C CALCIUM IN FOOD TISSUE AND FECES

#### Reference

- Corley, R C, and Denis W The Determination of Calcium in Tissue Feces and Milk *J Biol Chem* 66 601 608 (Dec ) 1925

#### Apparatus

- 1 A quantity of porcelain or platinum crucibles
- 2 Muffle oven set at 600 C
- 3 A quantity of good analytical filter paper (Whatman #40) and funnels
- 4 Rubber stoppers
- 5 Rest of apparatus same as for serum

#### Reagents

- 1 Concentrated nitric acid.
- 2 Hydrochloric acid one ml concentrated HCl to two parts of water
- 3 Alizarin red indicator in water
- 4 Litmus paper
- 5 Oxalic acid, 2.5 percent in water
- 6 Ammonium oxalate 3 percent in water
- 7 Ammonium oxalate, 0.5 percent in water
- 8 Sodium acetate 4.0 percent in water
- 9 Sulfuric acid (1 N)
- 10 Concentrated ammonia
- 11 The rest of the reagents are the same as for serum

#### Procedure

- 1 For food tissue and feces use samples weighing approximately 2 gm
- 2 Ignite in a porcelain or platinum crucible to a white ash with the aid of a few drops of conc HNO<sub>3</sub> and dissolve with HCl
- 3 Bring the volume to approximately 75 to 150 ml.
- 4 Make alkaline with concentrated ammonia adding it drop by drop Use litmus or alizarin red as an indicator
- 5 Then add concentrated HCl drop by drop until just acid to litmus
- 6 Add 10 drops of concentrated HCl (sp gr 1.20)
- 7 Add 10 ml of 2.5% oxalic acid
- 8 Stopper flask with a rubber stopper mix and stand for 12 hours
- 9 Add an excess of 3% ammonium oxalate and cool to room temperature
- 10 Then add 8 ml of 20% sodium acetate solution (For feces use 15 ml)
- 11 Stopper and shake vigorously for 10 minutes
- 12 Filter and wash precipitate with 0.5% ammonium oxalate (This is to wash free of chloride)
- 13 Then wash 3 times with cold distilled water The exact number of washings has to be determined in each laboratory
- 14 Transfer precipitate quantitatively to a 50 ml beaker washing last traces off filter paper with approximately 5 ml of 1 N sulfuric acid
- 15 Add approximately 25 ml of water and stir until precipitate is all dissolved. It may be necessary to add more sulfuric acid or to warm gently
- 16 Dilute with water to exactly 100 ml in a volumetric flask.
- 17 Titrate 10 ml aliquots with 0.01 N KMnO<sub>4</sub> as for urine

**Calculation**

$$\text{mg Ca per 100 gm sample} = \frac{(\text{ml titration} - \text{ml blank})}{(\text{gm sample} \times 10/100)} \times (\text{normality KMnO}_4) \times (0.04) \times (100)$$

**Example**

A sample of food weighed 1.9546 gm. Normality of  $\text{KMnO}_4$  was 0.0095. Sample was diluted to 100 ml and 10 ml aliquots were titrated. Titration figure was 50 ml. Blank was 0.04 ml.


$$\text{mg Ca/100 gm food} = \frac{(2.50 - 0.04)}{(1.9546 \times 10/100)} \times 0.0095 \times 0.04 \times 100 = 73.5$$

**6 Magnesium in Serum****Reference**

Briggs A. P. A Colorimetric Method for the Determination of Small Amounts of Magnesium. *J Biol Chem* 52: 349-355 (June) 1922

**Principle**

Magnesium is precipitated as magnesium ammonium phosphate and the phosphate is estimated by the colorimetric method of Tschopp and Tschopp or that of Fiske and Subbarow.

**Apparatus**

1. A quantity of 15 ml conical centrifuge tubes
2. An electric centrifuge
3. A quantity of rubber policeman
4. Refrigerator
5. Water bath set between 60-70°C
6. Coleman Jr Spectrophotometer Model 6
7. Cuvettes 19 x 150 mm

**Reagents**

1. Monobasic potassium phosphate ( $\text{KH}_2\text{PO}_4$ ) approximately 2% solution in water
2. Ammoniacal alcohol. This is prepared from 200 ml of 95%  $\text{EtOH}$  and 50 ml of concentrated  $\text{NH}_4\text{OH}$  diluted to one l with water
3. Trichloroacetic acid ( $\text{CCl}_3\text{COOH}$ ) 15% solution in water. Store in refrigerator
4. Tschopp and Tschopp phosphorus reagents (See phosphorus method for reagents)
5. Ammonium oxalate ( $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$ ) approximately 40% in water
6. Concentrated ammonium hydroxide ( $\text{NH}_4\text{OH}$ ) low in magnesium and phosphate. Hydrochloric acid approximately 1 N and N/5
7. Standard magnesium solution. 0.1413 gm of magnesium ammonium phosphate ( $\text{MgNH}_4\text{PO}_4 \cdot 6\text{H}_2\text{O}$ ) per liter in approximately 0.01 N sulfuric acid. Add 2 ml of chloroform for a preservative. Of this 5 ml = 0.07 mg Mg
8. Molybdic acid photorex and sodium bisulfite as described in the method on phosphorus

**Procedure**

- 1 Mix in a 15 ml conical centrifuge tube
  - 20 ml of serum
  - 20 ml of water
  - 10 ml of ammonium oxalate solution
 Mix by inversion and centrifuge  
 (The precipitate may be used for determination of calcium)
- 2 Into a conical centrifuge tube pipette 40 ml of the supernatant liquid
- 3 Add 60 ml of 1% trichloroacetic acid, mix and centrifuge
- 4 Transfer 50 ml of supernatant fluid to cone tipped 15 ml centrifuge tube
- 5 Add 1 ml of 2% potassium acid phosphate solution ( $\text{KH}_2\text{PO}_4$ )
- 6 Add 4 ml of concentrated ammonia ( $\text{NH}_4\text{OH}$ )
- 7 Scratch side of tube with rubber policeman until precipitate forms and rinse rod with a few ml of water
- 8 Let stand overnight in refrigerator
- 9 Centrifuge decant and wash with 5 ml of ammoniacal alcohol
- 10 Centrifuge and decant again
- 11 Dissolve the precipitate in 1 ml of 1 N HCl
- 12 Transfer quantitatively to a calibrated 25 ml colorimeter tube
- 13 Estimate phosphate (equivalent to magnesium) by adding the following phosphate reagents
  - 0.4 ml of  $\text{NaHSO}_4$
  - 0.4 ml of molybdic acid
  - 0.4 ml of Photorex
- 14 Dilute to exactly 25 ml
- 15 Heat for 20 minutes at 60-70°C in water bath then cool under running water
- 16 Set reagent blank at 100% T at wave length 660 mμ and read unknowns

**Calculation**

- 1 A standard curve is based on samples containing 0.01 to 0.05 mg of magnesium per milliliter and treated as above
- 2 mg of magnesium per 100 milliliters of serum =

$$\begin{aligned}
 & (\text{mg magnesium in cuvette}) \times \frac{2}{(r)} \times \frac{4}{(10)} \times 100 = \\
 & (\text{mg magnesium in cuvette}) \times 16
 \end{aligned}$$

**Example**

Magnesium in cuvette was determined to be 0.04 mg  
 $\text{mg magnesium}/100 \text{ ml serum} = 0.04 \times 16 = 2.4$

**Precautions**

- 1 All reagents must be as low in magnesium as it is possible to procure them

**7 Inorganic Iron in Serum****Reference**

Kitzes G, Elvehjem C A and Schuette H A The Determination of Blood Plasma Iron *J Biol. Chem* 155 653-660 (Oct) 1944.

**Principle**

Serum or plasma is boiled in water to form a gelatinous solution the proteins are precipitated out with trichloroacetic acid and the filtrate is treated with alpha, alpha bipyridine to produce a pink color which is measured photo colorimetrically

## Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 x 150 mm
- 3 A number of 15 ml conical graduated centrifuge tubes
- 4 A number of 1, 3, 4 and 5 ml volumetric pipettes
- 5 A quantity of glass stirring rods
- 6 A constant temperature bath set at 90-95 C
- 7 An electric centrifuge

## Reagents

- 1 Trichloroacetic acid 5% in water If not free of iron this acid should be redistilled.
- Indicator p-nitrophenol 0.1% solution in water
- 3 Ammonium hydroxide solution approximately 6 N in water
- 4 Buffer solution (pH 4.58) To prepare this 27.2 ml of glacial acetic acid and 33.4 gm of anhydrous sodium acetate are dissolved in water and made up to 250 ml
- 5 Thioglycolic acid solution (Eastman Practical Grade)
- 6 Alpha-alpha-bipyridine reagent Dissolve 0.5 gm of reagent in 5 ml of glacial acetic acid and dilute to 100 ml with water
- 7 Redistilled water (iron free) Redistill in a glass still
- 8 Standard iron solution Dissolve iron wire (99.8%) in a mixture of nitric acid and hydrochloric acid and dilute with water to give a concentration of 1 mg/ml. From this solution standard solutions of 100 and 1 mcg per ml are made

## Procedure

- 1 Pipette exactly 3 to 5 ml of plasma or serum (containing approximately 3 to 9 mcg of inorganic iron) into a 15 ml conical centrifuge tube
- 2 Add 3 ml of redistilled water and mix with plasma
- 3 A blank should be prepared at the same time using all the reagents
- 4 Place the tubes in boiling water for 3 minutes or until the solution becomes opaque. This must be done very carefully so as not to get a gelatinous coagulation from which it will be difficult to extract the iron
- 5 Remove the tubes while the solution is still liquid and opaque and allow to cool by placing them in a beaker of cold water
- 6 Add 1 ml of trichloroacetic acid (5%) stir carefully and thoroughly so that the acid is mixed intimately with the solution. The precipitate that is formed becomes fine and flocculant. A small blunt end stirring rod is preferred so as not to run the end of it through the centrifuge tube
- 7 Place the tubes in a water bath set between 90-95 C for about three minutes stirring the solution once or twice
- 8 Remove the tubes from the bath and cool again in cold water
- 9 Then centrifuge the tubes at 2000 to 3000 rpm for about five minutes
- 10 Remove the tubes and carefully decant the supernatant fluid into a 15 ml graduated tube
- 11 To the precipitate remaining in the tubes add 4 ml of redistilled water and 1 ml trichloroacetic acid. Then break up the precipitate and stir it well
- 12 Again place in the water bath (90-95 C) for 3 minutes stirring each solution once or twice. Remove and cool the tubes in cold water
- 13 Centrifuge the tubes for 5 minutes and then decant the supernatant fluid into the tube containing the first filtrate
- 14 Add 1 drop of p-nitrophenol indicator and then add ammonium hydroxide solution drop by drop until the solution becomes yellow
- 15 Add 1 ml of the buffer solution and dilute to exactly 15 ml. Mix thoroughly and centrifuge

- 16 Pipette 10 ml into a cuvette and then add 2 drops of thioglycolic solution. (*Do not get this solution on fingers or skin as it is very toxic*) Mix by tapping the tube or by a slight rotary motion
- 17 Add 1 ml alpha alpha-bipyridine reagent and a faint pink color will develop Mix carefully but thoroughly
- 18 Set the blank at 100% T at 510 mμ and read the unknowns and blank

#### Calculations

A standardization curve is set up on the Coleman Jr Spectrophotometer with standard solutions of iron ranging from one gamma to twenty gamma. The calibration constant *K* is computed as in the section on spectrophotometry

Then mcg inorganic iron/100 ml serum =

$$\frac{K(2 - \log \% T) \times 15 \times 100}{\text{ml aliquot} \times \text{ml serum}}$$

#### Example

Five ml of serum were analyzed *K* was 45.32 %T'' was 8.5 Aliquots of 10 ml were used.

mcg inorganic iron/100 ml serum =

$$\frac{45.32 (2 - \log 8.5) \times 15 \times 100}{10 \times 5} = 113.5$$

#### Precautions

Thioglycolic acid is very toxic Be careful with it

## 8 Copper in Serum

### A. METHOD I

#### Reference

Robinson J C A Simple Method for the Determination of Serum Copper *J Biol Chem* 179 1103 1109 (July) 1949

#### Principle

The addition of sodium diethyldithiocarbamate to a solution containing copper produces a golden brown color that can be measured colorimetrically. The copper salt can be rapidly and quantitatively extracted from an aqueous solution by isomyl alcohol which also intensifies the color. This color is stable for at least two hours at a pH between 5.7 and 9.2

Iron interferes with the color but the addition of sodium pyrophosphate converts it to iron pyrophosphate which does not react with the carbamate

#### Apparatus

- 1 Test tubes 13 x 100 mm with ground glass stoppers
- 2 Pipettes. Serological 0.2 ml graduated to 0.01 ml
- 3 Mechanical shaker Kahn
- 4 A water bath set between 50 to 60 C
- 5 Spectrophotometer In this laboratory a Coleman Model 6 equipped with 1 x 75 mm cuvettes is used
- 6 Syringe pipettes exactly 0.1 ml 0.2 ml 1.0 ml and 3.0 ml.
- 7 Glass still for metal free water
- 8 Electric centrifuge

## Reagents

- 1 Sodium pyrophosphate Saturate red distilled water with the salt
- 2 Ammonium hydroxide 10 to 12% Distill concentrated ammonia (30% reagent) into red distilled water to saturation Distill to required concentration with red distilled water
- 3 Sodium diethylthiocarbamate 30% aqueous solution Remove copper by shaking with carbon tetrachloride which has been distilled over calcium oxide The solution is stable for some weeks if stored in the dark. (Sodium diethylthiocarbamate the reagent used to react with copper in this procedure gives colored complexes with a number of heavy metal cations In biological material however on is the chief interfering substance This interference of iron is satisfactorily eliminated by developing the copper carbamate complex in an ammoniacal solution containing pyrophosphate and having pH 9 or greater The color intensities of other interfering metals such as nickel and cobalt are weight for weight only 1/20 1/30 as great as that given by copper Because the concentration of these ions in serum is of such relatively low order the interference is negligible)
- 4 Metal free water Distilled water redistilled from resistant glass apparatus should be used in the preparation of all reagents and in the final running of equipment
- 5 Standard copper solution 0.100 mg Cu/ml Prepare stock solution by dissolving 0.98 gm of uneffloresced  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in red distilled water and diluting to one liter From this stock solution on the working standards may be prepared It is best to use a solution containing 1 mcg per ml because the final volume of the aqueous phase should be kept at a minimum To prepare this the dilution is 1 ml of stock solution to 100 ml with water
- 6 Isoamyl alcohol Distill from all glass apparatus saturate with water and store at 4°C

## Procedure

- 1 Into each of glass stoppered test tubes accurately deliver 10 ml serum Add to each of the tubes 0.2 ml of saturated solution of sodium pyrophosphate Mix by tapping gently
- 2 To each of the tubes add exactly 0.04 ml of 10.1% ammonium hydroxide Again mix by tapping
- 3 Add 0.2 ml sodium diethylthiocarbamate to one of the tubes and 0.1 ml water to the other which is the serum blank Mix and allow to stand for one hour
- 4 Pipette exactly 3.0 ml of isoamyl alcohol into each tube
- 5 Place tubes in lab rack and securely fasten into a horizontal position on a mechanical shaker (It is convenient to place a sponge against one end of the shaker and clamp the rack so that the stoppered ends of the tubes are firmly pressed against the sponge) Shake for 15 minutes
- 6 Place tubes in centrifuge cups and cool until water used in balancing the cups begins to form crystals of ice
- 7 Centrifuge at 4000 rpm for 10 minutes If the minimum required volume (15 ml for 1 × 75 mm cuvette) of extract is not obtained gently shake the precipitate free into the isoamyl alcohol and repeat centrifugation
- 8 Carefully pipette the isoamyl alcohol fractions into 1 × 5 mm cuvettes stopper and measure optical densities at 440 mμ (If cloudy extracts are obtained place cuvettes for a few seconds in water warmed to 50-60°C)
- 9 It is necessary to prepare a reagent blank substituting distilled water for serum for each series of determinations Very low readings about 99-100 percent transmittance are obtained.

**Calculations**

Let  $O_c$  represent the optical density due to the serum copper,  $O$  the optical density of the carbamate extract from serum,  $O$  the optical density of the serum extract and  $O$  the optical density of the reagent blank

Then  $O_c = O - (O + O)$

and Copper mcg per 100 ml serum  $\approx$

$O_c \times \text{calibration constant} \times \text{dilution factor} \times 100$




Fig 16—Example of calibration curve for copper in serum. Ordinate optical density abscissa copper mcg per ml of water

**Example**

$O$  was 0.09  $O$  was 0.06  $O$  was 0.01, and  $K$  (mcg copper per ml per unit optical density) was 17.5

In the analysis, the copper of 1 ml serum was extracted into 3 ml of isoamyl alcohol therefore copper mcg per 100 ml serum  $\approx$

$$(0.09 - 0.06 - 0.01) \times 17.5 \times \frac{3}{1} \times 100 = 105$$

## Experimental Considerations

As previously mentioned the final color is due to three components serum copper, copper as an impurity in the reagents, and substances other than copper carbamate. A large percentage of the density of the final color is due to materials extracted from the serum other than copper carbamate. This point is demonstrated in Figs 16 and 17 which show that the maximum light absorption for both an isoamyl alcohol extract of copper carbamate from aqueous solution and an isoamyl alcohol extract of the serum blank occurs in the same general region of the spectrum.




Fig 17—Spectral transmission curve of copper carbamate and serum. Ordinate optical density, abscissa wave length in millimicrons.

The ratio of color to the concentration of copper carbamate in isoamyl alcohol is linear up to a concentration of about 7 mg. Above this value however there is an appreciable deviation from Beer's law. This should be borne in mind particularly when doing recovery experiments and in such experiments concentration figures should be obtained from a calibration curve rather than use of the calibration constant.

## Precautions

1. Be certain the tubes are cool before centrifuging the specimens.
2. All reagents must be nearly free of metals of the iron, cobalt, nickel series as well as of copper.



**B METHOD II****Reference**

Cartwright E Jones J and Wintrobe M A Method for the Determination of Copper in Blood Serum *J Biol Chem* 160 593 600 (Oct ) 1945

**Principle**

The addition of sodium diethyldithiocarbamate to a solution containing copper produces a golden brown color that can be measured colorimetrically The copper salt can be rapidly and quantitatively extracted from an aqueous solution by isoamyl alcohol which also intensifies the color This color is stable for at least 2 hours at pH 5.7 to 9.2

Iron interferes with the color but the addition of sodium pyrophosphate converts it to iron pyrophosphate which does not react with the carbamate

**Apparatus**

- 1 Coleman Jr Spectrophotometer Model #6
- 2 Cuvettes 19 x 150 mm
- 3 A quantity of 15 ml centrifuge tubes graduated
- 4 Electric centrifuge
- 5 A water bath set at 90-95 C
- 6 A quantity of glass stirring rods
- 7 Glass water still

**Reagents**

- 1 Redistilled water is used in all the reagents
- 2 Trichloroacetic acid 20% in water
- 3 Sodium pyrophosphate saturated solution
- 4 Ammonium hydroxide solution 10 to 12 percent
- 5 Sodium diethyldithiocarbamate solution 1%
- 6 Standard copper solution as used in the method of Robinson (1949)
- 7 Isoamyl alcohol

**Procedure**

All glassware used must be carefully cleaned and then rinsed three times with redistilled copper free water

- 1 Pipette 3 to 5 ml samples preferably 5 ml of serum or plasma into ungraduated 15 ml Pyrex centrifuge tubes It is convenient to set up that number of tubes which can be centrifuged simultaneously
- 2 Add 3 ml of redistilled water to each mix with glass rods and place in boiling water until they become opaque This must be done very carefully or else the solutions will become very gelatinous and coagulation of this type will hinder the extractions which are to follow
- 3 When the tubes are opaque and still liquid, remove them and cool in a beaker containing cold water
- 4 Add 2 ml of 20% trichloroacetic acid mix thoroughly and break up the precipitate that forms so that it is finely flocculant Then place in water bath at 90-95 C for five minutes
- 5 Stir this solution carefully and frequently
- 6 Cool as before, remove the rods and centrifuge the tubes at 3000 rpm for about five to ten minutes
- 7 Remove and decant the supernatant liquid into a clean dry 15 ml graduated centrifuge tube
- 8 To the original tubes add 2 ml of water and 1 ml of trichloroacetic acid being sure to break up the precipitate well by stirring the solution carefully and thoroughly

- 9 Heat again at 90-95 °C for five minutes by carefully stirring the solution to enhance the extraction of the copper
- 10 Cool again in the cold water centrifuge and decant the solution into the same tubes containing the first filtrate
- 11 Repeat this operation once again using the same quantities of the water and trichloroacetic acid.
- 12 Make the filtrate up to 15 ml. It will contain some precipitate which came over with each decantation
- 13 Centrifuge this down draw off a 10 ml aliquot and place this aliquot in a Coleman Jr 150 x 19 mm spectrophotometer tube
- 14 A blank of all the reagents is prepared in the same way. This may be prepared directly in the Coleman tube
- 15 To the filtrate and blank in the Coleman tubes add 1 ml of a saturated solution of sodium pyrophosphate
- 16 Add 4 ml of ammonium hydroxide solution
- 17 Add 1 ml of the 0.1% sodium diethylthiocarbamate solution
- 18 Add 1 ml of distilled water to bring the total volume up to 15 ml. Mix the solutions thoroughly and carefully
- 19 Set the blank at 100% T (Optical density) and read the unknowns as soon as possible

#### Calculation

With the blank set at 100% T copper mcg per 100 ml serum = (concentration in cuvette)  $\times \frac{(\text{final dilution})}{(\text{ml serum})} \times \frac{(\text{total dilution})}{(\text{aliquot})} \times 100$

#### Example

5 ml of serum were diluted finally to 15 ml in trichloroacetic acid. The photometer reading was 9% T in this case equivalent to 0.3 mcg of copper per ml in the cuvette

Therefore copper mcg per 100 ml serum =  $0.3 \times \frac{15}{5} \times \frac{15}{10} \times 100 = 144$

#### Precautions

- 1 Be careful that the solutions do not become gelatinous as this will vitiate results
- 2 All reagents must be nearly free of heavy metals

### 9 Chloride in Serum and Urine




#### A. CHLORIDE IN SERUM

##### Reference

Keys A. The Microdetermination of Chlorides in Biological Materials J Biol Chem 119: 39-403 (July) 1937

##### Principle

The chloride is precipitated as silver chloride with standard silver nitrate and the excess silver is titrated against thiocyanate with ferric alum as indicator



##### Apparatus

- 1 A number of round bottom test or centrifuge tubes of 15 ml capacity
- A variety of volumetric pipettes including 1 ml, 2 ml and 10 ml graduated in 0.1 ml divisions

- 3 One dropping bottle for storing indicator
- 4 Water bath (boiling)
- 5 Assorted bottles, beakers and flasks for reagents and samples
- 6 Automatic syringe pipettes to deliver 0.1 and 0.2 ml
- 7 Micro burette, 3 ml

### Reagents

For field use it is convenient to carry all reagents except nitric acid weighed out dry in small bottles and to dilute them with distilled water as needed. When many analyses are to be done the best procedure is to have the facilities to make up 'strong standards' and to prepare 'working standards' by diluting the 'strong standard'.

- 1 'Strong'  $\text{AgNO}_3$ . In a mixing cylinder dissolve 23.785 gm  $\text{AgNO}_3$  and dilute to exactly 1 liter with distilled water. Tap water may be used if it is chloride free.
- 2 'Working standard'  $\text{AgNO}_3$ . Dilute exactly 100 ml of strong standard to 1 liter with distilled or chloride free tap water.
- 3 'Strong'  $\text{KCNS}$ . Dilute 18.434 gm dry crystals to 1 liter with distilled or chloride free tap water.
- 4 'Working standard'  $\text{KCNS}$ . Dilute exactly 100 ml of 'strong  $\text{KCNS}$ ' to 1 liter with distilled or chloride free tap water.
- 5  $\text{HNO}_3$ , concentrated, chloride free, reagent grade.
- 6 Ferric alum indicator. Prepare a saturated solution of ferric ammonium sulphate in distilled or chloride free tap water using approximately 40 gm to every 100 ml of water. Just before use, add conc  $\text{HNO}_3$  drop by drop until the color is yellow not red. This converts ferrous ions to ferric. For use, the indicator is stored in a dropping bottle.
- 7 Standard solution of  $\text{NaCl}$ , 100 mEq/liter. Dilute 5.845 gm pure  $\text{NaCl}$  to exactly 1 liter. This reagent is used solely for standardizing the silver nitrate and potassium thiocyanate and keeps indefinitely.
- 8 Distilled or chloride free tap water.

### Procedure

- 1 Prepare and store samples of serum as described under the section on handling blood.
- 2 Deliver exactly 0.2 ml of serum by means of an automatic pipette into a 15 ml test tube.
- 3 Deliver into the same tube 5 to 8 successive washings of the pipette with distilled or chloride free tap water.
- 4 Add exactly 2.0 ml of 'working standard'  $\text{AgNO}_3$  to each tube.
- 5 Add approximately  $\frac{1}{4}$  ml of conc  $\text{HNO}_3$  to each tube. Digest on a boiling water bath not less than a half hour. More does not matter.
- 6 Cool tubes in ice or cold water bath. This step is essential for a sharp end point.
- 7 Add 3 or more drops of ferric alum indicator.
- 8 Fill the microburette with 'working standard'  $\text{KCNS}$  and titrate to the first permanent pink color lasting a half minute or more.

### Standardization

If the quality of the chemicals is good and weights are accurate the 'strong'  $\text{AgNO}_3$  will be 0.1400 N and the 'working standard' 0.0140 N. The 'strong'  $\text{KCNS}$  will be 0.2000 N and the 'working standard' 0.0200 N. This ideal situation rarely occurs and the solutions must be standardized as follows.

- 1 Pipette exactly 5 ml of 'working standard'  $\text{AgNO}_3$  into a small Erlenmeyer flask.
- 2 Add 10 drops of ferric alum indicator.

- 3 Add approximately 1 ml conc HNO
- 4 With working standard KCNS in the burette titrate to the first pink color lasting 1 minute
- 5 Repeat until satisfactory triplicate estimations are obtained
- 6 Run the standard solution of sodium chloride exactly as for serum until satisfactory triplicates are obtained
- 7 The normalities of the working standards of AgNO and KCNS may be calculated from the two equations

$$(a) \text{ Normality AgNO} = \frac{\text{Titration in step 5}}{5} \times \text{Normality KCNS}$$

$$(b) \text{ Normality NaCl} =$$

$$\frac{(\text{Normality AgNO}) - (\text{ml step 6} \times \text{Normality KCNS})}{0.00}$$

- 8 With the known values substituted equations (a) and (b) may be solved as simultaneous equations
- 9 An example of standardization is

Concentration of NaCl standard = 0.100 N  
 Calibration of syringe pipette = 0.00 ml  
 Titration in Step 5 = 3.66 ml  
 Titration in Step 6 = 0.440 ml  
 Let normality of AgNO = (AgNO) and  
 normality of KCNS = (KCNS)

$$(a) (\text{AgNO}) = \frac{3.66}{5} \times (\text{KCNS}) = 0.732 (\text{KCNS})$$

$$(b) 0.100 = \frac{2 (\text{AgNO}) - 0.440 (\text{KCNS})}{0.00}$$

$$0.00 = 2 (\text{AgNO}) - 0.440 (\text{KCNS})$$

Substituting from (a) above

$$0.00 = 2 \times 0.732 (\text{KCNS}) - 0.440 (\text{KCNS})$$

$$0.00 = 1.024 (\text{KCNS})$$

$$\text{and } (\text{KCNS}) = \frac{0.000}{1.024} = 0.0195$$

$$\text{From (a) above } (\text{AgNO}) = (0.732) \times (0.0195) = 0.0143$$

#### Calculation

- 1 In the case of serum mEq chloride/liter serum =

$$\frac{(2 \times \text{Normality AgNO}) - (\text{ml titration} \times \text{Normality KCNS})}{\text{Calibration of Syringe Pipette in ml}} \times 1000$$

- 2 Line charts are useful for routine calculations when reagents and pipettes are always the same. The line chart for serum is constructed by marking off the vertical axis according to the titration figure and the horizontal axis for serum chloride. Two points only are required with a straight line drawn accurately between. These two points are obtained by substituting all of the figures which are constant in the equation above and then calculating the serum value for titration on figures of 0.200 and also 0.500 ml.

For urine set up a chart with the vertical axis titration figure from 0 to 1.500 ml and urinary chloride from 0 to 2 gm sodium chloride per 100 ml as the

horizontal axis The two necessary points are obtained by substituting all of the constant values in the equation above and then calculating for titration figures of 0 and 1.350 ml

**Example**

Working standard ' AgNO = 0.0139 N

Working standard' KCNS = 0.0234 N

Titration figure = 0.370 ml of KCNS

$$\text{mEq chloride/liter serum} = \frac{(2 \times 0.0139) - (0.370 \times 0.0234)}{0.200} \times 1000 = 95 \text{ mEq/l}$$

**Precautions**

- 1 Be sure to titrate with a good daylight or daylight bulb against a white back ground.
- 2 Do not expose to bright light the silver standards or the test tubes containing unknowns to which silver nitrate has been added.
- 3 Chloride in serum and urine is of course stable indefinitely at any temperature. Therefore samples may be kept for months before analysis provided there is no evaporation.
- 4 The operator may leave chloride analyses for hours at any stage provided the tubes containing silver nitrate are not left in the light.
- 5 The titration must be carried out in the cold. Otherwise the results will be erratic.

**B CHLORIDE IN URINE****Principle Apparatus Reagents**

Same as for serum

**Procedure**

- 1 Collect and store samples as described previously.
- 2 Deliver exactly 0.1 ml of urine by means of an automatic syringe pipette into a 15 ml test tube.
- 3 4 5 Steps 3, 4 and 5 are identical for urine and serum except urine need not be digested in a boiling water bath unless there is good reason to suspect that the urine contains protein.
- 6 7 8 Steps 6, 7 and 8 are identical for urine and serum.
- 9 When the urine contains a high concentration of chloride 2 ml of working standard AgNO may not suffice to precipitate all the chloride. When this is the case the first drop of thiocyanate produces a red color. In such cases add 1 ml more of working standard AgNO, and proceed with the titration. The titration figure is then the final reading minus the reading before the first drop of KCNS.

**Calculation**

- 1 Grams NaCl/100 ml urine =

$$\frac{(\text{ml AgNO} \times \text{Normality AgNO}) - (\text{ml KCNS} \times \text{Normality KCNS})}{(\text{calibration of syringe pipette in ml})} \times 5.85$$

- 2 The use of line charts was described under serum.

**Example**

Working standard AgNO = 0.0139 N

Working standard KCNS = 0.005 N

Syringe pipette contains 0.1 ml

3 ml AgNO had to be used

Titration figure 0.550 ml of KCNS

Therefore gm NaCl/100 ml =

$$\frac{(3 \times 0.0139) - (0.550 \times 0.0055)}{0.100} \times 5.85 = 1.78 \text{ gm/100 ml}$$

#### Precautions

Same as for serum

### Inorganic Phosphorus in Serum Urine Food and Feces

#### N SERUM

#### Reference

Fiske C H and Subbarow Y The Colorimetric Determination of Phosphorus  
J Biol Chem 66 375-400 (Dec) 1925

#### Principle

Inorganic phosphate reacts with molybdate and a blue color is formed by reduction with para-aminonaphthol sulfonic acid

#### Apparatus

- 1 Coleman J Spectrophotometer Model #6
- Electric centrifuge
- 3 A series of pipettes including 0.5 ml, 5.0 ml and two 10 ml graduated in 0.1 ml divisions
- 4 Round bottomed 15 ml centrifuge tubes with #1 rubber stoppers
- 5 Coleman Jr cuvettes 19 x 150 mm
- 6 Pan of water and source of heat for water bath
- 7 An accurately calibrated syringe pipette 1.0 ml
- 8 A quantity of 50 ml volumetric flasks

#### Reagents

- 1 Trichloroacetic acid approximately 10% in water
- 2 Sulfuric acid approximately 10 N In a mixing cylinder carefully pour 100 ml of conc H<sub>2</sub>SO<sub>4</sub> into 900 ml of water After cooling dilute to exactly 300 ml with water
- 3 Molybdate I (for standards and urine) 0.5% ammonium molybdate in 5 N H<sub>2</sub>SO<sub>4</sub> In a one liter mixing cylinder dissolve 0.5 gm of salt in 500 ml of 10 N H<sub>2</sub>SO<sub>4</sub> and dilute to one liter with water
- 4 Molybdate II (for serum) 2.5% ammonium molybdate in 3.3 N H<sub>2</sub>SO<sub>4</sub> In a one liter mixing cylinder dissolve 0.5 gm of the salt in 100 ml of 10 N H<sub>2</sub>SO<sub>4</sub> and dilute to one liter with water
- 5 Para-aminonaphtholsulfonic acid (ANSA) To a 250 ml mixing cylinder add 0.5 gm ANSA and 195 ml of 15% sodium bisulfite (metabisulfite gives the same results) Add 5 ml of 0% sodium sulfite Mix If the ANSA does not dissolve within a half hour add more sulfite 5 ml at a time until solution is complete Store in a brown bottle It keeps for about 2 weeks  
(If the dry ANSA is more than slightly colored brown or purple it must be recrystallized as follows and sealed airtight if necessary)

In approximately one liter of water at about 90 C dissolve 150 gm of sodium bisulfite and 10 gm of sodium sulfite Add about 30 gm ANSA and stand 15 minutes with intermittent shaking Filter by suction while still hot Cool the filtrate and add 10 ml of concentrated HCl Let stand for 10 minutes Filter off the crystals wash once with 300 ml of cold water and twice with 100 ml portions of 95% alcohol Dry in air in the dark and store in a brown bottle.)

## 110 *Biochemical Procedures*

- 6 **Standard phosphate** In a two liter mixing cylinder place exactly 0.351 gm pure  $\text{KH}_2\text{PO}_4$  approximately 20 ml 10 N  $\text{H}_2\text{SO}_4$  and approximately one liter distilled water. When dissolved dilute to exactly two liters. This solution keeps indefinitely in refrigerator. 1 ml equals 0.04 mg P.

### Procedure

- 1 To a 15 ml centrifuge tube add exactly 9.0 ml of 10% trichloroacetic acid.
- 2 Add exactly 1.0 ml serum mix and centrifuge till clear.
- 3 To a cuvette add exactly 5.0 ml of filtrate.
- 4 Add exactly 1.0 ml Molybdate II.
- 5 Then add 0.5 ml ANSA and 3.5 ml of distilled water. Stopper and mix. Read in 10 minutes.
- 6 Prepare a blank using Molybdate I and set at 100% T at 660  $\text{m}\mu$ . Read unknowns.
- 7 Prepare a standard curve using Molybdate I with 0.01 to 0.05 mg P at 0.01 mg intervals.

### Calculation

$$\text{mg P/100 ml serum} = \frac{(\text{mg P in cuvette}) \times 10 \times 100}{(\text{ml aliquot}) \times (\text{ml serum})}$$

### Example

Sample of serum was 1 ml. Aliquot was 5 ml of a 10 ml trichloroacetic acid filtrate. % T was 41.0 equivalent on the calibration curve to 0.02 mg P in the cuvette.

$$\text{mg inorganic P/100 ml serum} = \frac{(0.02) \times 10 \times 100}{(5) \times (1)} = 4.0 \text{ mg P}$$

### Precautions

1. The quality of the ANSA is critical. If there is any doubt recrystallize it before using.

## B IN URINE

### Procedure

- 1 Dilute exactly 1.0 ml of urine to 50 ml with distilled water.
- 2 To a cuvette add 5.0 ml of diluted urine.
- 3 Add 1.0 ml of Molybdate I.
- 4 0.5 ml of ANSA and 3.5 ml of water. Stopper and mix.
- 5 Place blank at 100% T and read the unknowns at 660  $\text{m}\mu$ .

## C IN FOOD AND FECES

### Reference

Official A O A C Method

Corley R C and Denis W. The Determination of Calcium in Tissues, Feces and Milk. *J Biol Chem* 66: 601-608 (Dec) 1925.

### Apparatus

- 1 Steam bath
- 2 Crucibles (porcelain or quartz)
- 3 Muffle furnace at 600 C
- 4 Coleman Spectrophotometer Model #6
- 5 Cuvettes 19 x 150 mm

**Reagents**

- 1 Standard potassium dihydrogen phosphate solution  
Dissolve 0.4394 gm of pure dry  $\text{KH}_2\text{PO}_4$  in water and make up to a liter. 50 ml of this solution is diluted to 200 ml and gives a standard of which 2 ml equals 0.05 mg of P.
- 2 Ammonium molybdate solution  
Dissolve 5 gm of ammonium molybdate in 300 ml of water. Dilute 75 ml of  $\text{H}_2\text{SO}_4$  to 200 ml and add to the ammonium molybdate solution.
- 3 Magnesium nitrate solution  
Dissolve 160 gm of  $\text{MgO}$  in  $\text{HNO}_3$  (1 volume conc  $\text{HNO}_3$  plus 1 volume water), avoiding excess of the acid. Add a little  $\text{MgO}$  in excess. Boil, filter from the excess  $\text{MgO}$ ,  $\text{Fe}_2\text{O}_3$ , etc. and dilute to 1 liter.
- 4 Hydrochloric acid 1 part plus 4 parts of water.

**Preparation of Solution**

- 1 To 1 gm of sample in a small crucible add 1 ml of the magnesium nitrate solution and place on steam bath.
- 2 After a few minutes cautiously add a few drops of  $\text{HCl}$  taking care that formation of gas bubbles does not push portions of sample over edge of crucible.
- 3 Make or 3 further additions of a few drops of  $\text{HCl}$  while sample is on bath so that as it approaches dryness there is a tendency for it to char.
- 4 If contents of crucible become so viscous that no further drying may be obtained on bath, complete drying on a hot plate. Then put on crucible cover, transfer to cold muffle and ignite at dull red heat for at least 6 hours (600°C) or until an even gray ash is obtained. (It may be necessary to cool crucible, dissolve ash in a little water or alcoholic glycerol, evaporate to dryness and return uncovered to muffle for 4-5 hours longer.)
- 5 Cool, take up with  $\text{HCl}$  (1 plus 4) and transfer to 100 ml beaker.
- 6 Add 5 ml of  $\text{HCl}$  and evaporate to dryness on steam bath to dehydrate  $\text{SiO}_2$ .
- 7 Moisten residue with 1 ml of  $\text{HCl}$ , add about 50 ml of water, heat a few minutes on bath, transfer to 100 ml volumetric flask, cool, immediately make to volume and filter, discarding first portion of filtrate.
- 8 From this point on for feces take 3 ml of dilution and for food take 10 ml of dilution.
- 9 The procedure is the same as for phosphorus in urine, steps 2-5.

**D. ALTERNATIVE METHOD FOR SERUM****Reference**

Tschopp E. and Tschopp E. *Helv Chim Acta* 15: 793-809, 1932.

**Principle**

A blue complex forms upon the interaction of phosphates and para-methylamino phenol sulfonic acid. The intensity of the color is measured colorimetrically.

**Apparatus**

- 1 A quantity of 15 ml conical centrifuge tubes.
- 2 A quantity of analytical filter paper and funnels to match.
- 3 Coleman Jr Spectrophotometer Model 6.
- 4 Cuvettes 19 x 150 mm.
- 5 Centrifuge.
- 6 A water bath set at approximately 60°C.
- 7 Syringe pipette 0.1, 1.0 ml and 4.7 ml.

**Reagents**

- 1 Trichloroacetic acid ( $\text{CCl}_3\text{COOH}$ ) 10% filtered. Keep in refrigerator.
- 2 Sodium bisulfite ( $\text{NaHSO}_3$ ) 0% solution in water.



## 112 Biochemical Procedures

- 3 Ammonium molybdate (Merck Reagent Grade 81% MoO<sub>3</sub>) in 1 N sulfuric acid. Dissolve 25 gm phosphate free salt in water and 50 ml 10 N H<sub>2</sub>SO<sub>4</sub>. Dilute to 500 ml. Do not use if blue tinge develops. Keep in wax lined soft glass bottle to guard against decomposition.
- 4 Photorex  
Dissolve in water and dilute to 1 liter  
20 gm p methyl amino phenol sulfonic acid  
2000 gm sodium bisulfite  
100 gm sodium sulfite
- 5 Phosphate standard solution  
Dissolve in water and dilute to 1 liter  
0.35 gm potassium phosphate monobasic  
(pure dry salt) and 100 ml of 10 N sulfuric acid  
1 ml is equivalent to 0.03 mg P

### Procedure

- 1 To a 10 ml cone tipped centrifuge tube add exactly 1 ml serum and 9 ml tri chloroacetic acid
- 2 Mix by inversion and centrifuge for 5 minutes. Filter supernatant fluid through quantitative filter paper
- 3 Pipette 5 ml of filtrate into a cuvette (preferably graduated at 10 ml)
- 4 Add 0.1 ml sodium bisulfite and mix
- 5 Add 0.1 ml Molybdate and mix
- 6 Add 0.1 ml photorex and mix.
- 7 Dilute to exactly 10 ml and mix or add 4.7 ml of water
- 8 Immediately heat for 20 minutes in hot water bath at 60 C. After 20 minutes cool the tube in running water to room temperature
- 9 For a blank make up with water and all reagents as for unknowns. Heat similarly. Set at 100% T at 660 mμ (or 635 mμ). Then read the unknowns

### Standard Curve

Dilute standard in such a way that 5 ml of solution contains 0.05 mg through 0.4 mg of phosphorus

### Calculation

From the phosphate standard curve determine value of galvanometer readings for unknowns in mg of phosphate. As 5 ml of filtrate is equivalent to 0.5 ml serum multiply by 2. To calculate mg of phosphorus per 100 ml multiply by 100

mg P/100 ml serum = (mg P in cuvette obtained from standard table)  $\times 2 \times 100$

### Example

A reading of 41.0% T was equivalent to 0.02 mg of phosphorus on the calibration chart. In the cuvette was the equivalent of 0.5 ml serum (5 ml of a 1 to 10 dilution). Therefore mg P/100 ml serum = (0.02)  $\times 200 = 4.0$

### Precaution

- 1 Be sure to heat all tubes immediately at 60 C and for exactly 20 minutes.

## 11 Sulfur in Urine

### A. TOTAL SULFATE IN URINE

#### Reference

- Folin O. On Sulphur and Sulphur Determinations J Biol Chem 1 131 159 (Jan) 1906

**Principle**

The sulfuric acid of the conjugated sulfates is set free by boiling with acid. The total sulfates are precipitated as barium sulfate and weighed

**Apparatus**

- 1 Porcelain crucibles
- Analytical balance
- 3 Muffle oven set at 600 C
- 4 A desiccator for drying crucibles over calcium oxide
- 5 250 ml Erlenmeyer flasks
- 6 Analytical grade filter paper and funnels

**Reagents**

- 1 Barium chloride ( $\text{BaCl}_2$ ) 5% in distilled water
- 2 Hydrochloric acid ( $\text{HCl}$ ) concentrated.

**Procedure**

- 1 Pipette 25 or 50 ml of urine into a 250 ml Erlenmeyer flask (Less than 25 ml gives erratic results)
- 2 Add 5 ml of conc  $\text{HCl}$  and 10 ml water
- 3 Cover with watch glass boil gently for 10-30 minutes and cool
- 4 Dilute to approximately 150 ml with water
- 5 Add 5 ml of  $\text{BaCl}_2$  (5%) along the side of the flask
- 6 Do not shake or disturb but mix gently by tilting the flask back and forth
- 7 Let stand for 1 hour and filter through analytical grade filter paper
- 8 Wash filter with 50 ml of distilled water and place filter paper plus precipitate into a weighed crucible
- 9 Dry and ignite in muffle oven at 600 C
- 10 Cool in a desiccator and weigh

**Calculation**

- 1 gm S/liter urine =

$$(\text{gm ppt}) \times \frac{\text{atomic wgt S}}{\text{molecular wgt BaSO}_4} \times \frac{1000}{\text{ml urine}}$$

$$= (\text{gm ppt}) \times \frac{32.06}{330.43} \times \frac{1000}{\text{ml urine}}$$

$$= \frac{\text{gm ppt}}{\text{ml urine}} \times 137.5$$

**Example**

- 1 Precipitate of  $\text{BaSO}_4$  1.00 mg
- Urine analyzed 5 ml

$$\text{gm S/liter urine} = \frac{0.100}{5} \times 137.5 = 0.688$$

This computation is facilitated by Table G

**Precautions**

- 1 At least 5 ml of urine have to be used in this method. Otherwise recovery of added sulfate is erratic and usually low

TABLE 6  
SULFUR IN URINEWhen 25 ml of urine are analyzed, gm S/liter urine = (gm ppt)  $\times$  55 or (mg ppt)  $\times$  0.0055

| mg<br>ppt                       | 0    | 2    | 4    | 6    | 8    |
|---------------------------------|------|------|------|------|------|
| GRAMS SULFUR PER LITER OF URINE |      |      |      |      |      |
| 0                               | 0    | 01   | 02   | 03   | 04   |
| 10                              | 06   | 07   | 08   | 09   | 10   |
| 20                              | 11   | 12   | 13   | 14   | 15   |
| 30                              | 17   | 18   | 19   | 20   | 21   |
| 40                              | 22   | 23   | 24   | 25   | 26   |
| 50                              | 28   | 29   | 30   | 31   | 32   |
| 60                              | 33   | 34   | 35   | 36   | 37   |
| 70                              | 39   | 40   | 41   | 42   | 43   |
| 80                              | 44   | 45   | 46   | 47   | 48   |
| 90                              | 50   | 51   | 52   | 53   | 54   |
| 100                             | 55   | 56   | 57   | 58   | 59   |
| 110                             | 61   | 62   | 63   | 64   | 65   |
| 120                             | 66   | 67   | 68   | 69   | 70   |
| 130                             | 72   | 73   | 74   | 75   | 76   |
| 140                             | 77   | 78   | 79   | 80   | 81   |
| 150                             | 83   | 84   | 85   | 86   | 87   |
| 160                             | 88   | 89   | 90   | 91   | 92   |
| 170                             | 94   | 95   | 96   | 97   | 98   |
| 180                             | 99   | 1 00 | 1 01 | 1 02 | 1 03 |
| 190                             | 1 05 | 1 06 | 1 07 | 1 08 | 1 09 |
| 200                             | 1 10 | 1 11 | 1 12 | 1 13 | 1 14 |
| 210                             | 1 16 | 1 17 | 1 18 | 1 19 | 1 20 |
| 220                             | 1 21 | 1 22 | 1 23 | 1 24 | 1 25 |
| 230                             | 1 27 | 1 28 | 1 29 | 1 30 | 1 31 |
| 240                             | 1 32 | 1 33 | 1 34 | 1 35 | 1 36 |
| 250                             | 1 38 | 1 39 | 1 40 | 1 41 | 1 42 |
| 260                             | 1 43 | 1 44 | 1 45 | 1 46 | 1 47 |
| 270                             | 1 49 | 1 50 | 1 51 | 1 52 | 1 53 |
| 280                             | 1 54 | 1 55 | 1 56 | 1 57 | 1 58 |
| 290                             | 1 60 | 1 61 | 1 62 | 1 63 | 1 64 |
| 300                             | 1 65 | 1 66 | 1 67 | 1 68 | 1 69 |
| 310                             | 1 71 | 1 72 | 1 73 | 1 74 | 1 75 |
| 320                             | 1 76 | 1 77 | 1 78 | 1 79 | 1 80 |
| 330                             | 1 82 | 1 83 | 1 84 | 1 85 | 1 86 |
| 340                             | 1 87 | 1 88 | 1 89 | 1 90 | 1 91 |

## B INORGANIC SULFATE IN URINE

## Reference

Fohn O On Sulphate and Sulphur Determination *J Biol Chem.* 1 131 159 (Jan.)  
1906

## Principle

The inorganic sulfates are precipitated at once as barium sulfate

## Apparatus

- 1 A variety of 250 ml Erlenmeyer flasks
- 2 Crucibles (porcelain)
- 3 Desiccator
- 4 Analytical balance
- 5 Analytical filter paper
- 6 Muffle oven—600 C

## Reagents

- 1 Hydrochloric acid dilute One volume to 4 volumes of distilled water
- 2 Barium chloride—5% solution

**Procedure**

- 1 Pipette 5 ml of urine into a 250 ml Erlenmeyer flask
- 2 Add approximately 100 ml of water
- 3 Then acidify with 10 ml of hydrochloric acid.
- 4 Add 10 ml of barium chloride solution slowly drop by drop to the cold solution
- 5 Allow to stand for one hour and filter
- 6 Wash the filter paper with 50 ml of water and then place the filter paper plus the precipitate into a weighed crucible
- 7 Dry and ignite in muffle oven at 500-600°C (overnight)
- 8 Cool in a desiccator and weigh

**Calculation**

1 gm S per liter urine =

$$(\text{gm ppt}) \times \frac{\text{atomic wgt S}}{\text{molecular wgt BaSO}_4} \times \frac{1000}{\text{ml urine}}$$

$$= (\text{gm ppt}) \times \frac{3.06}{33.43} \times \frac{1000}{\text{ml urine}} = \frac{\text{gm ppt}}{\text{ml urine}} \times 137.5$$

**Example**

Precipitate of BaSO<sub>4</sub> = 105.0 mg  
 Urine analyzed = 25 ml  
 gm inorganic S per liter urine =  $\frac{0.1050}{25} \times 137.5 = 0.578 \text{ gm}$

(This calculation is facilitated by Table 6)

**Precautions**

- 1 At least 25 ml of urine have to be used in this method otherwise recovery of added sulfate is erratic and usually low

**ETHEREAL SULFATES****Reference**

Folin O On Sulphate and Sulphur Determination J Biol Chem 1:131-159 (Jan) 1906

**Calculation**

Total S - Inorganic S = Etheral sulfates

**Example**

0.686 gm - 0.578 = 0.108 gm

**2 Minerals Miscellaneous References****Alkaline Earths**

allmann S Quantitative Separation of Calcium Barium and Strontium Anal Chem 20:449-451 (May) 1949

**Aluminum**

trafford N and Wyatt P F The Separation and Determination of Very Small Amounts of Aluminum and Iron in Water Analyst 68:319-344 (Nov) 1943

Golff H Spektrochemische Untersuchungen über den Aluminiumgehalt des Blutes Biochem Ztschr 319:1-8 No 1 1943

**Antimony**

Jaren T H Colorimetric Micro-Determination of Antimony With Rhodamine B Anal Chem 19(7):487-491 (July) 1947

## 116 *Biochemical Procedures*

### Bromine

- Wuth A A Rational Bromide Treatment, *J A M A* 88 2013 2017 (June) 19 7  
Tsyvina, B S Determination of Small Quantities of Calcium Zavodskaya Lab 15  
142 144 1949  
Tyner E H Determining Small Amounts of Calcium in Plant Materials *Anal Chem* 20  
76 80 (Jan) 1948

### Chlorine

- Binkley F A Colorimetric Reaction of Chloride Ion *J Biol. Chem* 173 (1) 403 405  
(March) 1948  
Ellms J W and Hauser S J Ortho Toluidine as a Reagent for the Colorimetric Estima-  
tion of Small Quantities of Free Chlorine *J Indust Engin Chem* 5 915 917 (Nov)  
1913  
Hoffman W S Photoelectric Clinical Chemistry New York, 1941 W Morrow and Com-  
pany pp 129 137  
Schales O and Schales S S A Simple and Accurate Method for the Determination of  
Chloride in Biological Fluids *J Biol Chem* 140 (3) 819 884 (Sept) 1941  
Taylor W A Lamotte Wuth Bromide Comparator *J Lab & Clin Med* 13 49 496  
(Oct March) 19 7 1948  
Stiff, H A Jr The Colorimetric Determination of Blood Chloride by the Iodometric  
Method *J Biol Chem* 172 (2) 695 698 (Feb) 1948  
Swift E H Arcand G M Lutwack R and Meier, D J Determination of Chloride  
*Anal Chem* 22 (2) 306 308 (Feb) 1950  
Whitehorn J C Simplified Method for the Determination of Chlorides in Blood or  
Plasma *J Biol Chem* 40 419-460 (Feb) 1921  
Zimmerman W J and Layton E M Jr Polarographic Micromethod for the Determina-  
tion of Blood Chloride *J Biol Chem* 181 (1) 141 147 (Nov) 1949

### Calcium

- Alten F Weiland H and Knippenberg E Kolorimetrische Calciumbestimmung uber  
das Ickrolonat *Biochem Ztschr* 265 266 85 89 1933  
DeLuca H A Photometric Microdetermination of Calcium *Canad. J Research* 25  
B 449 454 (Sept) 1947  
Klement R Rapid and Simple Determination of Small Quantities of Calcium *Z Anal  
Chem* 128 431 435 1948  
Korenman I M and Gutnik G B Determination of Small Quantities of Calcium  
Zavodskaya Lab 15 136 138 1949  
Murnaghan M F Micro Diffusion Method for the Determination of Calcium in Blood  
*Serum Analyst* 74 546 549 (Oct) 1949  
Ringbom A and Merikanto B A Titrimetric Method for the Determination of Calcium  
*Acta Chem Scandinav* 3 99 31 No 1 1949  
Roe J H and Kahn B S The Colorimetric Determination of Blood Calcium *J Biol  
Chem* 81 18 (Jan) 1929  
Shohl A T The Effect of Hydrogen Ion Concentration Upon the Determination of  
Calcium *J Biol Chem* 50 527 536 (Feb) 1929  
Shohl A T, and Pedley F G A Rapid and Accurate Method for Calcium in Urine  
*J Biol Chem* 50 337 344 (Feb) 1929  
Sobel A E and Sklarsky S A Direct Acidimetric Microtitration Method for Calcium  
*J Biol Chem* 122 665 672 (Feb) 1938  
Trevan J W and Bainbridge H W The Estimation of Calcium in Blood Serum *Bio-  
chem. J* 20 4 3 426 No 6 1926

### Cobalt

- Claassen A and Westerveld W Photometric Determination of Cobalt With Nitroso E  
Salt *Rec Trav Chim* 67 720 724 1948  
McNaught K J The Determination of Cobalt in Animal Tissue Modification of the  
Method *Analyst* 67 97 98 (March) 1949  
Van Klooster H S Nitroso R Salt a New Reagent for the Detection of Cobalt *J Am  
Chem Soc* 43 746 749 (Jan) 1921  
Willard H H and Kaufman S Colorimetric Determination of Cobalt Using Nitroso R  
Salt *Anal Chem* 19 505 (July) 1947

### Copper

- Oglethorpe C C and Smith C G Determination of Copper Volumetrically by the  
Iodine Thiocyanate Method *Analyst* 68 325 3 8 (Nov) 1943

## Iodine

- Barker S B Determination of Protein Bound Iodine J Biol Chem 173 715 74 (April) 1948
- Connor A C Curtis G M and Swenson R E A Simplified Method for the Determination of the Protein Bound Blood Iodine and Its Clinical Application J Clin Endocrinol 9 (11) 1185 1189 (Nov) 1949
- Dobrouz P Monge Hedde M F and Fondarai J Colorimetric Microdetermination of Iodine Bull Soc Chim Biol France pp 898 899 1947
- Fine S D Report on Iodine J Assoc Offic Agric Chemists 32 (3) 555 558 1949
- Gomberg H J Beierwaltes W H and Lampe I Radio Iodine Uptake in Humans V Absolute Determination of Radio Iodine in the Thyroid Proc Soc Exper Biol & Med 73 40 408 (March) 1950
- Gross W G Wood L K and McHargue J S Spectrophotometric Determination of Iodine Anal Chem 20 900 901 (Oct) 1948
- Lacourt A and Timmermans A M Contributions to the Volumetric Micro Method for the Determination of Organic Iodine Anal Chim Acta 1 (2) 140 148 1947
- McCullagh D R A New Method for the Determination of Iodine J Biol Chem 107 (1) 35 44 (Oct) 1934
- Salter W T and McKay E A Iodine in Blood and Thyroid of Man and Small Animals Endocrinol 35 90 99 (Nov) 1944
- Stimmel B F and McCullagh D R A Note Concerning the Determination of Iodine J Biol Chem 116 (1) 1 4 (Nov) 1936

## Iron

- Cooke W D Hazel F and McNabb W M Volumetric Determination of Small Amounts of Iron Anal Chem 21 1011 (Aug) 1949
- Hahn P F Bale W F Hettig R A Hamen M D and Whipple G H Radioactive Iron and Its Excretion in Urine Bile and Feces J Exper Med 70 443 451 (Nov) 1939
- Worwood A J The Determination of Iron in Biological Material Biochem J 41 (1) 39 41 1947

## Magnesium

- Flaschka H Rapid Determination of Calcium and Magnesium Monatsch 80 506 509 1949
- Garner R J Colorimetric Determination of Magnesium in Plasma or Serum by Means of Titan Yellow Biochem J 40 (5/6) 833 831 1946
- Kunkel H O Pearson P R and Schweigert B S The Photoelectric Determination of Magnesium in Body Fluid J Lab & Clin Med 32 (8) 107 1033 194
- Michaels G D Anderson C T Margen S and Kinsell L W Colorimetric Determination of Ca and Magnesium in Small Amounts of Urine Stool and Food J Biol Chem 180 175 180 1949

## Manganese

- Gate, I M and Ellis G H A Micro Colorimetric Method for the Determination of Manganese in Biological Materials With 4,4-Tetramethyldiaminotriphenylmethane J Biol Chem 168 3 544 (May) 1947
- Kun L Microdetermination of Manganese in Biological Material by Means of Catalysis J Biol Chem 170 509 514 (Oct) 1947

## Mixtures

- Boyle A J Whitehead T, Bird E J Batchelor T M, Iseri T, Jacobson S D and Myers G B The Use of the Emission Spectrograph for the Quantitative Determination of Na K Ca Mg and Fe in Plasma and Urine J Lab & Clin Med 31 (5) 63 630 1949
- Polley J P The Microdetermination of Sodium Potassium Calcium and Chloride in a Single Plasma Sample Canad J Research 26 E 183 195 (April) 1948

## Phosphorus

- Brigg A P Some Applications of the Colorimetric Phosphate Method J Biol Chem 59 3 65 (March) 1944
- Kuttner T and Cohen H R Micro Colorimetric Studies I. Molybdenic Acid Stannous Chloride Reagent The Micro Estimation of Phosphate and Calcium in Pus Plasma and Spinal Fluid J Biol Chem 75 517 53 (Nov) 197
- Soyenkov B A Micromethod of Phosphate Determination J Biol Chem 147 447 457 (May) 1947

**Potassium**

- Breh F, and Gaebler O H The Determination of Potassium in Blood Serum *J Biol Chem* 87 (1) 81 90 (May) 1930
- Chen Y C, and Shen S S Determination of Potassium by Zinc Cobaltinitrite *J Chinese Chem Soc* 8 (1) 12 14 1941
- Dahr E Determination of Potassium in Urine by Means of Periodic Acid *Acta physiol Scandinav* 12 229 234, 1946
- Folch J, and Lauren M Estimation of Potassium in Biological Materials as Potassium Phosphotungstate *J Biol Chem* 169 (3) 539 549 (Aug) 1947
- Harris J E A Modified Silver Cobaltinitrite Method for Potassium Determination *J Biol Chem* 136 619 627 (Dec) 1940
- Looney J M and Dyer C G A Photoelectric Method for the Determination of Potassium in Blood Serum *J Lab & Clin Med* 28 (3) 355 363 (Dec) 1947
- Shohl A T and Bennett H B A Micro Method for the Determination of Potassium as Iodoplatinate *J Biol Chem* 78 643 657 (Aug) 1928
- Van Rysselberge P J Quantitative Determination of Potassium by Sodium Cobaltinitrite Method *Indust & Engin Chem Anal Ed* 3 (1) 3 4 (Jan) 1931

**Sodium**

- Albanese, A A and Lein M The Microcolorimetric Determination of Sodium in Human Biologic Fluids *J Lab & Clin Med* 33 246 250 (Feb) 1948
- Bills C E McDonald F G Neidermeyer W and Schwartz M C Sodium and Potassium in Foods and Waters Determination by the Flame Photometer, *J Am Dietet A* 25 (4) 304 314 (April) 1949
- Butler A M and Tuthill E An Application of the Uranyl Zinc Acetate Method for Determination of Sodium in Biological Material *J Biol Chem* 93 171 180 (Sept) 1931
- Hald P M Heinsen A J and Peters J P The Estimation of Serum Sodium From Bicarbonate Plus Chloride *J Clin Investigation* 26 (5) 983 990 (Sept) 1947
- Marinis T P Muirhead F L Jones F and Hill J M Sodium and Potassium Determinations in Health and Disease *J Lab & Clin Med* 32 (10) 1 03 1 16 (Oct) 1947
- McCance R A and Shipp H L A Colorimetric Determination of Sodium *Biochem J* 25 449 456, No 2 1931
- Stone, G C H and Goldzieher J W A Rapid Colorimetric Method for the Determination of Sodium in Biological Fluids and Particularly in Serum *J Biol Chem* 181 511 522 (Dec) 1949
- Weinbach A P A Micromethod for the Determination of Sodium *J Biol Chem* 110 (1) 95 100 (June) 1935
- Willebrands A F Jr Determination of Na and K in Blood Serum and Urine by Means of the Flame Photometer *Nederl tijdschr v geneesk* 93 III 2351 2356 1949

**Sulfur**

- Erdos J Procedure for the Micro Determination of Sulphur *Mikrochemie ver Mikrochim Acta* 34 286 288 1949
- Kahn B S and Leiboff S L Colorimetric Determination of Inorganic Sulphate in Small Amounts of Urine *J Biol Chem* 80 623 630 (Dec) 1928
- Khimenko V G Microdetermination of Sulphur in Biological Materials *Biokhimiya* 14 14, 1949
- Leclerc R Volumetric Determinations of Sulphate *Mém Services Chim état* 37 15 19 1945
- Letonoff, T V and Reinhold J G A Colorimetric Method for the Determination of Inorganic Sulfate in Serum and Urine *J Biol Chem* 114 (1) 147 156 (May) 1936
- Power, M H and Wakefield E G A Volumetric Benzidine Method for the Determination of Inorganic and Etheral Sulphate in Serum *J Biol Chem* 123 665 678 (May) 1938
- Shu Chuan Liang Y W and Tzu Cheng C Gravimetric Determination of Sulphur as Barium Sulphate II. Elimination of Chromic Ion Interference by Means of Tannin *J Chinese Chem Soc* 14 9 16 1946
- Sperber I A Direct Turbidimetric Method for Determining Etheral Sulfates in Urine *J Biol Chem* 172 (2) 441 444 (Feb) 1948
- Steyermark A Bass E and Littman B Microdetermination of Sulfur in Organic Compounds *Anal Chem* 20 587 (June) 1948
- Stockholm M and Koch F C A Quantitative Method for the Determination of Total Sulfur in Biological Material *J Am. Chem Soc* 45 (9) 1953 1959 (Aug) 19 3
- Wakefield E G The Colorimetric Determination of Total and Inorganic Sulphates in Blood Serum Urine and Other Body Fluids *J Biol Chem* 81 713 711 (March) 19 9

**Total Base**

Fiske C H. A Method for the Estimation of Total Base in Urine J Biol Chem 51 (1) 55 61 (March) 1922

Stadie W C and Ross E C A Micro Method for the Determination of Base in Blood and Serum and Other Biological Materials J Biol Chem. 65 735 54 (Oct) 19 5

**Water**

Sapirstein L A A Simple Method for the Densitometric Determination of Heavy Water J Lab & Clin Med 35 793 794 (May) 1950

**Zinc**

Hoch F L and Vallee B L Precipitation by Trichloroacetic Acid as a Simplification in Determination of Zinc in Blood and Its Components J Biol. Chem 181 (1) 295 306 (Nov) 1949

Vallee B I and Gibson J G 2nd An Improved Dithizone Method for the Determinations of Small Quantities of Zinc in Blood and Tissue Samples J Biol Chem 176 435 443 (Oct) 1948



## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### B PROTEIN AND NITROGEN COMPOUNDS

##### 1 Total Nitrogen (Micro Kjeldahl) in Serum Urine, Food and Feces

###### A. IN SERUM AND URINE

###### References

- Keys A     A Rapid Micro Kjeldahl Method, *J Biol Chem* 132 181 187 (Jan.) 1940
- Ma T S and Zuazaga G     Micro Kjeldahl Determination of Nitrogen A New Indicator and an Improved Rapid Method, *J Ind & Eng Chem (Anal Ed)* 14 280 282 (March) 1942

###### Principle

Serum and urine are digested with a phosphate catalyst plus hot sulfuric acid and converted to ammonium sulfate. The ammonia is liberated by the addition of alkali, distilled into a boric acid solution and titrated directly with standard acid.

###### Apparatus

- 1 Calibrated syringe pipette, 0.2 ml (Krogh type)
- 2 A quantity of 100 ml Kjeldahl (micro) flasks with standard tapered joints
- 3 A manifold of six micro burners with a good hood
- 4 Micro Kjeldahl ammonia apparatus (Procured at Macalaster Bixnell Co Cambridge Mass. See Fig 18)
- 5 A 2 ml burette (Krogh syringe)
- 6 A number of 15 ml round bottom test tubes (Pohberg type)
- 7 A number of 100 ml volumetric flasks (for food and fecal analyses)
- 8 Luer lok syringe 5 ml

###### Reagents

- 1 Mixed indicator     Prepare in 95% alcohol 10 ml bromocresol green (0.1%) and 2 ml methyl red (0.1%). This solution keeps well at room temperature
- 2 Boric acid 2% (approximately)     Stored in a glass stoppered bottle it keeps indefinitely
- 3 Hydrochloric acid or sulfuric acid 0.3 N accurately titrated
- 4 Sodium hydroxide 30% (approximately)     Store in bottle with rubber stopper. Keeps indefinitely
- 5 Catalyst     Powder and mix potassium phosphate, dibasic ( $K_2HPO_4$ ) and hydrated copper sulfate ( $CuSO_4$ ) in the proportion (approximately) 4:1
- 6 Concentrated sulfuric acid C.P. (A.C.S.)
- 7 Grease for glass joints     Use pure glycerol or a mixture of 50 ml glycerol with 5 gm NaOH
- 8 Standard ammonium sulfate     Dilute exactly - 30 gm of ammonium sulfate to 1000 ml of water. This is equivalent to 0.5 mg N/ml




Fig 18.—Still for distillation of ammonia in a partial vacuum, showing complete apparatus in use during distillation. (It produced with permission of the manufacturers.)

**Procedure**

- 1 Pipette exactly 0.2 ml of serum or urine with a syringe pipette into digestion flask, taking care to deliver to the very bottom of the flask. Wash out pipette with distilled water 6 times
- 2 Add roughly 2 gm of catalyst (phosphate mixture)
- 3 Add approximately 1.5 ml of concentrated sulfuric acid.
- 4 Digest over burners with exhaust or in hood until solution becomes a clear green. Let flasks cool gradually to room temperature
- 5 Pipette 2 ml of 2% boric acid into collecting tubes (15 ml)
- 6 Add one drop of mixed indicator to each tube
- 7 Add approximately 10 °0 ml distilled water down side of digestion flask.
- 8 Distill contents of flask into boric acid contained in Rehberg tube as follows
  - a Lubricate mouth of flask with glycerol—NaOH and attach to still
  - b Attach Rehberg tube containing boric acid to still so that delivery tube is just below surface of boric acid
  - c Turn on vacuum slowly
  - d When pressure is uniform throughout still add approximately 6 ml of sodium hydroxide to flask by opening cock under burette
  - e Wait 1 to 2 minutes until the evolution of gas decreases (when the rate of bubbling in the acid tube will decrease) and turn on full vacuum
  - f Light micro burner under flask with very low flame for the first few minutes gradually increasing the flame until it touches the flask
  - g Collect distillate. At intervals lower boric acid tube so that delivery tip remains only ½ inch below surface of distillate
  - h Collect 9-10 ml distillate
    - i. Pull acid tube below tip of delivery tube
    - j Distill 5-6 more drops and turn off gas flame
    - k In 15-20 seconds turn off vacuum
    - l Wait for bubbling to quiet down and open capillary cock slowly
  - m Allow time for pressure within and without still to equalize
  - n Remove collecting Rehberg tube and flask.
  - o Remove digestion flask by holding its bulb with a towel. Be sure to drain off excess alkaline solution on tip of glass tube which inserts into the flask.
- 9 Titrate contents of micro tube with 0.3 N  $H_2SO_4$  until solution first turns a distinct orange. Use micro burette (2 ml) and air pressure for stirring

**Blank**

Digest and titrate as above a solution of distilled water with 1 gram digest mixture and 1 ml of acid.

**Standard**

Place 5 ml of standard ammonium sulfate  $(NH_4)_2SO_4$  in a digestion flask and proceed as usual

**Calculation**

- 1 Grams N per 100 ml of serum or urine =

$$\frac{0.014 \times (\text{ml titration} - \text{ml blank}) \times (\text{Normality of acid}) \times 100}{\text{Volume of sample in ml}}$$

- 2 To convert nitrogen to protein in serum multiply by the factor of 6.25

**Example**

- 1 Normality of standard acid = 0.335 N

$$(\text{Equivalent to } 14 \times \frac{1}{1000} \times 0.335 \text{ gm nitrogen per ml})$$

- 2 Amount of serum used in analyses was 0.231 ml
- 3 ml of standard acid used in titration was 0.560 ml
- 4 Blank analysis = 0.02 ml acid.
- 5 Therefore

$$\frac{0.14 \times (0.560 - 0.02) \times (0.5) \times 100}{0.231 \text{ ml}} = 1.09 \text{ gm nitrogen per 100 ml}$$

- 6 Then
- 1.09 gm nitrogen/100 ml  $\times 6.5 = 6.8$  gm protein/100 ml

#### Precautions

- 1 Never allow the distillate tap to be above the boric acid level
- 2 In digesting the sample do not allow the sample to dry. If necessary add more concentrated sulfuric acid
- 3 Rinse out the 0.2 ml syringe pipette at least 6 times with water adding all washings to the flask
- 4 Be sure that enough alkali is added to liberate the ammonia. This can be determined by a change in color from green to blue or brown.
- 5 The catalyst recommended above serves for most biological fluids. If unusual compounds might be present it is necessary to check the efficiency of the catalyst and perhaps to use selenium or mercury as catalyst

### B NITROGEN IN FECES AND FOOD

- 1 Once the food and feces are diluted with water and homogenized in a Waring blender it is convenient to pipette out 0.2 ml for digesting
- 2 We use a calibrated 5 ml Luer Lok syringe for measuring as it has a larger bore than the regular syringe
- 3 Add 0.2 ml of concentrated  $\text{H}_2\text{SO}_4$  and 4 gm of the catalyst
- 4 After digestion to a clear green the flask is cooled and diluted very slowly to 100 ml (actually only concentrated  $\text{H}_2\text{SO}_4$  is left in the flask so be very cautious)
- 5 Then distill either 10 or 20 ml of solution

### 2 Serum Proteins (Biuret Method)

#### Reference

Wolfson W Q, Cohn C, Calvary E and Ichiba F. Studies in Serum Proteins. V. A Rapid Procedure for the Estimation of Total Protein, True Albumin, Total Globulin, Alpha Globulin, Beta Globulin and Gamma Globulin in 1.0 ml of Serum. *Am J Clin Path* 18, 3730 (Sept.) 1943

#### Principle

The albumin and the three major globulin fractions are determined after precipitation and measured colorimetrically in a spectrophotometer by the addition of an alkaline copper sulfate reagent (biuret) which produces a bluish purple color.

#### Apparatus

- 1 Cuvettes 19  $\times$  150 mm
- 2 Coleman Jr Spectrophotometer Model #6
- 3 Centrifuge tubes—15 ml round bottom
- 4 Centrifuge tubes—15 ml conical
- 5 Centrifuge—electric
- 6 Incubator set at 37°C
- 7 Calibrated 0.1, 0.2, 0.4, 0.6 and 0.8 ml syringe pipettes.
- 8 Rubber stoppers #1 solid.

**Reagents**

- 1 23% sodium sulfate Exactly 23 grams of anhydrous sodium sulfate in distilled water at 37 C Make up to 100 ml and keep in an incubator at 37 C
- 2 28% sodium sulfite solution Exactly 28 gm in water at 28 C Shake and dissolve
- 3 Saline ammonium sulfate In a liter volumetric flask dissolve 193 gm of ammonium sulfate in approximately 500 ml of water Add 40 gm of sodium chloride (NaCl) dissolve and make up to 1000 ml
- 4 Biuret reagent (Weichselbaum T E *Am J Clin Path 16 Tech. Sec 10* 40 49 1946) Dissolve 90 gm of Rochelle salts in about 400 ml of 0.2 N NaOH. Add 10 gm of  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  When dissolved add 10 gm of potassium iodide and make up to 2000 ml with 0.2 N NaOH Store in rubber stoppered waxed glass bottle
- 5 Span ether reagent (We are indebted to the Atlas Powder Co of Wilmington Delaware for Span )
 

|                   |   |   |
|-------------------|---|---|
| 1 ml Span ~ 20    | } | Filter through fast paper and dilute to 100 ml with ether |
| 99 ml Ether (USP) |   |   |
- 6 Ether USP

**Procedure for Total Protein**

For a standard curve use saline as the diluent and use protein concentrations of 10 to 80 gm per 100 ml at intervals of one gram Use serum albumin concentration of 25 gm per 100 ml

- 1 Pipette exactly 0.1 ml of serum into a 15 ml round bottom tube and dilute to 5 ml with 4.9 ml of distilled water Mix well by inversion
- 2 Transfer exactly 3 ml to a cuvette
- 3 Add exactly 3 ml of biuret reagent and mix by inverting the tube
- 4 Prepare a blank with 3 ml of water and 3 ml of biuret reagent Save the same blank for determination of albumin plus alpha globulin and for gamma globulin.
- 5 Let stand for at least 30 minutes and read at 540  $m\mu$  placing the blank at 100% T The color is stable for 24 hours

**Procedure for Serum Albumin Plus Alpha Globulin**

- 1 Pipette exactly 2.4 ml of sodium sulfate solution into a test tube This solution is always kept in an incubator at 37 C
- 2 Add exactly 0.1 ml of serum and mix by inversion *Do not shake*
- 3 Add approximately 1 ml of ether stopper and mix by inversion for 30 seconds Centrifuge for 5 to 10 minutes at 1500 to 2000 r p m
- 4 Hold tube on a slant and withdraw exactly 1.5 ml of clear centrifugate and transfer to a cuvette
- 5 Add exactly 1.5 ml of distilled water and 3.0 ml of biuret reagent Mix well.
- 6 Read after 30 minutes at 540  $m\mu$ , setting blank at 100% T

**Procedure for Serum Albumin**

- 1 Pipette exactly 4.9 ml of sodium sulfite solution into a 15 ml round bottom tube.
- 2 Add exactly 0.1 ml of serum and mix thoroughly by inversion
- 3 Add about 1 ml Span ether reagent stopper and invert gently for 30 seconds *Do not shake* In some cases up to 30% of the albumin is lost by vigorous shaking Centrifuge as above
- 4 Pipette exactly 3.0 ml of the clear centrifugate into a cuvette
- 5 Add 3 ml of biuret reagent and mix well by inversion Prepare a blank by using 3 ml sodium sulfite plus 3 ml biuret reagent
- 6 Let stand 30 minutes and read at 540  $m\mu$  setting the blank at 100% T

## Procedure for Serum Gamma Globulin

- 1 Pipette exactly 9.8 ml of saline ammonium sulfate solution into 15 ml conical tube
- 2 Layer exactly 0.2 ml of serum on top of the solution. Stopper and mix by careful slow repeated inversion for about a minute or two until the turbidity has reached a maximum
- 3 Remove exactly 1.0 ml of the mixture and discard.
- 4 Cork and centrifuge at 250 to 450 rpm for 0 minutes. If turbid cool under tap water for a few minutes and centrifuge again. Supernatant fluid must be clear
- 5 Decant supernatant fluid and centrifuge for 5 minutes more. Invert tubes on a layer of paper toweling or filter paper for a few minutes
- 6 First pipette exactly 3.0 ml of biuret reagent to the precipitate and 3.0 ml water and shake briskly for 30 seconds
- 7 Allow to stand for 15 minutes centrifuge down any slight residual turbidity remaining and decant into a cuvette
- 8 For the blank use 3.0 ml biuret reagent plus 3.0 ml water. Allow to stand 30 minutes and read at 540 m $\mu$  with the blank at 100% T
- 9 Divide the protein value obtained by 3 to obtain the serum gamma globulin concentration. This is the reason for discarding the 1 ml in step 3

## Calculation

Serum total globulin = Total protein - Serum albumin

A/G Ratio = Serum albumin  $\div$  Serum globulin

Serum alpha globulin = (Serum albumin + alpha globulin) - Serum albumin

Serum beta + gamma globulin = Total protein - (albumin + alpha globulin)

Serum beta globulin = (Serum beta + gamma globulin) - Serum gamma globulin

## Example

- 1 Total protein. A reading of 5% T on the calibration chart is equivalent to 4.16 mg of protein in actual amount. In the final analyses we actually used 0.06 ml of serum (3 ml of a 0.1 to 5.0 dilution)  
Therefore 4.16 mg = 0.06 ml of serum  
or 69 gm of protein per 100 ml of serum.
- Albumin. A reading of 66% T is equivalent to 2.60 mg of albumin in actual amounts. We used 0.06 ml of serum (3 ml of a 0.1 to 5 ml dilution)  
Therefore 2.60 mg = 0.06 ml of serum  
or 43 gm albumin per 100 ml of serum
- 3 Alpha globulin plus albumin. A reading of 61% T is equivalent to 3.1 mg of albumin + alpha globulin in actual amounts. We used 0.06 ml of serum (15 ml of a 0.1 to 5 ml dilution)  
Therefore 3.1 mg = 0.06 ml of serum  
or 5.1 gm of albumin plus alpha globulin per 100 ml of serum
- 4 Serum gamma globulin. A reading of 77% T is equivalent to 1.56 mg gamma globulin in actual amounts. We used 0.18 ml of serum (9 ml of a 0.2 to 10 ml dilution)  
Therefore 1.56 = 0.18 ml serum  
or 0.87 gm of gamma globulin per 100 ml of serum.
- 5 Total protein (69 gm) - albumin (43 gm) = total globulin (26 gm)  
Therefore albumin/globulin ratio = 1.6
- 6 Alpha globulin = (Serum albumin + alpha globulin 5.4 gm) - albumin (4.3) = 1.1 gm.
- 7 Beta + gamma globulin = Total protein (69) - albumin + alpha globulin (5.4) = 23.6 gm
- 8 Beta globulin = Beta + gamma globulin (23.6) - gamma globulin (0.87) = 22.7 gm

**Precautions**

- 1 In the albumin determination, do not shake the tube vigorously. In some cases up to 30 per cent of the albumin is lost in this way.
- 2 This method depends on exact pipetting. Therefore we have used syringe pipettes throughout the procedure.
- 3 In the serum gamma globulin add the biuret reagent *first* before adding the water.
- 4 The serum albumin must be analyzed by the Kjeldahl method for exact analysis.

**3 Albumin and Globulin in Serum****Reference**

Howe P E The Use of Sodium Sulfate as the Globulin Precipitant in the Determination of Proteins in Blood J Biol Chem 49 93 107 (Nov) 19 1

**Principle**

Serum proteins are fractionated by precipitation of globulin in sodium sulfate. The albumin in solution is analyzed by one of the total nitrogen methods.

**Apparatus**

- 1 An incubator set at 37 C
- 2 Filter paper analytical grade (Whatman No 40)
- 3 Ammonia still and standard tapered flasks (see total nitrogen analyses)
- 4 Nitrogen digestion unit
- 5 An accurately calibrated 0.5 ml syringe pipette
- 6 A quantity of 50 ml Erlenmeyer flasks
- 7 A supply of funnels

**Reagents**

- 1 Sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) anhydrous a 22% solution in water. Keep in 37 C oven.
- 2 Thymol crystals
- 3 Other reagents are the same as those for *Nonprotein Nitrogen* and *Total Nitrogen*.

**Procedure**

- 1 Pipette exactly 0.5 ml of serum and 15 ml of sodium sulfate solution into a 50 ml Erlenmeyer flask.
- 2 Add 1 crystal of thymol.
- 3 Place in incubator at 37 C for 3 hours.
- 4 Heat a funnel, a 50 ml Erlenmeyer flask and filter paper (Whatman No 40) by placing in the same oven.
- 5 Filter solution through heated filter into heated flask to prevent precipitation of sodium sulfate. The filtrate contains albumin and nonprotein nitrogen. If the filtrate is not clear repeat filtration using the same filter paper.
- 6 Digest exactly 5 ml of the filtrate for nitrogen determination using *Keys* method.
- 7 Estimate nonprotein nitrogen in the original sample of serum following *Daly's* method as described in the description of nonprotein nitrogen.

**Calculation**

$$\begin{aligned} \text{Total protein} &= (\text{Total N} - \text{NPN}) \times 6.25 \\ \text{Albumin} &= (\text{N in filtrate} - \text{NPN}) \times 6.25 \\ \text{Globulin} &= \text{Total protein} - \text{albumin} \\ \text{A/G ratio} &= (\text{Albumin} \div \text{Globulin}) \end{aligned}$$

**Example**

Total serum nitrogen = 1.05 gm. per 100 ml serum

Nonprotein nitrogen = 40.0 mg per 100 ml serum

Total nitrogen in filtrate = 0.67 gm per 100 ml serum

Total protein =  $(1.05 - 0.04) \times 6.25 = 6.32$  gm/100 ml

Albumin =  $(0.67 - 0.04) \times 6.25 = 3.94$  gm/100 ml

Globulin =  $6.32 - 3.94 = 2.38$  gm/100 ml

A/G ratio =  $3.94 \div 2.38 = 1.66$

**Precautions**

- 1 Time and temperature are critical. Be sure to incubate for exactly 3 hours at exactly 37 C
- 2 Be sure to have filter paper funnel and flask at exactly 37 C

**Whole Blood Hemoglobin and Serum Protein (Copper Sulfate Specific Gravity)****References**

- Phillips R A Van Slyke D D Hamilton P B Dole V P Emerson K. Jr and Archibald R M Measurement of the Specific Gravities of Whole Blood and Plasma by Standard Copper Sulphate Solutions J Biol Chem. 183 305 330 (March) 1950
- Van Slyke D D Hüller A Phillips R A Hamilton P B, Dole V P Archibald R M and Eder H A The Estimation of Plasma Protein Concentration from Plasma Specific Gravity J Biol Chem 183 331 343 (March) 1950
- Van Slyke D D Phillips R A Dole V P Hamilton P B Archibald R M and Plazin J Calculation of Hemoglobin from Blood Specific Gravities J Biol Chem 183 349 360 (March) 1950

**Principle**

The concentration of protein in serum and hemoglobin in cells determines in large part their specific gravities the specific gravity of whole blood is accounted for by that of serum and cells. The present method estimates specific gravity by the behavior of drops of serum and whole blood in solutions of copper sulfate of various specific gravity. The serum protein and whole blood hemoglobin are then calculated or more simply read off from a line chart.

**Apparatus**

- 1 30 square oz bottles of copper sulfate solution from sp gr 1.021 through 1.033 and 1.05 through 1.068. These bottles should have their sp gr etched clearly on them.
- Four 0.5 ml tuberculin syringes.
- 3 Four hypodermic needles #22 or finer
- 4 Four stylets for cleaning needles.
- 5 Heparinized vials (8 ml) for collecting whole blood.

**Reagents**

- 1 Copper sulfate solid  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  finely granular or powdered
  - a. Copper sulfate solution sp gr 1.100. Allow 1 liter for every 500 samples of blood. It is most conveniently carried in large amounts from the base laboratory into the field but can be made up as follows. Four pounds of copper sulfate are placed in a gallon jug. About 2.5 liters distilled or rain water are added. (Tap water may be used if it is not more than sp gr 1.003 compared with distilled water as 1.000 at the same temperature. Use urinometer for this comparison.) The bottle is stoppered and shaken vigorously for 5 minutes at



the end of which time the temperature is recorded to the nearest  $\frac{1}{2}$  degree C. The saturated supernatant solution is immediately decanted and is filtered through cotton into a clean dry jug. The solution of sp gr 1.100 is now made at once by diluting accurately to 1 liter the amount of saturated solution indicated in Table 7. Once made the solution of sp gr 1.100 keeps in definitely. Prepare 4 liters at a time.

TABLE 7

ML OF SATURATED COPPER SULFATE SOLUTION TO BE DILUTED TO 1 LITER TO MAKE STOCK SOLUTION OF SPECIFIC GRAVITY 1.100

Temperature in C or F refers to the temperature of the saturated solution at the time of saturation (end of shaking for 5 minutes)

| TEMPERATURE |      |     | TEMPERATURE |      |     | TEMPERATURE |       |     |
|-------------|------|-----|-------------|------|-----|-------------|-------|-----|
| C           | F    | ml  | C           | F    | ml  | C           | F     | ml  |
| 10.0        | 50.0 | 578 | 20.0        | 68.0 | 483 | 30.0        | 86.0  | 455 |
| 10.5        | 50.9 | 573 | 20.5        | 68.9 | 484 | 30.5        | 86.9  | 453 |
| 11.0        | 51.8 | 568 | 21.0        | 69.8 | 480 | 31.0        | 87.8  | 450 |
| 11.5        | 52.7 | 563 | 21.5        | 70.7 | 477 | 31.5        | 88.7  | 447 |
| 12.0        | 53.6 | 558 | 22.0        | 71.6 | 473 | 32.0        | 89.6  | 444 |
| 12.5        | 54.5 | 553 | 22.5        | 72.5 | 469 | 32.5        | 90.5  | 441 |
| 13.0        | 55.4 | 548 | 23.0        | 73.4 | 466 | 33.0        | 91.4  | 439 |
| 13.5        | 56.3 | 543 | 23.5        | 74.3 | 463 | 33.5        | 92.3  | 436 |
| 14.0        | 57.2 | 539 | 24.0        | 75.2 | 460 | 34.0        | 93.2  | 433 |
| 14.5        | 58.1 | 534 | 24.5        | 76.1 | 456 | 34.5        | 94.1  | 431 |
| 15.0        | 59.0 | 529 | 25.0        | 77.0 | 453 | 35.0        | 95.0  | 428 |
| 15.5        | 59.9 | 525 | 25.5        | 77.9 | 450 | 35.5        | 95.9  | 425 |
| 16.0        | 60.8 | 521 | 26.0        | 78.8 | 447 | 36.0        | 96.8  | 423 |
| 16.5        | 61.7 | 516 | 26.5        | 79.7 | 445 | 36.5        | 97.7  | 420 |
| 17.0        | 62.6 | 512 | 27.0        | 80.6 | 442 | 37.0        | 98.6  | 418 |
| 17.5        | 63.5 | 508 | 27.5        | 81.5 | 439 | 37.5        | 99.5  | 415 |
| 18.0        | 64.4 | 504 | 28.0        | 82.4 | 436 | 38.0        | 100.4 | 413 |
| 18.5        | 65.3 | 500 | 28.5        | 83.3 | 434 | 38.5        | 101.3 | 410 |
| 19.0        | 66.2 | 496 | 29.0        | 84.2 | 431 | 39.0        | 102.2 | 408 |
| 19.5        | 67.1 | 492 | 29.5        | 85.1 | 428 | 39.5        | 103.1 | 405 |

b Copper sulfate sp gr 1.100. An alternative to method 2 a is to dissolve exactly 170 gm of pure  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in water. Amounts of water recommended for preparing the solution sp gr 1.100

| Temperature C | ml of Water |
|---------------|-------------|
| 10.17 C       | 1004        |
| 17.22 C       | 1005        |
| 26.26 C       | 1006        |
| 26.29 C       | 1007        |
| 29.33 C       | 1008        |
| 33.36 C       | 1009        |
| 36.38 C       | 1010        |
| 38.40 C       | 1011        |

3 Copper sulfate solutions. These solutions are most conveniently made according to Table 8 by running the stock solution of sp gr 1.100 from a burette into a 50 ml graduated mixing cylinder, diluting to the mark with water, mixing, adding to the appropriate bottle, rinsing the cylinder and starting over. These solutions should be changed when 50 drops of whole blood or serum have been added. It is therefore advisable to make up in the field about 500 ml each of those most commonly used.

## Procedure

- 1 The amounts needed are about 0.5 ml each of whole blood and serum prepared as described in a previous section. The whole blood must be mixed thoroughly by inversion of its container immediately before use. The whole blood serum and solutions of copper sulfate must all be at approximately the same temperature.

TABLE 8

TABLE OF STOCK COPPER SULFATE SOLUTION OF GRAVITY 1.0000 TO BE DILUTED TO 50 ML TO PREPARE STANDARD SOLUTIONS OF GRAVITY G TO WITHIN 0.0001

| CuSO <sub>4</sub> IN MIXING CYLINDER |       | CuSO <sub>4</sub> IN MIXING CYLINDER |       |
|--------------------------------------|-------|--------------------------------------|-------|
| SP GR                                | ML    | SP GR                                | ML    |
| 1.009                                | 3.67  | 1.043                                | 1.00  |
| 9                                    | 4.16  | 44                                   | 21.50 |
| 10                                   | 4.66  | 45                                   | 00    |
| 11                                   | 5.15  | 46                                   | 0.50  |
| 1                                    | 5.6   | 47                                   | 03.00 |
| 13                                   | 6.14  | 48                                   | 23.50 |
| 14                                   | 6.64  | 49                                   | 24.00 |
| 15                                   | 7.13  | 50                                   | 24.50 |
| 16                                   | 7.63  | 51                                   | 5.00  |
| 17                                   | 8.1   | 52                                   | 5.50  |
| 18                                   | 8.6   | 53                                   | 6.00  |
| 19                                   | 9.11  | 54                                   | 6.50  |
| 0                                    | 9.61  | 55                                   | 07.00 |
| 1                                    | 10.10 | 56                                   | 7.50  |
|                                      | 10.60 | 57                                   | 8.00  |
| 3                                    | 11.09 | 58                                   | 28.50 |
| 04                                   | 11.9  | 59                                   | 9.00  |
| 5                                    | 1.0   | 60                                   | 9.50  |
| 6                                    | 1.5   | 61                                   | 30.00 |
| 7                                    | 13.00 | 6                                    | 30.50 |
| 08                                   | 13.54 | 63                                   | 31.00 |
| 09                                   | 14.03 | 64                                   | 31.50 |
| 30                                   | 14.5  | 65                                   | 3.00  |
| 31                                   | 1.01  | 66                                   | 3.50  |
| 3                                    | 15.50 | 67                                   | 33.00 |
| 33                                   | 16.00 | 68                                   | 33.5  |
| 34                                   | 16.50 | 69                                   | 34.04 |
| 35                                   | 17.00 | 0                                    | 34.6  |
| 36                                   | 1.50  | 1                                    | 35.08 |
| 37                                   | 18.00 |                                      | 35.60 |
| 38                                   | 18.50 | 73                                   | 36.1  |
| 39                                   | 19.00 | 4                                    | 36.64 |
| 40                                   | 19.50 | 75                                   | 37.15 |
| 41                                   | 0.00  |                                      |       |
| 4                                    | 0.50  |                                      |       |

- 2 For whole blood the range 1.05 through 1.009 is used for serum the range 1.01 through 1.033. Fill a syringe with its needle attached to about the 0.2 ml mark. If the syringe has just been used for a preceding estimation, empty it so far as practicable by pulling the plunger back and forth several times before finally filling the syringe with the new sample.
- 3 With the copper sulfate rack in a convenient position, remove the stopper from sp gr 1.060 for whole blood and 1.08 for serum. Hold the syringe in an almost vertical position with the tip of the needle about one-half inch above the solution.

# LINE CHART FOR CALCULATING PLASMA PROTEINS, HEMOGLOBIN AND HEMATOCRIT FROM GRAVITIES OF PLASMA AND BLOOD




Fig. 19—(From Copper Sulfate Method for Measuring Specific Gravities of Whole Blood and Plasma by Robert A. Phillips, Donald D. Van Slyke, Inc. at P. Dole, Kendall Emerson, Jr. Paul B. Hamilton, Reginald M. Archibald, (Reproduced by permission of the author))

4. Deliver a drop from the syringe keeping steady pressure on the plunger. For satisfactory estimation the drop must break the surface and fall at least an inch into the solution without leaving a streamer attached to the surface.
5. Watch the progress of the drop in the solution. When its specific gravity is exactly the same as that of the solution the drop will fall for a few seconds, come to a complete halt without moving up or down for about 10 seconds, and will then slowly fall to the bottom because it will take up copper and become heavy. When its specific gravity is greater than that of the solution the drop will continue to fall to the bottom. If the drop is lighter than the solution it will fall for a few seconds, halt for a second or two, and begin to rise. It will continue to rise until it takes up enough copper to fall to the bottom. If a drop rises on one bottle and falls in the next lightest the specific gravity lies between those of the two bottles and must be estimated to the nearest 0.000 by noting carefully the relative rates of rise in one bottle and fall in the next.
6. Record the specific gravity to the nearest 0.000 and proceed to the next sample.

### Calculations

Calculate serum protein in grams per 100 ml and hemoglobin in grams per 100 ml by referring to the line chart Fig 19.

To find plasma protein use the two scales on the left hand line. (Note that for plasma the scales read from bottom to top and for the others from top to bottom.) To find hemoglobin or hematocrit from  $G_a$  and  $G$  stretch thread between the numbers on the two outside lines and read the answer on the inside line. To approximate hemoglobin or hematocrit from  $G$  alone stretch the thread between observed  $G$  and mean normal  $G$  and read from the middle line. The ranges for normal values are as follows:

|                              | Men             | Women           |
|------------------------------|-----------------|-----------------|
| Whole blood specific gravity | 1.055 to 1.0635 | 1.053 to 1.0603 |
| Plasma specific gravity      | 1.040 to 1.050  | 1.040 to 1.050  |

When oxalate is used for an anticoagulant a correction must be made. For each mg of 3 ammonium oxalate potassium oxalate added per ml of blood deduct 0.0004 from both  $G_a$  and  $G$ .

### Example

Serum sp gr = 1.064

Whole blood sp gr = 1.059

From line chart serum protein = 7.35 gm per 100 ml and  
hemoglobin = 15.9 gm per 100 ml

### Precautions

1. The samples and bottles of copper sulfate must be at approximately the same temperature otherwise very serious errors can be made.  
Thorough mixing of the whole blood is essential.
2. When particles remain on the surface or when the drops fail to break the surface cleanly the surface should be cleaned with a wooden applicator.
3. In the absence of tuberculin syringes ordinary medicine droppers can be used.
4. The ranges of specific gravity suggested above are for normal soldiers in whom hypoproteinemia and anemia are extremely rare. For other populations other ranges might be necessary.
5. The calculations for hemoglobin and hematocrit are valid only for normochromic normocytic blood.

## EQUATIONS ON WHICH THE LINE CHARTS ARE BASED

(This is copied from the revised method, Don Baxter, Inc, Glendale California whose permission has been given to reprint it here)

The equations given below may be used in place of the line charts to calculate results

*Total plasma protein concentration from plasma specific gravity* The equation used is

$$(1) P = 369 (G_p - 1.0065)$$

P is the plasma protein concentration in grams per 100 ml  $G_p$  is the specific gravity of the plasma 1.0065 is the gravity of protein free plasma ultrafiltrate. In this equation as formerly used in the literature and in previous publications of the copper sulfate method the constant 343 was used instead of the present 369. The value 343 was originally derived by Moore and Van Slyke from observations on nephritic patients for use with such patients

*Hematocrit values* We have used the formulation of Ashworth and Tigertt for calculating hematocrit values from gravities. Let  $G_b$ ,  $G_p$  and  $G_c$  represent the specific gravities of blood, plasma and centrifuged cells, respectively and H the hematocrit number in terms of ml of cells per 100 ml of blood. Then the weight of 100 ml of blood 100  $G_b$  can be equated to the sum of weight of its cells and plasma as follows

$$(2) 100 G_b = H G_c + (100 - H) G_p$$

Rearrangement of Equation 2 to calculate H gives

$$(3) H = 100 \times \frac{G_b - G_p}{G_c - G_p}$$

Ashworth and Tigertt pointed out that if cell gravity in different bloods is nearly constant one can substitute its mean value for  $G_c$  in Equation 3 and use the latter to estimate hematocrit values from observed  $G_b$  and  $G_p$  values

To estimate the value of  $G_c$  from observed  $G_b$ ,  $G_p$  and H Equation 3 is rearranged to give Equation 4

$$(4) G_c = \frac{100 (G_b - G_p) + H G_p}{H}$$

Ashworth and Adams measured in a series of normal and pathological bloods the values of  $G_b$ ,  $G_p$  and H, and calculated  $G_c$  by Equation 4. They thus found the mean value of  $G_c$  to be 1.0971. We have made similar observations in a series of 20 normal bloods and have obtained a mean  $G_c$  of 1.0970 identical with the value found by Ashworth and Adams. Substituting this value for  $G_c$  in Equation 3 gives

$$(5) H = 100 \times \frac{G_b - G_p}{1.0970 - G_p}$$

Equation 5 was used to calculate the hematocrit scale

In the experiments both of Ashworth and Adams and of the writers the hematocrit values used in Equation 4 to calculate  $G_c$  were obtained by the usual conventional technique for clinical hematocrits using Wintrobe tubes. Heparinized blood was centrifuged at 3000 rpm for 60 minutes the center of the centrifuge tube being 18 centimeters from the axis.

When a centrifuge was used in which the centers of the tubes were only 9 cm. from the axis 3000 rpm for an hour gave 51 instead of 4 as the mean



## 134 Biochemical Procedures

### Apparatus

- 1 Special graduated centrifuge tube (65 ml)
- 2 Electric centrifuge
- 3 Accurate stopwatch

### Reagents

- 1 Tsuchiya's solution Mix 15 gm of phosphotungstic acid 50 ml of concentrated hydrochloric acid (HCl) and 1000 ml of 95% ethyl alcohol

### Procedure

- 1 Pipette exactly 4.0 ml urine into a special centrifuge tube that is calibrated to 6.5 ml
- 2 Add 2.5 ml of phosphotungstic reagent
- 3 Mix several times by inversion and allow to stand *exactly* 10 minutes
- 4 Centrifuge for exactly 10 minutes at 1800 r p m and read amount of precipitate in terms of ml

### Calculation

gram protein/liter  $\approx$  ml of precipitate  $\times 7.2 \times$  dilution (if any)

### Example

Urine diluted 1 to 10

Reading 0.11 ml

$$\text{gram protein/liter} \approx (0.11) \times (7.2) \times \frac{(10)}{1} = 7.9$$

### Precautions

- 1 The timing is very important exactly 10 minutes being required
- 2 With urine from patients with nephritis it may be necessary to dilute the urine before analysis. A dilution of 1 to 10 is usually satisfactory

## 6 Fibrinogen in Plasma

### Reference

Butler, A. M. and Montgomery, H. The Solubility of the Plasma Proteins *J. Biol. Chem.* 99: 1, 3, 195 (Dec) 1932

### Principle

Fibrinogen is precipitated in a phosphate buffer and the supernatant fluid is analyzed for nitrogen. Fibrinogen is calculated by difference.

### Apparatus

- 1 A quantity of 25 ml flasks or test tubes
- 2 Rubber stoppers
- 3 Analytical grade filter paper
- 4 Funnels
- 5 Apparatus for micro Kjeldahl determinations (see total nitrogen)

### Reagents

- 1 Phosphate buffer 3.0 molar pH 6.5 Into a one liter flask put 408.5 gm of monobasic potassium phosphate ( $\text{KH}_2\text{PO}_4$ ) add 375 ml of accurate 4 N potassium hydroxide (KOH) and dilute to 1 liter with water
- 2 Other reagents are the same as those for the total nitrogen method.

### Procedure

- 1 To exactly 1.0 ml of phosphate buffer add 1 ml of plasma from citrated blood
- 2 Mix stopper. Let stand overnight at room temperature

- 3 Filter through good analytical filter paper
- 4 Analyze a 3 ml al quot for nitrogen by the micro Kjeldahl method
- 5 Determ ne total n trogen of the plasma by the m cro Kjeldahl method

#### Inculation

Fibrinogen = (total nitrogen - soluble nitrogen)  $\times$  6.25

#### Sample

Total nitrogen = 1.05 gm per 100 ml plasma  
 Soluble nitrogen = 0.98 gm per 100 ml plasma  
 grams fibrinogen per 100 ml plasma =  $(1.0 - 0.98) \times 6.25 = 0.44$

#### Precautions

- 1 The pH of the precipitating fluid is critical. It should be 6.5
- 2 Be sure that the blood is e treated (0.5 ml of saturated sodium citrate per 100 ml of blood)

### Nonprotein Nitrogen in Serum and Whole Blood

#### Reference

Daly C A. The Determination of Non Protein Nitrogen With Special Preference to the Koch McNeekin Method. J Lab Clin Med 18: 179-1285 (Sept) 1933

#### Principle

Combined nitrogen is converted to ammonia which forms a yellow color in the presence of Nessler's solution

#### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6  
 A variety of 15 ml round bottomed centrifuge tubes with rubber stoppers (#1 solid)
- 3 Electric centrifuge
- 4 Micro burners in the laboratory or gasoline stove in the field
- 5 Fohn digestion tubes graduated at 1.5 and 5 ml
- 6 Glass beads for digestion. These should be about 1 mm in diameter
- 7 A variety of pipettes including 0.5 ml, 1.0 ml, 5.0 ml graduated in 0.1 ml and 10 ml graduated in 0.1 ml
- 8 Cuvettes 19  $\times$  150 mm

#### Reagents

- 1 Nitrogen standard. Exactly 0.236 gm of dry ammonium sulfate is made up to exactly one liter with distilled water. One ml of this contains 0.05 mg N. It keeps for weeks in the cold.  
 Trichloroacetic acid approximately 8 gm in 100 ml water. It keeps for at least 2 weeks in the cold.
- 3 Digestion mixture. A stock concentrated solution is prepared by pouring slowly and carefully with cooling 2.5 ml of concentrated nitrogen free sulfuric acid into 45 ml of water. The working mixture is prepared by diluting exactly 31 ml of the stock solution to 100 ml with distilled water.
- 4 Hydrogen peroxide. Superoxol (30% H<sub>2</sub>O<sub>2</sub>) is diluted with distilled water in the proportion 1 to 10. This 3% H<sub>2</sub>O<sub>2</sub> keeps for a week in the cold.
- 5 Gum ghatti. Approximately 1 gm is dissolved in 100 ml of water and the solution is filtered through cotton.
- 6 Nessler's solution (after Koch). In a 2 liter mixing cylinder place 40 gm of mercuric iodide, 20 ml of water and 20 gm of potassium iodide. Mix until dissolved and dilute to 200 ml. Add with mixing 27 ml of 10% sodium citrate followed by



975 ml of exactly 10% sodium hydroxide. After the solution has cooled to room temperature dilute to exactly 1600 ml with water and store in a Pyrex vessel. It keeps indefinitely but does throw down a deposit.

### Procedure

- 1 Precipitate protein by adding exactly 0.5 ml of blood or serum to exactly 4.5 ml of 8% trichloroacetic acid in a 15 ml centrifuge tube. Stopper and mix by inversion. Do not shake or foaming will be annoying.
- 2 Centrifuge at high speed for 5 minutes. If necessary filter to take out particles.
- 3 To a digestion tube add exactly 1.0 ml filtrate, 1.0 ml of digestion mixture and 2 small glass beads.
- 4 Place the digestion tube in a slanting position directly over a free flame. When boiling commences add exactly 2 drops of 3%  $\text{H}_2\text{O}_2$  (Superoxol usually has some combined nitrogen in it and therefore the exact amount of  $\text{H}_2\text{O}_2$  is important). During digestion the fluid should turn brown and then become clear. At the end dense white fumes of  $\text{SO}_2$  will appear in the middle part of the tube.
- 5 Cool in the air. Inspect the digest to be sure there is no brown in it, and if there is, redigest.
- 6 At the same time as steps 3, 4 and 5 set up a blank tube containing 1 ml of water and a standard tube containing exactly 1 ml (= 0.050 mg N) of standard and treat them exactly like the unknowns.
- 7 Add approximately 10 ml of water and the 0.5 ml of gum ghatti solution.
- 8 Dilute to exactly 19.5 ml with water.
- 9 Add exactly 5 ml of Nessler's solution with continuous mixing.
- 10 Set spectrophotometer at 5.0 m $\mu$  with blank at 100% T and read unknowns.
- 11 Make up a standard curve using a range of 0.01 mg to 0.07 mg N at intervals of 0.01 mg preferably from standards run through the whole method.

### Calculations

- 1 mg nonprotein N/100 ml serum =

$$\frac{(\text{mg N in cuvette}) \times 5 \times 100}{(\text{ml filtrate}) \times (\text{ml serum})}$$

### Example

- 1 Unknown read 74.5% T equivalent to 0.03 mg N in cuvette. Serum was diluted 0.5 to 5.0 with trichloroacetic acid. Of the filtrate 1 ml was used for analysis.

$$\text{mg nonprotein N/100 ml serum} = \frac{(0.03) \times 5 \times 100}{(1.0) \times 0.5} = 30.0 \text{ mg N}$$

- 2 Since all samples of serum and blood are treated always alike, it is a saving of time to set up the standard curve in the manner shown in Fig. 20 and to read the final concentration directly.

### Precautions


- 1 Both Nessler's solution and gum ghatti deposit sediments on standing for several days. Hence before each day's run it is desirable to decant off enough of the clear supernatant fluid from each for the analyses in prospect.
- 2 The commonest source of error is splattering during digestion. This can be avoided by keeping the tubes in gentle motion until boiling starts.

## 8 Creatinine in Urine

### Reference

- Folin O. and Wu H. A System of Blood Analysis. *J. Biol. Chem.* 38: 93-100 (March) 1919.

# NON-PROTEIN NITROGEN IN SERUM AND BLOOD



## %T VS NPN — DIRECT CONVERSION

| %<br>T | 0                            | 1    | 2  | 3    | 4  | 5    | 6    | 7    | 8    | 9    |
|--------|------------------------------|------|----|------|----|------|------|------|------|------|
|        | Mg NPN/100 ml Blood or Serum |      |    |      |    |      |      |      |      |      |
| 50     | 70                           | 68   | 66 | 64   | 62 | 60   | 58.5 | 57   | 55   | 53   |
| 60     | 51.5                         | 50   | 48 | 46.5 | 45 | 43   | 42   | 40.5 | 39   | 37   |
| 70     | 38                           | 34.5 | 33 | 31.5 | 30 | 28.5 | 27.5 | 26   | 25   | 23.5 |
| 80     | 22                           | 21   | 20 | 18.5 | 17 | 16   | 15   | 13.5 | 12.5 | 12   |
| 90     | 10                           | 9    | 8  | 7    | 6  | 5    | 4    | 2.5  | 1.5  | 0.5  |
|        | NORMAL. 25 to 35 mg/100 ml   |      |    |      |    |      |      |      |      |      |

G.E. BORDON

Fig. 3—Use of calibration curve and conversion table for facilitating computation of results in estimating nonprotein nitrogen. The same format is useful for all kinds of colorimetric procedures.

**Principle**

Creatinine reacts with picric acid in alkaline solutions at room temperature and in a few minutes a stable intense orange color is produced which can be measured in a photocolormeter

**Apparatus**

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 x 150 mm
- 3 100 ml volumetric flasks
- 4 Calibrated syringe pipettes 0.1 ml and 10 ml
- 5 Volumetric pipettes of 20 ml capacity

**Reagents**

- 1 Sodium hydroxide, a 10% solution in water
- 2 Picric acid a 1% solution in water Gentle heating may be required to dissolve the acid Mix 100 ml NaOH and 100 ml picric acid and dilute to 1000 ml. This solution is called the alkaline picrate solution
- 3 Crystalline creatinine zinc chloride Dissolve and dilute 1.6106 gm to 1000 ml in 0.1 N hydrochloric acid This solution contains 1.0 mg of creatinine per ml

**Procedure**

- 1 Into a 100 ml volumetric flask measure accurately 0.1 ml of urine with a syringe pipette
- 2 Dilute to approximately 90 ml with distilled water and add exactly 20 ml of the mixed sodium hydroxide and picric acid mixture
- 3 Allow to stand for exactly 15 minutes after a gentle mixing and then dilute to 100 ml mark
- 4 A blank using exactly 90 ml of the alkaline picrate mixture is set up and diluted to 100 ml
- 5 An aliquot (10 ml) of the blank is pipetted into a cuvette and set at 100% T at a wave length of 5.0 mμ
- 6 Then read aliquots of the unknowns

**Calculation**

- 1 Prepare a standard curve from creatinine zinc chloride solution
  - a Dilute exactly 5 ml of standard to 100 ml This gives a solution of 0.05 mg creatinine per ml
  - b Into 100 ml volumetric flasks pipette 0.1 2 3 4 and 5 ml of the dilute standard
  - c Add 20 19 18 17 16 and 15 ml of water respectively
  - d Add exactly 20 ml of the same alkaline picrate used for analysis of unknowns
  - e Stand exactly 15 minutes and then dilute to 100 ml
  - f Set blank at 100% T and read all standards at 5.0 mμ
  - g Plot the %T curve on semi log paper

$$\text{gm creatinine/liter urine} = \frac{(\text{mg creatinine/ml of final solution}) \times 100 \times 1000}{(1000) \times (0.1)} = (\text{mg creatinine/ml final solution}) \times 1000$$

**Example**

A reading of 81.0% T was equivalent to 0.0009 mg/ml of final solution Therefore  
 gm creatinine/l urine = 0.0009 x 1000 = 0.9 gm/l

### Precautions

- 1 Many variables affect the intensity of the creatinine picrate. Important among these are time of standing, temperature, concentration of alkali, and presence of noncreatinine color producing substances. Therefore it is essential to establish standard conditions and to adhere strictly to these.

## 9 Creatine in Urine

### Reference

- Folin O and Wu H. A System of Blood Analysis. J Biol Chem 38 98 100 (March) 1919

### Principle

Creatine is converted to creatinine when boiled or autoclaved with acid. The creatinine is determined before and after the boiling treatment and creatine originally present in the urine is calculated by difference.

### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- Cuvettes 10 150 mm
- 3 A quantity of 100 ml volumetric flasks
- 4 Calibrated single pipettes 0.1 ml and 10 ml
- 5 Autoclave set at 15 pounds pressure
- 6 A quantity of volumetric pipettes of 10 and 50 ml capacity

### Reagents

- 1 Sodium hydroxide a 10% solution
- Picric acid a 1% solution in water. It may be necessary to warm this to obtain a clear solution.
- 3 Crytalline creatinine zinc chloride. The same concentration as for creatinine (1 mg/ml)
- 4 Hydrochloric acid 1 N. 10 ml of concentrated HCl plus 100 ml of distilled water

### Procedure

- 1 Into a 100 ml volumetric flask measure accurately 0.1 ml of urine with a syringe pipette. Add exactly 10 ml of picric acid solution, approximately 50 ml of water and 10 ml of HCl solution.
- 3 Autoclave at 15 pounds pressure at a temperature of 115 to 120 C for 30 minutes.
- 4 Allow to cool to room temperature.
- 5 Add exactly 10 ml of sodium hydroxide and dilute to 100 ml after exactly 15 minutes of standing.
- 6 Run a blank using 10 ml of picric acid solution and 50 ml of water. Autoclave as above and run exactly as unknowns.
- 7 Set blank at 100% T at 50 mμ.

### Calculation

- 1 (Total creatinine - preformed creatinine) = Creatine (as creatinine)  
To convert creatine expressed as creatinine into the amount of creatine itself multiply by 1.16
- 3 (Total creatinine - preformed creatinine)  $\times 1.16$  = Creatine

### Example

- 1 Total creatine = 11.0 gm per liter urine
- 2 Preformed creatine = 0.00 gm per liter
- 3 (11.0 - 0.00)  $\times 1.16$  = 0.3 gm of creatine per liter urine

**Precautions**

- 1 The precautions are the same as for creatinine with one important addition the hydrolysis of creatine may be affected markedly by temperature and acidity. Hence conditions during hydrolysis must be controlled rigidly.

**10 Urea Nitrogen****A. UREA NITROGEN IN BLOOD AND SERUM****Reference**

- Getzkyow C J and Masen J M An Accurate Method for the Determination of Blood Urea Nitrogen by Direct Nesslerization *J Biol Chem* 143 531-544 (Feb) 1942

**Principle**

The enzyme urease hydrolyzes urea into ammonium carbonate quantitatively at ordinary temperatures.

**Apparatus**

- 1 Coleman Jr Spectrophotometer Model #6
- 2 Cuvettes 19 x 150 mm
- 3 15 ml round bottom centrifuge tubes
- 4 Water bath set at 55 C
- 5 Electric centrifuge
- 6 Syringe pipettes exactly 0.2 ml and 0.3 ml

**Reagents**

- 1 Urease 10% solution Mix 4 gm of urease (Squibbs Powder) 27 ml of distilled water and 13 ml of 90% ethyl alcohol. This solution keeps only for one working day.
- 2 Sodium tungstate a 10% solution in water
- 3 Sulfuric acid 2% N Exactly 6.6 ml concentrated sulfuric acid ( $\text{H}_2\text{SO}_4$ ) is made up to 360 ml with distilled water.
- 4 Dilute Nessler's reagent Place 500 ml of 10% NaOH solution into a liter flask add 150 ml of Nessler's (for NPN) and dilute to 1000 ml with water. Mix and use on same day.
- 5 Standard solution of nitrogen Exactly 0.236 gm of dry ammonium sulfate is diluted to exactly one liter with distilled water. Store in refrigerator. 1.0 ml = 0.05 mg of nitrogen.
- 6 Permutit prepared according to Folin (Folin O and Bell R D Applications of a New Reagent for the Separation of Ammonia *J Biol Chem* 29 399 March 1917). Permutit is an insoluble sodium aluminum silicate (zeolite) that is used to absorb ammonia from solutions. It is best prepared (freeing the ammonia) by washing with 10% NaOH and then with 2% acetic acid. Decant wash with distilled water and dry in the air without heat. The recommended powder to use is 60 mesh.

**Procedure**

- 1 Add exactly 0.0 ml distilled water to a 15 ml centrifuge tube
- 2 Add exactly 0.1 ml serum with a syringe pipette
- 3 Add exactly 0.2 ml urease solution and incubate at 55 C for 15 minutes
- 4 Add exactly 0.3 ml of sodium tungstate solution
- 5 Add exactly 0.3 ml of 2/3 N sulfuric acid and allow to flocculate
- 6 Add exactly 1.0 ml of distilled water. Mix by inversion. Centrifuge until clear.

- 7 Into a cuvette add exactly 5 ml of supernatant fluid and 5 ml of water
- 8 Prepare a blank by adding 5 ml of dilute Nessler's solution to 10 ml of water
- 9 Set blank at 100% T at a wave length of 500 mμ
- 10 Add 2 ml of Nessler's to the unknowns and read immediately
- 11 Prepare a standard curve using a range of 0.1 mg to 0.07 mg

**Calculation**

$$\text{mg urea N/100 ml serum} = \frac{(\text{mg N in cuvette}) \times 10 \times 100}{(\text{ml aliquot}) \times (\text{ml serum})}$$

**Example**

A reading of 65% T was obtained equal to 0.0 mg N in the cuvette. Of serum 0.5 ml was used and a 5 ml aliquot of the supernatant fluid.

$$\text{mg urea N/100 ml serum} = \frac{(0.0) \times 10 \times 100}{(5) \times (0.5)} = 0.0 \text{ mg urea N}$$

**Precautions**

- 1 The samples must be read immediately as they cloud up in 1 to 2 minutes. Incubation time and temperature are important.
- 3 The urease solution is unstable and must be replaced daily.

**B UREA NITROGEN IN URINE****Procedure**

- 1 To a 50 ml flask add approximately 5 ml of urine
- 2 Add approximately 5 gm of permittit. Stopper and shake for 5 minutes. (The permittit absorbs the ammonia from the urine)
- 3 Allow to settle and decant off clear supernatant fluid
- 4 Dilute exactly 1 ml of ammonia free urine to exactly 50 ml with water
- 5 Pipette exactly 10 ml of diluted urine into a 15 ml centrifuge tube
- 6 Add exactly 10 ml distilled water
- 7 Add exactly 0.5 ml of urease solution and incubate for 15 minutes at 55°C
- 8 Add 0.4 ml of sodium tungstate solution
- 9 Add 0.4 ml of 1/3 N sulfuric acid and dilute to 10 ml (5 ml water)
- 10 Allow to flocculate and centrifuge until clear
- 11 Use 0.5 ml of supernatant fluid in a cuvette
- 12 Add 8.0 ml distilled water
- 13 Add 2.0 ml dilute Nessler's reagent and read immediately at 500 mμ
- 14 Set blank consisting of 10 ml of water and 0.5 ml of Nessler's solution at 100% T
- 15 Use same standard curve as for serum

**Calculation**

$$\text{gm urea N/liter urine} = \frac{(\text{mg N in cuvette}) \times (\text{ml dilution}) \times (1000)}{(1000) \times (\text{ml aliquot}) \times (\text{ml urine})}$$

**Example**

A reading of 63.0% T was equivalent to 0.01 mg N in the cuvette. Dilution was 10. One ml of urine was diluted to 50. Of this 1 ml was diluted to 10. An aliquot of 2.0 ml was taken for analysis.

$$\text{gm urea N/liter urine} = \frac{(0.01) \times (50) \times (10) \times (1000)}{(1000) \times (2) \times (1)} = 5.0 \text{ gm N}$$

**Precautions**

Same as for urea nitrogen in serum and blood

**Precautions**

- 1 The precautions are the same as for creatinine with one important addition the hydrolysis of creatine may be affected markedly by temperature and acidity. Hence conditions during hydrolysis must be controlled rigidly.

**10 Urea Nitrogen****A. UREA NITROGEN IN BLOOD AND SERUM****Reference**

- Gentzkow C J and Masen J M. An Accurate Method for the Determination of Blood Urea Nitrogen by Direct Nesslerization. *J Biol Chem* 143: 531-544 (Feb) 1942

**Principle**

The enzyme urease hydrolyzes urea into ammonium carbonate quantitatively at ordinary temperatures.

**Apparatus**

- 1 Coleman Jr Spectrophotometer Model #6
- 2 Cuvettes 10 x 150 mm
- 3 15 ml round bottom centrifuge tubes
- 4 Water bath set at 55 °C
- 5 Electric centrifuge
- 6 Syringe pipettes exactly 0.2 ml and 0.3 ml

**Reagents**

- 1 Urease 10% solution. Mix 4 gm of urease (Squibbs Powder) 27 ml of distilled water and 13 ml of 95% ethyl alcohol. This solution keeps only for one working day.
- 2 Sodium tungstate a 10% solution in water.
- 3 Sulfuric acid 2/3 N. Exactly 66 ml concentrated sulfuric acid ( $H_2SO_4$ ) is made up to 80 ml with distilled water.
- 4 Dilute Nessler's reagent. Place 500 ml of 10% NaOH solution into a liter flask add 150 ml of Nessler's (for NPN) and dilute to 1000 ml with water. Mix and use on same day.
- 5 Standard solution of nitrogen. Exactly 0.236 gm of dry ammonium sulfate is diluted to exactly one liter with distilled water. Store in refrigerator. 1.0 ml = 0.05 mg of nitrogen.
- 6 Permunit prepared according to Folin (Folin O and Bell R D. Applications of a New Reagent for the Separation of Ammonia. *J Biol Chem.* 29: 3-9 March 1917). Permunit is an insoluble sodium aluminum silicate (zeolite) that is used to absorb ammonia from solutions. It is best prepared (freeing the ammonia) by washing with 10% NaOH and then with 2% acetic acid. Decant wash with distilled water and dry in the air without heat. The recommended powder to use is 60 mesh.

**Procedure**

- 1 Add exactly 2.0 ml distilled water to a 15 ml centrifuge tube.
- 2 Add exactly 0.2 ml serum with a syringe pipette.
- 3 Add exactly 0.2 ml urease solution and incubate at 55 °C for 15 minutes.
- 4 Add exactly 0.3 ml of sodium tungstate solution.
- 5 Add exactly 0.3 ml of 2/3 N sulfuric acid and allow to flocculate.
- 6 Add exactly 7 ml of distilled water. Mix by inversion. Centrifuge until clear.

- Into a cuvette add exactly 5 ml of supernatant fluid and 5 ml of water
- 8 Prepare a blank by adding . ml of dilute Nessler's solution to 10 ml of water
  - 9 Set blank at 100% T at a wave length of 500 mμ
  - 10 Add . ml of Nessler's to the unknowns and read immediately
  - 11 Prepare a standard curve using a range of 0.1 mg to 0.07 mg

#### Calculation

$$\text{mg urea N/100 ml serum} = \frac{(\text{mg N in cuvette}) \times 10 \times 100}{(\text{ml aliquot}) \times (\text{ml serum})}$$

#### Example

A reading of 65% T was obtained equal to 0.0 mg N in the cuvette. Of serum 0.2 ml was used and a 5 ml aliquot of the supernatant fluid

$$\text{mg urea N/100 ml serum} = \frac{(0.0) \times 10 \times 100}{(5) \times (0.2)} = 0.0 \text{ mg urea N}$$

#### Precautions

- 1 The samples must be read immediately as they cloud up in 1 to 2 minutes. Incubation time and temperature are important
- 2 The urease solution is unstable and must be replaced daily

### B UREA NITROGEN IN URINE

#### Procedure

- 1 To a 50 ml flask add approximately 5 ml of urine
- 2 Add approximately . gm of permitt. Stopper and shake for 5 minutes (The permitt. absorbs the ammonia from the urine)
- 3 Allow to settle and decant off clear supernatant fluid
- 4 Dilute exactly 1 ml of ammonia free urine to exactly 50 ml with water
- 5 Pipette exactly 1.0 ml of diluted urine into a 15 ml centrifuge tube
- 6 Add exactly 1.0 ml of distilled water
- 7 Add exactly 0.2 ml of urease solution and incubate for 15 minutes at 55 C
- 8 Add 0.4 ml of sodium tungstate solution
- 9 Add 0.4 ml of 2% V sulfuric acid and dilute to 10 ml (7 ml water)
- 10 Allow to flocculate and centrifuge until clear
- 11 Use .0 ml of supernatant fluid in a cuvette
- 12 Add 8.0 ml of distilled water
- 13 Add .0 ml dilute Nessler's reagent and read immediately at 500 mμ
- 14 Set blank consisting of 10 ml of water and . ml of Nessler's solution at 100% T
- 15 Use same standard curve as for serum

#### Calculation

$$\text{gm urea N/liter urine} = \frac{(\text{mg N in cuvette}) \times (\text{ml dilution}) \times (1000)}{(1000) \times (\text{ml aliquot}) \times (\text{ml urine})}$$

#### Example

A reading of 63.0% T was equivalent to 0.01 mg N in the cuvette. Dilution was .0 ml of urine was diluted to 50 of this 1 ml was diluted to 10. An aliquot of .0 ml was taken for analysis

$$\text{gm urea N/liter urine} = \frac{(0.01) \times (50) \times (10) \times (1000)}{(1000) \times (.0) \times (1)} = 5.0 \text{ gm N}$$

#### Precautions

Same as for urea nitrogen in serum and blood



## 11 Amino Acid Nitrogen in Urine

## Reference

Folin O A Colorimetric Determination of the Amino Acid Nitrogen in Normal Urine *J Biol Chem.* 51 393 394 (April) 1922

## Principle

The 1 to 2 percent of the total nitrogen which is combined as amino acids in the urine is estimated colorimetrically

## Apparatus

- 1 Cuvettes 25 x 105 mm  
Coleman Jr Spectrophotometer Model #6
- 3 Tubes graduated to 25 ml
- 4 Dark room or equivalent
- 5 Accurately calibrated 5 ml average pipette

## Reagents

- 1 Permunit
- 2 1.25% sodium carbonate in water
- 3 0.1 N hydrochloric acid (HCl)
- 4 Sodium beta naphthoquinone-4 sulfonate 0.5% in water
- 5 Phenolphthalein indicator
- 6 Sodium acetate 2.5% in a 2.5% aqueous solution of acetic acid
- 7 Sodium thiosulfate 4% in water
- 8 Standard solution Amino acetic acid 0.536 gm plus 2.0 gm sodium benzoate up to 100 ml 0.1 N HCl (0.1 mg of amino acid nitrogen per ml)

## Procedure

- 1 Dilute exactly 5.0 ml of urine to 25 ml add 1.0 gm of permunit and shake (This step removes free ammonia)
- 2 Decant into another 3 gm of permunit Shake and decant again
- 3 Pipette 5.0 ml of ammonia free diluted urine into a 5 ml cuvette
- 4 Add 1.0 ml of 0.1 N HCl
- 5 Add 1.0 ml of 1.25% sodium carbonate and dilute to 10 ml
- 6 Then add 5 ml of sodium beta naphthoquinone-4 sulfonate solution and mix.
- 7 Stand overnight in dark and make acid to phenolphthalein
- 8 Add 5.0 ml of sodium thiosulfate and dilute to 25 ml
- 9 Set blank at 100% T and read at 480 mμ
- 10 Set up standard curve using 100 to 400 mg amino acid nitrogen

## Calculation

$$\text{mg amino N/liter urine} = \frac{(\text{mg amino N in cuvette}) \times (\text{dilution}) \times 1000}{(\text{ml aliquot in cuvette}) \times (\text{ml urine in original sample})}$$

## Example

A reading of 29.0% T was equivalent to 0.20 mg amino nitrogen in cuvette. Of urine 5 ml samples were diluted to 25 ml of which 5 ml was analyzed.

$$\text{mg amino N/liter urine} = \frac{(0.20) \times (25) \times (1000)}{(5) \times (5)} = .400 \text{ mg amino nitrogen}$$

## Precautions

- 1 Depending on the intensity of color developed at Step 8 more or less diluted urine may be required at Step 3

## 12 Uric Acid

## A. URIC ACID IN SERUM

## Reference

Benedict S R and Franke E A Method for the Direct Determination of Uric Acid in Urine J Biol Chem 52 38 391 (June) 19 2

## Principle

Uric acid reduces phosphotungstic acid and cyanide with the formation of a blue color which is measured color metrically

## Apparatus

- 1 Coleman Jr Spectrophotometer Model #6
- 2 Cuvettes 19 x 150 mm
- 3 A quantity of 15 ml round bottom centrifuge tubes
- 4 Electric centrifuge
- 5 Water bath set at 5 C
- 6 Calibrated syringe pipette 0.5 ml
- 7 Burette

## Reagents

- 1 Sodium tungstate ( $\text{Na}_2\text{WO}_4 \cdot \text{H}_2\text{O}$ ) 10% solution in water
- 2 Sulfuric acid ( $\text{H}_2\text{SO}_4$ ) /3 N  
6.6 ml concentrated sulfuric acid made up to 360 ml with distilled water
- 3 Lithium chloride ( $\text{LiCl}$ )  
3 gm  $\text{LiCl}$  and 1 ml hydrochloric acid (conc) made up to 1 liter with water
- 4 Silver nitrate ( $\text{AgNO}_3$ )  
11.6 gm silver nitrate per liter distilled water
- 5 Sodium cyanide ( $\text{NaCN}$ ) a 10% solution in water freshly made *POISON*
- 6 Color reagent Arsenophosphotungstic acid *POISON*  
Stock. Dissolve 100 gm sodium tungstate ( $\text{Na}_2\text{WO}_4 \cdot \text{H}_2\text{O}$ ) in about 600 ml of water in a 1 liter flask. Add 50 gm pure arsenic pentoxide followed by 5 ml of syrupy phosphoric acid (85%) and 0 ml concentrated hydrochloric acid. Bring to a boil and add 10 ml of bromine water (The stock solution is 0 ml water saturated with bromine). Boil until color disappears and continue to boil until all traces of bromine are driven off. Dilute to 1000 ml. This color reagent keeps indefinitely.  
For use dilute stock color reagent with equal portions of water before use
- 7 Benedict Hitchcock Stock Standard Solution  
Dissolve 9 gm of pure crystallized sodium phosphate ( $\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$ ) and 1 gm monobasic sodium phosphate ( $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ ) in about 50 ml hot water. Filter if not clear. Make up to 500 ml with hot water and pour this solution on exactly 0.0 mg of pure uric acid suspended in a few ml of water in a 1 liter volumetric flask. Shake gently until uric acid is in solution. Cool and add exactly 1.4 ml of glacial acetic acid dilute to 1 liter and mix well. Add 5 ml chloroform to prevent bacterial growth. Keep in refrigerator. Good for two months or more. 1 ml = 0.2 mg uric acid
- 8 Dilute uric acid standard made by diluting stock uric acid standard 1:10 so that 1 ml contains 0.01 mg uric acid

## Procedure (Serum)

- 1 Into a 15 ml round bottom centrifuge tube pipette exactly 0.5 ml serum from an accurately calibrated syringe pipette
- 2 Add exactly 8.5 ml water
- 3 Add exactly 0.5 ml sodium tungstate solution ( $\text{Na}_2\text{WO}_4 \cdot \text{H}_2\text{O}$ )

- 4 Then add exactly 0.5 ml sulfuric acid ( $\text{H}_2\text{SO}_4$ )
- 5 Stopper shake vigorously for approximately 15 seconds and centrifuge for 5 minutes
- 6 Transfer exactly 5 ml of the above filtrate into another 15 ml centrifuge tube
- 7 Add exactly 2.5 ml lithium chloride solution ( $\text{LiCl}$ ) and mix by inversion
- 8 Then add exactly 2.5 ml silver nitrate solution (Steps 7 and 8 are said to remove non uric acid interfering substances)
- 9 Stopper mix and centrifuge for 5-10 minutes. Disregard cloudiness of solution
- 10 Drain quantitatively into a cuvette
- 11 Add exactly 0.5 ml of color reagent *POISON*
- 12 Add 1 ml fresh  $\text{NaCN}$ . Use burette or syringe pipette for both steps 11 and 12
- 13 Let stand for 15 minutes at  $25^\circ\text{C}$  then read in a spectrophotometer with wave length at 540  $\text{m}\mu$ . The instrument should be set at 100% T against a blank of 10 ml water 0.5 ml color reagent and 1.0 ml cyanide

**Calculation**

- 1 Make up a standard curve using from 5 meg to 35 meg of uric acid per sample at intervals of 5 gamma
- 2  $\text{mg uric acid}/100 \text{ ml serum} =$

$$(\text{mg uric acid in cuvette}) \times \frac{(\text{dilution})}{(\text{ml aliquot})} \times \frac{100}{(\text{ml serum})}$$

**Example**

A reading of 78% T was equivalent to 0.008 mg uric acid in the cuvette. Of serum 0.5 ml was diluted to 10 ml of which 5 were taken for analysis

$$\text{mg uric acid}/100 \text{ ml serum} = (0.008) \times \frac{10}{5} \times \frac{100}{0.5} = 3.2 \text{ mg}$$

**Precautions**

- 1 Use a burette for both the cyanide and color reagent
- 2 The cyanide should be made up every two weeks and kept in the refrigerator
- 3 Time and temperature are critical

**B URIC ACID IN URINE****Principle Apparatus Reagents**

Same as for uric acid in serum

**Procedure**

- 1 Dilute the urine routinely 1 to 20 ml being satisfactory. Occasionally 1 to 50 may be necessary
- 2 Pipette exactly 1 ml of the dilute urine into a cuvette
- 3 Add exactly 9 ml of water
- 4 Proceed as for Serum steps 11, 12 and 13

**Calculation**

- 1 Use the same standard curve as for Serum
- 2  $\text{mg uric acid}/\text{liter urine} =$

$$(\text{mg uric acid in cuvette}) \times \frac{(\text{dilution})}{(\text{ml aliquot})} \times \frac{(1000)}{(\text{ml urine used})}$$

**Example**

A reading of 67% T was equivalent to 0.012 mg uric acid in the cuvette. One ml of urine was diluted to 10 of which one ml was used for analysis

$$\text{mg uric acid}/\text{l urine} = 0.012 \times \frac{20}{1} \times \frac{1000}{1} = 240 \text{ mg uric acid}/\text{l}$$

**Precautions**

Same as for Serum

## 13 Uric Acid (Uricase Method)

## A. URIC ACID IN URINE

## Reference

- Buchanan O H Block W D and Christman A A The Metabolism of the Methylated Purines I The Enzymatic Determination of Urinary Uric Acid *J Biol Chem* 157 181 '01 (Jan ) 1945

## Principle

The enzyme uricase when incubated with urine specifically destroys all the uric acid. Therefore by measuring the color intensity of a specimen before and after incubation one can calculate the true uric acid in the urine

## Apparatus

- 1 Coleman Jr Spectrophotometer Model No 6
- 2 Cuvettes 19 x 150 mm
- 3 A quantity of 15 ml round bottom centrifuge tubes
- 4 An electric centrifuge
- 5 A water bath set at 5 C
- 6 A water bath set at 45 C
- 7 Calibrated syringe pipettes 1 ml ( ) and 0.5 ml
- 8 Calibrated ml burettes ( )
- 9 A good food chopper or Waring Blender
- 10 A refrigerator set at 4 C
- 11 A quantity of 50 ml volumetric flasks
- 12 A quantity of filter paper
- 13 An interval timer

## Reagents

- 1 Color reagent arsenophosphotungstic acid (for preparation see uric acid in serum) *POISON*
- 2 Urea cyanide solution Dissolve 10 gm of sodium cyanide and 50 gm of anhydrous sodium carbonate in 400 ml of water Cool add 5 gm of urea and dilute to 500 ml with water This solution is stable for several months even though a slight precipitate will settle out on standing *POISON*
- 3 Uric acid standard solution 0.01 per cent This solution is made up in a 0.15 per cent solution of lithium carbonate (15 gm dissolved in 1000 ml of distilled water)
- 4 Sodium carbonate 0.1 M solution Dissolve 10 gm  $\text{Na}_2\text{CO}_3$  in 1000 ml of water
- 5 Toluene C.P. (reagent grade)
- 6 Sodium tungstate a 10% solution in water
- 7 Sulfuric acid 0.3 N
- 8 Acetone C.P. (reagent grade)
- 9 Benzene C.P. (reagent grade)
- 10 Sodium hydroxide an 0.1 N solution
- 11 Borate buffer (pH 9) Dissolve 1.4 gm of boric acid in 1000 ml of 0.1 N sodium hydroxide
- 12 Thymol blue indicator 0.04% solution
- 13 Uricase preparation of (M B Blanch and F C Koch A New Method for the Determination of Uric Acid in Blood With Uricase *J Biol Chem* 130 443 454 Oct 1939)
  - a One pound of fresh beef kidney is mechanically defatted and then ground finely in a food chopper or Waring Blender
  - b Dehydrate with 4 washings of acetone (500 ml each) over a period of 48 hours at 4 C 12 hours being allowed for each extraction

- c Dry the tissue quickly with the aid of a fan at room temperature
  - d Treat the dry material twice (at room temperature) with 100 ml of benzene the first treatment for 6 hours followed by filtering and the second for 19 hours
  - e The material is dried at room temperature and then ground to 60 mesh
  - f Store at 4 C
- 14 *Measurement of uricolytic activity*
- a To 10 ml of a standard uric acid solution (0.01% in 0.15% lithium carbonate) add 1 ml of 0.125% sodium carbonate
  - b Add 1 drop of toluene and increasing amounts of the kidney pulp
  - c Incubate for two hours at 45 C with occasional stirring
  - d Add 30 ml of distilled water
  - e Add 1 ml of a 10% sodium tungstate solution
  - f Add 2 ml of 2/3 N sulfuric acid
  - g Shake and filter
  - h Then determine the uric acid content in the clear filtrate
  - i The amount of uricase to be used per analysis should be determined as the point where all the uric acid is destroyed 200 mg of the powder should destroy 3.4 mg of uric acid after incubation

**Procedure Total Color**

- 1 Pipette exactly 1 ml of urine (using a syringe pipette) into a 50 ml flask, dilute to the mark with distilled water and mix
- 2 To another 50 ml flask transfer exactly 10 ml of the diluted urine and add 15 ml of water
- 3 To a second flask add only 20 ml of water This sample will be the blank
- 4 To each flask add 2.5 ml of urea cyanide solution using a syringe pipette
- 5 Then add with a syringe pipette 1 ml of arsenophosphotungstic acid
- 6 Immediately after the addition of the color reagent dilute to the 50 ml mark and mix thoroughly
- 7 Allow to stand at room temperature for exactly 30 minutes and read the samples at 690 mμ, setting the blank at 100% T

**Procedure Residual Color**

- 1 Pipette exactly 5 ml of undiluted urine into a small flask add 10 ml of water and a few drops of 0.04 per cent thymol blue solution as indicator
- 2 Titrate with 0.1 N sodium hydroxide to a definite blue tint and record the amount of alkali required Discard the mixture
- 3 Into a 50 ml flask pipette exactly 5 ml of undiluted urine
- 4 Add the predetermined amount of alkali (without indicator)
- 5 Add 250 mg of uricase powder
- 6 Add 5 ml of borate buffer and wash down the sides of the flask with 10 ml of water
- 7 Place in a water bath at 45 C for 2 hours
- 8 Then add 1 ml of sodium tungstate solution and 1.5 ml of 2/3 N sulfuric acid.
- 9 Dilute to 50 ml mark with water mix well and pour into a dry filter
- 10 Pipette exactly 10 ml of the filtrate into a 50 ml volumetric flask and add 15 ml of water
- 11 Continue with the addition of the color reagents and measure the color in the same way as described above for the total color
- 12 The same blank may be used or a 5 ml portion may be treated with uricase and all the other reagents

**Calculation**

- 1 (Total color) minus (Residual color) = true uric acid content
- 2 mg true uric acid per 100 ml urine =

$$\left[ \text{mg (total color) uric acid in cuvette} \times \frac{\text{dilution}}{\text{ml al quot}} \times \frac{100}{\text{ml urine used}} \right] \\ \text{minus} \\ \left[ \text{mg (re dual color) uric acid in cuvette} \times \frac{\text{dilution}}{\text{ml al quot}} \times \frac{100}{\text{ml urine used}} \right]$$

**Example**

- 1 The reading for the total true color was equivalent to 0.006 mg uric acid in the cuvette. One ml of urine was diluted to 50 ml of which 10 ml was used and diluted to 50 ml. 10 ml was read in the cuvette.
- 2 The reading for the residual color was equivalent to 0.004 mg uric acid. 5 ml of urine was diluted to 50 ml of which 10 ml was used in the final sample and diluted to 50 ml. 10 ml was read in the cuvette.

$$\text{Total color mg/100 ml} = 0.006 \times \frac{50}{10} \times 10 \times \frac{100}{1} = 30 \text{ mg/100 ml}$$

$$\text{Residual color mg/100 ml} = 0.004 \times \frac{50}{10} \times 10 \times \frac{100}{5} = 4 \text{ mg}$$

$$\text{Therefore True uric acid mg/100 ml} = 30 - 4 = 26 \text{ mg}$$

**Precautions**

- 1 Use a bulbette or a syringe pipette for the sodium cyanide and urea cyanide. They are strong poisons.
- 2 Be sure to measure the uricolytic activity of each new batch of uricase.
- 3 Uricase is available commercially. However for most laboratories the price is prohibitive.

**14 Paper Chromatography**

This technique has afforded a new tool whereby many biological factors may be studied. One may investigate the amino acids in plant and animal tissue with a minimum amount of preparation and destruction. Only small quantities are required to demonstrate the presence of amino acids. Nor is its use confined solely to the estimation of amino acids. Its simplicity has permitted its ready adaptation to the analysis of such substances as sugars and organic acids. With proper development it holds more promise to the chemistry of proteins, carbohydrates, organic acids, and to the isotopic derivatives thereof than any other technique in the last decade.

During the past few years filter paper chromatography has been applied to many different categories of compounds. The bibliography below is by no means complete but will serve to introduce the prospective user to some of the various phases of this versatile method.

**Methodology**

- Brimley R C. Quantitative Paper Chromatography. *Nature* 163: 215-16 (Feb.) 1949.  
 Bull H B. *Advances in Protein Chemistry* by M L Anson and J T Edsall. New York City 1947. Academic Press vol 3 p 97.

- Bull H B Hahn J W and Baptist, V H Filter Paper Chromatography *J Am. Chem. Soc* 71 550 553 (Feb) 1949
- Chargaff E Levine C and Green C Techniques for the Demonstration by Chromatography of Nitrogenous Lipide Constituents Sulfur Containing Amino Acids and Reducing Sugars *J Biol Chem* 175 67 71 (Aug) 1948
- Consden R and Gordon A H Effects of Salt on Partition Chromatograms, *Nature* 167 180 181 (July) 1948
- Dent C E Detection of Amino acids in Urine and Other Fluids *Lancet* 2 13 637 (Nov) 1948
- Dunabin J E Mason H Seyfang A P and Woodman F J Spectrographic Examination of Chromatographic Columns *Nature* 164 916 (Nov) 1949
- Fink R M Dent C E and Fink K Application of Filter Paper Partition Chromatography to Radioactive Tracer Studies *Nature* 160 801 803 (Dec) 1947
- Fink R M and Fink K Radiocarbon and Filter Paper Partition Chromatography *Science* 107 253 254 (March) 1949
- Fisher, R B Parsons D S and Morrison G A Quantitative Paper Chromatography *Nature* 161 764 765 (May) 1948
- Fisher, R B Parsons D S and Holmes R Quantitative Paper Chromatography *Nature* 164 183 (July) 1949
- Keston A S Udenfriend S and Levy M Paper Chromatography Applied to the Isotopic Derivative Method of Analysis *J Am Chem. Soc* 69 3151 3156 (Dec) 1947
- Ja Cour L F and Drew R Partition Chromatography and Living Cells *Nature* 159 4035 307 (March) 1947
- Longnecker W H Glass Trough for Filter Paper Partition Chromatography *Science* 107 93 24 (Jan) 1948
- Lugg J W H and Overell B T Partition Chromatography of Organic Acids on a Paper Sheet Support *Nature* 160 4055 87 88 (July) 1947
- Martin A J P Partition Chromatography *Ann N Y Acad Sci* 49 249 264 (Feb) 1948
- Meinhard J E and Hall N F Surface Chromatography *Anal Chem* 21 185 188 (Jan) 1949
- Meinhard J E and Hall N F Surface Chromatography *Anal Chem* 22 (2) 344 351 (Feb) 1950
- Miettinen J K and Virtanen A I A New Technique in Paper Chromatography *Acta Chem Scandinav* 3 459 464 No 5 1949
- Nicholson D E Modified Technique for the Development of Paper Chromatograms *Nature* 163 954 (June) 1949
- Phillips D M P Use of Ultra Violet Fluorescence in Paper Chromatography *Nature* 161 53 (Jan) 1948
- Pollard F H McOmie I F W and Elbeih I I M Inorganic Paper Chromatography and Detection of Cations by Fluorescence *Nature* 163 216 219 (Feb) 1949
- Rutter L A Modified Technique in Filter Paper Chromatography, *Nature* 161 435 436 (March) 1948
- Schroeder W A Some Experiments in Systematic Quantitative Chromatography *Ann N Y Acad Sci* 49 24 217 (Feb) 1948
- Williams R J and Kirby H Paper Chromatography Using Capillary Ascent, *Science* 107 481 483 (May) 1948
- Winsten W A A Simplified Apparatus for One Dimensional Paper Partition Chromatography, *Science* 107 600 (June) 1948
- Wynn V Peptide like Contaminant of Filter Paper *Nature* 164 445 1949
- Zechmeister L History Scope and Methods of Chromatography *Ann. N Y Acad. Sci.* 49 145 160 (Feb) 1948
- Amino Acids**
- Agren G A Note on the Amino Acid Content of Bence Jones Protein *Acta Chem Scand* 3 (3) 301 302 1949
- Agren G and Nilsson, T Paper Chromatographic Analysis of Amino Acids and Other Ninhydrin Reacting Substances in Deproteinized Human Plasma *Acta Chem Scand* 3 525 538 No 5 1949
- Ames B R and Risley H A Aminoaciduria in Progressive Muscular Dystrophy *Fed. Proc* 7 1 142 (March) 1948
- Awapara J Application of Paper Chromatography to the Estimation of Some Free Amino Acids in Tissues of the Rat *J Biol Chem* 178 113 116 (March) 1949
- Berry H K and Cain L Biochemical Individuality IV Paper Chromatographic Technique for Determining Excretion of Amino Acids in the Presence of Interfering Substances *Arch Biochem* 24 (1) 1 9 189 (Nov) 1949
- Block R J Quantitative Estimation of Amino Acids on Paper Chromatograms *Science* 103 608 609 (Nov) 1948

- Blois R J Quantitative Paper Chromatography A Simplified Procedure *Proc Soc Exper Biol & Med* 72 (6) 337-341 (Nov) 1949
- Consden P Gordon A H and Martin A J Qualitative Analysis of Proteins A Partition Chromatographic Method Using Paper *Biochem J* 39 224-23 1944
- Consden R Gordon A H and Martin A J The Identification of Amino Acids Derived From Cystine in Chemically Modified Wool, *Biochem J* 40 390-597 1946
- Consden P Gordon A H and Martin A J Gramicidin S The Sequence of Amino Acid Residues *Biochem J* 41 596-607 1947
- Dent C E Partition Chromatography on Paper as Applied to the Investigation of Amino Acids and Peptides in Normal and Pathological Urine *Biochem J (Proc)* 40 (4) xlv-xlv 1946
- Dent C E Partition Chromatography on Paper A Qualitative Method for Identifying Amino Acids in Urine and Its Application to the Fanconi Syndrome *Trans of Fourteenth Meeting Conference on Metabolic Aspects of Convalescence New York (Nov)* 1946
- Dent C E The Amino Aciduria in Fanconi Syndrome A Study Making Extensive Use of Techniques Based on Paper Partition Chromatography *Biochem J* 41 (2) 240-253 1947
- Dent C E Methionine Metabolism and Alpha Amino Acid Science 105 37 336 (March) 1947
- Dent C E Stepka W and Steward F C Detection of Free Amino Acids of Plant Cells by Partition Chromatography *Nature* 160 40 2 687 (Nov) 1947
- Goldberg H J V Gilda J E and Tshkoff G H Paper Partition Chromatography Free Amino Acids in Saliva *J Dent Research* 27 493-496 (Aug) 1948
- Gordon A H Martin A J and Syge R L Partition Chromatography of Free Amino Acids and Peptides *Biochem J* 37 Proc xlv 1943
- Heston A S Ulfstrand G and Lott M Quantitative Analysis of Protein Hydrolysates on Paper Chromatograms by Means of the Isotopic Derivative Method *Fed Proc* 107 1 164 (March) 1948
- Landau A J and Awapara J Use of Modified Ninhydrin Reagent in Quantitative Determination of Amino Acids by Paper Chromatography *Science* 109 384 (April) 1949
- Martin A J P and Mittelmann R Quantitative Microanalysis of Amino Acid Mixtures on Paper Partition Chromatograms *Biochem J* 43 53 58 1949
- Moore S and Stein W H Partition Chromatography of Amino Acids on Starch *Ann NY Acad Sci* 49 8 (Feb) 1948
- Moore S and Stein W H Photometric Ninhydrin Method for Use in the Chromatography of Amino Acids *J Biol Chem* 176 367-394 (Oct) 1948
- Moore S and Stein W H Chromatography of Amino Acids on Starch Columns Solvent Mixtures for the Fractionation of Protein Hydrolysates *J Biol Chem* 178 53-78 (March) 1949
- Nastahn L Quantitative Chromatographic Estimation of  $\alpha$ -Amino Acids *Nature* 161 763-764 (May) 1948
- Poison A Moyle V J and Wyckoff R G W The Quantitative Chromatography of Silk Hydrolyzates *Science* 105 603-604 (June) 1948
- Pope C G and Steen M F The Determination of Amino Nitrogen Using a Copper Method *Biochem J* 43 10 0 10 9
- Rockland L B and Dunn M S Quantitative Determination of Amino Acids on Filter Paper Chromatograms by Direct Photometric *J Am Chem Soc* 71 (1) 41-44 (Dec) 1949
- Sanger F The Free Amino Groups of Insulin *Biochem J* 39 56 315 1944
- Stein W H and Moore S Chromatography of Amino Acids on Starch Columns Separation of Phenylalanine, Leucine, Isoleucine, Methionine, Tyrosine and Valine *J Biol Chem* 178 337-365 (Oct) 1948
- de Verdier C H and Agren G Paper Chromatographic Analysis of Amino Acids and Peptides in Tissue Extracts and Enzyme Hydrolyzed Proteins *Acta Chem Scandinav* 2 783-790 1948
- Winegard H M and Toennies G Detection of Sulfur Containing Amino Acids on Paper Chromatograms *Science* 108 506-507 (Nov) 1949
- Winwood A J Micro Estimation of Amino Nitrogen and Its Application to Paper Partition Chromatography *Nature* 161 169 (Jan) 1948
- Winwood A J Method for the Estimation of Micro Amounts of Amino Nitrogen and Its Application to Paper Partition Chromatography *Biochem J* 45 417-417 1949
- Young N F and Homburger F The Application of Paper Chromatography to the Study of Amino aciduria in Patients With Liver Disease *Led. Pro* 1 31 (March) 1948



**Carbohydrates**

- Brown F Hirst E L Hough L Jones J K N and Wadman H Separation and Identification of Methylated Sugars on the Paper Chromatogram, *Nature* 181 770 (May) 1948
- Flood A E Hirst E L and Jones J K N Quantitative Analyses of Mixtures of Sugars by the Method of Partition Chromatography Part I. Standardization of Procedure *J Chem Soc* 128 1679 1683 1948
- Flood A E Hirst E L and Jones J K Quantitative Estimation of Mixtures of Sugars by the Paper Chromatogram Method *Nature* 160 405 86 (July) 1947
- Horrocks R H and Manning G B Partition Chromatography on Paper Identification of Reducing Substances in Urine *Lancet* 256 1042 1045 (June) 1949
- Hough L Jones J K N and Wadman W H Application of Paper Partition Chromatography to the Separation of the Sugars and Their Methylated Derivatives on a Column of Powdered Cellulose *Nature* 162 448 (Sept) 1948
- Partridge S M with a note by Westall F G Filter Paper Partition Chromatography of Sugars 1 General Description and Application to the Qualitative Analysis of Sugars in Apple Juice Egg White and Foetal Blood of Sheep *Biochem J* 4 38 750 No 2 1948
- Partridge S M Application of Paper Partition Chromatogram to the Qualitative Analysis of Reducing Sugars *Nature* 158 4008 270 (Aug) 1946

**Miscellaneous**

- Allsopp A Chromatographical Study of Meristematic Plant Tissue *Nature* 181 4100 833 835 (May) 1948
- Ames S and Risley H A Determination of Creatine Creatinine and Related Compounds in Urine by Means of Paper Chromatography *Proc Soc Exper Biol & Med* 69 267 269 (Oct Dec) 1948
- Bate-Smith E C Paper Chromatography of Anthocyanins and Related Substances in Petal Extracts *Nature* 161 4100 83 833 (May) 1948
- Bentley H R and Whitehead J K Use of Furan Derivatives in Paper Chromatography *Nature* 164 182 183 (July) 1949
- Cavallini D Frontali N and Toschi G Determination of Keto Acids by Partition Chromatography on Filter Paper *Nature* 163 568 569 (April) 1949
- Christ C L Burton C J and Botty M C Use of X Ray and Electron Diffraction as Methods of Analysis in Biochemical Chromatography, *Science* 108 91 92 (July) 1943
- Crammer J L Paper Chromatography of Flavin Nucleotides *Nature* 161 4088 350 (March) 1949
- Dingemans E Huizendt Velt L G and de Laat B M Clinical Method for the Chromatographic Colorimetric Determination of Urinary 17 Keto Steroids *J Clin Endocrinol* 6 (8) 535 548 (Aug) 1946
- Edman P Hammarsten E Low B and Reichard P Partition Chromatographic Separation of Adenine and Guanine *J Biol Chem* 178 395 398 (March) 1949
- Ekman B Paper Chromatography of Primary Aromatic Amines *Acta Chem. Scand* 2 383 384 No 4 1948
- Franklin A F and Quastel J H Paper Chromatography of Proteins and Enzymes *Science* 110 447 451 (Oct) 1949
- Hais I M and Pecáková L Paper Partition Chromatography of Riboflavin Decomposition Products *Nature* 163 768 (May) 1949
- Hess W G Chromatographic Separation of Cholesterol and Cholesterol Esters in Blood 1 *Lab & Clin Med* 32 1163 1168 (Aug) 1947
- Holiday E R and Johnson E A Location of Paper Chromatogram Spots of Purine and Pyrimidine Derivatives in Ultra violet Light *Nature* 163 216 17 (Feb) 1949
- Klaeener R G A Paper Chromatographic Method for the Quantitative Estimation of Penicillin Entities *J Bact* 57 101 109 (Jan) 1949
- Lacourt A Sommeveves G Degeyndt E Baruh J and Gillard J Quantitative Inorganic Paper Chromatography Sub Micro Separation and Determination of Aluminum Iron and Titanium *Nature* 163 999 1000 (June) 1949
- Lederer M Chromatographic Separation of Antimony *Anal Chim Acta* 2 61 6 1948
- Lederer M Paper Chromatography of the Noble Metals *Nature* 162 776 777 (Nov) 1948
- Markham R and Smith J D Chromatography of Nucleic Acid Derivatives *Nature* 163 250 251 (Feb) 1949
- Maw G A The Detection of Creatine and Creatinine by Partition Chromatography *Biochem J* 43 139 142 1948
- Mitchell H A Gordon M and Haskins F A Separation of Enzymes on Filter Paper Chromatopile *J Biol Chem* 180 1071 1076 (Oct) 1949

- Polson A New Amino Acids in Bacterial Hydrolysates *Nature* 161 4088 151 15 (March) 1948
- Ramsey L L Separation of *n* Butyric and Isobutyric Acids by Partition Chromatography *J Assoc Offic Agr Chemists* 31 164 1948
- Ramsey L L and Patterson, W I Separation and Determination of the Straight Chain Saturated Fatty Acids C to C by Partition Chromatography *J Assoc Offic Agr Chemists* 31 139 150 1948
- Tabone J Robert D and Troestler J Application of Paper Strip Chromatography to the Analysis of Biological Arylamides *Bull Soc Chim Biol* 30 547 557 1948
- Taugro A Long W and Chaikoff I L Fractionation of Thyroid Iodine by Means of Filter Paper Partition Chromatography *Nature* 164 181 18 (July) 1949
- Tishkoff G H, Zaffaroni A and Tesluk H Purified Liver Extract Chemical Nature as Determined by Paper Partition Chromatography *J Biol Chem* 175 867 86 (Sept) 1949
- Tomarelli R M and Florey K Use of Papergrams in the Study of the Urinary Excretion of Radioactive Sulfur Compounds *Science* 107 630 631 (June) 1948
- Urbach K F and Giscafre L Identification of Histamine in Blood by Paper Chromatography *Proc Soc Exper Biol & Med* 68 430 431 (May) 1948
- Vischer E and Chargaff E The Separation and Characterization of Purines in Minute Amounts of Nucleic Acid Hydrolysates *J Biol Chem* 168 781 8 (May) 1947
- Vischer E and Chargaff E The Separation and Quantitative Estimation of Purines and Pyrimidines in Minute Amounts *J Biol Chem* 176 703 714 (Nov) 1948
- White J W Jr Chromatographic Separation of Aliphatic 4 Dinitrophenylhydrazones *Anal Chem* 20 796 7 8 (Aug) 1948
- White J W Jr and Dryden E C Separation of Aliphatic Alcohols by Chromatographic Adsorption of Their 35 Dinitrobenzoates *Anal Chem* 20 853 855 (Sept) 1948
- Winsten W A and Egen E Paper Chromatography of Vitamin B and Related Bacteriological Growth Factors *J Biol Chem* 181 (1) 109 1 0 (Nov) 1949
- Woodruff H B and Foster J C Analysis of Vitamin B and Vitamin B by Paper Strip Chromatography *J Biol Chem* 183 569 576 (Apr) 1950
- Zaffaroni A Burton E B and Keutmann E H The Application of Paper Partition Chromatography to Steroid Analysis *J Biol Chem* 177 109 116 (Jan) 1949
- Zaffaroni A, Burton E B and Keutmann E H Adrenal Cortical Hormones Analysis by Paper Partition Chromatography and Occurrence in the Urine of Normal Persons *Science* 111 68 (Jan) 1950

## 15 Proteins and Nitrogenous Compounds Miscellaneous References

### Alanine

- Alexander B, and Seligman A M A Colorimetric Method for the Microdetermination of  $\alpha$  Alanine in Blood *J Biol Chem* 159 9 19 (June) 1949
- Wiss O Microdetermination of Alanine in Blood *Helv Chim Acta* 31 6 1948

### Albumin and Globulin

- Chow B F The Correlation Between the Albumin and Alpha Globulin Contents of Plasma *J Clin Investigation* 26 (5) 883 886 (Sept) 1947
- Cohn C and Wolfson W Q Studies in Serum Proteins I The Chemical Estimation of Albumin and of the Globulin Fractions in Serum *J Lab & Clin Med* 32 (10) 1 03 1 07 (Oct) 1947
- Cohn C and Wolfson W Q Studies in Serum Proteins II A Rapid Clinical Method for the Accurate Determination of Albumin and Globulin in Serum or Plasma *J Lab & Clin Med* 33 (3) 367 370 (March) 1948
- Ishister J The Estimation of Serum Albumin by the Precipitation of Globulin with NaHSO *Australian J Exper Biol & Med* 27 61 64 (Jan) 1949
- Kingsley G P A Rapid Method for the Separation of Serum Albumin and Globulin *J Biol Chem* 133 731 735 (May) 1940
- Petermann M L, Young N F and Hogness K R A Comparison of the Howe and the Electrophoretic Methods for the Determination of Plasma Albumin *J Biol Chem* 169 (2) 379 387 (July) 1947

### Allantoin

- Iarson H W A Colorimetric Method for the Determination of Allantoin *J Biol Chem* 94 (3) 727 7 8 (Jan) 1933

### Amino Acids

- Moubasher R Estimation of  $\alpha$  Amino Acids in Pure Solutions in Blood and in Urine with Perinaphthendian 34 Trione Hydrazide *J Biol Chem* 175 187 193 (Aug) 1948
- Sahyun M Relationship of Amino Acids to the Nutritive Value of Proteins *Am J Digestive Disease* 16 (8) 301 06 (Aug) 1949

**Amino Nitrogen**

Danielson I S Amino Acid Nitrogen in Blood and Its Determination *J Biol Chem* 101 505 502 (July) 1933

**Creatine**

Fisher R B and Wilhelm A E The Metabolism of Creatine I A Micro Method for the Determination of Creatine and Creatinine *Biochem J* 31 ( ) 1131 1135 1937

Rose W C Helmar O M and Chanutin A A Modified Method for the Estimation of Total Creatinine in Small Amounts of Tissues *J Biol Chem* 75 (4) 543 548 (Nov) 1907

Addis T Barrett E and Menzies J T A Clinical Method for the Approximate Determination of Serum Creatinine Concentration *J Clin Investigation* 26 (5) 879 88 (Sept) 1947

Barclay J A and Kenney R A A Method for the Estimation of Creatinine *Biochem J* 41 (4) 586 589 1947

**Cysteine**

Nakamura K and Binkley F Colorimetric Estimation of Cysteine *J Biol Chem* 173 407 410 (March) 1948

**Histidine**

Langley W D Urinary Histidine Determination of Histidine in Urine Normal and Pregnancy Urines *J Biol Chem* 137 955 966 (Jan) 1941

**Non protein Nitrogen**

Fee D A Cruger D and Collier, H B A Photometric Modification of the Hypobromite Method for Non protein Nitrogen *J Lab & Clin Med* 34 (6) 873 876 (June) 1949

Rappaport F and Lichhorn F Rapid Titrimetric Micromethod for the Determination of Non protein Nitrogen *J Lab & Clin Med* 32 (8) 1034 1036 1947

Wong S Y The Use of Persulphate in the Estimation of Nitrogen by Folin's Direct Nesslerization Method *J Biol Chem* 85 431 436 (March) 19 3

**Nucleic Acid**

Steele R Sfortunato T and Ottolenghi L A Micromethod for the Determination of the Nucleic Acid *J Biol Chem* 177 321 325 (Jan) 1949

**Plasma Protein**

Albanese A A Saur B and Irby V The Microcolorimetric Estimation of Plasma Protein *J Lab & Clin Med* 52 (2) 296 299 1947

Barbour H G and Hamilton W F The Falling Drop Method for Determining Specific Gravity *J Biol Chem* 69 (2) 6 5 640 (Aug) 1906

Cohn E J Strong L E Hughes W L Jr Mulford D J Ashworth J N Mehn M and Taylor H L Preparation and Properties of Serum and Plasma Proteins IV A System for the Separation into Fractions of the Protein and Lipoprotein Components of Biological Tissues and Fluids *J Am Chem Soc* 68 459 475 (March) 1946

Hunter D Nesslerization Technique and Its Effect on Serum Protein Values *Am J Clin Path* 17 (8) 650 653 (Aug) 1947

Kagan B M A Simple Method for the Estimation of Total Protein Content of Plasma and Serum I A Falling Drop Method for the Determination of Specific Gravity II The Estimation of Total Protein Content of Human Plasma and Serum by the Use of the Falling Drop Method *J Clin Investigation* 17 369 376 (July) 1938

Lowry O H and Hunter T H The Determination of Serum Protein Concentration With a Gradient Tube *J Biol Chem* 159 465 464 (July) 1945

Mulford D J Derivatives of Blood Plasma *Ann Rev Physiol* 9 3 1 356 1947

Roberts S and White A Studies on the Origin of Serum Proteins *J Biol Chem* 180 (4) 505 516 (Sept) 1949

**Protein**

Hiller A Greif R L and Beckman W W Determination of Protein in Urine by the Biuret Method *J Biol Chem* 176 1491 1499 (Dec) 1948

**Purines**

Hitchings G H and Fiske C H The Determination of the Purines *J Biol Chem* 140 (9) 491 499 (Aug) 1941

Hunter D and Givens H S The Metabolism of Endogenous and Exogenous Purines in the Monkey *J Biol Chem* 17 37 53 (Feb) 1914

Soodak M Piroo A and Cerecedo L R A Colorimetric Method for the Estimation of Uracil and Cytosine *J Biol Chem* 181 713 718 (Dec) 1949

# Total Nitrogen

- Carpenter A M Dilution of Capillary Blood in the Microanalysis of Protein Am J Clin Path 17 (7) 565-68 (July) 1947  
Kirk P L Kjeldahl Method for Total Nitrogen Anal Chem 22 ( ) 354-358 (Feb) 1950

# Tryptophan

- Albanese A A Davis V I and Levin M T Utilization of  $\alpha$ -Amino Acids by Man VIII Tryptophan and Acetyltryptophan J Biol Chem. 172 (1) 29-44 (Jan) 1948

# Tryptophane

- Albanese A. A. and Frankston J E Estimation of Tryptophane in Human Urine J Biol Chem. 157 39-63 (Jan) 1945  
Berg C P and Fohle W G The Tryptophane Content of Normal Human Urine J Biol Chem 170 : , 51-59 (Oct) 1944  
Carpenter D C Spectrophotometric Estimation of Tryptophan Anal Chem 20 33C-33 (June) 1948  
Gordon M and Mitchell H A Fluorometric Method for the Estimation of Tryptophan J Biol Chem. 180 (2) 1065-1070 (Oct) 1944  
Spies J R and Chambers D C Chemical Determination of Tryptophane Anal. Chem. 20 30-3 (Jan) 1948  
Steer E and Savage M. G. Micromethod for Determination of Tryptophane in Bacteria and Protein Anal Chem 21 641-64 (May) 1944

# Tyrosine

- Goodwin T V and Morton R A The Spectrophotometric Determination of Tyrosine and Tryptophane in Proteins Biochim J 40 648-657 1946

# Urea

- Fohn O and Svedberg A An Improved Distillation Method for the Determination of Urea in Blood J Biol Chem 88 7-83 (Aug) 1930  
Obermer L and Milton R Further Notes on a New Colorimetric Technique for the Estimation of Urea in Urine Analyst 63 434-44 (May) 1934  
Van Slyke D D and Cullen G E A Permanent Preparation of Urease and Its Use in the Determination of Urea J Biol Chem 19 211-8 (Sept) 1914  
Van Slyke D D and Cullen G E The Determination of Urea by the Urease Method J Biol Chem 24 117-127 (Jan) 1916

# Uric Acid

- Block W D and Geis N C An Enzymatic Method for the Determination of Uric Acid in Whole Blood J Biol Chem. 165 (2) 74-76 (May) 1947  
Dulger H A and Johns H F The Determination of Plasma Uric Acid J Biol Chem 140 407-440 (Aug) 1941  
Fohn O Standardized Methods for the Determination of Uric Acid in Unclotted Blood and in Urine J Biol Chem 101 111-125 (June) 1930  
Fohn O The Preparation of Sodium Tungstate Free from Molybdate Together With a Simplified Process for the Preparation of a Correct Uric Acid Reagent J Biol Chem 108 311-314 (Aug) 1934  
Schaffer N E. The Determination of Uric Acid in Urine With Crude Uricase J Biol. Chem 153 163-166 (April) 1944  
Wolfson W Q Huddleston B and Levine P The Transport and Excretion of Uric Acid in Man II The Endogenous Uric Acid in Chromatography of Biological Fluids: J Clin. Investigation 26 (5) 995-1001 (Sept) 1947

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### C CARBOHYDRATES AND DERIVATIVES

##### 1 Blood Glucose (Colorimetric)

###### Reference

Folin O and Malmros H. An Improved Form of Folin's Micro Method for Blood Glucose Determination J Biol Chem 83 115 1<sup>st</sup> (July) 19 9

###### Principle

The glucose is oxidized with alkaline ferricyanide and the ferrocyanide produced is measured colorimetrically as Prussian blue. Gum ghatti is added to keep the Prussian blue in colloidal suspension. So intense is the resulting color that two separate determinations can be made on the extract from 0.1 ml of normal blood. The final color is a mixture of the yellows of the excess ferricyanide and ferric sulfate with the blue of the Prussian blue.



###### Prussian blue

###### Apparatus

- 1 Coleman Junior Spectrophotometer Model 6
- 2 Cuvettes 19 x 150 mm
- 3 Tubes blood sugar type graduated at 1<sup>st</sup> 5 and 25 ml
- 4 Centrifuged tubes 15 ml round bottom
- 5 Micro pipettes, 0.1 ml. These pipettes are calibrated to be washed out
- 6 A water bath with boiling water
- 7 An electric centrifuge
- 8 A wire rack to hold the tubes while boiling

###### Reagents

- 1 Sulfuric acid. A 2/3 N solution made by diluting 67 ml of concentrated H<sub>2</sub>SO<sub>4</sub> to 360 ml with water
- 2 Sodium tungstate. A 10% solution in water
- 3 Dilute tungstate solution. Transfer 20 ml of 10% sodium tungstate solution into a liter volumetric flask. Dilute to about 800 ml and while shaking add 10 ml 2/3 N H<sub>2</sub>SO<sub>4</sub>. Dilute to 1000 ml with water
- 4 Potassium ferricyanide solution K<sub>3</sub>Fe(CN)<sub>6</sub>. Dissolve 2 gm of pure potassium ferricyanide (recrystallized) in 500 ml flask and dilute to mark. Keep in dark bottle away from light
- 5 Sodium carbonate. Dissolve 8 gm of sodium carbonate in 500 ml of water
- 6 Ferric iron solution. Fill a liter cylinder to the mark with cold tap water. Push into this cylinder a circular piece of copper screening large enough to form a bowl well below the surface of the water. Place in it 20 gm gum ghatti # 1

Cover cylinder and set as do for 18 to 24 hours. At the end of this time remove screen and undissolved residue. Strain through double thickness of towelling to remove dirt. By aid of heat dissolve 5 gm of anhydrous ferric sulfate in 75 ml of 8.7% phosphoric acid and 100 ml distilled water. Cool and mix with gum ghatti. To the final mixture add 1% potassium permanganate solution at first 5 ml at a time and later 3 ml at a time until the pink color remains perceptible for at least three minutes. This step is essential for oxidation of certain materials in the gum ghatti which reduces ferricyanide. Precaution: Avoid excess of  $\text{KMnO}_4$  which leads to high and erratic results.

- 7 Standard stock glucose solution. Dissolve gm of benzoic acid in about 500 ml of hot water in a 1 liter flask. Weigh out exactly 20 gm of anhydrous glucose (Bureau of Standards) and place in the flask. Dilute to about 900 ml and cool. Then dilute to volume. The concentration of this solution is 2 gm/liter or 2 mg/ml.
- 8 Working glucose solution standard. Dissolve 2 gm of benzoic acid in about 500 ml of hot water. Add 10 ml of stock standard to a liter flask. Then add the hot benzoic acid solution and dilute to about 1500 ml. When cool dilute to mark. This solution contains 0.01 mg of glucose per ml. Store in icebox.

#### Procedure

- 1 Measure accurately 0.1 ml blood in a centrifuge tube containing 10 ml of dilute tungstic acid. Rin out pipette 8 times. Stopper and shake tube vigorously.
- 2 Centrifuge 35 minutes at 1000 r.p.m.
- 3 Pipette out exactly 4 ml of filtrate into a tube graduated at 1.5 and 5 ml.
- 4 Add ml of 0.4% potassium ferricyanide.
- 5 Add 1 ml of sodium carbonate solution just before boiling.
- 6 At the same time set up tubes containing 4 ml of water for a blank and 4 ml of working glucose solution as a check on the standard curve. Treat as in Steps 4 and 5.
- 7 Heat in boiling water for exactly 8 minutes and cool as soon as possible in running cold water 1 or 2 minutes.
- 8 Add 5 ml of ferric iron alongside tube to prevent foaming.
- 9 Mix and after 5 minutes dilute to 1 ml mark.
- 10 Set blank at 100% T at 5.0  $\mu$  wave length and read standard and unknowns.

#### Calculation

- 1 Set up a calibration curve with glucose at intervals in the range 0.01 to 10 mg glucose in cuvette setting blank at 100% T wave length 5.0  $\mu$ .
- mg glucose/100 ml blood =

$$\frac{(\text{mg glucose in cuvette}) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml blood sample})}$$

#### Example

A reading of 50.5% T was equivalent to 0.04 mg glucose in the cuvette. Samples of blood were 0.1 ml diluted to 10.0 ml of which 4 ml were used for analysis. (The precipitate is about 0.1 ml containing very little glucose.)

$$\text{mg glucose/100 ml blood} = \frac{(0.04) \times (10) \times (100)}{(4) \times (0.1)} = 100 \text{ mg/100 ml}$$

#### Precautions



- 1 Be sure the water is boiling before adding the sodium carbonate.
- 2 The boiling time is exactly 8 minutes.
- 3 The micro pipettes must be washed out at least 8 times to get out all the blood.

**2 Glucose (Volumetric)****References**

- Shaffer P A and Hartman A F The Iodometric Determination of Copper and Its Use in Sugar Analysis *J Biol Chem.* 45 349-390 (Jan) 1921  
 Somogyi M Notes on Sugar Determinations *J Biol Chem* 70 599-619 (Nov) 1928  
 Somogyi M A New Copper Reagent *J Biol Chem* 160 61-63 (Sept) 1945

**Principle**

The chief advantage of this method is its usefulness when samples contain widely varying amounts of glucose. Being titrimetric it covers a very wide range without necessitating extra dilutions. Cupric ions are reduced to cuprous in a hot medium. In acid solution iodine is released stoichiometrically from an iodide iodate mixture and is titrated with thiosulfate.

**Apparatus**

- 1 Filter paper quantitative
- 2 A number of special sugar tubes (see references)
- 3 Boiling water bath
- 4 A calibrated burette 20 ml in capacity
- 5 A quantity of 1.5 ml centrifuge tubes (round bottom) with #1 solid rubber stoppers

**Reagents**

- 1 Phosphate tartrate copper iodide solution. Dissolve 28 gm disodium phosphate in 600 ml water. Add 40 gm Rochelle salts dissolve and add 100 ml of 1 N sodium hydroxide. Add 80 ml of a 10% solution of cupric sulfate and with constant stirring 180 gm anhydrous sodium sulfate. Dilute to one liter with distilled water and add 8 gm potassium iodide and 25 ml of 3.5% potassium iodate.
- 2 Sodium thiosulfate. Prepare and standardize 0.1000 N thiosulfate as described in the section on volumetric analysis. As needed prepare 0.003 N thiosulfate by diluting 25 ml of 0.1 N to 500 ml.
- 3 Starch indicator. To 50 ml of boiling water add 1 gm of soluble starch and stir until dissolved. Then dilute to 100 ml with cold water. To prevent molds dissolve 1 gm of salicylic acid in the water before using it in the starch solution.
- 4 Sulfuric acid. Dilute 20 ml of concentrated sulfuric acid to 360 ml with distilled water. This is approximately 2 N.
- 5 Zinc sulfate. 1.5% solution in water.
- 6 Barium hydroxide. An 0.3 N solution is prepared by dissolving 5.0 gm  $\text{Ba}(\text{OH})_2$  in 1000 ml warm water. Titrate against zinc sulfate. These two solutions must be equal in strength.

**Procedure**

- 1 To a 15 ml centrifuge tube add 4.5 ml of barium hydroxide solution.
  - 2 Add exactly 1 ml of blood.
  - 3 Then add 4.5 ml of zinc hydroxide stopper and shake vigorously.
  - 4 Filter.
  - 5 Pipette exactly 2 ml of the filtrate into a special sugar tube.
  - 6 Add 2 ml of copper reagent.
- The blank consists of 2 ml of distilled water plus 2 ml of copper reagent.

- 8 Cover both tubes and heat in a boiling water bath for exactly 10 minutes. Cool under running water.
- 9 Add 1 ml of sulfuric acid rapidly to each tube while mixing gently.
- 10 Add 1 to 2 drops of starch indicator and titrate with 0.005 N sodium thiosulfate. The end point is the disappearance of the blue color.

#### Calculation

Addition of iodide to the copper reagent increases the speed of this method but cuts down its accuracy in the titration range 0 to 1 ml. Above a titration figure of 1 ml 1 ml of 0.005 N thiosulfate is equivalent to 0.135 mg glucose. However, to save time in computation it is best to use a table derived empirically for the first ml (Table 9).

#### Example

The blank titration was 0.12 ml the unknown titration was 1.6 ml. One ml of blood was diluted to 10 ml of which 2 ml aliquots were analyzed. Reference to Table 9 shows that (unknown - blank) i.e. (1.6 - 0.12) is equivalent to 110 mg glucose per 100 ml blood.

#### Precautions

- 1 Strict attention to every detail of this method is required. Especially important are:
  - a Timing
  - b Avoidance of reoxidation of cuprous compounds once they are formed
  - c Temperature

TABLE 9  
BLOOD GLUCOSE SHAFER HASTMAN METHOD

(Direct conversion of titration ml 0.005 N thiosulfate into blood glucose when 2 ml aliquots of a 1:10 dilution are analyzed)

| ml CVK<br>MINUS<br>ml BLK | 0                             | 0.1 | 0.2 | 0.3 | 0.4 | 0.5 | 0.6 | 0.7 | 0.8 | 0.9 |
|---------------------------|-------------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
|                           | mg GLUCOSE/100 ml WHOLE BLOOD |     |     |     |     |     |     |     |     |     |
| 0                         | 0                             | 11  | 16  | 20  | 3   | 40  | 47  | 54  | 61  | 64  |
| 1                         | 5                             | 8   | 89  | 96  | 103 | 110 | 117 | 124 | 131 | 137 |
| 2                         | 143                           | 149 | 155 | 161 | 169 | 175 | 181 | 189 | 195 | 20  |
| 3                         | 09                            | 15  | 21  | 8   | 34  | 41  | 48  | 54  | 61  | 69  |
| 4                         | 75                            | 81  | 288 | 0   | 36  | 309 | 315 | 32  | 39  | 336 |
| 5                         | 34                            | 34  | 356 | 363 | 370 | 376 | 383 | 389 | 396 | 403 |
| 6                         | 400                           | 416 | 4   | 49  | 46  | 443 | 449 | 456 | 46  | 463 |
| 7                         | 476                           | 48  | 489 | 496 | 503 | 510 | 516 | 523 | 530 | 536 |
| 8                         | 543                           | 549 | 556 | 563 | 570 | 576 | 583 | 590 | 596 | 603 |
| 9                         | 610                           | 616 | 623 | 630 | 637 | 644 | 650 | 657 | 663 | 670 |

### 3 Lactic Acid


#### A. IN BLOOD AND SERUM

##### Reference

Edwards H. T. A Simplified Estimation of Lactate in Normal Human Blood  
J Biol Chem 125 571 583 (Oct) 1938

##### Principle

Lactic acid is oxidized to an aldehyde which is volatile. The acetaldehyde is absorbed by bisulfite and the bound bisulfite is estimated iodometrically. The steps in the procedure are:







**2 Glucose (Volumetric)****References**

- Shaffer P. A. and Hartman A. F. The Iodometric Determination of Copper and Its Use in Sugar Analysis *J Biol Chem* 45 349 390 (Jan) 1921  
 Somogyi M. Notes on Sugar Determinations *J Biol Chem* 70 599 612 (Nov) 1926  
 Somogyi M. A New Copper Reagent *J Biol Chem* 100 61 68 (Sept) 1935

**Principle**

The chief advantage of this method is its usefulness when samples contain widely varying amounts of glucose. Being titrimetric it covers a very wide range without necessitating extra dilutions. Cupric ions are reduced to cuprous in a hot medium. In acid solution iodine is released stoichiometrically from an iodide iodate mixture and is titrated with thiosulfate.

**Apparatus**

- 1 Filter paper quantitative
- 2 A number of special sugar tubes (see references)
- 3 Boiling water bath
- 4 A calibrated burette 10 ml in capacity
- 5 A quantity of 15 ml centrifuge tubes (round bottom) with #1 solid rubber stoppers

**Reagents**

- 1 Phosphate tartrate copper iodide solution. Dissolve 28 gm disodium phosphate in 600 ml water. Add 40 gm Rochelle salts dissolve and add 100 ml of 1 N sodium hydroxide. Add 80 ml of a 10% solution of cupric sulfate and with constant stirring 180 gm anhydrous sodium sulfate. Dilute to one liter with distilled water and add 8 gm potassium iodide and 25 ml of 3.5% potassium iodate.
- 2 Sodium thiosulfate. Prepare and standardize 0.005 N thiosulfate as described in the section on volumetric analysis. As needed prepare 0.005 N thiosulfate by diluting 50 ml of 0.1 N to 500 ml.
- 3 Starch indicator. To 50 ml of boiling water add 1 gm of soluble starch and stir until dissolved. Then dilute to 100 ml with cold water. To prevent molds dissolve 1 gm of salicylic acid in the water before using it in the starch solution.
- 4 Sulfuric acid. Dilute 10 ml of concentrated sulfuric acid to 260 ml with distilled water. This is approximately 2 N.
- 5 Zinc sulfate. A 5% solution in water.
- 6 Barium hydroxide. An 0.3 N solution is prepared by dissolving 36.0 gm  $\text{Ba}(\text{OH})_2$  in 1000 ml warm water. Titrate against zinc sulfate. These two solutions must be equal in strength.

**Procedure**

- 1 To a 15 ml centrifuge tube add 4.5 ml of barium hydroxide solution.
- 2 Add exactly 1 ml of blood.
- 3 Then add 4.5 ml of zinc hydroxide stopper and shake vigorously.
- 4 Filter.
- 5 Pipette exactly 2 ml of the filtrate into a special sugar tube.
- 6 Add 2 ml of copper reagent.
- 7 The blank consists of 2 ml of distilled water plus 2 ml of copper reagent.

- 8 Cover both tubes and heat in a boiling water bath for exactly 10 minutes. Cool under running water.
- 9 Add 1 ml of sulfuric acid rapidly to each tube while mixing gently.
- 10 Add 1 to 2 drops of starch indicator and titrate with 0.005 N sodium metabisulfite. The endpoint is the disappearance of the blue color.

#### Calculation

Addition of iodide to the copper reagent increases the speed of this method but cuts down its accuracy in the titration range 0 to 1 ml. Above a titration figure of 1 ml 1 ml of 0.005 N metabisulfite is equivalent to 0.135 mg glucose. However, to save time in computation it is best to use a table derived empirically for the first ml (Table 9).

#### Example

The blank titration was 0.1 ml the unknown titration was 1.6 ml. One ml of blood was diluted to 10 ml of which 1 ml aliquots were analyzed. Reference to Table 9 shows that (unknown - blank) = (1.6 - 0.1) is equivalent to 110 mg glucose per 100 ml blood.

#### Precautions

- 1 Strict attention to every detail of this method is required. Especially important are:
  - a Timing
  - b Avoidance of oxidation of cuprous compounds once they are formed
  - c Temperature

TABLE 9  
BLOOD GLUCOSE SHAWVER-HARTMAN METHOD

(Direct conversion of titration ml 0.005 N metabisulfite into blood glucose when 1 ml aliquots of a 1:10 dilution are analyzed.)

| ml LNK<br>ml N S<br>ml BLK | 0                 | 0.1 | 0.2 | 0.3 | 0.4 | 0.5         | 0.6 | 0.7 | 0.8 | 0.9 |
|----------------------------|-------------------|-----|-----|-----|-----|-------------|-----|-----|-----|-----|
|                            | mg GLUCOSE/100 ml |     |     |     |     | WHOLE BLOOD |     |     |     |     |
| 0                          | 0                 | 11  | 18  | 25  | 3   | 40          | 47  | 54  | 61  | 68  |
| 1                          | 75                | 8   | 89  | 96  | 103 | 110         | 117 | 124 | 131 | 137 |
| 2                          | 143               | 149 | 155 | 16  | 169 | 175         | 18  | 189 | 195 | 0   |
| 3                          | 09                | 215 | 1   | 28  | 234 | 41          | 48  | 54  | 61  | 68  |
| 4                          | 25                | 81  | 88  | 95  | 30  | 309         | 315 | 3   | 39  | 336 |
| 5                          | 342               | 349 | 356 | 363 | 370 | 376         | 383 | 389 | 396 | 403 |
| 6                          | 403               | 418 | 4   | 49  | 436 | 443         | 449 | 456 | 46  | 469 |
| 7                          | 476               | 482 | 489 | 496 | 503 | 510         | 516 | 523 | 529 | 536 |
| 8                          | 543               | 549 | 556 | 563 | 570 | 577         | 583 | 590 | 596 | 603 |
| 9                          | 610               | 616 | 623 | 630 | 637 | 644         | 650 | 657 | 663 | 670 |

### 3 Lactic Acid

#### A. IN BLOOD AND SERUM

##### Reference


Edwards, H. T. A Simplified Estimation of Lactate in Normal Human Blood  
J Biol Chem 125: 571-583 (Oct.) 1938

##### Principle

Lactic acid is oxidized to acetaldehyde which is volatile. The acetaldehyde is absorbed by bisulfite and the bound bisulfite is estimated iodometrically. The steps in the procedure are:

- 1  $\text{CH}_3\text{CH}(\text{OH})(\text{COOH}) + \text{O} \rightleftharpoons \text{CH}_3\text{CHO} + \text{CO}_2 + \text{H}_2\text{O}$  (at 100°C in  $\text{KMnO}_4$ )
- 2  $\text{CH}_3\text{CHO} + \text{NaHSO}_3 \rightleftharpoons \text{CH}_3\text{CH}(\text{OH})\text{NaHSO}_3$  (at pH 3.5)

# LACTIC ACID UNIT



gram of the apparatus to scale. The reaction flask is 15 cm its capacity 5 ml.

- 3  $\text{HSO}_3^- + 2\text{I}^- + \text{H}^+ = \text{HSO}_4^- + 2\text{HI}$
- 4  $\text{CH}_3\text{CHO} + \text{NaHSO}_3 = \text{CH}_3\text{CH}(\text{OH})\text{SO}_3\text{Na}$  (at pH 8.9)
- 5  $\text{HSO}_3^- + \text{I}^- + \text{H}^+ = \text{HSO}_4^- + 2\text{HI}$

#### Apparatus

- 1 Special lactic acid apparatus (See Fig. 91)
- 2 50 ml centrifuge tubes
- 3 Electric centrifuge
- 4 0.5 ml pipettes or accurately calibrated syringe pipette
- 5 Microburette 3 ml accurately calibrated

#### Reagents

- 1 Sodium tungstate ( $\text{Na}_2\text{WO}_4 \cdot \text{H}_2\text{O}$ ) approximately 10% Adjust to pH 8.9 using thymol blue as an indicator with dilute NaOH or  $\text{H}_2\text{SO}_4$ . (High blanks can sometimes be avoided by boiling 1 liter of this with 0.5 ml of 3%  $\text{H}_2\text{O}_2$ )
- 2 Sulfuric acid ( $\text{H}_2\text{SO}_4$ ) 2/3 N made from special concentrated Pregl acid 66 ml  $\text{H}_2\text{O}$  made up to 360 ml of water
- 3 Manganous sulfate ( $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ ) approximately 10%
- 4 Sulfuric acid approximately 10 N made from special Pregl sulfuric acid (to get a low blank) in distilled water 100 ml is made up to 360 ml with water (Concentrated acid is 36 N)
- 5 Sodium metabisulfite ( $\text{Na}_2\text{S}_2\text{O}_5$ ) approximately 5% Keep in glass stoppered bottle Do not use after two weeks
- 6 Potassium permanganate ( $\text{KMnO}_4$ ) approximately N/1000
- 7 Starch solution approximately 1%
- 8 Sodium metabisulfite ( $\text{Na}_2\text{S}_2\text{O}_5$ ) 0.05% Make fresh daily
- 9 Iodine ( $\text{I}_2$ ) approximately N/3 3 gm  $\text{I}_2$  and 48 gm of KI made up to 1 liter Use crude iodine and USP KI Add more KI if necessary to get the iodine in solution
- 10 Iodine approximately N/10 Dilute from above
- 11 Iodine standard N/100 made by diluting standard N/10 iodine N/200 is stable for 23 days. N/10 is stable for several months Make up and standardize N/10 iodine
- 12 Disodium phosphate ( $\text{H}_2\text{N}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$ ) ground to a powder Use Elmer and Amend sodium phosphate secondary USP granular
- 13 Standard lithium lactate  $\text{CH}_3\text{CHOHCOOH}$  mol wt 95.940 Prepare the salt according to Friedemann and Graessner (J. Biol. Chem. 100: 291-308 March 1933) Prepare an accurate 0.1 M solution in 0.2 N  $\text{H}_2\text{SO}_4$ . Dilute this stock standard 1 to 500 with 0.1 N  $\text{H}_2\text{SO}_4$ . The diluted standard keeps well in the refrigerator. Of this 15 ml is equivalent to 0.270 mg lactic acid or 1.0 ml N/100 iodine

#### Procedure

This method is also valid for trichloroacetic acid centrifugates from blood and muscle. Once the blood is hemolyzed it can be kept in the icebox at any stage of precipitation the value being constant for a week or more. It is best to remove the filtrate from precipitate if it will not be analyzed within a few days.

- 1 For whole blood pipette exactly 1 ml of whole blood into 45 ml of distilled water. For plasma pipette 1 ml plasma into 17 ml of distilled water. Mix. For cells pipette 1 ml of cells into 53 ml of distilled water. Mix. For cells it is best to use a pipette calibrated to contain and wash out. Add mixing between additions.  
For whole blood pipette 0.5 ml of 2/3 N sulfuric acid and 2 ml of 10% tungstate.  
For plasma 1 ml of 2/3 N sulfuric acid and 1 ml of 10% tungstate.  
For cells 3 ml of 2/3 N sulfuric acid and 3 ml of 10% tungstate.
- 2 Centrifuge 5-10 minutes at 3000 rpm. Note:  
a. If alcohol is present boil the filtrate for 5 minutes with new reagents in Steps 4 & 8 before adding permanganate.

- b If ether is present bubble cold air through the filtrate 15 minutes before using
- c In mild ketosis use the filtrate just as with normal blood
- 4 Clean absorption towers and test tubes with distilled water
- 5 Pipette 10 or 15 ml aliquots of the centrifugate into the reaction flasks (Lift up the  $KMnO$  reservoir 6 inches, insert pipette lower reservoir to hold it, and allow pipette to drain while you do the other pipettings. This saves time.) For standard lactate use about 0.2 mg of lactate
- 6 Add approximately 5 ml 10 N sulfuric acid approximately 10 ml of 10%  $MnSO$  and approximately 50 ml distilled water to each reaction flask. Mix these three in graduate and add all together. When first run is finished add fresh unknown solutions directly to the used solutions in the reaction flasks even when the old  $MnSO$   $H_2SO$  solutions are deep brown or hot, or both. As many as 9 subsequent runs can be made without addition of fresh reagents or cleaning out the flasks (Add 0.5 ml 10 N  $H_2SO$  every 4 runs)
- 7 Add approximately 15 ml 5%  $NaHSO$  to each absorption tube but not through tower adjust absorption tower fill air intake reservoir with  $KMnO$ , (N/200) (approximately 9 ml) with stopcocks shut and connect all the rubber stoppers
- 8 Make sure water in condenser is flowing. Adjust the air to flow briskly (The exact rate does not matter provided it is brisk.) Light the burners bring to boiling point and add  $KMnO$  by adjusting stopcocks to drip so that after 10 minutes a faint pink shows in each flask. About 1 drop every 5 seconds is needed. If  $KMnO$  collects at dripping point a vigorous shake after the run will dispel the collection (In subsequent runs a brown precipitate obscures the pink color but drop in the  $KMnO$  at the usual rate.) Be sure to add all of the  $KMnO$  in dropping reservoir in order to oxidize all unknown
- 9 Add  $KMnO$  for a total of exactly 15 minutes. Continue boiling for 5 additional minutes
- 10 At the end of the 20 minutes turn off the burners and turn off the air flow. Loosen one tube at a time and wash the tower while holding the tube with the tip of the tower near the top of the tube
- 11 Wash the towers with four 2 ml portions of distilled water from a fast running 2 ml pipette. Do this by washing with the first portion around the top of the absorption tower and with the next three so as to cover  $\frac{1}{4}$  of the tower each time. Wash the tip of the tower with about 1 ml of water

### **Titration**

Note If room temperature is above  $74^\circ$ , the tubes must be kept cold in a pan of ice

- 1 Add 4 drops of starch to the tube
- 2 Add strong iodine first then 0.01 N iodine until a blue color persists (Avoid adding a large excess of I because this raises the blank. There is no need to hurry with the titration at this point)
- 3 Add 0.05%  $NaHSO$  until the blue color is discharged (Clean this burette with distilled water and add fresh 0.05%  $NaHSO$  each day)
- 4 Add 0.005 N iodine until a faint blue color persists 15 seconds. When you think you have reached the end point tip the tube so that the contents rinse down the sides. This prevents drops of iodine or  $NaHSO$  on the sides of the tube giving bad results. Be sure the end point is still present after tipping
- 5 Add 2 spatulas full of  $Na_2HPO \cdot 12H_2O$  to discharge the blue color and mix thoroughly by swirling the tube
- 6 Read the burette containing the standard N/200 iodine
- 7 Add N/200 iodine until the first faint blue color that persists for 15 seconds agitate the tube constantly. It does not matter if this end point fades after 15 seconds
- 8 Read the burette again.
- 9 Subtract the readings

- 10 Subtract the blank from this result (The blank should never be over 0.04 ml of 0.005 N iodine) It is well if possible to run a distilled water blank and a standard lactate solution with each run. The reagent blank should be the same as a water blank. If it is not change the reagents. The commonest sources of high blanks are the  $\text{H}_2\text{SO}_4$ , the  $\text{Na}_2\text{WO}_4$  and the excess iodine during titration. The blank appears to be due to the distilled water used.
- 11 Duplicate titrations must agree to 0.0 ml of 0.005 N  $\text{I}_2$ .

#### Calculation

mg lactic acid/100 ml blood =

$$\frac{\left( \begin{array}{c} \text{ml titration} \\ \sim \text{ml blank} \end{array} \right) \times \left( \begin{array}{c} \text{normality of} \\ \text{iodine} \end{array} \right) \times (45) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml blood sample})}$$

#### Example

The standard iodine was 0.0049 N. One ml of whole blood was diluted to 50 ml of which 15 ml aliquots were used. The average titration was 1.09 ml and the blank 0.04 ml of iodine.

mg lactic acid/100 ml blood =

$$\frac{(1.09 - 0.04) \times 0.0049 \times 45 \times 50 \times 100}{15 \times 1} = 77.2 \text{ mg/100 ml}$$

#### Precautions

- 1 There are many steps in this method three of which (oxidation, binding of aldehyde, titration) are subject to many errors. Most technicians need a great deal of practice before the method runs smoothly.
- 2 Cleaning the apparatus is best done by standing in a mixture of 4 vols concentrated  $\text{HCl}$  + 1 vol of 5%  $\text{NaHSO}_3$  until the precipitated  $\text{MnO}_2$  is dissolved. Then wash thoroughly with cold tap water and rinse off with distilled water. Precipitates of tungstic oxide can be removed with dilute  $\text{NaOH}$ . The burette containing 0.005 N  $\text{NaHSO}_3$  and containing 0.005 N  $\text{I}_2$  must be rinsed out with distilled water each day. The burette containing more concentrated iodine may be left until its contents have been completely used.

### B. LACTIC ACID IN MUSCLE

#### Preparation of Specimen

- 1 Set the tube of a homogenizer tube in a pan of solid  $\text{CO}_2$ . Pipette in 2 ml of 10% trichloroacetic acid and allow to freeze.
- 3 The muscle or tissue (approximately 250 mg) is weighed on a Roller-Smith type balance and frozen in powdered  $\text{CO}_2$ .
- 4 It is then dropped into the test tube covered with 3 ml more of trichloroacetic acid and allowed to freeze solidly.
- 5 The tube is dipped in cold water until melting just begins. The pestle is inserted, the motor started and the muscle is ground to a fine pulp. This takes about 2 minutes.
- 6 The pestle is withdrawn and washed off quantitatively into the test tube with ice-cold trichloroacetic acid.
- 7 Dilute to the 25 ml mark, insert a clean rubber stopper and shake.
- 8 Centrifuge for 5 minutes at 1000 rpm.
- 9 Pipette out an aliquot and proceed as for whole blood.

### 4 Glycogen in Tissue

#### Reference

- Good C. A., Kramer H. and Somogyi M. The Determination of Glycogen. *J. Biol. Chem.* 100: 485-491 (April) 1933.
- Association of Official Agricultural Chemists. *Methods of Analysis*, ed. 5. Washington D. C. 1940. Assoc. of Official Agricultural Chemists.

**Principle**

The tissue is hydrolyzed with boiling KOH which does not destroy glycogen. The glycogen is precipitated with ethyl alcohol and is hydrolyzed with boiling HCl. Sugar is then estimated by any appropriate macro method such as that of Shaffer Hartman Somogyi.

**Apparatus**

- 1 Conical 15 ml graduated centrifuge tubes
- 2 Hollow glass balls which fit the necks of the centrifuge tubes. These act as reflux condensers

**Reagents**

- 1 Potassium hydroxide (KOH) approx 54 N
- 2 Potassium hydroxide approx 0.01 N
- 3 Hydrochloric acid (HCl) 0.55 N. For this solution take 31.4 ml of 35% acid and dilute it to 500 ml with distilled water
- 4 Thymol blue approx 1% aqueous solution
- 5 Ethyl alcohol (EtOH) 95%
- 6 Glycogen standard solution 80 mg/100 ml. Weigh out 0.08 gm glycogen and dissolve in 100 ml of 20% KOH solution. This solution is stable in icebox for at least one week
- 7 Shaffer Hartman Somogyi sugar reagents

**Procedure**

- 1 Place 1-2 ml (enough to cover the muscle) of either cold or hot 20% KOH in the tube fitted with the air condenser
- 2 After weighing the frozen muscle (about 0.5 gm) drop it into the tube immediately replace the air condenser and place the tube in a boiling water bath. Boil for 30 minutes occasionally removing the tube and shaking it to aid the hydrolysis
- 3 Cool the tube to room temperature and add 95% alcohol. The volume added should equal 11 to 12 times the volume of solution already in the tube. Stir the solution vigorously with a platinum wire until a complete mixing results. Any glycogen will begin to precipitate here. Place the tube in warm water. Slowly and carefully bring the alcoholic solution to the boiling point and immediately remove the tubes from the water bath. Cool to room temperature. This completes the coagulation of the glycogen
- 4 After the precipitate has been allowed to settle somewhat note carefully the amount present. This value will be used later in the procedure
- 5 Remove the condensers and centrifuge the tube at the maximum speed for 15 minutes. Decant the solution and drain the tubes for a few minutes. If the rest of the estimation must be delayed this is the best place to stop. Place the tubes inverted in the icebox until the estimation can be resumed.
- 6 Remove the remaining alcohol by placing the tubes in boiling water or in an oven at 100° for about 5 minutes or until no more alcohol can be detected by its odor. This is not necessary if the tubes have stood in the icebox overnight
- 7 Cool and add 22% HCl in amounts according to the amount of precipitated glycogen noted earlier

|                              |              |
|------------------------------|--------------|
| Use about 6 ml HCl for       | 0.304 ml ppt |
| Use about 5 to 5½ ml HCl for | 0.203 ml ppt |
| Use about 4 ml HCl for       | 0.102 ml ppt |

Never use less than 1.5 to 2 ml HCl. Place the tubes with attached air condensers in boiling water and hydrolyze for 2½ hours. Cool to room temperature

- 8 Add 1 drop of thymol blue to each tube. If the solution is not to be estimated for over 5 hours adjust the pH in the range of 2 to 6 by adding KOH from a dropper until the indicator is yellow in color and place the covered tubes in the icebox until the determination can be resumed. Then adjust the pH of the solution to 8 with 30% KOH. 0.5% KOH and 2% HCl. The end point is reached when the indicator has just reached a pale purple color.
- 9 Estimation of glucose by Shaffer Hartman Somogyi method. This method is recommended because it covers a wide range of glucose concentrations without necessitating many dilutions.
- Adjust the pH to 8 as in Step 8.
  - Dilute to a measured volume e.g. 5 ml and use aliquots containing 0.25-1.50 mg glucose. The volume of aliquot should be around 5 ml.
  - Measure glucose by the Shaffer Hartman Somogyi method.
- 10 Estimation of Glucose by Method of Folin and Wu. This method may be preferred by those who use it routinely for blood.
- Adjust pH to 8 as in Step 8.
  - Dilute the solution to a known volume with distilled water according to the amount of precipitated glycogen noted earlier. If the precipitate was greater than 0.4 ml dilute to 15 ml and use 1 ml of the solution and 1 ml of distilled water for the sugar determination.
- If ppt was 0.3-0.4 ml make up to 15 ml and use 2 ml  
 If ppt was 0.2-0.3 ml make up to 15 ml and use 2 ml  
 If ppt was 0.1-0.2 ml make up to 8-10 ml and use 2 ml  
 If ppt was 0.1 or less make up to 6-8 ml and use 2 ml
- Proceed as described under the method of Folin and Wu.

#### Calculation

$$1 \text{ mg glucose/100 gm tissue} = \frac{(\text{mg glucose/ml hydrolysate}) \times (\text{dilution}) \times (100)}{(\text{ml aliquot}) \times (\text{gm tissue})}$$

For the Shaffer Hartman Somogyi method an empirical calibration curve is required. Set up glycogen standards at intervals in the range 0.25 mg-1.50 mg/5 ml run through the hydrolysis technique previously described neutralize and estimate the glucose.

#### Example

A sample of muscle weighed 26.0 mg. After hydrolysis the solution was diluted to 25 ml and 5 ml aliquots were analyzed for glucose. The average concentration was found to be 0.16 mg glucose/ml of dilute hydrolysate.

$$\text{mg glycogen (as glucose)/100 gm tissue} = \frac{(0.16) \times 5 \times 100}{5 \times 0.6} = 30.0 \text{ mg/100 gm}$$

#### Precautions

- The chief source of error in this method is glycogenolysis in the tissue before the glycolytic system can be inhibited or inactivated. Extreme care is required in handling the tissue.

## 5 Alcohol in Blood and Urine

#### Reference

Gibson J. O. H. and Blodner H. The Determination of Ethyl Alcohol in Blood and Urine With the Photoelectric Colorimeter. *J. Biol. Chem.* 126: 5-1559 (Nov) 1938.

#### Principle

Alcohol is separated from blood filtrate by distillation. Alcohol is oxidized by potassium dichromate causing a loss of color.





**Apparatus**

- 1 An alcohol still (or a Kjeldahl ammonia still)
- 2 A quantity of 50 ml round bottom test tubes with a number of #1 solid rubber stoppers
- 3 An electric centrifuge
- 4 A quantity of 25 ml glass stoppered graduated cylinders
- 5 A Coleman Junior Spectrophotometer
- 6 Cuvettes 19 x 100 mm
- 7 An accurately calibrated 2 ml syringe pipette

**Reagents**

- 1 Potassium dichromate ( $K_2Cr_2O_7$ ) Dissolve 3.1943 gm pure salt in 1000 ml water
- 2 Sodium tungstate ( $Na_2WO_4$ ) Dissolve 100 gm to make 1 liter
- 3 Sulfuric acid ( $H_2SO_4$ ) Prepare 2/3 N  $H_2SO_4$  by diluting 6.6 ml of concentrated acid to 360 ml with water
- 4 Sulfuric acid ( $H_2SO_4$ ) Concentrated and 18 N both used
- 5 Potassium oxalate ( $K_2C_2O_4$ ) Prepare a saturated solution in water One drop will prevent clotting in 10 ml of blood
- 6 Standard anhydrous aldehyde free ethyl alcohol Mix together one liter absolute EtOH 0.5 gm NaOH and .5 gm  $AgNO_3$  dissolved in 3 ml water Reflux 2 to 3 hours Distill and store first 800 ml in brown bottle for use as standard

**Procedure**

- 1 Pipette exactly 2 ml of oxalated blood into 14 ml water in a 50 ml centrifuge tube
- 2 Add 2 ml of 2/3 N sulfuric acid mix
- 3 Add 2 ml of 10% sodium tungstate
- 4 Mix and centrifuge for 10 minutes
- 5 Place in a 100 ml Kjeldahl flask 10 ml of blood filtrate For urine use 1 to 10 ml depending on the concentration of alcohol suspected and add water to make 10 ml in the flask
- 6 Distill off 6 to 7 ml into 25 ml glass stoppered cylinder containing 5 ml of 18 N sulfuric acid and 2 ml of potassium dichromate The volume in the cylinder should be about 15 ml after distillation The cylinder is cooled in an ice bath during distillation
- 7 Mix and after sealing glass stopper with a drop of conc sulfuric acid heat cylinder and contents in a water bath for 20 minutes at 75-80 C
- 8 Then cool to room temperature and dilute to exactly 20 ml
- 9 Mix and decant into a colorimeter tube
- 10 At the same time prepare a reagent blank by mixing 5 ml 18 N  $H_2SO_4$  13 ml water and 2 ml of potassium dichromate
- 11 Set colorimeter at 440 m $\mu$  and a water tube at 100% T

**Calculation**

- 1 Prepare a calibration curve from standard alcohol solutions containing 0 to 150 mg per 100 ml at intervals of 15 mg per 100 ml and running them through the whole method Note that in this case the higher the % T the more concentrated is the alcohol
- 2 mg alcohol/100 ml fluid = 
$$\frac{(\text{mg alcohol in cuvette}) \times (\text{dilution}) \times (100)}{(\text{ml aliquot}) \times (\text{ml fluid sample})}$$

**Example**

A reading of 60% T was equivalent to 0.15 mg alcohol in the cuvette During the analysis 2 ml blood were diluted to 20 ml of which 10 ml aliquots were analyzed

$$\text{mg alcohol/100 ml blood} = \frac{(0.15) \times 20 \times 100}{10 \times 2} = 150 \text{ mg/100 ml}$$

**Precautions**

- 1 In this method the commonest sources of error are
  - a Loss of alcohol during distillation
  - b Erratic color from failure to observe conditions of oxidation exactly
- 2 The alcohol standard is critical it must be water and aldehyde free

**6 Total Bisulfite Binding Capacity****Reference**

- Thompson R H S and Johnson R E Blood Pyruvate in Vitamin B Deficiency  
 Biochem J 29 674 700 (May) 1935

**Principle**

On occasion as in exercise or in thiamine deficiency an increase in the keto acids of blood is observed. The present method gives a measure of the increase of total keto acids in the blood by virtue of their capacity to combine with bisulfite in acid solution. These complexes are broken down in alkaline solution and the bisulfite is titrated iodometrically. The present method is not useful in the presence of ketosis.

For specific compounds such as pyruvic acid and alpha ketoglutaric acid the reader should use methods recommended in the bibliography of the succeeding pages.

**Apparatus**

- 1 An electric centrifuge
- 2 An accurately calibrated 3 ml burette
- 3 A number of 5 ml round bottom tubes
- 4 An ice bath
- 5 A number of 25 ml volumetric flasks
- 6 A number of 15 ml round bottom centrifuge tubes

**Reagents**

- 1 Trichloroacetic acid a 4% solution in water
- 2 Thymol blue indicator
- 3 Sodium meta bisulphite a 5% solution in water
- 4 Starch indicator (see lactate method)
- 5 Standard iodine 0.005 N (see lactate method)

**Procedure**

- 1 Pipette 9 ml of 4% trichloroacetic acid into a 15 ml centrifuge tube
- 2 Add exactly 1 ml of blood and mix by inversion
- 3 Centrifuge for 5 minutes
- 4 Transfer 8 ml to 25 ml volumetric flask and adjust pH to 7.30 using thymol blue (orange color)
- 5 Dilute to 15 ml with water
- 6 Pipette out 10 ml samples into lactate titration tubes
- 7 Add 1 ml of sodium bisulfite (5%)
- 8 Allow to stand at room temperature for 15 minutes
- 9 Titrate in ice bath as in lactate. The end point is difficult
- 10 Set up a blank from trichloroacetic acid treated as blood filtrate

**Calculation**

$$\text{Bisulfite binding capacity expressed as mg pyruvic acid/100 ml blood} = \frac{(\text{ml titration} - \text{ml blank}) \times (45) \times (\text{normality of iodine}) \times \text{dilution} \times 100}{(\text{ml aliquot}) \times (\text{ml blood used})}$$

**Example**

A titration of 0.05 ml and a blank of 0.04 ml were obtained. Of blood 1 ml had been precipitated in 10 ml. 8 ml of this were diluted to 25 ml and 10 ml aliquots were titrated. Iodine was 0.005 N.

Bisulfite binding capacity expressed as mg pyruvic acid =

$$\frac{(0.25 - 0.04) \times 45 \times 0.005 \times 25 \times 10 \times 100}{10 \times 8 \times 1} = 147$$

**Precautions**

The end point is difficult and titration must be carried out in the cold.

**7 Carbohydrates Miscellaneous References****Aldehydes**

Siggia S and Maxcy W. Improved Procedure for Determination of Aldehydes. *Anal Chem* 19 10<sup>2</sup>3 10<sup>2</sup>5 (Dec) 1941.

**Citric Acid**

Goldberg A S and Beraheim A R. Citric Acid Determination. *J Biol Chem* 156 33 46 (Nov) 1944.

Natelson S, Lugovoy J H and Pincus J B. Determination of Micro Quantities of Citric Acid in Biological Fluids. *J Biol Chem* 170 597 606 (Oct) 1947.

Saffran M and Denstedt O F. A Rapid Method for the Determination of Citric Acid. *J Biol Chem* 175 849 855 (Sept) 1948.

Tausky H H. A Micro Colorimetric Method for the Determination of Citric Acid. II. A Note on the Substitution of Ferrous Sulphate for Hydrazine Sulphate as the Reducing Agent. *J Biol Chem* 181 195 198 (Nov) 1949.

Wolcott G H, Boyer P D. A Colorimetric Method for the Determination of Citric Acid in Blood and Plasma. *J Biol Chem* 172 729 736.

**Glucose**

Benedict, S R. The Analysis of Whole Blood. II. The Determination of Sugar and Saccharoids (Non fermentable Copper Reducing Substances). *J Biol Chem* 92 141 159 (June) 1931.

Fohn O and Wu H. A System of Blood Analysis. *J Biol Chem* 33 106 110 (May) 1919.

Hall D A. A Simplified Method of Blood Sugar Estimation. *Brit. M J* pages 351 353 (Feb) 1950.

Kingslev G R and Reinhold J G. The Determination of True Glucose in Blood by Reduction of Ferriyanaide. *J Lab & Clin Med* 34 (5) 713 719 (May) 1949.

Pohls B D, and Sortwell M. Rapid Photocolorimetric Micro Procedure for Blood Sugar Using Copper Reduction With Perchloric Acid Deproteinized Filtrates. *Arch Biochem* 11 (2) 233 (Oct) 1946.

Schaless O and Schaless S S. A Simple Method for Determination of Glucose in Blood. *Arch Biochem* 8 (2) 285 290 (Nov) 1945.

**Glycogen**

Wagner R. The Estimation of Glycogen in Whole Blood and White Blood Cells. *Arch. Biochem* 11 (2) 249 258 (Oct) 1946.

**Inulin**

Little J M. A Modified Diphenylamine Procedure for the Determination of Inulin. *J Biol Chem* 180 747 754 (Sept) 1949.

Robson J S, Ferguson M H, Olbrich O and Stewart C P. The Determination of the Renal Clearance of Inulin in Man. *Quart J Exper Physiol* 35 111 134 (June) 1949.

Roe J H, Epstein, J H and Goldstein N P. A Photometric Method for the Determination of Inulin in Plasma and Urine. *J Biol Chem* 178 839 845 (April) 1949.

**Lactic Acid**

Friedemann T E and Graesser J B. The Determination of Lactic Acid. *J Biol Chem* 100 291 303 (March) 1933.

Milton R. Notes on Mendel and Goldscheider's Method for Determining Lactic Acid in Blood. *Analyst* 61 91 96 (Feb) 1936.

**Mannitol**

Corcoran A C and Page I H A Method for the Determination of Mannitol in Plasma and Urine *J Biol Chem* 170 165 171 (Sept ) 1947

**Methylpentose**

Dische Z and Shettles L B A Specific Color Reaction of Methylpentoses and a Spectrophotometric Micro method for Their Determination *J Biol Chem* 175 595 603 (Sept ) 1948

**Pentose**

Dische Z Spectrophotometric Method for the Determination of Free Pentose and Pentose in Nucleotides *J Biol Chem* 181 319 329 (Nov ) 1949

Dunstan S and Gillam A E A Micro Method for the Determination of Pentoses by Photoelectric Spectrophotometry *J Chem Soc Suppl Issue No 1* 8140 144 1949

Roe J H and Rice E W A Photometric Method for the Determination of Free Pentoses in Animal Tissues *J Biol Chem* 173 50 51 (April) 1948

**Pyruvic Acid**

Klem, D The Determination of Pyruvic Acid in Blood in Presence of Acetoacetic Acid *J Biol Chem* 137 311 316 (Jan ) 1941

Tsao M and Brown S Pyruvic Acid Determination A Micro Method *J Lab & Clin Med* 35 (2) 303 303 (Feb ) 1950

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### D FATS AND DERIVATIVES

##### 1 Total Lipids in Serum

###### Reference

Bloor W R The Determination of Small Amounts of Lipid in Blood Plasma  
J Biol Chem 77 53 73 (April) 1928

###### Principle

The lipids are extracted with ether alcohol and saponified with sodium ethylate. Fatty acids and cholesterol are extracted with petroleum ether and oxidized with silver and potassium dichromate. The excess dichromate is titrated with sodium thio sulfate.

###### Apparatus

- 1 A quantity of 25 and 50 ml volumetric flasks
- 2 A steam bath
- 3 Special quantitative filter paper (Whatman 40)
- 4 A quantity of 125 ml glass stoppered Erlenmeyer flasks
- 5 An electric oven set at 124 C

###### Reagents

- 1 Ether alcohol mixture 3 volumes of 95% ethyl alcohol ( $C_2H_5OH$ ) and 1 volume of ethyl ether ( $C_2H_5O$ )
- 2 Sodium thiosulfate ( $Na_2S_2O_3$ ) 0.1000 N (See section on volumetric methods.)
- 3 Potassium iodide (KI) a 10% solution made freshly
- 4 Starch suspension (See lactate method)
- 5 Silver dichromate ( $Ag_2CrO_4$ ) Dissolve 5 gm of silver nitrate in 25 ml of water. Add 5 gm of  $K_2Cr_2O_7$  dissolved in 50 ml of water. Stir and centrifuge. Wash precipitate twice with water to remove  $NO_3^-$ . Dissolve in 500 ml of concentrated  $H_2SO_4$ .
- 6 Petroleum ether b.p. 40-60 C
- 7 Sodium ethylate ( $C_2H_5ONa$ ) Introduce 3 gm of fresh sodium ethylate into 100 ml of absolute  $EtOH$  containing 1 ml of 95% alcohol ( $EtOH$ ). Cool. Store in cool dark place.
- 8 Dilute sulfuric acid To one volume concentrated  $H_2SO_4$  add 3 volumes of water.
- 9 Potassium dichromate ( $K_2Cr_2O_7$ ) 1.0000 N (See section under standardization)
- 10 Methyl red indicator (0.1% in water)

###### Procedure

- 1 Into 40 ml of the ether alcohol mixture in a 50 ml volumetric flask pipette exactly 3 ml of serum or plasma.
- 2 Let stand 1 hour. Make up to 50 ml with ether alcohol.
- 3 Mix. Filter through filter paper.
- 4 Add 0 ml of filtrate to 1-5 ml Erlenmeyer flask.
- 5 Add 2 ml of sodium ethylate.
- 6 Evaporate on steam bath to disappearance of alcohol odor. Do not overheat.
- 7 Remove last traces of alcohol by stream of air.
- 8 To the residue add 1 ml of dilute sulfuric acid.
- 9 Heat mixture on steam bath for one minute.

- 10 Pour 10 ml petroleum ether onto hot mixture and allow to come to a boil
- 11 Rotate flask gently  $\approx$  3 minutes
- 12 Add 1 drop of methyl red to see the separation between ether-water layer
- 13 Pour only hot ether extract into 25 ml volumetric flask
- 14 Repeat extractions with 3 to 5 ml portions of petroleum ether combining the extracts
- 15 Cool and dilute to 5 ml with petroleum ether. Mix
- 16 Pipette 10 ml aliquots of the extract into 15 ml glass stoppered Erlenmeyer flasks
- 17 Evaporate remove last traces of ether with a stream of air
- 18 Add 5 ml of  $\text{Ag}_2\text{CrO}_4$  reagent and 3 ml of 10000 N  $\text{K}_2\text{CrO}_7$  rotating flask while adding reagents
- 19 Heat loosely stoppered flask at 140 C in an electric oven for 30 minutes a blank being treated with  $\text{Ag}_2\text{CrO}_4$  and  $\text{K}_2\text{CrO}_7$  in the same manner
- 20 Stir solutions and again leave at 124 C for 40 minutes more. If at end of 5 minutes the solution is green add more of the dichromates
- 21 To hot flasks add 75 ml of water
- 22 Add 10 ml of KI. Do not stir
- 23 Titrate with thiosulfate until the yellow color fades. Add starch indicator and titrate to the disappearance of blue color. The flask may be agitated vigorously as the iodine is taken up by the thiosulfate

#### Calculation

- 1 This method measures principally palmitic, oleic and stearic acids and cholesterol. A factor of 1 mg total fat = 0.360 ml 1 N dichromate has been determined empirically for the mixture typical of plasma

mg total lipids/100 ml plasma =

$$\frac{(\text{ml thiosulfate blank} - \text{ml unknown}) \times (\text{dilution}) \times (\text{dilution 1}) \times (100)}{(3.60) \times (\text{ml aliquot}) \times (\text{ml aliquot 1}) \times (\text{ml plasma used})}$$

#### Example

Three ml plasma were extracted in 50 ml. 20 ml aliquots were saponified and diluted to 5 ml and 10 ml aliquots were oxidized. Titration of blank was 33 ml 0.1 N thio-sulfate that of unknown was 26.0 ml

$$\text{mg total lipids/100 ml plasma} = \frac{(33 - 26) \times 5 \times 50 \times 100}{3.60 \times 10 \times 10 \times 3} = 406 \text{ mg/100 ml}$$

#### Precautions

Being rather empirical this method requires rigid adherence to time-temperature specifications

## 2 Total Cholesterol

#### Reference

Modification of Bloor W. R. The Fatty Acids of Blood Plasma J Biol Chem 56 711, 4 (July) 1933

#### Principle

For some clinical purposes the quick method recommended here for determining total cholesterol will prove satisfactory. Cholesterol is extracted into a mixture of ether and alcohol which is then evaporated. A greenish blue color is developed from the residue in the presence of chloroform, acetic anhydride and sulfuric acid.

#### Apparatus

- 1 A quantity of 15 ml round bottomed test tubes
- Calibrated 15 ml conical test tubes

- 3 0.2 ml accurately calibrated syringe pipette
- 4 A number of 10 ml pipettes
- 5 Graduated pipettes
- 6 Cork and rubber stoppers
- 7 Coleman Jr Spectrophotometer Model 6
- 8 Cuvettes 19 x 150 mm
- 9 An electric centrifuge
- 10 A variety of #4 cork stoppers
- 11 A boiling water bath

### Reagents

- 1 Mixture of ether and alcohol. 3 volumes of anhydrous ethyl alcohol ( $C_2H_5OH$ ) and 1 volume of ethyl ether ( $C_2H_5)_2O$
- 2 Chloroform ( $CHCl_3$ ) dry Merck's reagent
- 3 Acetic anhydride ( $CH_3CO$ )<sub>2</sub>O Merck's reagent
- 4 Sulfuric acid ( $H_2SO_4$ ) conc ACS
- 5 Stock standard cholesterol 400 mg cholesterol in 100 ml chloroform
- 6 Working standard cholesterol dilute stock standard 1 to 100 with chloroform  
10 ml = 0.4 mg cholesterol

### Procedure

- 1 To 9.8 ml of ether alcohol in a calibrated cone tipped centrifuge tube add exactly 0.2 ml serum. Stopper with a rubber stopper
- 2 Allow to rest on side 60 minutes with occasional shaking or place in a rotator for 30 minutes
- 3 Adjust volume to exactly 10 ml with ether alcohol mix, install tight cork stopper and centrifuge 5 to 10 minutes. Check the volume after centrifuging
- 4 Place an 8 ml aliquot of the ether alcohol supernatant fluid in a 15 ml round bottomed test tube
- 5 Evaporate to dryness in a water bath. AVOID OPEN FLAME. This takes about 10 minutes. Cool
- 6 Add 10 ml chloroform, 4 ml acetic anhydride and 0.2 ml sulfuric acid. Stopper with rubber stopper. Mix
- 7 Stand 10 minutes at  $25^\circ C \pm 1^\circ C$  preferably in water bath
- 8 At the same time set up a reagent blank starting with steps 1 and omitting serum.
- 9 Transfer unknowns and blanks to cuvettes and read at 660 m $\mu$  with the blank set at 100% T

### Calculation

- 1 Set up a calibration curve covering the range 0 to 0.5 mg cholesterol in the cuvette with samples of standard cholesterol run through the whole method
- 2  $\text{mg cholesterol/100 ml serum} = \frac{(\text{mg cholesterol in cuvette}) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml serum used})}$

### Example

A reading of 72% T was equivalent to 0.15 mg cholesterol in the cuvette. For analysis 0.2 ml of serum was extracted in 10 ml solvent, of which 8 ml was analyzed

$$\text{mg cholesterol/100 ml serum} = \frac{(0.15) \times 10 \times 100}{8 \times 0.2} = 93.7 \text{ mg/100 ml}$$

### Precautions

- 1 This method is recommended only for clinical purposes not for cases in which high precision and fractionation are required
- 2 Beware of open flames at all stages

### 3 Free and Combined Cholesterol in Blood

#### References

- Schoenheimer R and Sperry W M A Micro Method for the Determination of Free and Combined Cholesterol J Biol Chem 106 ,45 760 (Sept) 1934  
 Sperry W M and Brand F C The Colorimetric Determination of Cholesterol J Biol Chem 150 315 4 (Oct) 1943

#### Principle

Cholesterol exists in plasma partly free and partly esterified or combined. Differential precipitation with digitonin enables the two forms to be calculated by difference when measured colorimetrically in acetic anhydride and sulfuric acid.

#### Apparatus

- 1 Filter paper quantitative—extracted with hot alcohol until completely free of sterols
- 2 Preserving jars either pint or quart with rubber gaskets
- 3 A dark cabinet containing a water bath. A wooden packing box equipped with door or curtain. The cabinet is fitted with a pan or tray about 4 inches deep. The larger the water bath the easier it will be to control the temperature. Provisions should be made for a thermometer and a funnel for the addition of water by drilling holes in the top of the box.
- 4 Cvettes 19 × 150 mm
- 5 Coleman Jr Spectrophotometer Model 6
- 6 A quantity of 15 ml round bottom centrifuge tubes
- 7 A 110°C oven
- 8 An incubator set at 7°C
- 9 A quantity of 5 ml volumetric flasks
- 10 A quantity of stirring rods
- 11 A steam bath

#### Reagents

- 1 Acetone alcohol mixture. One volume of redistilled acetone is mixed with one volume of absolute alcohol.
- 2 Ether peroxide free
- 3 Acetone ether mixture. One volume of redistilled acetone to two volumes of peroxide free ether.
- 4 Digitonin solution. Dissolve 400 mg in 100 ml of distilled water by heating on a steam bath. The solution must be clear. A sediment will develop on standing but this can be removed by filtration or centrifuging before use.
- 5 Potassium hydroxide solution 33% in water
- 6 Phenolphthalein solution a 1% solution in alcohol
- 7 Acetic acid solution. Dilute 10 ml acetic acid up to 100 ml with distilled water.
- 8 Glacial acetic acid (ACS or reagent grade)
- 9 Acetic anhydride 99 to 100% chloride free
- 10 Sulfuric acid concentrated
- 11 Stock cholesterol solution in glacial acetic acid. 100 mg of cholesterol per 100 ml of glacial acetic acid. The cholesterol must be pure white and odorless and it must melt sharply at not below 147°C. If not recrystallize from an anhydrous solvent such as absolute alcohol, methyl alcohol or ethylene chloride.
- 12 Use a working standard of 0.1 mg/ml. Dilute 10 ml of the stock cholesterol to 100 ml with glacial acetic acid.
- 13 Acetic anhydride-sulfuric acid mixture. Add acetic anhydride (enough for a series of determination) into a glass stoppered flask chill it in an ice bath and add concentrated sulfuric acid in the proportion of 1 ml to 20 ml of anhydride. agitate the contents during the addition while the flask is kept in the ice bath. Insert the stopper remove from the bath shake vigorously for a few moments and return to the bath.



**Procedure****1 Preparation of Sample**

- a Pipette 7 to 10 ml of acetone alcohol solution into a 20 ml volumetric flask and pipette exactly 1 ml of serum into the solution. Mix after the addition of the serum by swirling the solution.
- b Heat the flask over a steam bath and rotate the flask to prevent bumping.
- c When the solution boils, cool the flask to room temperature and dilute to the mark with acetone alcohol solution.
- d Filter into a clean flask.
- e Pipette samples for free cholesterol and total cholesterol at once to avoid concentration changes due to evaporation.

**2 Isolation of Free Cholesterol**

- a Pipette 7 ml of the filtrate into a 15 ml centrifuge tube, add approximately 3.5 ml of digitonin solution and one drop of 10% acetic acid solution.
- b Stir thoroughly and leave rod in the tube. Place the tube in a preserving jar, cover tightly and leave it overnight at room temperature.
- c Transfer tube to a rack, stir to break up clumps of precipitate and to free particles which may adhere to the sides of the tube near the surface of the liquid. Remove the rod and place it in the rack, making certain that no precipitate on rod is lost.
- d Centrifuge until the precipitate is packed tightly enough to permit decanting with out any loss. Then drain the tube for a few moments by touching the lip to a clean towel.
- e Replace the stirring rod and wash down the wall of the tube and rod with approximately 4 ml of acetone ether solution.
- f Stir thoroughly, return rod to rack, centrifuge and decant again.
- g Wash the precipitate twice more with ether instead of the acetone ether mixture.
- h Replace the rod and set the tube aside for color development. (Be sure that the ether has been completely evaporated before proceeding.)

**3 Total Cholesterol Digitonide**

- a Add 3 drops of the 33% KOH solution to a 15 ml centrifuge tube, pipette 3 ml of the extract into the tube and stir at intervals with a vigorous up and down motion of a stirring rod until no droplets of the alkali can be seen at the tip of the tube.
- b Place a layer of sand about 5 cm deep in a preserving jar and heat it in a water bath until the temperature of the sand is about 40°C. Place the tube in the sand, cover the jar tightly and put it in an incubator 37 to 40°C for 30 minutes.
- c Remove the tube to a rack, allow to cool, raise the stirring rod and add alcohol acetone to the 6 ml mark.
- d Then add 1 drop of phenolphthalein solution and titrate with 10% acetic acid solution. About 0.4 ml should be required. Add one drop in excess.
- e Add 3 ml of the digitonin solution, stir thoroughly and place the tube in a preserving jar. Cover it tightly and leave it at room temperature for at least 3 hours, preferably overnight. Carry out the centrifuging and washing as described for free cholesterol, except that in this case only one ether washing is necessary.

**4 Development of Color in All Specimens**

- a Heat a shallow pan containing a layer of sand (about 3 cm deep) to 110 to 115°C in an oven. Place the tubes containing the precipitate in the sand in the order in which readings are to be taken and return the pans to the oven for 30 minutes.
- b While the precipitates are drying in the oven, adjust the water bath to 5°C and maintain it at this temperature during the rest of the procedure.
- c After 30 minutes remove the sand bath from the oven and pipette 2 ml of glacial acetic acid into the first tube while it is still in the hot sand.
- d Allow the acid to wash down the rod and the wall of the tube. Stir the contents vigorously and have the tube containing the rod in the sand while acid is being added to the next two or three tubes (2 or 3 minutes in all).

- e Stir again remove the tube from the sand allow to cool and place it in a rack in the water bath Continue in this manner until acet c acid has been added to all the centrifuge tubes in the series to be read and all have been placed in the water bath
- f In each series of determinations run at least two standards Pipette 1 ml of stand ard cholesterol (0.1 mg per ml) into centrifuge tubes insert stirring rods and place in the water bath along with the unknowns Place one at the beginning and one at the end of the series of unknowns
- g Start a stopwatch Then p ette 2 ml of glac al acetic acid and 4 ml of the mixed reagent into a cuvette mix thoroughly and use as a blank
- h After 9 minutes remove one of the standards from the water bath wipe it dry and pipette 4 ml of the mixed reagent into it
- i Stir vigorously for a few moments and pour into a cuvette Leave in water bath for 3 minutes
- j Then add 4 ml of mixed reagent to all the tubes shake vigorously and note time Pour into a cuvette and leave in a bath for 3 minutes
- k Set the blank at 100% T (or D = 0) at 6.5  $\mu$  Then remove each tube from the bath dry with a clean dry towel and read (be careful not to stir up the small quantity of sediment at the bottom of the tube)

#### Calculation

$$\text{mg } \propto \text{ total cholesterol} = C_T = \frac{166.7}{D} \times D_u$$

$$\text{mg } \propto \text{ free cholesterol} = C = \frac{71.4}{D} \times D_u$$

$$\text{mg } \propto \text{ combined cholesterol} = C = C_T - C$$

$$\text{Per cent cholesterol combined} = \frac{100 \cdot C}{C_T}$$

Where  $D_s$  is the mean optical density of the standards and  
 $D_u$  is the density of the unknown

#### Precautions

- 1 The centrifuging for free should be for 15 minutes at 800 rpm at the first stage For total 5 to 7 minutes is usually enough for the acetone ether and ether washings  
 The acetic anhydride reagent should be made up for each series of determinations Some heat is produced when mixed and a period of 10 minutes should be allowed before the reagent is used as an equilibrium reached at that time
- 2 The water bath should be maintained at 30°C
- 3 Be careful of water as it interferes with the reaction (avoid splashing)

## 4 Lipid Phosphorus in Blood and Serum

#### Reference

Youngburg G E and Youngburg M V Phosphorus Metabolism I A System of Blood Phosphorus Analysis J Lab & Clin Med 16 159 166 (Nov) 1930

#### Principle

The extracted lipids are oxidized with sulfuric acid and hydrogen peroxide and the phosphate present is measured by the method of Fiske and Subbarow

#### Apparatus

- 1 Large test tubes graduated at 10 ml
- 2 A boiling water bath
- 3 A quantity of filter paper (quantitative)
- 4 A quantity of broken sea ware
- 5 Wire racks

- 6 An electric hot plate or microburners
- 7 A quantity of 5 ml volumetric flasks
- 8 Other apparatus as in phosphorus method (Fiske and Subbarow)

**Reagents**

- 1 Alcohol ether mixture Mix 3 volumes of 95% redistilled alcohol and one volume of redistilled ether
- 2 Sulfuric acid 5 N Dilute 50 ml of concentrated sulfuric acid to 360 ml with distilled water
- 3 Hydrogen peroxide 30% (Superoxol)
- 4 Ammonium molybdate a 2.5 per cent solution
- 5 Amino-naphtholsulfonic acid solution (ANSA) (see phosphorus method)

**Procedure**

- 1 Into a graduated 20 ml test tube pipette 18 ml of the alcohol ether mixture
- 2 Add exactly 1 ml of plasma or serum drop by drop
- 3 Mix and place in a boiling water bath until the contents of the tube boil Allow to cool at room temperature
- 4 Dilute to 20 ml with alcohol ether mixture mix and filter through quantitative filter paper
- 5 Pipette exactly 8 ml of the filtrate into a large test tube Add a silica pebble and place in a wire rack
- 6 Heat the rack and tubes over an electric hot plate until evaporated to dryness
- 7 Then add 2.5 ml of 5 N sulfuric acid and digest over a hot plate or microburner (Place the tube in a slanting position and heat until the mixture turns brown or black)
- 8 Allow to cool and add 1 drop of hydrogen peroxide allowing it to fall directly into the mixture
- 9 Replace the tube and continue heating until colorless
- 10 When colorless cool the tube add some water heat again and transfer the contents to a 5 ml volumetric flask washing with water to about 15 ml volume
- 11 Then add 0.5 ml of ammonium molybdate solution
- 12 Add 1 ml of ANSA reagent mix and dilute to 5 ml with water
- 13 For a blank add 2 ml of water and digest with sulfuric acid and peroxide as above in Steps 7-12
- 14 Using the Fiske and Subbarow phosphorus curve set the blank at 100% T at 660 mμ and read the unknowns

**Calculation**

- 1 Use the phosphorus curve of the method of Fiske and Subbarow
- 2 mg lipid phosphorus (as P)/100 ml serum =

$$\frac{(\text{mg P/ml unknown}) \times (\text{dilution 1}) \times (\text{dilution 2})}{(\text{aliquot 1}) \times (\text{ml serum})} \times 100$$

**Example**

A reading of 85% T was equivalent to a concentration of 0.0015 mgP/ml in the cuvette In the analysis one ml of serum was extracted in 20 ml of alcohol ether of which 8 ml were finally diluted to 25 ml in the phosphorus analysis

mg lipid phosphorus (as P)/100 ml serum =

$$\frac{(0.0015) \times 20 \times 25 \times 100}{8 \times 1} = 9.4 \text{ mg P/100 ml}$$

**Precautions**

- 1 These in general are the same as for phosphorus estimations

## 5 Acetone Bodies

### A. TOTAL IN BLOOD

#### References

- Van Slyke D D Studies of Acidosis VII J Biol Chem 32 455-494 (Dec) 1911  
 Van Slyke D D and Fitz P Studies in Acidosis VIII The Determination of  $\beta$  Hydroxybutyric Acid Acetoacetic Acid and Acetone in Blood J Biol Chem 32 49-497 (Dec) 1917

#### Principle

Acetone in hot sulfuric acid solutions forms an insoluble precipitate with mercuric sulfate. Acetoacetic acid breaks down to acetone on heating and beta hydroxybutyric acid treated with potassium dichromate in hot acid solutions forms acetone. The three compounds may therefore be estimated by weighing the mercuric acetone sulfate.

#### Apparatus

1. Reflux condensers with interchangeable 500 ml Erlenmeyer flasks  
An oven set at 110°C
3. Sintered glass Pyrex crucibles with medium porosity
4. Volumetric pipettes 8 ml or an accurately calibrated syringe pipette
5. A desiccator
6. An accurate analytical balance
7. Quantitative filter paper
8. Graduated cylinders 50 ml capacity
9. Volumetric flasks 100 ml

#### Reagents

1. Copper sulfate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ) a 10% solution in water
- Calcium hydroxide ( $\text{CaOH}$ ) a 10% suspension in water
3. Mercuric sulfate ( $\text{HgSO}_4$ ) Dissolve 73 gm of red oxide of mercury in 1000 ml of 4 N sulfuric acid. Gentle warming may be required.
4. Sulfuric acid Make a 4 N solution by diluting 40 ml of concentrated  $\text{H}_2\text{SO}_4$  to 60 ml of water.
5. Potassium dichromate ( $\text{K}_2\text{Cr}_2\text{O}_7$ ) a 5% solution in water
6. Sulfuric acid ( $\text{H}_2\text{SO}_4$ ) 17 N 1.0 ml of concentrated sulfuric acid diluted to 360 ml of water.
7. A mixture of  $\text{HgSO}_4$  solution and 17 N sulfuric acid in the ratio 3 to 1 volumes  $\text{HgSO}_4$  to 10 volumes  $\text{H}_2\text{SO}_4$
8. Sodium hydroxide 5 N 40 gm of sodium hydroxide made up to 100 ml of water

#### Procedure

##### 1. Preparation of Blood

- a. Pipette exactly 8 ml of blood into 80 ml of water in 100 or 50 ml graduated cylinders. Mix.
- b. Add with mixing 16 ml of 10%  $\text{HgSO}_4$ . Allow to stand one-half hour or more.
- c. Dilute to 100 ml with distilled water. Mix and allow to stand 1 hour or overnight in refrigerator.
- d. Filter through quantitative filter paper.
- e. Use 100-125 ml of filtrate (= 5 ml blood) for analysis.

##### Preparation of Sintered Glass Filter

- a. Add 10-15 ml of aqua regia to cup and let stand overnight. Filter and wash out by running 100 ml of hot distilled water through filter. Dry 1 hour at 110°C. Cool for one-half hour in desiccator and weigh.

a *Precipitation of Mercuric Acetone Sulfate*

- a To a wide mouthed flask add 45 ml of  $\text{HgSO}_4 - \text{H}_2\text{SO}_4$  reagent
- b Add 125 ml of blood filtrate
- c Attach to a reflux condenser
- d Light a medium sized flame
- e Bring to a boil
- f Add 5 ml of potassium dichromate solution If the dichromate is added before solution comes to a boil the oxidation is to acetic acid instead of to acetone so that there is a significant loss
- g Boil for 90 minutes after adding dichromate keep boiling until just before filtering them through Pyrex crucibles
- h Filter while still very hot using suction Wash with 100-200 ml of water
- i Dry one hour at  $110^\circ\text{C}$  Weigh after cooling for one half hour

**Calculation**

- 1 In this method certain assumptions have to be made in the calculation
  - a That 1 mg acetone yields 20.0 mg precipitate, that 1 mg acetoacetic acid yields 11.4 mg precipitate and that both substances react quantitatively
  - b That beta hydroxybutyric acid yields only 75% under the conditions of this method and that 1 mg of beta hydroxybutyric acid yields 8.45 mg of precipitate
  - c That in most specimens of blood and urine the molecular proportion of beta hydroxybutyric acid is 75% of the total acetone bodies
  - d Under the conditions described the following factors are valid when the equivalent of 5 ml of blood is analyzed and calculated as acetone For total acetone bodies 1 mg of precipitate equals 1.28 mg of total acetone bodies/100 ml blood for beta hydroxybutyric acid 1 mg of precipitate equals 1.37 and for acetone + acetoacetic acid, 1 mg of precipitate equals 1.00
- 2  $\text{mg total acetone bodies (as acetone)}/100 \text{ ml blood} = (\text{mg ppt}) \times 1.28$

**Example**

The precipitate from the blood of a diabetic weighed 1.1 mg

$\text{mg total acetone bodies (as acetone)}/100 \text{ ml blood} = 1.1 \times 1.28 = 15.5 \text{ mg}/100 \text{ ml}$

**Precautions**

This method is rather empirical Hence great care must be taken to follow the instructions explicitly

**B TOTAL IN URINE****Procedure**1 *Preparation of Urine*

- a Pipette exactly 10 ml of urine into 100 ml volumetric flask Add about 30 ml water
- b Add .0 ml copper sulfate reagent
- c If not turbid add approximately 5 N NaOH drop by drop until turbidity persists
- d Add 20 ml calcium hydroxide reagent
- e Dilute to mark shake and keep in refrigerator at least one hour preferably over night
- f Filter

2 *Blank for Urine*

- a Add 25 ml filtrate to a 125 ml flask
- b Add 0.5 ml 17 N sulfuric acid
- c Boil not less than 5 nor more than 10 minutes
- d Filter and wash with 100 ml of water into 45 ml of mixed reagent ( $\text{HgSO}_4 - \text{H}_2\text{SO}_4$ )

- e Boil in condenser not more than 45 minutes Do not add dichromate
- f Remove and filter through sintered glass filter
- g Weigh

### 3 Preparation of Mercuric Acetone Sulfate

- a To a 500 ml flask add 45 ml of  $\text{HgSO}_4$  H<sub>2</sub>SO<sub>4</sub> reagent 25 ml of urine filtrate and 100 ml water
- b Proceed as under blood

### Calculations

- 1 For urine the same considerations hold as for blood For urine under the conditions of this method (i.e. the equivalent of 5 ml urine analyzed) the factors are total acetone bodies 1 mg ppt = 256 for beta hydroxybutyric acid 64 and for acetone 00
- $$\text{mg total acetone bodies (as acetone)/100 ml urine} = (\text{mg ppt} - \text{mg blank}) \times 256$$

### Example

The urine from a mild diabetic was analyzed The blank weighed 1 mg the unknown 56

$\text{mg total acetone bodies (as acetone)/100 ml urine} = (56 - 12) \times 256 = 112 \text{ mg/100 ml}$

### Precautions

Same as for blood

## C INDIVIDUAL ACETONE BODIES

### Principle

Acetone reacts directly with mercuric sulfate in hot acid solution acetoacetic acid breaks down quantitatively to carbon dioxide and acetone which reacts at once beta hydroxybutyric acid yields acetone only when an oxidizing agent (dichromate) is present

### Apparatus and Reagents

Same as for blood total acetone bodies

### Procedure

- 1 In either blood or urine determine acetone plus acetoacetic acid by running through the procedures described above but with the following alterations
  - a Do not add dichromate
  - b Boil for not less than 30 nor more than 45 minutes
- In either blood or urine determine beta hydroxybutyric acid by running through the whole procedure including dichromate after the following preliminary steps to remove acetone and acetoacetic acid
  - a The filtrates equivalent to 5 ml blood or 25 ml urine are mixed with 100 ml water and 2 ml 50% sulfuric acid.
  - b They are boiled 10 minutes in an open flask and cooled
  - c They are diluted to 150 ml with water

### Calculation

- 1  $\text{mg (acetone plus acetoacetic acid as acetone)/100 ml} = \text{mg ppt} \times 100 \text{ (for blood)}$   
 $(\text{mg ppt} - \text{mg blank}) \times 200 \text{ (for urine)}$
- 2  $\text{mg beta hydroxybutyric acid (as acetone)/100 ml} = \text{mg ppt} \times 130 \text{ (for blood)}$   
 $(\text{mg ppt} - \text{mg blank}) \times 64 \text{ (for urine)}$

**6 Fats and Derivatives, Miscellaneous References****Acetoacetic Acid**

Rosenthal S M A Colorimetric Method for the Estimation of Acetoacetic Acid in the Blood *J Biol Chem* 179 1 35 1<sup>st</sup> 44 (July) 1949

**Acetone**

Adams C A and Nicholls J R The Analysis of Mixtures Containing Acetone Ethyl Alcohol and Isopropyl Alcohol, *Analyst* 54 29 (Jan) 19 9

Behre J A A Modified Salicylaldehyde Method for the Determination of Acetone Bodies in Blood and Urine *J Biol Chem* 136 25 34 (Oct) 1940

Greenberg L A and Lester D A Micro Method for the Determination of Acetone and Ketone Bodies *J Biol Chem* 154 (1) 177 190 (June) 1944

Huckabay W B Newton C J and Mettler A V Optimum Conditions for Titrimetric Determinations Determination of Acetone, *Anal Chem* 19 838 841 (Nov) 1947

Rappaport, F and Baner B Microdetermination of Acetone and Diacetic Acid in Blood *J Lab & Clin Med* 28 1770 1772 (Nov) 1943

Rothera A C H Note on the Sodium Nitroprusside Reaction for Acetone *J Physiol* 37 491 494 (Dec) 1908

**Acetyl Acetone**

Witter R F Snyder J and Stotz E Colorimetric Determination of Acetylacetone and Related  $\beta$  Diketones *J Biol Chem* 176 493 500 (Nov) 1948

**Cholesterol**

Kingsley, G R and Schaffert R R Determination of Free and Total Cholesterol by Direct Chloroform Extraction *J Biol Chem* 180 (1) 315 3<sup>rd</sup> 3 (Aug) 1949

**Fatty Acids**

Bloor W R A Colorimetric Procedure for the Determination of Small Amounts of Fatty Acid *J Biol Chem* 170 671 6<sup>th</sup> 4 (Oct) 1947

**Formic Acid**

Bastrup J T Method for the Determination of Formic Acid in Urine *Acta Pharmacol et Toxicol* 3 303 311 1947

Grant W M Colorimetric Micromethod for Determination of Formic Acid *Anal Chem* 19 (3) 206 207 (March) 1947

**Hydroxybutyric Acid**

Shaffer P A and Marriott W M The Determination of Oxylbutyric Acid *J Biol Chem* 16 26, 80 (Sept Oct) 1914

**Phospholipid**

Egsgaard, J Determination of the Phosphatide Content of Serum *Acta Physiol Scandinav* 16 ( /3) 171 178 1948

Hack M H Estimation of the Phospholipides in Human Blood *J Biol Chem* 169 (1) 137 143 (June) 1947

Harnes A R The Colorimetric Determination of Lipoid Phosphorus in Blood *J Biol Chem* 77 403 408 (May) 19 3

**Total Fat**

Bauer F C Jr, and Hirsch E F Esterified Fatty Acid Levels of Normal Human Sera *Arch Biochem* 23 (1) 137 140 (Aug) 1949

Fowweather F S and Anderson W N A Method for the Determination of Fat in *Feces Biochem J* 40 (3) 350 351 1946

Stoddard J L and Drury I E A Titration Method for Blood Fat *J Biol Chem* 84 41 748 (Nov) 19 2

van de Kamer J H ten Bokkel Huinink H and Weyers H A Rapid Method for the Determination of Fat in *Feces J Biol Chem* 177 347 355 (Jan) 1949

Zuckerman J L Zymaris M C and Natelson S The Simple Method for the Determination of Fecal Fat and Fatty Acids *J Lab & Clin Med* 34 28, 286 (Feb) 1949

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### E VITAMINS

##### 1 Vitamin A and Carotenoids (Carr Price)

###### A. IN BLOOD AND SERUM

###### Reference

May C D Blackfan K D McCreary J F and Allen F H Jr Clinical Studies of Vitamin A in Infants and Children *Am J Dis Child* 59 1167 1184 (June) 1940

###### Principle

Materials containing vitamin A when treated by a chloroform solution of antimony trichloride form a blue color (Carr Price reaction). Experience has shown a definite association between the blue color and vitamin A content of a wide variety of materials as determined by biological assay.

###### Apparatus

- 1 A special container for the antimony trichloride. This reagent has a tendency to pick up water easily and as a result the blue color is impossible to read in the colorimeter.
- 2 Cuvettes—Coleman 10 x 75 mm.
- 3 Air pressure to evaporate the ether.
- 4 Centrifuge tubes—15 ml round bottom.
- 5 An electric centrifuge.
- 6 Coleman Jr Spectrophotometer Model 6 with the smallest cuvette holder (#6 108).  
A water bath set at 45°C.
- 8 Syringe pipettes 1 and 0.1 ml.

###### Reagents

- 1 Antimony trichloride 5% in chloroform.
- Ethyl alcohol 95%.
- Petroleum benzine (ether) USP.
- Vitamin A alcohol standard. Made up in petroleum ether one gamma per ml.
- Beta carotene standard 3 gamma/ml in petroleum ether (1 gm vitamin A alcohol is 4 500 000 International Units i.e. 1 IU is 0.6 micrograms. 1 IU is 0.6 micrograms beta carotene).

###### Procedure

- 1 To exactly 1 ml of serum add 2 ml of 95% ethyl alcohol.
- 2 Then add 1 ml of petroleum ether stopper tightly and shake vigorously and continuously for 10 minutes for complete extraction.



## 1) Biochemical Procedures

- 3 Centrifuge and pipette out accurately 1 ml of the ether layer into a small cuvette. Read immediately at 440 m $\mu$  the petroleum ether blank being set at 100% T.
- 4 The yellow color of the ether layer is due to the carotenoid pigments the bile pigments being unextracted under these conditions. The principal carotenoids in blood serum are carotene and xanthophyll. Estimation of the total amount of carotenoids is necessary in the determination of vitamin A, in order that the measurement of the total blue color formed by the Carr Price reaction may be corrected for that portion arising from the carotenoids.
- 5 Evaporate ether at 40 to 45 C to dryness by gently blowing dry air into the cuvette while it is partially immersed in water. Make sure that the ether odor is completely gone and that the cuvette is perfectly dry. Otherwise the Carr Price color will be a cloudy blue.
- 6 Place an antimony trichloride blank at 100% T at 690 m $\mu$ .
- 7 Add exactly 1 ml of antimony trichloride to the dry cuvette and read immediately. The galvanometer reading is taken at the point of maximum absorption of light that is maximum color formation. This is necessary as the blue color fades very fast ( $\approx$  40 seconds).

### Calculation

- 1 Run a calibration curve for carotenoids using beta carotene in the range 0 to 4 mcg/ml at 440 m $\mu$ .
- 2 Run a calibration curve for vitamin A alcohol in the range 0 to 15 mcg/ml at 690 m $\mu$  and following all the steps in the procedure.
- 3 
$$\text{mcg carotenoids/100 ml serum} = \frac{(\text{mcg in cuvette at 440 m}\mu) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml serum})}$$
- 4 
$$\text{mcg vitamin A alcohol/100 ml serum} = \frac{(\text{mcg in cuvette at 690 m}\mu) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml serum})} - (0.11 \times \text{mcg carotenoids})$$
- 5 The correction factor for blue color due to carotenoids is necessary in Equation 4 and was established empirically for a mixture of carotene and xanthophyll in the proportions usually present in serum.

### Example

At 440 m $\mu$  the reading was 85.0% T equivalent to 0.70 mcg beta carotene in the cuvette.

At 690 m $\mu$  the reading was 80.5% T equivalent to 0.40 mcg in the cuvette. One ml of serum was extracted into 2 ml of petroleum ether of which 1 ml was used for analysis.

$$1 \text{ mcg carotenoids/100 ml serum} = \frac{(0.70) \times (2) \times (100)}{(1) \times (1)} = 140 \text{ mcg/100 ml}$$

$$2 \text{ mcg vitamin A/100 ml} = \frac{(0.40) \times (2) \times (100)}{(1) \times (1)} - (0.11 \times 140) \\ = 80 - 15 = 65 \text{ mcg/100 ml}$$

### Precautions

- 1 Be sure the antimony trichloride contains no water otherwise the solution will be cloudy.
- 2 Be sure the sample is dry before the addition of antimony trichloride.

**B VITAMIN A IN FECES****Reference**

Friedemann T    Personal communication

**Apparatus**

- 1 Water bath set at 60 C
- 2 Electric centrifuge etc. as for serum
- 3 An analytical balance
- 4 A quantity of 100 ml volumetric flasks
- 5 A quantity of 100 ml graduated cylinders

**Reagents**

- 1 Phosphoric acid (syrupy)
- 2 Sodium sulfate (anhydrous)
- 3 Sodium hydroxide (NaOH) 10 N 40 gm/100 ml water
- 4 Caprylic alcohol
- 5 Ethyl alcohol 50% in water
- 6 Methyl alcohol 70% in water
- 7 Other reagents are the same as for serum vitamin A

**Procedure**

- 1 Weigh out 10 gm of mixed feces into a 100 ml volumetric flask
- 2 Add 50 ml of 95% alcohol 5 ml of 10 N NaOH and 1 ml of caprylic alcohol
- 3 Heat for 10 minutes at 60 C cool and dilute to 100 ml
- 4 Pipette 10 ml of diluted feces into a 100 ml cylinder
- 5 Add 10 ml of water and 10 ml of petroleum ether
- 6 Shake for 5 minutes and remove top layer to a new test tube
- 7 Add 5 ml of syrupy phosphoric acid and shake well. (The phosphoric acid takes away interfering material without removing any of the vitamin A)
- 8 Remove phosphoric acid layer and wash with 10 ml of 50% ethyl alcohol
- 9 Dry out water with sodium sulfate and decant solution to measure color
- 10 To remove carotene shake two times with 70% methyl alcohol. This step will take out the carotene and also a little of the vitamin A
- 11 Continue as in method for serum vitamin A

**Calculation and Precautions**

See vitamin A in serum

**2 Vitamin A and Carotenoids (Dichlorohydrin)****Reference**

Sobel A E and Snow S D    The Estimation of Serum Vitamin A With Activated Glycerol Dichlorohydrin *J Biol Chem* 171 211-53<sup>n</sup> (Dec) 194

**A. USING THE COLEMAN UNIVERSAL MODEL 11****Principle**

A solution of activated glycerol dichlorohydrin (GDH) when added to a chloroform solution containing vitamin A produces a violet color. Advantages over the commonly trichloride method are principally two: first the color is stable for two to 10 minutes; second the reagent is not affected by traces of moisture. The interference of vitamins D and D<sub>2</sub> ergosterol 7 dehydrocholesterol and cholesterol with the GDH reaction has been found to be negligible.

**Apparatus**

- 1 It is necessary to measure samples at 800 mμ or above. Hence the Coleman Universal Spectrophotometer Model 11 is recommended.

- 2 13 cm cuvettes and cuvette carrier obtained from the Coleman Electric Company
- 3 70 cm cuvettes obtained from the Pyrocell Manufacturing Company 07 East 84th Street New York 28 Carrier obtained from the Coleman Electric Company
- 4 555 mμ filter (filter combination 348° M608 10 mm and 5300, M954 40 mm) obtained from the Corning Glass Works Corning New York
- 5 An electric oven set at 60 C
- 6 A water bath set between 40-50 C
- 7 An electric centrifuge
- 8 Syringe pipette 1 ml accurately calibrated
- 9 A quantity of flat tipped glass stirring rods

#### Reagents

- 1 95 percent ethanol (USP)
- 2 1 N KOH in approximately 90 percent ethanol freshly prepared 5.6 gm KOH/100 ml
- 3 Petroleum ether b.p. 30-60 C
- 4 Nitrogen Ohio Chemical Company
- 5 Sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) anhydrous analytical reagent grade
- 6 Chloroform Analytical reagent grade is washed with water distilled the first and last fractions are discarded and the remainder is kept over dry  $\text{Na}_2\text{SO}_4$ .
- 7 Antimony trichloride, analytical reagent grade
- 8 Carotene 90 percent beta 10 percent alpha, from General Biochemicals Inc
- 9 Standard vitamin A solution About 40 to 50 mg of a vitamin A concentrate which comes in gelatin capsules control PC 3 from Distillation Products Inc were diluted in chloroform to give the desired concentration of vitamin A. The concentrate has  $E_{\frac{1\%}{1\text{ cm}}}$  value in absolute ethanol of 100.75. This value was supplied by the manufacturers and agreed with that found by the authors. Multiplying this value by 2000 the standard commercial conversion factor permits estimation of the vitamin A potency at 2,015,000 USP units per gm. If the extinction of the concentrate is divided by the extinction of crystalline vitamin A alcohol in ethanol 1780 the vitamin A content of the concentrate becomes by simple proportion 5.56 percent.
- 10 Activated GDH Glycerol dichlorohydrin (Eastman Kodak practical grade a mixture of 1, 3 and 2, 3 dichlorohydrin from the Shell Chemical Company) is vacuum distilled in the presence of approximately 1 percent by weight of  $\text{SbCl}_3$  at from 10 to 40 mm of pressure the first and last fractions being discarded. The distilled reagent should be colorless and clear (a sample of the reagent should not become even slightly cloudy on the addition of water). The presence of  $\text{SbCl}_3$  in the distillate (shown when a sample of GDH becomes cloudy upon the addition of water) can be remedied by a second vacuum distillation. The activated reagent may be obtained from the J. B. Shohan Laboratories 78 Wheeler Point Road Newark 5 New Jersey.

#### Procedure

- 1 Exactly 1 ml of serum is pipetted into a  $\frac{3}{4} \times 4$  inch test tube
- 2 Add 1 ml of 95 percent ethanol and mix by tapping the tube (or for saponification 1 ml of the 1 N KOH in 90 percent ethanol is added and the contents of the tube mixed and placed in a 60 C oven for 20 minutes)
- 3 Add 2 ml of analytical reagent petroleum ether and shake for 10 minutes
- 4 After shaking the tube is centrifuged for about 30 seconds
- 5 The supernatant petroleum ether is aspirated and placed in a  $\frac{1}{2} \times 4$  inch test tube. The aspirator is a fine tipped dropper.

- 6 With another 9 ml of the petroleum ether and shaking for only five minutes the extraction procedure is repeated
- 7 The extract is evaporated to dryness by placing the tube in a 40-50°C water bath and running a stream of nitrogen over it
- 8 Add 1 ml of analytical reagent grade chloroform to bring the dried extract into a solution
- 9 Then add 4 ml of GDH and mix with a flat tipped stirring rod.
- 10 Two minutes after mixture the solution is placed in a 10 mm cuvette. The absorption of the solution is measured first at 550 mμ (with the 555 mμ filter) against a blank consisting of 4 ml of GDH and 1 ml of chloroform. When the concentration of vitamin A is high read at 4 minutes
- 11 The wave length dial is then turned to 800 mμ the 555 mμ filter replaced with filter PC 5 and at 4 minutes after the initial mixture of the reagents the absorption is read at 800 mμ

#### Calculation

- 1 Prepare the calibration curve for vitamin A as follows
  - a Prepare solutions of vitamin A in chloroform solution containing standard solution as well as the GDH preheated to 25°C before use
  - b To 1 ml of the standard in a glass stoppered cylinder are added 4 ml of GDH
  - c The contents of the cylinder are mixed by inversion and are placed in a 25°C water bath for minutes. The solution is then poured into the 5 cm cuvette
  - d The absorption of the pink color produced is read against a blank consisting of 4 ml of GDH and 1 ml of chloroform. The instrument is used with the wave length scale set at 550 mμ in conjunction with a 555 mμ filter
- Prepare the calibration curve for carotene at 800 mμ as follows
  - a Prepare standards of carotene in chloroform in the range 0.10 mcg/ml
  - b Carry these through all steps of the procedure reading them against a chloroform blank at 800 mμ using filter PC 5
- 3 Prepare a curve for correcting for carotene interference at 550 mμ as follows
  - a Prepare carotene standards in chloroform in the range 0.10 mcg/ml.
  - b Run these through the vitamin A method reading them against the blank set at 100°C at 50 mμ with the 555 mμ filter
  - c Determine the interference factor of carotene at 550 mμ i.e. the apparent vitamin A in mcg due to 1 mcg of carotene/100 ml serum when measured at 550 mμ

$$4 \text{ mcg carotene/100 ml serum} = \frac{(\text{mcg in cuvette at } 800 \text{ m}\mu) \times 100}{(\text{ml serum})}$$

$$5 \text{ mcg vitamin A/100 ml serum} = \frac{(\text{mcg in cuvette at } 550 \text{ m}\mu) \times 100}{(\text{ml serum})} - (\text{Factor} \times \text{mcg carotene/100 ml serum})$$

#### Example

- 1 At 800 mμ the reading was 78°C equivalent to 16 mcg carotene in the cuvette. At 550 mμ the reading was 84°C equivalent to 1.2 mcg vitamin A in the cuvette. One ml of serum had been extracted into 4 ml of petroleum ether all of which was evaporated to dryness and treated in the cuvette with 1 ml chloroform and 4 ml of GDH

$$2 \text{ mcg carotene/100 ml serum} = \frac{(16) \times 100}{1} = 160 \text{ mcg/100 ml}$$

$$3 \text{ mcg vitamin A/100 ml serum} = \frac{(1.2) \times 100}{1} - (0.15 \times 160) = (12 - 24) = 96 \text{ mcg/100 ml}$$

**Precautions**

- 1 In this method exact timing is necessary for good results
- 2 When the concentration of vitamin A is high read at 4 minutes instead of "

**B VITAMIN A (DICHLOROHYDRIN METHOD) WITH BECKMAN SPECTROPHOTOMETER****Apparatus**

- 1 A Beckman Spectrophotometer Model DU
- 2 Four matched Corex cells with 10 mm light path absorption cells
- 3 A four place cell holder
- 4 Syringe pipettes 0.6 ml and 2.4 ml accurately calibrated
- 5 Other apparatus same as for Coleman Universal Model 11

**Procedure**

- 1 Same as procedure for Coleman Universal up to and including step 7
- 2 Add 0.6 ml of chloroform with a syringe pipette
- 3 Stir and add " 4 ml GDH
- 4 Read at 2 minutes but if the sample is high (a tolerance test) read between 3 and 4 minutes
- 5 Read at 555 m $\mu$  setting the blank at 100% T
- 6 Change wave length to 830 m $\mu$  and read again setting the blank at 100% T

**Calculation**

- 1 The calculation and calibration is the same as for the Coleman Model 11

**Precaution**

- 1 Use glass stoppered tubes for extraction
- 2 Saponification is a critical step complete saponification is essential

**3 Thiamine in Urine Food Blood and Feces****A. IN URINE****Reference**

- Hennessey D J and Cerecedo L R The Determination of Free and Phosphorylated Thiamin by a Modified Thiochrome Assay *J Am Chem Soc* 61 179 183 (Jan) 1939
- Connor R T and Straub G J Combined Determination of Riboflavin and Thiamine in Food Products *Ind and Eng Chem Anal Ed* 13 385 389 (June) 1941

**Principle**

Thiamine is adsorbed from urine in acid solution by means of activated zeolite. It is eluted with concentrated potassium chloride and is converted to thiochrome by ferricyanide in alkaline solution. Thiochrome is extracted into isobutyl alcohol and measured fluorometrically.

**Apparatus**

- 1 Photofluorometer Coleman Model 12 or 12B
- 2 Hennessey adsorption tubes or columns with long (5 to 5.5 inch) stems
- 3 Fine glass wool
- 4 Reaction vessels 25 ml glass stoppered (interchangeable stoppers) vessels with a conical bottom
- 5 Reading cuvettes Coleman 19 x 150 mm
- 6 Syringe pipettes—5ml 10 ml
- 7 Wooden applicators for inserting the glass wool into the columns

## Reagents

- 1 Potassium chloride 0.5% solution Dissolve 0.5 gm KCl in 100 ml 0.1 N HCl with the aid of heat Filter  
Bromocresol green indicator 0.4% in 70% alcohol.
- 3 Sodium acetate 2.5 M Dissolve 330 gm in 1000 ml distilled water
- 4 Decalco Wash 60 to 80 mesh Decalco 3 times with 3% acetic acid once with 0.5% KCl again with 3% acetic acid and then several times with distilled water Each washing consists of stirring the Decalco 15 minutes in the wash settling and decanting
- 5 Sodium hydroxide 1 N solution Dissolve 40 gm of sodium hydroxide in 1000 ml of distilled water
- 6 Sodium sulfate anhydrous
- 7 Acetic acid 0.5% solution in water
- 8 Isobutyl alcohol redistilled
- 9 Potassium ferricyanide a 1% solution in water This solution is stable for 6 months if kept in a dark bottle
- 10 Potassium hydroxide a 15% solution in water
- 11 Oxidizing reagent One part of 1% ferricyanide solution is mixed with 9 parts of potassium hydroxide (15%) Prepare fresh every day
- 12 Quinine sulfate standard  
Stock standard 100 mg USP quinine sulfate per liter of 0.1 N  $H_2SO_4$  Stable for 1 year  
Intermediate standard 0.5 ml stock standard dilute to 100 ml with 0.1 N  $H_2SO_4$   
Working standard 5 ml intermediate standard to 500 ml 0.1 N  $H_2SO_4$  Prepare daily
- 13 Thiamine standard (100 mcg/ml) Dissolve 50 mg USP thiamine standard in 200 ml of 95% C<sub>2</sub>H<sub>5</sub>OH Add 5 ml 0.1 N HCl and bring to 500 ml with distilled water Store in refrigerator It is stable for 6 months.

## Collection of Sample

Individual samples are preserved by adding 0.5 ml of 0.1 N  $H_2SO_4$  to each 100 ml of sample 24 hour samples are collected in amber colored 1 gallon jugs which contain 5 ml of 0.1 N  $H_2SO_4$  The vitamin content of the acidified sample will remain unchanged over several weeks if the sample is kept in the refrigerator A suitable aliquot is transferred to a 150 ml extraction flask the reaction is adjusted to pH 3.5 and the adsorption is carried out as described below

## Procedure

- 1 Place a small piece of glass wool in the tip of the aluminum adsorption column, fill the column with activated Decalco so that the bottom of the bottle is just covered
- 2 Wash the column with 3 ml 0.5% acetic acid and allow to drain
- 3 Take a 5 ml aliquot of urine in a 50 ml Erlenmeyer flask add 4 drops of bromocresol green and neutralize to a light greenish yellow color with 1 N NaOH
- 4 A blank consisting of 5 ml 0.5% acetic acid is run through the method with each set of samples Run the samples through the columns and allow to drain by gravity into a waste pan
- 5 Rinse the flasks with 3 ml 0.5% acetic acid and put this through the columns
- 6 After this has drained through the columns are completely drained by using a small amount of air pressure at the top of the column Recoveries of thiamine standard added to urine should also be run along with each batch of unknowns
- 7 Reaction vessels are placed under the adsorption columns and the columns are eluted with 5 ml 25% KCl After this has drained the columns are again blown out with a small amount of air pressure
- 8 To the 5 ml eluate 0.7 ml of oxidizing reagent is added with mixing

- 9 Then add 10 ml isobutanol stopper with the glass stopper and shake 50 times. Each sample is done individually and after the layers have settled the bottom layer is drawn off with suction.
- 10 A small amount of anhydrous  $\text{Na}_2\text{SO}_4$  is added and the vessel is agitated gently. The isobutanol layer must be clear after the sulfate has settled. If not repeat this step.
- 11 The sample is then decanted into cuvettes to be read on the Coleman Model 19 Photofluorometer.
- 12 The B filters are used and after the machine has warmed up and has been balanced the working quinine sulfate standard is used to set the machine at 60.

#### Preparation of Standard Curve

- 1 Dilute the stock thiamine solution to 1 mcg/ml in 0.1 N  $\text{H}_2\text{SO}_4$  to which has been added 10 ml 2.5 M sodium acetate/100 ml. Run a series of standards in the range 0.2 mcg per tube following all steps of the procedure.
- 2 Plot these values as described in the section on fluorometry.

#### Calculation

$$1 \text{ mcg thiamine/100 ml urine} = \frac{(\text{mcg in cuvette} - \text{blank}) \times (\text{total ml eluate}) \times 100}{(\text{ml eluate}) \times (\text{ml urine})}$$

- 2 The problem of the blank has never been solved completely in this procedure because urine contains many interfering substances. In particular N-methyl nicotinamide is adsorbed and eluted like thiamine fluoresces in alkaline solution and is only partially destroyed by ferricyanide. Najjar has recommended a blank consisting of urine treated in every way as for thiamine omitting the ferricyanide step and making an empirical correction. Mickelsen has recommended washing the final isobutyl layer with acid to remove non-thiochrome fluorescence. The reader will have to make his own choice of blank.

#### Example

The unknown samples read an average of 55 and the blank 5 equivalent to 0.68 mcg and 0.06 mcg of thiamine respectively in the cuvette. Five ml of urine had been used and eluted with 5 ml of KCl all of which was oxidized.

$$\text{mcg thiamine/100 ml urine} = \frac{(0.68 - 0.06) \times 5 \times 100}{5 \times 5} = 12.4 \text{ mcg/100 ml}$$

#### Precautions

- 1 The samples must be shaken immediately after the addition of the oxidizing solution.
- 2 Be sure the Decalco is activated before using it.
- 3 The isobutyl alcohol must have a low blank if it does not redistill it until it does read below 5 on the scale.

### B THIAMINE IN FOOD, BLOOD AND FECES

#### Reference

Same as for thiamine in urine.

#### Principle

Same as for thiamine in urine except that in food, blood and feces bound thiamine has to be hydrolyzed to liberate free thiamine.

#### Apparatus

- 1 Digestion tubes. The same tubes are used for digestion as for nonprotein nitrogen. They are sufficiently accurate for many types of samples to serve as receiving tubes for the eluate. Heavy walled tubes made from ignition

tubes should be specified. It is necessary that they be heavy walled since they are also used as centrifuge tubes in some analyses (blood). As purchased they are calibrated at 1, 5, 35 and 50 ml.

- 2 J shaped stirring rods. These are made from 1 inch lengths of 3 mm glass rod. The round bend at one end fits the bottom of the tube. With such a rod the contents of the tube can be violently agitated by rolling the rod back and forth between the thumb and index finger.
- 3 Hennessey adsorption tubes with the long (5 to 5.5 inch) stems are recommended.
- 4 Reaction vessels. Centrifuge tubes of about 50 ml capacity which are glass stoppered and pointed at the bottom (Wilkins Anderson #7105 GMR).
- 5 Waring Blender
- 6 Coleman Photofluorometer
- 7 Cuvettes 10 x 150 mm
- 8 An incubator set at 45 C

#### Reagents

These are the same as for thiamine in urine with a few additions.

- 1 Hydrochloric acid 0.1 N. Mix 1 ml of conc. HCl with 110 ml of water.
- Sodium bicarbonate. To prepare a molar solution dissolve 84 gm up to 1000 ml with water.
- 3 Metaphosphoric acid. A 10 per cent solution is prepared at frequent intervals preferably each day. Since the acid deteriorates rapidly at room temperature it should be dissolved in cold water and stored in a refrigerator.
- 4 Enzymes. A 10 per cent solution of Takadiastase or Clarase is prepared at frequent intervals preferably each day. The preparation should yield a pale yellow solution and it should dissolve readily with very little residue. The insoluble material contains considerable quantities of thiamine. The solution is centrifuged and the clear supernatant solution is kept in a refrigerator.

#### Collection and Preparation of Samples

- 1 Blood. Blood is collected and is prevented from clotting with a minimum of heparin or oxalate. The sample is cooled immediately and the analysis is made within a few hours of collection. Five ml are transferred to an NPN tube the volume is adjusted to about 5 ml and 0.0 ml of 0.1 N HCl are added.
- Feces. Samples are collected in large glass vessels and stored in the cold. For digestion 5 gm samples are weighed on filter papers and filter paper and samples are transferred to NPN tubes. Three ml of 0.1 N HCl are added and the volume is adjusted to 35 ml. (Because of many factors which affect digestion absorption and intestinal motility the composition of the stool may vary considerably from day to day even when the human subject is receiving a constant weighed diet. In order to be of any significance the collection should be carried out for at least 4 days. The storage bottles must be kept cold. Bacterial synthesis of thiamine due to growth of bacteria may greatly increase the vitamin content of the sample if it is not immediately analyzed or if it is not kept cold.)
- 3 Cereals and dry samples of food. The sample is finely ground by suitable means.
- 4 Fresh tissues animal and plant. The sample is rapidly cut into small strips or chunks. A weighed portion is homogenized with approximately 100 ml of water in the Waring Blender. It is then transferred quantitatively to a volumetric flask and the volume is brought to the mark. A suitable aliquot is transferred to an NPN tube.

#### Procedure -- Extraction

- 1 The finely ground or homogenized sample is weighed or measured by means of a syringe pipette. It is transferred to an NPN tube the volume is brought to the 35 ml mark and a measured volume of 0.1 N HCl is added. The quantity of acid to be added should be determined for each type of material. Enough acid should be



added so that the acidity after 10 or more minutes in the boiling water bath is approximately pH 3. The contents of the tube are thoroughly emulsified by means of a J shaped stirring rod. They are then heated 30 minutes with frequent stirring in the boiling water bath.

- 2 The tube is now cooled in a cold bath of running water after which 1 drop of caprylic alcohol and a volume of M NaHCO<sub>3</sub> which is equivalent to the N HCl is added. The contents after stirring should have a reaction of pH 4 to 5.
- 3 Several blank tubes which contain equal volumes of N HCl and M NaHCO<sub>3</sub> are prepared at this time. They should be carried through the subsequent operations along with the sample.

#### Procedure—Hydrolysis of Bound Thiamine

- 1 Add 2 ml of enzyme solution. More of enzyme preparation is added to samples which contain much starch or viscous material. The tube is incubated 2 hours in a water bath at 40 to 50 C. The mixture is stirred frequently. (Some samples like blood required only 1 hour of incubation. Long incubation may introduce serious errors due to bacterial growth. The temperature should be kept near 45 C in order to prevent such growth.)
- 2 The tube is cooled in a bath of running water.

#### Procedure—Precipitation of Protein

- 1 The contents of the tube are transferred quantitatively to a volumetric flask of 50 ml capacity. 20 ml of HPO are added and the volume is brought to the mark. Sufficient time should be given to allow diffusion of the vitamin from the particles into the solution.
- 2 Filter or centrifuge. (The analysis may be discontinued at this time. The flasks may stand overnight in the refrigerator.)
- 3 Blood and samples which contain less than 1 microgram of thiamine are treated somewhat differently. Ten ml (or less if the protein content is small) of HPO solution are added, the volume is brought to 50 ml and the contents are mixed thoroughly. The supernatant solution is poured off after centrifugation. The residue is emulsified with 25 ml HPO solution (using the J rod) and the volume is again brought to the 50 ml mark. The clear supernatant solution is added to the first supernatant.

#### Procedure—Adsorption, Elution, Oxidation, Extraction

- 1 The necessary adsorption columns are prepared as under thiamine in urine.
- 2 An aliquot of the extract sufficient to contain 1 to 5 mcg of thiamine is adjusted to pH 3.5 by means of 0.1 N NaOH with bromocresol green as indicator. The extract is then passed through the column. After passage of all of the extract the flask is rinsed with several successive portions of about 5 ml each of distilled water and the washings are sent through the column. The fluid is allowed to drain completely after each of the washings.
- 3 With the tip of the column in an NPN tube four 5 ml portions of 5% KCl are sent through the column. Finally enough 5% KCl is sent through to bring the volume to exactly 25 ml.
- 4 Five ml portions of eluate (10 for blood) are added to extraction vessels and the rest of the procedure is followed as under thiamine in urine.

#### Calculation

- 1 Prepare calibration curves with standard thiamine as described under urine.
- 2 mcg total thiamine/100 gm initial substance =

$$\frac{(\text{mcg in cuvette} - \text{blank}) \times (\text{total eluate}) \times (\text{total extract}) \times 100}{(\text{ml eluate used}) \times (\text{ml extract used}) \times (\text{gm sample})}$$

**Example**

A sample of feces weighed 4.701 gm. The fluorescence of the unknown was equivalent to 0.60 mcg thiamine in the tube that of the blank was equivalent to 0.07 mcg. The sample had been eluted to 50 ml of which 50 ml had been sent through the zeolite. The zeolite had been eluted with 25 ml of HCl of which 5 ml samples were analyzed.

$$\frac{\text{mcg total thiamine}/100 \text{ gm feces} = (0.60 - 0.07) \times (5) \times (250) \times 100}{(5) \times (50) \times (4.701)} = 8 \text{ mcg}/100 \text{ gm}$$

**Precautions**

Same as for thiamine in urine

**4 N<sup>1</sup> Methylnicotinamide in Urine****Reference**

Coulson R. A. Ellinger P. and Holden M. A Method for the Estimation of N<sup>1</sup> cotinamide Methochloride in Urine Biochem J 38 150-154 (1944)

**Principle**

N<sup>1</sup> methylnicotinamide a metabolic product of niacin breaks down in alkaline solution to a fluorescent compound which is soluble in isobutyl alcohol

**Apparatus**

Same as for thiamine in urine

**Reagents**

Same as for thiamine in urine with a few additions

- 1 Standard N<sup>1</sup> methylnicotinamide Dissolve 10 mg of the pure substance in 50 ml of 1% acetic acid (40 mcg/ml of solution)

**Procedure**

Same as for thiamine in urine down to the oxidation step. After eluting the column with HCl proceed as follows

- 1 To 5 ml of HCl eluate add 10 ml of isobutyl alcohol. Add 0.7 ml of 15% NaOH stopper and shake immediately for at least 50 times. Stand at least 10 minutes before reading.
- 2 Proceed as in thiamine analysis with quinine standard at 60 (maximum fluorescence develops slowly so that a period of standing is required before reading)

**Standard Curve**

Run a series of standard N<sup>1</sup> methylnicotinamide solutions. The range should be 40 gamma to 100 gamma per sample with intervals of 40 gamma.

**Calculation**

$$\text{mg N}^1 \text{ methylnicotinamide}/100 \text{ ml urine} = \frac{(\text{mcg in cuvette}) \times (\text{total eluate}) \times 100}{(1000) \times (\text{ml aliquot}) \times (\text{ml urine})}$$

**Example**

A reading of 30 on the gal. anometer was equivalent to 43 mcg in the cuvette. The blank was equivalent to 10 mcg. Five ml of urine were used for analysis and eluted with 5 ml of HCl all of which was used for analysis.

$$\text{mg N}^1 \text{ methylnicotinamide}/100 \text{ ml urine} = \frac{(43 - 10) \times 5 \times 100}{1000 \times 5 \times 5} = 0.66 \text{ mg}/100 \text{ ml}$$

**Precautions**

This is a method fraught with inaccuracies. At every step there are significant losses. Hence the results must be considered only semi-quantitative.

## 5 Riboflavin in Urine, Blood, Food and Feces

## A. IN URINE

## Reference

Connor, R T and Straub G J Combined Determination of Riboflavin and Thiamin in Food Products *Ind and Eng Chem Anal Ed* 13 385 389 (June) 1941

## Principle

Riboflavin is measured fluorometrically after interfering substances are destroyed. An internal standard is used and the blank is determined after reduction of riboflavin to the leuco form which is not fluorescent.

## Apparatus

- 1 Photofluorometer Coleman Model 12 or 12B
- 2 Cuvettes Coleman 19 x 150 mm
- 3 Syringe pipettes one and 2 ml
- 4 A number of 75 ml graduated cylinders with ground glass stoppers

## Reagents

- 1 Buffer solution pH 4.7 Dissolve 111 gm of sodium acetate plus 54.4 ml of glacial acetic acid in 1000 ml of water
- 2 Potassium permanganate a 4% solution in water
- 3 Hydrogen peroxide a 3% solution in water
- 4 Standard riboflavin solutions Prepare a stock solution (100 mcg/ml) by dissolving 10 mg of standard in exactly 100 ml of 3% acetic acid. Store in a dark bottle in a refrigerator. This solution is stable for six months. For a working standard (1 mcg/ml) dilute the stock 1 to 100 with water.

TABLE 10

## RIBOFLAVIN IN URINE

(For analyses in which 1 ml of urine is analyzed in 10 ml of solution. A is galvanometer reading with unknown. B with 1 mcg riboflavin internal standard and C with blank. This table gives approximately correct values in the range of normal urine.)

| A C (ACROSS) | 18             | 0  | 22 | 24  | 26  | 8    | 30   | 32   | 34   | 36   |
|--------------|----------------|----|----|-----|-----|------|------|------|------|------|
| B A (DOWN)   | mcg B/ml URINE |    |    |     |     |      |      |      |      |      |
| 20           | 82             | 89 | 98 | 107 | 115 | 1.3  | 1.9  | 1.35 | 1.45 | 1.54 |
| 1            | 78             | 85 | 94 | 10  | 110 | 1.18 | 1.24 | 1.32 | 1.39 | 1.46 |
| 22           | 75             | 82 | 90 | 98  | 105 | 1.13 | 1.19 | 1.27 | 1.33 | 1.41 |
| 23           | 1              | 79 | 87 | 94  | 101 | 1.09 | 1.13 | 1.22 | 1.28 | 1.35 |
| 4            | 69             | 75 | 84 | 90  | 98  | 1.04 | 1.10 | 1.17 | 1.24 | 1.31 |
| 5            | 66             | 72 | 81 | 88  | 94  | 1.00 | 1.07 | 1.13 | 1.20 | 1.26 |
| 26           | 64             | 70 | 8  | 84  | 90  | 97   | 1.03 | 1.09 | 1.15 | 1.22 |
| 27           | 6              | 68 | 75 | 81  | 88  | 94   | 1.00 | 1.06 | 1.12 | 1.19 |
| 8            | 60             | 65 | 73 | 79  | 85  | 90   | 97   | 1.02 | 1.08 | 1.14 |
| 29           | 58             | 63 | 70 | 76  | 82  | 88   | 93   | 99   | 1.05 | 1.10 |
| 30           | 56             | 61 | 68 | 74  | 79  | 85   | 90   | 96   | 1.0  | 1.07 |
| 31           | 51             | 56 | 66 | 72  | 77  | 82   | 88   | 93   | 99   | 1.04 |
| 32           | 52             | 58 | 65 | 70  | 75  | 79   | 85   | 91   | 96   | 1.01 |
| 33           | 51             | 56 | 63 | 68  | 73  |      | 82   | 88   | 93   | 99   |
| 34           | 50             | 55 | 61 | 66  | 71  | 75   | 80   | 85   | 91   | 96   |
| 35           | 48             | 53 | 59 | 64  | 69  | 73   | 78   | 82   | 88   | 93   |
| 36           | 47             | 52 | 57 | 62  | 67  | 71   | 75   | 80   | 86   | 91   |
| 37           | 46             | 51 | 56 | 61  | 65  | 69   | 73   | 78   | 84   | 88   |
| 38           | 45             | 49 | 54 | 59  | 64  | 67   | 71   | 75   | 8    | 86   |
| 39           | 44             | 48 | 53 | 58  | 63  | 65   | 69   | 74   | 80   | 84   |
| 40           | 43             | 47 | 52 | 56  | 61  | 64   | 67   | 7    | 77   | 8    |
| 41           | 42             | 46 | 51 | 55  | 59  | 63   | 66   | 70   | 75   | 80   |

- 5 Sodium fluorescein reference standard Dilute 10 mg of sodium fluorescein to 1000 ml in 3% acetic acid For a working standard (0.5 mcg/ml) dilute the stock 1 to 20 with water
- 6 Pure sodium hydrosulfite powdered

#### Procedure

- 1 Pipette exactly 1 ml urine plus 3 ml buffer solution into a 95 ml graduated stoppered cylinder
- 2 Add one or two drops of 4%  $\text{KMnO}_4$  until a pink or violet color persists for one minute
- 3 Then add two to four drops of  $\text{H}_2\text{O}$  (3%) until the color disappears
- 4 Make up to 10 ml with buffer solution mix and transfer to cuvettes
- 5 For a reagent blank use 1 ml distilled water instead of 1 ml of urine
- 6 Read samples within one hour (keeping in dark meanwhile) as follows
  - a Warm up instrument for 15 minutes before using Use B filter
  - b Balance at zero and set galvanometer at 25 with a sample of the dilute fluorescein reference in the cuvette
  - c Read the unknown (Reading A)
  - d To the sample add exactly 1 ml of standard B solution (10 mcg B)
  - e Mix thoroughly with stirring rod and read immediately (Reading B)
  - f To the same cuvette now add a small amount of sodium hydrosulfite on the end of a spatula at ring until dissolved
  - g Read immediately (Reading C)

#### Calculation

- 1 Reading A represents the galvanometer deflection due to riboflavin plus blank in 10 ml
- Reading B represents the galvanometer deflection due to added riboflavin in 11 ml
- Reading C represents the galvanometer reading due to the blank alone in 11 ml

TABLE 10 (CONT'D)

| A C (ACROSS) | 35                           | 40  | 4   | 44  | 46  | 48  | 50  | 5   | 54  | 56  |
|--------------|------------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| B A (DOWN)   | mcg B <sub>2</sub> /ml URINE |     |     |     |     |     |     |     |     |     |
| 20           | 139                          | 16  | 174 | 180 | 188 | 194 | 00  | 07  | 14  | 20  |
| 21           | 153                          | 160 | 167 | 173 | 180 | 189 | 19  | 00  | 0   | 21  |
| 3            | 147                          | 154 | 160 | 167 | 174 | 181 | 185 | 19  | 198 | 204 |
| 4            | 141                          | 148 | 154 | 161 | 167 | 174 | 178 | 185 | 191 | 19  |
| 5            | 136                          | 143 | 149 | 155 | 161 | 168 | 174 | 178 | 184 | 190 |
| 6            | 131                          | 138 | 143 | 150 | 156 | 16  | 166 | 173 | 178 | 184 |
| 7            | 1                            | 123 | 139 | 14  | 151 | 157 | 161 | 167 | 173 | 178 |
| 8            | 123                          | 129 | 135 | 140 | 146 | 15  | 156 | 16  | 167 | 173 |
| 9            | 119                          | 125 | 130 | 136 | 141 | 14  | 151 | 15  | 167 | 168 |
| 10           | 115                          | 121 | 126 | 132 | 137 | 143 | 147 | 15  | 158 | 163 |
| 11           | 11                           | 116 | 121 | 126 | 133 | 137 | 143 | 148 | 153 | 158 |
| 12           | 109                          | 114 | 119 | 124 | 130 | 131 | 135 | 144 | 149 | 154 |
| 13           | 106                          | 111 | 116 | 121 | 126 | 131 | 135 | 140 | 145 | 150 |
| 14           | 103                          | 108 | 113 | 118 | 123 | 128 | 131 | 136 | 141 | 146 |
| 15           | 100                          | 105 | 110 | 115 | 119 | 124 | 128 | 133 | 137 | 142 |
| 16           | 9                            | 103 | 107 | 111 | 116 | 12  | 125 | 130 | 134 | 138 |
| 17           | 95                           | 100 | 104 | 109 | 113 | 119 | 124 | 127 | 131 | 135 |
| 18           | 93                           | 98  | 101 | 106 | 111 | 116 | 119 | 124 | 128 | 133 |
| 19           | 91                           | 95  | 99  | 103 | 108 | 113 | 116 | 121 | 125 | 129 |
| 20           | 88                           | 93  | 97  | 101 | 106 | 110 | 113 | 118 | 122 | 126 |
| 21           | 86                           | 91  | 95  | 99  | 103 | 108 | 110 | 115 | 119 | 123 |
| 22           | 84                           | 89  | 93  | 97  | 101 | 106 | 108 | 113 | 117 | 120 |

2 mcg riboflavin/100 ml urine =

$$\frac{(A - C \times \frac{11}{10})}{(B \times \frac{11}{10} - A)} \times \text{mcg riboflavin added} \times \frac{\text{ml dilution}}{\text{ml urine}} \times 100$$

#### Example

- 1 One ml of urine was diluted to 10 and read 15.0 (Reading A). After 1 mcg of riboflavin was added the reading was 50.0 (Reading B). After hydrosulfite the reading was 5.0 (Reading C).

mcg riboflavin/100 ml urine =

$$\frac{15.0 - 5 \times \frac{11}{10}}{50.0 \times \frac{11}{10} - 15.0} \times 1 \times \frac{10}{1} \times 100 = 275 \text{ mcg/100 ml}$$

- 2 A table (Table 10) has been constructed which for most practical purposes gives a good approximation to the correct value over most of the normal range. The assumptions made are: (a) the blank is usually small so that  $\text{blank} \times \frac{11}{10}$  is equivalent to blank; (b) Reading B - Reading A is almost constant so that  $B \times \frac{11}{10} - A$  is also almost constant.

#### Precautions

- 1 Keep samples in dark as much as possible.
- 2 Because of the danger of reoxidation, read immediately after adding hydrosulfite.

### B IN BLOOD, FOOD AND FECES

#### Reference

Same as for riboflavin in urine.

#### Principle

The determination consists of

- (1) liberation of bound alloxazine nucleotides by means of heat and acid
- (2) clarification of the sample and conversion of alloxazine nucleotides into riboflavin by means of Clarase or Takadiastase
- (3) removal of many impurities by adsorption on Supersorb and subsequent elution with pyridine solution at pH 5
- (4) oxidation of impurities in the eluate by means of  $\text{KMnO}_4$  and
- (5) determination of fluorescence

#### Apparatus

Same as for thiamine in blood, food and feces.

#### Reagents

- 1 Hydrochloric acid approx N. Conc. HCl (sp. gr. 1.19) is diluted 1 L. times, or 83 ml are diluted to a volume of 1000 ml in distilled water.
- 2 Sodium bicarbonate approx N. Dilute 84 gm  $\text{NaHCO}_3$  to 1000 ml.
- 3 Isoamyl alcohol.
- 4 Enzyme solution. A 10% solution of Takadiastase (Parke Davis, Detroit) or Clarase (Takamine Laboratories, Clifton, N. J.) is prepared at frequent intervals, preferably each day. The preparation should yield a pale yellow solution and it

should dissolve readily with very little residue. The solution is centrifuged and the clear supernatant solution is kept in the refrigerator when not in use.

- 5 Supersorb (Florida Co., Warren, Pa.) The adsorbent as bought contains much fine material which must be removed. It also must be activated before it can be used. One pound or more is covered with a 1 per cent solution of acetic acid (in a large round bottom flask with short neck or 3 liter Erlenmeyer flask). The mixture is stirred very briefly by rotation allowed to stand until most of the adsorbent has settled and the supernatant solution is then poured off. This is repeated twice. It is similarly washed by decantation several times with distilled water after which it is transferred to a large Buchner funnel and washed further with distilled water. The Supersorb is then dried at room temperature in a large pan. In order to have the material in its maximum expanded state before use a small quantity is kept moistened in a wide mouthed bottle.
- 6 Pyridine solution. Approximately 500 ml of distilled water, 200 ml of redistilled pyridine and 0 ml of glacial acetic acid are transferred to a 1000 ml volumetric flask after which the volume is made to the mark with distilled water.
- 7 Potassium permanganate. Dissolve 31.6 gm  $KMnO_4$  in one l distilled water. Keep this approximately N solution in a brown dropping bottle.
- 8 Hydrogen peroxide USP 3%. Keep in a brown dropping bottle.
- 9 Riboflavin stock solution. The reference standard (Dr. E. F. Cook, 43rd St and Woodland Ave Philadelphia) is dried for 4 hours over phosphorus pentoxide. 40 mg are transferred to a flask of 1000 ml capacity. Water to about 800 ml and 100 ml of approximately N  $H_2SO_4$  are added. The solution is frequently shaken and is kept in the refrigerator or in a cool dark place until all of the vitamin is dissolved. The contents are warmed to a temperature of  $0^\circ C$  after which the volume is brought to the mark. The flask is protected from light in double brown paper bags. If kept in the refrigerator the stock standard solution is stable for at least 6 months (40 mcg/ml).
- 10 Riboflavin working solution. Dilute 2 ml of the stock solution to a volume of 50 ml or 10 ml to a volume of 250 ml. Five ml of this intermediate standard and 2.5 ml (graduated cylinder) of pyridine solution are transferred to a 50 ml volumetric flask and the volume is brought to the mark. The sensitivity is adjusted to a galvanometer reading of 50 (0.16 mcg/ml).
- 11 Sodium hydrosulfite solution (Lycopon). Since sodium hydrosulfite reacts rapidly with atmospheric oxygen solutions must be prepared immediately before use. They must not be unnecessarily shaken. The surfaces should be small and solutions should be used not longer than 2 hours. To 5 ml of ice cold water contained in an NPN tube are added 1 gm of sodium hydrosulfite and 1 gm of  $NaHCO_3$ .
- 12 Caprylic alcohol.

#### Collection and Preparation of Samples

The conditions for collection and preparation of samples are the same as described in the directions for the determination of thiamine by the thiochrome method.

Special precautions must be taken to avoid unnecessary exposure of samples to strong light especially to sunlight.

#### Procedure

- 1 The finely ground or homogenized sample is weighed or measured by means of a pipette. It is transferred to an NPN tube the volume is brought to the 35 ml mark and a measured volume of N  $HCl$  (from 3 to 5 ml) is added. Enough acid should be added to bring the acidity of the boiled sample to pH 1 to 2. The quantity of acid to be added should be determined for each type of material.

- 2 The contents of the tube are thoroughly emulsified by means of a J shaped stirring rod. They are then heated 30 minutes to 1 hour with frequent stirring especially during the first 10 minutes in the boiling water bath.
- 3 The tube is now cooled in a cold bath of running water, after which 1 drop of caprylic alcohol and a volume of M NaHCO<sub>3</sub>, which is equivalent to the N HCl are added. The contents after stirring, should have a reaction of pH 4 to 5.
- 4 Several blank tubes which contain equal volumes of N HCl and M NaHCO<sub>3</sub> are prepared at this time. They should be carried through the subsequent operations along with the sample.
- 5 Add 0.0 ml of Clarase or Takadiastase solution. More of the enzyme solution is added to samples which contain much starch or viscous material. The tube is incubated 2 hours in a water bath at 40 to 50 C the contents being stirred frequently.
- 6 The tube is cooled in a bath of running water after which the contents are transferred quantitatively to a volumetric flask of 250 ml capacity. The volume is brought to the mark. At this point the analysis may be discontinued and the flask may be stored in the refrigerator.
- 7 The solution is allowed to settle overnight or if analyzed immediately it is centrifuged. Filtration may cause large and variable errors due to adsorption of the vitamin by the paper. If filtered the first 25 ml of the solution should be discarded.
- 8 The columns are filled with Supersorb as described in directions for the determination of thiamine.
- 9 A suitable aliquot of the clear supernatant solution is transferred to a 150 ml extraction flask. Bromocresol green indicator and enough water to bring the volume to about 100 ml are added. The reaction is then adjusted to approximately pH 4.5 by means of 0.1 N HCl or NaOH.
- 10 The extract is passed through the column. After passage of all the extract the flask is rinsed with two successive portions of about 5 ml each of distilled water and the washings are sent through the column. The fluid is allowed to drain completely after each of the washings. The container below the column is removed and the liquid which has passed through the column is discarded. A calibrated NPN tube is put in its place below the column.
- 11 Small portions of pyridine solution are now added until exactly 25 ml of eluate have been collected at which point the elution is discontinued.
- 12 Oxidation of impurities in the pyridine eluate. Permanganate solution is added drop by drop until an excess is present as indicated by a violet color. After an interval of 1 to 5 minutes H<sub>2</sub>O is added drop by drop until the solution is completely clear. Usually 4 drops of KMnO<sub>4</sub> and 8 drops of H<sub>2</sub>O solutions are sufficient.
- 13 The volume is brought to the 50 ml mark and the solution is mixed by means of a J stirring rod.
- 14 Several tubes or cuvettes containing the standard should be prepared and they should be used alternately to prevent rise of temperature. It should be emphasized that the fluorescence changes with temperature. The standards should therefore be kept under the same conditions as the solutions which are to be examined. The standard tubes should be stoppered if they are to be used over a period of time.
- 15 The cuvette or tube is filled to the mark. The sensitivity of the apparatus is set by means of the standard after which the reading is taken with the unknown solution. This is reading A. An excess (0.2 ml) of sodium hydrosulfite solution is then added the contents are quickly mixed and the reading is made immediately. This is reading B. (The solution should be mixed with gentle rotation after addition of the hydrosulfite. Violent mixing produces bubbles which scatter the light and thus give high results. Riboflavin becomes nonfluorescent in

stantly. Biological materials contain many substances whose fluorescence is slowly destroyed by the hydrosulfite. The reading therefore should be made immediately after addition and mixing of the reagent.)

#### Standardization

1. Known quantities of the standard are adsorbed on Supersorb, eluted with pyridine and treated with  $\text{KMnO}_4$  as described above. Water blanks should be included in the series. The readings corrected by the blank are then plotted as mcg riboflavin/ml solution in the cuvette.

The presence of light adsorbing substances in the eluate may markedly alter the results even though the recovery of added riboflavin may be quantitative. Milk, urine and feces for example contain such substances. However their effect is almost entirely eliminated by the adsorption and the  $\text{KMnO}_4$  oxidation. The recovery should be determined for every type of sample. If repeated analyses yield quantitative recovery of added riboflavin it is not necessary to make such determinations in subsequent analyses. It is suggested however that the recovery be tested at frequent intervals as a check on the accuracy of the technique.

#### Calculation

1. mcg riboflavin/ml solution in cuvette =

$$\frac{(\text{unknown reading } A - B) - (\text{Blank reading } A - B)}{\% \text{ recovery}} \times 100 \times \text{calibration factor}$$

2. mcg riboflavin/100 ml original material =

$$\frac{(\text{mcg riboflavin/ml in cuvette}) \times (\text{ml eluate}) \times (\text{total ml extract}) \times 100}{(\text{ml extract through column}) \times (\text{gm sample})}$$

#### Example

Weight of food 4 gm Dilution after hydrolysis 50 ml Total volume adsorbed 50 ml Total volume elution 50 ml Reading A 55 Reading B 11 Blank reading A 9 Blank reading B 5 Recovery 94%

1. mcg riboflavin/ml solution in cuvette =

$$\frac{(55 - 11) - (9 - 5)}{94} \times 100 \times \frac{0.16}{50} = \frac{40}{94} \times 0.3 = 0.136 \text{ mcg/ml}$$

$$\text{mcg riboflavin/100 gm food} = \frac{0.136 \times 95 \times 50 \times 100}{50 \times 4} = 404 \text{ mcg/100 gm}$$

#### Precautions

Same as for riboflavin in urine

## 6 Ascorbic Acid in Whole Blood, Serum and Urine

### A. WHOLE BLOOD AND SERUM (MACROMETHOD)

#### Reference

Roe J. H. and Wether C. A. The Determination of Ascorbic Acid in Whole Blood and Urine Through the 4-dinitrophenylhydrazine Derivative of Dehydroascorbic Acid. *J. Biol. Chem.* 147: 393-401 (Feb.) 1943

#### Principle

Ascorbic acid is converted to dehydroascorbic acid. In the presence of 4-dinitrophenylhydrazine dehydroascorbic acid is coupled with the phenylhydrazine. In the presence of strong sulfuric acid a colored compound is formed from the hydrazone.



**Apparatus**

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 × 150 mm
- 3 Centrifuge tubes capacity of 50 ml
- 4 Incubator set at 37 C
- 5 A quantity of 125 ml Erlenmeyer flasks
- 6 An electric centrifuge
- 7 A quantity of 50 ml centrifuge tubes
- 8 Quantitative filter paper Whatman 40
- 9 An ice bath

**Reagents**

- 1 2,4-dinitrophenylhydrazine reagent Dissolve 2 gm of hydrazine in 100 ml of 9 N sulfuric acid (3 parts of water and one part of concentrated sulfuric acid) Filter
- 2 Acid washed Norit Place 200 gm Norit in a large flask Add 1 liter of 10% HCl and heat to boiling Filter with suction Remove Norit cake to a large beaker and add 1 liter distilled water Stir thoroughly and filter Repeat until washings give a negative or very faint test for ferric iron. (Test filtrate with 1% potassium ferrocyanide)  
Dry Norit cake in oven overnight at 110–120 C (Norit is an activated char coal)
- 3 Trichloroacetic acid a 6% solution in water
- 4 Sulfuric acid 85% solution To 100 ml of distilled water add 900 ml of concentrated sulfuric acid (sp gr 1.84)
- 5 Thiourea a 10% solution Dissolve 10 gm thiourea in 100 ml of 50% (by volume) aqueous ethyl alcohol This reagent keeps satisfactorily for two months but to check it see that it readily reduces HgCl<sub>2</sub> or KMnO<sub>4</sub>
- 6 Standard Vitamin C Dissolve exactly 100 mg in 100 ml of 5% acetic acid (1 mg/ml) Dilute 10 to 100 with 6% trichloroacetic acid (100 mcg/ml) Add 1 ml of thiourea before making to volume

**Procedure**

- 1 Pipette 15 ml of 6% trichloroacetic acid into a 50 ml Pyrex centrifuge tube Add 0.5 ml of whole blood (plasma) drop by drop and with a glass rod stir until fine suspension is produced
- 2 Let stand for at least 5 minutes and centrifuge
- 3 Transfer supernatant to Pyrex test tube Add 0.75 gm (3 scoops) acid washed Norit and shake vigorously Filter through 9 cm filter paper
- 4 Place 4 ml Norit filtrate in each of two matched photoelectric colorimeter tubes Add 1 drop of 10% thiourea to each
- 5 Hold one tube as reserve for blank
- 6 Add 1 ml hydrazine reagent to other
- 7 Place second tube in constant temperature bath at 37 C Keep tube immersed for exactly 3 hours
- 8 Remove tube and place this and the blank in an ice bath
- 9 To each tube add slowly 5 ml of 85% sulfuric acid Add the acid drop by drop from the burette taking at least one to one and a half minutes for total addition Mix thoroughly with tubes immersed
- 10 Add 1 ml hydrazine reagent to blank.
- 11 Remove and after 30 minutes at room temperature wipe the tubes clean and dry
- 12 Transfer to cuvettes
- 13 Set blank at 100% T at 540 mμ and read unknowns

**Standard Curve**

Prepare a standard curve by running a series of ascorbic acid standards through the method The range should be from 0 to 35 mcg per tube at intervals of 5 mcg Plot on semilog paper

**Calculation**

$$\text{mg total ascorbic acid/100 ml serum} = \frac{\text{mcg in cuvette} \times \text{total dilution} \times 100}{1000 \times \text{ml aliquot} \times \text{ml serum}}$$

**Example**

A reading of 77.5% T was equivalent to 10 mcg of ascorbic acid in the cuvette. Five ml of blood were diluted to 10 ml of which 4 ml were taken for analysis.

$$\text{mg total ascorbic acid/100 ml serum} = \frac{10 \times 20 \times 100}{1000 \times 4 \times 5} = 10 \text{ mg/100 ml}$$

**Precautions**

- 1 Norit must be free from ferric iron whose ions may irreversibly oxidize ascorbic acid. Ferric ions also react with hydrazine reagent producing interfering colors. Thiourea by producing a reducing environment protects against interference by ferric ions and other oxidizing substances. The coupling reaction is not appreciably slowed by thiourea.
- 3 The final color is stable for 40 minutes. After 18 hours the maximum fading is 0.5 points on the galvanometer scale.
- 4 Glucose, xylose, fructose and glucuronic acid do not interfere in ordinary physiologic and pathologic ranges.

**B. IN URINE (MACROMETHOD)****Reference Principle Apparatus and Reagents**

Same as for ascorbic acid in blood.

**Procedure**

- 1 To a 15 ml Erlenmeyer flask add exactly 5 ml of urine.
- 2 Add 5 ml of 6*N* trichloroacetic acid. (This stabilizes the ascorbic acid during extraction.)
- 3 Add 0.5 gm acid washed Norit and shake vigorously. (Norit oxidizes ascorbic acid to dehydro form.)
- 4 Filter.
- 5 Pipette 4 ml of filtrate into two colorimeter tubes.
- 6 Add 1 drop of thiourea to both.
- 7 Add 1 ml of 4-dinitrophenylhydrazine to one, excluding the blank.
- 8 Mix and incubate for 3 hours at 37°C.
- 9 Remove and place in ice bath.
- 10 Add 5 ml of 85% sulfuric acid slowly so as not to raise the temperature (1 to 1½ minutes).
- 11 Add 1 ml of hydrazine reagent to the blank.
- 12 Set blank at 100% T and read all specimens at 540 mμ.

**Standard Curve**

Same as for ascorbic acid in blood.

**Calculation**

$$\text{mg ascorbic acid/100 ml urine} = \frac{\text{mcg in cuvette} \times \text{total ml dilution} \times 100}{1000 \times \text{ml aliquot} \times \text{ml urine}}$$

**Example**

A reading of 7.5% T was equivalent to 10 mcg ascorbic acid in the cuvette. Five ml of urine had been diluted to 10 ml of which 4 ml aliquots were analyzed.

$$\text{mg total ascorbic acid/100 ml urine} = \frac{10 \times 10 \times 100}{1000 \times 4 \times 5} = 0.5 \text{ mg/100 ml}$$

**Precautions**

Same as for ascorbic acid in blood.

# WHOLE BLOOD AND SERUM (MICROMETHOD)

## Reference

- Bessey O A Vitamin Methods Vol 1, Academic Press Inc Publishers Edited by P Gyorgy New York, pages 303-304, 1950
- Lowry O H, Lopez J A and Bessey O A The Determination of Ascorbic Acid in Small Amounts of Blood Serum, J Biol Chem 160 609-615 (Oct) 1945

## Principle

Same as for the method of Roe and Kuether

## Apparatus

- 1 A Coleman Jr Spectrophotometer or a Beckman Spectrophotometer Model DU
- 2 A quantity of cuvettes 12 x 75 mm
- 3 1 cuvette holder for 12 x 75 mm cuvettes
- 4 A water bath or an incubator set at 38 C
- 5 An ice bath
- 6 Quantitative filter paper or an electric centrifuge
- 7 A quantity of 15 ml round bottom centrifuge tubes
- 8 A quantity of #1 solid rubber stoppers

## Reagents

- 1 A 2% solution of 2,4-dinitrophenylhydrazine made up in 10 N sulfuric acid.
- 2 Thiourea a 5% solution in water This solution keeps well for one month if left at 4 C
- 3 Copper sulfate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ) an 0.6% solution made up in 0.01 N sulfuric acid This solution is stable indefinitely
- 4 Sulfuric acid 65% solution Dilute 70 ml of concentrated sulfuric acid with 30 ml of distilled water Keep in the refrigerator
- 5 Trichloroacetic acid 4% solution
- 6 Trichloroacetic acid a 5% aqueous solution
- 7 For the combined reagent dilute 18 ml of the 2,4-dinitrophenylhydrazine solution 1 ml of copper sulfate solution and 1 ml of thiourea (This solution is stable for at least one week and should be centrifuged unless crystal clear)

## Procedure

- 1 Pipette 4 ml of 5% trichloroacetic acid into a 15 ml round bottom centrifuge tube
- 2 Add exactly 1 ml of whole blood or serum stopper and mix
- 3 Centrifuge at 3000 rpm for 10 minutes
- 4 Pipette exactly 1 ml of the aliquot into each of two small test tubes
- 5 Into one tube add exactly 0.25 ml of the mixed thiourea copper sulfate and 2,4-dinitrophenylhydrazine reagent
- 6 Mix by tapping and incubate for 4 hours at 38 C
- 7 Chill in ice water and add 1.5 ml of cold 65% sulfuric acid solution The cold acid may be added at once
- 8 Now add 0.25 ml of the mixed reagent to the blank tube and mix by shaking
- 9 Set the blank at 100% T at 590 m $\mu$  and read the unknowns
- 10 Treat 1 ml of standard ascorbic acid solution (in 4% trichloroacetic acid solution) in the same way These solutions should cover a range of 3 to 21 mcg/ml at intervals of 3 mcg

## Calculation

$$\text{mg total ascorbic acid/100 ml of serum} = \frac{\text{mcg in cuvette} \times \text{total dilution} \times 100}{1000 \times \text{ml aliquot} \times \text{ml serum}}$$

**Example**

One ml of blood was diluted to 5 ml of which one ml was used in the analyses. Ascorbic acid in the cuvette was 3.0 mcg.

Therefore

$$\text{mg total ascorbic acid/100 ml of blood} = \frac{3 \times 5 \times 100}{1000 \times 1 \times 1} = 1.5$$

**Precautions**

- 1 The final concentration of acid is critical in this method and must be watched carefully. Variations in the concentration of thiourea will cause large variations in the final color.
- 3 The age of the trichloroacetic acid appears to be immaterial.

**7 Ascorbic Acid in Tissue****Reference**

Roe J H, Mills M B, Oesterling M J and Damron C M. The Determination of Diketogulononic Acid, Dehydroascorbic Acid and Ascorbic acid in the Same Tissue Extract by the 2,4-Dinitrophenylhydrazine Method. *J Biol. Chem.* 174: 201-208 (May) 1948.

**Principle**

Same as for serum ascorbic acid.

**Apparatus**

The same as for serum ascorbic acid with the addition of

- 1 A mortar and pestle
- An H<sub>2</sub>S generator
- Carbon dioxide gas

**Reagents**

The same as for serum ascorbic acid with the addition of

- 1 Metaphosphoric acid 5% in water
- Stannous chloride solution. Dissolve 10 gm stannous chloride in 100 ml 5% metaphosphoric acid.
- 3 Thiourea, powdered
- 4 Bromine

**Preparation of Tissue**

- 1 Weigh accurately a sample of tissue of from 1 to 5 gm. This sample is ground in a mortar with 1 ml of 10% SnCl<sub>2</sub> solution. Add 100 ml of 5% metaphosphoric acid and mix thoroughly. (Do not use a Waring Blendor or homogenizing apparatus that would introduce increased amounts of oxygen into the slurry). The slurry is thoroughly mixed and filtered. (A solution with a dilution of 1 part of tissue in 100 parts of extract is desirable and in dilutions greater than this a concentration 1 to 10 per ml of total 4A, DHA and DKA is recommended. It is necessary to work rapidly to prevent changes in partition of the three compounds. A conversion of DHA to DKA of 7 percent in a solution containing 0.5 percent SnCl<sub>2</sub> - 5 percent HPO was observed in 1 hour at room temperature.

**Procedure for Diketogulononic Acid (DGA)**

- 1 Place 100 ml of the original SnCl<sub>2</sub>/HPO filtrate in a large tube. Introduce a gas filter tube with a sintered glass filter about 0.5 mm in diameter and pass H<sub>2</sub>S through for 15 minutes.

- 2 To 40 ml of the  $\text{H}_2\text{S}$  saturated extract add 0.4 gm of powdered thiourea shake until dissolved and filter
- 3 Bubble  $\text{CO}_2$  into the filtrate for 5 minutes
- 4 Pipette 4 ml aliquots into each of three colorimeter tubes
- 5 To two tubes add 1 ml of 2,4-dinitrophenylhydrazine The third tube will be a blank
- 6 Place in water bath at 37° C for 6 hours
- 7 From here on the analysis is the same as for serum ascorbic acid.

**Procedure for the Sum of Dehydro-1 Ascorbic Acid (DHA) and Diketo-1-Gulonic Acid (DGA)**

- 1 Pipette 4 ml aliquots of the original  $\text{SnCl}_2/\text{HPO}_4$  filtrate into each of three cuvettes
- 2 Add 1 ml of 2,4-dinitrophenylhydrazine to two of the tubes the third being used as a blank
- 3 Place in a water bath at 37° C for 6 hours
- 4 Continue as for serum ascorbic acid.

**Procedure for the Sum of 1 Ascorbic Acid (AA) Dehydro 1 Ascorbic Acid (DHA) and Diketo-1-Gulonic Acid (DGA)**

- 1 The remainder of the  $\text{H}_2\text{S}$  saturated solution is filtered to remove  $\text{SnS}$  (thiourea is not added) into a test tube suitable for bubbling gas Air is drawn through a water trap before passing through the  $\text{H}_2\text{S}$  solution to minimize volume changes Sufficient bromine is added to color the solution Excess bromine is removed by bubbling air through the solution
- 2 Add powdered thiourea in sufficient quantities to make a 1% solution
- 3 Pipette 4 ml aliquots into three colorimeter tubes
- 4 Add 1 ml of 2,4-dinitrophenylhydrazine to two tubes the third being used as a blank.
- 5 Place in a water bath at 37° C for 6 hours
- 6 Continue as for serum ascorbic acid

**Preparation of Calibration Curves**

- 1 Because the rate of coupling is influenced by the acid solution used and the antioxidant present it is important in preparing a calibration curve to use the same concentration of  $\text{HPO}_4$  and thiourea or  $\text{SnCl}_2$  as is used in the procedure followed for the unknown A calibration curve is made with standard dehydroascorbic acid solution prepared by bromine oxidation of ascorbic acid as directed by Roe and Oettinger The standards are made up in 5 percent  $\text{HPO}_4$  solution containing 1 percent of thiourea The tubes are incubated for 6 hours at 37° C The curve obtained with 1 percent thiourea is practically the same as that obtained with 0.5 percent  $\text{SnCl}_2$  Therefore calculations of DHA, DGA and AA are all made with the same standard curve

**Calculation**

- 1 All three substances are of practically the same molecular weight and are expressed in terms of mg ascorbic acid, according to the procedure outlined under ascorbic acid.
- 2 DGA is measured directly
- 3  $\text{DHA} = (\text{DHA} + \text{DGA}) - \text{DGA}$
- 4  $\text{AA} = (\text{AA} + \text{DHA} + \text{DGA}) - (\text{DHA} + \text{DGA})$

**Precautions**

- 1 Speed in handling the use and timing of reactions are the most important aspects of this method.

## 8 Rapid Field Methods for Ascorbic Acid, Riboflavin, Thiamine and N<sup>1</sup> Methylnicotinamide

### Reference

- 1 For ascorbic acid see Farmer C J and Abt A F Proc Soc Exper Biol & Med 34 146 150 1936
- 2 For riboflavin thiamine and N<sup>1</sup> methylnicotinamide see Johnson R E Sargent F H Robinson P F and Conclazio C F Estimation of Riboflavin Thiamine and N<sup>1</sup> methylnicotinamide Rapid Field Methods Ind and Eng Chem. Anal. Ed 17 384 387 (June) 1945

## A. ASCORBIC ACID IN SERUM AND URINE.

### Principle

The blue dye 6-dichlorophenolindophenol in acid solution is decolorized by ascorbic acid. An aqueous solution of the dye standardized against known amounts of ascorbic acid can therefore be used to titrate samples of serum and urine.

### Apparatus

- 1 Burette 5 ml capacity graduated to 0.5 ml  
Round bottom test tubes 15 ml capacity
- 2 Pipettes to deliver exactly 0.1 ml 0.5 ml and approx 3.5 ml Automatic syringe pipettes are the most satisfactory type
- 3 Suitable titration bench with good daylight or daylight electric shining against a white background.

### Reagents

- 1 Metaphosphoric acid approximately 3% and 6% in aqueous solution. This takes a long time to dissolve and one should start making it several hours before use. It keeps for several days.  
Buffered solution of 6-dichlorophenolindophenol. In the field this should be prepared from a dry mixture weighed out in the laboratory of 1.0 mg of the dye 1.8 gm KH<sub>2</sub>PO<sub>4</sub> and 8 gm of Na<sub>2</sub>HPO<sub>4</sub> all mixed together and stored dry in small envelopes. One envelope is dissolved in about 500 ml of water either distilled or tap if the latter's blank is known to be low. The aqueous solution of buffered dye will keep for several days but changes its standardization so that it must be standardized daily as described below.
- 2 Ascorbic acid 100 mg water soluble tablets prepared by reputable drug houses are satisfactory for making a standard solution but each fresh lot of tablets must be assayed before use. In a 100 ml mixing cylinder place about 0 ml of 3% metaphosphoric acid add 1 tablet ascorbic acid and fill accurately to the mark with water. The tablet is dissolved by inversion. After complete solution 1 ml of this strong standard is pipetted accurately into a 100 ml mixing cylinder containing about 0 ml of metaphosphoric acid. The cylinder is filled accurately to the mark and mixing is effected by inversion. This weak standard contains 1.0 mg of ascorbic acid per 100 ml and is used for standardizing the dye as described below.

### Procedure

- 1 Before each day's work the dye must be standardized against a fresh solution of ascorbic acid. 4 ml of the weak standard ascorbic acid are pipetted accurately into a 15 ml test tube. Dye is added from the burette until the first pink color lasting 30 sec is reached. Satisfactory duplicate estimations must be obtained before any further steps are taken. The titration figure should be about 3 ml. A reagent blank must be run before each day's work. 4 ml of the 3% metaphosphoric acid solution are pipetted accurately into a test tube and are titrated

against the dye. For satisfactory results the titration figure must be below 0.2 ml. If it is higher there is something wrong with the reagents most probably the water and this error must be corrected before any further steps are taken.

- 3 Estimation of vitamin C in urine. Urine is collected and stored with oxalic acid as a stabilizer as described in the section on urine. In the case of fasting specimens exactly 0.5 ml is pipetted into a test tube and approximately 3.5 ml of 3% metaphosphoric acid is added. In the case of loaded specimens exactly 0.1 ml of urine is used. In either case dye is added from the burette rather rapidly. The first definite pink lasting 15 sec is taken as the end point and no notice is taken of fading after this stage is reached. There are in urine reducing substances other than ascorbic acid which reduce the dye somewhat more slowly than does ascorbic acid. Rapid titration is therefore absolutely essential.
- 4 Estimation in serum. Exactly 2 ml of serum are added to exactly 6 ml of 6% metaphosphoric acid in a 15 ml centrifuge tube. The tube is stoppered with a clean rubber stopper and is gently inverted several times. Do not shake vigorously because of froth formation. The tube is centrifuged at high speed and exactly 4 ml of the supernatant fluid are pipetted into a test tube. Titration is carried out directly on the supernatant fluid, which already contains metaphosphoric acid.

#### Calculations

##### A. Standardization of the dye

$$\text{mg ascorbic acid equivalent to 1 ml of dye} = \frac{0.04}{\text{titration in ml} - \text{reagent blank in ml}}$$

##### B. Calculation for serum.

$$\text{mg ascorbic acid in 100 ml} = (\text{titration figure in ml} - \text{reagent blank in ml}) \times \text{equivalent of dye} \times 100$$

##### C. Calculation for fasting urine

$$\text{mg ascorbic acid in 100 ml} = (\text{titration figure in ml} - \text{reagent blank in ml}) \times \text{equivalent of dye} \times 200$$

##### D. Calculation for loaded urine

$$\text{mg ascorbic acid in 100 ml} = (\text{titration figure in ml} - \text{reagent blank in ml}) \times \text{equivalent of dye} \times 1000$$

#### Examples

- 1 Reagent blank = 0.10. 3% ml dye titrated 4 ml of weak standard therefore 1 ml dye =  $0.04 \div (3.04 - 0.10) = 0.0127$  mg ascorbic acid.
- 2 Titration figure for serum = 0.85 ml. Therefore serum contains  $(0.85 - 0.10) \times 0.0127 \times 100 = 1.0$  mg ascorbic acid per 100 ml.
- 3 Titration figure for fasting urine = 0.55 ml dye. Therefore fasting urine contains  $(0.55 - 0.10) \times 0.0127 \times 200 = 1.1$  mg ascorbic acid in 100 ml.
- 4 Titration figure for loaded urine = 1.00 ml dye. Therefore loaded urine contains  $(1.00 - 0.10) \times 0.0127 \times 1000 = 11.4$  mg ascorbic acid in 100 ml.
- 5 Use of constant factors. During a single day the standardization of the dye does not change. Hence calculation is facilitated by computing the Factor F representing the constant number by which the titration figure is multiplied to give the result directly in terms of mg ascorbic acid per 100 ml of serum, fasting urine or loaded urine. In example 1 above 1 ml of dye = 0.0127 mg ascorbic acid.

$$\text{Hence } F \text{ for serum} = (0.0127) \times (100) = 1.27$$

$$F \text{ for fasting urine} = (0.0127) \times (200) = 2.54$$

$$F \text{ for loaded urine} = (0.0127) \times (1000) = 12.7$$

# ASCORBIC ACID




Fig. 2 —Line chart for ascorbic acid. For use, see text.



- 6 Use of line chart The nomogram (Fig 22) is for serum fasting urine and loaded urine. A string is stretched between the left hand line (titration means blank) and the right hand line (titration of 0.04 mg ascorbic acid minus blank) and the concentration of ascorbic acid in 100 ml of serum fasting urine or loaded urine is read off on the middle line. Section I of the right hand line is for use with loaded urine. Section II for fasting urine and Section III of serum when the procedure is exactly as described above.

#### Precautions

- 1 The reagents must be standardized at least once daily and in hot weather twice a day.
- 2 Titration of urine must be carried on as rapidly as is consistent with accurate work. Serum need not be titrated so rapidly because it does not contain significant amounts of non-ascorbic acid reducing substances.
- 3 It is emphasized that even in individuals whose body stores of vitamin C are very low the apparent urine level by this method never reaches 0 because of the presence in urine of reducing substances which are not ascorbic acid.

### B RIBOFLAVIN IN URINE

#### Principle

Riboflavin in urine is extracted into isobutyl alcohol after a preliminary treatment in the presence of pyridine and acetic acid with potassium permanganate followed by hydrogen peroxide. The isobutyl layer is cleared with anhydrous sodium sulfate and is read fluorometrically.

#### Apparatus

- 1 Apparatus for collecting and storing urine as described in Section II.
- 2 Apparatus for visual fluorometry as described in the section on fluorometry. The comparator rack for reading solutions of riboflavin should have fixed on the side nearest the eye a suitable yellow glass filter with maximal transmission at 560 microns.
- 3 At least two dozen glass stoppered test tubes of 10 ml capacity. These are conveniently made from standard 10 ml glass stoppered footed mixing cylinders by removing the foot with a sharp tap of a tack hammer and rounding the bottom in the blast lamp.

#### Reagents (Assuming 100 specimens)

- 1 Reagents for storing urine as described in Section II.
- 2 A mixture of 1 part pyridine and 1 part glacial acetic acid. This should be made up fresh each day. 50 ml for 100 specimens.
- 3 5 ml aqueous 4% potassium permanganate.
- 4 10 ml aqueous 3% hydrogen peroxide. This keeps well and the potency of any batch should be such that two drops destroy 1 drop of 4% permanganate in the presence of acetic acid.
- 5 Isobutyl alcohol of low blank. Normal butyl alcohol is equally satisfactory. Each fresh batch of alcohol should be examined before the ultraviolet light and must be redistilled if there is perceptible fluorescence.
- 6 Anhydrous sodium sulfate 40 gm.
- 7 Water soluble tablets of riboflavin 5 mg for standards. Each batch must be assayed before use.

#### Procedure

- 1 All manipulations should be carried out in diffuse light. Only direct sunlight destroys riboflavin rapidly enough to make extreme precautions necessary against light.

- In the case of fasting urine exactly 0.5 ml of urine is pipetted into a 10 ml glass stoppered test tube. In the case of loaded urine exactly 0.1 ml is added and 0.4 ml of water to keep volume relations constant.
- 0.5 ml of pyridine acetic acid mixture is added and mixed by gentle shaking of the rack.
- 1 drop of potassium permanganate is added and mixed by a few gentle shakes of the rack. At this point the tubes may stand anywhere from  $\frac{1}{2}$  to 5 minutes without affecting results.
- Two drops of hydrogen peroxide are added and mixed by gentle shaking of the rack. The permanganate should be destroyed within a few seconds. If it is not another drop of hydrogen peroxide may be necessary. In temperatures near freezing gentle warming of the tube may be necessary. This is conveniently done in a water bath at about 70° F.
- Exactly 1.5 ml of isobutyl alcohol are added by means of either an ordinary pipette or preferably from an automatic syringe pipette which avoids getting unpleasant vapors in the mouth. The glass stopper is inserted and the tube is shaken vigorously by hand with an up and down motion 25 times. The distribution coefficient of riboflavin between the aqueous phase and the alcohol pyridine acetic acid phase is such that more shaking has been proved to be unnecessary.
- The test tube is allowed to stand in the rack until the aqueous phase is nearly separated. This usually takes from 1 to 5 minutes being faster the higher the temperature.
- A small wide mouthed funnel is placed in the neck of the tube after the glass stopper has been removed and a small spatula of sodium sulfate is added to the tube. A gentle rotary motion imparted to the tube assists clearing of the alcohol layer.
- The tube is allowed to stand in the rack until it is water clear which usually takes one to two minutes. However a few turns on a hand centrifuge will materially hasten this step.
- Approximately 1 ml of the alcohol layer is transferred to a standard thin walled Pyrex test tube of dimensions 10 x 75 mm. This is best achieved with an automatic syringe pipette as described under fluorometry but any pipette can be used.
- Reading is conducted before the ultraviolet lamp as described under fluorometry. The alcohol solutions are stable in the dark for at least hours even at high temperatures and therefore there is no hurry in reading them. Make sure that standards and unknowns are at approximately the same temperature.
- After this initial reading if one wishes to obtain true riboflavin the tubes must be exposed to strong ultraviolet light for an hour or more. Where sunshine can be counted on it may be used but usually an ultraviolet lamp is essential. If this step is omitted satisfactory comparative results are still obtained since the amount of fluorescence not destroyed by ultraviolet light is remarkably constant at about 50% of the initial value. Following loading tests because of the high concentrations of riboflavin exposure to ultraviolet light may be safely omitted and a constant factor used.

#### Standard Solutions

Two types of standards will be described. The first and theoretically best is made by running aqueous solutions through the whole method. Such standards tend to be unstable and when a large number of estimations have to be done the second type of standard is recommended. These consist of aqueous solutions of riboflavin which are considerably more stable but somewhat less desirable otherwise. Both types of standard must be checked at least once a day for deterioration.

**A Standards run through the whole method**

- 1 Prepare two liters of water containing 0.3 ml of conc H<sub>2</sub>SO<sub>4</sub> per liter. Acidulated water is used in all subsequent steps because riboflavin is most stable in 0.01 N H<sub>2</sub>SO<sub>4</sub>.
- 2 Dissolve one 5 mg tablet in exactly 1 liter of acidulated water. This strong standard contains 5 mcg riboflavin per ml.
- 3 In a glass stoppered mixing cylinder make standard solutions by adding from a pipette the amount of strong standard shown in the table below diluting to exactly 25 ml with acidulated water mixing transferring to a dry 1 oz amber bottle for storage rinsing the mixing cylinder with water and starting the next dilution in the same cylinder.

| CONCENTRATION RIBOFLAVIN<br>REQUIRED mcg/100 ml | STRONG STANDARD REQUIRED<br>ml |   |
|---|--------------------------------|---|
| 0   | 0                              | In each case dilute to 25 ml<br>with acidulated water |
| 25  | 1.25                           |   |
| 50  | 2.5                            |   |
| 75  | 3.75                           |   |
| 100   | 5.0                            |   |
| 125   | 6.5                            |   |
| 150   | 7.50                           |   |
| 200   | 10.0                           |   |
| 250   | 12.50                          |   |
| 300   | 15.0                           |   |

- 4 Run 1/2 ml samples of the standards through the method as described in steps 1 to 12 above. The final alcohol layers should be placed in clean dry glass stoppered standard test tubes the racks of which should be clearly marked with the value of the initial standard solutions.

**B Aqueous standards for use with large numbers of estimations**

- 1 Prepare two liters of water containing 0.3 ml of conc H<sub>2</sub>SO<sub>4</sub>.
- 2 Dissolve one 5 mg pill in exactly one liter of acidulated water. This strong standard contains 5 mcg riboflavin per ml.
- 3 Prepare a working standard by diluting exactly 30 ml of strong standard to exactly 300 ml with acidulated water.
- 4 Prepare a series of aqueous standards as in step A 3 but using the following table which is constructed on the assumption that the over all recovery of riboflavin in the method is 80% and that 0.5 ml samples of urine are to be used.

| STANDARD DESIRED<br>mcg RIBOFLAVIN PER 100 ml | WORKING STANDARD REQUIRED<br>ml |   |
|---|---------------------------------|---|
| 0   | 0                               | In each case dilute to exactly<br>25 ml with acidulated water |
| 25  | 2.0                             |   |
| 50  | 4.0                             |   |
| 75  | 6.5                             |   |
| 100   | 8.5                             |   |
| 125   | 10.5                            |   |
| 150   | 12.5                            |   |
| 200   | 16.5                            |   |
| 250   | 21.0                            |   |
| 300   | 25.0                            |   |

- 5 Store these in brown bottles in a cool place. They will keep for one week.
- 6 For use as standards the aqueous solutions as they stand should be placed directly in the standard tubes without running them through the method. The racks of the tubes should be labeled with the numbers in the left hand column. The values then represent mcg of riboflavin per 100 ml of urine when 0.5 ml samples are analysed.
- 7 Replace the standards in the tubes morning and afternoon when large numbers of estimations are being done.

**Calculations**

The calculations for fasting urines are as follows

$\text{mcg}/100 \text{ ml urine} = (\text{1st reading} - \text{reading after 1 hr in UV light})$  If the sample has not been exposed to sunlight or ultraviolet light  $\text{mcg}/100 \text{ ml urine} = (\text{1st reading}) \times (0.5)$

The calculations for loaded urines are as follows

$\text{mcg}/100 \text{ ml urine} = (\text{1st reading} - \text{reading after 1 hr in UV light}) \times 5$  If the sample has not been exposed to sunlight or ultraviolet light  $\text{mcg}/100 \text{ ml urine} = (\text{1st reading} \times 5) - (1.5)$

**Examples**

Fasting Urine First reading 50 Reading after sunshine 1.5 Therefore  $\text{mcg riboflavin}/100 \text{ ml} = (50 - 1.5) = 48.5$

Loaded Urine First reading 90 Reading after sunshine 5 Therefore  $\text{mcg riboflavin}/100 \text{ ml} = (90 - 5) \times 5 = 425$

**Precautions**

- 1 The commonest source of error in this method is the instability of the standards. These must be checked at least every day and when many readings are being made twice a day. Our experience has been that an ultraviolet bulb is far preferable to sunlight because of the extreme variability of the sunshine in climates where there is intermittent cloudiness.
- 2 Standards and unknowns must be at the same temperature when read.

**THIAMINE AND *N*' METHYLNICOTINAMIDE****Principle**

These two substances are estimated by methods so similar that they are best described together except for the last steps. Both compounds are absorbed from urine on activated zeolite and eluted with potassium chloride. Thiamine is converted to thiochrome by ferricyanide in alkaline solution. *N*'methylnicotinamide (hereafter referred to as NMN) is itself nonfluorescent but is converted to the fluorescent Factor F by mere addition of alkali. The thiochrome and Factor F are shaken into isobutyl alcohol and are estimated fluorometrically.

**Apparatus**

- 1 Apparatus for collection and storage of urine as described in Section II.
- 2 Apparatus for visual fluorometry as described in Section III F.
- 3 Glass stoppered 10 ml test tubes as under riboflavin.
- 4 Beakers, pipettes, funnels and if available automatic syringe pipettes as under riboflavin.

**Reagents**

- 1 Activated zeolite. This requires considerable time, patience and sometimes experimentation in order to get a satisfactory product which will settle out rapidly and also absorb thiamine well. Commercial zeolite is crushed in a stone crusher or suitable mortar (In our experience permutit and decalco are as good as zeolite.) That portion of the product which will pass a 100 mesh sieve is saved and suspended in 1% acetic acid in distilled water in large cylinders. The granules which settle to the bottom rapidly that is in about 3 minutes are separated from the lighter particles by decanting. These heavier particles are boiled three times with 1% acetic acid with settling and decanting between fresh additions of acetic acid. The product is washed with distilled water and dried in the oven at 110°C. It should be a white homogeneous powder and when about 200 mg are shaken with water in a 10 ml mixing cylinder all the particles should settle rapidly to the bottom in not more than 2 min. We have encountered great vari-

bility from batch to batch. When test runs with synthetic thiamine are made recovery is sometimes poor and the product may have to be reactivated by the technique of Najar and Wood (1940). Sometimes addition of ferricyanide results in the formation of Prussian blue. In such cases the product has to be washed with warm fairly strong solutions of hydrochloric or nitric acid. When a satisfactory product is finally obtained its activity remains constant for months if it is kept dry.

- 2 Aqueous potassium chloride approximately 5%. At temperatures near freezing the solubility properties necessitate lower concentrations.
- 3 Aqueous sodium hydroxide approximately 15%.
- 4 Bromocresol green, 1% in water.
- 5 Aqueous potassium ferricyanide 0.25%. This should be made up fresh every day, most conveniently from weighed samples stored dry in small stoppered test tubes.
- 6 Isobutyl alcohol. This must have a blank of essentially zero and each new batch has to be tested before use. It can be recovered by redistillation of the residue after analysis. Normal butyl alcohol can be used instead of isobutyl.
- 7 Acidulated water. This is most conveniently made up a gallon at a time. Approximately 10 ml of glacial acid are added to a gallon of distilled water. Tap water of low blank may be used but is usually inferior to distilled water.

#### Procedure

- 1 Urine is collected and stored as described in Section II. Thiamine and NMN are stable for at least a week even at high temperatures provided the pH of the sample is between 3 and 6. This pH is almost always attained with the oxalic acid in the storage bottle.
- 2 The acidity of the urine is checked soon after collection with bromocresol green as external indicator. The acidity is correct when the indicator turns any shade of green or yellow. If it is blue add glacial acetic acid drop by drop to the storage bottle until the acidity is correct.
- 3 Steps 4 through 8 which follow are identical for NMN and thiamine. It is therefore convenient to run both at the same time but with separate aliquots until the final stages which will be differentiated below. 2 ml of fasting urine are pipetted into a glass stoppered test tube or 0.5 ml of loaded urine.
- 4 About 100 mg of zeolite are delivered from a small spatula into the tube most easily through a small funnel and mixed by about ten rapid vigorous shakes of the rack. The time of standing at this stage is not critical up to 1/2 hr. No more than ten shakes are required.
- 5 Approximately 8 ml of acidulated water are added to the tube from a large syringe and the tube is mixed by inversion 10 times with the finger as a stopper. After a few seconds standing in the rack until the particles of zeolite have begun to settle the few particles adherent to the top and sides are washed down by closing the top with a clean finger and administering a single sharp upward jerk. The tube is replaced in the rack until the zeolite particles have settled to the bottom. In urine containing oxalic acid there is usually slight permanent turbidity which is disregarded in deciding when the zeolite has settled.
- 6 The supernatant fluid is withdrawn by means of suction through a long needle of stainless steel and is discarded. A mechanical vacuum pump is highly desirable but a 50 or 100 ml syringe may be used if necessary.
- 7 Approximately 8 ml of acidulated water are added and the processes of standing, inversion and removal of the supernatant fluid are repeated. This washing step is critical in both the estimation of NMN and thiamine but for different reasons. In the case of thiamine if washing is not effective the fluorescence obtained in the final step is off color with a silvery blue admixture to the true thiochrome mauve and reading is wholly unsatisfactory. In our experience tap water can be

used instead of distilled water but it may be necessary in this case to add more acetic acid to the water or even sometimes to wash three times before the thiochrome color is satisfactory. Thiamine is firmly bound to the zeolite and repeated washings remove the interfering substances without affecting the thiamine. Experimentation determines how many washings with the available water supply are required in order to get satisfactory results for thiamine. In the case of NMN in contrast to thiamine each washing diminishes the final fluorescence. Apparently NMN is relatively easily washed off the zeolite. In making routine surveys it is therefore essential to adopt and follow scrupulously a fixed technique of washing in order to be able to compare values for NMN from one place to another.

- 8 0.5 ml of potassium chloride is added to the tube and the rack is shaken gently so as to mix the zeolite with the potassium chloride but not to streak the zeolite far up the sides. The time necessary for elution is short not exceeding 30 seconds. From this point on the estimations of NMN and thiamine are different.

#### Final steps for thiamine

- 9 0.1 ml of potassium ferricyanide is added and mixed by a few gentle shakes of the rack.
- 10 0.5 ml of sodium hydroxide is added and mixed with a gentle shake of the rack. The time necessary for completing the oxidation from thiamine to thiochrome is only a few seconds. Once formed in this mixture it is stable for at least 1/2 hour.
- 11 Add 1 ml of isobutyl alcohol. The glass stopper is inserted and the tube is shaken vigorously 20 times up and down its long axis. This number of shakes insures complete distribution of thiochrome between the alcohol and the aqueous phases. More shaking is superfluous.
- 12 The tube is allowed to stand in the rack until separation of the phases is almost complete and is then centrifuged if necessary a few seconds in a hand centrifuge to insure complete clearing.
- 13 About 1 ml of the supernatant fluid is transferred to a small test tube by means of a pipette or a syringe and needle.
- 14 The fluorescence is matched against standards in the visual comparator as described in the section on fluorometry care being taken that standards and unknowns are at approximately the same temperature.

#### Final steps for *N*<sup>1</sup> methylnicotinamide

- 9 1 ml of isobutyl alcohol are delivered into the tube. No mixing is necessary at this point.
- 10 0.5 ml of sodium hydroxide is added the tube is stoppered with its glass stopper and is shaken 10 times as with thiamine. The addition of sodium hydroxide should be followed by shaking at once because Factor F in aqueous alkaline solution breaks down to nonfluorescent products within a very few minutes but is stable in the isobutyl layer for days.
- 11 The tube is replaced in the rack and when separation of the phases is more or less complete it is centrifuged if necessary for a few seconds in a hand centrifuge.
- 12 1 ml of the supernatant layer is transferred to small test tube as under thiamine.
- 13 Fluorescence is matched in the comparator against standard solutions as described in the section on fluorometry. The color of F<sub>2</sub> develops relatively slowly after step 10. It is absolutely necessary to allow at least 5 minutes between step 10 and the final reading. The development of color goes on at the same rate either at step 11 or at step 12 becomes maximal in 5 minutes and then remains constant for at least an hour. Care should be taken that standards and unknowns are at approximately the same temperature.

**Standards****A For thiamine**

In order to compensate for variability in the reagents standards are made by running solutions of thiamine hydrochloride through the method

- 1 Prepare 2 liters of acidulated distilled water by adding approximately 20 ml of glacial acetic to 2 liters of distilled water. This acidulated water is used in all subsequent steps since the vitamin is maximally stable at about pH 4 and is quickly destroyed above pH 7
- 2 One 5 mg tablet of thiamine hydrochloride is dissolved in exactly one liter of acidulated water. This strong standard contains 5 mcg of thiamine hydrochloride per ml
- 3 Dilute exactly 75 ml of strong standard to exactly 250 ml with acidulated water. This working standard contains 15 mcg of thiamine hydrochloride per ml
- 4 Using a mixing cylinder and graduated pipettes make up the following series of standards

CONCENTRATION OF THIAMINE HYDROCHLORIDE WORKING STANDARD REQUIRED  
DESIRED (mcg/100 ml)

|     | (ml) |  |
|-----|------|--|
| 0   | 0    | In each case dilute to exactly 25 ml with acidulated water |
| 5   | 0.83 |  |
| 10  | 1.67 |  |
| 15  | 2.5  |  |
| 20  | 3.33 |  |
| 30  | 5.0  |  |
| 50  | 8.3  |  |
| 75  | 12.5 |  |
| 100 | 16.7 |  |
| 150 | 25.0 |  |

- 5 The above standards are stored in 1 oz amber bottles in a dark cool place. They keep for at least 2 weeks
- 6 Two ml of each of the series are run through the thiamine method. One ml of the isobutyl layer is placed in the standard comparator tube the rack of which is clearly labeled with the concentration of the original solution. Standards must be checked daily for deterioration

**B For N<sup>1</sup> methylnicotinamide**

- 1 Dissolve 10.0 mg N<sup>1</sup> methylnicotinamide chloride in exactly 133 ml of 1% acetic acid in water
- 2 Make up the following for storage in 1 oz brown bottles

CONCENTRATION DESIRED  
(mg PER 100 ml)

|      | ml   |   |
|------|------|---|
| 0    | 0    | In each case dilute to exactly 25.0 ml with 1% acetic acid. |
| 0.25 | 0.8  |   |
| 0.50 | 1.7  |   |
| 0.75 | 2.5  |   |
| 1.0  | 3.3  |   |
| 1.5  | 5.0  |   |
| 2.5  | 8.3  |   |
| 3.75 | 12.5 |   |
| 5.0  | 16.6 |   |
| 7.5  | 25.0 |   |

- 3 Store these standards in a cool place. They keep for 12 weeks. N<sup>1</sup> methyl nicotinamide is not fluorescent but is converted to F<sup>+</sup> in alkaline isobutyl alcohol

**Calculations**

- 1 Fasting specimens of urine  
mcg of thiamine per 100 ml and mg of NMN per 100 ml = reading on rack of standard.

Loaded urine

mcg of thiamine per 100 ml and mcg of NMN per 100 ml = reading on rack of standard  $\times 4$ .

- 3 On theoretical grounds most of the standard laboratory methods subtract *F* from thiochrome in calculating thiamine. We do not make this correction in the present method because of the demonstration by Hajjar and Keiron (1944) that *F* is about 80% destroyed by the addition of ferricyanide. Therefore only when the *F* is large in comparison with *B* is correction for *F* necessary. This has very rarely happened in our field surveys.

#### Examples

- 1 Fasting urine

Readings *B* 15 *F* 15 Therefore mcg *B* per 100 ml urine = 30 and mcg NMN per 100 ml = 15

- 2 Loaded urine

Readings *B* 30 *F* 15 Therefore mcg *B* per 100 ml urine =  $4 \times 30 = 120$   
mcg NMN per 100 ml urine =  $4 \times 40 = 160$

#### Precautions

##### A For both thiamine and NMN

- 1 The chief sources of error in these methods are two—the difficulty of getting a satisfactory preparation of zeolite and the washing stages. Both of these difficulties have to be settled by experimentation.
- 2 Once prepared the standards for thiamine usually are stable for one to two weeks. They should be checked every day or two against freshly prepared standards.
- 3 The standards and unknowns should be at the same temperature when read because fluorescence is strongly affected by changes in temperature.

##### B For thiamine

- 1 For urine of subjects taking quinine or atabrine preliminary treatment of the specimen is essential.
  - (a) Pour approximately 7 ml of urine into a centrifuge tube.
  - (b) Add approximately 35 ml of isobutyl alcohol and about 10 mg Norite ( $\frac{1}{4}$  of a small spatula). The exact amount of charcoal is critical as too much will destroy the thiamine.
  - (c) Close with thumb and shake 5 times.
  - (d) Centrifuge and use 2 ml samples of the water layer below the alcohol layer.

##### C For *N*-methylnicotinamide

- 1 In the final step be sure to use the following order: 0.5 ml KCl mix, 2 ml isobutyl alcohol, 0.2 ml NaOH. Shake at once. *F* is unstable in watery alkaline solutions but stable in the isobutyl layer. After the shaking with NaOH *F* develops relatively slowly in the isobutyl layer; at least 5 minutes must elapse after the shaking before the specimen is read. You may either let the glass stoppered tubes sit for 5 minutes or else you may draw off the isobutyl layer into the reading tubes and read after 5 minutes. A minimum of 5 minutes waiting is absolutely essential. The color is stable up to one hour so that speed in reading is not necessary.
- 2 For urine of subjects taking quinine or atabrine preliminary treatment with charcoal is essential. The procedure is:
  - (a) Pour approx. 5 ml of urine into a centrifuge tube.
  - (b) Add approx. 10 mg Norite A (about  $\frac{1}{4}$  of a small spatula). The exact amount is not critical for *F* but enough is essential.
  - (c) Close with thumb and shake 25 times.
  - (d) Centrifuge and use 2 ml samples of supernatant fluid.



9 **Vitamins and Nutrition Miscellaneous References****Books**

- American Medical Association *The Vitamins A Symposium* Chicago 1939 American Medical Association
- Association of Official Agricultural Chemists *Official and Tentative Methods of Analysis of the Association of Official Agricultural Chemists* ed 6 Washington 1945 AOAC
- Association of Vitamin Chemists Inc *Methods of Vitamin Assay*, New York 1947 Interscience
- Bicknell F and Prescott F *The Vitamins in Medicine* ed 2 reprinted 1947 New York 1941 Grune and Stratton
- Coward K H *The Biological Standardization of the Vitamins* ed 2 London 1947 Bailliere Tindall and Cox
- Dann W J and Satterfield G H *Biological Symposia VII Estimation of the Vitamins* New York 1947 Jacques Cattell Press vol VII
- György P *Vitamin Methods* New York 1950 Academic Press Inc Publishers vol I
- McLester J S *Nutrition and Diet in Health and Disease* ed 5 Philadelphia and London 1949, W B Saunders Company
- Rosenberg H B *Chemistry and Physiology of the Vitamins* New York 1946 Interscience
- Sahyun M *Proteins and Amino Acids in Nutrition* New York 1948 Reinhold Publishers Vol xv
- Sherman H C *Chemistry of Food and Nutrition* ed 7 New York 1946 The Macmillan Company
- Williams P R and Spies T D *Vitamin B (Thiamin) and Its Use in Medicine* New York 1940 The Macmillan Company
- U S Pharmacopoeial Conventions *The Pharmacopoeia of the U S of America* 13th Revision (USP XIII) Washington 1947
- Younts J B *Nutritional Deficiencies Diagnosis and Treatment* ed 2 Philadelphia Montreal London 1943 J B Lippincott Company

**Ascorbic Acid**

- Adams J Acker M and Friedman H A *Iodimetric Titrimetric Determination of Ascorbic Acid* *J Am Pharmaceut Assoc Sci Ed* 36 (6) 170 172 (July) 1947
- Be sey O A *A Method for the Determination of Small Quantities of Ascorbic Acid and Dehydroascorbic Acid in Turbid and Colored Solutions in the Presence of Other Reducing Substances* *J Biol Chem* 126 (2) 771 784 (Dec) 1933
- Farmer C J and Abt A F *Determination of Reduced Ascorbic Acid in Small Amounts of Blood* *Proc Soc Exp Biol & Med* 34 146 150 (March) 1936
- Hochberg M Melnick D and Oser B L *Photometric Determination of Reduced and Total Ascorbic Acid* *Ind Eng Chem (Anal Ed)* 15 (3) 182 183 (March) 1943
- Huelin F L and Stephens I M *The Influence of Ferrous Iron in the Determination of Ascorbic Acid* *Australian J Exper Biol & M Sc* 25 (1) 17 23 (March) 1947
- King C G *Chemical Methods for Determination of Vitamin C* *Ind Eng Chem (Anal Ed)* 13 (4) 5 997 (April) 1941
- Mills M B and Roe J H *A Critical Study of Proposed Modifications of the Roe and Kuehner Method for the Determination of Ascorbic Acid With Further Contributions to the Chemistry of This Procedure* *J Biol Chem* 170 (1) 153 164 (Sept) 1947
- Mundlin F L and Butler A M *The Determination of Ascorbic Acid in Plasma A Macromethod and Micromethod* *J Biol Chem* 122 6 3 636 (Feb) 1933
- Pecover V *The Microdetermination of Ascorbic Acid in Blood Plasma* *Australian J Exper Biol & M Sc* 25 (2) 175 178 (June) 1947

**Carotene**

- Johnson R M and Baumann C A *Studies on the Reaction of Certain Carotenoids With Antimony Trichloride* *J Biol Chem* 169 (1) 83 90 (June) 1947

**Choline**

- Marenzi A D and Cardini C E *The Colorimetric Determination of Choline* *J Biol Chem* 147 (2) 363 3 0 (Feb) 1943
- Entenman C Taurig A and Chaskoff I L *The Determination of Choline in Phospholipids* *J Biol Chem* 155 (1) 13 18 (Sept) 1944

**Folic Acid**

- Allfrey V Topley L J Gessen C and King C G *A Fluorometric Method for the Determination of Pteroylglutamic Acid* *J Biol Chem* 178 465 481 (March) 1949
- Olazake A J and Wolf L M *Colorimetric Determination of Folic Acid and Adenine* *Arch. Biochem.* 21 241 243 (March) 1949

**Nicotinic Acid**

- Banerjee S Ghosh N C and Bhattacharya G A Chemical Method of the Estimation of Nicotinic Acid in the Presence of Sugar J Biol Chem 172 (2) 495 499 (Feb) 1945
- Harris L J and Raymond W D Assessment of the Level of Nutrition A Method for the Estimation of Nicotinic Acid in Urine Biochem J 33 (1) 537 551 (Dec) 1939
- Mueller A and Fox S H Chemical Determination of Niacin J Biol Chem 167 (1) 291 9 (Jan) 1947
- Perlzweig W A Levy F D and Sarett H P Nicotinic Acid Derivatives in Human Urine and Their Determination J Biol Chem 136 7 9745 (Dec) 1940
- Teeri A F and Shimer S R A Colorimetric Determination of Nicotinic Acid J Biol Chem 163 (1) 307 311 (April) 1944

**Pyridone of N<sup>1</sup> Methylnicotinamide**

- Rosen F Perlzweig W A and Handler P A Fluorometric Assay for N<sup>1</sup> Methyl 6 pyridone 3 carboxylamide Proc Fed Soc 7 181 (March) 1943
- Rosen F Perlzweig W A and Leder I G A Fluorometric Method for the Determination of the 6 Pyridone of N<sup>1</sup> Methylnicotinamide in Urine J Biol Chem. 179 157 168 (May) 1944

**4 Pyridoxic Acid**

- Huff J W and Perlzweig W A A Product of Oxidative Metabolism of Pyridoxine 2 Methyl 3 hydroxy 4 carboxy 5 hydroxymethylpyridine (4 pyridoxic acid) I Isolation From Urine Structure and Synthesis J Biol Chem 153 343 355 (Sept) 1944

**Pyridoxine**

- Hochberg H, Melnick D and Oser B L Chemical Determination of Pyridoxine in Biological Materials and Pharmaceutical Products The Multiple Nature of Vitamin B J Biol Chem 155 (1) 119 1 9 (Sept) 1944

**Riboflavin**

- Emmett A D Bird O D Brown R A Peacock G and Vandenbelt J M Determination of Vitamin B<sub>2</sub> (Riboflavin) Comparison of Bioassay Microbiological and Fluorometric Methods Ind. Eng. Chem. (Anal. Ed.) 13 (4) 192 1 (April) 1941
- Loy H W Jr Report on Pilo flavin (Microbiological)—Revised Method (Chemical)—Fluorometric Method Assoc. of Offic. Agric. Chemists 57 (3) 461 464 1949
- Morell D B and Slater F C The Fluorimetric Determination of Riboflavin in Urine Biochem. J 40 (5/6) C, 657 1946
- Rubin S H De Ritter F Schuman F L and Bauernfeind J C Determination of Riboflavin in Low Potency Foods and Feeds Ind. Eng. Chem. (Anal. Ed.) 17 (3) 1 6 140 (March) 1945
- Slater F C and Morell D B A Modification of the Fluorimetric Method of Determining Riboflavin in Biological Materials Biochem. J 40 644 65 1946

**Thiamine**

- Fgassa F and Meiklejohn A P The Estimation of Thiamine in Urine J Biol Chem 141 853 870 (Dec) 1941
- Friedemann T L and Kmiecik T C The Determination of Thiamine in Blood J Lab & Clin. Med. 28 1 60 1 63 (July) 1943
- Hennery D J Chemical Methods for the Determination of Vitamin B Ind. Eng. Chem. (Anal. Ed.) 13 16 13 (April) 1941
- Hennessey D J, and Cerecedo L R The Determination of Free and Phosphorylated Thiamin by a Modified Thiochrome Assay J Am. Chem. Soc. 61 (1) 1 9 183 (Jan) 1939
- Hochberg M Melnick D, and Oser B I Simplified Colorimetric Determination of Thiamine in Cereal Product Cereal Chem 22 ( ) 83 90 (March) 1945
- Mickelson O Condiff, H and Keys, A. Determination of Thiamine in Urine by Means of the Thiochrome Technique J Biol. Chem 160 (1) 361 370 (Sept) 1945

**Tocopherol**

- Bencze R The Determination of Tocopherol by the Color Reagent of Furter and Meyer Z Physiol Chem. 232 15 157 194
- Minot C The Determination of Tocopherol in Blood Serum J Lab & Clin. Med. 29 (7) 700 (July) 1944
- Qualsbein M L and Buchler R A Simplified Hydrogenation Technique for the Determination of Blood Plasma Tocopherols J Biol. Chem 159 (3) 663 665 (Aug) 1945

Quaife M L, and Harris P L The Chemical Estimation of Tocopherols in Blood Plasma *J Biol Chem* 156 (2) 499 505 (Dec) 1944

Quaife M L Scrimshaw, N S and Lowry O H A Micromethod for Assay of Total Tocopherols in Blood Serum *J Biol Chem* 180 1-9 1935 (Oct) 1949

#### Trigonelline

Kodicek E and Wang Y L A New Method for the Estimation of Trigonelline in Urine and Foodstuffs *Nature* 148 324 (July) 1941

#### Vitamin A and Carotene

Adlersberg D Kann S Maurer A P Newerly K Winternitz W and Sobotka, H Studies on Serum Carotene in Man *Am J Digest Dis* 16 (9) 333 337 (Sept) 1949

Clausen S W, and McCoord A B The Carotinoids and Vitamin A in the Blood *J Pediat* 13 (5) 635 650 (Nov) 1938

Kawasaki, C and Suetsugu S Colorimetric Estimation of Vitamin A by Glycero-dichlorohydrin Reagent *J Pharm Soc Japan* 69 460 463 1949

Koch W and Kaplan D A Simultaneous Carr Price Reaction for the Determination of Vitamin A *J Biol Chem* 173 (1) 363 369 (March) 1948

Thompson, S Y A Photoelectric Spectrophotometer Suitable for the Measurement of Vitamin A by the Antimony Trichloride Reaction, *British J Nut* 3 (1) 43 50 1949

#### Vitamin D

DeWitt J B and Sullivan M X Spectrophotometric Procedure for Quantitative Estimation of Vitamin D *Ind Eng Chem (Anal Ed)* 18 117 119 (Feb) 1946

Muller P B The Colorimetric Determination of Vitamins A D and  $\beta$  Carotene With Particular Regard to Vitamin D *Helv Chim Acta* 30 (5) 117-1190 1949

Sobel A E Mayer A M and Kramer B A New Colorimetric Reaction of Vitamin D and D and Their Provitamins *Ind Eng Chem (Anal Ed)* 17 (3) 160 165 (March) 1945

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### F HORMONES

##### 1 Ketosteroids in Urine

###### Reference

Consolazio W V and Talbott J H The Extraction and Determination of Ketosteroids in Urine *J Endocrinol* 27 335-39 (Sept.) 1940

###### Principle

The steroid linked to glucuronic acid is set free by hydrolyzing with hydrochloric acid and is extracted with carbon tetrachloride. Interfering substances are washed out and the resulting extract is equilibrated with alkaline meta-dinitrobenzene. The final color is measured.

###### Apparatus

- 1 One continuous extractor unit. (See reference Made by MacAlaster Dicknell Company Cambridge Mass See Fig 3)
- 2 Steam hydrolyzer (According to Consolazio and Talbott.)
- 3 All glass still to fit round bottom extraction flask
- 4 Automatic syringe pipette 0.1 ml capacity  
Six pear-shaped separatory funnels (50 ml capacity) with interchangeable ground glass stoppers.
- 6 Swivel separatory funnel rack
- 6 Coleman Jr Spectrophotometer
- 8 A quantity of 19 x 100 mm cuvettes.
- 9 A quantity of unglazed porcelain
- 10 A water bath set at 3 C.
- 11 An electric centrifuge
- 1 A steam bath

###### Reagents

- 1 Potassium hydroxide (KOH) exactly 2 N  
Ethyl alcohol (EtOH) aldehyde free
- 3 Meta-dinitrobenzene  $C_6H_3(NO_2)_2$  2% Dissolve 20 gm in 100 ml of EtOH. Heat to 40 C and add 100 ml of 2 N NaOH After 5 minutes cool and add 500 ml of water Filter and wash thoroughly with water Dry by suction. Recrystallize twice from 10 ml and 50 ml of absolute aldehyde free EtOH using Norite Melting point 90.5-91.0 C. A mixture of 1% solution with equal portions of aqueous 2 N NaOH gives no color after 60 minutes. A 2% solution in absolute aldehyde free EtOH stored in the dark may be used for 10-14 days.
- 4 Hydrochloric acid (HCl) concentrated  
Carbon tetrachloride ( $CCl_4$ ) C.I.
- 6 Ethyl alcohol ( $C_2H_5OH$ ) absolute
- 7 Sodium hydroxide (NaOH) approximately 2 N Dissolve 80 gm in 1000 ml water

- Quaife M L, and Harris P L. The Chemical Estimation of Tocopherols in Blood Plasma, *J Biol Chem* 156 (2) 499-505 (Dec) 1944
- Quaife, M L., Scrimshaw, N S and Lowry O H. A Micromethod for Assay of Total Tocopherols in Blood Serum *J Biol Chem* 180 1-9 1235 (Oct) 1949

**Trigonelline**

- Kodicek E, and Vang, I L. A New Method for the Estimation of Trigonelline in Urine and Food stuffs, *Nature* 148 23-24 (July) 1941

**Vitamin A and Carotene**

- Adlersberg D, Hann, S, Maurer A P, Nevery K, Winternitz, W, and Sobotka H. Studies on Serum Carotene in Man. *Am. J Digest Dis* 16 (9) 333-337 (Sept) 1949
- Clausen S W and McCoord A B. The Carotinoids and Vitamin A in the Blood, *J Pediatr* 13 (5) 635-640 (Nov) 1938
- Kawasaki C, and Suenaga S. Colorimetric Estimation of Vitamin A by Glycerodichlorohydrin Reagent *J Pharm Soc Japan* 69 460-461 1949
- Koch W and Kaplan, D. A Simultaneous Carr Price Reaction for the Determination of Vitamin A *J Biol Chem* 173 (1) 363-369 (March) 1948
- Thompson S I. A Photoelectric Spectrophotometer Suitable for the Measurement of Vitamin A by the Antimony Trichloride Reaction *British J Nut* 3 (1) 43-50, 1949

**Vitamin D**

- DeWitt, J B and Sullivan M Y. Spectrophotometric Procedure for Quantitative Estimation of Vitamin D. *Ind Eng Chem (Anal Ed)* 18 11, 119 (Feb) 1946
- Muller P B. The Colorimetric Determination of Vitamins A, D and  $\beta$  Carotene With Particular Regard to Vitamin D. *Helv Chim Acta* 30 (5) 1112-1190 1949
- Sobel A E, Mayer A M and Kramer B. A New Colorimetric Reaction of Vitamin D<sub>2</sub> and D and Their Provitamin. *Ind Eng Chem. (Anal Ed)* 17 (3) 160-165 (March) 1945

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### F HORMONES

##### 1 Ketosteroids in Urine

###### Reference

Consolazio W V and Talbott J H The Extraction and Determination of 17 Ketosteroids in Urine J Endocrinol 27 355 9 (Sept) 1940

###### Principle

The steroid linked to glucuronic acid is set free by hydrolyzing with hydrochloric acid and is extracted with carbon tetrachloride. Interfering substances are washed out and the resulting extract is equilibrated with alkaline meta-dinitrobenzene. The final color is measured.

###### Apparatus

- 1 One continuous extractor unit (See reference. Made by MacAlaster Bicknell Company, Cambridge, Mass. See Fig. 3)
- 2 Steam hydrolyzer (According to Consolazio and Talbott)
- 3 All glass set to fit round bottom extraction flask
- 4 Automatic syringe pipette 0.1 ml capacity
- 5 Six pear-shaped separatory funnels (250 ml capacity) with interchangeable ground glass stoppers
- 6 Swivel separatory funnel rack
- 7 Coleman Jr Spectrophotometer
- 8 A quantity of 19 x 150 mm cuvettes
- 9 A quantity of unglazed porcelain
- 10 A water bath set at 25 C
- 11 An electric centrifuge
- 12 A steam bath

###### Reagents

- 1 Potassium hydroxide (KOH) exactly 2 N  
Ethyl alcohol (EtOH) aldehyde free
- 3 Meta-dinitrobenzene  $\text{CH}(\text{NO}_2)_2$  2% Dissolve 0.2 gm in 750 ml of EtOH. Heat to 40 C and add 100 ml of 2 N NaOH. After 5 minutes cool and add 2500 ml of water. Filter and wash thoroughly with water. Dry by suction. Recrystallize twice from 120 ml and 10 ml of absolute aldehyde-free EtOH using Norite. Melting point 90.5-91.0 C. A mixture of 1% solution with equal parts of aqueous 0.1 N NaOH gives no color after 60 minutes. A 1% solution in absolute aldehyde-free EtOH stored in the dark may be used for 10-14 days.
- 4 Hydrochloric acid (HCl) concentrated
- 5 Carbon tetrachloride ( $\text{CCl}_4$ ) C.P.
- 6 Ethyl alcohol ( $\text{C}_2\text{H}_5\text{OH}$ ) absolute
- 7 Sodium hydroxide (NaOH) approximately 2 N Dissolve 80 gm in 1000 ml water




Fig. 23—Extraction unit (Reproduced by permission of the manufacturers)

- 8 Sodium hydrosulfite (Lycopen) (Obtained from Central Scientific Company)  
To 100 ml of 1 N NaOH add 10 gm of Lycopen Mix quickly pour into a separatory funnel and cover with light mineral oil
- 9 And osterone standard. Dissolve 250 mg androsterone in 1000 ml aldehyde free absolute EtOH One ml = 0.25 mg
- 10 Mineral oil.
- 11 Norite

#### Collection and Storage of Urine

Collect 3 overnight specimens on successive nights keeping the specimens cold until brought to the laboratory Use no preservative Pool the two specimens and if they measure not more than 1000 ml extract the total pool No sample should be assayed which shows signs of putrefaction If necessary to keep the urine as dry as possible and keep it in a refrigerator

#### Procedure

- 1 Into a steam hydrolyzer measure 500 ml of urine and 50 ml of concentrated HCl Hydrolyze with steam for 30 minutes
- 2 Cool with running water
- 3 Add approximately 50 ml of CCl<sub>4</sub> to extractor and approximately 75 ml of CCl<sub>4</sub> to a round bottomed flask
- 4 Run hydrolyzed urine into extractor make proper connections and extract for 60 minutes
- 5 Into 4 separatory funnels pipette 10 ml of 2 N sodium hydroxide
- 6 Suck CCl<sub>4</sub> extract into first funnel and wash extract flask with two 10 ml portions of CCl<sub>4</sub>
- 7 Wash 17 ketosteroids (a) 4 times with 10 ml of 2 N sodium hydroxide (b) 1 time with 10 ml of 10% Lycopen in 1 N NaOH until the pigment disappears Shake for 30 seconds longer Pigments disappear when all oxygen has been absorbed Do not open separatory funnel until complete separation has occurred Then deliver rapidly into a separatory funnel containing 550 ml of water If this step is not carried out as described the red pigment appears again (c) Wash with water
- 8 After washing with water let settle for 10 minutes and deliver directly into distilling flask Add a piece of unglazed porcelain
- 9 Distill off CCl<sub>4</sub>
- 10 Dry thoroughly with suction and dilute to 50 ml or 100 ml with aldehyde free EtOH
- 11 Into dry centrifuge tube pipette 0.2 ml of alcohol extract
- 12 Add 0.2 ml of 2% meta dinitrobenzene
- 13 Add 0.2 ml of 2 N potassium hydroxide Stopper with new clean cork
- 14 Equilibrate for 5 minutes at 25°C
- 15 Add 100 ml absolute EtOH Mix.
- 16 Centrifuge Transfer to a cuvette (Caution Centrifuge must be cold.)
- 17 Read within 15 minutes at 500 mμ setting the blank at 100% T

#### Standard Curve

- 1 0.2 ml of androsterone in EtOH.
- 2 0.2 ml of meta dinitrobenzene in alcohol
- 3 0.2 ml of 2 N KOH
- 4 Equilibrate 5 minutes at 25°C
- 5 Add 100 ml of absolute EtOH (If androsterone is dissolved in more than 0.2 ml of EtOH add androsterone to tube and evaporate EtOH with air current Then add 0.2 ml aldehyde free EtOH 0.2 ml meta dinitrobenzene and proceed) The points in the standard curve range from 0.05 mg through 0.5 mg of androsterone in the cuvette



**Calculation**

$$\text{mg 17 Ketosteroids/1000 ml} = \frac{(\text{mcg in cuvette}) \times (\text{total ml extract}) \times 1000}{1000 \times (\text{ml aliquot}) \times (\text{ml urine sample})}$$

**Example**

Volume of urine extracted 500 ml    Volume of extract 10 ml    Volume of extract assayed 0.1 ml    Reading 0.67 equivalent to 180 mcg in cuvette

$$\text{mg 17 Ketosteroids/1000 ml urine} = \frac{180 \times 10 \times 1000}{1000 \times 0.1 \times 500} = 18 \text{ mg/1000 ml}$$

**Precautions**

- 1 All glassware must be chemically clean especially after separation of urine. The glassware used after the ether phase must be dry. A useful additional precaution is to heat glassware in the oven to destroy any steroids. Stop-cock grease should not come in contact with solutions if possible. It is necessary to grease stopcocks each time they are used.
- 2 Cork stoppers must never come in contact with solutions. Use tinfoil over them or use ground glass junctions.

**2 Hormones Miscellaneous References****Acetylcholine**

- Tower D B and McEachern, D. Experiences With the 'Venus Heart Method for Determining Acetylcholine. *Canad. J Research* 26B 183 187 (April) 1948
- Wait R B. The Action of Acetylcholine on the Isolated Heart of *Venus Mercenaria*. *Biol Bull* 85 79 85 (Aug) 1943

**Adrenaline**

- Annersten S, Gronwall A and Koiv E. The Fluorimetric Determination of Adrenaline in Blood Plasma. *Scandinav J Clin & Lab Invest* 1 60 69 1940
- Bloor W P and Bullen S S. The Determination of Adrenalin in Blood. *J Biol Chem* 138 (2) 27 739 (April) 1941
- Jorgen en, K. S. Determination of Adrenaline Content of Blood With the Fluorescence Method. *Acta Pharmacol Toxicol* 1 27 739 1941
- Kobro M. Investigations of Adrenalin Concentration in the Blood. I. A Quantitative Chemical Method for Determination of the Concentration of Adrenalin in the Blood. *Acta med. Scandinav* 124 511 521 1946
- Whitehorn J C. A Chemical Method for Estimating Epinephrine in Blood. *J Biol Chem* 103 (3) 633 643 (March) 1935

**Androgens**

- Callow N H, Callow R K and Fimmen C W. Colorimetric Determination of Substances Containing the Grouping CH<sub>2</sub>CO in Urine Extracts as an Indication of Androgen Content. *Biochem J* 32 II 1312 1331 (June) 1938
- Hulmer P E and Hess W C. Spectrophotometric Determination of Androsterone and Testosterone. *Anal Chem* 21 82 83 (July) 1949
- Langstroth G O, Talbot N B, and Fineman A. Spectrochemical Assay of Androsterone and Dehydroisandrosterone in Simple Solutions. *J Biol Chem* 130 (2) 585 591 (Oct) 1939
- Munson P J, Jones M E, McCall P J and Gallagher T F. A Colorimetric Method for Estimations of Dehydroisandrosterone and Its Applications to Urine Extracts. *J Biol Chem* 176 73 82 (Oct) 1948
- Nielsen, A T. Colorimetric Assay of Testosterone. *Acta Endocrinol* 1 367 374 1944

**Cortical Steroids**

- Schneider J J. Studies on the Excretion of Adrenocortical Compounds. I. Isolation of the 17 Hydroxy 11 Dehydrocorticosterone in Normal Males. *J Biol Chem* 183 363 376 (March) 1950

**Corticosteroids**

- Corcoran A C and Page I H. Methods for the Chemical Determination of Corticosteroids in Urine and Plasma. *J Lab & Clin Med* 33 176 183 (Oct) 1948
- Daughaday W H, Jaffe H and Williams F H. Chemical Assay for 'Cortin'. Determination of Formaldehyde Liberated on Oxidation With Periodic Acid. *J Clin Endocrinol* 8 166 174 (Feb) 1948

- Furchgott R F Rovenkrantz H and Shorr E Infra red Absorption Spectra of Steroids IV Adrenal Cortical Hormones and Related Steroids *J Biol Chem* 171 ( ) 523 59 (Dec) 1947
- Talbot N B Saltzman A H Wom R L and Wolfe J K The Colorimetric Assay of Urinary Corticosteroid like Substances *J Biol Chem* 160 ( ) 335 346 (Oct) 1945
- Venning E H Kazmin V E and Bell J C Biological Assay of Adrenal Corticoids *Endocrinol* 38 ( ) 989 (Feb) 1946

### Estrogens

- Cohen H and Bates R W A Simple Quantitative Colorimetric Method for Estrogenic Steroids *J Clin Endocrinol* 7 (10) 701 707 (Oct) 1947
- Dunford R A The Colorimetric Determination of Some Synthetic Oestrogens *Canad J Research* 7 B 646 65 (July) 1949
- Evans J S Varney R F and Koch F C The Mouse Uterine Weight Method for the Assay of Estrogens *Endocrinol* 28 47 752 (May) 1941
- Finkelstein M Microdetermination of Steroid Estrogens in Urine by Fluorometry *Proc Soc Exper Biol & Med* 69 181 184 (Nov Dec) 1945
- Jailer J W A Fluorometric Method for the Determination of Estrogens *Endocrinol* 41 198 201 (Aug) 1947
- Malpress F H The Colorimetric Estimation of Stilboestrol Hexoestrol and Their Glucuronides in Urine *Biochem J* 43 13 16 1948
- Stevenson M F and Marran G F The Determination of Oestrogens in Human Pregnancy Urine A New Method of Correcting for the Brown Colour Developed in the Kober Reaction by Non oestrogenic Substances *Biochem J* 41 (4) 507 511 1947
- Sego C M and Roberts S The Determination of Protein Bound Blood Estrogen *Endocrinol* 41 3 3 4 (Oct) 1947
- Tubis M and Bloom A Colorimetric Determination of Diethylstilbestrol *Ind Eng Chem (Anal Ed)* 14 309 194
- Van Bruggen J T A Comparison of Methods Used for the Hydrolysis of Conjugated Urinary Estrogens *J Lab & Clin Med* 33 215 (Feb) 1949

### Gonadotrophins

- Dekanski J The Kaolin Adsorption Method for the Quantitative Assay of Urinary Gonadotrophins *Brit J Exper Pathol* 30 (4) 37 39 (Aug) 1949
- Haskins A I Jr and Sherman A I Quantitative Bio Assay of Chorionic Gonadotrophin With the Male Frog *Endocrinol* 44 54 545 (June) 1949
- Smith I H Albright F and Dodge E Modifications in Method for the Precipitation and Assay of Increased Amounts of Pituitary Gonadotropic Substances in the Urine *J Lab & Clin Med* 28 1761 1766 (Nov) 1943
- Varney R F and Koch F C A Method for the Assay of the Gonadotropin Content of Normal Human Urine *Endocrinology* 30 399 407 (March) 1942

### Histamine

- Mintire F C Poth L W and Shaw J L The Purification of Histamine for Bioassay *J Biol Chem* 170 (2) 53 544 (Oct) 1947

### Parathormone

- Opińska Blauth J Preparation and Determination of Parathyroid Hormones *Ann Univ Maria Curie Skłodowska Lublin Polonia Sec D* 2 13 38 1947

### Pituitary Adrenocorticotropin

- Fosham P H Thorn G W Prunty F T G and Hills A G Clinical Studies With Pituitary Adrenocorticotropin *J Clin Endocrinol* 8 15 66 (Jan) 1948
- Locke W Albert A and Kepler E J Adrenocortical Activity of Human Urine *Proc Soc Exper Biol & Med* 72 ( ) 470 474 (Nov) 1949
- Sayers M A Sayers G and Woodbury L A The Assay of Adrenocorticotrophic Hormone by the Adrenal Isotopic Acid Depletion Method *Endocrinol* 42 379 383 (May) 1948

### Protein Bound Iodine

- Ba sett A M Coons A H Salter W T and Simmons S M Protein Bound Iodine in Blood V Naturally Occurring Iodine Fractions and Their Chemical Behavior *Am J Med Sci* 202 516 5 1941

### Steroids

- Iarnett J Henly A A and Moss C J O R The Polarographic Estimation of Steroid Hormones I Polarography of Neutral 17 Ketosteroids in Urinary Extracts *Biochem J* 40 445 449 1946

- Clark L C Jr and Thompson H A New Series of Reagents for the Colorimetric Determination of Steroids *Science* 107 499 431 (April) 1948
- Dobriner K Lieberman S and Rhoads C P Studies in Steroid Metabolism I. Methods for the Isolation and Qualitative Estimation of Neutral Steroids Present in Human Urine *J Biol Chem* 172 241 261 (Jan) 1948
- Drekter I J Pearson S, Bartczak E and McGavack T H A Rapid Method for the Determination of Total Urinary 17 Ketosteroids *J Clin Endocrinol* 7 (1) 795 800 (Dec) 1947
- Engel L L Patterson H R Wilson H and Schinkel M The Colorimetric Estimation of Urinary Neutral Steroid Alcohols *J Biol Chem* 183 47 72 (March) 1950
- Fraser R W, Forbes A P Albright F Sulkowitch H, and Reifensstein E Jr Colorimetric Assay of 17 Ketosteroids in Urine *J Clin Endocrinol* 1 234 36 (March) 1941
- Heard R D H and Sobel H VIII A Colorimetric Method for the Estimation of Reducing Steroids *J Biol Chem* 165 (2) 687 698 (Oct) 1946
- Holtorf, A F and Koch F C The Colorimetric Determination of 17 Ketosteroids and Their Application to Urine Extracts *J Biol Chem* 135 377 392 (Sept) 1940
- Lieberman S Dobriner K, Hill B R Fieser L F and Rhoads C P Studies in Steroid Metabolism II Identification and Characterization of Ketosteroids Isolated From Urine of Healthy and Diseased Persons *J Biol Chem* 172 263 295 (Jan) 1948
- Talbot N B Butler A M and McLachlan E The Colorimetric Assay of Total  $\alpha$  and  $\beta$  17 Ketosteroids in Extracts of Human Urine *J Biol Chem* 132 (2) 595 603 (Feb) 1940
- Tompsett S L and Oastler E G The Excretion of Corticosteroids The Determination of the Total Free Reducing Ketosteroids in Urine *Glasgow M J* 28 349 365 (Nov) 1947
- Tompsett S L The Determination of Steroids in Urine *The Analyst* 74 6 17 (Jan.) 1949
- Wilson H and Carter P Stabilization of the Alcoholic Potassium Hydroxide in Colorimetric 17 Ketosteroids Determinations *Endocrinol* 41 (5) 417 421 (Nov) 1947
- Thyroxine**
- Pitt Rivers R and Lerman J Physiological Activity of the Optically Active Isomers of Thyroxine *J Endocrinol* 5 (5) 223 228 (June) 1948

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### G ENZYMES

##### 1 Lipase in Serum

###### Reference

Johnson T A and Bockus H L Diagnostic Significance of Determinations of Serum Lipase Arch Int Med 66 6 7 (July) 1940

###### Principle

This method depends on the estimation of fatty acids produced by the hydrolysis of an olive oil emulsion on which lipase in the blood serum has been acting during a 4 hour period under optimum conditions of pH and temperature

###### Apparatus

- 1 A water bath set at 100 C
- 2 An incubator set at 37.5 C
- 3 Accurately calibrated 0.5 ml and 1 ml syringe pipette
- 4 A 10 ml burette
- 5 A Waring Blendor

###### Reagents

- 1 Phenolphthalein solution 1% in alcohol
- 2 Disodium hydrogen phosphate ( $\text{Na}_2\text{HPO}_4 \cdot \text{H}_2\text{O}$ ) Dissolve 59.35 gm in 1000 ml of distilled water to make an M/3 solution
- 3 Potassium dihydrogen phosphate ( $\text{KH}_2\text{PO}_4$ ) Dissolve 45.38 gm in 1000 ml of distilled water to make M/3 solution
- 4 Mixed buffer solution Mix 10 parts of  $\text{Na}_2\text{HPO}_4 \cdot \text{H}_2\text{O}$  with 3 parts of  $\text{KH}_2\text{PO}_4$
- 5 Olive oil emulsion a 50% solution Prepare by using a 5% solution of acacia as the emulsifying agent Add 0.5 percent of sodium benzoate as a preservative Mix in Waring Blendor
- 6 Ethyl alcohol 95%
- 7 Sodium hydroxide 0.05 N See under volumetric analysis for making up 1 N sodium hydroxide and dilute exactly 20 times

###### Procedure

- 1 Pipette 3 ml of distilled water into each of two 15 ml centrifuge tubes A and B
- 2 Add exactly 1 ml serum to each and shake
- 3 Place tube A (control) in a 100 C water bath for 5 minutes to inactivate the lipase Be sure to cool before proceeding
- 4 Add 0.5 ml of buffer solution to each tube
- 5 Then add 0.5 ml of olive oil emulsion to each and shake
- 6 Incubate at 37.5 C for 4 hours
- 7 At the end of that time add 3 ml of 95% alcohol and 2 drops of phenolphthalein indicator
- 8 Titrate each tube to a faint pink color with 0.05 N sodium hydroxide

3 **Cholinesterase in Serum and Cells****Reference**

Michel H O Electrometric Method for the Determination of Red Blood Cell and Plasma Cholinesterase Activity, *Fed Proc* 8 99 (March) 1949

**Principle**

Cholinesterase hydrolyses acetylcholine bromide with the production of acetic acid. The resultant drop in pH is measured in a barbital phosphate buffer when serum is mixed with a known excess of acetylcholine bromide. Results are expressed as decrease in pH units per hour.

**Apparatus**

- 1 Constant temperature bath at 25 C
- 2 Beckman pH meter Model G
- 3 A number of 15 ml thin wall glass test tubes
- 4 Accurately calibrated 1.0 and 10.0 ml syringe pipettes
- 5 0.1 ml volumetric serology pipettes
- 6 Small glass beakers 20 ml capacity

**Reagents**1 **Preparation of the Stock Buffer**

A stock buffer, which will keep approximately 1 month in a glass stoppered bottle at room temperature and which contains 0.1 M sodium diethylbarbiturate (Barbital Sodium Merck) and 0.0167 M  $\text{KH}_2\text{PO}_4$  in distilled water is prepared as follows:

a 90.618 gm of barbital sodium and 2.9092 gm of  $\text{KH}_2\text{PO}_4$  are weighed accurately and placed in a clean one liter beaker. About 50 ml of distilled water are added. This solution will have a pH of about 9.9. Dilute acid (approximately 0.1 N HCl, 1 ml of concentrated HCl to 1.0 ml of water) is added slowly and with constant stirring to bring the pH of the buffer down to between 8.2 and 8.4. If strong acid is used, if addition of the acid is rapid or if stirring is inadequate white crystals of barbital may appear. These may be redissolved by heating and they will remain in solution if the solution is cooled slowly. The buffer is then transferred with washing into a volumetric flask and brought to volume with distilled water.

2 **Preparation of Working Buffer**

A working buffer must be made up fresh each day. It is made from the stock buffer with a salt solution containing 1.5 N NaCl and 0.30 N  $\text{MgSO}_4$  in distilled water (87.68 gm NaCl and 18.06 gm of  $\text{MgSO}_4$  dissolved in water up to one liter) in the proportions: 3 ml stock buffer, 10 ml salt solution and 87 ml distilled water. This gives 100 ml of working buffer. In the preparation of the working buffer exactly 1 ml of stock buffer and 10 ml of salt solution are mixed with about 80 ml of distilled water. The pH is then adjusted to exactly 8.00 by addition of approximately 0.1 N HCl if its pH is over 8.00 or approximately 0.1 N NaOH (4 gm NaOH per liter) if its pH is under 8.00. The buffer is then transferred with washing into a volumetric flask and brought to 100 ml with distilled water.

3 **Preparation of Acetylcholine Solution**

A solution of acetylcholine bromide is used 3.73% in distilled as obtained commercially is impure and must be purified. Commercially obtained acetylcholine is dissolved in a glass stoppered flask. This is cooled in ice water. Cold ether is then added until turbid. The mixture is placed in a refrigerator for 48 hours. The mixture is then filtered by filtration through a fine filter. The filtrate is dried over  $\text{NaOH}$  and the solution will keep.

The 3.73%

The material is purified by dissolving in 95% alcohol for 2 hours, placed in a separated funnel and dried over  $\text{NaOH}$  and

## Procedure

- 1 Pipette 10 ml of working buffer into a thin glass test tube and place in a water bath at 5 C for at least 10 minutes
- 2 Add 0.1 ml of unknown serum (or plasma) and 10 ml of acetylcholine solution to the buffer
- 3 The exact time of addition of acetylcholine solution is noted. Mixing with a glass stirring rod is performed on the addition of the acetylcholine solution.
- 4 The test tube is stoppered (rubber) and placed back in the 5 C water bath
- 5 Readings are conveniently taken at exactly 1 hour but if this is impracticable they may be taken at any time after an hour and the decrease in pH corrected to one hour. For sera with very low cholinesterase values incubation for longer periods than one hour (e.g. up to 5 hours) may be required to give reliable readings. The reading is noted and the corrected value for pH drop per hour determined from the accompanying table (Table 11) which makes correction for non enzymatic hydrolysis of acetylcholine and change in rate of pH drop with varying pH's
- 6 It should be noted that all glassware must be scrupulously clean as small changes in pH may be caused by dirt and this will affect readings. It is suggested that all readings be done in duplicate as a check and that with every set of unknown sera on which determinations are being run a serum whose cholinesterase activity is known be run simultaneously as a check on the method

## Calculation

- 1 Two significant corrections have to be made in the computation of units  $\Delta$  pH/hr
  - a. Non enzymatic hydrolysis of the substrate
  - b. The effect of pH on the specific activity of cholinesterase
- 2 The complete calculation is corrected  $\Delta$  pH units/hour =
 
$$\frac{(\text{8.00} - \text{final pH} + \text{blank} \pm \text{specific effect of pH})}{(\text{Total hours})}$$
- 3 Table 11 incorporates the two necessary corrections and in addition the subtraction from pH 8.00

TABLE 11

MICHEL'S METHOD FOR ESTIMATING CHOLINESTERASE ACTIVITY CORRECTED VALUES FOR pH DECREMENT PER HOUR FOR pH 6.3 TO pH 7.99

(This table incorporates into the observed  $\Delta$  pH corrections for nonenzymatic hydrolysis of acetylcholine bromide and for effect of pH on specific activity and it subtracts the corrected pH from 8.00)

| pH Units<br>as Read | 0.00                        | 0.0  | 0.04 | 0.06 | 0.08 |
|---------------------|-----------------------------|------|------|------|------|
|                     | Corrected $\Delta$ pH Units |      |      |      |      |
| 6                   | 1.74                        | 1.7  | 1.70 | 1.67 | 1.64 |
| 6.3                 | 1.63                        | 1.60 | 1.58 | 1.56 | 1.54 |
| 6.4                 | 1.5                         | 1.50 | 1.48 | 1.46 | 1.44 |
| 6.5                 | 1.42                        | 1.40 | 1.38 | 1.36 | 1.34 |
| 6.6                 | 1.3                         | 1.30 | 1.29 | 1.27 | 1.26 |
| 6.7                 | 1.2                         | 1.2  | 1.2  | 1.19 | 1.17 |
| 6.8                 | 1.12                        | 1.12 | 1.11 | 1.10 | 1.08 |
| 6.9                 | 1.06                        | 1.04 | 1.03 | 1.01 | 1.00 |
| 7.0                 | 0.98                        | 0.96 | 0.94 | 0.93 | 0.91 |
| 7.1                 | 0.89                        | 0.87 | 0.85 | 0.83 | 0.81 |
| 7.2                 | 0.8                         | 0.77 | 0.75 | 0.73 | 0.71 |
| 7.3                 | 0.69                        | 0.67 | 0.65 | 0.63 | 0.61 |
| 7.4                 | 0.59                        | 0.56 | 0.54 | 0.52 | 0.50 |
| 7.5                 | 0.47                        | 0.45 | 0.43 | 0.40 | 0.38 |
| 6                   | 0.36                        | 0.34 | 0.32 | 0.30 | 0.27 |
| 7.7                 | 0.2                         | 0.2  | 0.20 | 0.18 | 0.16 |
| 7.8                 | 0.14                        | 0.12 | 0.10 | 0.08 | 0.06 |
| 7.9                 | 0.04                        | 0.03 | 0.0  | 0.01 | 0.00 |

### Chymotrypsin

Ishii B M Huang H T, and Niemann C A Colorimetric Method for the Determination of Chymotrypsin Activity, *J Biol Chem* 183 403 408 (April) 1950

### Cytochrome Oxidase

Smith F, and Stotz, E A Colorimetric Method for the Determination of Cytochrome Oxidase, *J Biol Chem* 179 891 902 (June) 1949

### Hyaluronidase

Kass E H, and Seastone C V Role of the Mucoid Polysaccharide (Hyaluronic Acid) in the Virulence of Group A Hemolytic Streptococci, *J Exper Med* 79 319 330 (March) 1944

Tolksdorf S McCready M H McCullagh, D R and Schwenk E The Turbidimetric Assay of Hyaluronidase *J Lab & Clin Med* 34 74 89 (Jan) 1949

Warren G H Durso J G and Levin N R A Modified Turbidimetric Method for the Assay of Hyaluronidase *Endocrinology* 43 48 51 (July) 1948

Dorfman A and Ott M L A Turbidimetric Method for the Assay of Hyaluronidase *J Biol Chem* 172 367 370 (Feb) 1948

### Lipase

Seligman A L and Nachlas M M The Colorimetric Determination of Lipase and Esterase in Human Serum *J Clin Investigation* 29 (1) 31 36 (Jan) 1950

### Lipoxidase

Sumner J B and Smith G N The Estimation of Lipoxidase Activity *Arch Biochem* 14 87 92 (July) 1947

### Peptic Activity

Fodor P J Price V E and Greenstein J P Enzymatic Hydrolysis of Saturated and Unsaturated Tripeptides *J Biol Chem* 180 (1) 193 208 (Aug) 1949

Tomarelli, R M Charney J and Harding M L The Use of Azoalbumin as a Substrate in the Colorimetric Determination of Peptic and Tryptic Activity, *J Lab & Clin Med* 34 478 483 (March) 1949

Wesselman, H J and Hiltz W W The Turbidimetric Determination of Peptic Activity, *J Am Pharmaceut Assoc Sci Ed* 37 357 359 (Sept) 1948

### Peroxidase

Eitner, J The Estimation of Peroxidase Activity *Biochem J* 44 30 38 (Jan) 1949

### Phosphatase

Andersch M A and Szczypinski A J Use of p Nitrophenylphosphate as the Substrate in Determination of Serum Acid Phosphatase *Am J Clin Path* 17 (7) 571 574 (July) 1947

Hudson P B Brendler H and Scott W W A Simple Method for the Determination of Serum Acid Phosphatase *J Urol* 58 (1) 89 92 (July) 1947

King E J and Armstrong A R A Convenient Method for Determining Serum and Bile Phosphatase Activity *Canad M A J* 31 (4) 376 381 (Oct) 1934

Sherlock S P V and Walshe V Hepatic Alkaline Phosphatase Histological and Micro Chemical Studies on Liver Tissue in Normal Subjects and in Liver and Bone Disease *J Path & Bact* 59 (4) 610 630 (Oct) 1947

### Proteolytic Activity

Ishii B M and Niemann C A Procedure for the Determination of Proteolytic Activity *J Biol Chem* 182 821 8 9 (Feb) 1950

Schwartz T B and Engel F L A Photometric Ninhydrin Method for the Measurement of Proteolysis *J Biol Chem* 184 197 200 (May) 1950

### Trypsin

Lineweaver H Fraenkel Conrat H and Bean, B S Determination of Trypsin in the Presence of Egg White Trypsin Inhibitor and Demonstration of Absence of Trypsin From Egg White *J Biol Chem* 177 205 207 (Jan) 1949

Spencer W H Gastrointestinal Studies VIII A Method for the Quantitative Estimation of Trypsin in Gastric Contents *J Biol Chem* 21 165 167 (June) 1915

White F D and Bowman J M The Estimation of Trypsin *Canad J Research* 25 E 153 161 (Aug) 1947

## SECTION IV

### BIOCHEMICAL PROCEDURES (Continued)

---

#### H PIGMENTS

##### 1 Hemoglobin (Sanford Sheard)

###### Reference

Sanford A H and Sheard C The Determination of Hemoglobin With the Photoelectrometer J Lab & Clin Med 15 483 489 (Feb) 1930

###### Principle

Blood is diluted with enough water to insure that all of the hemoglobin is converted to oxyhemoglobin. Complete solution is assisted by the addition of a small quantity of sodium carbonate.

###### Apparatus

- 1 A Coleman Jr Spectrophotometer
- Cuvettes 10 x 150 mm
- 3 Micro pipettes 0.1 ml accurately calibrated
- 4 A quantity of 50 ml volumetric flasks

###### Reagents

- 1 Sodium carbonate a 5% solution in water

###### Procedure

- 1 Into a 50 ml volumetric flask containing approximately 5 ml of distilled water pipette exactly 0.1 ml of blood & using the pipette out at least 6 times. Add 0.2 ml carbonate dilute to 50 ml & transfer to cuvette set water blank at 100% T at 550 mμ and read unknowns.

###### Calibration Curve

Determine the exact oxygen capacity of blood specimens equilibrated 0.30 minutes with room air using the Van Slyke technique for analysis. The values in volumes percent oxygen are first corrected for dissolved oxygen and then divided by a factor of 1.34 for conversion to grams of hemoglobin. Set up dilutions of the blood using the Sanford Sheard method in the range 5 to 0% gm hemoglobin/100 ml blood.

###### Calculation

$$\text{gm hemoglobin/100 ml blood} = \frac{(\text{mg Hb/ml in cuvette}) \times (\text{total dilution}) \times 100}{1000 \times \text{ml blood}}$$

###### Example

ml blood 0.1    Total dilution 50 ml  
mg Hb/ml in cuvette 0.30

$$\text{gm hemoglobin/100 ml blood} = \frac{(0.30) \times 50 \times 100}{1000 \times 0.1} = 15.0 \text{ gm/100 ml}$$



**Precautions**

- 1 The color is stable for at least 30 minutes at usual room temperatures, but is unstable in the heat
- 2 This method is good for clinical purposes. For more precision and for fractionating the heme pigment fractions the method of Evelyn and Malloy is recommended.

**Normal Hemoglobin Standards**

Various authors have proposed different values for 'normal' or 100% hemoglobin. For instance Haldane used 13.8 gm/100 ml and Hales 15.6 gm/100 ml. The following distribution curve was derived from analyses of the oxygen capacity of 791 healthy males

| OXYGEN CAPACITY<br>(ml/100 ml BLOOD) | NUMBER OF SUBJECTS<br>(% OF TOTAL) |
|--------------------------------------|------------------------------------|
| From 15.6 to 17.0                    | 1.2                                |
| From 17.1 to 18.0                    | 3.9                                |
| From 18.1 to 19.0                    | 15.4                               |
| From 19.1 to 20.0                    | 5.3                                |
| From 20.1 to 21.0                    | 30.6                               |
| From 21.1 to 22.0                    | 16.2                               |
| From 22.1 to 23.0                    | 4.6                                |
| From 23.1 to 23.6                    | 0.8                                |

The average of all analyses was 20.11 ml oxygen/100 ml or 15.0 gm Hb/100 ml

**EXAMPLE OF HEMOGLOBIN CALIBRATION CHART**

Direct conversion gm Hb/100 ml blood when 0.1 ml of blood is diluted to 50 ml

| °C | 0                  | 1    | 2    | 3    | 4    | 5    | 6    | 7    | 8    | 9    |
|----|--------------------|------|------|------|------|------|------|------|------|------|
|    | gm Hb/100 ml BLOOD |      |      |      |      |      |      |      |      |      |
| 40 | 19.6               | 19.1 | 18.6 | 18.1 | 17.5 | 17.0 | 16.6 | 16.1 | 15.7 | 15.2 |
| 50 | 14.8               | 14.4 | 14.1 | 13.6 | 13.2 | 12.8 | 12.5 | 12.0 | 11.6 | 11.3 |
| 60 | 10.9               | 10.6 | 10.1 | 9.8  | 9.5  | 9.2  | 8.8  | 8.4  | 8.0  | 7.9  |
| 70 | 7.6                | 7.3  | 6.9  | 6.6  | 6.4  | 6.1  | 5.8  | 5.5  | 5.0  | 4.9  |
| 80 | 4.7                | 4.5  | 4.3  | 3.9  | 3.6  | 3.4  | 3.3  | 3.1  | 2.9  | 2.7  |

**2 Total Hemoglobin (Evelyn Malloy)****Reference**

Evelyn, R. A. and Malloy H. T. (Modified) Microdetermination of Oxyhemoglobin, Methemoglobin and Sulfhemoglobin in a Single Sample of Blood, *J Biol Chem* 126: 655-662 (Dec) 1938

**Principle**

Total hemoglobin is determined by converting all the hemoglobin present to cyanmethemoglobin and measuring the light absorption at 440 mμ

**Apparatus**

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Micro pipettes 0.1 ml accurately calibrated
- 3 A variety of 50 ml volumetric flasks
- 4 Cuvettes 19 x 150 mm

**Reagents**

- 1 Phosphate buffer (pH 6.6) 0.1 M (See section on buffers)
- 2 Potassium ferricyanide a 1.0% solution in water
- 3 Sodium cyanide or potassium cyanide a 5% solution in water **POISON!**
- 4 Concentrated ammonium hydroxide

### Procedure

- 1 Introduce exactly 0.1 ml of whole blood into a 50 ml volumetric flask containing approximately 25 ml of distilled water. Rinse pipette at least 6 times.
- 2 Add 10 ml of 1/10 M phosphate buffer pH 6.6 and mix.
- 3 Add 3 drops of 20% potassium ferricyanide mix and let stand for 10 minutes.
- 4 Add 3 drops of 5% NaCN or KCN mix and let stand 2 minutes.
- 5 Add 1 or 2 drops of concentrated NH<sub>4</sub>OH and mix.
- 6 Dilute to 50 ml with distilled water.
- 7 Prepare a blank tube in same manner as above being careful to add exactly the same quantity of potassium ferricyanide.
- 8 Set blank at 100% T at 540 mμ and read the unknowns.

### Calibration Curve

- 1 Determine the oxygen capacity of samples of blood as described under the method of Sanford and Sheard.  
Run appropriate dilutions through the Evelyn Malloy procedure in the range 5 to 25 g n hemoglobin/100 ml blood.

### Calculation

$$\text{gm hemoglobin/100 ml blood} = \frac{(\text{mg/ml in cuvette}) \times (\text{dilution}) \times (100)}{(1000) \times (\text{ml blood sample})}$$

### Example

62.5% T was equivalent to 0.30 mg Hb/ml in cuvette  
For analysis 0.1 ml blood was finally diluted to 50 ml

$$\text{gm hemoglobin/100 ml blood} = \frac{(0.30) \times (50) \times (100)}{(1000) \times (0.1)} = 15.0 \text{ gm/100 ml}$$

### Precautions

- 1 Be sure to wash out the micro pipette at least 6 times in the water.
- 2 The color is stable for at least 30 minutes.

## Oxyhemoglobin, Methemoglobin, Sulfhemoglobin

### Reference

Evelyn K. A. and Malloy H. T. Microdetermination of Oxyhemoglobin, Methemoglobin and Sulfhemoglobin in a Single Sample of Blood. *J. Biol. Chem.* 126: 655-666 (Dec.) 1938.

### Principle

The methemoglobin is measured photometrically after being converted to cyanmethemoglobin. The oxyhemoglobin is converted first to methemoglobin with ferricyanide and then to cyanmethemoglobin. Sulfhemoglobin is estimated by direct observation at 60 mμ with a correction for the absorption by oxyhemoglobin and cyanmethemoglobin at this wavelength.

### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 × 150 mm
- 3 Micro pipettes 0.1 ml accurately calibrated
- 4 Dropping bottles

### Reagents

- 1 Phosphate buffer M/15 pH 6.6 (See Sørensen's tables in the section on buffers).  
Phosphate buffer M/60 diluted from above when needed.
- 2 Potassium ferricyanide, K<sub>3</sub>Fe(CN)<sub>6</sub> 20%.

- 4 Sodium cyanide (NaCN), a 10% solution in water POISON!
- 5 Acetic acid (CH<sub>3</sub>COOH) a 12% solution in water
- 6 Solution of 10% NaCN mixed within ONE HOUR of use with an equal volume of 12% acetic acid.
- 7 Ammonium hydroxide (NH<sub>4</sub>OH) concentrated

**Procedure for Methemoglobin**

- 1 Pipette exactly 0.1 ml of fresh whole blood into 9.9 ml of M/60 phosphate buffer in cuvette. Wash out the pipette at least 6 times
- 2 Let stand for 2 minutes
- 3 Make up a blank with 10 ml of M/60 buffer and set at 0 density at a wave length of 635 mμ.
- 4 Read the unknown sample (Reading D)
- 5 Add 1 drop of neutralized NaCN to convert MHb to MHbCN
- 6 After 2 minutes read again at 635 mμ with the same blank. This is D<sub>2</sub>. The difference (D - D<sub>2</sub>) is proportional to the concentration of MHb

**Procedure for Sulfhemoglobin**

- 1 Clear the solution for SHb determination by addition of 1 drop of concentrated ammonium hydroxide
- 2 Read in colorimeter at 690 mμ with water blank at 0 density. This reading is D

**Procedure for Total Hemoglobin**

- 1 Pipette 2 ml of this solution into a second cuvette
- 2 Add 8 ml of M/15 phosphate buffer of pH 6.6 and 0.05 ml of 20% potassium ferricyanide
- 3 Let stand for 2 minutes to allow HbO to be converted to MHb
- 4 Add 1 drop of NaCN reagent to convert MHb to MHbCN for total hemoglobin
- 5 Read D<sub>2</sub> after 2 minutes at 540 mμ, with a blank consisting of 10 ml water 1 drop cyanide reagent and 0.05 ml 20% potassium ferricyanide

**Calibration Curves**

- 1 Estimate the total hemoglobin of samples of blood by means of the oxygen capacity as described under the Sanford Sheard method for hemoglobin
- 2 Following the procedure described above run samples through the method for total hemoglobin and compute K<sub>th</sub>, the calibration constant for total hemoglobin
- 3 At the same time convert HbO to MHb as under total hemoglobin but run through as for methemoglobin and compute the calibration constant K<sub>mh</sub>
- 4 The calibration constant for sulfhemoglobin K<sub>sh</sub> is assumed to be 10 unless a sample of pure sulfhemoglobin can be prepared

**Calculation**

- 1  $T = \text{gm total hemoglobin}/100 \text{ ml blood} = \frac{100 \times D}{K_{th}}$
- 2  $M = \text{gm methemoglobin}/100 \text{ ml blood} = \frac{100 (D - D_2)}{K_{mh}}$
- 3  $S = \text{gm sulfhemoglobin}/100 \text{ ml blood} = \frac{100 \times D - (8.5 \times M) - (4.4 \times T)}{10}$
- 4 Corrected value of T = T from equation 1 - 0.2 × S
- 5 HbO = Corrected T - M - S

**Precautions**

- 1 The least accurate part of this method is that for sulfhemoglobin. Fortunately this compound is rarely present in amounts over 0.5 gm/100 ml blood.

#### 4 Pigments Miscellaneous References

##### Carboxyhemoglobin

Klenshøj, N. C., Feldstein, M. and Sprague, A. L. The Spectrophotometric Determination of Carbon Monoxide. *J Biol Chem* 183: 303 (March) 1950

##### Chlorophyll

Goodwin, R. H. Fluorometric Method of Estimating Small Amounts of Chlorophyll a. *Anal Chem* 19: 789 (Oct) 1947

Petering, H. G., Benne, E. J. and Morgal, P. W. Simplification of the Petering-Wolman-Hibbard Method for Determination of Chlorophyll and Carotene. *Ind Eng Chem (Anal Ed)* 13: 236 (April) 1941

##### Chromoproteins

Drabkin, D. L. The Distribution of the Chromoproteins Hb, Myoglobin and Cytochrome C in the Tissues of Different Species and the Relationship of the Total Content of Each Chromoprotein to Body Mass. *J Biol Chem* 182: 317-34 (Jan) 1950

##### Cytochrome C

Prader, A. and Gonella, A. The Spectrophotometric Determination of Cytochrome C in Tissues. *Experientia* 3: 46-464 (1947)

Rosenthal, O. and Drabkin, D. L. Spectrophotometric Studies. XI. The Direct Micro Spectrophotometric Determination of Cytochrome C. *J Biol Chem* 149: 437-450 (Aug) 1943

##### Hemoglobin

Holden, H. F. On the Colorimetric Determination of Haemoglobin. *Australian J Exper Biol & Med Sci* 25 (1): 57-60 (March) 1947

Pennell, R. B., Smith, W. E. and Werkheiser, W. C. Preparation of Hemoglobin Solutions Containing Hemoglobin Reducing Enzymes. *Proc Soc Exper Biol & Med* 65 (2): 35-36 (June) 1947

Peters, J. T. A New Micro-iron Method for the Determination of Hemoglobin in Blood From One Drop of Fingertip Blood. *South M J* 40 (11): 949-6 (Nov) 1947

##### Porphyrins

Grinstein, M. and Wintrobe, M. M. Spectrophotometric Micromethod for the Quantitative Determination of the Free Erythrocyte Protoporphyrin. *J Biol Chem* 172 (2): 459-467 (Feb) 1948

Nicholas, R. E. H. and Rimington, C. Qualitative Analysis of Porphyrins by Partition Chromatography. *Scandinav J Clin & Lab Investigation* 1: 1-18 (1949)

Sveinsson, S. L., Rimington, C. and Barnes, H. D. Complete Porphyrin Analysis of Pathological Urines. *Scandinav J Clin & Lab Investigation* 1: 1-11 (1949)

Watson, C. J., Schwartz, S. and Hawkinson, V. Studies of the Uroporphyrins. II. Further Studies of the Porphyrins of the Urine, Feces, Bile and Liver in Cases of Porphyrin. With Particular Reference to a Waldenström Type Porphyrin Behaving as an Entity on the Tswett Column. *J Biol Chem* 157: 345-361 (Jan) 1945

# BIOCHEMICAL PROCEDURES (Continued)

## SECTION IV

---

### I HYDROGEN ION CONCENTRATION

#### 1 Indicators

Indicators are weak organic acids or bases whose salts are different in color from the acid. They are buffers with buffer ranges characteristic for each indicator. As the ratio of acid to salt form of the indicator changes during a titration, a gradual change in color ensues, and the range of pH over which this change occurs determines what indicator may be used for any particular titration.

When a strong mineral acid is titrated with a strong base, almost any indicator may be used since the pH changes from acid to alkaline with a very small amount of alkali at the end point.

When a weak organic acid such as acetic is titrated with a strong base, the indicator must be one that changes color at alkaline pH. This is because the weak organic acids are buffers and in order to titrate them completely, one must titrate well beyond their buffer range. Phenolphthalein is often used for this kind of titration, but cannot be used in the presence of ammonium salts or of carbonic acid.

When a weak base is titrated with a strong acid, one must use an indicator whose color change is in the acid range. Bromphenol blue is useful for these titrations.

LaMotte, American Aniline Company, Eastman and Elmer and Amend supply either the dry indicator from which solutions can be prepared or standard solutions of indicators prepared according to Clark. (See Table 13.)

A useful universal indicator is prepared and used as follows:

5 mg thymol blue  
25 mg methyl red  
60 mg bromthymol blue  
60 mg phenolphthalein

Dissolve and make up to 100 ml of 5% ethyl alcohol. Titrate with 0.01 N NaOH to a green color.

The range of this solution is from pH 4 to pH 10. The color changes are:

|        |      |        |       |
|--------|------|--------|-------|
| Red    | pH 4 | Green  | pH 7  |
| Orange | pH 5 | Blue   | pH 8  |
| Yellow | pH 6 | Indigo | pH 9  |
|        |      | Violet | pH 10 |

#### 2 Buffer Solutions

##### Reference

Clark, W. Mansfield. *The Determination of Hydrogen Ions*, ed. 3. Baltimore 1928.  
Williams and Wilkins Company.

TABLE 13  
INDICATORS SELECTED BY CLARK AND LUBS (1917)  
SUPPLEMENTED BY COHEN (1917)

| COMMON NAME        | MOL WT | A   | pK   | RANGE | COLOR CHANGE |        | B        |    | C       |           | ABSORPTION MAXIMUM |           |
|--------------------|--------|-----|------|-------|--------------|--------|----------|----|---------|-----------|--------------------|-----------|
|                    |        |     |      |       | ACID         | ALK.   | pH       | pH | ACID mM | ALKALI mM | ACID mM            | ALKALI mM |
| Metacresol purple  | 38     | 6   | 1.1  | 12.28 | red          | yellow | conc HCl | 6  | 533     |           |                    |           |
| Thymol blue        | 466    | 215 | 1    | 1° 8  | red          | yellow | conc HCl | 6  | 541     |           |                    |           |
| Bromphenol blue    | 699    | 149 | 3.98 | 30.46 | yellow       | blue   | 0        | 7  |         |           | 50                 |           |
| Bromocresol green  | 638    | 143 | 4.67 | 38.54 | yellow       | blue   | 1        | 8  |         |           | 617                |           |
| Chlorophenol       | 43     | 36  | 5.98 | 46.64 | yellow       | red    |          | 9  |         |           | 573                |           |
| Bromphenol red     | 512    | 195 | 6.16 | 5.68  | yellow       | red    | 3        | 10 |         |           | 574                |           |
| Bromocresol purple | 540    | 185 | 6.3  | 5.68  | yellow       | purple | 3        | 10 |         |           | 591                |           |
| Bromthymol blue    | 64     | 160 | 7.0  | 5.70  | yellow       | blue   | 4        | 10 |         |           | 617                |           |
| Phenol red         | 354    | 982 | 7.9  | 68.84 | yellow       | red    | 5        | 11 |         |           | 58                 |           |
| Cresol red         | 38     | 62  | 8.3  | 72.88 | yellow       | red    | 5        | 11 |         |           | 7                  |           |
| Metacresol purple  | 38     | 26° | 8.3  | 74.90 | yellow       | purple | 5        | 11 |         |           | 580                |           |
| Thymol blue        | 466    | 215 | 8.9  | 80.96 | yellow       | blue   | 6        | 1  |         |           | 596                |           |
| Cresol phthalein   |        |     | 9.4  | 82.98 | colorless    | red    | 6        | 1° |         |           |                    |           |

A = ml of 0.01 N NaOH required per 0.10 gm indicator to form mono sodium salt

Dilute to 250 ml for 0.04% solution

B = approx. pH value of solution required for full acid color

C = approx. pH value of solution required for full alkaline color

Clark W M The Determination of Hydrogen Ions ed 3 Baltimore 1928 Williams & Wilkins Company

#### Definitions

- 1 A buffer is any substance or combination of substances which when dissolved in water produces a solution which resists a change in its hydrogen ion concentration upon the addition of an acid or an alkali
- 2 pH = - log hydrogen ion concentration expressed in moles/liter (the pH of 0.01 N HCl in 0.09 N KCl at 20°C is taken as 2.038)
- 3 pK = pH at which a buffer is half ionized and half unionized.

#### Buffer Solutions

The various portions of Tables 14A to 14E give the directions for preparing buffer solutions of known hydrogen ion concentration as recommended by different authors

### 3 Glass Electrodes

(Measurement of pH of Blood Plasma and Red Cell Under Anaerobic Conditions)

#### References

- Dill D B Daly C A and Forbes W H pH of Serum and Red Cells J Biol Chem. 117 569-579 (Feb) 1937
- Clark W M The Determination of Hydrogen Ions ed 3 Baltimore 1928 Williams and Wilkins Company p 10

#### Principle

Glass acts like a fluid in which hydrogen ions dissolve. It is therefore possible to measure the boundary potentials when fluids of different hydrogen ion concentrations are on either side of a thin glass membrane. In the present method fluids may be handled anaerobically within a thin glass tube which acts as the glass electrode.

#### Apparatus

- 1 Glass electrode with 0.10 N HCl AgCl half cell and water jacket for temperature control (See Fig 24). The electrode proper of Corning Glass No 015 is a

TABLE 14A

SØRENSEN'S PHOSPHATE MIXTURE FOR pH 5.8-8.2  
(*Biochem. Z.* 21:131 1909 and 22:352 1909)

| 1 M/15 Primary phosphate 9.08 gm $\text{KH}_2\text{PO}_4$ /liter           |  |                                       |                    |
|--|--|---------------------------------------|--------------------|
| 2 M/15 Secondary phosphate 9.41 gm $\text{Na}_2\text{HPO}_4$ (anhyd)/liter |  |                                       |                    |
| pH at 20 C   | M/15 $\text{Na}_2\text{HPO}_4$<br>(ml) | M/15 $\text{KH}_2\text{PO}_4$<br>(ml) | Phosphate<br>mEq/l |
| 5.8  | 8.0                                    | 92.0                                  | 108.0              |
| 5.9  | 9.0                                    | 90.1                                  | 108.2              |
| 6.0  | 12.2                                   | 87.8                                  | 111                |
| 6.1  | 15.3                                   | 84.7                                  | 115.3              |
| 6.2  | 18.6                                   | 81.4                                  | 118.6              |
| 6.3  | 22.4                                   | 77.6                                  | 122.4              |
| 6.4  | 26.7                                   | 73.3                                  | 126.7              |
| 6.5  | 31.8                                   | 68.2                                  | 131.8              |
| 6.6  | 37.5                                   | 62.5                                  | 137.5              |
| 6.7  | 43.5                                   | 56.5                                  | 143.5              |
| 6.809  | 49.6                                   | 50.4                                  | 149.6              |
| 6.862  | 52.5                                   | 47.5                                  | 152.5              |
| 6.909  | 55.4                                   | 44.6                                  | 155.4              |
| 6.958  | 58.2                                   | 41.8                                  | 158.2              |
| 7.005  | 61.1                                   | 38.9                                  | 161.1              |
| 7.052  | 63.9                                   | 36.1                                  | 163.9              |
| 7.103  | 66.6                                   | 33.4                                  | 166.6              |
| 7.154  | 69.2                                   | 30.8                                  | 169.2              |
| 7.212  | 72.0                                   | 28.0                                  | 172.0              |
| 7.261  | 74.4                                   | 25.6                                  | 174.4              |
| 7.313  | 76.8                                   | 23.2                                  | 176.8              |
| 7.364  | 78.9                                   | 21.1                                  | 178.9              |
| 7.412  | 80.8                                   | 19.2                                  | 180.8              |
| 7.462  | 82.5                                   | 17.5                                  | 182.5              |
| 7.504  | 84.1                                   | 15.9                                  | 184.1              |
| 7.561  | 85.7                                   | 14.3                                  | 185.7              |
| 7.610  | 87.0                                   | 13.0                                  | 187.0              |
| 7.655  | 88.2                                   | 11.8                                  | 188.2              |
| 7.705  | 89.4                                   | 10.6                                  | 189.4              |
| 7.754  | 90.5                                   | 9.5                                   | 190.5              |
| 7.806  | 91.5                                   | 8.5                                   | 191.5              |
| 7.909  | 93.2                                   | 6.8                                   | 193.2              |
| 8.018  | 94.7                                   | 5.3                                   | 194.7              |
| 8.10   | 95.8                                   | 4.2                                   | 195.8              |
| 8.20   | 97.0                                   | 3.0                                   | 197.0              |

U-shaped tube with a wall less than 0.1 mm in thickness. With a length of 200 mm and inside diameter of 1.5 mm the volume is about 0.10 ml. The tube is bent by its own weight while a positive pressure is maintained within the tube and the mid portion heated gently with a small luminous flame. Then the wall of the tube is thinned by thoroughly cleaning and dipping in hydrofluoric acid. The process is controlled by frequently removing the electrode running in water and then observing the change in the amount of pressure required to bend the tube. A satisfactory electrode will give a potential of approximately 0.0615 volt per unit change in pH at 37 C. The electrode is held in position between the halves of a rubber stopper which is sealed and fastened together with thick rubber cement and a wire band. 0.1 N HCl surrounds the electrode in the inner chamber and fills the side tube F. An Ag/AgCl electrode connected to an amplifier at E, is placed within a paraffined tube at D and a large glass bead is placed above the constriction in order to minimize movements of solution around the Ag/AgCl electrode. Water from a constant temperature bath maintains the temperature by circulating in the outer chamber. We have used the electrode successfully at 0.5, 25, and 37 C. A small copper wire is grounded at G and passed

TABLE 14B

WALPOLE'S ACETATE BUFFER MIXTURES FOR pH 3.6 TO 5.8  
(*J Chem Soc* 105 2001 1914)

| 1 Sodium acetate 0.1 N 7 gm $\text{NaCH}_3\text{CO}_2 \cdot \text{H}_2\text{O}$ diluted to 1000 ml with water   |                           |                        |
|---|---------------------------|------------------------|
| 2 Acetic acid 0.1 N approximately 113 ml of 99% glacial acetic acid diluted to 1000 ml with water Standardize against 0.1 N sodium hydroxide with phenolphthalein as an indicator |                           |                        |
| pH  | 0.1 N Sodium acetate (ml) | 0.1 N Acetic Acid (ml) |
| 3.6   | 75                        | 92.5                   |
| 3.8   | 100                       | 89.0                   |
| 4.0   | 130                       | 86.0                   |
| 4.2   | 165                       | 83.5                   |
| 4.4   | 200                       | 80.0                   |
| 4.6   | 230                       | 77.0                   |
| 4.8   | 260                       | 74.0                   |
| 5.0   | 290                       | 71.0                   |
| 5.2   | 320                       | 68.0                   |
| 5.4   | 350                       | 65.0                   |
| 5.6   | 380                       | 62.0                   |
| 5.8   | 410                       | 59.0                   |

TABLE 14C

CLARK'S 0.05 M BORATE BUFFER MIXTURES FOR pH 7.8 TO 10  
(The Determination of Hydrogen Ions 1938)

| 1 0.2 M boric acid 1.405 gm $\text{H}_3\text{BO}_3$ /l      |                                    |                          |
|---|------------------------------------|--------------------------|
| 0.1 M potassium chloride 14.91 gm $\text{KCl}$ /l           |                                    |                          |
| 0.1 M sodium hydroxide from conc. $\text{CO}_2$ free alkali |                                    |                          |
| pH  | 0.1 M $\text{H}_3\text{BO}_3$ (ml) | 0.1 M $\text{NaOH}$ (ml) |
| 7.8   | 50                                 | 6.5                      |
| 8.0   | 50                                 | 4.00                     |
| 8.2   | 50                                 | 5.00                     |
| 8.4   | 50                                 | 8.55                     |
| 8.6   | 50                                 | 12.00                    |
| 8.8   | 50                                 | 18.40                    |
| 9.0   | 50                                 | 14.00                    |
| 9.2   | 50                                 | 16.70                    |
| 9.4   | 50                                 | 20.00                    |
| 9.6   | 50                                 | 28.85                    |
| 9.8   | 50                                 | 40.80                    |
| 10.0  | 50                                 | 43.90                    |

Dilute to 200 ml with distilled water in each case

through the rubber tubing and connection into the outer chamber in order to remove stray electric currents which might result from the circulating water

- Calomel half cell with saturated  $\text{KCl}$  and  $\text{KCl}$  reservoir
- Apparatus for amplification and measurement of potentials Any circuit which will measure accurately to 0.1 mv will serve
- A constant temperature water bath set at 31°C
- A number of 2 ml syringes

#### Reagents

- Potassium chloride ( $\text{KCl}$ ) a saturated solution in water
- Sodium chloride ( $\text{NaCl}$ ) isotonic 0.9% solution in water
- Buffers on each side of the range of the unknowns
- For blood work use M/L phosphate buffers These should be kept in an ice box as two separate solutions (a) and (b) and mixed as needed in the proportions shown in the accompanying Table 15



TABLE 14D

Sørensen's GLYCINE HYDROCHLORIC ACID MIXTURE FOR pH 1.1-3.5  
(*Ergeb. Physiol.* 12: 393, 1912)

| 1 Glycine sodium chloride solution 0.1 M<br>Dissolve 7.505 gm of glycine and 5.85 gm of sodium chloride in 1000 ml of water |             |       |
|---|-------------|-------|
| 2 Hydrochloric acid an 0.1 N solution   |             |       |
| Temperature = 18 C  |             |       |
| Glycine-NaCl<br>(ml)  | HCl<br>(ml) | pH    |
| 0.0   | 10.0        | 1.038 |
| 1.0   | 9.0         | 1.146 |
| 2.0   | 8.0         | 1.251 |
| 3.0   | 7.0         | 1.419 |
| 4.0   | 6.0         | 1.645 |
| 5.0   | 5.0         | 1.932 |
| 6.0   | 4.0         | 2.279 |
| 7.0   | 3.0         | 2.607 |
| 8.0   | 2.0         | 2.92  |
| 9.0   | 1.0         | 3.341 |
| 9.5   | 0.5         | 3.679 |

TABLE 14E

Sørensen's GLYCINE SODIUM HYDROXIDE MIXTURES FOR pH 8.5-12.5  
(*Ergeb. Physiol.* 12: 393, 1912)

## Reagents

- 1 Dissolve 7.505 gm of glycine and 5.85 gm of sodium chloride in 1000 ml of water (0.1 M glycine plus 0.1 M NaCl)
- 2 Sodium hydroxide an 0.1 N solution

| Glycine<br>(ml) | Sodium Hydroxide<br>(ml) | 18 C<br>pH | 20 C<br>pH | 30 C<br>pH |
|-----------------|--------------------------|------------|------------|------------|
| 9.5             | 0.5                      | 8.58       | 8.53       | 8.3        |
| 9.0             | 1.0                      | 8.93       | 8.88       | 8.67       |
| 8.0             | 2.0                      | 9.36       | 9.31       | 9.08       |
| 7.0             | 3.0                      | 9.71       | 9.66       | 9.42       |
| 6.0             | 4.0                      | 10.14      | 10.09      | 9.83       |
| 5.5             | 4.5                      | 10.48      | 10.42      | 10.17      |
| 5.1             | 4.9                      | 11.07      | 11.01      | 10.74      |
| 5.0             | 5.0                      | 11.31      | 11.25      | 10.97      |
| 4.9             | 5.1                      | 11.57      | 11.51      | 11.22      |
| 4.5             | 5.5                      | 12.10      | 12.04      | 11.74      |
| 4.0             | 6.0                      | 12.40      | 12.33      | 12.03      |
| 3.0             | 7.0                      | 12.67      | 12.60      | 12.9       |
| 2.0             | 8.0                      | 12.86      | 12.79      | 12.47      |
| 1.0             | 9.0                      | 12.97      | 12.90      | 12.51      |

(a) Monobasic potassium phosphate ( $\text{KH}_2\text{PO}_4$ ) 9.078 gm/1000 ml

(b) Dibasic sodium phosphate ( $\text{Na}_2\text{HPO}_4$ ), 9.469 gm/1000 ml

The salts must be dried at 95-100 C. Pure salts may be obtained from La Motte and Co

- 5 (a) For buffer of pH 6.21 to 37 C 5 ml of  $\text{Na}_2\text{HPO}_4$  are mixed with 0 ml of  $\text{KH}_2\text{PO}_4$

(b) For buffer of pH 7.35 at 37 C, 20 ml of  $\text{Na}_2\text{HPO}_4$  are mixed with 5 ml of  $\text{KH}_2\text{PO}_4$

## Procedure—Testing System

- 1 Rinse electrode twice with saline
- 2 Rinse twice with buffer
- 3 Put in sample of buffer




Fig 24—Glass electrode for electrochemical work.

TABLE 15

SØRENSEN'S PHOSPHATE BUFFER SOLUTIONS AT 20°C AND 37°C FOR pH RANGE 5.8 TO 8.2  
 (M/15 acid phosphate 9.08 gm dry  $\text{KH}_2\text{PO}_4$ /liter,  
 M/15 basic phosphate 9.469 gm anhydrous  $\text{Na}_2\text{HPO}_4$ /liter)

| pH      |         | M/15 Phosphate Solution        |                               |
|---------|---------|--------------------------------|-------------------------------|
| At 20°C | At 37°C | $\text{Na}_2\text{HPO}_4$ (ml) | $\text{KH}_2\text{PO}_4$ (ml) |
| 5.80    | 5.77    | 8.0                            | 92.0                          |
| 5.90    | 5.87    | 9.9                            | 90.1                          |
| 6.00    | 5.97    | 12.2                           | 87.8                          |
| 6.10    | 6.07    | 15.3                           | 84.7                          |
| 6.20    | 6.17    | 18.6                           | 81.4                          |
| 6.24    | 6.21    | 20.0                           | 80.0                          |
| 6.30    | 6.27    | 22.4                           | 77.6                          |
| 6.40    | 6.37    | 26.7                           | 73.3                          |
| 6.50    | 6.47    | 31.8                           | 68.2                          |
| 6.60    | 6.57    | 37.5                           | 62.5                          |
| 6.70    | 6.67    | 43.5                           | 56.5                          |
| 6.80    | 6.77    | 49.6                           | 50.4                          |
| 6.90    | 6.87    | 55.4                           | 44.6                          |
| 7.00    | 6.97    | 61.1                           | 38.9                          |
| 7.10    | 7.07    | 66.6                           | 33.4                          |
| 7.20    | 7.17    | 72.0                           | 28.0                          |
| 7.30    | 7.27    | 76.8                           | 23.2                          |
| 7.35    | 7.35    | 80.0                           | 20.0                          |
| 7.40    | 7.37    | 80.8                           | 19.2                          |
| 7.50    | 7.47    | 84.1                           | 15.9                          |
| 7.60    | 7.57    | 87.0                           | 13.0                          |
| 7.70    | 7.67    | 89.4                           | 10.6                          |
| 7.80    | 7.77    | 91.5                           | 8.5                           |
| 7.90    | 7.87    | 93.2                           | 6.8                           |
| 8.00    | 7.97    | 94.7                           | 5.3                           |
| 8.10    | 8.07    | 95.8                           | 4.2                           |
| 8.20    | 8.17    | 97.0                           | 3.0                           |

4. Flush stopcock with KCl.

5. Turn stopcock to connect KCl with buffer and read at once.

6. Read again in one minute.

- The two readings should be within 0.5 millivolts. Temperature is probably drifting if greater divergence is noted. This drift must be eliminated. When a stable reading is obtained repeat steps 1 to 6 with second buffer. Then repeat 3 & 6 with a second sample.
- These 4 readings on the second buffer should be within 1 mv of each other and should differ from the readings of the other buffer by 70 mv (i.e. 61.5 mv per pH unit at 37°C or 58.1 mv at 20°C). There is a linear relation between temperature and EMF.
- If the difference between these two buffers (pH 6.21 and 7.35) exceeds 70 mv the buffers are not correct. If it is less than 70 the buffers may be wrong but it is more probable that the trouble lies in the electrode. It is either too thick, too old, cracked, or perhaps there is an electrical leak around it. However if the difference is constant and not less than 60 mv the electrode is usable.

#### Procedure—Measuring Unknown

- Wash electrode four times with warm isotonic salt solution.
- Put in blood sample by means of small syringe. Sample should be fresh and must be handled as anaerobically as possible, avoid getting oil into the electrode. Push at least a quarter of a ml slowly through the electrode so that the blood in contact with the membrane is not the first blood layer to come from the syringe.

- 3 Rinse stopcock with KCl
- 4 Turn stopcock to connect KCl with sample and read at once
- 5 Read again in 1 minute
- 6 Follow E M F for 5 minutes by occasional readings. There will probably be a greater change with the blood than with the buffers during the first minute but it should not be more than 1 mv. The change between 1 and 5 minutes should be less than 1 mv. Take the 1 minute reading.
- 7 Repeat 2, 3, 4 and 5 if there is enough of the sample.
- 8 Repeat steps 1 through 6 with the buffer which is nearest in pH to the unknown.
- 9 Repeat steps through 6 with the other buffer. If the temperature is constant and the electrode is working satisfactorily it is safe to use as many as 3 samples between rechecks against a buffer and the check against the two buffers may be left to the end.

#### Calculation

$$\begin{aligned} \text{pH of unknown} = & \\ & \text{pH acid buffer} + [\text{E M F unknown} - \text{E M F acid buffer}] + \\ & \left[ \frac{\text{E M F alkaline buffer} - \text{E M F acid buffer}}{\text{pH alkaline buffer} - \text{pH acid buffer}} \right] \end{aligned}$$

#### Example

|                     |        |
|---------------------|--------|
| pH alkaline buffer  | 7.35   |
| pH acid buffer      | 6.01   |
| EMF alkaline buffer | 301 mv |
| EMF acid buffer     | 251 mv |
| EMF unknown         | 314 mv |

$$\begin{aligned} \text{pH of unknown} &= 6.01 + \left[ (314 - 251) + \frac{(301 - 251)}{(7.35 - 6.01)} \right] \\ &= 6.01 + (59.3 + 59.6) \\ &= 7.00 \end{aligned}$$

#### Precautions

- 1 All solutions put into the electrode should be within 1 or 2°C of the temperature of the solution to be tested. For blood this is 37°C. Electrode must be constant within 0.1 and should be within 1°C of the temperature required. It is not necessary to have the calomel cell and KCl at the same temperature as the electrode but they must not change in temperature during the experiment.
- 2 Electrode should be left filled with water if in constant use and if not with 0.01 N HCl.
- 3 Saturated KCl should not come in contact with electrode.
- 4 The solutions are introduced with small syringes and should be at about the temperature of the electrode.
- 5 The main difficulties with the glass electrode arise from the fragility of the electrode and its very high electrical resistance. One is measuring a small potential in series with a resistance of several megohms and great care must be taken to avoid electrical leaks over moist surfaces. The N/10 HCl, the KCl and the NaCl if split promiscuously are likely to lead to leaks especially in damp weather and if the electrode is used at low temperatures condensation will do the same.
- 6 The temperature coefficient for drawn human blood under the conditions of measurement is approximately 0.01 pH units per 1°C. As the temperature rises

the pH falls. The same is true to a much smaller extent of the phosphate buffers. However the EMF of the whole system changes greatly with temperature and one should work as near as possible to the desired temperature and be particularly careful to avoid fluctuations.

## 4 Measurement of pH in Solutions

(Conditions Not Anaerobic Beckman pH Meter)

### Principle

A glass electrode, when immersed in an unknown solution, develops a potential difference the magnitude depending on the hydrogen ion concentration of the solution. This potential difference is measured by combining the glass electrode with some standard half cell such as the saturated calomel electrode and measuring the voltage of the system.

### Apparatus

- 1 A Beckman pH meter Model M

### Reagents

- 1 Potassium chloride a saturated solution approximately 28 gm diluted up to 100 ml with water
- 2 Buffer solution pH 7.0 (see section on glass electrode under anaerobic conditions)

### Procedure

- 1 Remove rubber cap from the end of the calomel electrode. Mount the electrodes in their holder so that the calomel electrode projects slightly lower than the glass electrode and be sure they cannot touch the bottom of the beaker when the holder is lowered against the stop. Remove the stopper from the calomel electrode during use.
- 2 With control now in off position the meter should read exactly pH 7.0. Adjust if necessary by turning the screw below the meter's ale.
- 3 With push button in lower right hand corner in up position, turn switch to on and leave in this position for a few seconds.
- 4 Turn switch to 7.14 range and adjust meter to pH 7 with amplifier control knob in lower left corner.
- 5 Turn switch to 0.7 range and adjust to pH 7 with standard knob. Check the amplifier adjustment occasionally during use.
- 6 Immerse electrodes in a buffer solution of pH 7.0.
- 7 Place switch at proper range, press button and adjust meter to buffer pH by means of the AP knob.
- 8 To determine pH of test solutions, clean electrodes, immerse in test solution and press button. The meter will indicate the pH automatically. If the needle goes off the scale, turn the switch to the other range.

### Calculation

- 1 pH is read directly.

### Precautions

- 1 The standard knob should require only infrequent adjustment.
- 2 The electrodes should be immersed a full inch in the solution to avoid slight residual effects.
- 3 When not in use the electrodes should be immersed in a beaker of water and the calomel electrode stopper replaced.
- 4 Keep the electrodes clean at all times and do not allow the test solution to dry on them.

- 5 Keep the solutions at approximately 25 C as the apparatus has been calibrated at that temperature
- 6 When testing solutions of pH 9.0 or higher containing alkali metal ions readings should be corrected for alkali ion error. This correction may be determined experimentally

## 5 Hydrogen Ion Concentration, Miscellaneous References

### Plasma pH

- Van Slyke D D Weisiger J R and Van Slyke K K. Photometric Measurement of Plasma pH *J Biol Chem* 178 743-756 (June) 1949

### Urine pH

- Van Slyke D D Weisiger J R and Van Slyke K K. Photometric Measurement of Urine pH *J Biol Chem* 178 777-781 (June) 1949

## SECTION V

# MICROBIOLOGICAL PROCEDURES

---

### 1 Pteroylglutamic Acid (Folic Acid)

#### Reference

Modified Microbiological Method of Tepley, L. J. and Elvehjem, C. A. The Titrimetric Determination of Lactobacillus Casei Factor and Folic Acid J Biol Chem 157 303-309 (Jan.) 1945

#### Principle

The response (acid production) of *Lactobacillus casei* E to graded increments of the sample and of a standard solution of pteroylglutamic acid added to media furnishing all the nutrients required by this microorganism except the nutrient under assay is measured by titration with sodium hydroxide. The acid produced is a direct function of the original concentration of pteroylglutamic acid.

#### Apparatus

- 1 Autoclave horizontal type
- 2 Beakers assorted 25-1000 ml
- 3 Burette 50 ml
- 4 Centrifuge
- 5 Cotton nonabsorbent
- 6 Graduated cylinders glass stoppered assorted 25-1000 ml
- 7 Filter paper qualitative
- 8 Flasks Erlenmeyer 50 and 125 ml
- 9 Flasks volumetric glass stoppered assorted 25-2000 ml
- 10 Funnels 2-3 inch diameter
- 11 Incubator water jacketed electric ( $\pm 0.1^\circ\text{C}$ )
- 12 Inoculating needle platinum or nichrome steel
- 13 Automatic pipette (Fisher Volumstat)
- 14 Graduated pipettes 0.5, 1.0, 2.0, 5.0, 10.0 and 25.0 ml
- 15 Volumetric transfer pipettes 0.5, 1.0, 2.0, 5.0, 10.0, 25.0 and 50.0 ml
- 16 Refrigerator
- 17 Arnold sterilizer
- 18 Hot air sterilizer
- 19 Sterilizing cans
- 20 Hypodermic syringes 5 to 10 ml
- 21 Interval timers
- 22 Test tubes uniform lipless Pyrex
- 23 Test tube covers aluminum rectangular
- 24 Test tube supports metal
- 25 Water bath constant temperature

#### Reagents

- 1 Water, glass distilled.
- 2 0.02 N NaOH and 0.1 N NaOH
- 3 0.1% bromthymol blue
- 4 1 N NaOH (approx.)

- 5 0.1 N HCl (approx)
- 6 Chloroform USP
- 7 Toluene USP
- 8 Vitamin free acid hydrolyzed casein 10<sup>6</sup>
- 9 Adenine Guanine and Uracil These are dissolved in 75 ml distilled water and 5 ml or more conc HCl Heat and dilute to 1 m<sup>3</sup> of each per ml of final solution
- 10 Asparagine Dissolve 10 gm asparagine in 1 liter of water Add a few drops of chloroform and toluene
- 11 Alanine Dissolve 5 gm dl alanine in 50 ml of water Add a few drops of chloroform and toluene
- 12 Cystine and Tryptophane Each of these is dissolved in a small amount of 10 N HCl and diluted to 4 mg per ml and 2 m<sup>3</sup> per ml respectively
- 13 Salt Solution A Dissolve 2.0 gm KH<sub>2</sub>PO<sub>4</sub> and 2.0 gm K<sub>2</sub>HPO<sub>4</sub> in water to make 50 ml Keep under toluene and chloroform
- 14 Salt Solution B Dissolve 10 gm of MgSO<sub>4</sub> 7H<sub>2</sub>O 0.5 gm NaCl 0.5 gm FeSO<sub>4</sub> 7H<sub>2</sub>O and 0.4 gm MnSO<sub>4</sub> H<sub>2</sub>O in water and dilute to 50 ml Add 10 drops of HCl to clear solution Add chloroform and toluene
- 15 Vitamin Solution 10 mg of thiamin hydrochloride 60 mg pyridoxine hydrochloride 30 mg nicotinic acid 20 mg calcium pantothenate 0.5 mg p aminobenzoic acid and 1 ml salt solution A are added to 100 ml water Add chloroform and toluene and store in refrigerator
- 16 Pteroylglutamic Acid Standard Solution Weigh 10 mg pteroylglutamic acid and transfer to 1000 ml volumetric flask using 0.0 N NaOH to make up to volume Using 1 ml of this solution plus 1 ml of salt solution A dilute to 100 ml with 0.0 N NaOH for a standard stock solution of 100 mcg per ml Dilute to 0.1 milli micrograms per ml for the working standard
- 17 Biotin 25 mcg of crystalline (free acid) biotin are transferred to a 50 ml volumetric flask Add 7 ml of salt solution A and adjust to volume with water
- 18 Riboflavin Weigh 50 mg USP reference standard riboflavin (dried in vacuum desiccator) and transfer to 2 liter volumetric flask Add 1500 ml H<sub>2</sub>O 24 ml glacial acetic acid and warm Cool and make to volume This solution contains 25 mcg riboflavin per ml of 0.0 N acetic acid
- 19 Yeast Dextrose Agar 15 gm agar (Difco Bacto) are dissolved in about 150 ml of hot water After agar is in solution add 5 gm glucose 20 gm yeast extract and dilute to 1 liter Tube in 7.8 ml portions Plug tubes with nonabsorbent cotton and autoclave at 15 lbs pressure for 15 minutes Cool and store in refrigerator
- 20 Yeast Dextrose Broth Dissolve 5 gm Bacto yeast extract 10 gm NaC<sub>2</sub>H<sub>3</sub>O<sub>2</sub> 3H<sub>2</sub>O 5 gm glucose 5 ml salt A and 5 ml salt B in water Adjust to pH 6.8-7.0 and dilute to 1 liter Tube in 10 ml portions plug with nonabsorbent cotton and autoclave at 15 lbs pressure for 15 minutes Cool and store in refrigerator Maximum amount of pteroylglutamic acid used in preparation of standard curve may be incorporated with each tube as a supplement to the broth
- 21 Isotonic Salt Solution Tube in 10 ml portions plug with cotton and sterilize in autoclave
- 22 Lactobacillus Stock Culture The test organism is a pure culture obtained originally from the American Type Culture Collections Georgetown University Medical School Washington D C as number 7469 The organism is carried as a stab culture in a yeast dextrose agar medium and transferred at weekly intervals
- 23 Takadiastase Dissolve 2 gm of commercial takadiastase in 100 ml of water

#### Procedure

- 1 Preparation of the Inoculum On the day prior to use transfer cells from the stock culture to a sterile tube of inoculum culture medium (Peagent "0") Inoculum



bate this culture for 16-24 hours at 37.5 C. After incubation, the cells are centrifuged out, resuspended in saline and 1 drop added to 10 ml of 0.9 per cent saline. For the inoculum one drop of this cell suspension is added per tube.

## 2. Preparation of Samples

a. Food Weigh accurately about 1 gm food. Add 25 ml acetate buffer pH 4.5, 1 ml takadiastase enzyme (2%) solution, 2 ml of toluene, plug with cotton and incubate 37.5 C for 18-24 hours. Adjust pH to 6.6-6.8, steam for 10 minutes in autoclave, dilute 1:20 for normal diet and filter.

b. Feces About 0.8 gm feces is treated by enzyme method above, final dilution is 1:25,000.

c. Urine Dilute 2 ml urine to 100 ml with water in a dilution of 1 to 50.

d. Plasma and Blood Add 8 ml acetate buffer to 25 ml Erlenmeyer flask and add 0.1 ml blood or plasma. Wash pipette several times with acetate buffer. Add 0.1 ml 2% takadiastase solution. Gently mix and then add 1 drop of toluol to each Erlenmeyer flask. Always run blank consisting of all components except sample. Incubate at 37.5 C for 24 hours  $\pm$  2 hours. Steam samples in autoclave for 15 minutes. Cool and adjust to pH 6.6-6.8. Dilute blood 1:500, plasma 1:50.

## 3. Preparation of Basal Medium Prepare the basal medium having the composition shown below. About 500 ml is required for 100 tubes.

20 gm anhydrous glucose

33.4 gm sodium acetate

1.25 gm  $K_2HPO_4$

50 ml casein hydrolysate 10%

50 ml l-tryptophan solution (2 mg per ml)

50 ml l-cystine (4 mg per ml)

10 ml adenine, guanine, uracil (1 mg each per ml)

10 ml asparagine (10 mg per ml)

10 ml D-alanine (20 mg per ml)

0.5 ml salt B

4 ml vitamin solution

50 mcg thiamin

300 mcg pyridoxine

150 mcg niacin

100 mcg pantothenic acid

2.5 mcg p-aminobenzoic acid

8 ml biotin 0.1 mcg per ml

8 ml riboflavin 0.5 mcg per ml

} Per ml

Adjust pH to 6.6-6.8 and dilute to 500 ml.

4. Preparation of Standard Tubes To duplicate tubes 13 x 100 mm, add 0.0, 0.15, 0.25, 0.5, 0.75, 1.0, 1.25, 1.5, 1.75 and 2.0 ml of standard pteroylglutamic acid solution (0.1 millimicrogram per ml). Add sufficient water to bring the volume in each tube to 2 ml. To each of these tubes add 1 ml of basal media. An automatic pipetting machine is a valuable asset in this step.

5. Preparation of Assay Tubes To duplicate tubes 13 x 100 mm add 0.5, 1.0, 1.5 and 2.0 ml of the specimen. Add sufficient water to bring the volume to 2.0 ml. To each of these tubes add 1 ml of basal medium.

6. Sterilization Place aluminum covers over the assay tubes in racks. Autoclave at 15 lbs pressure for 10 minutes.

7. Inoculation and Incubation Cool all the tubes to the incubation temperature. Aseptically inoculate each tube with one drop of inoculum. Incubate for 72 hours at 37.5 C.

8. Titration Transfer the contents of a tube to a 50 ml Erlenmeyer flask rinsing once with 10 ml of water. Titrate the contents of the flask with 0.02 N NaOH to a pH of 6.8 using bromthymol blue as an internal indicator.

**Calculation**

- 1 Draw a standard curve for the assay by plotting millimicrograms of pteroyl glutamic acid (abscissa) in the standard series against ml 0.0 N NaOH (ordinate)
- 2 Determine the pteroylglutamic acid content of the tubes in the unknown series by interpolation of the titer values on the standard curve. Use only values which do not differ from the average by more than 10%.

$$\text{millimicrograms folic acid/ml specimen} = \frac{\text{average/tube}}{\text{aliquot volume}} \times \text{dilution factor}$$

**Example**

The method of calculation is shown in the table below

| ASSAY TUBE<br>ml EXTRACT | TITRATION AFTER<br>INCUBATION<br>ml | PTEROYLGLUTAMIC<br>ACID EVALUATED<br>FROM STD. CURVE<br>mμ | PTEROYLGLUTAMIC<br>ACID PER ml OF<br>TEST EXTRACT<br>mμ |             |
|--------------------------|-------------------------------------|--|---|-------------|
| 0.5                      | 4.9                                 | 0.019  | 0.038   | } Av. 0.031 |
| 1.0                      | 5.9                                 | 0.030  | 0.060   |             |
| 1.5                      | 6.5                                 | 0.041  | 0.077   |             |
| 2.0                      | 8.0                                 | 0.068  | 0.034   |             |

Dilution was 1 to 40. Final answer is therefore  $40 \times 0.031 = 1.24$

**Precautions**

- 1 In order to avoid erratic results only glass distilled water and suitably pure re search products should be used.
- 2 Digestion with takad stage and papain does not give maximum values for the folic acid content of many natural materials. Each kind of material needs individual research to determine the best means for hydrolysis.

**2 Nicotinic Acid****References**

- 1 Modified Microbiological Method of Snell E. E. and Wright L. D. A Micro biological Method for the Determination of Nicotinic Acid J. Biol. Chem. 139 65-685 (May) 1941
- Krehl W. A. Strong F. M. and Flvehjem C. A. Determination of Nicotinic Acid Ind. Eng. Chem. Anal. Ed. 15 471-475 (July) 1943

**Principle**

The microbiological determination of nicotinic acid is based upon the growth stimulation of *Lactobacillus arabinosus*. In a basal medium that is niacin free but otherwise complete growth responses of the organism are measured quantitatively in standard and unknown solutions by titration of the acid produced.

**Apparatus**

- 1 Same as in pteroylglutamic acid assay
- 2 **Vitamin Solution** Incorporate 5 mg. thiamin hydrochloride 5 mg. calcium pantothenate 5 mg. p-aminobenzoic acid 10 mg. pyridoxine 1 ml. salt solution A and dilute to 100 ml. Add chloroform and toluene and store in refrigerator.
- 3 **Standard Nicotinic Acid Solution** A solution of nicotinic acid 0.5 mg. per ml. (500 micrograms per ml.) is made by dissolving 1.15 mg. of crystalline vitamin in the proper amount of distilled water. 0 ml. of this solution plus 1 ml. salt solution A

are made up to 100 ml in a volumetric flask to give a concentration of 100 mcg per ml. Dilute to 50 millimicrograms per ml for working standard for micro-assay

- 4 *Stock Culture* Stock cultures are carried as stab cultures in a yeast dextrose agar medium. These cultures are transferred at weekly intervals incubated 24 hours at 37.5 C and then stored in the refrigerator until needed. The organism used is *Lactobacillus arabinosus* 175 and may be secured from the American Type Culture Collection, where it is listed as No. 8014.

#### Procedure

- 1 *Preparation of the Inoculum* On the day prior to use transfer cells from the stock culture to a sterile tube of yeast dextrose broth supplemented with the maximum quantity of nicotinic acid used for preparation of standard curve. For the inoculum the cells are centrifuged out, the supernatant liquid decanted and the cells re-suspended in 10 ml isotonic saline solution. One drop of this cell suspension is added to 10 ml of saline to provide a working cell suspension. Use 1 drop per tube for inoculation.
- 2 *Preparation of Samples*
  - a. Food—Weigh accurately about 1 gm food, add 20 ml 1 N HCl, autoclave 30 minutes at 15 lbs., neutralize to pH 6.6-6.8, dilute 1 to 100 in water and filter.
  - b. Blood—Dilute 0.5 ml to 50 ml with water.
  - c. Urine—To 10 ml urine add 10 ml 1 N HCl, autoclave 30 minutes at 15 lbs., cool and adjust pH to 6.6-6.8. For assay dilute this solution 1 to 4 with water.
  - d. Weigh accurately about 1 gm feces. Add 20 ml 1 N HCl, autoclave 30 minutes at 15 lbs., cool, adjust to pH 6.6-6.8, dilute to 100 ml with water and filter. For use dilute filtrate 1 to 10 with water (i.e. final dilution 1:1000).
- 3 *Preparation of Basal Medium* Prepare a basal medium having the composition shown below:
  - a. Mix the following:
    - 20 gm anhydrous glucose
    - 33.4 gm Na acetate 3H<sub>2</sub>O (20 gm anhydrous)
    - 50 ml casein hydrolysate 10%
    - 50 ml tryptophane (2 mg l tryptophane per ml)
    - 50 ml l cystine (4 mg per ml)
    - 10 ml adenine guanine uracil (1 mg of each per ml)
    - 2 ml vitamin solution (see reagents)
    - 5 ml salt mixture A
    - 5 ml salt mixture B
    - 4 ml riboflavin (25 mcg per ml)
    - 2 ml biotin (0.1 mcg per ml)
  - b. Adjust the pH to 6.6-6.8. Dilute to 500 ml with water.
- 4 *Preparation of Standard Tubes* Same as for pteroylglutamic acid assay using standard nicotinic acid solution (50 millimicrograms per ml) for working standard.
- 5 *Preparation of Assay Tubes* *Sterilization, Inoculation and Incubation.* Titration. These follow the same general microbiological techniques used in pteroylglutamic acid assay.

#### Calculation

Same as for pteroylglutamic acid.

#### Precautions

- 1 *Lactobacillus arabinosus* responds equally to niacin and its amide. Differential microbiological assays have to be used to distinguish between the several metabolites of this vitamin.

## 3 Biotin

## References

- 1 Modified Microbiological Method of Shull G M Hutchings B L and Peterson V H. A Microbiological Assay for Biotin J Biol Chem 142 913-910 (Feb) 1944
- 2 Tepley L J and Elvehjem C A. The Titrimetric Determination of Lactobacillus Casei Factor and Folic Acid J Biol Chem 157 303-309 (Jan) 1945

## Principle

Biotin may be determined microbiologically by measurement of its effect upon the growth of *Lactobacillus casei* E. The titratable acidity produced by *L. casei* E is a function of the biotin concentration of the medium.

## Apparatus

- 1 Same as for pteroylglutamic acid assay

## Reagents

Reagents are the same as for the assay of pteroylglutamic acid with a few additions.

- 1 Manganese dioxide (MnO USP)
- 2 Norite or activated charcoal
- 3 Superoxol (30% H<sub>2</sub>O ACS)
- 4 Peptone (Difco Bacto)
- 5 Peptone Supplement. Dissolve 5 gm of peptone in 100 ml of warm distilled water. Adjust to pH 3. Stir with 25 gm Norite for 1 hour. Filter. Dilute to 90 ml. Add 10 ml Superoxol (30% H<sub>2</sub>O). Do not pipette!! Allow to stand at room temperature 24 hours. Adjust to pH 7. Add about 15 gm MnO<sub>2</sub> slowly with stirring. When evolution of oxygen ceases filter solution and dilute to 50 ml with water. Add chloroform and toluene and store in refrigerator.

## Stock Culture

Same as for pteroylglutamic acid.

## Procedure

- 1 Preparation of the Inoculum. Same as for pteroylglutamic acid assay except instead of pteroylglutamic acid the yeast dextrose broth should be supplemented with the minimum quantity of biotin used in preparation of the standard error.
- 2 Preparation of Samples
  - a Feces. Weigh accurately about 0.7 gm feces into a 50 ml Erlenmeyer flask. Add 5 ml acetate buffer, 1 ml 1% papain (emulsify 0.5 gm papain in 0.5 ml glycerine and dilute with 4 ml water), 1 ml 2% alkaline phosphatase and 2 ml toluene. Plug with cotton, incubate 18-24 hours at 37°C, adjust to 6.6-6.8 steam 10 minutes and cool. Dilute 1 to 100 with water, filter and dilute the filtrate to a final dilution of 1 to 10,000.
  - b Food. Weigh accurately a sample of about 1 gm. Treat like feces through the stage of steaming. Dilute 1 to 50, filter and use directly. An alternative method is sometimes required: add 1 gm food to 5 ml 1 N HCl, autoclave 30 minutes at 15 lbs, adjust to pH 6.6-6.8, dilute 1 to 50, filter and use.
  - c Urine. Dilute one ml to 50 with water.
  - d Blood and Plasma. Mix 10 ml blood or plasma with 10 ml 1 N HCl, autoclave at 15 lbs for 30 minutes and cool. Adjust to pH 6.6-6.8, dilute 1 to 5, filter and use.

Preparation of Basal Medium. Prepare the basal medium having the composition shown below.

- a Mix together

0 gm glucose

33.4 gm Na<sub>2</sub> acetate 3H<sub>2</sub>O (0 gm anhydrous)

1.25 gm  $K_2HPO_4$

50 ml casein hydrolysate 10%

50 ml l-tryptophane (2 mg per ml)

50 ml l-cystine (4 mg per ml)

10 ml adenine guanine uracil solution (1 mg each per ml)

10 ml asparagine (10 mg per ml)

10 ml d,l alpha alanine (20 mg per ml)

10 ml peptone supplement (norite and peroxide treated)

40 ml vitamin solution (same as for pteroylglutamic acid)

10 ml folic acid (1 mcg per ml)

80 ml riboflavin (25 mcg per ml)

50 ml salt mixture B

b Adjust to pH 6.6-6.8 and dilute to 500 ml with water

4 *Preparation of Standard Tubes* Same as for pteroylglutamic acid assay using standard biotin solution (0.2 millimicrogram per ml) for working standard.

5 *Preparation of assay tube* sterilization inoculation and incubation titration and calculation are all according to the same general microbiological techniques as used for pteroylglutamic acid assay

## 4 Pantothenic Acid

### References

1. Modified Microbiological Method of Cheldelin V. H. Hoag, E. H. and Sarett H. P. The Pantothenic Acid Requirements of Lactic Acid Bacteria *J. Bact.* 49: 41-45 (Jan.) 1945
2. Hoag E. H. Sarett H. P. and Cheldelin V. H. Use of *Lactobacillus Arabinosus* 175 for Microassay of Pantothenic Acid *Ind. and Eng. Chem. Anal. Ed.* 17: 60-62 (Jan.) 1945

### Principle

Pantothenic acid is estimated from its growth stimulation of *Lactobacillus arabinosus* 175. Lactic acid formation is estimated titrimetrically. The bound vitamin is first liberated by means of digestive enzymes.

### Apparatus

- 1 Same as for pteroylglutamic acid assay

### Reagents

Reagents are the same as for pteroylglutamic acid assay with a few additions.

- 1 *Alkali-Treated Peptone* To a solution of 50 gm peptone in 300 ml of distilled water is added 25 gm of NaOH also dissolved in 300 ml of distilled water. The solution is placed under ultraviolet light and allowed to stand 24 hours. The pH is adjusted to 6.6-6.8 with glacial acetic acid and 14 gm of sodium acetate are added and the volume made up to 1 liter. This solution contains 5% peptones and 6% sodium acetate. Store in a refrigerator. Toluene and chloroform are added.
- 2 *Vitamin Solution* Incorporate 5 mg thiamine hydrochloride 5 mg p-aminobenzoic acid 10 mg nicotinic acid 10 mg pyridoxine hydrochloride 1 ml of salt solution A and dilute to 100 ml. Add chloroform and toluene and store in refrigerator.
- 3 *Standard Solution of Pantothenic Acid* A solution of calcium pantothenate 1 mg per ml, is made up by dissolving 10-20 mg of the crystalline vitamin in the proper amount of water. Of this solution 54 ml plus 1 ml salt solution A are made up to 100 ml in a volumetric flask to give a concentration of 50 mcg per ml. Dilute to 10 millimicrograms per ml for the working standard.
- 4 *Stock Culture* The organism used is *Lactobacillus arabinosus* 175 (American Type Culture Collection No. 8014). It is maintained on yeast agar slabs.

## Procedure

- 1 *Preparation of the Inoculum* Prepare the inoculum for the assay tubes by transferring from the stock culture to a sterile tube of yeast dextrose broth supplemented with 10 millimicrograms of calcium pantothenate. Incubate for 24 hours at 37.5 C. Centrifuge and discard supernatant liquid. Resuspend the cells in 10 ml of physiological saline. One drop of this cell suspension is added to 10 ml saline and one drop of this resuspended cell suspension is added per assay tube.
- 2 *Preparation of Samples*
  - a. Food About 1 gm food is treated with 1 ml 2% sakadiastase solution, 25 ml acetate buffer pH 4.5 and 1 ml toluene. Plug with cotton, incubate at 37 C for 18-24 hours. Adjust pH to 6.6-8.0, steam 10 minutes in autoclave, dilute and filter. For use dilute 1 to 100.
  - b. Feces About 0.8 gm feces are treated with enzyme as above. final dilution is 1:2000.
  - c. Whole Blood, 10 ml blood is treated as above. final dilution is 1:40.
  - d. Plasma 10 ml plasma enzyme is treated as above. final dilution is 1:20.
  - e. Urine Dilute 1 ml to 100 with water.
- 3 *Preparation of Basal Medium* Prepare the basal medium having the composition below.
  - a. Mix together
    - 0 gm glucose
    - 10.6 gm sodium acetate
    - 50 gm l-tryptophan (2 mg per ml)
    - 10 ml l-cystine (4 mg per ml)
    - 100 ml peptone alkali treated (5% peptone 6% NaAc)
    - 10 ml adenine guanine uracil (1 mg each per ml)
    - 50 ml casein hydrolysate 10%
    - 5 ml salt solution A
    - 5 ml salt solution B
    - 3 ml vitamin solution (thiamin and PABA 50 mcg per ml, niacin and pyridoxine 100 mcg per ml)
    - 3 ml biotin (0.1 mcg per ml)
    - 8 ml riboflavin (2 mcg per ml)
  - b. Adjust pH to 6.7 and dilute to 500 ml.
- 4 *Preparation of Standard Tubes* Same as for pteroylglutamic acid assay using standard calcium pantothenate solution (10 millimicrograms per ml) for the working standard.
- 5 *Preparation of assay tubes* sterilization, inoculation and incubation treatment and calculations are by the same general microbiological techniques as were used in pteroylglutamic acid assay.

## Precautions

- 1 *Laetobacillus arabinosus* does not respond to bound pantothenic acid and since the vitamin is labile to heat in acid or alkaline solution digestion of the sample must be accomplished at neutral pH by means of digestive enzymes.

5 Pyridoxine (Vitamin B<sub>6</sub>)

## Reference

Modified Microbiological Method of Atkin L. Schultz, A. S. Williams, W. L. and Frey, C. N. Yeast Microbiological Methods for Determination of Vitamins Ind and Eng Chem Anal Ed 15 141-144 (Feb.) 1943

**Principle**

The growth stimulation of a strain of the yeast *Saccharomyces carlsbergensis* is used for assaying pyridoxine microbiologically. Growth is estimated turbidimetrically.

**Apparatus**

- 1 Same as for pteroylglutamic acid assay
- 2 Photoelectric spectrophotometer Coleman Jr
- 3 Mechanical shaker

**Reagents**

Same as for pteroylglutamic acid assay with some additions

- 1 *Sugar and Salt Solution* 200 gm CP dextrose 2.2 gm monopotassium phosphate 17 gm potassium chloride 0.5 gm calcium chloride 0.5 gm magnesium sulfate 0.01 gm ferric chloride and 0.01 gm manganese sulfate made up to 1 liter with water
- 2 *Potassium Citrate Buffer* 100 gm potassium citrate and 20 gm citric acid made up to one liter with water
- 3 *Casein Hydrolysate Solution* 10% Adjust the pH to 4.6
- 4 *Thiamine Solution* 10 micrograms per ml
- 5 *Inositol Solution* 1 mg per ml
- 6 *Biotin Solution* 0.1 microgram per ml
- 7 *Calcium Pantothenate* 400 micrograms per ml (0.2 gm per liter with 10 ml salt solution A per liter)
- 8 *Bacto Malt Extract* (Difco Bacto) and fresh moist baker's yeast
- 9 0.44 N and 0.55 N  $\text{H}_2\text{SO}_4$
- 10 *Standard Pyridoxine Hydrochloride* A solution of pyridoxine hydrochloride 1 mg per ml is made up by dissolving 15.0 mg of crystalline vitamin in water. 1 ml of this solution plus 1 ml salt solution A are made up to 100 ml in a volumetric flask to give a concentration of 100 mg of pyridoxine per ml. Two working standards are used 0.01 mcg per ml and 0.005 mcg per ml
- 11 *Stock Culture*  
*Saccharomyces carlsbergensis* American Type Culture Collection No 4008, is carried on Difco malt agar slants. Fresh slants are incubated for 24 hours at 30 C and stored in a refrigerator for no more than two weeks.

**Procedure**

- 1 *Preparation of the Yeast Inoculum* Prepare a fresh slant of *Saccharomyces carlsbergensis* on Difco malt agar and incubate 24 hours at 30 C. Remove a quantity of fresh growth with a sterile wire loop and suspend in 10 ml of sterile isotonic saline in a colorimeter tube. With the aid of a photoelectric colorimeter adjust the concentration of the yeast to an equivalent of 1 mg of moist baker's yeast per ml by adding sterile saline. Dilute 10 ml of the adjusted suspension with 90 ml of saline in a sterile Erlenmeyer flask providing a yeast suspension with final concentration of 0.1 mg per ml.
- 2 *Preparation of Samples*
  - a *Food*—Extract 1 gm of the specimen with 90 ml 0.44 N  $\text{H}_2\text{SO}_4$  by autoclaving at 15 pounds pressure for one hour. Cool, neutralize to pH 5.52, dilute to 100 ml and filter.
  - b *Feces*—Extract 0.5 gm of specimen with 90 ml 0.44 N  $\text{H}_2\text{SO}_4$  by autoclaving at 15 lbs pressure for one hour. Cool, neutralize to pH 5.52, dilute to 100 ml and filter.
  - c *Urine*—Extract 5 ml of specimen with 80 ml 0.05 N  $\text{H}_2\text{SO}_4$  by autoclaving at 15 lbs. pressure for 1 hour. Cool, neutralize at pH 5.52, dilute to 100 ml and filter.
- 3 *Preparation of Basal Medium* Prepare the basal medium having the composition below.
 

|                                   |        |
|-----------------------------------|--------|
| a. Mix together                   |        |
| Sugar salt solution               | 250 ml |
| Potassium citrate buffer solution | 50 ml  |

|                               |                                |
|-------------------------------|--------------------------------|
| Casein hydrolysate            | 50 ml                          |
| Thiamine solution             | 5 ml                           |
| Inositol solution             | 5 ml                           |
| Biotin solution               | 8 ml                           |
| Calcium pantothenate solution | 125 ml or 50 ml P. A. standard |

b Dilute to 500 ml

- 4 *Preparation of Standard Tubes* Place 5 ml of basal medium in each of a series of 18 mm Pyrex test tubes. Into duplicate tubes pipette 0.0, 0.5, 1.0, 1.5, 2.0, 2.5, 3.0, 3.5 and 4.0 ml of diluted standard. Adjust volume to 9.0 ml with water.
- 5 *Preparation of Assay Tubes* Pipette into successive tubes 0.5, 1.0, 2.0, 3.0 and 4.0 ml of the extract of specimen. Adjust the volume to 4.0 ml and add 5.0 ml of the basal medium.
- 6 *Sterilization* The tubes are covered with aluminum covers, steamed for 10 minutes and cooled.
- 7 *Inoculation and Incubation* Under aseptic conditions introduce into each tube 1 ml of the yeast inoculum. Place the tubes in a mechanical shaker for 16-18 hours at 30° C (rate of about 100 strokes a minute).
- 8 *Estimation* Immediately after incubation estimate the yeast growth turbidimetrically in a photoelectric spectrophotometer. Set the wave length at 680 mμ and 100% T with distilled water. The yeast growth inoculum should read between 80 and 90% T.

#### Calculation

- 1 On ordinary graph paper plot the results of the reference series in per cent transmittance against concentration of pyridoxine hydrochloride. Estimate the values for the unknowns from the graph. Calculate the pyridoxine content of the specimen for each tube and average all values which agree within 10% of their mean following the same general formula as for pteroylglutamic acid.

#### Precautions

- 1 The vitamin B group consists of pyridoxine, pyridoxal and one or more labile factors. All members of the complex show comparable activity for *Saccharomyces carlsbergensis*. Hence this procedure is not specific for pyridoxine.
- 2 Careful attention should be given to the hydrolytic procedures employed to liberate the bound forms of the vitamin; it is extracted from most plant and animal tissues with difficulty.
- 3 The chief disadvantage of this method is its low degree of precision. On the average values may be regarded as reproducible to within 10% but frequently much greater variations are encountered for unexplained reasons.

### Riboflavin

#### References

- 1 Modified Microbiological Method of Snell E. E. and Strong F. M. A Microbiological Assay for Riboflavin. *Ind. and Eng. Chem. Anal. Ed.* 11: 346-350 (June) 1939.
- 2 Scott M. L., Randall F. E. and Hessel F. H. A Modification of the Snell and Strong Microbiological Method for Determining Riboflavin. *J. Biol. Chem.* 141: 53-56 (Oct.) 1941.
- 3 Strong F. M. and Carpenter L. F. Preparation of Samples for Microbiological Determination of Riboflavin. *Ind. and Eng. Chem. Anal. Ed.* 14: 909-913 (Nov) 1942.

#### Principle

Riboflavin is determined by measurement of the growth stimulation of *Lactobacillus casei*. The acid produced by the microorganism is determined by titration with sodium hydroxide.



## 254 Microbiological Procedures

### Apparatus

- 1 Same as pteroylglutamic acid assay

### Reagents

Reagents are the same as for pteroylglutamic acid assay with a few additions

- 1 *Alkali Treated Peptone Solution* Dissolve 40 gm of peptone (Difco Bacto or Wilson's) in 250 ml of hot distilled water. Add 40 gm of NaOH in 250 ml of distilled water. Mix the two solutions. Allow to stand for 18<sup>24</sup> hours then neutralize with glacial acetic acid. Add 14 gm of dry  $\text{Na}_2\text{C}_2\text{H}_3\text{O}_2$  (or 3.2 gm  $\text{Na}_2\text{C}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$ ) and sufficient water to make up to 800 ml. Preserve under chloroform and toluene in a refrigerator.
- 2 *Yeast Supplement Solution* Dissolve 100 gm of yeast extract (Difco Bacto) in 500 ml of water and 150 gm of basic lead acetate in 500 ml of water. Mix the two solutions. Adjust the pH to 10 with concentrated  $\text{NH}_4\text{OH}$  and filter or centrifuge. Adjust the filtrate to pH 6.5 with glacial acetic acid. Precipitate the excess lead with  $\text{H}_2\text{S}$  and add distilled water to make 1000 ml. Preserve under toluene in the refrigerator.
- 3 *Cystine Solution* Dissolve 1 gm l-cystine in 20 ml of 10% hydrochloric acid and dilute to 1000 ml with water. Store under toluene in refrigerator not below 10° C.
- 4 *Reference Standard Solution* Dissolve 50 mg of U.S.P. riboflavin reference standard in 1500 ml of distilled water containing 4 ml of glacial acetic acid and warm to aid solution. After cooling to room temperature make to 2000 ml volume with water. Preserve in the dark under toluene in a refrigerator. Dilute to 0.1 mcg of riboflavin per ml immediately before use for reference working standard.

### Stock Culture

- 1 Same as for pteroylglutamic acid assay

### Procedure

- 1 *Preparation of the Inoculum* On the day prior to use transfer cells from the stock culture to a sterile tube of yeast dextrose broth supplemented with 0.5 mcg of riboflavin/ml and incubate 24 hours at 37.5° C. For the inoculum the cells are centrifuged out and the supernatant liquid decanted. Suspend the cells in 10 ml of sterile isotonic saline solution and use 1 drop of this cell suspension.
- 2 *Preparation of Samples*
  - a. *Foods* To sample (10 or more micrograms riboflavin) add 50 ml 0.1 N HCl and autoclave for 15 minutes at 15 lbs pressure. Add 1 ml of freshly prepared papain-takadiastase solution (20 mg/ml). Incubate at 37.5° C for 4 hours. Heat with flowing steam for 10 minutes. Adjust to pH 4.5 make to 100 ml volume and filter.
  - b. *Urine* Dilute with water so that the specimen contains between 0.05 and 0.5 mcg riboflavin per ml.
- 3 *Preparation of Basal Medium* Prepare a basal medium having the composition below
  - a. *Mix together*

|                        |       |
|------------------------|-------|
| Alkali treated peptone | 50 ml |
| 0.1% cystine solution  | 50 ml |
| Yeast supplement       | 5 ml  |
| Salt solution A        | 5 ml  |
| Salt solution B        | 25 ml |
| Dextrose anhydrous     | 10 gm |
  - b. Adjust to pH 6.6-6.8. Dilute to 50 ml

- 4 *Preparation of Standard Tubes* To duplicate tubes 16 by 150 mm in size add 0.05, 0.1, 0.15, 0.2, 0.3 and 0.5 ml of standard riboflavin solution (0.1 mcg per ml). Add sufficient water to bring the volume in each tube to 5.0 ml. To each of these tubes add 5 ml of basal medium.
- 5 *Preparation of Assay Tubes* To duplicate tubes add 0.5, 1.0, 1.5 and 2.0 ml of the extract of the test material. Add sufficient water to bring the volume in each tube to 5.0 ml and add 5.0 ml of basal medium.
- 6 *Sterilization, inoculation and incubation* are performed as in the pteroylglutamic acid assay.
- 7 *Titration* Transfer the contents of each tube to an Erlenmeyer flask using 10 ml of distilled water for rinsing. Titrate with 0.1 N NaOH using bromthymol blue as an indicator.

#### Calculation

- 1 The same type of calculation is used as for pteroylglutamic acid assay.

## 7 Tryptophane

#### References

- 1 Modified Microbiological Method of Snell E. E. and Wright L. D. A Microbiological Method for the Determination of Nicotinic Acid *J Biol Chem* 139: 675-685 (June) 1941.
- 2 Greene R. D. and Black A. The Microbiological Assay of Tryptophane in Protein and Foods *J Biol Chem* 155: 18 (Sept.) 1944.
- 3 Wooley J. G. and Sebrell W. H. Two Microbiological Methods for the Determination of Tryptophane in Proteins and Other Complex Substances *J Biol Chem* 157: 141-151 (Jan.) 1945.

#### Principle

The growth stimulation of *Lactobacillus arabinosus* 175 in relation to the tryptophane content of a basal medium is employed in the microbiological assay for this amino acid. The lactic acid produced is titrated.

#### Apparatus

Same as for pteroylglutamic acid.

#### Reagents

These are the same as for pteroylglutamic acid with some additions:

- 1 *Standard Tryptophane Solution* Dissolve 0.40 mg of L-tryptophane in the proper amount of hot water to provide a solution containing 2 mg tryptophane per ml. Incorporate 10 ml of the solution with 1 ml salt solution A and make up to 100 ml volume to give 0.02 mcg tryptophane per ml. Dilute to 1 mcg per ml for working standard.
- Vitamin Solution* Same as for pantothenic acid assay.
- 3 *Stock Culture*  
Same as for nicotinic acid assay.

#### Procedure

- 1 *Preparation of Inoculum* Same as for nicotinic acid assay.
- 2 *Preparation of Samples* Either of two methods may be employed.
  - a. *Use of Enzymes* Place sample (5 gm food, 0.5 ml blood, 0.5 ml plasma, 1 ml urine) in 100 ml beaker. Add 5 ml 0.1 N  $H_2SO_4$ , 0.5 ml 1% pepsin solution and incubate 18-4 hours at 37.5°C with occasional mixing. Add 0.6 gm  $K_2HPO_4$  and adjust pH to 8.4. Add 2.5 ml 1% trypsin solution, 10 ml toluene and a few drops of chloroform. Incubate at 18°C 4 hours at 37°C in 15 ml Erlenmeyer flask with occasional mixing. Adjust pH to 8.8, transfer to a 50 ml Erlenmeyer flask and add 5 ml 1% chymotrypsin solution, a few drops of chloroform and 10 ml toluene.

Incubate 72 hours at 37 C with occasional mixing. Replace toluene as necessary. Adjust pH to 6.0-6.5, including blank, dilute to 200 ml and filter until clear using same filter paper for repeated filtrations. Make final dilution as follows: blood 1 to 3000, plasma 1 to 1000, diet 1 to 400, enzyme blank, 1 to 1000 on basis of pepsin, urine 1 to 20.

- b. Use of Barium Hydroxide. Place sample containing about 250 mcg tryptophane in a 250 ml Erlenmeyer flask (diet 1 to 2 gm, blood 0.1 ml, plasma 0 ml, urine 0.2 ml). Add 6.5 gm barium hydroxide  $\text{Ba}(\text{OH})_2 \cdot 8\text{H}_2\text{O}$  and 0 ml water mix thoroughly and autoclave 16 hours at 15 lbs. Autoclaved samples are transferred quantitatively while hot to 250 ml beakers. After cooling the pH is made acid to 2.5 with 2 N  $\text{H}_2\text{SO}_4$ . Samples are neutralized to pH 6.0-6.5 with NaOH and diluted to volume plus 1 ml to allow for  $\text{BaSO}_4$  formed. Final dilutions are as follows: food 1 to 50 or 1 to 150 depending on tryptophane content, blood 1 to 1000, plasma 1 to 500, urine 1 to 100. Heat to 60-70 C, filter and assay. *Note:* Because of racemization when barium hydroxide is used multiply all assay values by 2.

3. Preparation of Basal Medium. The basal medium recommended here differs from that of Snell and Wright's for niacin only in the omission of tryptophane and the inclusion of niacin.

a. Mix together

- 20 gm glucose
- 33.4 gm sodium acetate ( $\text{Na}_2\text{C}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$ )
- 50 ml casein hydrolysate 10%
- 50 ml L-cystine (4 mg per ml)
- 10 ml adenine, guanine, uracil (1 mg each per ml)
- 50 ml salt solution A
- 50 ml salt solution B
- 30 ml vitamin solution (thiamin and PABA 50 mcg per ml, niacin and pyridoxine 100 mcg per ml)
- 30 ml biotin (0.1 mcg per ml)
- 80 ml riboflavin (25 mcg per ml)
- 30 ml pantothenic acid solution (50 mcg  $\text{P}_4$  per ml)

- b. Adjust to pH 6.6-6.8 and dilute to 500 ml

4. Preparation of Standard Tubes. Same as for riboflavin as above except for using standard tryptophane (1 mcg per ml) for the working standard.
5. Preparation of assay tubes: sterilization, inoculation and incubation, titration and calculation are by the same general procedures used in riboflavin assay.

#### Precautions

1. *Lactobacillus arabinosus* does not respond to D-tryptophane and complete racemization of hydrolysates prepared with base is assumed. Therefore direct analytical results are obtained with a D standard. If L-tryptophane is used as a standard the analytical results should be multiplied by 2. However recent investigations indicate that there is a good reason to believe that there is destruction of tryptophane after treatment with alkali. Cystine added as a stabilizing agent affords good results.

## 8. Agar Plate Assay for Biotin

#### References

1. Genghof, D. S., Partridge, C. W. H. and Carpenter, F. H. An Agar Plate Assay for Biotin. *Arch. Biochem.* 17: 413-420 (June) 1948.
2. Snell, E. E. and Wright, L. D. A Microbiological Method for the Determination of Nicotinic Acid. *J. Biol. Chem.* 139: 615-686 (June) 1941.

# Principle

Biotin solutions of increasing concentrations are added to filter paper discs on a biotin free but otherwise complete medium seeded with *Lactobacillus arabinosus* 17 5. Graded zones of growth are produced. A linear relationship exists between the diameter of the zone of growth and the logarithm of the dose of biotin.

# Apparatus

- 1 Pen cillin reader (Fisher)
- 2 Petri dishes (100 mm diameter) with unglazed porcelain tops (Coors)
- 3 Filter paper discs (Schleicher and Schuell No 740 E)
- 4 Incubator set at 30 C
- 5 Autoclave
- 6 Refrigerator
- 7 Constant temperature bath set at 45 C
- 8 Pipettes (0.1 ml graduated in hundredths)
- 9 Centrifuge
- 10 Brewer automatic pipette to deliver 16 ml amounts (Manual delivery of 16 ml is time consuming)
- 11 Graph paper semi log

# Reagents

- 1 Water glass distilled
- 2 Biotin ampoules or crystalline
- 3 Casein hydrolysate 10%
- 4 L (-) Cystine
- 5 Tryptophane
- 6 Sodium acetate
- 7 Glucose
- 8 Potassium monohydrogen phosphate
- 9 Potassium dihydrogen phosphate
- 10 Sodium chloride
- 11 Manganese sulfate (tetrahydrate)
- 12 Magnesium sulfate (heptahydrate)
- 13 Ferrous sulfate (heptahydrate)
- 14 Guanine hydrochloride
- 15 Adenine sulfate
- 16 Uracil
- 17 Riboflavin
- 18 Nicotinic acid
- 19 Thiamine
- 20 Pantothenic acid
- 21 Pyridoxine
- 22 Para amino benzoic acid (PABA)
- 23 Agar
- 24 Stock Culture

AGU (1 mg each per ml)

The organism used is *Lactobacillus arabinosus* 17 5 (American Type Culture Collection Number 6014). It is maintained on yeast agar slabs.

# Procedure

- 1 Sterilize Petri dishes at 165 C for at least 90 minutes in oven
- Inoculate a culture of *L. arabinosus* 17 5 into 10 ml of sterile basal medium plus 0.0002 gamma biotin per tube. After incubation at 30 C for 18-20 hours the organisms are centrifuged and washed twice with 10 ml of sterile saline and resuspended in 5 ml of saline.

Incubate 72 hours at 37 C with occasional mixing. Replace toluene as necessary. Adjust pH to 6.0-6.5, including blank dilute to 200 ml, and filter until clear using same filter paper for repeated filtrations. Make final dilution as follows: blood 1 to 3000, plasma 1 to 2000, diet 1 to 400, enzyme blank 1 to 1000 on basis of pepsin, urine 1 to 20.

- b *Use of Barium Hydroxide* Place sample containing about .50 mcg tryptophane in a .50 ml Erlenmeyer flask (diet 1 to 2 gm, blood 0.1 ml, plasma 0.1 ml, urine, 0.2 ml). Add 6.5 gm barium hydroxide  $\text{Ba}(\text{OH})_2 \cdot 8\text{H}_2\text{O}$  and 20 ml water, mix thoroughly and autoclave 16 hours at 15 lbs. Autoclaved samples are transferred quantitatively while hot to 250 ml beakers. After cooling the pH is made acid to 2.5 with 2 N  $\text{H}_2\text{SO}_4$ . Samples are neutralized to pH 6.0-6.5 with NaOH and diluted to volume plus 1 ml to allow for  $\text{BaSO}_4$  formed. Final dilutions are as follows: food 1 to 50 or 1 to 150 depending on tryptophane content, blood 1 to 1000, plasma 1 to 500, urine 1 to 100. Heat to 60-70 C, filter and assay. *Note* Because of racemization when barium hydroxide is used multiply all assay values by 2.

- 3 *Preparation of Pabal Medium* The basal medium recommended here differs from that of Snell and Wright's for niacin only in the omission of tryptophane and the inclusion of niacin.

a *Mix together*

20 gm glucose

33.4 gm sodium acetate ( $\text{NaC}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$ )

50 ml casein hydrolysate 10%

50 ml L-cystine (4 mg per ml)

10 ml adenine, guanine, uracil (1 mg each per ml)

50 ml salt solution A

50 ml salt solution B

30 ml vitamin solution (thiamin and PABA 50 mcg per ml, niacin and pyridoxine 100 mcg per ml)

30 ml biotin (0.1 mcg per ml)

80 ml riboflavin (25 mcg per ml)

30 ml pantothenic acid solution (50 mcg PA per ml)

b *Adjust to pH 6.6-6.8 and dilute to 100 ml*

- 4 *Preparation of Standard Tubes* Same as for riboflavin assay except for using standard tryptophane (1 mg per ml) for the working standard.

- 5 *Preparation of assay tubes* sterilization, inoculation and incubation, titration and calculation are by the same general procedures used in riboflavin assay.

**Precautions**

- 1 *Lactobacillus arabinosus* does not respond to D-tryptophane and complete racemization of hydrolysates prepared with base is assumed. Therefore direct analytical results are obtained with a DL standard. If L-tryptophane is used as a standard the analytical results should be multiplied by 2. However recent investigations indicate that there is a good reason to believe that there is destruction of tryptophane after treatment with alkali. Cystine added as a stabilizing agent affords good results.

## 8 Agar Plate Assay for Biotin

**References**

- 1 Genghof D S, Partridge C W H and Carpenter F H. An Agar Plate Assay for Biotin. *Arch Biochem* 17: 413-40 (June) 1949.
- 2 Snell E E and Wright L D. A Microbiological Method for the Determination of Nicotinic Acid. *J Biol Chem* 139: 675-686 (June) 1941.

## Principle

Biotin solutions of increasing concentrations are added to filter paper discs on a biotin free but otherwise complete medium seeded with *Lactobacillus arabinosus* 175. Graded zones of growth are produced. A linear relationship exists between the diameter of the zone of growth and the logarithm of the dose of biotin.

## Apparatus

- 1 Penicillin reader (Fisher)
- 2 Petri dishes (100 mm diameter) with unglazed porcelain tops (Coors)
- 3 Filter paper discs (Schleicher and Schuell No. 740 E)
- 4 Incubator set at 30 C
- 5 Autoclave
- 6 Refrigerator
- 7 Constant temperature bath set at 45 C
- 8 Pipettes (0.1 ml graduated in hundredths)
- 9 Centrifuge
- 10 Brewer automatic pipette to deliver 16 ml amounts (Manual delivery of 16 ml is time consuming)
- 11 Graph paper semi log

## Reagents

- 1 Water glass distilled
- 2 Biotin ampoules or crystalline
- 3 Casein hydrolysate 10%
- 4 1 (-) Cystine
- 5 Tryptophane
- 6 Sodium acetate
- 7 Glucose
- 8 Potassium monohydrogen phosphate
- 9 Potassium dihydrogen phosphate
- 10 Sodium chloride
- 11 Manganese sulfate (tetrahydrate)
- 12 Magnesium sulfate (heptahydrate)
- 13 Ferrous sulfate (heptahydrate)
- 14 Guanine hydrochloride
- 15 Adenine sulfate
- 16 Uracil
- 17 Riboflavin
- 18 Niacin
- 19 Thiamine
- 20 Pantothenic acid
- 21 Pyridoxine
- 22 Para amino benzoic acid (PABA)
- 23 Agar
- 24 Stock Culture

} AGU (1 mg each per ml)

The organism used is *Lactobacillus arabinosus* 175 (American Type Culture Collection Number 8014). It is maintained on yeast agar slabs.

## Procedure

- 1 Sterilize Petri dishes at 165 C for at least 90 minutes in an oven.
- Inoculate a culture of *L. arabinosus* 175 into 10 ml of sterile basal medium plus 0.000 gamma biotin per tube. After incubation at 30 C for 18-20 hours the organisms are centrifuged and washed twice with 10 ml of sterile saline and resuspended in 5 ml of saline.

Incubate 72 hours at 37 C with occasional mixing. Replace toluene as necessary. Adjust pH to 6.0-6.5 including blank, dilute to 200 ml and filter until clear, using same filter paper for repeated filtrations. Make final dilution as follows: blood, 1 to 3000; plasma, 1 to 2000; diet, 1 to 400; enzyme blank, 1 to 1000 on basis of pepsin; urine, 1 to 20.

- b. Use of Barium Hydroxide. Place sample containing about 250 mcg tryptophane in a 250 ml Erlenmeyer flask (diet, 1 to 2 gm; blood, 0.1 ml; plasma, 0.2 ml; urine, 0.2 ml). Add 6.5 gm barium hydroxide  $\text{Ba}(\text{OH})_2 \cdot 8\text{H}_2\text{O}$  and 70 ml water, mix thoroughly and autoclave 16 hours at 15 lbs. Autoclaved samples are transferred quantitatively while hot to 250 ml beakers. After cooling the pH is made acid to 2.5 with 2 N  $\text{H}_2\text{SO}_4$ . Samples are neutralized to pH 6.0-6.5 with NaOH and diluted to volume plus 1 ml to allow for  $\text{BaSO}_4$  formed. Final dilutions are as follows: food, 1 to 50 or 1 to 100 depending on tryptophane content; blood, 1 to 1000; plasma, 1 to 500; urine, 1 to 100. Heat to 60.0 C, filter and assay. *Note:* Because of racemization when barium hydroxide is used multiply all assay values by 2.

3. *Preparation of Basal Medium.* The basal medium recommended here differs from that of Snell and Wright's for niacin only in the omission of tryptophane and the inclusion of niacin.

a. Mix together

- 20 gm glucose
- 33.4 gm sodium acetate ( $\text{NaC}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$ )
- 50 ml casein hydrolysate 10%
- 50 ml L-cystine (4 mg per ml)
- 10 ml adenine, guanine, uracil (1 mg each per ml)
- 50 ml salt solution A
- 50 ml salt solution B
- 30 ml vitamin solution (thiamin and PABA 50 mcg per ml; niacin and pyridoxine 100 mcg per ml)
- 30 ml biotin (0.1 mcg per ml)
- 80 ml riboflavin (25 mcg per ml)
- 30 ml pantothenic acid solution (50 mcg PA per ml)

- b. Adjust to pH 6.6-6.8 and dilute to 500 ml

4. *Preparation of Standard Tubes.* Same as for riboflavin assay except for using standard tryptophane (1 mcg per ml) for the working standard.
5. *Preparation of assay tubes.* Sterilization, inoculation and incubation, titration and calculation are by the same general procedures used in riboflavin assay.

### Precautions

1. *Lactobacillus arabinosus* does not respond to D-tryptophane and complete racemization of hydrolyzates prepared with base is assumed. Therefore direct analytical results are obtained with a DL standard. If L-tryptophane is used as a standard the analytical results should be multiplied by 2. However, recent investigations indicate that there is a good reason to believe that there is destruction of tryptophane after treatment with alkali. Cystine added as a stabilizing agent affords good results.

## 8. Agar Plate Assay for Biotin

### References

1. Genghof, D. S., Partridge, C. W. H. and Carpenter, F. H. An Agar Plate Assay for Biotin. *Arch. Biochem.* 17: 413-420 (June) 1948.
2. Snell, E. E. and Wright, L. D. A Microbiological Method for the Determination of Nicotinic Acid. *J. Biol. Chem.* 139: 675-686 (June) 1941.

**Principle**

Biotin solutions of increasing concentrations are added to filter paper discs on a biotin free but otherwise complete medium seeded with *Lactobacillus arabinosus* 175. Graded zones of growth are produced. A linear relationship exists between the diameter of the zone of growth and the logarithm of the dose of biotin.

**Apparatus**

- 1 Penicillin reader (Fisher)
- 2 Petri dishes (100 mm diameter) with unglazed porcelain tops (Coors)
- 3 Filter paper discs (Schleicher and Schnell No 740 E)
- 4 Incubator set at 30 C
- 5 Autoclave
- 6 Refrigerator
- 7 Constant temperature bath set at 45 C
- 8 Pipettes (0.1 ml graduated in hundredths)
- 9 Centrifuge
- 10 Brewer automatic pipette to deliver 16 ml amounts (Manual delivery of 16 ml is time consuming)
- 11 Graph paper semi-log

**Reagents**

- 1 Water glass distilled.
- 2 Biotin ampoules or crystalline
- 3 Casein hydrolysate 10%
- 4 L (-) Cystine
- 5 Tryptophane
- 6 Sodium acetate
- 7 Glucose
- 8 Potassium monohydrogen phosphate
- 9 Potassium dihydrogen phosphate
- 10 Sodium chloride
- 11 Manganese sulfate (tetrahydrate)
- 12 Magnesium sulfate (heptahydrate)
- 13 Ferrous sulfate (heptahydrate)
- 14 Guanine hydrochloride
- 15 Adenine sulfate
- 16 Uracil
- 17 Riboflavin
- 18 Nicotin
- 19 Thiamine
- 20 Pantothenic acid
- 21 Pyridoxine
- 22 Para amino benzoic acid (PABA)
- 23 Agar
- 24 Stock Culture

} AGU (1 mg each per ml)

The organism used is *Lactobacillus arabinosus* 175 (American Type Culture Collection Number 6014). It is maintained on yeast agar slabs.

**Procedure**

- 1 Sterilize Petri dishes at 163 C for at least 90 minutes in oven.
- 2 Inoculate a culture of *L. arabinosus* 175 into 10 ml of sterile basal medium plus 0.0006 gamma biotin per tube. After incubation at 30 C for 18-20 hours the organisms are centrifuged and washed twice with 10 ml of sterile saline and resuspended in 5 ml of saline.



Incubate 72 hours at 37 C with occasional mixing. Replace toluene as necessary. Adjust pH to 6.0-6.5 including blank dilute to 200 ml and filter until clear using same filter paper for repeated filtrations. Make final dilution as follows: blood 1 to 3000, plasma, 1 to 2000, diet, 1 to 400, enzyme blank 1 to 1000 on basis of pepsin urine 1 to 20.

- b. Use of Barium Hydroxide. Place sample containing about 250 mcg tryptophane in a 250 ml Erlenmeyer flask (diet 1 to 2 gm, blood 0.1 ml, plasma 0.2 ml, urine 0.2 ml). Add 6.5 gm barium hydroxide,  $\text{Ba}(\text{OH})_2 \cdot 8\text{H}_2\text{O}$  and 70 ml water mix thoroughly and autoclave 16 hours at 15 lbs. Autoclaved samples are transferred quantitatively while hot to 200 ml beakers. After cooling the pH is made acid to 2.5 with 2 N  $\text{H}_2\text{SO}_4$ . Samples are neutralized to pH 6.0-6.5 with NaOH and diluted to volume plus 1 ml to allow for  $\text{BaSO}_4$  formed. Final dilutions are as follows: food 1 to 50 or 1 to 150 depending on tryptophane content, blood 1 to 1000, plasma 1 to 500, urine 1 to 100. Heat to 60-70 C, filter and assay. *Note:* Because of racemization when barium hydroxide is used multiply all assay values by 2.

3. *Preparation of Basal Medium.* The basal medium recommended here differs from that of Snell and Wright's for niacin only in the omission of tryptophane and the inclusion of niacin.

a. Mix together

20 gm glucose

33.4 gm sodium acetate ( $\text{NaCH}_3\text{CO}_2 \cdot 3\text{H}_2\text{O}$ )

50 ml casein hydrolysate 10%

50 ml L-cystine (4 mg per ml)

10 ml adenine, guanine, uracil (1 mg each per ml)

50 ml salt solution A

50 ml salt solution B

30 ml vitamin solution (thiamin and PABA 50 mcg per ml, niacin and pyridoxine 100 mcg per ml)

30 ml biotin (0.1 mcg per ml)

80 ml riboflavin (0.5 mcg per ml)

30 ml pantothenic acid solution (50 mcg PA per ml)

- b. Adjust to pH 6.6-6.8 and dilute to 500 ml

4. *Preparation of Standard Tubes.* Same as for riboflavin assay except for using standard tryptophane (1 mcg per ml) for the working standard.

5. *Preparation of assay tubes:* sterilization, inoculation and incubation, titration and calculation are by the same general procedures used in riboflavin assay.

### Precautions

1. *Lactobacillus arabinosus* does not respond to D-tryptophane and complete racemization of hydrolysates prepared with base is assumed. Therefore diet analytical results are obtained with a D standard. If L-tryptophane is used as a standard the analytical results should be multiplied by 2. However, recent investigations indicate that there is a good reason to believe that there is destruction of tryptophane after treatment with alkali. Cystine added as a stabilizing agent affords good results.

## 8. Agar Plate Assay for Biotin

### References

1. Genghof, D. S., Partridge, C. W. H. and Carpenter, F. H. An Agar Plate Assay for Biotin. *Arch. Biochem.* 17: 413-40 (June) 1948.
2. Snell, E. E. and Wright, L. D. A Microbiological Method for the Determination of Nicotinic Acid. *J. Biol. Chem.* 139: 675-686 (June) 1941.

## Principle

Biotin solutions of increasing concentrations are added to filter paper discs on a biotin free but otherwise complete medium seeded with *Lactobacillus arabinosus* 175. Graded zones of growth are produced. A linear relationship exists between the diameter of the zone of growth and the logarithm of the dose of biotin.

## Apparatus

- 1 Petri dish reader (Fisher)
- 2 Petri dishes (100 mm diameter) with unglazed porcelain tops (Coors)
- 3 Filter paper discs (Schleicher and Schuell No. 740 E)
- 4 Incubator set at 30 C
- 5 Autoclave
- 6 Refrigerator
- 7 Constant temperature bath set at 45 C
- 8 Pipettes (0.1 ml graduated in hundredths)
- 9 Centrifuge
- 10 Bower automatic pipette to deliver 16 ml amounts (Manual delivery of 16 ml is time consuming)
- 11 Graph paper semi log

## Reagents

- 1 Water glass distilled
- 2 Biotin ampoules or crystalline
- 3 Casein hydrolysate 10%
- 4 L (-) Cystine
- 5 Tryptophane
- 6 Sodium acetate
- 7 Glucose
- 8 Potassium monohydrogen phosphate
- 9 Potassium dihydrogen phosphate
- 10 Sodium chloride
- 11 Manganese sulfate (tetrahydrate)
- 12 Magnesium sulfate (heptahydrate)
- 13 Ferrous sulfate (heptahydrate)
- 14 Guanine hydrochloride
- 15 Adenine sulfate
- 16 Uracil
- 17 Riboflavin
- 18 Nicotin
- 19 Thiamine
- 20 Pantothenic acid
- 1 Pyridoxine
- Para amino benzoic acid (PABA)
- 3 Agar
- 4 Stock Culture

} AGU (1 mg each per ml)

The organism used is *Lactobacillus arabinosus* 175 (American Type Culture Collection Number 3014). It is maintained on yeast agar slabs.

## Procedure

- 1 Sterilize Petri dishes at 160 C for at least 90 minutes in oven
- Inoculate a culture of *L. arabinosus* 175 into 10 ml of sterile basal medium plus 0.0001 gamma biotin per tube. After incubation at 30 C for 18-20 hours the organisms are centrifuged and washed twice with 10 ml of sterile saline and resuspended in 5 ml of saline.

Incubate 72 hours at 37 C with occasional mixing. Replace toluene as necessary. Adjust pH to 6.0-6.5, including blank, dilute to 200 ml and filter until clear, using same filter paper for repeated filtrations. Make final dilution as follows: blood, 1 to 3000; plasma, 1 to 2000; diet, 1 to 400; enzyme blank, 1 to 1000 on basis of pepsin; urine, 1 to 20.

- b **Use of Barium Hydroxide** Place sample containing about 250 mcg tryptophane in a 250 ml Erlenmeyer flask (diet 1 to 2 gm, blood, 0.1 ml, plasma 0.0 ml, urine 0.2 ml). Add 6.5 gm barium hydroxide  $\text{Ba}(\text{OH})_2 \cdot 8\text{H}_2\text{O}$  and 90 ml water mix thoroughly and autoclave 16 hours at 15 lbs. Autoclaved samples are transferred quantitatively while hot to 250 ml beakers. After cooling the pH is made acid to 2.5 with 2 N  $\text{H}_2\text{SO}_4$ . Samples are neutralized to pH 6.0-6.5 with NaOH and diluted to volume plus 1 ml to allow for  $\text{BaSO}_4$  formed. Final dilutions are as follows: food 1 to 100 or 1 to 150 depending on tryptophane content; blood 1 to 1000; plasma 1 to 500; urine 1 to 100. Heat to 60-70 C, filter and assay. *Note:* Because of racemization when barium hydroxide is used multiply all assay values by 2.

- 3 **Preparation of Basal Medium** The basal medium recommended here differs from that of Snell and Wright's for niacin only in the omission of tryptophane and the inclusion of niacin.

a **Mix together**

- 20 gm glucose
- 33.4 gm sodium acetate ( $\text{NaC}_2\text{H}_3\text{O}_2 \cdot 3\text{H}_2\text{O}$ )
- 50 ml casein hydrolysate 10%
- 50 ml L-cystine (4 mg per ml)
- 10 ml adenine, guanine, uracil (1 mg each per ml)
- 50 ml salt solution A
- 50 ml salt solution B
- 30 ml vitamin solution (thiamin and PABA 50 mcg per ml, niacin and pyridoxine 100 mcg per ml)
- 30 ml biotin (0.1 mcg per ml)
- 80 ml riboflavin (25 mcg per ml)
- 30 ml pantothenic acid solution (50 mcg PA per ml)

- b Adjust to pH 6.6-6.8 and dilute to 500 ml

- 4 **Preparation of Standard Tubes** Same as for riboflavin assay except for using standard tryptophane (1 mcg per ml) for the working standard.
- 5 **Preparation of assay tubes** sterilization, inoculation and incubation, titration and calculation are by the same general procedures used in riboflavin assay.

### Precautions

- 1 *Lactobacillus arabinosus* does not respond to *d*-tryptophane and complete racemization of hydrolysates prepared with base is assumed. Therefore direct analytical results are obtained with a *d*-standard. If *l*-tryptophane is used as a standard the analytical results should be multiplied by 2. However recent investigations indicate that there is a good reason to believe that there is destruction of tryptophane after treatment with alkali. Cystine added as a stabilizing agent affords good results.

## 8 Agar Plate Assay for Biotin

### References

- 1 Genghof, D. S., Partridge, C. W. H. and Carpenter, F. H. An Agar Plate Assay for Biotin. *Arch. Biochem.* 17: 413-420 (June) 1948.
- 2 Snell, E. E. and Wright, L. D. A Microbiological Method for the Determination of Nicotinic Acid. *J. Biol. Chem.* 139: 675-686 (June) 1941.



## 3 Prepare the following medium

|  |           |
|--|-----------|
| Casein hydrolysate 10%                                       | 50 ml     |
| 1 (-) Tryptophane  | 100 mg    |
| 1 (-) Cystine  | 100 mg    |
| Glucose  | 10 gm     |
| Na <sub>2</sub> C <sub>2</sub> H <sub>3</sub> O <sub>2</sub> | 6 gm      |
| AGU  | 10 mg     |
| Salt A   | 5 ml      |
| Salt B   | 5 ml      |
| Niacin   | 500 gamma |
| PABA   | 100 gamma |
| Calcium pantothenate   | 100 gamma |
| Pyridoxine hydrochloride                                     | 150 gamma |
| Thiamine   | 100 gamma |
| Riboflavin   | 200 gamma |
| Bacto agar   | 15 gm     |
| Adjust to pH 6.6-6.8 and dilute up to 1 liter                |           |

4 Autoclave the medium for 15 minutes at 120 C cool to 45 C and inoculate with resuspended cells (20 ml of the *Lactobacillus arabinosus* suspension per liter of medium)

## 5 To each sterile Petri dish add 16 ml of the above inoculated medium. Place plates in refrigerator as soon as the agar is hardened. These plates can be stored at 5 C for at least one week before use.

## 6 Standard Curve. Range 1 milligram to 1000 milligram of biotin. Biotin solutions of the following concentration are made

|                        |
|------------------------|
| 10 milligram per ml    |
| 25 milligram per ml    |
| 50 milligram per ml    |
| 100 milligram per ml   |
| 250 milligram per ml   |
| 500 milligram per ml   |
| 1000 milligram per ml  |
| 2500 milligram per ml  |
| 5000 milligram per ml  |
| 10000 milligram per ml |

## 7 Filter paper discs are placed on the agar and 0.08 ml of the biotin solutions delivered upon them. Run all dilutions in duplicate with 2 filter pads per plate through 10 milligram for other dilutions only 1 pad per plate.

## 8 Invert plates and incubate for 18-20 hours at 30 C.

## 9 After incubation the diameters of growth about the filter pads are measured to the closest 0.5 mm with the penicillin reader.

## 10 A standard curve is plotted on semi log graph paper. The mean diameter of the circle of growth is the abscissa the log of the concentration of biotin is the ordinate. A linear relationship exists between the diameter of the growth zone in mm and the log of the dose.

## 11 After urine samples are adjusted to pH 6.6-6.8 with 2 N HCl or 4 N NaOH and the volume change noted if any the samples are carried through the same procedure as the standard solutions of biotin.

**Calculation**

- 1 From the standard curve read directly the concentration of biotin in milligrams per ml. It is preferable to have no dilution of urine samples. If there is a volume change due to adjustment of pH multiply the direct reading by the dilution factor.

**Example**

- 1 On the standard curve 25 milligram of biotin per ml gave a reading of 77 mm on the Penicillin reader. A specimen of urine, undiluted produced a growth zone of 27 mm. Thus the sample had 250 milligram biotin per ml. If the sample had been diluted 9 ml of urine plus 1 ml of acid or base when adjusting pH and the zone of growth read 27 mm the sample tested would contain

$$250 \times \frac{10}{9} = 278 \text{ milligram/ml}$$

**Precautions**

- 1 Since the filter pads rapidly absorb water from the moist agar the solutions must be added to the pads as soon as possible after they have been placed upon the agar
- 2 The same volume of fluid must be delivered to each pad. Care is required at this step in the procedure
- 3 After inoculation of agar with organisms do not flame the flask each time before or after adding the 16 ml amounts to Petri dishes
- 4 Agar for inoculation must not be over 45 C nor cool enough that the agar starts to harden

**9 Microbiology Miscellaneous References****Microbiology—General**

Snell E E. Microbiological Methods in Vitamin Research. Look Vitamin Methods by Paul Gyorgy. Academic Press Inc. New York pages 3 303 1950

**Amino Acids**

Dunn M. S. Determination of Amino Acids by Microbiological Assay. Physiol. Rev. 29 219 250 (July) 1949

Guirard B M, Snell E E., and Williams R J. Microbiological Determination of Amino Acids. IV. Lysine, Histidine, Arginine and Valine. Proc. Soc. Exper. Biol. & Med. 61 (2) 153 161 (Feb.) 1946

Schwegerl B S. and Snell E E. Microbiological Methods for the Estimation of Amino Acids. Nutrition Abstracts & Rev. 16 (3) 437 510 (Jan.) 1946

Woodson H W, Herber S W, Solomon J D. and Bergman, O. Urinary Excretion of Amino Acids by Human Subjects on Normal Diets. J. Biol. Chem. 172 (2) 613 618 (Feb) 1948

**Arginine**

Horn M J, Jones D B. and Blum A E. Microbiological Determination of Arginine in Proteins and Foods. J. Biol. Chem. 178 59 63 (Oct.) 1949

**Choline**

Bandelin F J. The Microbiological Determination of Choline. J. Am. Pharm. A. (Scient. Ed.) 38 (6) 304 307 (June) 1949

Horowitz N H., and Beadle G W. A Microbiological Method for the Determination of Choline by Use of a Mutant of *Neurospora*. J. Biol. Chem. 150 (2) 3 5 333 (Oct) 1943

Luecke R. W. and Pearson P B. The Microbiological Determination of Free Choline in Plasma and Urine. J. Biol. Chem. 153 (1) 50 63 (April) 1944

**Glutamic Acid**

Dunn M S, Camien M N, Shankman S. and Block H. The Determination of d (-) and l (+) Glutamic Acid in Protein Hydrolysates. J. Biol. Chem. 163 (1) 43 49 (April) 1947

**Histamine**

Chen, O, Anson, C R. and Clarke I G. The Biological Assay of Histamine and Diphenhydramine Hydrochloride (Benadryl Hydrochloride). J. Pharmacol. and Exp. Therapeutics 97 90 9 (Jan.) 1948

**Isoleucine**

Horn M J., Jones, D B. and Blum, A. E. Microbiological Determination of Isoleucine in Proteins and Food. J. Biol. Chem. 180 69 91 (Sept.) 1949

## 260 Microbiological Procedures

### Lysine

Horn M. J. Jones D. B., and Blum A. E. Microbiological Determination of Lysine in Proteins and Foods, *J Biol. Chem* 169 (1) 71-76 (June) 1947

### Manganese

Bentley O. G., Snell E. E. and Phillips P. H. A Microbiological Method for the Determination of Manganese *J Biol Chem* 170 (1) 343-350 (Sept) 1947

### Pantothenic Acid

Bushirk H. H. Bergdahl, A. M., and Delor P. A. Enzymatic Digestion of Samples for Microbiological Assay of Pantothenic Acid *J Biol Chem* 172 (2) 671-675 (Feb) 1948

### Para-Amino Benzoic Acid

Landy M. and Dicken D. M. A Microbiological Method for the Determination of p-Amino Benzoic Acid, *J Biol. Chem* 146 109-114 (Nov) 1942

Lewis J. C. A Lactobacillus Assay Method for p-Amino Benzoic Acid, *J Biol. Chem* 146 441-450 (Dec) 1942

Thompson R. C. Isbell F. F. and Mitchell H. K. A Microbiological Assay Method for p-Amino Benzoic Acid *J Biol Chem* 148 281-287 (May) 1943

### Phenylalanine

Prescott B. A. Borek, E. Brecher A. and Waelisch H. Studies on Oxyphosphorylation of Phenylpyruvate. I. Microbiological Determination of L- and D-Phenylalanine and of Phenyllactic Acid *J Biol Chem* 181 273-280 (Nov) 1949

### Progesterone

Hooker, C. W. and Forbes T. R. A Bio Assay for Minute Amounts of Progesterone, *Endocrinol* 41 158-169 (Aug) 1947

### Pteric Acid

Franklin, A. L. Stokstad E. L. R. and Jukes T. H. Urinary Excretion Studies Following the Administration of Pteric Acid *Proc Soc Exper Biol & Med* 66 (3) 516-518 (Dec) 1947

### Riboflavin

Roberts E. C. and Snell E. E. An Improved Medium for Microbiological Assays With Lactobacillus Casei *J Bio Chem* 163 (2) 499-509 (May) 1946

### Sterols

Krámlí A. and Horvath, J. Microbiological Oxidation of Sterols, *Nature* 163 19 (Feb) 1949

### Tryptophane

Dunn M. S., Schoff H. F., Frankl W. and Rockland L. B. Investigation of Amino Acids, Peptides and Proteins. XX The Determination of Apparent Free Tryptophan in Blood by a Microbiological Method *J Biol Chem* 157 387-394 (Jan) 1945

### Vitamin B<sub>12</sub>

Capps B. F. Hobbs N. L., and Fox S. H. A Method for the Microbiological Assay of Vitamin B<sub>12</sub>, *J Biol Chem* 178 517-518 (March) 1949

Caswell M. C. Roditschek L. K. and Hendlin D. The Microbiological Estimation of Lactobacillus Lactis Dornier Activity With Vitamin B<sub>12</sub> as a Standard *J Biol. Chem* 180 125-132 (Aug) 1949

Hoffmann, C. E., Stokstad E. L. R. Hutchings B. L., Doranbush, A. C. and Jukes T. H. The Microbiological Assay of Vitamin B<sub>12</sub> With Lactobacillus Leichmannii *J Biol Chem* 181 635-644 (Dec) 1949

Rickes, E. L., Brink N. G., Koniusy F. R., Wood T. R. Folkers, K. Vitamin B<sub>12</sub> A Cobalt Complex *Science* 108 134 (Aug) 1949

Shorb M. S. Activity of Vitamin B<sub>12</sub> for the Growth of Lactobacillus Lactis *Science* 107 397-398 (April) 1948

## SECTION VI

# PHYSIOLOGICAL MEASUREMENTS

---

### A BLOOD GASES

#### 1 Calculation of Gas Volumes

For many physiological purposes it becomes necessary to compute gas volumes dry at 0°C and 760 mm Hg barometric pressure from observations made on gas saturated with water vapor at any temperature and barometric pressure

The general equation for this calculation is

$$\text{ml gas dry at 0°C and 760 mm Hg} = \text{ml observed} \times \frac{(273.2)}{T + 273.2} \times \frac{(B - V)}{(760)}$$

Where T is ambient temperature in °C B is ambient pressure in mm Hg and V is water vapor tension at the ambient temperature a value which must be obtained from empirical tables in standard handbooks and 273.2 is the absolute temperature for 0°C

The line chart (Fig 2a) was prepared by Dr Robert C Darling it gives at once the factor by which observed gas volume is multiplied to reduce it to volumes dry at 0°C and 760 mm Hg

#### 2 Techniques for Obtaining Arterial Blood and for Handling Blood Anaerobically

##### Usual Method for Obtaining Arterial Blood

Arterial blood is obtained by a procedure quite similar to that for venous blood. A # 0 gauge needle 1½ inches long with a short concave bevel is preferable. Particular attention is given to oiling the syringe. Care should be taken that the oil is free from air bubbles and air bubbles should not be introduced during the aspiration of blood into the syringe. Under no circumstances is a tourniquet to be applied. The brachial artery just proximal to the antecubital fossa is usually selected for puncture. It may be helpful to hyperextend the arm slightly during an arterial puncture by inserting a small pad under the elbow. The position of the artery is palpated with the left index finger. Arterial punctures are performed mainly by palpation rather than by observation of the vessel. Attention should be given to avoid penetrating a superficial vein before entering the artery. The skin is cleansed with alcohol. The syringe is held in the right hand at an angle of 45° and the skin and subcutaneous tissue entered. The needle and syringe are continued along a plane parallel to the artery. As soon as the artery is entered blood should flow into the syringe so that little or no decompression of the plunger is needed. Once the needle is in the artery its position should not be altered until the desired quantity of blood is removed. A large sponge is then applied forcefully over the assumed site of entrance into the artery and the syringe is removed. A firm bandage is immediately applied over the sponge with the pressure on and above the puncture never below



LINE CHART FOR DETERMINING FACTORS TO REDUCE SATURATED GAS VOLUMES TO  
 DRY VOLUMES AT 0°C AND 760 MM HG




Fig 25—Line chart for converting wet gas volumes to STP and dry (This chart was prepared by Dr Robert C Darling)

The blood is next transferred to a blood tube which contains an anticoagulant and  $\frac{1}{2}$  inch layer of mineral oil. During the initial stages of the transfer the syringe is held in a horizontal position and the plunger is not touched. The needle is removed and replaced by a rubber connection to a 6 inch capillary glass tube. The capillary tube is filled with blood by gentle compression of the plunger. The tip of the tube is next placed at the bottom of the blood tube below the oil and the remainder of the blood is gently delivered into the bottom of the tube. The blood may be mixed with the anticoagulant by repeated exchange of blood between syringe and tube meanwhile avoiding all air bubbles. An alternate procedure is gently stirring the blood in the tube with a glass stirring rod taking care to avoid emulsifying the blood with the oil.

**Technique of Goldschmidt and Light for Arterializing Venous Blood** (*J Biol. Chem.* 64 53 19 5)

The hand is immersed in hot water at 45 C for 20 minutes. Blood is drawn anaerobically from any convenient vein on the back of the hand.

**Technique of Lillenthal and Riley for Arterializing Venous Blood** (*J Clin. Investigation* 23 904 1944)

The ear is heated by radiant heat for 10 minutes at an air temperature of 45 C measured next to the pinna. The lobe is pricked and drops of blood are caught in oxalate or heparin in a specially constructed funnel from which samples may be drawn with a micropipette.

**Mook's Technique for Handling Blood Anaerobically** (*Biochem Z.* 242 3 9 1931)

One uses small glass tubes of length about 35 mm internal diameter about 3 mm and total content about .50 cmm. Each tube has a slight flare to the bottom.

The warmed finger is pricked and a tube is firmly applied to the hole. The diffusion of  $\text{CO}_2$  out of the blood and  $\text{O}_2$  into it under these conditions is such that when the tube is filled samples may be drawn with micropipettes from the center of the tube and the results are essentially the same as if the samples were drawn under oil.

**Scholander's Technique for Handling Blood Anaerobically** (*J Biol. Chem.* 148 541 1943 See Figs. 46 and 7)

The plunger of an all glass 3 or 5 ml syringe is lubricated with 1 or 2 drops of heavy oil. A little heparin solution is drawn into the syringe and the walls of the barrel moistened with the solution by drawing the plunger all the way down and pushing it up again the dead space at the top and the nozzle being left full of heparin. The solution is freed from air bubbles by rapid twisting of the syringe in the vertical position the syringe being held by the nozzle or needle the plunger being arrested by a clip. This procedure centrifuges the bubbles to the center so that they can be easily expelled. The nozzle is attached to a needle and the blood vessel punctured. The air bubble from the needle is cautiously expelled the needle removed and the nozzle dipped into a dish containing mercury a little of which is drawn into the syringe. The air bubble from the needle can be avoided by inserting the needle into the blood vessel with another syringe drawing the needle full of blood and then exchanging the first syringe for the syringe containing heparin. A  $\frac{1}{2}$  inch length of rubber tube is joined to the nozzle filled with blood and closed with a glass plug without trapping an air bubble. The mercury is used for mixing the blood in the syringe which is immersed in a beaker of ice and water. For high accuracy it may be necessary to correct for the slight dilution with the heparin solution.

For transfer to the pipette the syringe is removed from the beaker dried the blood well mixed and the fluid in the rubber tube pumped a few times into the main bulk of the blood by working the glass plug up and down. The plug and rubber tube are then removed and the pipette with the tip protruding about an inch over the edge of the beaker is held with its conical tip pressed against the opening of the syringe nozzle. If difficulty




Fig 26—Anaerobic storage of blood by Scholander's technique




Fig 27—Transfer of blood from storage syringe into pipette.

is experienced in this step the pipette tip may be covered with a rubber plug with a capillary bore. The pipette is then loaded by screwing in the plunger. The syringe is stoppered without trapping air and returned to the ice bath.

### 3 The Van Slyke-Neill Manometric Apparatus

#### Reference

Van Slyke D D and Neill J M The Determination of Gases in Blood and Other Solutions by Vacuum Extraction and Manometric Measurements *J Biol Chem* 61 5 3 538 (Sept) 19 4

#### Principle

In the Van Slyke apparatus advantage is taken of a vacuum. The gas is liberated chemically extracted by vacuum and its partial pressure measured in a partial vacuum on a mercury manometer calibrated in centimeters. All reagents employed must be free from physically dissolved gases in significant amounts.

#### Apparatus

(See Fig 28)

#### Preparation of Air Free Solutions

Dissolved gases are removed from the solution by adding a suitable volume of reagent to the chamber sealing the stopcock to the chamber with mercury and shaking the evacuated chamber. This is done by lowering the leveling bulb and opening the leveling stopcock. The liquid is allowed to fall to the point where the mercury is approximately at the 50 ml mark in the chamber. Shaking is allowed to proceed for three minutes then the machine is stopped and the mercury stopcock is opened allowing the solution to come to atmospheric pressure. The leveling bulb is now raised to the higher bracket on the machine and the stopcock to the chamber is opened allowing only the gas to escape. Repeat the procedure again.

About one ml of mercury is now added to the cup and a storing device usually some form of separatory funnel containing mineral oil or a burette seated in the chamber cup. The mercury acts as a gasket. The stopcock to the storage apparatus is opened and this is followed by opening the chamber stopcock forcing the fluid up into the storage tower where it is stored under oil.


Most reagents are evacuated and stored in this manner. Under the heading *Reagents* in each method the reagents that need evacuating and storage in this manner will be mentioned.

#### Machine Correction (C Corrections)

In all gasometric determinations a machine correction must be established (C). This is done by using the same procedure as for blood except that where blood is used a substituted volume of air free water is used or in the case of the determination of oxygen 1 ml extra oxygen reagent is substituted for blood. The C correction is really a correction that compensates for the pressure change due to the addition of the absorption reagents. For oxygen the C correction is about 10 mm for carbon dioxide about 15.0 mm of mercury.

#### Moisture in the Manometer

Small amounts of moisture find their way into the manometer as the mercury flows back and forth from the analysis chamber. Some error from variation of the vapor pressure might occur if it were left in the manometer tube. To absorb the water vapor a few drops of ethylene glycol (E. K. Co.) are admitted through the stopcock at the top of the manometer and permitted to flow down the tube for 10 cm. Allow the fluid to run up and down the manometer by raising and lowering the mercury leveling bulb. Finally the fluid is forced through the top stopcock and fluid is left behind but not enough to interfere with the readings of the mercury meniscus. Renew occasionally in the same manner.



**Precautions**

- 1 Always remember two things
  - a. Whenever taking in samples or solutions the leveling bulb must always be on the lower support
  - b. Whenever expelling any mixture from the chamber always place the leveling bulb in the upper support
- 2 When the chamber has been evacuated never open the main stopcock as the vacuum will send the mixture in the chamber back into the manometer

**1 Oxygen in Whole Blood (Van Slyke)****Reference**

- Van Slyke D D and Neill J M The Determination of Gases in Blood and Other Solutions by Vacuum Extraction and Manometric Measurements VII. The Determination of Oxygen in Blood *J Biol Chem* 61 55 554 (Sept) 19 4

**Principle**

Oxygen in blood is released from oxyhemoglobin by shaking in a vacuum with potassium ferric oxide. The liberated oxygen is measured manometrically by absorbing the oxygen with sodium hydrosulfite.

**Apparatus**

- 1 A Van Slyke manometric apparatus with an electric motor  
Van Slyke Neill stopcock pipettes calibrated between marks capacities of 0.5 ml 10 ml and 10 ml
- 2 Two 50 ml burettes for storing air free sodium hydroxide (NaOH) and sodium hydrosulfite ( $\text{Na}_2\text{S}_2\text{O}_4$ ) solutions. Keep two to three inches of oil in the burettes to prevent the solution from coming into contact with the outside air
- 3 Rubber tips for delivering sample into evacuating chamber

**Reagents**

- 1 Oxygen reagent 3 gm of pure white saponin 6 gm of potassium ferricyanide and 3 ml caprylic alcohol Dissolve and dilute to 1000 ml with water Shake well before using
- 2 Sodium hydroxide Approximately 1 N Dissolve 40 gm NaOH in water and dilute to 1000 ml. Make about 530 ml air free by shaking in a vacuum two to three times and expel gas between shakings Store in 50 ml burette under oil Make up air free NaOH each day
- 3 Sodium hydrosulfite ( $\text{Na}_2\text{S}_2\text{O}_4$  Lycopon) Mix 5 gm quickly in 0.5 ml of 1 N NaOH and make air free as soon as possible Store in the other 50 ml burette under paraffin oil. Do not try to dissolve any more than will go into solution by merely rotating flask for a few seconds Make fresh daily
- 4 Fieser's modified reagent (may be used for absorbing the oxygen) Dissolve in the following order in 100 ml of distilled water 14 gm of potassium hydroxide 16 gm of sodium hydrosulfite 3 gm sodium anthraquinone beta sulfonate and add 3 drops of 10% ferric chloride solution Store in rubber stoppered flask in the cold Make fresh every week
- 5 Mercury Although a cleaning tower may be used a simpler device is to use an airating bottle By bubbling air through the mercury in the latter case the mercury is brought into contact with the diluted  $\text{HNO}_3$  Bubble air through the Hg for several hours Remove  $\text{HNO}_3$  and repeat with several washes using distilled water Repeat once with glass-distilled water Drain off mercury and absorb visible water with filter paper
- 6 Stopcock grease (Lubri-Seal A. H. Thomas Co)

- 7 Lactic acid solution approximately 1 N Use C P lactic acid solution and dilute 92 ml to 1000 ml with distilled water
- 8 Caprylic alcohol (Eastman Kodak Co)
- 9 Oil (heavy paraffin)

#### Procedure

- 1 Rinse out apparatus by shaking 2 ml of 1 N lactic acid solution plus 10 ml of distilled water in a vacuum Expel into waste jar
- 2 The waste bore of the chamber stopcock is sealed with mercury Vacuum chamber and main bore of main cock are also filled with mercury Main cock is closed Leveling cock is open Mercury bulb is at atmospheric level (lower level)
- 3 Pour 6 ml oxygen reagent and 2 drops of caprylic alcohol into entry cup Open main cock and take in reagent to top of capillary between cup and chamber Close cock
- 4 Seal chamber and cock Pour a few drops of mercury into cup Open and close cock so that main bore fills with mercury from end to end. Suck off excess from cup
- 5 Evacuate solution Drop mercury level to 50 ml mark in chamber (Open leveling cock and lower bulb until sufficient mercury has left chamber Then close cock.) Shake gently for 3 minutes
- 6 Expel air from chamber Raise bulb to high level. Open leveling cock slowly (so that mercury will not strike top of manometer tube suddenly) Open main cock slowly and release air Let out reagent to 55 ml mark in cup Close main cock and support bulb at atmospheric level.
- 7 Admit blood Using main cock take in exactly 1 ml blood from a Van Slyke pipette and take oxygen reagent into cup to the 35 ml mark Close cock and seal Suck off excess fluid and mercury
- 8 Evacuate Drop mercury level to 50 ml mark as above Shake for 3 minutes
- 9 Absorb carbon dioxide Reduce gas phase to about 5 ml by opening and closing leveling cock Add 1 ml NaOH (air free) to cup Take in 0.5 ml drop by drop by manipulation of main cock Seal and suck off mercury and remains of the reagent Let reagent drain for a minute
- 10 Read P Adjust gas phase to 2 ml mark at a constant rate using leveling cock. P is the pressure of gaseous oxygen and nitrogen Read manometer to the nearest 1/10 mm.
- 11 Absorb oxygen Add 1 ml hydrosulfite solution to cup and take in 0.5 ml drop by drop over a period of 2 minutes When half has been added open leveling cock part way and allow contents to rise in chamber while remainder is added. Seal and suck off excess fluid cup There should be only a tiny bubble of nitrogen at top of vacuum chamber when oxygen is completely absorbed.
- 12 Read P Adjust gas phase to 2 ml mark using leveling cock and read manometer to nearest 1/10 mm
- 13 Record the temperature of the chamber to the nearest 1/10 C
- 14 Release pressure and expel the mixture into the waste jar
- 15 Rinse out with lactic acid and water again and expel

#### Calculation

- 1 mm partial pressure total O at T of water jacket =  $P_t - P - C$
- 2 ml total  $O_2$ /100 ml blood =  
(mm partial pressure O)  $\times$  (temperature factor from Table 16)
- 3 ml  $O_2$  bound by hemoglobin/100 ml blood =  
(ml total O - ml dissolved O at T and pO observed)
- a For venous blood dissolved O is 0.10 ml/100 ml
- b For arterial blood dissolved O is 0.25 ml/100 ml.

TABLE 16

## OXYGEN FACTORS (VAN SLYKE)

To obtain ml  $O_2$ /100 ml blood mult ply pO at jacket temperature by factor in table (One ml samples analyzed)

| JACKET<br>TEMP<br>C | $O_2$<br>FACTOR | JACKET<br>TEMP<br>C | $O$<br>FACTOR | JACKET<br>TEMP<br>C | $O$<br>FACTOR | JACKET<br>TEMP<br>C | $O$<br>FACTOR |
|---------------------|-----------------|---------------------|---------------|---------------------|---------------|---------------------|---------------|
| 10.0                | 493             | 0.0                 | 350           | 5.0                 | 408           | 30.0                | 360           |
| 1                   | 9               | 1                   | 49            | 1                   | 0.5           | 1                   | 65            |
| 2                   | 91              | 2                   | 48            | 2                   | 01            | 2                   | 64            |
| 3                   | 91              | 3                   | 47            | 3                   | 03            | 3                   | 64            |
| 4                   | 90              | 4                   | 46            | 4                   | 03            | 4                   | 63            |
| 5                   | 89              | 5                   | 45            | 5                   | 0             | 5                   | 6             |
| 6                   | 88              | 6                   | 45            | 6                   | 01            | 6                   | 61            |
| 7                   | 87              | 7                   | 44            | 7                   | 100           | 7                   | 60            |
| 8                   | 87              | 8                   | 43            | 8                   | 100           | 8                   | 60            |
| 9                   | 86              | 9                   | 4             | 9                   | 399           | 9                   | 9             |
| 10.0                | 485             | 1.0                 | 441           | 6.0                 | 398           | 31.0                | 358           |
| 1                   | 84              | 1                   | 40            | 1                   | 97            | 1                   | 57            |
| 2                   | 84              | 2                   | 39            | 2                   | 96            | 2                   | 56            |
| 3                   | 83              | 3                   | 38            | 3                   | 96            | 3                   | 56            |
| 4                   | 8               | 4                   | 37            | 4                   | 9             | 4                   | 55            |
| 5                   | 81              | 5                   | 36            | 5                   | 94            | 5                   | 54            |
| 6                   | 81              | 6                   | 36            | 6                   | 93            | 6                   | 53            |
| 7                   | 80              | 7                   | 3             | 7                   | 9             | 7                   |               |
| 8                   | 79              | 8                   | 34            | 8                   | 9             | 8                   |               |
| 9                   | 78              | 9                   | 33            | 9                   | 91            | 9                   | 51            |
| 11.0                | 48              | 0                   | 4             | 0                   | 399           | 4.0                 | 350           |
| 1                   | 77              | 1                   | 31            | 1                   | 89            | 1                   | 49            |
| 2                   | 76              | 2                   | 0             | 2                   | 88            | 2                   | 48            |
| 3                   | 7               | 3                   | 9             | 3                   | 88            | 3                   | 48            |
| 4                   | 4               | 4                   | 98            | 4                   | 8             | 4                   | 47            |
| 5                   | 73              | 5                   | 7             | 5                   | 86            | 5                   | 46            |
| 6                   | 7               | 6                   | 7             | 6                   | 85            | 6                   | 4             |
| 7                   | 71              | 7                   | 96            | 7                   | 84            | 7                   | 44            |
| 8                   | 70              | 8                   | 5             | 8                   | 81            | 8                   | 44            |
| 9                   | 69              | 9                   | 94            | 9                   | 83            | 9                   | 43            |
| 12.0                | 468             | 1.0                 | 411           | 8.0                 | 383           | 33.0                | 351           |
| 1                   | 67              | 1                   | 31            | 1                   | 81            | 1                   | 41            |
| 2                   | 66              | 2                   | 31            | 2                   | 80            | 2                   | 40            |
| 3                   | 66              | 3                   | 0             | 3                   | 9             | 3                   | 39            |
| 4                   | 65              | 4                   | 19            | 4                   | 9             | 4                   | 38            |
| 5                   | 64              | 5                   | 18            | 5                   | 8             | 5                   | 3             |
| 6                   | 63              | 6                   | 18            | 6                   | 7             | 6                   | 3             |
| 7                   | 6               | 7                   | 1             | 7                   | 6             | 7                   | 36            |
| 8                   | 61              | 8                   | 16            | 8                   | 6             | 8                   | 35            |
| 9                   | 60              | 9                   | 15            | 9                   | 5             | 9                   | 34            |
| 13.0                | 459             | 1.0                 | 414           | 9.0                 | 384           | 34.0                | 351           |
| 1                   | 58              | 1                   | 33            | 1                   | 73            |                     |               |
| 2                   | 57              | 2                   | 32            | 2                   |               |                     |               |
| 3                   | 56              | 3                   | 32            | 3                   |               |                     |               |
| 4                   | 55              | 4                   | 11            | 4                   | 1             |                     |               |
| 5                   | 54              | 5                   | 10            | 5                   | 0             |                     |               |
| 6                   | 51              | 6                   | 99            | 6                   | 69            |                     |               |
| 7                   | 51              | 7                   | 98            | 7                   | 68            |                     |               |
| 8                   | 5               | 8                   | 98            | 8                   | 68            |                     |               |
| 9                   | 51              | 9                   | 0             | 9                   | 6             |                     |               |



c For blood equilibrated with air at 760 mm and 25 C dissolved O is 0.6 ml/100 ml

$$4 \text{ gm hemoglobin/100 ml blood} = \frac{(\text{ml bound oxygen})}{1.34}$$

#### Example

P = 153.3 mm P = 72.5 mm C = 0.7 mm T = 25 C

1 ml total O<sub>2</sub>/100 ml = (153.2 - 72.5 - 0.7) × (0.2406) = 19.25

2 ml dissolved O /100 ml = 0.62

3 ml bound O = 18.63/100 ml blood

4 gm Hb/100 ml blood = 18.63/1.34 = 13.9

#### Precautions

See section on apparatus

### 5 Carbon Dioxide in Whole Blood and Serum (Van Slyke)

#### Reference

Van Slyke D D and Neill J M The Determination of Gases in Blood by Vacuum Extraction and Manometric Measurements V Determination of Carbon Dioxide  
J Biol Chem 61 543-553 (Sept) 1924

#### Principle

Both bound carbonate and free carbon dioxide are liberated from the blood or serum with lactic acid and the liberated carbon dioxide is measured manometrically by absorption with sodium hydroxide

#### Apparatus

The same as for Oxygen.

#### Reagents

- 1 Caprylic alcohol (Eastman Kodak Co)
- 2 Distilled water Make air free and store under oil as for oxygen reagents.
- 3 Sodium hydroxide approximately 1 N 40 gm dissolved in 1000 ml distilled water Make air free
- 4 Lactic acid approximately 1 N 92 ml USP lactic acid diluted to 1000 ml with water
- 5 Very clean mercury (See oxygen method)
- 6 Paraffin oil heavy

#### Procedure

- 1 Clean out apparatus by shaking in a vacuum with 2 ml lactic acid solution and approximately 10 ml of distilled water Shake approximately one minute and expel into waste jar
- 2 Seal off the main stopcock at both ends with mercury by raising the leveling bulb to the upper position
- 3 Place leveling bulb at lower position and admit one to 2 drops of caprylic alcohol into the capillary of the main stopcock
- 4 Place 2.5 ml of air free water into the entry cup
- 5 Pipette exactly 1 ml of serum or blood (when using blood make certain that the blood specimen has been stirred well as too much of either the red cells or plasma will cause very erroneous results) and introduce into the chamber through a rubber tip
- 6 Then wash in the sample with the remainder of the water
- 7 Add 6 drops (0.3 ml) of lactic acid solution and admit into the chamber
- 8 Seal the main stopcock with mercury and draw off the excess mercury in the cup by means of a vacuum

TABLE 17

## CO FACTORS FOR VAN SLYKE APPARATUS

(For 1 ml samples of blood or serum. To obtain ml total CO per 100 ml sample multiply pCO by factor for jacket temperature.)

| JACKET<br>TEMP<br>C | CO<br>FACTOR | JACKET<br>TEMP<br>C | CO<br>FACTOR | JACKET<br>TEMP<br>C | CO <sub>2</sub><br>FACTOR | JACKET<br>TEMP<br>C | CO<br>FACTOR |
|---------------------|--------------|---------------------|--------------|---------------------|---------------------------|---------------------|--------------|
| 15.0                | .35          | 0.0                 | .66          | 5.0                 | .94                       | 30.0                | .53          |
| 1                   | 3            | 1                   | .61          | 1                   | .93                       | 1                   | .5           |
|                     | 30           |                     | .60          | "                   | .91                       | 2                   | .31          |
| 3                   | 8            | 3                   | .59          | 3                   | .90                       | 3                   | .30          |
| 4                   | 6            | 4                   | .57          | 4                   | .89                       | 4                   | .29          |
| "                   | 4            | 5                   | .5           | 5                   | .88                       | 5                   |              |
| 6                   | .23          | 6                   | .54          | 6                   | .86                       | 6                   | .6           |
| 7                   | .2           | 7                   | .53          | 7                   | .8                        |                     | .5           |
| 8                   | .21          | 8                   | .52          | 8                   | .84                       | 8                   | .4           |
| 9                   | .20          | 9                   | .50          | 9                   | .83                       | 9                   | .3           |
| 16.0                | .719         | 21.0                | .648         | 6.0                 | .581                      | 31.0                | .5           |
| 1                   | .17          | 1                   | .47          | 1                   | .80                       | 1                   | .1           |
|                     | .16          |                     | .46          | "                   | .79                       | "                   | .0           |
| 3                   | .15          | 3                   | .45          | 3                   | .77                       | 3                   | .19          |
| 4                   | .13          | 4                   | .44          | 4                   | .76                       | 4                   | .18          |
| 5                   | .11          | 5                   | .43          | 5                   | .5                        | .5                  | .16          |
| 6                   | .10          | 6                   | .41          | 6                   | .74                       | 6                   | .15          |
| 7                   | .8           | 7                   | .39          | 7                   | .73                       | 7                   | .14          |
| 8                   | .7           | 8                   | .37          | 8                   | .71                       | 8                   | .13          |
| 9                   | .5           | 9                   | .36          | 9                   | .70                       | 9                   | .12          |
| 17.0                | .04          | 0                   | .634         | 7.0                 | .569                      | 32.0                | .511         |
| 1                   | .03          | 1                   | .33          | 1                   | .68                       | 1                   | .10          |
|                     | 0            |                     | .3           | "                   | .67                       | "                   | .09          |
| 3                   | .00          | 3                   | .31          | 3                   | .65                       | 3                   | .08          |
| 4                   | .699         | 4                   | .30          | 4                   | .64                       | 4                   | .07          |
|                     | .69          | 5                   | .9           | 5                   | .63                       | 5                   | .05          |
| 6                   | .96          | 6                   | .6           | 6                   | .6                        | 6                   | .04          |
|                     | .94          | 7                   | .4           |                     | .61                       | 7                   | .03          |
| 8                   | .93          | 8                   | "            | 8                   | .59                       | 8                   | .02          |
| 9                   | .91          | 9                   | .1           | 9                   | .58                       | 9                   | .01          |
| 18.0                | .090         | 3.0                 | .60          | 8.0                 | .55                       | 33.0                | .2500        |
| 1                   | .89          | 1                   | .19          | 1                   | .56                       | 1                   | .499         |
|                     | .87          | "                   | .18          | .2                  | .55                       | 2                   | .98          |
| 3                   | .86          | .3                  | .16          | 3                   | .53                       | .3                  | .97          |
| 4                   | .84          | 4                   | .1           | 4                   | .5                        | 4                   | .96          |
|                     | .82          | 5                   | .14          | 5                   | .51                       | 5                   | .94          |
| 6                   | .81          | 6                   | .1           | 6                   | .50                       | 6                   | .93          |
|                     | .9           | 7                   | .11          | 7                   | .49                       | 7                   | .9           |
| 8                   | .78          | 8                   | .10          | 8                   | .4                        | 8                   | .91          |
| 9                   | .76          | 9                   | .09          | 9                   | .46                       | 9                   | .90          |
| 19.0                | .675         | 4.0                 | .60          | 9.0                 | .54                       | 34.0                | .159         |
| 1                   | .74          | 1                   | .06          | 1                   | .44                       |                     |              |
|                     | .73          | "                   | .04          | .2                  | .43                       |                     |              |
| 3                   | .7           | 3                   | .03          | 3                   | .41                       |                     |              |
| 4                   | .71          | 4                   | .0           | 4                   | .40                       |                     |              |
| 5                   | .69          | .5                  | .01          | .5                  | .39                       |                     |              |
| 6                   | .68          | 6                   | .99          | 6                   | .38                       |                     |              |
| 7                   | .67          | 7                   | .98          |                     | .37                       |                     |              |
| 8                   | .66          | 8                   | .97          | .8                  | .35                       |                     |              |
| 9                   | .64          | 9                   | .95          | 9                   | .34                       |                     |              |

- 9 Take leveling bulb from lower support and create a vacuum in the mixing chamber by lowering the bulb below the lower support. When the blood mixture reaches the 50 ml mark close the vacuum controlling stopcock. Place leveling bulb on lower support and shake for three minutes.
- 10 Raise the contents of the chamber by opening the vacuum controlling stopcock at a fairly rapid rate to prevent the reabsorption of the carbon dioxide. Stop when the gas phase reaches the 2 ml mark and read immediately (P). If you overshoot the 2 ml mark draw back to the 50 ml mark and shake again for three minutes. Read again on the 2 ml mark.
- 11 Drop gas phase to approximately 5 ml region of the chamber by creating a vacuum again.
- 12 Add approx 15 ml of air free NaOH to the cup. By opening the main stopcock admit drop by drop one ml of NaOH.
- 13 Seal the main stopcock again with mercury and draw off the excess NaOH in the cup by means of a vacuum.
- 14 Allow to stand one minute for complete absorption of the carbon dioxide and adjust the gas phase to the 2 ml mark. This is done by opening the vacuum controlling stopcock slowly. If the mark is overshoot just lower the gas phase below the 2 ml mark allow to drain again and read (P). Release pressure and place leveling bulb on upper support. Expel the blood mixture into the waste jar. Place the leveling bulb at the lower support and clean the apparatus again with lactic acid and water.
- 15 Record temperature of water jacket.

**Calculation**

ml total  $\text{CO}_2$ /100 ml blood =

$(P_1 - P - C) \times \text{temperature factor for CO from Table 17}$

**Example**

P 271.3 mm P 100 mm C 1.3 mm

Water jacket temperature 25 C

ml total  $\text{CO}_2$ /100 ml blood =

$(271.3 - 70.0 - 1.3) \times 0.0594 = 51.9 \text{ ml}/100 \text{ ml}$

**Precautions**

See section on apparatus

## 6 Oxygen and Carbon Dioxide in Single Samples of Blood (Van Slyke)

**Reference**

Van Slyke D D and Neill J M The Determination of Gases in Blood and Other Solutions by Vacuum Extraction and Manometric Measurements VIII. Combined Determination of Oxygen and Carbon Dioxide in Blood *J Biol. Chem.* 61: 561-562 (Sept.) 1924

**Principle**

By combining the technique for oxygen alone and carbon dioxide alone a single sample suffices for both substances. In our experience the most consistent results are obtained when the size of samples is 0.5 ml.

**Apparatus**

The Van Slyke manometric apparatus

**Reagents**

1 Caprylic alcohol (Eastman Kodak Co.)

2 Combined reagent Mix equal volumes of the following

a Lactic acid 12% in water

- b. Ferricyanide solution containing in 1000 ml water 33 gm potassium ferricyanide  $[K_3Fe(CN)_6]$  8 gm pure sodium and 3 ml caprylic alcohol. Shake before mixing with lactic acid solution
- 3 Sodium hydroxide 1 N Dissolve 40 gm in 1000 ml of water
- 4 Sodium hydrosulfite See the Oxygen method.

#### Procedure

- 1 Rinse out apparatus with lactic acid and water before proceeding with the analysis. Admit drops of caprylic alcohol into capillary of main (entry) cock
- 2 Admit 6 ml combined reagent to vacuum chamber. Seal with mercury. Evacuate. Shake for 3 minutes
- 3 Expel air and let out reagent to 5 ml mark in cup
- 4 Take in exactly 0.5 ml blood and lower reagent in cup to 3.0 ml mark. If using 1.0 ml blood lower reagent to 3.5 ml mark
- 5 Seal chamber and suck off excess fluid from cup. Evacuate to 50 ml mark and shake for three minutes
- 7 Adjust gas phase to 0 ml mark in chamber
- 8 Read P value on mercury manometer. If mark is overshot evacuate and shake for three minutes more
- 9 Lower contents to 5 ml region. Add 1.0 ml NaOH to cup and admit 0.5 ml drop by drop into chamber
- 10 Seal with mercury. Wait a few seconds while CO is absorbed.
- 11 Adjust gas phase to 2 ml mark and read P value. Recheck P by lowering to approx. 5 ml and adjusting to the 2 ml mark again. The reading should be constant before continuing with the next step
- 12 From here on proceed as in oxygen method. The final reading P = pressure of gaseous N
- 13 Record the temperature in the water jacket

#### Calculation (0.5 ml samples)

- 1 ml  $CO_2/100$  ml blood =  $(P - P - C_{CO_2}) \times (\text{temp factor for } CO_2 \text{ Table 17}) \times "$
- ml total  $O_2/100$  ml blood =  $(P - P - C_O) \times (\text{temp factor for } O_2 \text{ Table 16}) \times 2$

#### Example

Size of blood sample 0.5 ml  
 P 193.0 mm P 108 mm.  $P_a$  0.0 mm  
 Water Jacket 5 C C for  $CO_2$  1.3 C for O 0.7  
 ml  $CO_2/100$  ml blood =  $(193.0 - 108 - 1.3) \times 0.534 \times 2 = 45.9$   
 ml  $O_2/100$  ml blood =  $(108.2 - 0.0 - 0.7) \times 0.2406 \times 2 = 18.0$

#### Precautions

See section on apparatus. We have obtained more reproducible results using 0.5 ml samples than when using 1 ml samples

#### Oxygen Capacity (Van Slyke)

##### Reference

Boughton F J W Darling H C and Root W S Factors Affecting the Determination of Oxygen Capacity Content and Pressure in Human Arterial Blood  
 Am J Physiol 142 33-50 (Dec) 1944

- 9 Take leveling bulb from lower support and create a vacuum in the mixing chamber by lowering the bulb below the lower support. When the blood mixture reaches the 50 ml mark close the vacuum controlling stopcock. Place leveling bulb on lower support and shake for three minutes.
- 10 Raise the contents of the chamber by opening the vacuum controlling stopcock at a fairly rapid rate to prevent the reabsorption of the carbon dioxide. Stop when the gas phase reaches the 2 ml mark and read immediately (P). If you overshoot the 2 ml mark draw back to the 50 ml mark and shake again for three minutes. Read again on the 2 ml mark.
- 11 Drop gas phase to approximately 5 ml region of the chamber by creating a vacuum again.
- 12 Add approx 15 ml of air free NaOH to the cup. By opening the main stopcock admit drop by drop one ml of NaOH.
- 13 Seal the main stopcock again with mercury and draw off the excess NaOH in the cup by means of a vacuum.
- 14 Allow to stand one minute for complete absorption of the carbon dioxide and adjust the gas phase to the 2 ml mark. This is done by opening the vacuum controlling stopcock slowly. If the mark is overshoot just lower the gas phase below the 2 ml mark allow to drain again and read (P). Release pressure and place leveling bulb on upper support. Expel the blood mixture into the waste jar. Place the leveling bulb at the lower support and clean the apparatus again with lactic acid and water.
- 15 Record temperature of water jacket.

**Calculation**

ml total  $\text{CO}_2$ /100 ml blood =

$(P_1 - P - C) \times \text{temperature factor for CO from Table 17}$

**Example**

P 271.3 mm P 70.0 mm C 1.3 mm

Water jacket temperature 25°C

ml total  $\text{CO}_2$ /100 ml blood =

$(271.3 - 70.0 - 1.3) \times 0.2594 = 51.9 \text{ ml}/100 \text{ ml}$

**Precautions**

See section on apparatus

## 6 Oxygen and Carbon Dioxide in Single Samples of Blood (Van Slyke)

**Reference**

Van Slyke D D and Neill J M The Determination of Gases in Blood and Other Solutions by Vacuum Extraction and Manometric Measurements VIII. Combined Determination of Oxygen and Carbon Dioxide in Blood J Biol Chem 61 561 562 (Sept) 1924

**Principle**

By combining the technique for oxygen alone and carbon dioxide alone a single sample suffices for both substances. In our experience the most consistent results are obtained when the size of samples is 0.5 ml.

**Apparatus**

The Van Slyke manometric apparatus

**Reagents**

1 Caprylic alcohol (Eastman Kodak Co)

2 Combined reagent Mix equal volumes of the following

a. Lactic acid 1.2% in water

- b Ferricyanide solution containing in 1000 ml water 3 gm potassium ferricyanide  $[K_3Fe(CN)_6]$  6 gm pure saponin and 3 ml caprylic alcohol. Shake before mixing with lactic acid solution.
- 3 Sodium hydroxide 1 N Dissolve 40 gm in 1000 ml of water
- 4 Sodium hydrosulfite See the Oxygen method.

#### Procedure

- 1 Rinse out apparatus with lactic acid and water before proceeding with the analysis. Admit 2 drops of caprylic alcohol into capillary of main (entry) cock.
- 2 Admit 6 ml combined reagent to vacuum chamber. Seal with mercury. Evacuate. Shake for 3 minutes.
- 3 Expel air and let out reagent to 5 ml mark in cup.
- 4 Take in exactly 0.5 ml blood and lower reagent in cup to 3.0 ml mark.
- 5 If using 1.0 ml blood lower reagent to 3.5 ml mark.
- 6 Seal chamber and suck off excess fluid from cup. Evacuate to 50 ml mark and shake for three minutes.
- 7 Adjust gas phase to 1 ml mark in chamber.
- 8 Read P value on mercury manometer. If mark is overshot evacuate and shake for three minutes more.
- 9 Lower contents to 5 ml reg on. Add 1.0 ml NaOH to cup and admit 0.5 ml drop by drop into chamber.
- 10 Seal with mercury. Wait a few seconds while CO is absorbed.
- 11 Adjust gas phase to 2 ml mark and read P value. Recbeck P by lowering to approx 5 ml and adjust ng to the 1 ml mark again. The reading should be constant before continuing with the next step.
- 12 From here on proceed as in oxygen method. The final reading  $P_2$  = pressure of gaseous  $N_2$ .
- 13 Record the temperature in the water jacket.

#### Calculation (0.5 ml samples)

- 1 ml  $CO_2/100$  ml blood =  
 $(P - P - C_{CO}) \times (\text{temp factor for CO Table 17}) \times 2$
- 2 ml total  $O_2/100$  ml blood =  
 $(P - P - C_O) \times (\text{temp factor for O Table 16}) \times 2$

#### Example

Size of blood sample 0.5 ml.  
 $P_0$  193.0 mm  $P$  108.2 mm.  $P_2$  90 mm  
 Water Jacket  $^{\circ}C$   $C$  for CO 13  $C$  for O 0  
 ml  $CO_2/100$  ml blood =  
 $(193.0 - 108.2 - 13) \times 0.94 \times 2 = 45.9$   
 ml  $O_2/100$  ml blood =  
 $(108.2 - 90 - 0.7) \times 0.406 \times 2 = 18.05$

#### Precautions

See section on apparatus. We have obtained more reproducible results using 0.5 ml samples than when using 1 ml samples.

## 7 Oxygen Capacity (Van Slyke)

#### Reference

Poughton F J W Darling L C and Foot, W S Factors Affecting the Determination of Oxygen Capacity Content and Pressure in Human Arterial Blood  
 Am J Physiol 14 93 0 (Dec) 1944

- 9 Take leveling bulb from lower support and create a vacuum in the mixing chamber by lowering the bulb below the lower support. When the blood mixture reaches the 50 ml mark close the vacuum controlling stopcock. Place leveling bulb on lower support and shake for three minutes.
- 10 Raise the contents of the chamber by opening the vacuum controlling stopcock at a fairly rapid rate to prevent the reabsorption of the carbon dioxide. Stop when the gas phase reaches the 2 ml mark and read immediately (P). If you overshoot the 2 ml mark draw back to the 50 ml mark and shake again for three minutes. Read again on the 2 ml mark.
- 11 Drop gas phase to approximately 5 ml region of the chamber by creating a vacuum again.
- 12 Add approx 15 ml of air free NaOH to the cup. By opening the main stopcock admit drop by drop one ml of NaOH.
- 13 Seal the main stopcock again with mercury and draw off the excess NaOH in the cup by means of a vacuum.
- 14 Allow to stand one minute for complete absorption of the carbon dioxide and adjust the gas phase to the 2 ml mark. This is done by opening the vacuum controlling stopcock slowly. If the mark is overshoot just lower the gas phase below the 2 ml mark allow to drain again and read (P). Release pressure and place leveling bulb on upper support. Expel the blood mixture into the waste jar. Place the leveling bulb at the lower support and clean the apparatus again with lactic acid and water.
- 15 Record temperature of water jacket.

**Calculation**

ml total  $\text{CO}_2$ /100 ml blood =

$(P - P - C) \times \text{temperature factor for CO from Table 17}$

**Example**

P 2713 mm P 700 mm C 13 mm

Water jacket temperature  $^{\circ}\text{C}$

ml total  $\text{CO}_2$ /100 ml blood =

$(2713 - 700 - 13) \times 0.2594 = 51.9 \text{ ml}/100 \text{ ml}$

**Precautions**

See section on apparatus

## 6 Oxygen and Carbon Dioxide in Single Samples of Blood (Van Slyke)

**Reference**

Van Slyke D. D. and Neill, J. M. The Determination of Gases in Blood and Other Solutions by Vacuum Extraction and Manometric Measurements. VIII. Combined Determination of Oxygen and Carbon Dioxide in Blood. *J. Biol. Chem.* 61: 561-567 (Sept.) 1944.

**Principle**

By combining the technique for oxygen alone and carbon dioxide alone a single sample suffices for both substances. In our experience the most consistent results are obtained when the size of samples is 0.5 ml.

**Apparatus**

The Van Slyke manometric apparatus

**Reagents**

- 1 Caprylic alcohol (Eastman Kodak Co.)
- 2 Combined reagent. Mix equal volumes of the following
  - a. Lactic acid 12% in water

- b. Ferricyanide solution containing in 1000 ml water 3 gm potassium ferricyanide [ $K_3Fe(CN)_6$ ] 8 gm pure saponin and 3 ml caprylic alcohol. Shake before mixing with lactic acid solution.
3. Sodium hydroxide 1 N. Dissolve 40 gm in 1000 ml of water.
4. Sodium hydrosulfite. See the Oxygen method.

#### Procedure

1. Rinse out apparatus with lactic acid and water before proceeding with the analysis. Admit 2 drops of caprylic alcohol into capillary of main (entry) cock.
2. Admit 6 ml combined reagent to vacuum chamber. Seal with mercury. Evacuate. Shake for 3 minutes.
3. Expel air and let out reagent to 5 ml mark in cup.
4. Take in exactly 0.5 ml blood and lower reagent in cup to 3.0 ml mark.
- ✓ If using 1.0 ml blood lower reagent to 3.5 ml mark.
6. Seal chamber and suck off excess fluid from cup. Evacuate to 50 ml mark and shake for three minutes.
7. Adjust gas phase to 2 ml mark in chamber.
8. Read P value on mercury manometer. If mark is overshoot, evacuate and shake for three minutes more.
9. Lower contents to 5 ml region. Add 1.0 ml NaOH to cup and admit 0.5 ml drop by drop into chamber.
10. Seal with mercury. Wait a few seconds while CO is absorbed.
11. Adjust gas phase to 2 ml mark and read P value. Recheck P by lowering to approx. 5 ml and adjusting to the 2 ml mark again. The reading should be constant before continuing with the next step.
12. From here on proceed as in oxygen method. The final reading P = pressure of gaseous  $N_2$ .
13. Record the temperature in the water jacket.

#### Calculation (0.5 ml samples)

- 1 ml  $CO_2/100$  ml blood =  
 $(P - P - C_{CO}) \times (\text{temp factor for } CO_2, \text{ Table 17}) \times 2$
- ml total  $O_2/100$  ml blood =  
 $(P - P - C_O) \times (\text{temp factor for } O_2, \text{ Table 16}) \times 2$

#### Example

Size of blood sample 0.5 ml.  
 $P_b$  198.0 mm I. 108.2 mm  $P_w$  0.0 mm  
 Water Jacket 5 C. C for  $CO_2$  13 C for O 0  
 ml  $CO_2/100$  ml blood =  
 $(198.0 - 108.2 - 13) \times 0.594 \times 2 = 43.9$   
 ml  $O_2/100$  ml blood =  
 $(108.2 - 0.0 - 0.7) \times 0.406 \times 2 = 88.0$

#### Precautions

See section on apparatus. We have obtained more reproducible results using 0.5 ml samples than when using 1 ml samples.

## 7 Oxygen Capacity (Van Slyke)

#### Reference

- Foughton F. J. W., Darling R. C. and Root W. S. Factors Affecting the Determination of Oxygen Capacity Content and Pressure in Human Arterial Blood. *Am J Physiol* 142: 35-50 (Dec.) 1944.



**Principle**

The reference cited discusses significant errors in the common technique for determining oxygen capacity, i.e. equilibrating blood with air in a tonometer. The authors recommend equilibration directly in the Van Slyke extraction chamber.

**Apparatus**

- 1 Van Slyke manometric apparatus
- 2 Roughton Scholander syringe pipette

**Reagents**

- 1 Sodium chloride solution, 0.9%, equilibrated in a tonometer with 100 per cent oxygen
- 2 A solution of 3 % potassium ferricyanide and 8% of saponin deaerated in a tonometer and transferred to a burette (no oil) with no foam
- 3 Sodium hydroxide 0.1 N Deaerate and store under oil
- 4 Pyrogalllic acid 15% in sodium hydroxide solution (20%)
- 5 Caprylic alcohol
- 6 Clean mercury

**Procedure**

- 1 Clean chamber as in previous analyses
- 2 Admit two drops of caprylic alcohol and 1 ml of oxygenated saline into the chamber
- 3 Add exactly 1 ml blood into the cup and wash into the chamber with 5 ml of oxygenated saline
- 4 Mix contents of chamber and let in room air until mercury level is at the 50 ml mark.
- 5 Shake for two minutes and expel all the air quantitatively. Note that this is the time at which the blood is saturated with oxygen.
- 6 Add 0.25 ml of ferricyanide saponin solution and seal the chamber with mercury
- 7 Lower the mercury level to the 50 ml mark and shake the chamber for three minutes
- 8 Release the pressure of the gas phase to approximately 5 to 10 ml
- 9 Add 15 ml of 0.1 N NaOH in the cup and admit 1 ml into the chamber slowly. Seal with mercury and read at the 2 ml mark (P). When reading is constant proceed with the analyses
- 10 Add 15 ml of pyrogalllic acid solution to the cup and admit 1 ml very slowly keeping liquid level near the 2 ml mark. Allow two to three minutes for the absorption
- 11 Before sealing mix contents in the chamber open lower stopcock allowing gas to contract to a small bubble seal slowly so that last oxygen is absorbed by the last few drops of pyrogallol
- 12 Lower just below the two ml mark allow to drain and read (P) at the 0 ml mark

**Calculation**

See Oxygen in blood (Van Slyke)

**Example**

See Oxygen in blood (Van Slyke)

**Precautions**

- 1 The oxygen capacity of arterial blood as obtained by 15 minutes saturation with air in a tonometer followed by analysis in the Van Slyke Neill manometric apparatus is considered to be about 2 per cent higher than the value in the blood at the moment of withdrawal from the blood vessel. Three factors contribute to this discrepancy

- a Drainage errors in the saturator These cause the blood sample removed for analysis to be unduly rich in red cells (average effect equals about 0.7 per cent)
- b Gradual reversion of some of the inactive pigment to the active gas combining form. This effect is very variable but on the average equals 1.0 to 1.5 per cent of the total capacity Owing to the smallness of the effect the nature of the inactive reverting pigment has not been established. It may be in part methemoglobin
- c Presence of traces of COHb in the blood.
- \* Correction for these factors raises the average percent O Hb of normal man at sea level from 95 to 97 the arterial pO<sub>2</sub> as calculated from these figures and the blood dissociation curves comes out to 100 mm in agreement with alveolar pO<sub>2</sub> rather than to the average figure of 80 mm given in recent papers The latter figure is believed to be false

## 8 Carbon Monoxide in Blood (Total Heme Pigments)

### Reference

Horvath S M and Roughton F J W Improvements in the Gasometric Estimation of Carbon Monoxide in Blood *J Biol Chem* 144 747 755 (June) 194

### Principle

Dissolved carbon monoxide dissolved nitrogen and dissolved carbon dioxide are liberated into the gas phase of the Van Slyke manometric apparatus Bound carbon monoxide and carbon dioxide are released at pH 6 by shaking with a hydro sulfite ferri cyanide solution Carbon dioxide is absorbed with sodium hydroxide and the residual gas is solely carbon monoxide

### Apparatus

- 1 Van Slyke Neill manometric apparatus
- \* Roughton Scholander syringe pipette

### Reagents

- 1 Caprylic alcohol (Eastman Kodak Co)
- 2 Saponin a 1% solution of pure white saponin
- 3 Sodium hydrosulfite ( $\text{Na}_2\text{S}_2\text{O}_4$ ) a 1% solution in saturated sodium borate A stock of 3% sodium borate ( $\text{Na}_2\text{B}_4\text{O}_7 \cdot 2\text{H}_2\text{O}$ ) is prepared. At the beginning of the day 50 ml are placed in a 50 ml Erlenmeyer flask 1 gm of solid  $\text{Na}_2\text{S}_2\text{O}_4$  is dropped in the flask is quickly corked and shaken with only a minute bubble of air therein. The solution is at once transferred to a 50 ml burette
- 4. Potassium ferri cyanide [ $\text{K}_3\text{Fe}(\text{CN})_6$ ] Prepare a 30% air free solution by shaking in an evacuated tonometer and transferring to a burette with rubber outlet and pinchcock. If it darkens a fresh supply must be made.
- 5 Acetate buffer 70 gm of sodium acetate plus 100 ml water and 15 ml of glacial acetic acid.
- 6 Mixed acetate buffer and  $\text{K}_3\text{Fe}(\text{CN})_6$  reagent (air free) 50 ml of  $\text{K}_3\text{Fe}(\text{CN})_6$  solution plus 1.5 ml of the acetate buffer
- 7 Sodium hydroxide ( $\text{NaOH}$ ) a 10% solution (air free)
- 8 Water saturated with carbon monoxide at one atmosphere. Keep under paraffin oil in a 50 ml burette
- 9 Clean mercury (See oxygen method)

### Procedure

- 1. Admit 5 ml of water into a clean chamber seal with mercury evacuate and shake for 2 minutes.

## 276 *Physiological Measurements*

- 2 Expel all the water into the waste jar and seal off both ends of the main stop cock with mercury
- 3 5 drops of caprylic alcohol are drawn into the extracting chamber
- 4 Place 15 ml of 1% saponin solution into the cup
- 5 Exactly 1 or 2 ml of blood (according to whether the per cent COHb is greater or less than 50) are pipetted directly into the chamber followed by the remainder of the saponin solution that is in the cup
- 6 Mix in chamber by raising and lowering the leveling bulb, and allow one minute for complete laking
- 7 Place 2 ml of the hydrosulfite borate solution into the cup and draw 15 ml into the chamber
- 8 Mix again and allow to stand for one minute
- 9 Draw off and discard the excess hydrosulfite borate from the cup and add 18 ml of the carbon monoxide water into the chamber
- 10 Seal off with mercury and lower the mercury level in the chamber to the 50 ml mark
- 11 Cover the chamber with black paper and shake the solution for two minutes.
- 12 The evolved gases are quantitatively ejected. Seal again with mercury and shake again for two minutes more. Expel gases again ( $O_2$  and  $N_2$ ) and remove paper from chamber
- 13 Place 30 ml of air free ferricyanide buffer solution into the cup and take 10 ml into the chamber
- 14 Seal with mercury and discard remaining buffer &  $Fe(CN)_6^{3-}$  solution
- 15 Lower the mercury level to approximately 1 cm below the 50 ml mark and shake the mixture for a total of 5 minutes
- 16 Once during the shaking the motor is temporarily stopped and the mercury raised to the 50 ml mark, so as to mix the solution in the stem of the chamber with the main body of the solution. In this way destruction of the ferricyanide by the mercury is minimized
- 17 At the end of 5 minutes the gas is then compressed nearly to atmospheric pressure and 2.5 ml of the air free NaOH placed in the cup
- 18 Take 15 ml of the NaOH into the chamber and allow about one minute for complete carbon dioxide absorption. Try not to let the solution get above the 2 ml mark
- 19 Adjust the solution meniscus to the 2 ml mark and read the pressure (P)
- 20 Using a Poughton-Scholander syringe pipette quantitatively eject the gas sample into it for CO analysis by the technique described in the sections on Poughton-Scholander methods
- 21 Lower the solution below the 2 ml mark allow to drain and read again at the 2 ml mark (P)
- 22 Record temperature of the water jacket

### Calculation

- 1 ml total CO/100 ml blood =  
 $(P_1 - P - C) \times \text{temp factor for } O_2 \text{ Table 16}$
- 2 The correction for CO is calculated from the analysis of the sample in the syringe pipette

### Example

One ml of blood was analyzed from a tonometer in which the blood was saturated with CO. P was 134.2 P was 47.8 C was 0.5 Water jacket was at 21.6 C

ml total CO/100 ml blood =  
 $(134.2 - 47.8 - 0.5) \times 0.2445 = 20.84 \text{ ml CO/100 ml blood}$

### Precautions

- 1 Be sure to cover the chamber with black paper as light will destroy carbon monoxide heme complexes

## Oxygen in Blood (Roughton Scholander)

## Reference

Roughton F J W and Scholander, P F Micro Gasometric Estimation of the Blood Gases I Oxygen *J Biol Chem.* 143 541 550 (June) 1943

## Principle


By mixing of blood in a special syringe with ferricyanide containing potassium bicarbonate and saponin and adding an acetate buffer a gas phase is provided for the extraction of the blood gases in the same way the vacuum acts in the Van Slyke method. The CO thus produced is then absorbed with 10 percent NaOH. The residual gas bubble is driven into the syringe's capillary and its volume measured before and after absorption of O with alkaline pyrogallol. The difference in these two volumes represents the O content of the blood and reagents. Then a blank is run for the O content of the reagents alone this figure is subtracted from the former and the remainder when multiplied by the usual correction factor for temperature and pressure gives the O content of the blood.

## Apparatus

- 1 A one ml Pyrex tuberculin syringe with arresting clip and with standard bore precision 0.5 mm Pyrex capillary fused to its nozzle. The top of the capillary is expanded to a cylindrical cup of about 2.5 mm bore and 1.5 mm length. The capillary (7 to 8 cm length) is graduated into 30 divisions each of 2 mm length (Fig. 9).
- 2 Blood pipette made of thin walled glass tubing (1 to 1.5 mm bore) with a tip to fit into the bottom of the glass cup. The volume delivered from the mark to the tip is equal to 100 divisions of the capillary and is actually 39.3 cmm.
- 3 A detachable rubber cup of about 1 ml capacity is fitted to the top of the glass cup when required.
- 4 Several 2 ml and 5 ml syringes with fine tipped rubber connected glass nozzle to store the reagents.
- 5 Metal or glass beakers one of them containing water at room temperature to place the syringe for equilibration between readings.
- 6 Suction line and fine tipped glass point.
- 7 When suction system is not available a regular syringe or dropper can be used.
- 8 Lens paper or a piece of cloth to wipe the capillary between readings.
- 9 A piece of fine wire to remove air bubbles from the cup.
- 10 Thermometer for water bath.

## Reagents

- 1 Distilled water.
- 2 Caprylic alcohol.
- 3 Ferricyanide solution. 1.5 gm of  $K_3Fe(CN)_6$ , 3 gm of KIICO and 15 gm of saponin are dissolved in water and made up to 50 ml. It is not advisable to make larger amounts as the solution should be renewed every 3 days.
- 4 Acetate buffer. 70 gm of sodium acetate  $NaC_2H_3O_2 \cdot 3H_2O$  are dissolved in 100 ml of water and 15 ml of glacial acetic acid are added.
- 5 Urea 45 percent. This is used as a cleaning solution because of its protein dissolving property.
- 6 Sodium hydroxide 10 per cent.
- 7 Pyrogallol solution. 15 gm of powdered pyrogallol are added to 100 ml of 20 percent NaOH in a rubber stoppered bottle and covered with a layer of oil 1 cm thick. The pyrogallol is dissolved under the oil by stirring with a glass rod.
- 8 The distilled water is placed in a bottle above the apparatus and is delivered through rubber tubing with a glass nozzle. It can also be delivered by a 20 ml syringe with a needle or a glass nozzle.



Figs 29 31

Fig 29 (A B C) —Technique for absorption of excess of  $\text{CO}_2$ .

Fig 30 (D) —Technique for absorption of  $\text{CO}$  (or other absorbable gas) in glass cup.

Fig 31 (E) —Temperature equilibration of gas bubble in capillary.

- 9 Reagents 2 3 4 5 and 7 are conveniently stored in 1 or 5 ml syringes and Reagent 6 in a 10 or 20 ml syringe. These syringes have fine tipped glass nozzles attached.

# Procedure

- 1 The syringe is held vertically and any liquid in the cup of the syringe is withdrawn by suction. The plunger is pushed up and the cup filled with ferricyanide solution. The solution is drawn down to the bottom of the syringe and expelled through the cup and removed. This procedure is repeated twice with fresh lots of ferricyanide without trapping air bubbles. No grease or oil is used in the syringe.
- 2 The glass cup is filled to the mark with ferricyanide and the latter drawn down to the bottom of the cup.
- 3 A drop of caprylic alcohol is deposited on the bottom of the cup.
- 4 The pipette is filled with blood to the mark, wiped and cautiously introduced into the glass cup and its tip pressed snugly but not too vigorously against the bottom of the cup. The pipette should be held at a slight angle to the horizontal so that the blood does not run out when both ends of the pipette are open to the air (see Figs 30 31 32 and 33).
- 5 By pulling out the plunger gradually the blood is slowly and evenly drawn down into the capillary followed by a bubble of air of about 1 mm length. No appreciable amount of caprylic alcohol should be drawn in during this step. The bubble of air prevents any blood being sucked back when the pipette is removed. (See Fig 33 B).
- 6 The pipette is quickly removed and the bubble of air is then expelled through the caprylic alcohol with the aid of the wire.
- 7 A trace of caprylic alcohol is about two divisions length of the capillary is drawn down onto the top of the blood and the rest of the caprylic alcohol is sucked out from the cup.
- 8 The cup is filled to the mark with acetate buffer and the latter drawn down to the bottom of the cup.
- 9 The cup is then immediately filled to the top with 45 percent urea solution and then closed firmly with the finger.
- 10 The arresting clip is then removed and the closed apparatus vigorously shaken in the horizontal position. The plunger will be gradually drawn out by itself as the CO and other gases are evolved the gas pressure in the syringe being kept roughly atmospheric. The total volume evolved is usually about 0.5 ml. Shaking is continued for a total of two minutes (see Figs 33 and 34).
- 11 The arresting clip is put back in place.
- 12 The finger is cautiously released the gas meniscus being kept in the capillary by manipulating the plunger.
- 13 A small amount of urea is allowed to run down into the capillary and left there until the walls are perfectly clean.
- 14 Three quarters of the urea solution in the glass cup is removed by gentle suction.
- 15 The rubber cup is adjusted to the cup and filled with 20 per cent NaOH without trapping air bubble. Use the fine wire if necessary.
- 16 A little NaOH is drawn into the syringe. This absorbs some CO causing a partial vacuum which quickly sucks in more NaOH until only a small bubble consisting of O<sub>2</sub>, N<sub>2</sub> and CO (if any was originally present in the blood) is left at the top of the syringe. The absorption takes a few seconds and just before it is complete the residual bubble is screwed slowly and carefully up into the capillary by manipulation of the plunger.
- 17 The rubber cup is removed. As it still contains some NaOH it is advisable to remove the cup under a stream of running water. The glass cup is then emptied by suction.




Fig 3 —Transfer of blood direct into capillary




Fig 33 (A B C) —Further details of transfer of blood into capillary

- 18 The capillary is placed for half a minute in a beaker of water at room temperature and the temperature read (Fig 31)
- 19 It is then removed dried by light wiping care being taken that the capillary is not handled and the volume of the bubble read  $= V_1$  divisions.
- 20 The glass cup is then filled with pyrogallol solution and the O of the bubble is absorbed by pulling the gas bubble down to the bottom of the capillary and back again a few times. Finally the bubble is moved very slowly up into the top part of the capillary and after a further temperature equilibration its volume is read again  $V_2$  division. If  $V_2$  is only a few divisions the second temperature equilibration can be omitted.
- 21 A blank is obtained by following exactly the same steps except the adding of the blood sample.

### Calculations

The oxygen content of the blood equals

$$(V_1 - V_2 - c) \times f \quad \text{where}$$

- $V_1$  First reading after absorption of the CO with NaOH  
 $V_2$  Second reading after absorption of the O with pyrogallol  
 $c$  Blank correction for O content of reagents  
 $f$  Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from the line chart at the beginning of the section on blood gases (Fig 3)




Fig 34—Shaking of syringe and extraction of gas.

### Example

|  |       |   |
|--|-------|---|
| First reading  | 27.8  | ( $V_1$ )   |
| Second reading   | 5.5   | ( $V_2$ )   |
| Blank figure   | 1.1   | ( $c$ )   |
| Correction factor  | 0.000 | ( $f$ ) for $T = 2^\circ \text{C}$ and bar. pres. $= 60 \text{ mm}$ |
| $27.8 - 5.5 - 1.1$ multiplied by $0.000 = 21.2 \times 0.0 = 19.1 \text{ ml of O per 100 ml blood}$ |       |   |

### Precautions

- 1 When blood is drawn from the pipette into the syringe care should be taken that the pipette be pressed snugly but not too vigorously against the bottom of the cup in order to prevent breakage of the tip.
- 2 The arresting clip taken off at the beginning of the shaking should be put back in place before releasing finger.
- 3 Speedy adjustment of the bubble into the capillary at the end of the CO absorption helps to eliminate reabsorption of O at the critical moment at which the partial pressure of the  $O_2$  becomes high (i.e. nearly 1 atmosphere).
- 4 While there is gas in the syringe or bubbles in the capillary the hands should not touch either one so as to prevent any effect of the hands' warmth on the volume of the gases.



- 5 Slipping back of the plunger is prevented by a proper adjustment of the arresting clip
- 6 Temperature equilibration of the bubble before the measurement of  $V$  has proved necessary for with long bubbles a contraction of nearly a division may be observed as a result of immersion of the capillary in the water bath
- 7 To wash the instrument the plunger is pulled out under a stream of running water and the blood mixture poured out. The syringe is filled and emptied several times with water before the plunger is restored. The plunger should never be forced inwards if there is a resistance due to precipitates inside it. Occasionally the whole syringe should be rinsed with dichromate cleaning solution.

## 10 Dissolved Nitrogen in Blood (Edwards Scholander Roughton)

### Reference

Edwards G A Scholander P F and Roughton F J W *Micro Gasometric Estimation of the Blood Gases III Nitrogen* *J Biol Chem* 148 565 571 (June) 1943

### Principle

A  $\text{CO}_2$  phase is evolved by mixing the blood sample with bicarbonate and acid phosphate and thus the whole of the dissolved nitrogen of the blood and reagents is extracted as well as some oxygen. No  $\text{CO}$  is extracted. The  $\text{CO}$  and  $\text{O}$  are then absorbed with alkaline hydrosulphite and the bubble of  $\text{N}_2$  which is left is measured in the syringe capillary. The reading less a small blank from the reagents is multiplied by the usual correction factor for temperature and pressure and the figure so obtained gives the  $\text{N}_2$  content of the blood in volumes percent.

### Apparatus

- 1 One ml Pyrex tuberculin syringe with arresting clip and with a standard bore precision 0.5 mm Pyrex capillary fused to its nozzle. The top of the capillary is expanded to a cylindrical cup of about 2.5 mm bore and 1.5 cm length. The capillary (7 to 8 cm length) is graduated into 30 divisions each of 2 mm length.
- 2 Blood pipette made of thin walled glass tubing (1 to 1.5 mm bore) with a tip to fit into the bottom of the glass cup. The pipette has two marks one near the tip where it is calibrated to deliver an amount equal to 100 divisions of the capillary and which actually is 39.3 cmm. The other mark is calibrated to contain thrice as much as the former amount. (See Figs 29 34)
- 3 A detachable rubber cup of about 1 ml capacity is fitted to the top of the glass cup when required.
- 4 Several 2 ml and 5 ml syringes with fine-tipped rubber connected glass nozzles to store the reagents.
- 5 Metal or glass beakers one of them containing water at room temperature to place the syringe for equilibration between readings.
- 6 Suction line and fine tipped glass point.
- 7 A regular syringe or a dropper when suction system is not available.
- 8 Lens paper or a piece of cloth to wipe the capillary before readings.
- 9 A piece of fine wire to remove air bubbles from the cup.
- 10 Thermometer for water bath.

### Reagents

- 1 Aerated distilled water
- 2 Caprylic alcohol
- 3 Bicarbonate solution 11 gm of  $\text{KHCO}_3$  are dissolved in 100 gm of water
- 4 Acid phosphate buffer 95 gm of  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  are dissolved in 100 gm of warm water. The resultant volume is about 142 ml at 0°C and the strength about 5 M.
- 5 Urea 45 percent in water

- 6 Hydrosulfite solution 15 gm of a mixture of sodium hydrosulfite 10 parts and sodium anthraquinone beta sulfonate 1 part are added to 50 ml of 20 percent KOH. This solution should be made with the least possible contamination from air and stored in a rubber stoppered bottle.
- 7 Reagents <sup>a</sup> to 4 are stored in well greased ml syringes provided with glass tips. Reagents 5 and 6 may be stored in 10 ml syringes provided with glass tips. The distilled water may best be stored in a bottle with rubber tubing and glass nozzle hanging down over the apparatus. A syringe may serve for this purpose as well.

#### Procedure

- 1 The syringe is held vertically and any liquid in the cup of the syringe is withdrawn. The plunger is pushed up and the syringe cup filled with bicarbonate solution. The solution is then drawn down into the syringe barrel the plunger being pulled practically the whole way out. The plunger is then pushed up without any bubbles being trapped and the bicarbonate withdrawn from the syringe cup. The step is repeated once more. The dead space in the syringe normally contains enough bicarbonate to give 0.75 ml of CO<sub>2</sub> on subsequent treatment with the acid buffer.
- <sup>a</sup> The walls of the syringe cup are dried with a roll of cotton or cleaning paper and a drop of caprylic alcohol placed in the bottom of the cup without trapping any air bubbles.
- 3 The blood is measured from the second mark of the pipette into the capillary under the caprylic alcohol. With the syringe at the same angle the pipette is cautiously introduced into the glass cup and its tip pressed snugly but not too vigorously against the bottom of the cup.
- 4 The blood is drawn down into the syringe until the upper end of the column of blood is 7 mm below the top of the syringe capillary with a small bubble of air above it. Care must be taken that no caprylic alcohol is drawn down with the blood. The pipette is then removed the air bubble is expelled through the alcohol by gently screwing up the plunger and 2 divisions of caprylic alcohol are left on top of the blood in the capillary.
- 5 The cup is filled to the mark with the acid buffer and the latter drawn down to the bottom of the cup.
- 6 The cup is filled to the top with urea solution and is closed firmly with the finger.
- 7 The arresting clip is taken off the syringe barrel.
- 8 The syringe is shaken vigorously for 2 minutes in the horizontal position the plunger being withdrawn as the CO<sub>2</sub> is evolved. The gas pressure in the syringe must be kept at atmospheric pressure. The volume of the gas finally evolved must be between 0.6 and 0.5 ml.
- 9 The arresting clip is put back in place and then the finger tip is cautiously released the plunger being manipulated so as to keep the gas meniscus in the capillary.
- 10 The urea solution is drawn down to the bottom of the capillary without being allowed to enter the syringe barrel. This prevents subsequent contamination of the capillary wall with precipitates from the blood.
- 11 The syringe is held in the vertical position the rubber cup attached and about 1 ml of the hydrosulfite solution added without trapping air bubbles in the syringe cup.
- 12 A little of the hydrosulfite is gently drawn down into the syringe creating a vacuum which sucks in the rest of the solution needed for absorption of the CO<sub>2</sub> and such traces of O<sub>2</sub> as are present.
- 13 The bubble is pushed up into the lower part of the capillary. The hydrosulfite is sucked out of the rubber cup and the latter detached under a stream of running

water The syringe cup is then filled with water three fourths of which is drawn down over the bubble layering itself on top of the blood mixture

- 14 The bubble is pushed up into the clean capillary very gently and the capillary is placed in a beaker of water at room temperature to bring about temperature equilibrium The temperature is read.
- 15 The syringe is removed from the beaker, dried by light wiping care being taken that the capillary is not handled and the volume of N read V divisions
- 16 A blank is run to get the blank correction figure It is exactly the same procedure but no blood sample is used. Some difficulty may be encountered in first determining the blank since in the blank analysis the blood has a damping effect upon the reaction between the bicarbonate and the acid so that the  $\text{CO}_2$  is not evolved for 10 or more seconds after the acid is drawn down into the syringe In the blank analysis on the other hand the reaction begins immediately when the acid comes in contact with the bicarbonate and hence some bubbles of gas belonging in the syringe may be lost to the outside To avoid this difficulty one must
  - a Make sure that the 2 divisions of caprylic alcohol are not washed out when the acid is put in the cup
  - b Draw the acid down into the syringe barrel quickly
  - c Fill the cup with urea solution immediately after the acid is drawn down and then clamp the finger over the glass cup at once
 The blank figure has been found to be constant and of the order of 1.3 to 1.5 units on the capillary depending upon the instrument used

#### Calculations

The nitrogen content of the blood equals

$$\frac{(V - c) \times f}{3}$$

V = Uncorrected N<sub>2</sub> reading

c = Blank figure

f = Correction factor for temperature and barometric pressure (See STP chart at beginning of section on blood gases)

The Factor 3 is inserted because a volume of blood equal to 300 scale divisions is delivered in Step 3

#### Example

N reading = 52

Blank = 14

Corr factor 0.900 (for T = 20°C, bar press = 760 mm)

That is ml N<sub>2</sub>/100 ml blood =  $\frac{(52 - 14) \times 0.900}{3} = 11.4$

#### Precautions

- 1 When blood is drawn from the pipette into the syringe care should be taken that the pipette be pressed snugly but not too vigorously against the bottom of the cup in order to prevent the breakage of the tip
- 2 The arresting clip taken off at the beginning of the shaking should be put back in place before releasing finger
- 3 Slipping back of the plunger is prevented by a proper adjustment of the arresting clip
- 4 Temperature equilibration of the bubble before the measurement has proved necessary for with long bubbles a contraction of nearly a division may be observed as a result of the immersion of the capillary in the water bath
- 5 Great caution must be exercised in handling blood of low N content as contamination with atmospheric nitrogen occurs rapidly due to the high gradient of nearly 4/5 of an atmosphere

6 To wash the instrument the plunger is pulled out under a stream of running water and the blood mixture is poured out. The syringe is filled and emptied several times with water before the plunger is restored. The plunger should never be forced inward if there is a resistance due to precipitate inside it. Occasionally the whole syringe should be run with dichromate cleaning solution.

## 11 Carbon Monoxide in Blood (Scholander Roughton)

### Reference

Scholander P F and Poulton F J W. Micro Gasometric Estimation of the Blood Gases. II Carbon Monoxide. *J Biol Chem* 148 551-563 (June) 1943

(This paper describes how for the microdetermination of carbon monoxide in the blood the syringe can be used in three different ways: i.e.

- A. General method for saturation ranging from 0 to 100 per cent COHb.
- B. Method for combined determination of O and CO on one sample of blood.
- C. Special method precise enough for blood volume determination in which the CO content is kept below 0.5 volumes per cent.)

### A. GENERAL METHOD FOR ESTIMATION OF CO IN BLOOD

#### Principle

The O, CO and N of the blood and reagents are extracted by excess of CO evolved by the mixing of ferricyanide solution containing potassium bicarbonate and saponin with acetate buffer. The CO and O are subsequently absorbed together by alkaline pyrogallol. The volume of the residual gas is measured in the capillary of the instrument before and after the absorption of CO by Winkler's solution. The difference in volume multiplied by the correction factor for temperature and pressure gives the CO content of the blood in volumes percent. Since the blood itself is the only source of CO and the reagents contain none a blank correction is not necessary.

#### Apparatus

- 1 One ml. Pyrex tuberculin syringe with arresting clip and with a standard bore precision 0.5 mm. Pyrex capillary fused to its nozzle. The top of the capillary is expanded to a cylindrical cup of about 0.5 mm bore and 1.5 cm length. The capillary (1 to 8 cm length) is graduated into 30 divisions each of 0.2 mm length.
- 2 Blood pipette made of thin walled glass tubing (1 to 1.5 mm bore) with a tip to fit into the bottom of the glass cup. The volume delivered from the mark to the tip is equal to 100 divisions of the capillary and is actually 39.3 cmm. (See Figs. 79-84.)
- 3 A detachable rubber cup of about 1 ml capacity is fitted to the top of the glass cup when required.
- 4 Several 2 ml and 0.5 ml syringes with fine-tipped rubber connected glass nozzles to store the reagents.
- 5 Metal or glass testers, one of them containing water at room temperature to cover the syringe for equilibrations between readings.
- 6 Suction line and fine-tipped glass point.  
A rubber syringe or a dropper when suction system is not available.
- 7 Cleaning paper or a piece of cloth to wipe the capillary before readings.
- 8 A piece of fine wire to remove air bubbles from the cup.
- 9 A thermometer for the water bath.

#### Reagents

- 1 Aerated distilled water
- 2 Caprylic alcohol.

- 3 Ferrieyanide solution 12.5 gm of  $K_2Fe(CN)_6$  3 gm of  $KHCO_3$  and 1.5 gm of saponin are ground in a mortar and dissolved in water and made up to 100 ml. The solution should be renewed every three days
- 4 Acetate buffer 70 gm of sodium acetate ( $NaC_2H_3O_2 \cdot 3H_2O$ ) are dissolved in 100 gm of water and 15 ml glacial acetic acid added
- 5 45 per cent urea.
- 6 Winkler's solution 20 gm of cuprous chloride 25 gm of ammonium chloride and 75 gm of water are placed in a bottle just large enough to contain them. The bottle is corked, shaken with as little air as possible and the precipitate then allowed to settle. A coil of copper wire is placed in the solution which is then covered with a layer of heavy oil. After some time the reagent becomes almost colorless.
- 7 Pyrogallol solution 15 gm of powdered pyrogallol are added to 100 ml of 10 per cent  $NaOH$  in a rubber stoppered bottle and covered with a layer of oil. The pyrogallol is dissolved under the oil by stirring with a glass rod. The pyrogallol solution is stored in a 10 ml syringe with the rear of the plunger well greased and a fine tipped glass nozzle attached to the syringe. The other reagents are similarly stored in smaller syringes. The water is placed in a bottle above the apparatus and is delivered through rubber tubing with a glass nozzle attached. A 20 ml syringe with a needle or glass nozzle attached can be used for the delivery of the water.

#### Procedure

- 1 The syringe is held vertically and any liquid in the cup of the syringe is with drawn by suction. The plunger is pushed up and the cup filled with ferrieyanide solution. The solution is drawn down to the bottom of the syringe and expelled through the cup and removed. This procedure is repeated twice with fresh lots of ferrieyanide the dead space finally being left full of ferrieyanide without trapping any air bubbles. No grease or oil is used in the syringe.
- 2 The glass cup is filled to the mark with ferrieyanide and the latter drawn down to the bottom of the cup.
- 3 A drop of caprylic alcohol is deposited on the bottom of the cup.
- 4 The pipette is filled with blood to the mark wiped and cautiously introduced into the glass cup and its tip pressed snugly but not too vigorously against the bottom of the cup. The pipette should be held at a slight angle to the horizontal so that the blood does not run out when both ends of the pipette are open to the air.
- 5 By pulling out the plunger gradually the blood is slowly and evenly drawn into the capillary followed by a bubble of air of about 1 mm length. If the tip is properly ground and the right amount of pressure applied no appreciable caprylic alcohol is drawn in during this step. The bubble of air prevents any blood from being sucked back into the tip when the pipette is removed.
- 6 The pipette is quickly removed and the bubble of air is then expelled through the caprylic alcohol with the aid of the wire.
- 7 A trace of caprylic alcohol is about 2 divisions length of the capillary is drawn down onto the top of the blood and the rest of the caprylic alcohol is removed from the cup.
- 8 The cup is filled to the mark with acetate buffer and the latter drawn down to the bottom of the cup.
- 9 The cup is then immediately filled to the top with 45 per cent urea solution and then closed firmly with the finger.
- 10 The arresting clip is then removed and the closed apparatus vigorously shaken in the horizontal position. The plunger will be gradually drawn out by itself as the  $CO_2$  and other gases are evolved the gas pressure in the syringe being kept roughly atmospheric. The total volume evolved is usually about 0.5 ml. Shaking is continued for a total of 2 minutes.

- 11 The arresting clip is put back in place
- 12 The finger is cautiously released the gas meniscus being kept in the capillary by manipulating the plunger
- 13 A small amount of urea solution is allowed to run down into the capillary and left there until the walls are perfectly clean. Then the rest of the urea solution (three quarters) is removed.
- 14 The rubber cup is adusted and filled with pyrogallol solution. Use the fine wire to avoid trapping of air bubble
- 15 A little of the pyrogallol solution is drawn into the syringe. This absorbs some CO and O causing a partial vacuum which quickly sucks in more pyrogallol until only a small bubble consisting of N and CO (if any was originally present in the blood) is left at the top of the syringe. The absorption takes a few seconds and just before it is complete the residual bubble is screwed slowly and carefully up into the capillary by manipulation of the plunger
- 16 The rubber cup is removed. As it still contains a fair amount of pyrogallol it is advisable to remove cup under a stream of running water. The glass cup is then emptied and filled again with water
- 17 The capillary is placed for half a minute in a beaker of water at room temperature and the temperature is read.
- 18 It is then removed dried by light wiping care being taken that the capillary is not handled and the volume of the bubble read  $\frac{1}{4}$  divisions
- 19 The glass cup is flushed clean with water and left filled. About three-fourths of this water is pulled quickly down into the syringe forming a layer on top of the heavier blood mixture. The bubble with clean water below it is then at once run up to the top of the capillary
- 20 The glass cup is emptied of water and filled with Winkler's solution.
- 21 The syringe is then pointed with the cup downward the capillary making a slight angle with the horizontal
- 22 By cautiously screwing in the plunger the gas bubble is driven out into the glass cup where it rests near the junction of the capillary and the cup. As soon as the bubble is free in the cup Winkler's solution is sucked behind it so as to half fill the capillary. Gentle rotation for a few seconds completes the absorption of the CO
- 23 The syringe is then turned to the vertical position with the cup downward the gas bubble is sucked back into the capillary and its volume  $\frac{1}{4}$  measured.

#### Calculations

The carbon monoxide content of the blood in volumes per cent is

$$(\frac{1}{4} - \frac{2}{4}) \times f \text{ where}$$

$\frac{1}{4}$  = First reading after absorption of the CO, and O<sub>2</sub>

$\frac{2}{4}$  = Second reading after absorption of the CO

$f$  = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig. 23 at the beginning of the section on blood gases

#### Example

First reading  $\frac{1}{4}$  = 5.5

Second reading  $\frac{2}{4}$  = 4.2

Correction factor  $f$  = 0.900 at  $T = 22^\circ \text{C}$ , bar pres. = 760 mm.

ml CO/100 ml blood =  $(5.5 - 4.2) \times 0.900 = 1.2$

#### Precautions

1. When the blood is drawn from the pipette into the syringe care should be taken that the pipette is pressed snugly but not too vigorously against the bottom of the cup in order to prevent breakage of the tip.

- 3 Ferricyanide solution 12.5 gm of  $K_3Fe(CN)_6$  3 gm of  $KHCO_3$  and 1.5 gm of saponin are ground in a mortar and dissolved in water and made up to 50 ml. The solution should be renewed every three days
- 4 Acetate buffer 70 gm of sodium acetate ( $NaC_2H_3O_2 \cdot 3H_2O$ ) are dissolved in 100 gm of water and 15 ml glacial acetic acid added
- 5 45 per cent urea.
- 6 Winkler's solution 20 gm of cuprous chloride 25 gm of ammonium chloride and 75 gm of water are placed in a bottle just large enough to contain them. The bottle is corked shaken with as little air as possible and the precipitate then allowed to settle. A coil of copper wire is placed in the solution which is then covered with a layer of heavy oil. After some time the reagent becomes almost colorless
- 7 Pyrogallol solution 15 gm of powdered pyrogallol are added to 100 ml of 0 per cent NaOH in a rubber stoppered bottle and covered with a layer of oil. The pyrogallol is dissolved under the oil by stirring with a glass rod. The pyrogallol solution is stored in a 10 ml syringe with the rear of the plunger well greased and a fine tipped glass nozzle attached to the syringe. The other reagents are similarly stored in smaller syringes. The water is placed in a bottle above the apparatus and is delivered through rubber tubing with a glass nozzle attached. A 20 ml syringe with a needle or glass nozzle attached can be used for the delivery of the water

#### Procedure

- 1 The syringe is held vertically and any liquid in the cup of the syringe is withdrawn by suction. The plunger is pushed up and the cup filled with ferricyanide solution. The solution is drawn down to the bottom of the syringe and expelled through the cup and removed. This procedure is repeated twice with fresh lots of ferricyanide, the dead space finally being left full of ferricyanide without trapping any air bubbles. No grease or oil is used in the syringe
- 2 The glass cup is filled to the mark with ferricyanide and the latter drawn down to the bottom of the cup
- 3 A drop of caprylic alcohol is deposited on the bottom of the cup
- 4 The pipette is filled with blood to the mark, wiped and cautiously introduced into the glass cup and its tip pressed snugly but not too vigorously against the bottom of the cup. The pipette should be held at a slight angle to the horizontal so that the blood does not run out when both ends of the pipette are open to the air
- 5 By pulling out the plunger gradually the blood is slowly and evenly drawn into the capillary followed by a bubble of air of about 1 mm length. If the tip is properly ground and the right amount of pressure applied no appreciable caprylic alcohol is drawn in during this step. The bubble of air prevents any blood from being sucked back into the tip when the pipette is removed.
- 6 The pipette is quickly removed and the bubble of air is then expelled through the caprylic alcohol with the aid of the wire
- 7 A trace of caprylic alcohol i.e. about 2 divisions length of the capillary is drawn down onto the top of the blood and the rest of the caprylic alcohol is removed from the cup
- 8 The cup is filled to the mark with acetate buffer and the latter drawn down to the bottom of the cup
- 9 The cup is then immediately filled to the top with 45 per cent urea solution and then closed firmly with the finger
- 10 The arresting clip is then removed and the closed apparatus vigorously shaken in the horizontal position. The plunger will be gradually drawn out by itself as the  $CO$  and other gases are evolved the gas pressure in the syringe being kept roughly atmospheric. The total volume evolved is usually about 0.75 ml. Shaking is continued for a total of 2 minutes

- 11 The arresting clip is put back in place.
- 12 The finger is cautiously released the gas meniscus being kept in the capillary by manipulating the plunger
- 13 A small amount of urea solution is allowed to run down into the capillary and left there until the walls are perfectly clean. Then the rest of the urea solution (three quarters) is removed.
- 14 The rubber cup is adjusted and filled with pyrogallol solution. Use the fine wire to avoid trapping of air bubbles.
- 15 A little of the pyrogallol solution is drawn into the syringe. This absorbs some CO and O<sub>2</sub> causing a partial vacuum which quickly sucks in more pyrogallol until only a small bubble consisting of N<sub>2</sub> and CO (if any was originally present in the blood) is left at the top of the syringe. The absorption takes a few seconds and just before it is complete the residual bubble is screwed slowly and carefully up into the capillary by manipulation of the plunger.
- 16 The rubber cup is removed. As it still contains a fair amount of pyrogallol it is advisable to remove cup under a stream of running water. The glass cup is then emptied and filled again with water.
- 17 The capillary is placed for half a minute in a beaker of water at room temperature and the temperature is read.
- 18 It is then removed dried by light wiping care being taken that the capillary is not handled and the volume of the bubble read in divisions.
- 19 The glass cup is flushed clean with water and left filled. About three fourths of this water is pulled quickly down into the syringe forming a layer on top of the heavier blood mixture. The bubble with clean water below it is then at once run up to the top of the capillary.
- 20 The glass cup is emptied of water and filled with Winkler's solution.
- 21 The syringe is then pointed with the cup downward, the capillary making a slight angle with the horizontal.
- 22 By cautiously screwing in the plunger the gas bubble is driven out into the glass cup where it rests near the junction of the capillary and the cup. As soon as the bubble is free in the cup Winkler's solution is sucked behind it so as to half fill the capillary. Gentle rotation for a few seconds completes the absorption of the CO.
- 23 The syringe is then turned to the vertical position with the cup downward the gas bubble is sucked back into the capillary and its volume is measured.

#### Calculations

The carbon monoxide content of the blood in volumes per cent is

$$(\text{V} - \text{V}') \times f \text{ where}$$

V = First reading after absorption of the CO and O<sub>2</sub>.

V' = Second reading after absorption of the CO

f = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig. 4.5 at the beginning of the section on blood gases.

#### Example

First reading V = 5.5

Second reading V' = 4.

Correction factor f = 0.900 at T = 37°C bar press = 760 mm

ml CO/100 ml blood = (5.5 - 4.2) × 0.900 = 1.2

#### Precautions

- 1 When the blood is drawn from the pipette into the syringe care should be taken that the pipette is pressed snugly but not too vigorously against the bottom of the cup in order to prevent breakage of the tip.



- 2 The arresting clip taken off at the beginning of the shaking should be put back in place before the releasing of the finger
- 3 Speedy adjustment of the bubble into the capillary at the end of the CO and O absorption with the pyrogallol helps to eliminate reabsorption of CO at the critical moment at which the partial pressure of the CO becomes high (i.e. nearly 1 atmosphere)
- 4 While there is gas in the syringe or bubbles in the capillary the hands should not touch either one so as to prevent the effect of the hands' warmth on the volume of the gases
- 5 Slipping back of the plunger is prevented by a proper adjustment of the arresting clip
- 6 Temperature equilibration of the bubble before the measurement of V has proved necessary for with long bubbles a contraction of nearly a division may be observed as a result of immersion of the capillary in the water bath
- 7 To wash the instrument the plunger is pulled out under a stream of running water and the blood mixture poured out. The syringe is filled and emptied several times with water before the plunger is restored. The plunger should never be forced inward if there is a resistance due to precipitates inside it. Occasionally the whole syringe should be rinsed with dichromate cleaning solution.

## **B SPECIAL METHOD FOR COMBINED ESTIMATION OF CO AND O IN SINGLE SAMPLES OF BLOOD**

### **Principle**

When it is desired to measure both oxygen and carbon monoxide a single sample of blood can be used.

### **Procedure**

- 1 Follow Steps 1 through 10 in the procedure for oxygen and when  $V$  and  $V_1$  are obtained
- 2 The glass cup and capillary are thoroughly washed with water no traces of the pyrogallol solution in the capillary in the cup being left
- 3 About three fourths of the clean water left in the cup is quickly pulled down into the syringe forming a layer on top of the heavier blood mixture. The bubble with clean water below it is then at once run up to the top of the capillary
- 4 The glass cup is emptied of water and filled with Winkler's solution
- 5 The syringe is then held with the cup downward, the capillary making a slight angle with the horizontal
- 6 By cautiously screwing in the plunger the gas bubble is driven out into the glass cup where it rests near the junction of the capillary and the cup. As soon as the bubble is free in the cup Winkler's solution is sucked behind it so as to half fill the capillary. Gentle rotation for a few seconds completes the absorption of the CO
- 7 The syringe is then turned to the vertical position with the cup downward and the gas bubble is sucked back into the capillary and its volume  $V$  is measured.

### **Calculations**

$$(V_1 - V_2 - c) \times f = \text{O content}$$

$$(V - V_1) \times f = \text{CO content}$$

### **Where**

$V$  = First reading after absorption of the CO with NaOH

$V_1$  = Second reading after absorption of the O with pyrogallol

$V_2$  = Third reading after absorption of the CO with Winkler's

$f$  = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig 25

$c$  = Blank correction which in this case can be taken as 1.1

**Example**

First reading = 68 (V)  
 Second reading = 55 (V)  
 Third reading = 40 (V)  
 Blank figure = 11 (c)  
 Corr factor = 0.000 (f) for Temp °C bar press 760 mm  
 ml O<sub>2</sub> per 100 ml blood = (68 - 55 - 11) (0.000) = 18  
 ml CO per 100 ml blood = (55 - 40) (0.000) = 14

**Precautions**

All precautions given for the general method are applied here

### C. SPECIAL METHOD FOR BLOOD CONTAINING SMALL AMOUNTS OF CARBON MONOXIDE

When the CO content of the blood is suspected to be small i.e. that it does not exceed 0.0 volumes percent this special method is used which requires three times (10 cmm) the amount of the blood sample used in the general method. Thus an accuracy of 0.03 to 0.05 volume percent is obtained.

**Principle**

The gases of the blood are extracted by means of the CO gas phase produced by adding an acetate buffer to a ferricyanide solution containing sodium bicarbonate and saponin in the usual way described in the General Method. A special method is necessary because three times more blood is used than in the regular method and the relatively large amount of O<sub>2</sub> liberated may cause an appreciable formation of CO during the pyrogallol absorption whereas when only 40 cmm are used as in the regular method no such CO formation has been detected. To avoid this extra formation of CO a weakly alkaline hydrosulfite is employed which will absorb the O<sub>2</sub> of the blood.

**Apparatus**

The same apparatus described in the general method is used here. However the pipette should have a mark where it has been calibrated to deliver thrice the usual blood volume.

**Reagents**

Those reagents described under the general method for estimation of CO in blood are used together with the following:

- 1 Hydrosulfite borate solution. A stock solution of 3 per cent borax ( $\text{Na}_2\text{B}_4\text{O}_{10} \cdot 10\text{H}_2\text{O}$ ) is prepared. 1 gm of sodium hydrosulfite ( $\text{Na}_2\text{S}_2\text{O}_4$ ) is dissolved in 50 ml of the borax solution and stored in a glass-tipped syringe. The final solution should not be kept more than 2 days.

**Procedure**

- 1 The glass cup of the syringe is filled with hydrosulfite borate solution which is drawn down to the bottom of the syringe and then expelled leaving the dead space and the capillary filled with the solution.
- 2 The blood is drawn down the capillary directly on top of the borate solution from a pipette calibrated to deliver 3 times the usual quantity of blood. No caprylic alcohol is needed at this stage.
- 3 The blood is drawn down to the bottom of the capillary and the syringe rotated horizontally round its axis so as to mix the hydrosulfite with the blood.
- 4 The syringe is then restored to the vertical position and the blood solution brought to the top of the capillary.
- 5 The cup is filled to the top with ferricyanide solution and the whole of this drawn down.

- 6 Two or more divisions of caprylic alcohol are placed on top of the ferricyanide solution in the capillary so as to protect it from contact with the acetate buffer used in the next step
- 7 The glass cup is rinsed with water and filled to the mark with acetate buffer
- 8 The acetate buffer is then drawn down a few divisions below the bottom of the cup. The cup is then filled with water and closed at the top with the tip of the finger as soon as the few divisions of air have been expelled from the top of the capillary
- 9 The arresting clip is taken off and the gases are extracted in the usual way by 3 minutes shaking
- 10 The arresting clip is put back in place
- 11 The syringe is then turned to the vertical position with the cup downward and the plunger pushed in so as to give a slight positive pressure
- 12 With the syringe held over a dish the finger tip is removed gradually from the top of the cup and fluid is expelled until the meniscus is halfway down the capillary. The gas containing part of the syringe chamber must not be touched with the warm hand
- 13 The glass cup is filled to the top with water and then slanted downward and the gas bubble pushed to the end of the capillary. The blood mixture passes into the cup and settles through the water
- 14 With the instrument still slanting water is sucked into the capillary
- 15 The instrument is then restored to the usual position (vertical) with the cup upward and the latter filled with urea solution
- 16 The capillary is filled with urea solution and is left until the walls are perfectly clean
- 17 Three quarters of the urea in the glass cup is removed and the rubber cup adjusted and filled with pyrogallol solution without trapping air bubbles
- 18 A little pyrogallol is drawn down into the syringe. This absorbs some CO and O causing a partial vacuum which quickly sucks in more pyrogallol until only a small bubble consisting of N and CO is left at the top of the syringe. The absorption takes a few seconds and just before it is complete the residual bubble is screwed slowly and carefully up into the capillary by manipulation of the plunger
- 19 The rubber cup is removed. As it still contains a fair amount of pyrogallol it is advisable to remove the cup under a stream of running water. The glass cup is then emptied and refilled with water
- 20 The capillary is placed for half a minute in a beaker containing water at room temperature and the temperature is read
- 21 It is then removed dried by light wiping care being taken that the capillary is not handled and the volume of the bubble read *V* divisions
- 22 The glass cup is flushed clean with water and left filled. About three quarters of this water is pulled quickly down into the syringe forming a layer on top of the heavier blood mixture. The bubble with clean water below it is then at once run up to the top of the capillary
- 23 The glass cup is emptied of water and filled with Winkler's solution
- 24 The syringe is then held with the cup downward the capillary making a slight angle with the horizontal
- 25 By cautiously screwing in the plunger the bubble is drawn into the glass cup where it rests near the junction of the capillary and the cup. As soon as the bubble is free in the cup Winkler's solution is sucked behind it so as to fill the capillary. Gentle rotation for a few seconds completes the absorption of the CO
- 6 The syringe is then turned to the vertical position with the cup downward and the gas bubble is sucked back into the capillary and its volume *V* measured.

## Calculations

$$\frac{(V - v) \times f}{3} = \text{CO content of blood}$$

## Where

$V$  = First reading after absorption of the CO and O

$v$  = Second reading after absorption of the CO

$f$  = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig 25 at the beginning of the section on blood gases

(As thrice the usual amount of blood has been the usual result has to be divided by 3)

## Example

First reading 65 ( $V$ )

Second reading 50 ( $v$ )

Corr factor 0.900 ( $f$ ) When Temp = C Bar press 760 mm

$$\frac{(65 - 50) \times 0.900}{3} = 0.45 \text{ ml CO per 100 ml blood}$$

## Precautions

All those given for the general CO method are observed here

## 12 Carbon Dioxide in Blood (Scholander Roughton)

## Reference

Scholander P F and Roughton F J W Micro Gasometric Estimation of the Blood Gases IV Carbon Dioxide J Biol. Chem. 148 573 580 (June) 1943

## Principle

While carbon dioxide is by comparison ideal as an extractant for the microestimation of the other gases and there is no practicable method using a gas other than CO for blood a modification of the technique was necessary in order to use the convenient syringe method for the estimation of CO. This has been accomplished by the introduction of the vacuum extraction to the general method.

The blood sample and an acid buffer are vacuum extracted in the syringe by closing the capillary with the especially constructed rubber tipped plug and drawing the plunger out to the position fixed by a spacer (see Apparatus) and then shaking for 2 minutes. The vacuum is released by letting in the plunger and the volume of gas is measured before and after absorption with alkali. The reabsorption of CO during the compression is one of the special difficulties but it is possible to check it by (a) using an acid buffer reagent of minimal CO solubility and (b) clearing the capillary of liquid before the shaking begins. When the plunger moves up the whole of the gas therefore goes directly into the well drained capillary so that during the final stages of the contraction when the pressure of CO is approaching atmospheric the gas is only in close contact with the minute amount of liquid contained in the drainage films on the walls of the capillary.

## Apparatus

1 One ml Pyrex tuberculin syringe with arresting clip and with a standard bore precision 0.1 mm Pyrex capillary fused to its nozzle. The top of the capillary is expanded to a cylindrical cup of about 5 mm bore and 1.5 cm length. The capillary (7 to 8 cm length) is graduated into 30 divisions each of .1 mm length. An additional mark scratched at 33.3 has proved convenient when in an alternative procedure is desired to use exactly one third of the normal pettite load.

- 2 Blood pipette made of thin walled glass tubing (1 to 1.5 mm bore) with a ground tip to fit into the bottom of the glass cup. The volume delivered from the mark to the tip is equal to 100 divisions of the capillary and is actually 39.3 mm.
- 3 Rubber tipped wooden plug to seal the capillary during the vacuum extraction. The end of a round toothpick is dipped in rubber latex and a small drop left on the tip. It is dried tip downward at a moderate temperature in a drying oven. About 1 inch of the coated plug is used as in Fig 35, A.
- 4 A spacer for holding out the syringe plunger in a fixed position during the vacuum extraction. A piece of light sheet metal about 1.5 cm wide and 5.5 cm long is folded into a U shaped channel. To apply the vacuum one end of the spacer is pressed against the expanded lower end (button) of the plunger and the latter is then drawn out until the other end of the spacer presses against the butt end of the syringe barrel. The length of the spacer is so adjusted as to give a gas phase of about 0.75 ml in the barrel. (See Figs 35A and B).
- 5 Several 2 ml and 5 ml syringes with fine tipped rubber connected glass nozzles to store the reagents.
- 6 Metal or glass beakers one of them containing water at room temperature to hold the syringe for equilibration between readings.
- 7 Suction line with a very fine tipped glass point.
- 8 Thermometer for water bath.




Fig 35—Carbon Dioxide in blood

A Rubberized wooden plug for vacuum sealing glass cup

B Spacer for keeping plunger extended in fixed position during vacuum extraction.

#### Reagents

- 1 CO free distilled water (Ordinary water can be freed of CO by boiling with a drop of sulfuric acid.)
- 2 Caprylic alcohol.
- 3 Acid buffer 95 gm of acid sodium phosphate  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  are dissolved in 100 gm of warm water to make an almost saturated solution. The strength is about 2 M and the solubility of CO therein is only about one twelfth that in water.
- 4 10 percent NaOH
- 5 Glycerol

Reagents 3, 3 and 4 may best be stored in 2 or 5 ml syringes with glass tips attached. The glycerol is best kept in a 2 ml syringe having a No. 20 needle the point of which has been ground off. The distilled water may be stored in a bottle with rubber tubing and glass tip attached hanging above the bench.

## Procedure

- 1 The clean syringe is taken apart and the water shaken out of the barrel and the plunger dried. The rear part of the plunger is lubricated with a few streaks of glycerol and returned to the moist barrel
- 2 The cup is filled with distilled water and the latter is drawn one fourth down the barrel and expelled through the cup leaving the dead space of the syringe full of water without trapping any air bubbles
- 3 The plunger is pulled out very slightly so that the water meniscus at the bottom of the cup is lowered 1 or 2 mm down into the capillary
- 4 The blood pipette is held against the opening of the capillary trapping a small air bubble. The blood is drawn down very slowly into the capillary the air bubble separating it from the water
- 5 If the blood supply is ample the quickest procedure is to dry the cup with a little cotton on a toothpick and then half fill it with blood trapping the small air bubble in the capillary. The blood is then drawn down into the capillary the top layer being used as a protective against CO loss. If the blood is stored in a syringe it is convenient to attach a capillary glass tip to the nozzle and then after good mixing and discarding of a few drops through the tip the cup is filled to the mark directly from the syringe tip
- 6 With the blood meniscus at the 30 to 35 mark the pipette is detached and the blood in the cup sucked off. The upper blood meniscus is slowly moved to the zero mark and the amount of blood read in divisions
- 7 (As an alternative procedure the blood is moved down to a special mark scratched at 33.3 divisions. With a fine suction tip the upper meniscus is adjusted exactly to the zero mark the lower meniscus of the blood being simultaneously kept at the 33.3 mark. In this way exactly one third of the normal pipette load is used and the volume of CO found here after the usual correction for temperature etc simply to be multiplied by 3 to give the CO content in volumes percent.)
- 8 A drop of caprylic alcohol is deposited on the bottom of the cup and the bubble of air above the blood ejected through the capillary with the aid if necessary of a piece of fine wire. Two divisions of caprylic alcohol are drawn down onto the top of the blood in the capillary. The remainder of the caprylic alcohol in the cup is removed.
- 9 The cup is filled to the mark with acid phosphate which is then pulled down very slowly into the syringe until the upper meniscus is 1 mm below the bottom of the cup
- 10 The rubber end of the wooden plug is moistened with the phosphate buffer and with a few drops adhering to it is inserted in the bottom of the cup trapping a small air bubble
- 11 With the plug resting loosely against the bottom of the cup the air bubble is gently screwed up until it touches the rubber tip. The plug is then pressed against the bottom of the cup leaving the air bubble in direct contact with the rubber. The free lower end of the plug is covered by the drop of phosphate. The capillary is kept closed in this way during steps (12) to (18) with the left hand
- 12 Fresh glycerol is introduced into the plunger bearing with the right hand
- 13 One end of the metal spacer is put in place around the plunger under the plunger head the right hand being used and is held there with the other end sticking out at an angle. The cup end of the syringe is pointed slightly upward and kept upward during steps (14) to (17)
- 14 With the spacer and the plunger button held the plunger is slowly moved out in such a way that the fluid meniscus under the stopper moves down the capillary very slowly and evenly just as for a reading. When the capillary and its barrel

opening are cleared and drained of the fluid, the plunger is drawn slowly out the free end of the spacer being simultaneously moved in until it rests against the butt of the syringe barrel

- 15 Glycerol is added to the plunger bearing
- 16 The syringe is shaken with the cup end upward so as to prevent any fluid from blocking the entrance to the capillary. If nevertheless the capillary should become bridged over it can be cleared by warming up the capillary with the hand. Shaking is continued for 2 minutes. In case of foaming the plunger is released so that the fluid and extracted gas are compressed into the top part of the syringe barrel where there are normally residual traces of caprylic alcohol. When the plunger is drawn out again the foam generally disappears.
- 17 With the entire capillary free from fluid the plunger is pulled out slightly so as to allow the spacer to fall out. The plunger is then allowed to rise up within the barrel at a controlled but rather rapid speed until the lower end of the meniscus is inside the capillary and the pressure of the gas is at atmospheric pressure. If the entrance to the capillary should get bridged by fluid while the plunger is being let in the speed of the plunger must be so adjusted that the bridge moves very slowly up the capillary thus enabling proper drainage to take place.
- 18 When the gas bubble is at atmospheric pressure the plug is removed and the upper meniscus moved slowly and evenly to the zero mark.
- 19 The capillary is placed for  $\frac{1}{4}$  minute in a beaker of water at room temperature and the temperature is read. It is then removed and dried by light wiping care being taken that the capillary is not handled.
- 20 The volume of the bubble is read  $\nabla$  divisions without any unnecessary adjustments as this may lead to appreciable reabsorption (the capillary should have been properly drained during step [14]).
- 21 The cup is filled to the top with water and three fourths of the latter is pulled down into the syringe forming a layer on top of the blood mixture. The bubble is returned to the capillary with water beneath it. Owing to absorption of CO by this washing the bubble is much shorter.
- 22 The cup is filled with 10 percent NaOH the cup pointed downward and the bubble expelled into the alkali some of which is drawn into the capillary as soon as the bubble is free.
- 23 The cup is pointed slightly down and rotated a few times so as to complete the absorption of CO. The syringe is then returned to the vertical position with the cup downward and then a gas bubble is sucked back into the capillary.
- 24 The temperature of the capillary is again adjusted and the bubble volume read  $\nabla$ .

#### Calculations

The original formula to obtain the CO content of the blood by this method was as follows

$$\text{CO content} = (\nabla - \nabla') \times f \times \frac{100}{b}$$

Where

$\nabla$  = First reading

$\nabla'$  = Second reading after absorption of the CO

$f$  = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig. 25 at the beginning of the section on blood gases

$b$  = Volume of blood expressed as the number of capillary divisions

No allowance was made however for the incomplete extraction of CO during shaking and for reabsorption during compression. It was necessary therefore to run a series of determinations of CO by the syringe method and compare them with the results obtained by more precise a method like the Van Slyke. By dividing the





opening are cleared and drained of the fluid the plunger is drawn slowly out, the free end of the spacer being simultaneously moved in until it rests against the butt of the syringe barrel

- 15 Glycerol is added to the plunger bearing
- 16 The syringe is shaken with the cup end upward so as to prevent any fluid from blocking the entrance to the capillary. If nevertheless the capillary should become bridged over it can be cleared by warming up the capillary with the hand. Shaking is continued for 2 minutes. In case of foaming the plunger is released so that the fluid and extracted gas are compressed into the top part of the syringe barrel where there are normally residual traces of caprylic alcohol. When the plunger is drawn out again the foam generally disappears.
- 17 With the entire capillary free from fluid the plunger is pulled out slightly so as to allow the spacer to fall out. The plunger is then allowed to rise up within the barrel at a controlled but rather rapid speed until the lower end of the meniscus is inside the capillary and the pressure of the gas is at atmospheric pressure. If the entrance to the capillary should get bridged by fluid while the plunger is being let in the speed of the plunger must be so adjusted that the bridge moves very slowly up the capillary thus enabling proper drainage to take place.
- 18 When the gas bubble is at atmospheric pressure the plug is removed and the upper meniscus moved slowly and evenly to the zero mark.
- 19 The capillary is placed for  $\frac{1}{2}$  minute in a beaker of water at room temperature and the temperature is read. It is then removed and dried by light wiping care being taken that the capillary is not handled.
- 20 The volume of the bubble is read  $\frac{1}{2}$  divisions without any unnecessary adjustments as this may lead to appreciable reabsorption (the capillary should have been properly drained during step [14]).
- 21 The cup is filled to the top with water and three fourths of the latter is pulled down into the syringe forming a layer on top of the blood mixture. The bubble is returned to the capillary with water beneath it. Owing to absorption of CO by this washing the bubble is much shorter.
- 22 The cup is filled with 10 percent NaOH, the cup pointed downward, and the bubble expelled into the alkali some of which is drawn into the capillary as soon as the bubble is free.
- 23 The cup is pointed slightly down and rotated a few times so as to complete the absorption of CO. The syringe is then returned to the vertical position with the cup downward and then a gas bubble is sucked back into the capillary.
- 24 The temperature of the capillary is again adjusted and the bubble volume read  $\frac{1}{2}$ .

#### Calculations

The original formula to obtain the CO<sub>2</sub> content of the blood by this method was as follows

$$\text{CO content} = (V - V') \times f \times \frac{100}{b}$$

Where

$V$  = First reading

$V'$  = Second reading after absorption of the CO

$f$  = Correction factor for temperature aqueous vapor pressure and barometric pressure as derived from Fig. 25 at the beginning of the section on blood gases

$b$  = Volume of blood expressed as the number of capillary divisions

No allowance was made however for the incomplete extraction of CO during shaking and for reabsorption during compression. It was necessary therefore to run a series of determinations of CO by the syringe method and compare them with the results obtained by more precise a method like the Van Slyke. By dividing the

## Reagents

- 1 Glycerol aquaresein mixture This mixture is made by mixing equal parts of the two reagents and storing it in a two ml syringe with a # 0 needle the tip of which has been ground off square (Aquaresein is obtainable from the Glyco Products Co Inc New York.)
  - 2 Ferricyanide solution. Dissolve 1 gm of potassium ferricyanide and 0.5 gm of saponin in 10 ml of distilled water. Then add 0.1 ml of the acetate buffer to the solution.
  - 3 Caprylic alcohol Reagent grade
  - 4 Acid phosphate buffer Dissolve 95 gm of acid sodium phosphate ( $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ ) in 100 ml of distilled water. This is almost a saturated solution.
  - 5 Sodium hydroxide a 10% solution in water
  - 6 Distilled water
  - 7 Pyrogallol solution Add 15 gm of powdered pyrogallol to 100 ml of 20% sodium hydroxide solution in a rubber stoppered bottle and cover the solution with a layer of oil 2 cm thick. The pyrogallol is dissolved by stirring it with a glass rod passing it through the oil.
- Reagents 1 and 7 are kept for current use in 4 or 5 ml syringes with glass tips attached through the rubber tubing.

## Procedure

- 1 Remove the plunger from the syringe. Rinse it well with hot water and dry it by suction. This is easily done by attaching a small glass tube to the vacuum line and inserting the tube to the bottom of the barrel so that air is pulled through the warm capillary and barrel.
- 2 Add a drop of glycerol aquaresein mixture into the barrel opening. Lubricate the dry plunger with a few streaks of the glycerol solution and return it to the barrel adding enough of the solution so that no air is carried in with the plunger. Unless the plunger is moved gently in the barrel the latter may break because of the viscosity of the glycerol solution. When the plunger strikes the bottom the glycerol solution should rise only a few mm in the capillary.
- 3 Fill the cup with the ferricyanide solution draw it down into the barrel and expel it without trapping any air bubbles. Remove the excess ferricyanide leaving only the dead space filled with the ferricyanide solution.
- 4 Place a drop of caprylic alcohol on top of the ferricyanide solution draw down a scale divisions of the alcohol and suck out the remainder. Draw out the plunger until the fluid meniscus is lowered 1 to 2 mm into the capillary.
- 5 Seat the plastic tipped wire plug in the cup trapping a small air bubble on top of the caprylic alcohol. The rubber band must be strong enough to keep the capillary firmly closed by the plastic tip. Add a little of the ferricyanide solution into the cup to complete the seal.
- 6 Pull the plunger out slowly until the spring clip engages the notch in the plunger. Then add glycerol aquaresein solution around the plunger bearing. If the caprylic alcohol stays in the capillary it may be moved down the barrel by warming the capillary behind it in the hand. Allow a few seconds for extraction and release the vacuum by letting the plunger in. Expel the bubble and repeat the extraction.
- 7 Lower the meniscus of the extracted fluid 1 to 2 mm into the capillary. Place the tip of the blood pipette firmly against the bottom of the cup trapping a small (0.5 to 1 division) air bubble between the blood and the extracted fluid. Draw the blood slowly in until the length (A) of the air bubble can be measured. Continue to draw the blood in very slowly and evenly until 13 to 20 divisions have entered. Remove the blood pipette suck off the excess blood adjust the up

Oxygen is released from chemical combination by ferricyanide and carbon dioxide is evolved by the addition of acid phosphate buffer. The gases are vacuum extracted in the syringe selectively absorbed and measured in the capillary.

#### Apparatus

- 1 The same syringe analyzer as used for the Scholander Roughton Blood Gas analysis method with modifications shown in Fig 36
- 2 The syringe analyzer described by Scholander and Roughton has been modified to make it easy to hold the plunger out during the vacuum extraction by providing a notch in the plunger which engages a spring steel clip. To cut the notch draw the plunger out 0.8 ml and mark it 3 to 4 mm past the syringe barrel. Cut the notch in the plunger with a triangular medium or fine India oil stone straight below the number of the plunger slanting the stone so that the spring clip will not disengage the notch during the vacuum extraction.

#### PILLARY SYRINGE AND SLING FOR THE COMBINED CARBON DIOXIDE AND OXYGEN ME




Fig 36

- 3 A plastic tipped wire plug to seal the capillary during the vacuum extraction. The tip is fashioned from a piece of 15 to 2 mm plastic cord by rotating and gently heating it in front of a micro flame and pulling it apart when the plastic becomes soft. The pointed tip so formed is cut off to a length of 2 to 3 mm pierced straight with a fine wire which is bent. A rubber band attaches the wire plug to the syringe barrel.
- 4 A wire clip to give the plunger a suitable friction and to hold the plunger out during vacuum extraction. A 6 to 7 mm wide ring is cut off from a piece of thick walled rubber tubing and is slipped firmly onto the syringe. A straight angle bow of 0.4 to 0.5 mm of music wire wide enough to slip over the syringe collar is stuck straight through the rubber tubing. With the clip engaged in the plunger notch the end of the plunger should read 80 (0.8 ml).
- 5 A hand sling for the syringe for whirling down liquid bridges in the capillary and for breaking up bubbles in the syringe. The sling is 4 cm wide by 50 cm long made from 15 mm sheet aluminum. It can be whirled around a handle. The other end is bent out at a straight angle. It is furnished with a plastic cradle which carries a screw hook padded with rubber tubing for fastening the syringe.

## Reagents

- 1 Glycerol aquaresin mixture This mixture is made by mixing equal parts of the two reagents and storing it in a two ml syringe with a # 0 needle the tip of which has been ground off square (Aquaresin is obtainable from the Glyco Products Co Inc. New York)
  - 2 Ferricyanide solution Dissolve 1 gm of potassium ferricyanide and 0.5 gm of saponin in 10 ml of distilled water. Then add 0.1 ml of the acetate buffer to the solution
  - 3 Caprylic alcohol Reagent grade
  - 4 Acid phosphate buffer Dissolve 95 gm of acid sodium phosphate ( $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ ) in 100 ml of distilled water. This is almost a saturated solution
  - 5 Sodium hydroxide a 10% solution in water
  - 6 Distilled water  
Pyrogallol solution Add 15 gm of powdered pyrogallol to 100 ml of 20% sodium hydroxide solution in a rubber stoppered bottle and cover the solution with a layer of oil 1 cm thick. The pyrogallol is dissolved by stirring it with a glass rod passing it through the oil.
- Reagents 2 and 5 are kept for current use in 2 or 5 ml syringes with glass tips attached through the rubber tubing

## Procedure

- 1 Remove the plunger from the syringe, rinse it well with hot water and dry it by suction. This is easily done by attaching a small glass tube to the vacuum line and inserting this tube to the bottom of the barrel so that air is pulled through the warm capillary and barrel.
- 2 Add a drop of glycerol aquaresin mixture into the barrel opening. Lubricate the dry plunger with a few streaks of the glycerol solution and return it to the barrel adding enough of the solution so that no air is carried in with the plunger. Unless the plunger is moved gently in the barrel the latter may break because of the viscosity of the glycerol solution. When the plunger strikes the bottom the glycerol solution should rise only a few mm in the capillary.
- 3 Fill the cup with the ferricyanide solution draw it down into the barrel and expel it without trapping any air bubbles. Remove the excess ferricyanide leaving only the dead space filled with the ferricyanide solution.
- 4 Place a drop of caprylic alcohol on top of the ferricyanide solution draw down a second volume of the alcohol and suck out the remainder. Draw out the plunger until the fluid meniscus is lowered 1 to 2 mm into the capillary.
- 5 Seat the plastic tipped wire plug in the cup trapping a small air bubble on top of the caprylic alcohol. The rubber band must be strong enough to keep the capillary firmly closed by the plastic tip. Add a little of the ferricyanide solution into the cup to complete the seal.
- 6 Pull the plunger out slowly until the spring clip engages the notch in the plunger. Then add glycerol aquaresin solution around the plunger bearing. If the caprylic alcohol stays in the capillary it may be moved down the barrel by warming the capillary behind it in the hand. Allow a few seconds for extraction and release the vacuum by letting the plunger in. Expel the bubble and repeat the extraction.
- 7 Lower the meniscus of the extracted fluid 1 to 2 mm into the capillary. Place the tip of the blood pipette firmly against the bottom of the cup trapping a small (0.5 to 1 division) air bubble between the blood and the extracted fluid. Draw the blood slowly until the length (A) of the air bubble can be measured. Continue to draw the blood in very slowly and evenly until 18 to 20 divisions have entered. Remove the blood pipette, suck off the excess blood, adjust the up-

- per blood meniscus to the zero mark, and read the volume of the blood sample in divisions (With practice the blood volume can be made very close to 30 divisions)
- 8 Fill the cup to the mark with the acid phosphate solution, leaving a small (0.5 to 1 division) air bubble between the buffer and the blood. Draw it slowly in until the length ( $A_1$ ) of the air bubble can be measured. Turn the syringe cup down and continue to draw in the phosphate buffer slowly and evenly until the bubble reaches the junction of the capillary with the syringe barrel, preventing the acid from entering the syringe barrel
  - 9 Rotate the syringe on an oblique axis with the cup down until the blood is completely laked to a uniform dark brown color by the ferricyanide-saponin solution. If more than a trace of acid buffer enters the blood will not lake properly and the oxygen will not be completely released from the hemoglobin.
  - 10 Hold the syringe cup up. Draw in the remainder of the acid phosphate solution and close the capillary with the plastic tipped wire plug trapping a small air bubble (0.5 division or less) underneath. Estimate the length of this bubble ( $A_2$ ). Then add a drop of acid buffer to the cup to complete the vacuum seal.
  - 11 Pull the plunger out slowly so that the capillary drains properly when the small bubble underneath the plug expands. Continue until the plunger notch engages the clip. Again add the glycerol solution around the plunger bearing.
  - 12 Attach the syringe to the sling. It is necessary that the plunger head should rest against the bent end of the sling. Whirl the blood solution down to the bottom of the syringe and tilt back and forth several times to have the liquids well mixed. Whirl the sling again for about 30 seconds. The capillary should be clear of bridges and the blood solution free of bubbles.
  - 13 Remove the syringe from the sling without touching the capillary, and hold the cup up. Release the plunger stop and move the fluid quickly up into the capillary until the atmospheric pressure is almost reached.
  - 14 Without delay remove the plug and adjust the upper gas meniscus to the zero mark and read the volume ( $V_1$ ) of the gas bubble. Since the capillary has not been handled at any time there is no necessity for temperature equilibration in a water bath. (Delay and fumbling in Steps 13 and 14 result in some reabsorption of  $\text{CO}_2$ .)
  - 15 Fill the cup with water and with the syringe held cup up, pull most of the water down into the syringe where it will form a layer on top of the heavy blood solution. Immediately return the bubble to the capillary leaving water behind it on both sides of the gas in the capillary. Owing to the absorption of  $\text{CO}$  by this washing the bubble is now much shorter.
  - 16 Remove the excess water and fill the cup with  $\text{NaOH}$  solution. Invert the syringe run the gas bubble out into the alkali, and rotate a few times to facilitate absorption. Draw some alkali into the capillary and then the gas bubble. Measure its volume ( $V$ ).
  - 17 Remove the  $\text{NaOH}$  solution from the cup and fill it with pyrogallol solution. Run the bubble out into the cup and absorb as in the previous step. Return the bubble to the capillary and read the volume record and designate it  $V_2$  in calculations.
  - 18 Remove the plunger gently from the barrel under the water faucet. Shake out the contents of the barrel and wash with running tap water several times before the plunger is again inserted. It is recommended that dichromate cleaning solution be used after each analysis to insure that no protein remains sticking to the capillary. While the syringe is still hot dry it by suction as described in Step 1.

# Calculations

## Carbon dioxide

Volumes percent of CO in blood =  $(V - V - c) \times f \times 100/b \times i$ .

1  $c$  is the blank correction for the CO content of the reagents (It is determined by running an analysis with all the reagents but no blood. It is usually about 0.1 capillary divisions)

$f$  is the factor from Fig. 5 for adjusting gas volumes from the observed temperature, barometric pressure and water vapor to their dry volumes at 0 and 760 mm.

3  $b$  is the blood volume in terms of division of the capillary.

4  $i$  is the combined correction factor for the incomplete extraction of CO and for the CO reabsorbed by the solutions. It was found to average 1.015 which is the same figure as that determined in the syringe method for CO in blood.

## Oxygen

Volumes percent of O in blood =  $(V - V - c - a) \times f \times 100/b \times i$ .

1  $c$  is the blank correction for the O contents of the reagents used. (Usually about 0 capillary division)

2  $a$  is the O content of the air bubbles introduced during analysis in order to separate reagents.  $a = 21$  per cent  $(A + A + A)$

3  $f$  is the factor for the adjustment of moist gas at observed temperature and barometric pressure to the dry condition at 0 and 760 mm (see Fig. 5)

4  $i$  is the combined correction factor for the incomplete extraction of O and for the reabsorption of O by the solutions. It was found to average 1.011 (The factor was obtained as the average of the factors by which the results with the syringe [in other respects corrected] must be multiplied to equal the results of the Van Slyke analysis)

## 14 Manometric Methods Miscellaneous References

### Amino Nitrogen

Chinard F P and Van Slyke D D. Comparison of a Modified Folin Photometric Procedure and the Ninydrin Manometric Method for the Determination of Amino Acid Nitrogen in Plasma. *J Biol Chem* 169 (3) 515-521 (Aug) 1944

Van Slyke D D. Improved Methods in the Gasometric Determination of Free and Conjugated Amino Acid Nitrogen in the Urine. *J Biol Chem* 161 134 (1951)

Van Slyke D D. Manometric Determination of Primary Amino Nitrogen and Its Application to Blood Analysis. *J Biol Chem* 83 454-457 (May) 1929

Van Slyke D D, MacFadyen D A and Hamilton P B. The Gasometric Determination of Amino Acids in Urine by the Ninydrin-Carbon Dioxide Method. *J Biol Chem* 150 51-53 (Sept) 1943

### Bicarbonate

Van Slyke D D and Cullen G F. Studies of Acidosis. I. The Bicarbonate Concentration of the Blood Plasma: Its Significance and Its Determination as a Measure of Acidosis. *J Biol Chem* 30 39-54 (May) 1917

### Carbon Dioxide

Van Slyke D D, Sendroy J Jr and Lavy S K. Manometric Analysis of Gas Mixtures. III. Manometric Determination of CO Tension and pH of Blood. *J Biol Chem* 95 45-58 (March) 1932

### Carbon Monoxide

Van Slyke D D and Robchest Robbins F S. The Gasometric Determination of Small Amounts of Carbon Monoxide in Blood and Its Application to Blood Volume Studies. *J Biol Chem* 72 39-50 (March) 1927

### Hemoglobin

Van Slyke D D and Hiller A. Gasometric Determination of Hemoglobin by the Carbon Monoxide Capacity Method. *J Biol Chem* 78 807-810 (Aug) 1923

**Hydrogen**

- Van Slyke D D, and Hanke M E *Manometric Analyses of Gas Mixtures IV Manometric Determination of Hydrogen and Oxygen by Combustion*, J Biol Chem 95 569 586 (March) 1932
- Van Slyke D D and Hanke M E *Manometric Analysis of Gas Mixtures V Manometric Analysis of Gas Mixtures Hydrogen Absorption With Paul's Picrate Palladium Solution* J Biol Chem 95 587 593 (March) 1932

**Iodates**

- Van Slyke D D Hiller A and Berthelsen K C *A Gasometric Micromethod for Determination of Iodates and Sulphates and Its Application to the Estimation of Total Base in Blood Serum* J Biol Chem 74 659 676 (Sept) 1927

**Lactic Acid**

- Avery B F and Hastings A B *A Gasometric Method for the Determination of Lactic Acid in Blood* J Biol Chem 94 273 280 (Nov) 1931

**Methemoglobin**

- Van Slyke D D and Hiller A *Gasometric Determination of Methemoglobin* J Biol Chem 84 205 216 (Oct) 1929

**Nonprotein Nitrogen**

- Stehle, R L *Gasometric Determination of Nitrogen and Its Application to the Estimation of the Nonprotein Nitrogen in the Blood* J Biol Chem 45 223 228 (Dec) 1920

**Oxalic Acid**

- Van Slyke D D and Sendroy J Jr *Gasometric Determination of Oxalic Acid and Calcium and Its Application to Serum Analyses* J Biol Chem 84 217 3 (Oct) 1929

**Potassium**

- Kramer B and Gittleman I *The Gasometric Determination of Potassium* Proc Soc Exper Biol & Med. 24 241 243 (Oct June) 1926 1927

**Sugar**

- Van Slyke D D and Hawkins J A *A Gasometric Method for the Determination of Reducing Sugars and Its Application to the Analysis of Blood and Urine* J Biol Chem 79 739 763 (Jan) 1927
- Van Slyke D D and Hawkins J A *Gasometric Determination of Fermentable Sugar in Blood and Urine* J Biol Chem 83 51 70 (May) 1929

**Total Nitrogen**

- Van Slyke D D *Gasometric Micro Kjeldahl Determination of Nitrogen* J Biol Chem 71 235 248 (Dec) 1926
- Van Slyke D D and Kugel V H *Improvements in Manometric Micro Kjeldahl and Blood Urea Methods* J Biol Chem 102 489 498 (Oct) 1933

**Urea**

- Mirkin A *A Gasometric Method for the Determination of Urea Nitrogen in the Blood* J Lab & Clin Med. 8 (1) 50 52 (Oct) 1922
- Van Slyke D D *Determination of Urea by Gasometric Measurement of the Carbon Dioxide Formed by the Action of Urease* J Biol Chem 73 695 704 (June) 1922

## SECTION VI

### PHYSIOLOGICAL MEASUREMENTS (Continued)

---

#### B PHYSICAL CHEMISTRY OF BLOOD GASES

##### 1 Technique for Equilibrating Arterial and Venous Blood Against Gas Mixtures of Known Composition

###### Reference

Henderson, L. J. Blood New Haven Conn 1938, Yale University Press Appendix page 388 by Dill D B

###### Principle

Blood is equilibrated at known tensions of carbon dioxide and oxygen in a constant temperature bath for a known period of time. The blood is then analyzed in the Van Slyke apparatus and the gas phase is analyzed in the Haldane gas analysis apparatus. This technique was used extensively by Henderson, Bock, Dill, Edwards, Talbott and associates in their studies of blood as a physicochemical system.

###### Apparatus

1. Tonometer Barcroft type (see Fig 37)
2. A number of syringes 5 and 10 ml
3. A burette 100 ml (see Fig 37)
4. A number of 50 ml round bottom centrifuge tubes
5. A number of  $\pm$  0 gauge needles (sterilized)
6. A water bath constant temperature set at 37 C. This bath contains a device for rotating blood tonometers
7. A number of blood sampling tubes for the storage of the blood
8. A Haldane gas analysis apparatus as in Fig 41
9. A Van Slyke manometric apparatus as in Fig 8
10. A rubber bag for the storage of gases
11. A number of glass tubes with capillary bore of 2 mm


###### Reagents

1. A tank of oxygen
2. A tank of carbon dioxide
3. A tank of nitrogen
4. Clean mercury
5. Stopcock grease (Lubri-Seal A H Thomas Co)
6. Heparin solution. One drop should inhibit coagulation of 10 ml blood for 24 hours

###### Procedure—Preparation of Gas Mixture

1. The tonometer in which the blood is to be equilibrated is ordinarily prepared in advance of the blood drawing. These tonometers are of the Barcroft type (see Fig 37) and have a capacity between 280 to 300 ml. The capacity is etched on each and serves as an identifying mark. The dry tonometer with the stopcock well lubricated and held in by rubber bands, is fitted with a rubber stopper





## PREPARATION OF GAS MIXTURES

Fig 37—Apparatus used by Henderson and associates for equilibrating blood with gas mixtures of known composition

Since in the calculation of arterial pH we are primarily interested in a partial pressure of 40 mm CO<sub>2</sub> we shall describe specifically the preparation of a mixture with pCO<sub>2</sub> 40 mm Hg and pO<sub>2</sub> 180 to 200 mm Hg. The same general technique is used for preparing tonometers for any other partial pressures.

- The general formula for computing the amount of gas to be added is

ml gas to be added to tonometer =

$$\frac{(p \text{ Gas}) \times (\text{Volume of tonometer})}{(\text{barometric pressure})} \quad (\text{Equation 1})$$

For example pCO<sub>2</sub> wanted is 40 mm pO<sub>2</sub> wanted is 200 mm capacity of tonometer is 300 ml and barometric pressure is 760 mm Hg

$$\text{Therefore ml CO}_2 = \frac{40 \times 300}{760} = 15.8 \text{ ml and}$$

$$\text{ml O}_2 = \frac{200 \times 300}{760} = 79 \text{ and}$$

$$\text{Sum of CO}_2 + \text{O}_2 = (15.8 + 79) = 94.8$$

### 3 Filling the tonometer is done as follows

- The tonometer is first flushed out with nitrogen and then connected to the 100 ml burette
- With the tonometer stopcock open to the outside air and with stopcock of the burette open to the tonometer stopcock the gas burette is filled with mercury by raising the mercury leveling bulb
- The tonometer stopcock is turned and 94.8 ml of gas is removed from the tonometer
- The tonometer stopcock is again turned and the gas is discharged through the side outlet
- The burette stopcock is turned to the rubber bag containing oxygen and the burette connections are flushed out with the oxygen twice
- Then measure out accurately 79.0 ml of oxygen in the burette and deliver it into the tonometer
- The bag of oxygen is replaced with a bag of carbon dioxide and after flushing out the burette and connections through the side outlet 15.8 ml of carbon dioxide are measured into the burette and then introduced into the tonometer
- When more than 15 ml of gas is removed from the tonometer account should be taken of the lowered pressure. This may be done approximately by suitable lowering of the leveling bulb or the volume of gas withdrawn may be measured by closing the burette stopcock.

### Procedure Drawing Preparation and Analysis of Arterial Blood

- For the collection of arterial blood round bottomed centrifuge tubes of 50 ml capacity are used. In the bottom are placed 4 drops of heparin solution and paraffin oil up to a height of 10 cm
- After establishment of local anesthesia with one per cent Novocain blood is withdrawn from the radial or brachial artery
- Particular care is taken in the case of arterial blood to avoid contact with the air. The syringe is turned up the needle replaced with a glass capillary tube 2 x 100 mm and blood is pushed through to the tip of the tube. The syringe is inverted and with the tube leading into the oil in the 50 ml centrifuge tube the blood is rapidly transferred. It is then stirred carefully until the heparin is in solution with the blood. (Any specimen of arterial blood normally represents the blood leaving the lungs in respect to its oxygen and carbon dioxide contents whereas arm vein blood may differ greatly from vena cava blood. Consequently it is of great interest to know the exact carbon dioxide and oxygen content of

arterial blood and of but little interest to secure such information with regard to a sample of venous blood) At this point one man proceeds promptly to the determination of oxygen and carbon dioxide contents of the arterial blood by the Van Slyke techniques

#### Procedure—Equilibration of Arterial Blood

- 1 With a 5 ml syringe 5 ml of arterial blood is introduced into the tonometer through its stopcock
- 2 With another syringe mercury is forced through the stopcock until the capillary bores are exactly filled
- 3 The time is noted and the tonometer is placed in a constant temperature bath set at 37 C where it is held in a mechanical rotator
- 4 Twenty minutes equilibration insures practically complete equilibrium under the above defined conditions
- 5 At the end of 20 minutes, the tonometer is clamped in an upright position in the water bath and a blood sampling tube is attached with a short piece of small bore tubing The sampling tube is like the gas sampling tube except that the capacity is only 6 ml
- 6 The tonometer stopcock is turned to the side arm the sampling tube stopcock is opened and by raising the leveling bulb some mercury is forced through the connection
- 7 The tonometer stopcock is then turned and by lowering of the leveling bulb the blood is drawn into the sampling tube
- 8 After closing of both stopcocks the rubber connection is removed and sampling tube and tonometer are taken from the bath
- 9 The sampling tube is placed in a rack and sealed off with mercury
- 10 The sample is kept at 5 C until it is analyzed for oxygen and carbon dioxide This is done after mixing the blood in the sampling tube thoroughly by inversion. The stopcock is carefully opened until all the mercury is expelled from the capillary bore A number 0 rubber stopper with a small bore makes a satisfactory connection with a Van Slyke Neill pipette The leveling bulb is supported in the rack slightly above the level of the sampling tube stopcock so that there is a positive pressure to fill the pipette when the stopcock is opened After the blood sample has been introduced into the Van Slyke apparatus the sampling tube is again sealed off with mercury

#### Procedure—Analysis of the Gas Phase

- 1 In general the partial pressure of gas in the tonometer at 37 C is calculated from the gas pressure and the temperature of the tonometer at the time of analysis, according to the following equation

pGas at 37 C =

$$\frac{\% \text{ Gas}}{100} \times (B - p\text{H}_2\text{O} + p) \times \frac{(310)}{(273 + T)} \times \frac{(V + H - 1)}{(V - b)} \quad (\text{Equation } )$$

Where % Gas is percentage as analyzed B is barometric pressure mm Hg  $p\text{H}_2\text{O}$  is partial pressure of water at temperature of analysis p is the difference in pressure (tonometer minus barometer) at time of analysis V is volume of tonometer H is the volume of gas drawn into the Haldane apparatus for analysis and b is the volume of blood during equilibration One ml of blood always stays on the walls of the tonometer

- 2 In analyzing the gas phase the steps are as follows
  - a Dry the tonometer and allow it to come to room temperature Expel blood from the exterior ends of the tonometer cock
  - b Connect to the Haldane apparatus
  - c Expel residual gas from the Haldane apparatus

- d. Connect to tonometer and use Haldane leveling bulb to draw gently the residual drop of blood from the interior bore of the tonometer cock and to expel it into the side bore of the exterior cock
- e. Draw 7 ml of gas into the Haldane
- f. Read difference in mercury levels to the nearest mm; this figure  $p$  is used in the calculation
- g. Continue to flush out the dead space and analyze for CO and O. Note that for fully oxygenated blood the  $pO$  is always very high and the value can be assumed without serious error

#### Example of Calculation

1. CO was 5.00 B was 760 mm T at analysis was 5 C  $p$  which is tonometer minus barometer was minus 6 mm Volume of gas drawn into Haldane was 7 ml Tonometer 300 ml blood 7 ml  $pH O$  at 37 C is 74 mm Hg

Therefore  $pCO$  mm Hg =

$$\frac{5}{100} \times (760 - 74 - 6) \times \frac{310}{273 + 37} \times \frac{(300 + 7 - 1)}{(300 - 7)} = 40.0$$

2. The values at different temperatures for  $pH O$  and the factor  $310/(273 + T)$  are consolidated into Table 18

TABLE 18

PARTIAL PRESSURE OF WATER AND TEMPERATURE FACTOR FOR EQUILIBRATED BLOOD AT DIFFERENT ROOM TEMPERATURES

| ROOM TEMPERATURE<br>C | $pH O$<br>mm Hg | FACTOR<br>$310/(273 + T)$<br>f |
|-----------------------|-----------------|--------------------------------|
| 0                     | 17              | 1.058                          |
| 21                    | 19              | 1.054                          |
|                       | 9               | 1.051                          |
| 23                    | 21              | 1.047                          |
| 24                    | 22              | 1.044                          |
| 5                     | 24              | 1.040                          |
| 26                    | 5               | 1.037                          |
| 27                    | 27              | 1.033                          |
| 28                    | 8               | 1.030                          |

## 2 Dissociation Curves of Blood

### References

1. Henderson L. J. Blood New Haven Conn. 1918 Yale University Press
2. Dill, D. B. Graybiel A. Hurtado A. and Taquini A. C. Der Gasaustausch in den Lungen im Alter. Sonderaus Zeits. Alters II (1) 033 1910

### Calculations Involving Oxygen and Hemoglobin

1. Total Oxygen of Blood. Total oxygen of blood (total O)<sub>b</sub> is calculated from results of analysis according to the method of Van Slyke and Neill (1944) as described previously
2. Physically Dissolved Oxygen of Blood (O) According to the solubility studies of Seidory Dillon and Van Slyke (1934) one liter of normal blood will dissolve 0.00138 millimole (mM) of O at 37 C and at 1/760 atmosphere. Hence
 
$$(O)_b \text{ in mM/l} = (0.00138) \times (pO) \quad \text{or} \quad \text{(Equation 1)}$$

$$(O) \text{ in ml } O_2/100 \text{ ml} = (0.0031 \times (pO)) \quad \text{(Equation 2)}$$
3. Oxygen Bound by Hemoglobin This is calculated as follows
 
$$(HbO) = (\text{total } O)_b - (O) \quad \text{(Equation 3)}$$

4 **Oxygen Capacity** This figure ( $\text{HbO}_2$  capacity), is ( $\text{HbO}_2$ )<sub>s</sub> of blood saturated with oxygen at  $p\text{O}_2$  180 to 200 mm Hg at 37 °C

5 **Percentage Oxygen Saturation**

$$\% \text{ Saturation} = \frac{100 \times (\text{HbO}_2)_a}{(\text{HbO}_2 \text{ capacity})_s} \quad (\text{Equation 4})$$

6 **Hemoglobin in Blood.** The hemoglobin is calculated from the oxygen capacity as follows

$$\text{gm Hb/100 ml blood} = \frac{\text{Oxygen capacity ml O}_2/100 \text{ ml blood}}{1.34} \quad (\text{Equation 5})$$

#### Calculations Involving Carbon Dioxide

1 **Total CO of Blood** This value ( $\text{total CO}_2$ )<sub>s</sub> is calculated after Van Slyke and Sendroy (1927)

2 **Glycolysis During Equilibration** Bock and associates (1929) calculated that ( $\text{total CO}_2$ )<sub>s</sub> is reduced about 0.4 vols % in 15 mins at 37 °C. This quantity is added to ( $\text{total CO}_2$ )<sub>a</sub> of equilibrated blood

3 **Carbon Dioxide Content of Blood Equilibrated at 37 °C and  $p\text{CO}_2$  40** Corrected for Glycolysis This value  $T_{40}$  may be found from the curve of Henderson (1930) in Fig. 67 at the back of this book. The data required are ( $\text{HbO}_2$  capacity)<sub>s</sub> and ( $\text{total CO}_2$ )<sub>s</sub> at known  $p\text{CO}_2$ . From these data ( $\text{total CO}_2$ )<sub>a</sub> at any other  $p\text{CO}_2$  may be derived. Corrected for glycolysis this value at  $p\text{CO}_2$  40 is called  $T_{40}$ .

4 **Total Carbon Dioxide in Plasma** If the equilibrated blood is handled under anaerobic conditions direct analysis of the plasma by the technique of Van Slyke and Neill gives the total carbon dioxide ( $\text{total CO}_2$ )<sub>a</sub>. As an alternative, Henderson (1928) developed Fig. 38 for computing  $\Delta(\text{CO}_2)$  difference between ( $\text{total CO}_2$ )<sub>a</sub> and ( $\text{total CO}_2$ )<sub>s</sub> when  $T_{40}$  and ( $\text{HbO}_2$  capacity)<sub>s</sub> are known

$$(\text{Total CO}_2)_a = (\text{total CO}_2)_s + \Delta(\text{CO}_2) \quad (\text{Equation 6})$$

5 **Free Carbonic Acid of Serum ( $\text{HCO}_2$ )** From solubility coefficients of Van Slyke (1928) one liter of normal serum will dissolve 0.0313 mM of  $\text{CO}_2$  at 37 °C and at 1/760 atmosphere. Hence

$$(\text{H}_2\text{CO}_3) \text{ in mM/l} = (0.0313) \times (p\text{CO}_2) \quad (\text{Equation 7}) \text{ or}$$

$$(\text{HCO}_2) \text{ in ml CO}_2/100 \text{ ml} = (0.0696) \times (p\text{CO}_2) \quad (\text{Equation 8})$$

6 **Bicarbonate in Serum ( $\text{BHCO}_2$ )** This value is computed from the equation  $(\text{BHCO}_2) \text{ in mM/l} = (\text{total CO}_2) - (\text{HCO}_2)$  (Equation 9)

7 **Serum pH** This value, pH, is calculated from the Henderson Hasselbalch equation

$$\text{pH} = \text{pK} + \log (\text{BHCO}_2) - \log (\text{HCO}_2) \quad (\text{Equation 10})$$

Dall, Daly and Forbes (1937) determined that the pK of normal human serum at 37 °C is 6.11

#### Calculation of Arterial $p\text{CO}_2$

1 **Knowing  $\text{CO}_2$  combining capacity ( $T_{40}$ ) and the  $\text{CO}_2$  content of arterial blood as drawn and the oxygen combining capacity one calculates  $p\text{CO}_2$  of the blood as drawn from Fig. 67 at the back of the book as follows**

a Place string at  $T_{40}$  on the left and 40 mm on P scale at right


b Fix a pin on A scale corresponding to the oxygen capacity in ml/100 ml blood

c Now swing the string about the pin setting the the left hand figure on the T scale at the  $\text{CO}_2$  content as analyzed.

d Keeping string straight read  $p\text{CO}_2$  on right hand scale

2 **Under ordinary conditions the  $\text{CO}_2$  combining capacity of arterial blood may be taken as 0.1 mM greater than that of oxygenated blood. Under unusual conditions the correction must be arrived at more precisely. After correction**

# HEMOGLOBIN



—C. BORDEN—

Fig 38—Hand some chart for computing whole blood and serum anion difference for blood with different concentration of hemoglobin.

is made for dissolved oxygen the proportion of reduced Hb is calculated. The  $\text{CO}_2$  combining capacity of arterial blood then is

$\text{CO combining capacity} \approx$

$$T_{50} + \frac{\Delta (\text{CO})_b \times (\text{Hb})_b}{(\text{HbO}_2 \text{ capacity})_b} \quad (\text{Equation 11})$$

Where  $T_{50}$ ,  $(\text{Hb})_b$ , and  $(\text{HbO}_2 \text{ capacity})_b$  are as defined above and  $\Delta (\text{CO})_b$  is the reduction of  $\text{CO}_2$  combining capacity for increasing concentrations of oxyhemoglobin. To determine the value of  $\Delta (\text{CO}_2)_b$  from the  $T_{50}$  of oxygenated blood and  $(\text{HbO}_2 \text{ capacity})_b$  use Fig. 39.





Fig. 39—Henderson's chart for determining the effect of oxygenation on the  $T_{50}$  of whole blood.

#### Dissociation Curves

- 1 For both oxygen and carbon dioxide, dissociation curves are plotted with some function of gas content as ordinate against some function of gas pressure as abscissa.



## OXYGEN DISSOCIATION CURVES, HUMAN BLOOD, AT DIFFERENT TEMPERATURES

Fig 40—Effects of temperature and pH on oxygen dissociation curves (according to Dill and associates)



- 2 Points on the oxygen dissociation curves are obtained by preparing gas mixtures of suitable composition followed by precise determination of both  $p\text{CO}_2$  and  $p\text{O}_2$  in the gas phase. The oxygen dissociation curve is usually plotted to show %  $\text{HbO}_2$  as a function of  $p\text{O}_2$ , the  $p\text{CO}_2$  being constant. An example of this kind of curve is shown in Fig. 40.
- 3 The logarithmic oxygen dissociation curve plots  $\log \frac{100(\text{Hb})_s}{(\text{HbO})_s}$  against  $\log p\text{O}_2$ , where  $(\text{Hb})_s/(\text{HbO})_s$  is the ratio of reduced to oxygenated hemoglobin. This curve has two advantages over the usual oxygen curve. It is nearly straight and the effect of a given change in pH (or in temperature) on the  $\log p\text{O}_2$  is the same for all values of  $100(\text{Hb})/(\text{HbO})$ .

### 3 Miscellaneous Calculations of the Physicochemical Properties of Blood

#### Calculation of the Donnan Ratio

##### Reference

Henderson L. J. Blood New Haven Conn. 1928 Yale University Press

##### Calculation

- 1 As applied to blood the Donnan ratio is the ratio of the concentration of an ion in the cells divided by its concentration in the plasma, concentration being best expressed as mEq of electrolyte per liter of water in the cells or serum.
- 2 The direct determination of the Donnan ratio is accomplished in the following steps:
  - a. Separate cells and plasma by centrifuging
  - b. Determine the concentration of the ion in the cells and plasma
  - c. Determine the water concentration in the cells and plasma by drying weighed volumes at  $110^\circ\text{C}$  (see estimations of water and solids)
  - d. Calculate the Donnan ratio according to the equation

$$r \text{ the Donnan ratio} = \frac{\text{mEq/l cells}}{\text{mEq/l plasma}} \times \frac{\text{ml water/1000 ml plasma}}{\text{ml water/1000 ml cells}}$$

- 3 The indirect determination of the Donnan ratio is accomplished with measurements on whole blood and plasma:
  - a. Determine the hematocrit i.e. ml packed red cells/100 ml blood.
  - b. Determine the concentration of the ion in whole blood and in plasma
  - c. Determine the water content of whole blood and plasma by drying weighed known volumes at  $110^\circ\text{C}$  (see estimation of water and solids)
  - d. Estimate the concentration of ion in plasma mEq/l
  - e. Estimate the concentration in whole blood mEq/l
  - f. Compute the concentration in cells according to the equation

$$\text{mEq/l cells} =$$

$$\frac{\text{mEq/l whole blood}}{\text{Hematocrit/100}} - \frac{\text{mEq/l plasma} \times (100 - \text{Hematocrit})}{\text{Hematocrit}}$$

- g. Compute the Donnan ratio

$$r \text{ the Donnan ratio} = \frac{\text{mEq/l cells}}{\text{mEq/l plasma}} \times \frac{\text{ml water/1000 ml plasma}}{\text{ml water/1000 ml cells}}$$

#### Composition of Normal Arterial Human Blood

For some purposes a table of some of the properties of normal human arterial blood is useful. Table 19 has been compiled from data accumulated at the Harvard Fatigue Laboratory.

## Temperature Effects on Base Bound by Plasma Protein

## Reference

- 1 Stadie W C Austin J H and Robinson H W The Effect of Temperature on the Acid Base Protein Equilibrium and Its Influence on the CO Absorption Curve of Whole Blood True and Separated Serum *J Biol Chem* 66 801-90 (Dec) 1955
- 2 Van Slyke D D Hastings A B Hiller A and Sendroy J Jr Studies of Gas and Electrolyte Equilibria in Blood XIV The Amounts of Alkali Bound by Serum Albumin and Globulin *J Biol Chem* 79 69-80 (Oct) 1928

## Calculation

The base bound by plasma protein, expressed as BP in mEq/l is a function of

- a Concentration of plasma protein P expressed in gm/l
- b pH
- c pI the pH of minimal base binding
- d Buffer value expressed by a constant 0.104 in the case of normal man

The equation at 33 C according to Van Slyke Hastings Hiller and Sendroy is

$$BP = (0.104) \times (P) \times (pH - 5.08)$$

At this temperature pI has the value of 5.03

The effect of temperature on pI according to Stadie Austin and Robinson is about -0.0 per 1 C that is

$$\Delta pI = -0.0 \text{ (Observed } T - 38 \text{ C)}$$

Hence the general equation for evaluating BP taking temperature into account is

$$BP = (0.104) \times (P) \times pH - (5.08 - \Delta pI)$$

By substitution this becomes

$$BP = [(0.104) \times (P)] - [0.08 - 0.0(t - 38 \text{ C})]$$

where t = observed temperature

TABLE 19

## COMPOSITION OF NORMAL ARTERIAL HUMAN BLOOD

Averages based on over 100 analyses at the Harvard Fatigue Laboratory with healthy young men as subjects)

| MEASUREMENT                | CELLS | PLASMA | WHOLE BLOOD |
|----------------------------|-------|--------|-------------|
| pH                         | 7.39  | 7.40   |             |
| Volume ml/l                | 448   | 55     | 1000        |
| Sp Gravity                 | 1.100 | 1.027  | 1.060       |
| Water ml/l                 | 720   | 938    | 840         |
| Total Solids gm/l          | 390   | 89     | 20          |
| Protein gm/l               | 356   | 70     | 130         |
| Glucose gm/l               | 0.9   | 1.1    | 1.0         |
| Sodium mEq/l               | 149   | 140    | 83.7        |
| Potassium mEq/l            | 93.7  | 5.3    | 44.9        |
| Calcium mEq/l              | 0.0   | 5.0    | 2.6         |
| Magnesium mEq/l            | 0     | 0.0    | 0.0         |
| Total cations mEq/l        | 110.6 | 15.3   | 133.2       |
| Chloride mEq/l             | 54    | 104.6  | 82.0        |
| Bicarbonate mEq/l          | 15.3  | 5.1    | 0.7         |
| Lactate mEq/l              | 0.7   | 1.4    | 1.1         |
| Inorganic Phosphorus mEq/l | 1.1   | 0      | 1.7         |
| Inorganic Sulfate mEq/l    | 0.6   | 0.6    | 0.6         |
| Protein mEq/l              | 41.7  | 17.4   | 78.2        |
| Total anions mEq/l         | 113.3 | 151.3  | 133.9       |

## SECTION VI

### PHYSIOLOGICAL MEASUREMENTS (Continued)

---

#### C ESTIMATION OF RESPIRATORY GASES

##### 1 The Haldane Apparatus for Measuring Carbon Dioxide and Oxygen

###### Reference

Haldane, J S and Priestley, J G *Respiration*, Oxford 1935, Clarendon Press

###### Principle

The  $\text{CO}_2$  in a known volume of air is absorbed with  $\text{KOH}$ , and the  $\text{O}_2$  with a strong reducing solution. The remainder of gas is considered to be all  $\text{N}_2$ .

###### Apparatus


- 1 Haldane gas analyzer (See Fig 41)
- 2 Gas sampling tubes. Oiled syringes may be used if storage time is only a few hours
- 3 Pipe cleaners for cleaning apparatus
- 4 50 ml syringes for handling reagents and cleaning apparatus
- 5 Dropping bottle for dilute  $\text{H}_2\text{SO}_4$
- 6 Small lengths of pressure tubing for cleaning apparatus

###### Reagents

- 1 Dilute  $\text{H}_2\text{SO}_4$  approximately 1 ml conc  $\text{H}_2\text{SO}_4$  diluted to 100 ml with water
- 2  $\text{KOH}$  approximately 14 gm dissolved in 100 ml water
- 3 Litmus paper for detecting acid or alkali in the wrong places
- 4 Oxygen absorbent. In a screw topped bottle place about 8 gm of sodium hydro sulfite ( $\text{Lycopon Na}_2\text{SO}_3$ ) 15 gm sodium anthraquinone beta sulfonate and a small crystal of ferric chloride. Just before the apparatus is to be made ready for use add approximately 50 ml of 14%  $\text{KOH}$  which has been warmed to about 50 C. Solution is complete in a few minutes. Throughout all manipulation of the absorbent unnecessary exposure to air must be avoided. It is wise to weigh out about 6 bottles of dry reagents at a time since 50 ml of absorbent will be enough for about 50 analyses.
- 5 Stopcock grease mercury and copper sulfate crystals

###### Procedure—Calibrating Burette

- 1 This is necessary before a new burette is used since the bore is not necessarily uniform. It need not be repeated unless structural alterations are made in the burette.
- 2 Fuse a one way stopcock with a fine tip on the lower end of the burette. (Most companies ship burettes with stopcock attached.)
- 3 Clean the burette with sulfuric acid dichromate mixture



## HALDANE APPARATUS

G. E. GORDON

Fig. 41.—Haldane's apparatus for measuring oxygen, carbon dioxide, and nitrogen in samples of air.

- 4 Rinse thoroughly with water and finally with dilute sulfuric acid and allow to drain for 10 minutes
- 5 With the short limb of Stopcock 1 to the left and the long limb open to air fill the burette with mercury exactly to the top of the stopcock bore Hang a thermometer right next to the burette
- 6 Into a weighed 25 ml Erlenmeyer flask deliver mercury through the fused bottom stopcock until the mark 7 00 is exactly reached by the upper part of the mercury meniscus Measure the temperature
- 7 Reweigh the flask
- 8 Now deliver mercury from the bottom stopcock into the flask until exactly 7 100 is reached measure the temperature and weigh again
- 9 Repeat this for 7 200 7 300 7 400 and so on to 10 000
- 10 The calibration is now calculated by the following steps
  - a List the burette readings (Column 1)
  - b List the corresponding temperatures and weights of Hg (Columns 2 and 3)
  - c From 7 0 on calculate the increment in weight for each 0 1 ml division on the burette (Column 4)
  - d List the density of mercury i e, gm of Hg per ml at given temperatures (Column 5) from table 20
  - e Calculate the increment in volume for every 0 1 increase in burette reading by dividing every number in Column 4 by the corresponding number in Column 5
  - f Calculate the actual volume of the burette at every 0 1 division by determining the actual volume at 7 0 (i e weight of Hg at 7 0 divided by density of Hg at the temperature of the burette when 7 0 was weighed) and adding to it successively each number in Column 6
  - g The final answer in a Haldane analysis is expressed as ml of gas per 100 ml of initial sample Hence it is now necessary to calculate all burette readings on the basis that burette reading 10 0 represents exactly 10 00 ml. This is done by multiplying all readings in Column 7 by the factor

$$\frac{10\ 000}{\text{Column 7 reading for } 10\ 000}$$

- h. The correction factor for each burette reading is now calculated (Column 9) This correction represents the number which has to be added to any given burette reading in order to give its value correctly on the basis of 10 000 = exactly 10 ml
- i The table of correction factors is placed on a card which is mounted for use on the Haldane apparatus which has been calibrated
- j Steps 10a through 10h are illustrated in Table 21 which presents the actual calibration of a new burette

TABLE 20  
DENSITY OF MERCURY AT VARIOUS TEMPERATURES

| C                                      | 0 0     | 0 1     | 0 2     | 0 3     | 0 4     | 0 5     |
|--|---------|---------|---------|---------|---------|---------|
| GM Hg PER ml AT DIFFERENT TEMPERATURES |         |         |         |         |         |         |
| 20                                     | 13 5462 | 13 5457 | 13 5452 | 13 5448 | 13 5443 | 13 5438 |
| 21                                     | 5438    | 5433    | 5428    | 5423    | 5418    | 5413    |
| 22                                     | 5413    | 5408    | 5403    | 5399    | 5394    | 5389    |
| 23                                     | 5389    | 5384    | 5379    | 5374    | 5369    | 5364    |
| 24                                     | 5364    | 5359    | 5354    | 5350    | 5345    | 5340    |
| 25                                     | 13 5340 | 13 5335 | 13 5330 | 13 5325 | 13 5320 | 13 5315 |
| 26                                     | 5315    | 5310    | 5305    | 5301    | 5296    | 5291    |
| 27                                     | 5 91    | 5 86    | 5 81    | 5 76    | 5 71    | 5 66    |
| 28                                     | 5 66    | 5 61    | 5 56    | 5 52    | 5 47    | 5 42    |
| 29                                     | 5 42    | 5 37    | 5 32    | 5 27    | 5 22    | 5 17    |

TABLE 21  
SAMPLE CALIBRATION OF A HALDANE BURETTE

| 1                           | 2                       | 3                     | 4                                | 5                         | 6                                | 7                      | 8                                     | 9                                 |
|-----------------------------|-------------------------|-----------------------|----------------------------------|---------------------------|----------------------------------|------------------------|---------------------------------------|-----------------------------------|
| BURETTE<br>READING<br>UNITS | TEMP OF<br>BURETTE<br>C | WEIGHT<br>OF Hg<br>gm | INCRE<br>MENT IN<br>WEIGHT<br>gm | DENSITY<br>OF Hg<br>gm/ml | INCRE<br>MENT IN<br>VOLUME<br>ml | ACTUAL<br>VOLUME<br>ml | VOLUME<br>ON BASIS<br>10.0 =<br>10 ml | CORREC<br>TION<br>FACTOR<br>UNITS |
| 70                          | 45                      | 98.0339               |                                  | 13.535                    |                                  | 7.43                   | 7.004                                 | +4                                |
| 71                          | 35                      | 99.4474               | 1.4135                           | 13.5376                   | 0.1044                           | 7.347                  | 7.105                                 | +5                                |
| 80                          | 46                      | 111.9599              | 1.3815                           | 13.5350                   | 0.100                            | 8.7                    | 8.000                                 | 0                                 |
| 81                          | 24.6                    | 113.364               | 1.4035                           | 13.5350                   | 0.1037                           | 8.375                  | 8.099                                 | -1                                |
| 89                          | 47                      | 124.5315              | 1.3935                           | 13.5347                   | 0.109                            | 9.00                   | 8.897                                 | -3                                |
| 90                          | 48                      | 15.9389               | 1.4074                           | 13.5345                   | 0.1040                           | 9.304                  | 8.998                                 | -2                                |
| 99                          | 25.1                    | 138.5453              | 1.3808                           | 13.5337                   | 0.100                            | 10.236                 | 9.900                                 | 0                                 |
| 100                         | 51                      | 139.9560              | 1.4112                           | 13.5337                   | 0.1043                           | 10.340                 | 10.000                                | 0                                 |

The correction factor is the number to be added in the third decimal place to the actual burette reading in order to make it represent the actual volume on the basis that the burette reading 10.0 is exactly 10 ml. Example: Burette reading 8.957 Correction Factor from Column 9 = 2 Corrected reading 8.955

#### Procedure—Setting Up and Preparing Apparatus for Use

- 1 Install the burette and dummy burette in the water jacket put a few crystals of  $\text{CuSO}_4$  in the bottom and fill with water to  $\frac{1}{2}$  inch of the top. The tube of an air bubble either hand driven or from a pressure pump should be installed at the bottom of the jacket. The dummy burette should contain 1%  $\text{H}_2\text{SO}_4$  up to about 8 ml.  
Remove all stopcocks and clean out their tubes and cores.
- 3 With stopcocks still out fill the mercury reservoir with mercury so that the hand bulb is about  $\frac{1}{3}$  full when the levels of mercury in the bulb and burette are the same and the mercury in the burette is almost up to the Stopcock 1. Suck 3 drops of 1%  $\text{H}_2\text{SO}_4$  in on top of the mercury column. Now hang the hand bulb on the right hand rack and pump low down so that the mercury is at about 10 ml in the burette.
- 4 Grease and insert Stopcock 1 leaving it in the full off position (not connected with any part of the machine or outside air).
- 5 With Stopcocks 2, 3 and 4 taken out and the KOH bulb in the position shown add 14% KOH to the bulb until the KOH level is at A and B and the bulb is about  $\frac{1}{4}$  full.
- 6 Grease and insert Stopcocks 2 and 3 leaving them in the full off position.
- 7 With Stopcock 4 taken out add 14% KOH to the far right hand glass tube until its left hand bulb is full and the KOH is almost to the bottom of its right hand bulb. This is necessary to protect oxygen absorbent from atmospheric oxygen.
- 8 Fill a 50 ml syringe with oxygen absorbent up to about the 40 ml mark. Unscrew the pinch clamp on the rubber nipple at the bottom of the middle U tube and insert the nozzle of the syringe. Force oxygen absorbent slowly up into the U tube until it is about 2 cm above C. Air will now be trapped in the right hand bulb of this U tube and oxygen will be absorbed. During this period lasting about 1 minute keep the level of absorbent at about 2 cm above C by slowly pushing it in as necessary from the syringe. When the rapid absorption of oxygen has ceased clamp off the rubber nipple at the bottom and remove the syringe. The absorbent should now be approximately 2 cm above level C. Oxygen absorption is not complete for  $\frac{1}{4}$  hour.
- 9 Grease and insert Stopcock 4 leaving it in the full off position.
- 10 Before an analysis can be performed the whole apparatus from Stopcock 1 to levels A and C must be full of nitrogen. This is accomplished by carrying out one analysis as described below but without any measurements.

**Procedure—Taking a Sample of Gas**

- 1 The gas sampling tube is fixed near the upper left hand corner of the machine and is connected as shown to the burette by means of rubber tubing through a 3 way stopcock of small bore. For field use, syringes are best and the following description applies only to them.
- 2 Stopcocks 1 and 5 are opened so as to connect the burette with outside air but not with the rest of the apparatus or the gas sampling syringe.
- 3 With the right hand mercury leveling bulb is detached from the rack and is carefully raised until the mercury in the burette almost reaches Stopcock 1. It must never enter Stopcock 1.
- 4 Stopcock 5 is turned with the left hand so as to connect the sampling tube with the burette but not outside air and the mercury bulb is lowered 5 or 10 cm so as to exert a negative pressure on the sampling tube and prevent contamination of the sample in Step 5.
- 5 The pinch clamp on the nozzle of the syringe is loosened with the left hand and with the right the mercury bulb is lowered until the burette reads about 10. Stopcock 5 is turned so as to connect the burette with outside air but not the sampling tube. (Keep the left hand on Stopcock 5 for smooth and accurate sampling. No pressure is required on the barrel of the syringe unless the sampling tube has been poorly oiled.)
- 6 Smoothly raise the mercury bulb. When the column is approaching Stopcock 1 turn Stopcock 5 to connect the burette with the sampling tube but not with outside air and simultaneously lower the mercury bulb until the column stands at about 10.  
Repeat Step 6 once more and hang the mercury bulb on the rack and pinion so that the mercury level is 10.0 or slightly less.
- 8 Turn Stopcock 1 to the full off position and close the pinch clamp on the sampling tube.

**Procedure—Analysis of Sample**

- 1 For analysis the machine must have nitrogen in it from Stopcock 1 to A and C and the fluid levels should be at A, B and C. Both of these requirements will have been fulfilled by a previous analysis.
- 2 Turn Stopcock 3 to be open to dummy burette, KOH and outside air.
- 3 Turn Stopcock 2 to be open to burette and KOH but not outside air.
- 4 Turn Stopcock 1 to be open to burette and KOH but not outside air.
- 5 Rapidly adjust the KOH bulb up and down until the fluid level is at B and turn Stopcock 3 to connect with dummy burette and KOH but not outside air.
- 6 Adjust fluid level to be exactly at A by moving the rack and pinion. Bubble air through the water jacket. (The purpose of the dummy burette is to compensate for possible changes in temperature or barometric pressure during the course of an analysis. This can be seen from Fig. 41. If the fluid levels are at A and B and the samples of gas in the burette and the dummy burette are at external barometric pressure as they will be if all the steps above are followed carefully any change of temperature in the apparatus or any change of pressure will affect burette and dummy burette equally. The gas in the burette is thus put in a condition which makes any change in its volume due solely to the removal of  $\text{CO}_2$  or  $\text{O}_2$ .)
- 7 Read the total volume,  $V$ , at the top of the mercury meniscus. Check reading after readjusting levels A and B. Correct  $V$  from the correction chart.
- 8 In order to absorb  $\text{CO}_2$  raise the mercury bulb until the mercury is in the enlarged portions of the burette. Raise and lower the bulb 6-8 times keeping the mercury oscillating between about 2 and 6 ml.

- 9 Carefully hold the mercury bulb next to the rack and pinion with the level at about 6 ml and slowly lower it until the KOH is about 1 cm below level A. Hang the bulb on the rack and pin on.
- 10 Adjust levels A and B, bubble air through the water jacket and read the burette V at the top of the mercury meniscus.
- 11 Repeat Steps 8 and 9 with 4 flushes into KOH.
- 12 Repeat Step 10.
- 13 If the second reading of V is not within 0.001 of the first, repeat Steps 8, 9 and 10 until checks are obtained. (If V becomes progressively smaller, there is probably a leak in Stopcock 1 or 2.) Correct V from the correction chart.
- 14 To absorb O<sub>2</sub>, turn Stopcock 2 to connect the burette with absorber but not with KOH.
- 15 Using the same technique as in Step 8, rinse the gas into the absorbent about 8 times.
- 16 With the right hand holding the mercury bulb, bring the absorber carefully up to approximately C and hold it there while with the left hand Stopcock 2 is turned to connect the burette with KOH but not absorbent.
- 17 In order to remove O<sub>2</sub> from the gas between Stopcock 2 and A, rinse gas into KOH 3 times.
- 18 Bring fluid level to approximately A and turn Stopcock 2 to connect burette with absorbent but not KOH.
- 19 Repeat Steps 15, 16 and 17 once.
- 20 Rinse into absorbent finally about 3 times, hang the bulb on the rack and pin on and bring the level exactly to C.
- 1 Turn Stopcock 2 to connect burette with KOH but not absorbent and adjust KOH levels exactly to A and B and bubble air through the water jacket. Read the burette V at the top of the mercury meniscus.
- 2 Repeat Steps 14, 15, 16, 17, 18 and 21. V must check to at least 0.001 ml. (If it drops progressively, there are probably leaks in the stopcocks or else the absorbent is inefficient.) Correct V from the correction chart.
- 23 Prepare for the next analysis by opening Stopcock 3 to dummy burette and KOH but not outside air, readjusting KOH levels to exactly A and B, closing Stopcock 3 to connect with dummy burette and KOH but not outside air, readjusting levels A and B and opening Stopcock 1 to connect with outside air but not to the machine.
4. When leaving the apparatus after a series of analyses, turn Stopcocks 1, 2 and 4 to full off position and turn Stopcock 3 so that it connects dummy burette outside air and KOH. The above precautions are necessary to prevent changes in temperature or barometric pressure sucking over the solutions into the wrong parts of the apparatus.

#### Procedure—Cleaning the Apparatus

- 1 When KOH or oxygen absorbent has been sucked into Stopcocks 4, 1 or some times even further, cleaning with water and finally dilute sulfuric acid is necessary. If there is a suspicion of such a misfortune, remove the stopcock and insert a piece of red litmus paper. If it turns blue, cleaning is necessary because alkali is present.  
Open taps 1 and 4 to outside air but not to burette, KOH or absorbent and disconnect Stopcock 5.
- 3 Attach a piece of rubber tubing to the glass tubing to the right of Stopcock 4.
- 4 Using a large syringe, force 100 ml of water through the whole line of glass tubing into a beaker held under the outside air connection of Stopcock 1.
- 5 Follow this with 50 ml of dilute sulfuric acid and remove excess fluid by pumping air through with the empty syringe.



- 6 If it is suspected that alkali has reached the burette remove Stopcock 1, run the mercury right up to the stopcock seat and test with litmus paper. If there is alkali, replace Stopcock 1 turn it to connect burette with outside air but not machine and rinse twice by sucking dilute sulfuric acid all the way beyond 10 and discarding it through the stopcock. Leave 2-3 drops of dilute sulfuric acid in the burette.
- 7 Remove all stopcocks wash them in hot water clean them with a pipe cleaner dry them, grease them and put them back in the apparatus.
- 8 Run an analysis without measurements in order to fill the apparatus with nitrogen.

**Calculations**

$$\text{ml CO}_2 \text{ per 100 ml sample} = \frac{V_{\text{corrected}} - V_1 \text{ corrected}}{V_{\text{corrected}}} \times 100$$

$$\text{ml O}_2 \text{ per 100 ml sample} = \frac{V_1 \text{ corrected} - V_{\text{corrected}}}{V_{\text{corrected}}} \times 100$$

$$\text{ml N}_2 \text{ per 100 ml sample} = \frac{V_{\text{corrected}}}{V_{\text{corrected}}} \times 100$$

**Example**

|   | READING        | BURETTE<br>CORRECTION | CORRECTED<br>READING |
|---|----------------|-----------------------|----------------------|
| V | 9.897<br>9.897 | 0                     | 9.897                |
| V | 9.643<br>9.640 | +2                    | 9.644                |
| V | 7.901<br>7.900 | +5                    | 7.905                |

$$\text{ml CO}_2 \text{ per 100 ml sample} = \frac{9.897 - 9.644}{9.897} \times 100 = \frac{.253}{9.897} = 2.56$$

$$\text{ml O}_2 \text{ per 100 ml sample} = \frac{9.644 - 7.905}{9.897} \times 100 = \frac{1.739}{9.897} = 1.76$$

$$\text{ml N}_2 \text{ per 100 ml sample} = \frac{7.905}{9.897} \times 100 = 79.87$$

$$\text{CO}_2 + \text{O}_2 + \text{N}_2 \text{ should} = 100$$
**Precautions—Reagents**

1. One filling of oxygen absorbent suffices for about 50 estimations. When it is becoming weak scum forms at level C and V checks are hard to obtain. It is then time to change the reagent.
2. One filling of KOH will last for about 1000 estimations.

**Precautions—Calibrating Burette**

1. There are no short cuts in accurate calibration. Each burette should be calibrated twice or until good agreement is obtained.

**Precautions—Setting up Apparatus**

1. A chief source of systematic errors is improper lubrication of stopcocks. When properly seated they show no air bubbles in the grease layer between the core and the shell of the stopcock and the capillaries are free of grease.

**Precautions—Taking Gas Samples**

1. It is essential to have the gas in the burette at external barometric pressure before connecting burette with KOH through Stopcocks 1 and 2. Otherwise KOH will inevitably be sucked over or else will be pushed far enough down to absorb

CO before  $V$  can be read. A properly oiled syringe will avoid this but if there is any binding during the taking of samples set the mercury at about 10 for the final step leave Stopcock 1 as it is and open Stopcock 5 to burette and outside air before closing Stopcock 1.

#### Precautions—Analysis of Sample

- 1 The commonest error is sucking KOH oxygen absorbent or both through Stopcocks 2 and 4. This can be avoided by
  - a Avoiding jerky abrupt movements
  - b Slowing down whenever mercury KOH or absorbent approach or leave constrictions in the apparatus
  - c Never moving the mercury bulb unless one is sure that the stopcocks are turned in the correct manner
  - d During absorption of CO or O making sure that the mercury in the burette does not fall below 70
 All readings must check to within 0.001 ml. When successive readings fall progressively
  - a In the case of  $V_0$ , there is either a leak in one or more stopcocks or else CO is being absorbed from the sample by KOH somewhere in the system
  - b In the case of  $V$  there is either a leak in one or more stopcocks or else oxygen absorbent has reached to or beyond Stopcock 2
  - c In the case of  $V$  there is either a leak in one or more stopcocks or else the oxygen absorbent has become inefficient
- 3 Erratic readings usually come (a) from errors in technique (b) if the burette does not contain acidulated water (c) if grease or mercury act as a ball valve somewhere in the apparatus
- 4 Be sure to follow Steps 23 and 24 carefully in order to avoid trouble in subsequent analyses

#### Precautions—Cleaning the Apparatus

- 1 Avoid using alcohol ether or other organic solvents. They are hard to remove completely from the system are in part dissolved in the reagents have high vapor pressure and give very erratic results.

#### Precautions—Calculations

- 1 If the volume of dilute sulfuric acid in the burette amounts to more than 0.005 ml it must be taken into account in the calculations since it takes up volume that would otherwise be occupied by gas. It is wise to adopt the following routine
  - a After the apparatus has stood undisturbed for at least an hour read the top of the mercury meniscus and the bottom of the dilute sulfuric acid meniscus
  - b  $V_0$ ,  $V$  and  $V$  are read as usual
  - c In all calculations the differences between  $V_0$ ,  $V$  and  $V$  are calculated as usual. However instead of dividing by  $V$  one divides by ( $V$  - amount of dilute sulfuric acid). If dilute sulfuric acid is less than 0.005 it is neglected.
- 2 Before proceeding to a series of analyses it is wise to analyze a sample of pure outdoor air. The percentage of CO should be  $0.03 \pm 0.01$  and the percentage of O should be  $20.93 \pm 0.0$ . If one cannot achieve this accuracy there is no use proceeding to analysis of unknown samples. More practice or else cleaning of the apparatus will be necessary.
- 3 In satisfactory duplicate samples of expired air analysis for CO must agree to within 0.03 percent and analyses for O must agree to within 0.04 percent. Until the operator can achieve these limits he is unreliable.

## 2 The Scholander Micrometer Gas Analyzer

*Reference*

Scholander, P. F. Analyzer for Accurate Estimation of Respiratory Gases in One Half Cubic Centimeter Samples *J Biol Chem.* 167 235-250 (Jan.) 1947

*Principle*

The present analyzer permits the determination of carbon dioxide oxygen and nitrogen in 0.5 ml samples with an accuracy of  $\pm 0.015$  volume percent. It will handle directly samples containing from zero to over 99 percent absorbable gases. The analysis requires 6 to 8 minutes. A gas sample is introduced into a reaction chamber, connected to a micrometer burette and is balanced by means of an indicator drop in a capillary against a compensating chamber. Absorbing fluids for carbon dioxide and oxygen can be tilted into the reaction chamber without causing any change in the total liquid content of the system. During absorption of gas mercury is delivered into the reaction chamber from the micrometer burette so as to maintain the balance of the gas against the compensating chamber. Volumes are read in terms of micrometer divisions. The rising fluid and absorbents are accurately adjusted to have the same vapor tension.

*Apparatus (See Fig. 42)*

1. **Reaction Chamber Unit** This part of the apparatus consists of a compensating chamber (A), a reaction chamber (B) and two side chambers (C and D) for storing absorbents for carbon dioxide and oxygen respectively.

The compensating chamber or thermobarometer (A) can be closed with a stop cock (S1), the lower end of which is ground to fit into the neck of A. Between analyses S1 is kept in a ground glass receptacle (G). A 1 mm bore capillary, with the upper end ground to a wide angled funnel, intrudes into the bottom of the compensating chamber. The capillary is ring marked in the middle. The lower end of the capillary bore is coarsely surface ground for 2 mm and is covered with a Clarite coating. The upper part of the compensating chamber is also made hydrophobic with a Clarite coating. The somewhat unconventional scheme of running the indicator drop vertically is correct as long as the weight of the drop remains constant during the analysis. The intruded capillary assures that no acid can trickle down and add to the weight of the drop. The lower end of the reaction chamber unit is ground to fit the micrometer burette unit. The side arms C and D of the reaction chamber are furnished with capillaries for filling. The ends are stoppered with tight fitting solid vaccine bottle rubber stoppers (E and F).

2. **Micrometer Burette Unit** A micrometer burette H is attached to a 3 way stop cock (S2) the upper limb of which is ground to fit the lower end of the reaction chamber unit. The side limb is connected with a mercury leveling bulb (I) through rubber or plastic tubing. A 5 mm thick round wooden dowel 10 cm long is fitted to the knurled ratchet on the handle of the micrometer by means of a 5 cm length of rubber tubing.

3. **Water Bath** The reaction chamber unit is immersed in a water bath which can be made of a glass bottle with the bottom cut off. A one hole rubber stopper is inserted in the neck of the bottle and the upper limb of the micrometer burette is inserted through the hole. The rubber stopper should be short enough to permit rocking the apparatus by tapping the micrometer burette unit. The water bath rests on a ring which can be tilted on a friction bushing. The position of the ring is controlled by means of an attached handle (J). The water level reaches to the neck of (A) so that the gas in the thermobarometer is totally immersed. The water bath is continuously stirred by bubbling air from the bottom. A 50 ml centrifuge tube (K) is suspended from the edge into the water bath.

## SCHOLANDER GAS ANALYZER




Fig. 1—Analyze for accurate estimation of respiratory gases. 0.5 ml sample. A, Gas inlet chamber (in room-barometer). B, reaction chamber. C, side arm for carbon dioxide absorbent. D, side arm for oxygen absorbent. E and F, solid vacuum bottle stoppers. G, receptacles for top-cock. S-1, H, micrometer bulb. J, leveling bulb. K, handle for tilting of apparatus. L, tub for setting a rinsing solution. M, pipette for rinsing acid. N, transfer pipette. The heavily drawn lower end of the capillary at B indicates where the capillary is surfaced ground and coated with Clarite.

and holds the acid rinsing solution described below. A 5 ml pipette (L) furnished with a rubber bulb is suspended on the outside of the water bath. A suction line with a 15 cm long fine tipped capillary glass tube is used for sucking off waste fluids. It is desirable to furnish the apparatus with a device to shake or vibrate it during the absorption of oxygen. A small motor with a crankshaft imparts motion to the micrometer burette through a rubber band attached by means of a metal ring to the crank shaft. The rate should be about 10 rps and the amplitude about 1 cm at the crank shaft.

### Reagents

- 1 Acid Rinsing Solution. In 400 gm distilled water dissolve 7 gm anhydrous sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) 1 ml conc  $\text{H}_2\text{SO}_4$  (delivered from a calibrated tuberculin syringe) and 21 ml glycerina (delivered from a calibrated syringe). Before use dissolve 50 mg potassium dichromate ( $\text{K}_2\text{Cr}_2\text{O}_7$ ) in 100 ml of the mixture and store in water bath of apparatus.
- 2 Carbon Dioxide Absorbent. In 100 gm of water dissolve 11 gm potassium hydroxide ( $\text{KOH}$ ) and 40 mg potassium dichromate ( $\text{K}_2\text{Cr}_2\text{O}_7$ ). Store in a syringe fitted with a No. 24 needle.
- 3 Oxygen Absorbent. Prepare (a) a solution of 100 gm water and 8 gm potassium hydroxide and (b) a powder mixture of 10 gm sodium hydrosulfite ( $\text{Na}_2\text{S}_2\text{O}_4$  Lycopen) and 0.1 gm sodium anthraquinone beta sulfonate. For use place 0.6 gm of powder (b) in an 8 ml vial add 5 ml of solution (a) close with finger dissolve anaerobically and transfer immediately to a 5 ml syringe fitted with a No. 24 needle.
- 4 Clarite Coating Solutions. Clarite X is a thermoplastic resin commonly used for mounting histologic sections. Two strengths of solution are required (a) dissolve 3 gm Clarite X in 2 ml of toluol (b) dissolve 1 gm Clarite X in 4 ml of toluol.
- 5 Petrolatum. White commercial petrolatum is used.

### Procedure—Assembling Apparatus

- 1 First insert the rubber stoppers and indicator band (115 mm binding tape) on the analysis capillary.
- 2 Then with the lower joint dried and the female part very lightly greased with petrolatum fix the upper part of apparatus back in place and jam it on.
- 3 Put water (room temperature) into the jar up to the neck of the compensating vessel. Connect to the leveling bulb and let the mercury up into the analysis capillary.
- 4 Leaving the apparatus connected with the leveling bulb insert an empty  $\text{CO}_2$  absorbing syringe via the needle through the middle of the rubber cork and remove the air.  $\text{Hg}$  will follow up into the syringe.
- 5 Remove the needle from the stopper at the same time pressing down the whole apparatus so that it does not disconnect at the ground point.
- 6 Repeat on the other side.

Then charge the respective sides with the corresponding absorbents (all bubbles must be excluded).

- 8 If the lower joint traps any air bubble remove the same by working the rinsing acid through it.

### Procedure—Disassembling and Coating of Capillary

- 1 Suck out all the water from the water bath and remove any rinsing acid from the apparatus. Lower the mercury below the ground point by lowering the reservoir. Lift off the upper part of the apparatus from the ground point.

2. Holding it over a beaker remove the plugs so that the Hg and absorbent empty into a beaker. Wash thoroughly under a water tap and remove most of the water by suction.
3. Pour toluol into the compensating chamber holding the finger in the lower end and rinse through the middle with toluol. Repeat wash thoroughly with acetone and dry by suction. Coating solution is 10 to 50% Clarite in toluol.
4. Draw out a glass tube into a fine capillary which will reach the ground ends inside of the capillary of the analysis apparatus both from above and below. Attach a rubber tube to the glass tube and fold the rubber tube then squeeze it together.
5. Place the end of the glass capillary in Clarite and relax the squeezing so that the Clarite fills the capillary. Then let go the fold of the rubber tube to prevent air from being sucked in when capillary is taken out of Clarite.
6. Wipe outside of the capillary. Fold rubber tube over and insert glass capillary into upper ground end of analysis capillary. Squeeze out drop (about 2 mm) of Clarite and withdraw capillary.
7. Move the drop of Clarite down the analysis capillary until it reaches the end of the ground part (i.e. about 1 cm) by placing finger on the opening of the compensating chamber and pressing. Similarly move the drop back and withdraw it by suction through its own capillary pipette rinsed with toluol to remove Clarite still there. Dry with vacuum.
8. The other end of the capillary is coated in the same way from the other end. The middle third of the capillary must be untouched by the coating solution.
9. A slight suction is applied to the apparatus overnight.
10. The upper half of the compensating chamber is also coated with Clarite so as to protect the joint above from being wetted by water. Clarite must be removed from this joint by toluol swabbing if it accidentally gets on it. The Clarite is swabbed on to the upper half by a plug of cotton on the end of a match.

#### Procedure—Washing Out Apparatus

1. Fill the thermobarometer up to the neck with acid rinsing solution.
2. With Stopcock in position III draw the rinsing solution up and down past the side arms by moving the leveling bulb (not the micrometer).
3. Suck out the rinsing solution and repeat the process until the solution remains yellow. This takes two to four times.
4. Connect the micrometer with B (position II) and adjust mercury to the top of the capillary.
5. Leave the manometer half full of acid rinsing solution.

#### Procedure—Analysis of Samples

1. Gas is contained in the usual sampling tube (Baccott type) with tap and Hg reservoir. From a syringe squirt a small drop of Hg into the bottom of the tube just above the tap. Hold a dry transfer pipette (M) vertically and press the rubber tip firmly against the top of the tube. Open the tap cautiously so that the Hg drop is gently blown up into the pipette bulb and allow the gas to bubble past it for a few seconds. Withdraw the pipette. The Hg drop starts to fall slowly down the tube. Before it has had time to get far fit the rubber tip on to the top of the capillary in the compensating vessel (A) there being a layer of rinsing solution a few mm thick to protect it from the outside air.
- Screw the micrometer (H) out until the gas meniscus reaches the lower mark, turn the stopcock (position II) to connect the micrometer and bulb and adjust the micrometer to zero. This is the reading M.

- 3 Turn the stopcock to connect the micrometer with the analysis chamber (position I) and screw the micrometer out to about 15 to 20 on the scale. Screw micrometer back a few mm so that the Hg drop moves up thereby indicating slight positive pressure in chamber. In this way no acid will leak in when the transfer pipette is removed.
- 4 Screw in about 2 mm of the rinsing fluid into the capillary and suck away all the rest of the rinsing fluid from the compensating chamber. (A) Screw the acid drop down until the lower meniscus of the drop coincides with the lower mark on the capillary.
- 5 Insert the tap of the compensating chamber and close the tap (B 1).
- 6 Readjust the drop meniscus accurately and repeat until the micrometer reading is constant (M).
- 7 Tilt the apparatus downward to the carbon dioxide absorbent side (C) admitting a small amount of the latter into the analysis chamber at the same time screwing in the micrometer to compensate for the downward movement of the acid drop in the capillary occasioned by the carbon dioxide absorption (tap with finger to absorb carbon dioxide). When the acid drop is stationary readjust to mark and read micrometer again (M).
- 8 Again tilt downwards on the oxygen absorbent side (D) drawing in similarly rather more of the latter again screwing in the micrometer to take care of the oxygen absorption.
- 9 Then switch on the motor shaker until the absorption is finished the acid drop being kept in a safe part of the capillary by the continuous screwing in of the micrometer.
- 10 Allow about 2 minutes more for the temperature rise (due to the oxygen +  $\text{Na}_2\text{S}_2\text{O}_4$  reaction) to disappear.
- 11 Then readjust the drop to the mark and read micrometer (M).
- 12 Open the compensating chamber to the air by turning the key of the tap (B 1) and after doing that transfer the tap to its socket (G).
- 13 Screw micrometer in until the top of the absorbing solution in the chamber coincides with the bottom mark in the capillary. Read micrometer again (M<sub>2</sub>). It should be zero.
- 14 Turn the 3 way tap to connect the leveling bulb with the absorption chamber (position III with S 2) and eject the absorbing solution into a vacuum junction tip not allowing the absorbing solution to spill into the compensating chamber.
- 15 Rinse twice with the acid solution as described at the beginning of the procedure.
- 16 To recharge the side arms (E and F) with respective absorbents puncture the rubber stoppers through their central hole taking great care not to introduce any gas bubbles. The side arms should not be filled more than  $\frac{1}{2}$  full with the absorbents.

#### Calculation

- 1 Carbon dioxide ml per 100 ml original sample dry =

$$\frac{M - M_1}{M_1} \times 100$$

- 2 Oxygen ml per 100 ml original sample dry =

$$\frac{M - M_2}{M_2} \times 100$$

- 3 Nitrogen ml per 100 ml original sample dry =

$$\frac{M_2}{M_1} \times 100$$

**Example**

$M$  was 9888  $M$  was 955  $M$  was 8000

$M_1$  and  $M_2$  were 000 which showed a satisfactory analysis

Carbon dioxide ml/100 ml dry =

$$\frac{9888 - 955}{9888} \times 100 = 3.671$$

Oxygen ml/100 ml dry =

$$\frac{955 - 8000}{9888} \times 100 = 15.423$$

Nitrogen ml/100 ml dry =

$$\frac{8000}{9888} \times 100 = 80.906$$

**Precautions**

- 1 *Apparatus* The micrometer must be all steel the plunger uniform from one end to the other as measured by another standard micrometer the screw without any play and the lap fit at the plunger bearing perfect as determined by visual inspection with the plunger screwed level. The glass portion of the glass metal junction in the micrometer burette unit has to be ground perfectly flat in order to make a mercury tight junction.
- 2 *Transfer Pipette* The tip of this is critical and must conform to the directions given. If the tip becomes plugged clear it with a fine wire or hold the rubber tip of the pipette against the nozzle of a 5 ml syringe and press air through the tip. Between samples the mercury drop and the drop of the acid rinsing solution are sucked out of the pipette by application of the suction tip to the rubber tip on the pipette.
- 3 *Analysis* The most critical part of the analysis is correct delivery of the gas sample avoiding suction on the micrometer burette. A proper line up between the bore of the transfer pipette and the bore of the reaction chamber capillary is easily accomplished by means of the funnel ground capillary and conical rubber tip to match a most valuable improvement of the original model contributed by Mr. Warner Love of the Department of Zoology Swarthmore. The transfer pipette must be kept clean. The easiest way to avoid clogging is to make use of a long and relatively wide bore for the air resistance in the tip instead of a short and very fine one.
  - a. The reason for not setting the micrometer at zero exactly at the beginning of the analysis is as follows. The rubber stoppers E and F are elastic and a trace of elasticity is likewise derived from unavoidably trapped minute air bubbles. Hence gross changes in mercury levels between the reaction chamber and the side arms will slightly affect the total volume of the apparatus. When the mercury level is lowered from the mark on the capillary to half way down the reaction chamber during introduction of the gas sample there is consequently a minute decrease in the total liquid capacity of the apparatus. The magnitude of this shrinkage is easily determined by reading the gas meniscus at the mark with the mercury almost to the mark then letting in all absorbing fluids from both side arms drawing down the mercury from the capillary and readjusting the gas meniscus to the mark. The reading will then be close to 0.5 division higher than before. Setting the micrometer at minus 0.5 instead of zero corrects empirically for this minute elasticity error.
  - b. If the concentration of absorbable gas in the sample is high larger than usual amounts of absorbent must be used. The CO absorbent will absorb 10 times its volume of CO without significantly altering vapor tension and the O absorbent will absorb 5 to 6 times its own volume of O before turning brown and giving poor results.



- c Carefully avoid sucking the rinsing solution down through S<sub>4</sub>. The rinsing solution will gradually remove the grease coating allowing the trapped traces of acid to move in the stopcock. This is usually first detected by the fact that the micrometer misses the zero after the end of the analysis. To check the tightness of S<sub>2</sub>, fill B with the acid solution and adjust the gas meniscus to the mark. With S<sub>2</sub> in Position I hang the leveling bulb low on the ring or handle J. If S<sub>2</sub> leaks the meniscus will move slowly down. If the bulb is hung high the meniscus will rise. Unless the meniscus stays at the mark, dry and regrease S<sub>2</sub> and tie the leveling bulb to the ring with a piece of string so that it cannot be lowered below the glass joint.

## SECTION VI

### PHYSIOLOGICAL MEASUREMENTS (Continued)

---

#### D RESPIRATORY METABOLISM

##### 1 Collection of Expired Air—The Tissot Gasometer

###### Reference

- Bock A. V. Dill D. B. and Talbott J. H. Studies in Muscular Activity I Determination of the Rate of Circulation of Blood in Man at Work J. Physiol. 61: 121-133 (Oct.) 1923

###### Principle

The subject's nose is clamped and through a mouthpiece he passes all his expired air into a calibrated counterpoised metal cylinder which is seated by means of a water seal in a second cylinder and which is free to move with respiration. From the calibration constant and the time the pulmonary ventilation may be calculated.

###### Apparatus

- 1 Tissot gasometer of 100 or 250 liter capacity (see Fig. 43) of use during walking. The inner bell must be counterpoised accurately so as not to impede respiration. The gasometer should be calibrated against a standardized wet gasometer.
- 2 One mouthpiece-nose clip assembly. Depending on whether the subject is to exercise voluntarily or not the valves should be flutter valves or butterfly valves. For all exercise flutter valves are required. From the expiratory valve corrugated 1 to 2 inch gas mask tubing leads to the intake of the Tissot gasometer. One comfortable screw type nose clip is needed.
- 3 One barometer either mercury or aneroid sensitive to 0.5 mm Hg.
- 4 One watch with a second hand preferably a stopwatch.
- 5 Gas sampling tubes Barcroft design in racks of four filled with mercury.
- 6 A three way large bore aluminum stopcock inserted between the exit valve of the mouthpiece assembly and the gasometer. Two positions are required: (a) subject to gasometer (b) subject to outside.

###### Reagents

- 1 Alcohol 80% for sterilizing the mouthpiece after use.

###### Procedure

- 1 Check the gasometer for counterpoise and for water level which should be 1 inch below the rim of the outer cylinder. If more water is needed be sure that its temperature is close to room temperature. Check that the hanging scale is not scraping against the sides of the guides.
- 2 Before giving the mouthpiece to the subject wash it in 80% alcohol and dry it.
- 3 Check leakage of valves by blowing through them and check the tightness of fit by blowing backward through them. There should be no leakage.
- 4 Place the nose clip firmly on the subject's nose and have him test for leaks past the clip by trying to blow air through his nose.
- 5 Be sure that the three way tap connects subject to outside air.
- 6 Start the subject performing the task whose energy requirements are to be measured.
- 7 Turn three way tap to connect the subject with the Tissot gasometer.



Fig 43—Use of Tissot gasometer for collecting specimens of expired air (Reproduced by courtesy of Dr H. S. Belding)

- 8 Allow bell to fill about one third full with air and turn the tap to connect the subject with air
- 9 Remove rubber stopper from thermometer hole and by manual pressure expel all air from the gasometer
- 10 Repeat steps 7 8 and 9 twice more This insures that there will be no mixing error from dead space air already in the gasometer
- 11 Now read the gasometer scale (Reading 1) to the nearest mm
- 12 Start stopwatch and turn tap to connect subject with gasometer
- 13 At exactly 1 minute and every minute desired thereafter read the scale This is a check on the consistency of the subject and apparatus
- 14 At a desired time usually 5 minutes for low rates of pulmonary ventilation the tap is turned to connect subject with outside air and the final scale reading is made (Reading 2)
- 15 Turn on mixing fan and read the temperature in the bell to the nearest 0.1 C
- 16 Read the barometer to the nearest 0.5 mm
- 17 Take specimens of the contents of the bell from the thermometer hole into Barcroft gas sampling tubes making sure that
  - a The contents of the bell have been mixed and
  - b The Barcroft gas sampling tube is flushed three times in the following order contents of sampling tube to outside air contents of bell into sampling tube contents of sampling tube to outside air

#### Calculation

Pulmonary Ventilation =

$$\frac{(\text{Reading 2} - \text{Reading 1})}{(\text{minutes})} \times \text{Tissot factor} \times \text{STP factor}$$

Where

Pulmonary ventilation is expressed in liters of gas expired per minute reduced to 0 C and 60 mm Hg Readings 1 and 2 are the last and first readings expressed in cm to the nearest mm Tissot factor is the liters of gas equivalent to a rise of 1 cm in the Tissot bell and STP factor is the factor for reducing any volume of wet gas at any temperature and pressure to 0 C and 60 mm Hg dry (see Fig. 45 in section on metabolic calculations)

#### Example

A subject was lying quietly

Reading 1 was 51 cm Reading 2 was 60.5 cm T was 4 C barometer was 750 mm and the time was 5 minutes The Tissot factor for this particular gasometer was 1.3 l/cm. (Note that every Tissot has its own calibration figure to be determined at purchase)

Pulmonary ventilation l/min dry STP =

$$\frac{(60.5 - 51)}{5} \times 1.3 \times 0.89 = 12.1 \text{ l/min}$$

#### Precautions

- 1 The chief sources of error in this method are
  - a Leakage around nose clip or in valves
  - b Inaccurate reading of Tissot scale
  - c Dilution of contents of bell with outside air which will not vitiate calculation of pulmonary ventilation but which will destroy the validity of all metabolic calculations
  - d Inaccurate use of Barcroft gas sampling tubes
- 2 If room air is significantly different in composition from outdoor air it may be necessary to pipe outdoor air into the intake valve of the subject's valve assembly All metabolic calculations assume outdoor air composition of the inspired air




Fig 44 —Valve assembly for collection of expired air when rate of pulmonary ventilation is moderate.

—C.E. GORDON

## 2 Collection of Expired Air—The Douglas Bag Technique

### Reference

Haldane J S and Priestley J G *Respiration* Oxford 1935 Clarendon Press

### Principle

The subject's nose is clamped and through a mouthpiece he passes all of his expired air into an air tight bag. Two one way valves allow him to breathe in only outside air and to expire only into the bag. The total pulmonary ventilation in a given time is measured by passing the contents of the bag through a gas meter and samples are taken for analysis.

### Apparatus

- 1 One barometer either mercury or aneroid sensitive to 0.5 mm Hg
- 2 One watch with a second hand preferably a stopwatch
- 3 One mouthpiece assembly (Fig 44). This should be usable even when the subject is exercising violently. A rubber mouthpiece is mounted in the expiratory hole of a gas mask of which the goggles are removed and the nose is cut out. The inspiratory cannula is removed and the head straps left intact. The rubber mouthpiece is attached to a pair of aluminum 1 inch tubes containing rubber flutter valves one inspiratory the other expiratory which lie snugly along the outside of the face of the mask. From the expiratory valve a piece of corrugated 1 inch gas mask tubing leads to the head of the subject to the bag as assembly described below.
- 4 One comfortable nose clip screw type
- 5 One Douglas bag assembly. This consists of one 1.5 l Douglas bag (Davol Rubber Company Providence Rhode Island) and a light pack board for the subject to wear on his back. The bag consists of two layers of canvas with rubber sheeting cemented between and must be absolutely leak proof when the exit tube is clamped off. The upper corners of the bag can be hooked to the upper corners of the pack board. The exit tube of the bag should then reach the level of the subject's neck. A strong flat brass clamp should be placed on the exit tube of the bag which is then attached to a 1 inch way aluminum tap. One limb of the tap is attached to the bag and one to the gas mask tubing from the subject's expiratory valve. The third limb is to outside air. The three way tap should allow either of two circuits for the subject's expired air i.e. mouth—valve—gas mask tubing—way tap—and thence either (a) outside air or (b) into Douglas bag.
- 6 One dry gas meter of at least 75 cu. ft. per minute capacity (Tufts Meter Company Albany N.Y.). The meter should be mounted in the intake side of the gasmeter.
- 7 Rack of gas sampling tubes. Orthodox mercury filled tubes must be avoided in the field when temperatures are extreme. A satisfactory sampling tube consists of a 50 ml syringe oiled with heavy oil and fitted with a 1 inch length of pressure tubing over the nozzle. The tubing is clamped off with a screw clamp. The 3 way aluminum tap must have a side arm for taking gas samples. This is closed with a rubber nipple all the time except when samples are taken.

### Reagents

- 1 Heavy oil for the sampling syringes
- 80% alcohol to sterilize the mouthpiece after use

### Procedure

- 1 Check the Douglas bag for leaks the night before using it. This is done by hooking the bag to the outlet of the gas meter running in a measured volume of air clamping off the neck of the bag placing 10 lbs or more on the bag for an

## 332 *Physiological Measurements*

- hour and running the contents back through the gasmeter. The readings should be nearly the same. If there are leaks they can be patched temporarily with ordinary tire patching.
- 2 Before asking the subject to take the mouthpiece wash it with 80% alcohol and dry it.
  - 3 Place the nose clip firmly on the subject's nose. Have him test for leaks past the clip by trying to blow air through his nose. This step should never be omitted.
  - 4 Place the mouthpiece assembly on the subject's head. Fit the mouthpiece so that it is comfortable and tie the headstraps snugly.
  - 5 Empty the Douglas bag completely of air by laying it flat on the ground with the exit open, rolling it up snugly and clamping off the exit tube with the brass clamp.
  - 6 Hang the Douglas bag on the pack board, connect the exit tube to the 3 way tap and finally connect the corrugated tube from the mouthpiece to the 3 way tap.
  - 7 Be sure that the 3 way tap is adjusted so that the subject's expired air is going to the outside.
  - 8 Start the subject performing the task whose energy requirement is to be measured. Let him continue this long enough to 'warm up' before Step 9.
  - 9 Take off the brass clamp. Simultaneously start the watch at 0 and turn the 3 way tap so that all of the expired air goes into the bag.
  - 10 Allow the subject to breathe into the bag while continuing his task until the bag seems about  $\frac{3}{4}$  full. During rest this takes 10 minutes; during light activity about 5; and during moderate activity about 3 minutes. Vigorous activity is difficult to control or measure. The bag may fill up in about a minute.
  - 11 Simultaneously note the time to the nearest second and turn the 3 way tap so that the subject is breathing to the outside. He may now stop work.
  - 12 Clamp off the neck of the bag, disconnect it along with the 3 way tap and remove it from the pack board to the neighborhood of the gas meter.
  - 13 Mix the bag's contents by laying it on the floor and kneading it a few times.
  - 14 Connect the 3 way tap with bag attached to one limb to the inlet of the gas meter.
  - 15 Unscrew the clamp from a sampling syringe and connect the pressure tubing to the sampling nipple on the 3 way tap.
  - 16 Take a reading on the dial of the gas meter  $V$ .
  - 17 Open the 3 way tap so that the contents of the bag flow into the gas meter *not* to the outside air.
  - 18 Put steady pressure on the bag and when it is about  $\frac{1}{4}$  empty take a sample by drawing the syringe full and empty 3 times while the gas is flowing steady through the 3 way tap. Finally fill the syringe, clamp off the rubber tube and after taking a duplicate sample close the sampling nipple of the 3 way tap.
  - 19 Express all of the bag's contents through the gas meter. Toward the end lay the bag flat on the ground and roll it up tightly from the bottom.
  - 20 Clamp off the neck of the bag and take the final reading on the gas meter  $V$ .
  - 21 Read the temperature in the intake of the gas meter to the nearest  $\frac{1}{2}$  °C and read the barometer to the nearest  $\frac{1}{2}$  mm Hg.

### Calculations

$$\text{Pulmonary ventilation} = \frac{(V - V_1)}{(\text{time})} \times \text{STP factor}$$

Where Pulmonary ventilation is expressed in liters of gas expired per minute reduced to 0 °C and 760 mm Hg dry.

$V$  is the first reading of the gas meter in liters.

$V_1$  is the second reading of the gas meter in liters.

STP factor is the factor for reducing gas at any temperature and pressure wet to 0 °C and 760 mm Hg dry as shown in Fig. 45.

If the gas meter reads cubic feet the conversion factor by which readings are multiplied to give liters is 28.3. Time is expressed in minutes and decimal fractions.

**Example**

V 91.5 l 185.3 Temp C Barometer 760 mm  
 Time of collection 4 minutes %G seconds = 4 and %G/60 = 4.43 mins

$$\text{Pulmonary ventilation} = \frac{(185.3 - 91.5)}{(4.43)} (0.90) = 19.0 \text{ l per min STP}$$

**Precautions**

- 1 Test the bag for leaks every few days  
 Be sure to have the subject test for leaks through the nose or around the mouth piece
- 3 Timing must coincide with turning the taps.
- 4 Every time a tap is turned be sure that the gas flow will be in the direction desired
- 5 Do not forget to read the temperature and barometric pressure every time samples are measured
- 6 In vigorous activity when a significant oxygen debt is accumulated it is necessary to have the subject exercise for a minute or so then stop and collect expired air while he is resting after work for about 5 minutes. If this is not done the estimate of caloric expenditure during work will be erroneously low

**3 Respiratory Metabolic Calculations****Reference**

Haldane J S and Priestley J G Respiration Oxford 1935 Clarendon Press

**Principle**

From the pulmonary ventilation and the carbon dioxide oxygen and nitrogen content of the expired air compared with the inspired air the consumption of oxygen and production of carbon dioxide are calculated. From these data energy expenditure may be computed.

**Pulmonary Ventilation**

This is expressed as liters of air expired per minute the volume of gas being reduced to the standard temperature and pressure 0 C and 760 mm Hg dry. The formula is

$$\frac{B - W}{60 (1 + 0.00367 t)}$$

Where B is ambient barometric pressure W is the vapor tension of water mm Hg at the temperature of the gasometer and t is the temperature of the gasometer in C. To save computation one uses the line chart which was devised by Dr R C Darling (Fig 45)

**Respiratory Quotient (R Q)**

This is the ratio of CO<sub>2</sub> expired divided by O<sub>2</sub> consumed i.e.

$$\frac{\text{CO}_2 \text{ production}}{\text{O}_2 \text{ consumption}}$$

and is used in the calculations which follow

**True Oxygen, True Carbon Dioxide and R.Q**

True Oxygen represents the number of ml of oxygen consumed for every 100 ml of air expired. It is based on several considerations as follows. One desires to know how much oxygen is removed from the air breathed in but the only measurements made are the volume of air expired and its oxygen carbon dioxide and nitrogen contents. The volume of inspired air is usually not the same as that of expired air. This is because the R.Q. has to be exactly 1.00 for inspired to equal expired air. If the R.Q. is less than one as is usually the case in rest or moderate



LINE CHART FOR DETERMINING FACTORS TO REDUCE SATURATED GAS VOLUMES TO  
 DRY VOLUMES AT 0°C AND 760 MM HG




Fig 45—Line chart for converting wet gas volumes to STP and dry as drawn by Dr H. C. Darling

exercise then the oxygen removed from inspired air is only partly replaced by carbon dioxide and if one l has been inspired less than one l will be expired. The nitrogen concentration in this case will be higher in expired than in inspired air. One knows the nitrogen concentration in outdoor air to be 79.04% and the concentration in expired air is determined by Haldane analysis. The volume of inspired air may be calculated from the expired air by the formula

$$\text{Vol inspired} = \text{Vol expired} \times \frac{\%N \text{ in expired air}}{79.04} \quad (\text{Equation 1})$$

The total volume of oxygen breathed in (not all consumed) is then

$$\text{Vol O inspired} = \text{Vol air inspired} \times \frac{\%O \text{ of inspired air}}{100} \quad (\text{Equation 2})$$

The percentage of oxygen in practically all outdoor air is 0.93. Hence

$$\text{Vol O inspired} = \text{Vol air inspired} \times \frac{0.93}{100} \quad (\text{Equation 3})$$

The volume of oxygen expired (amount not consumed) is

$$\text{Vol O expired} = \frac{\%O \text{ in expired air}}{100} \times \text{vol air expired} \quad (\text{Equation 4})$$

The amount of oxygen consumed is

$$O \text{ consumed} = \text{Vol O inspired} - \text{Vol O expired} \quad (\text{Equation 5})$$

Substituting values from equations (3) and (4)

$$O \text{ consumed} = \text{Vol air inspired} \times \frac{0.93}{100} - \text{Vol air expired} \times \frac{\%O \text{ in expired air}}{100}$$

Substituting values from equation (1)

$$O \text{ consumed} = \text{Vol air expired} \times \frac{N \text{ in expired air}}{79.04} \times \frac{0.93}{100} - \text{Vol air expired} \times \frac{\%O \text{ in expired air}}{100}$$

Simplifying

$$\text{Vol O consumed} = \frac{\text{Vol air expired}}{100} (\%N \text{ in expired air} \times 0.85 - \%O \text{ in expired air}) \quad (\text{Equation 6})$$

The factor  $(\%N \text{ in expired air} \times 0.85 - \%O \text{ in expired air})$  is the True Oxygen the number by which the volume of expired air (divided by 100) is multiplied to give the oxygen consumption.

By steps similar to the derivation of oxygen consumption the carbon dioxide expired becomes

$$\text{Vol CO expired} = \frac{\text{Vol air expired}}{100} (\%CO \text{ in expired air} - \frac{\%N \text{ in expired air}}{79.04} \times \%CO \text{ in inspired air})$$

In the case of ordinary outdoor air this becomes

$$\text{Vol CO expired} = \frac{\text{Vol air expired}}{100} (\%CO \text{ in expired air} - 0.03) \quad (\text{Equation 7})$$

The factor  $(\%CO \text{ in expired air} - 0.03)$  is the True Carbon Dioxide the number by which the volume of expired air (divided by 100) is multiplied to give the carbon dioxide production.

The respiratory quotient R Q is obtained by dividing equation (7) by equation (6) and is

$$R Q = \frac{(\%CO \text{ in expired air} - 0.03)}{(\%N \text{ in expired air} \times 0.85 - \%O \text{ in expired air})} \quad (\text{Equation 8})$$

To facilitate the computation of True O and RQ a line chart has been constructed. A string stretched between the %O<sub>2</sub> and %CO<sub>2</sub> in expired air passes across lines giving RQ and True O' directly. (See Fig 66 at end of this book)

### Oxygen Consumption and CO Production

The oxygen consumption in liters per minute is

$$O_2 \text{ l/min} = \frac{\text{Pulmonary ventilation l/min}}{100} \times \text{True O} \quad (\text{Equation 9})$$

The carbon dioxide production in liters per minute is

$$CO \text{ l/min} = \frac{\text{Pulmonary ventilation l/min}}{100} \times \text{True CO} \quad (\text{Equation 10})$$

### Energy Expenditure

It is conventional to express the metabolic cost of a given task in terms of Calories per hour. The Calorie is the amount of heat necessary to raise the temperature of 1 kilogram of water 1 C and is the basic unit of all nutritional energy calculations. To compute the cost of a given task the formula is used

$$\text{Cals/Hr} = (O \text{ consumption l/min}) \times 60 \times 50 = O \text{ l/min} \times 300$$

The factor 50 is an average figure representing the Caloric equivalent of 1 liter of oxygen consumed during moderate activity.

Since most of the common forms of work require an energy expenditure which is directly proportional to the body weight, most tables of metabolic costs of activity are expressed on the basis of a standard man weighing 150 lbs or 70 kg. To compensate for variations in weight, energy expenditure in work may be expressed in terms of Cals/hour and kg body weight or converted to a standard body weight as follows:

$$\text{Cals/Hr (standard man)} = O \text{ l/min} \times 300 \times \frac{150}{\text{body weight lbs}}$$

For some special purposes energy expenditure is computed on the basis of body area, not body weight. This will be discussed in the section on basal metabolism.

### Example of Metabolic Calculations for a Walking Man

Company marched 3.5 miles in 50 mins on the level

Weight of subject stripped 145 lbs

Weight of clothes helmet pack rifle canteen 29 lbs

Weight of metabolic apparatus 13 lbs

#### A. Douglas Bag Calculations

Collection started 7 hr 30.0 ended 7 hr 37.0 = 7.00 mins

Gas meter readings—Start 448.0 L. End 667.9 L = 219.9

Temperature 21.9 C Barometer 757 mm

STP factor (from Fig 45) = 0.880

$$\text{Pulmonary ventilation l/min STP dry} = \frac{(219.9)}{(7.00)} \times 0.880 = 27.0$$

#### B. Haldane Analysis

$$\%CO_2 = 3.81 \quad \%O_2 = 16.80$$

#### C. Calculation of O<sub>2</sub> consumption

From line chart for True O<sub>2</sub> (See Fig 66 in back)

$$RQ = 0.89 \quad \text{True O} = 4.23$$

$$O_2 \text{ consumption l/min} = \frac{\text{Pulmonary Vent}}{100} \times \text{True O} = \frac{27.0}{100} \times 4.23 = 1.14$$

$$CO_2 \text{ production l/min} = \frac{\text{Pulmonary Vent}}{100} \times \text{True CO} = \frac{27.0}{100} \times (3.81 - 0.03) = 1.0$$

**D Calculation of energy expenditure**

$$\begin{aligned}\text{Cals/min} &= \text{O l/min.} \times 5 \\ &= 114 \times 5 = 57\end{aligned}$$

$$\text{Cals/50 min} = 57 \times 50 = 285$$

Cals for 10 min rest assumed to be 0

Energy expended in one hour that consists of 3 miles marching in 50 minutes and 10 minutes rest  $= 285 + 0 = 285$  Cals

**Protein Metabolism During Work**

The above calculations and examples neglect two factors which are not important in field measurements but are theoretically important. Some of the Calories are derived from protein and one liter of oxygen consumed represents different amounts of heat depending on the ratio of fat, carbohydrate and protein being consumed by the body at the time of measurement. Protein metabolism is determined from the urinary nitrogen excretion and the non protein RQ is determined from the over all RQ (obtained as described above) and the protein metabolism. For details of this type of calculation see Peters and Van Slyke *Quantitative Clinical Chemistry* pp 65-115 Vol II Baltimore 1931 Williams and Wilkins

**4 Basal Metabolic Rate****Reference**

Du Bois E F *Basal Metabolism in Health and Disease* ed 3 Philadelphia 1937  
Lea and Febiger

**Principle**

The basal metabolic rate (BMR) is the energy expenditure of a subject resting quietly after at least 8 hours sleep and at least 1 hour after the last meal. When expressed on the basis of Calorie production per square meter of body area it is remarkably constant for healthy people of the same age and sex.

**Apparatus and Reagents**

Same as under measurement of expired air. In addition an oral clinical thermometer, a clinical weighing scale and a measuring stick are needed.

**Procedure**

- 1 For an accurate test the subject must have had at least 8 hours sleep and no food for at least 1 hour. This also includes alcoholic beverages, coffee, tea and tobacco. He lies down comfortably for at least 1/2 hour in a quiet place before the test.
- 2 Metabolism is determined while he is still resting as described in the section on expired air. Count the respirations for 1 minute during the process. If he is too excited the test will be invalid.
- 3 Before he gets up and immediately after the mask and mouthpiece are removed measure and record the pulse and oral temperature. This is necessary to insure that one is not dealing with an infection. Fever invalidates the test.
- 4 Measure the subject's height in stocking feet and weight with a minimum of clothing.
- 5 Record his age to the nearest year.

**Calculation**

- 1 Calculate oxygen consumption and RQ as described previously. If the RQ is outside the range 0.80-0.90 the test is invalid. From his height and weight determine the subject's surface area from the line chart (Fig 46). This nomogram was computed from the DuBois-Meeh formula for surface area.

# SURFACE AREA FROM HEIGHT AND WEIGHT





Fig 46—Nomogram constructed from the Du Bois Meek formula for surface area.


# HEAT PRODUCTION (Cal/M<sup>2</sup>/hr) FROM SURFACE AREA AND OXYGEN CONSUMPTION (R Q 0.85)



C E GORDON—

Fig 47—Nomogram for computing heat production

# BASAL METABOLIC RATE FROM HEAT PRODUCTION, AGE AND SEX



BY G.E. GORDON -

Fig 48—Nomogram for computing BMR from heat production, age, and sex. The Boothby Sandiford standards were employed.

- 3 From oxygen consumption and surface area compute heat production (Cals per square meter per hour) The line chart in Fig 47 was drawn to facilitate this computation.
- 4 Complete the subject's basal metabolic rate (i.e. ° above or below normal) by the formula

$$\text{BMR} = 100 \times \frac{(\text{Observed Cals/m}^2/\text{hr} - \text{Average normal})}{\text{Average normal}}$$

(Where average normal refers to the Boothby Sandiford values)

- 5 The BMR line chart in Fig 48 facilitates the above calculation. It was computed from the Boothby Sandiford normal values

#### Example

Subject R D Sex Male  
 Age 26 years  
 Ht 5 ft 7½ inches  
 Wt 19 lbs  
 Pulse rate 69  
 Temp 98.4 F  
 Pulmonary Ventilation l/min at STP 823  
 Expired air %CO 35, %O 17.08  
                   RQ 0.89 True O 3.95  
 Oxygen Consumption l/min 03.5  
 Surface Area 1.11 sq.m. (Fig 46)  
 Heat production 445 Cals/m<sup>2</sup>/hr (Fig 47)  
 BMR +12 (Fig 48)

#### Precautions

- 1 The commonest source of error are
  - a The subject is not truly basal.
  - b The subject is not comfortable and at ease
- 2 The above calculations and examples neglect the two factors protein metabolism and the varying heat equivalent of one liter of oxygen depending upon the ratio of protein fat and carbohydrate being consumed by the body. These factors are not important for field use but are theoretically important. The matter is discussed at the end of Metabolic Calculations.
- 3 If the pulse rate is above 80 or the respirations above 25 the subject is either not basal or else has some pathological condition. If the temperature is elevated he is either not basal or else has some disease. (Exercise of moderate severity elevates the body temperature.)

## 5 Metabolism During Work—The Treadmill

#### Reference

Robinson S. Experimental Studies of Physical Fitness in Relation to Age. *Arbeitsphysiologie* 10: 51-33 (No 3) 1938

#### Principle

The motor driven treadmill is the best apparatus yet devised for general studies of the response of man to exercise. A typical protocol is described for studying respiratory and cardiovascular responses during moderate exercise (walking) and exhausting exercise (running).

#### Apparatus

- 1 The treadmill (see photograph Fig 49) The treadmill has a motor driven end less leather belt that can be operated at any speed up to 20 miles per hour. The front end can be raised up to make a grade of 20°. It is possible to make the



work so easy that an ordinary individual can walk for hours and so severe that the best runners can be exhausted in a very short time (seconds) (Purchased at A R Young Co, Indianapolis Indiana)

- 2 The stationary sides of the treadmill carry a superstructure on which are suspended the rubber hose connection (1 inch inside diameter), and a two way flutter valve system. The hose connection brings fresh air from the outside to the valve and empties the expired air into the gasometer
- 3 The gasometer is a 500 liter copper bell type that is immersed in water. Each side of the bell is weighed down so that it balances perfectly when all the outlets are open. The measuring scale is in centimeters. Installed at the top of the bell is an electric stirrer for mixing the gas before sampling. At the inlet of the gasometer is a baffle mixer of 10 liter capacity from which samples of gas may be drawn while the subject is working, without interfering with the total collection in the gasometer




Fig 49—Electrically driven treadmill manufactured by A. R. Young Indianapolis Indiana.  
(Photograph by courtesy of A. R. Young Company Indianapolis Ind)

- 4 An accurate sweep second clock and a stopwatch
- 5 Spare Douglas bags that have a capacity of 700 liters
- 6 Rack of gas sampling tubes (Barcroft type) with an approximate capacity of 30 ml each
- 7 A cardiometer for recording pulse rates with electrodes and electrode paste. The electrodes are able to pick up electrical currents that are generated at each contraction of the heart. The current is then amplified to operate a relay and a pen on a moving tape thus yielding a record of the heart rate as work is going on. The three electrodes are fastened on two lengths of elastic tape. The first electrode is placed just to the left of the sternum at the level of the second rib

- The second electrode is placed over the apex of the heart and the third or ground which is on the same elastic tape as the second is on the right side
- 8 Micro pipettes for lactic acid and blood sugar of 0.5 ml and 0.1 ml capacity
- 9 Watch glasses for collecting the blood during exercise
- 10 A lancet for pricking finger to get blood.
- 11 Centrifuge tubes 15 ml and 50 ml round bottom with rubber stoppers
- 12 Barometer
- 13 A Haldane gas analysis apparatus

#### Reagents

- 1 Clean mercury for gas sampling tubes  
Heparin solution (1 drop dried on a watch glass)
- 3 Dilute acid tungstate for blood sugar 20 ml of 10% sodium tungstate solution  
0 ml of  $\frac{2}{3}$  N sulfuric acid ( $\text{H}_2\text{SO}_4$ ) Dilute to 1000 ml with water
- 4 Distilled water for lactic acid tubes

#### Procedure—Moderate Exercise

- 1 The subject dressed in gymnasium clothes is connected to the cardiographometer with the three electrodes fastened around his chest
- 2 The subject sits quietly at the side of the treadmill while his resting pulse is recorded and his blood pressure is determined  
Subject straddles the treadmill while the treadmill is started (3.5 m.p.h. on 8.6% grade)
- 3 Subject starts walking on signal from experimenter Notation on heart rate tape is made
- 4 After approximately five minutes walking the subject applies mouthpiece and nose clip and the valve is turned so that expired air enters the gasometer
- 5 In the course of the next four minutes the gasometer is washed out with three portions of expired air and finally a gasometer reading is taken after expulsion of the last wash portion
- 6 At the ninth minute of the walk collection of expired air is started and continued for five minutes the gasometer readings being noted every minute
- 7 At approximately ten minutes of the walk finger blood is taken for sugar and lactate
- 8 After a total of 15 minutes walk the subject steps off the treadmill sits down and rests for ten minutes during which time the following procedures are carried out
  - a Duplicate gas samples taken from gasometer and temperature noted Change gas circuit to mixing chamber
  - b Subject's pulse rate and blood pressure recorded at frequent intervals

#### Procedure—Exhausting Exercise

- 1 The usual grade for exhausting work is the same as for the walk (8.6%) but with the speed of the treadmill increased to 7.0 m.p.h. Ordinarily the subject is not required to run longer than five minutes To obtain exhaustion within five minutes in well trained subjects the grade may need to be increased If the grade is to be changed it is changed during the half minute before the run starts but not until the expired air and blood have been sampled as above
- 2 In total gasometer readings taken again
- 3 Subject prepares for run straddles treadmill holding the sides and is instructed to signal with his hand if he desires to stop before five minutes
- 4 The treadmill is started at 7.0 m.p.h. and on signal the following procedures occur simultaneously
  - a The subject starts running
  - b Valve is turned into gasometer
  - c Cardiographometer is marked

- 5 Gasometer readings are taken every half minute during run
- 6 Expired air samples are taken over an entire minute for the second, third fourth and fifth minutes (the first minute is not of sufficient value), or for as long as the subject runs. Samples of less than 20 seconds are of doubtful value

#### Procedure—Recovery After Exhausting Exercise

- 1 When subject stops running by jumping to a straddling position (keeping the mouthpiece in and the nose clip on) the following must be done precisely at the time he stops running
    - a Valve at treadmill superstructure is turned to collect expired air in Douglas bags for three minutes
    - b Cardiotachometer tape is marked
    - c Time is noted
    - d Treadmill is stopped. Subject sits on tall stool in place on treadmill still breathing through mouthpiece
  - 2 As soon as possible the following additional steps should be taken
    - a Final gasometer reading is taken
    - b Temperature of expired air in gasometer is noted
    - c Gasometer is emptied. Connection removed from mixing chamber and placed on direct connection to gasometer. Valve turned for collection and reading noted
  - 3 Blood pressure is taken as soon as possible preferably during the first minute after subject stops running
  - 4 At the end of three minutes valve at treadmill superstructure is turned to gasometer
  - 5 Gasometer readings are taken every minute during the remainder of the 15 minutes
  - 6 Finger blood samples for blood sugar and lactate are taken at 5, 10 and 15 minutes during recovery
  - 7 Cardiotachometer is left running at least for the first five minutes of recovery and may be turned off then and used only to take heart rate during 10th and 15th minutes of recovery
  - 8 At 15 minutes gasometer valve is closed. temperature is recorded and final reading is noted
  - 9 Two expired air samples are taken from gasometer (3, 15 minutes of recovery)
  - 10 After gasometer is emptied valves are turned connecting Douglas bag to gasometer and initial reading is noted. Two samples of contents of bag are drawn from small side tube on bag as gas is forced to large gasometer for measurement of volume. Bag is finally rolled to expel all air, connection clamped and gasometer volume and temperature are noted. (Three minutes of recovery)
  - 11 Barometric pressure is recorded
- (The work experiment as carried out above required at least three and preferably four men to guarantee smooth operation.)

#### Calculations and Data Obtained from Work Experiments

- 1 Blood pressure at rest (sitting) and immediately following exhausting run
- 2 Continuous heart rate record counted to obtain following
  - a Rest (sitting)
  - b Maximum during walk (sixth minute)
  - c Maximum during run (usually last minute)
  - d First five minutes of recovery. May be counted by six second intervals for first minute of recovery and longer intervals up to five minutes
  - e At 10 minutes during recovery
  - f At 15 minutes during recovery

2 Blood sugar and blood lactate in mg per 100 ml.

a. During walk

b. At 5, 10 and 15 minutes during recovery

4 Metabolism during walk, run and recovery are calculated as described in previous sections

5 Oxygen debt. This measurement is the difference in oxygen consumption for fifteen minutes following the run and the basal requirement for that period. To be complete the period should be 60 minutes.

#### Example of a Protocol of a Typical Experiment

| SUBJECT B<br>Gals 86°C                    |                              | DATE<br>8/3/49 | BAR PRESSURE 61 mm Hg<br>GASOMETER 16MP 34 C |          |         | TISSOT FACTOR 5<br>× CM STP FACTOR 0.898<br>COMBINED FACTOR 4.68 |                |                                 |  |
|---|------------------------------|----------------|--|----------|---------|--|----------------|---------------------------------|--|
| WALK<br>3.5 mph                           | SPIROMETER<br>READINGS<br>cm | DIFF<br>CM     | PULMONARY<br>VENTIL<br>L MIN                 | CO<br>CC | O<br>CC | RQ   | TRK<br>O<br>CC | O <sub>2</sub><br>USED<br>L/MIN |  |
| 4 mins                                    | 46                           |                |  |          |         |  |                |                                 |  |
| 5   | 131                          | 85             | 30   | 459      | 16.8    | 0.97   | 4.0            | 18                              |  |
| 6   | 10                           | 9              | 36.9   | 4.1      | 15.5    | 0.91   | 5.1            | 19.1                            |  |
| Run 7.0 mph                               |                              |                |  |          |         |  |                |                                 |  |
| 0 mins                                    | 45                           |                |  |          |         |  |                |                                 |  |
| 1   | 169                          |                |  |          |         |  |                |                                 |  |
| 1½  | 36                           |                |  |          |         |  |                |                                 |  |
|   | 310                          | 1.1            | 80.0   | 4.0      | 16.48   | 1.06   | 4.40           | 25                              |  |
| ¾   | 440                          |                |  |          |         |  |                |                                 |  |
| 3   | 534                          | 19.4           | 90.8   | 4.9      | 16.64   | 1.14   | 4.16           | 37.8                            |  |
| 3½  | 631                          |                |  |          |         |  |                |                                 |  |
| 4   | 31                           | 19             | 9  | 4.89     | 16.64   | 1.17   | 4.17           | 38.4                            |  |
| 4½  | 830                          |                |  |          |         |  |                |                                 |  |
| 5   | 8.5                          | 19.4           | 90.8   | 4.97     | 16.61   | 1.19   | 4.17           | 37.9                            |  |
| Recovery                                  |                              |                |  |          |         |  |                |                                 |  |
| Bag                                       |                              |                |  |          |         |  |                |                                 |  |
| 0 min                                     | 4                            |                |  |          |         |  |                |                                 |  |
| 2   | 44.5                         | 39.8           | 6.0  | 4.01     | 18.09   | 1.56   | 9.00           | 15.8                            |  |
| Gasometer                                 |                              |                |  |          |         |  |                |                                 |  |
| 3 mins                                    | 48                           |                |  |          |         |  |                |                                 |  |
| 5   | 165                          |                |  |          |         |  |                |                                 |  |
| 7   | 247                          |                |  |          |         |  |                |                                 |  |
| 9   | 3.5                          |                |  |          |         |  |                |                                 |  |
| 11  | 38.8                         |                |  |          |         |  |                |                                 |  |
| 13  | 45.8                         |                |  |          |         |  |                |                                 |  |
| 15  | 51.6                         | 46.8           | 19.3   | 3.57     | 1.9     | 1.06   | 3.30           | 0.60                            |  |
| Partial O <sub>2</sub> consumption = 0.30 |                              |                |  |          |         |  |                |                                 |  |

Net O<sub>2</sub> debt in 15 minutes

$$3 \times (1.58 - 0.30) + 1 \times (0.60 - 0.30) = 7.48 \text{ liters}$$

#### Precautions

- 1 The whole method depends on good sampling so be very careful. Rinse out at least 2 times before drawing the final sample.  
*Accurate timing is almost as important as accurate sampling.*
- 2 It will be noted that the two values for the walk should and do check quite closely and similarly for the last three minutes of the run. Generally a man must run at least two minutes at the same grade of work before the maximum is obtained.
- 3 We have found that the oxygen debt for 15 minutes is only roughly proportional to the total oxygen debt. For precise measurement the period should be 60 minutes or until the O<sub>2</sub> consumption approaches the basal level.

- 5 Gasometer readings are taken every half minute during run
- 6 Expired air samples are taken over an entire minute for the second third fourth and fifth minutes (the first minute is not of sufficient value), or for as long as the subject runs. Samples of less than 20 seconds are of doubtful value

#### Procedure—Recovery After Exhausting Exercise

- 1 When subject stops running by jumping to a straddling position (keeping the mouthpiece in and the nose clip on) the following must be done precisely at the time he stops running
  - a Valve at treadmill superstructure is turned to collect expired air in Douglas bags for three minutes
  - b Cardi tachometer tape is marked.
  - c Time is noted.
  - d Treadmill is stopped. Subject sits on tall stool in place on treadmill still breathing through mouthpiece
- 2 As soon as possible the following additional steps should be taken
  - a. Final gasometer reading is taken.
  - b Temperature of expired air in gasometer is noted
  - c Gasometer is emptied. Connection removed from mixing chamber and placed on direct connection to gasometer. Valve turned for collection, and reading noted
- 3 Blood pressure is taken as soon as possible preferably during the first minute after subject stops running
- 4 At the end of three minutes, valve at treadmill superstructure is turned to gasometer
- 5 Gasometer readings are taken every minute during the remainder of the 15 minutes.
- 6 Finger blood samples for blood sugar and lactate are taken at 5 10 and 15 minutes during recovery
- 7 Cardi tachometer is left running at least for the first five minutes of recovery and may be turned off then and used only to take heart rate during 10th and 15th minutes of recovery
- 8 At 15 minutes gasometer valve is closed temperature is recorded and final reading is noted.
- 9 Two expired air samples are taken from gasometer (3 15 minutes of recovery)
- 10 After gasometer is emptied valves are turned connecting Douglas bag to gasometer and initial reading is noted. Two samples of contents of bag are drawn from small side tube on bag as gas is forced to large gasometer for measurement of volume. Bag is finally rolled to expel all air connection clamped and gasometer volume and temperature are noted. (Three minutes of recovery)
- 11 Barometric pressure is recorded.  
(The work experiment as carried out above required at least three and preferably four men to guarantee smooth operation.)

#### Calculations and Data Obtained from Work Experiments

- 1 Blood pressure at rest (sitting) and immediately following exhausting run.
- 2 Continuous heart rate record counted to obtain following
  - a. Rest (sitting)
  - b Maximum during walk (sixth minute)
  - c Maximum during run (usually last minute)
  - d. First five minutes of recovery. May be counted by six second intervals for first minute of recovery and longer intervals up to five minutes
  - e At 10 minutes during recovery
  - f At 15 minutes during recovery

- 3 Blood sugar and blood lactate in mg per 100 ml
  - a During walk.
  - b At 5, 10 and 15 minutes during recovery
- 4 Metabolism during walk, run and recovery are calculated as described in previous sections
- 5 Oxygen debt. This measurement is the difference in oxygen consumption for fifteen minutes following the run and the basal requirement for that period. To be complete the period should be 60 minutes

#### Example of a Protocol of a Typical Experiment

| SUBJECT W<br>Grade 860'            | DATE<br>8-3-49                 | BAROMETRIC PRESSURE 761 mm Hg<br>GASOMETER TEMP 23.4°C |                              |         | THERMISTOR FACTOR 5<br>× CM STP FACTOR 0.896<br>LOBBIN FACTOR 1.05 |      |           |                                  |
|------------------------------------|--------------------------------|--|------------------------------|---------|--|------|-----------|----------------------------------|
|                                    |                                | DIFF<br>CM   | PULMONARY<br>VENTIL<br>L MIN | CO<br>% | O  | RQ   | TRLE<br>% | O <sub>2</sub><br>L/STP<br>L/STP |
| WALK<br>3.5 mph                    | SP. P. METER<br>READINGS<br>CM |  |                              |         |  |      |           |                                  |
| 4 mins                             | 40                             |  |                              |         |  |      |           |                                  |
| 5                                  | 13.1                           | 8.0  | 39                           | 4.59    | 16.6   | 0.97 | 4.0       | 1.87                             |
| 6                                  | 10                             | 0  | 3.9                          | 4.1     | 15.87  | 0.91 | 5.17      | 1.91                             |
| Run 7.0 mph.                       |                                |  |                              |         |  |      |           |                                  |
| 0 mins                             | 4.0                            |  |                              |         |  |      |           |                                  |
| 1                                  | 16.9                           |  |                              |         |  |      |           |                                  |
| 1½                                 | 25.6                           |  |                              |         |  |      |           |                                  |
|                                    | 34.0                           | 17.1   | 80.0                         | 4.0     | 16.49  | 1.06 | 4.40      | 3.5                              |
| ½                                  | 44.0                           |  |                              |         |  |      |           |                                  |
| 3                                  | 5.4                            | 19.4   | 93.8                         | 4.0     | 16.60  | 1.14 | 4.16      | 3.78                             |
| 3½                                 | 63.1                           |  |                              |         |  |      |           |                                  |
| 4                                  | 3.1                            | 19   | 9                            | 4.89    | 16.6   | 1.17 | 4.17      | 3.84                             |
| 4½                                 | 87.0                           |  |                              |         |  |      |           |                                  |
| 5                                  | 8.5                            | 19.4   | 90.8                         | 4.9     | 16.61  | 1.19 | 4.17      | 3.9                              |
| Recovery                           |                                |  |                              |         |  |      |           |                                  |
| Bag                                |                                |  |                              |         |  |      |           |                                  |
| 0 min                              | 4.7                            |  |                              |         |  |      |           |                                  |
| 3                                  | 44.5                           | 39.8   | 8.0                          | 4.04    | 18.09  | 1.6  | 5.0       | 1.58                             |
| Gasometer                          |                                |  |                              |         |  |      |           |                                  |
| 3 min                              | 4.8                            |  |                              |         |  |      |           |                                  |
|                                    | 16.5                           |  |                              |         |  |      |           |                                  |
| 7                                  | 24.7                           |  |                              |         |  |      |           |                                  |
| 9                                  | 5                              |  |                              |         |  |      |           |                                  |
| 11                                 | 39.8                           |  |                              |         |  |      |           |                                  |
| 13                                 | 45.9                           |  |                              |         |  |      |           |                                  |
| 1                                  | 51.6                           | 46.8   | 19.3                         | 3.3     | 1.9  | 1.08 | 3.30      | 0.60                             |
| Total O <sub>2</sub> consumption — |                                |  |                              |         |  |      |           | 0.30                             |

Net O<sub>2</sub> debt in 15 minutes

$$3 \times (1.8 - 0.30) + 12 (0.60 - 0.30) = 4.8 \text{ l/STP}$$

#### Precautions

- 1 The whole method depends on good sampling so be very careful. Puncture at least 3 times before drawing the final sample. Accurate timing is almost as important as accurate sampling.
- 2 It will be noted that the two values for the walk should and do check quite closely and similarly for the last three minutes of the run. Generally a man must run at least two minutes at this grade of work before the maximum is obtained.
- 3 We have found that the oxygen debt for 15 minutes is only roughly proportional to the total oxygen debt. For precise measurement the period should be 60 minutes or until the O<sub>2</sub> consumption approaches the basal level.

- 5 Gasometer readings are taken every half minute during run
- 6 Expired air samples are taken over an entire minute for the second third, fourth and fifth minutes (the first minute is not of sufficient value) or for as long as the subject runs. Samples of less than 20 seconds are of doubtful value

#### Procedure—Recovery After Exhausting Exercise

- 1 When subject stops running by jumping to a straddling position (keeping the mouthpiece in and the nose clip on) the following must be done precisely at the time he stops running.
  - a Valve at treadmill superstructure is turned to collect expired air in Douglas bags for three minutes
  - b Cardiotachometer tape is marked
  - c Time is noted
  - d Treadmill is stopped. Subject sits on tall stool in place on treadmill still breathing through mouthpiece
- 2 As soon as possible the following additional steps should be taken
  - a. Final gasometer reading is taken
  - b Temperature of expired air in gasometer is noted
  - c Gasometer is emptied. Connection removed from mixing chamber and placed on direct connection to gasometer. Valve turned for collection and reading noted
- 3 Blood pressure is taken as soon as possible preferably during the first minute after subject stops running
- 4 At the end of three minutes valve at treadmill superstructure is turned to gasometer
- 5 Gasometer readings are taken every minute during the remainder of the 15 minutes
- 6 Finger blood samples for blood sugar and lactate are taken at 5 10 and 15 minutes during recovery
- 7 Cardiotachometer is left running at least for the first five minutes of recovery and may be turned off then and used only to take heart rate during 10th and 15th minutes of recovery
- 8 At 15 minutes gasometer valve is closed temperature is recorded and final reading is noted
- 9 Two expired air samples are taken from gasometer (3 15 minutes of recovery)
- 10 After gasometer is emptied valves are turned connecting Douglas bag to gasometer and initial reading is noted. Two samples of contents of bag are drawn from small side tube on bag as gas is forced to large gasometer for measurement of volume. Bag is finally rolled to expel all air connection clamped and gasometer volume and temperature are noted. (Three minutes of recovery)
- 11 Barometric pressure is recorded  
(The work experiment as carried out above required at least three and preferably four men to guarantee smooth operation)

#### Calculations and Data Obtained from Work Experiments

1. Blood pressure at rest (sitting) and immediately following exhausting run
- 2 Continuous heart rate record counted to obtain following
  - a. Rest (sitting)
  - b Maximum during walk (sixth minute)
  - c Maximum during run (usually last minute)
  - d First five minutes of recovery. May be counted by six second intervals for first minute of recovery and longer intervals up to five minutes
  - e At 10 minutes during recovery
  - f At 15 minutes during recovery

- 3 Blood sugar and blood lactate in mg per 100 ml.
  - a. During walk
  - b. At 5, 10 and 15 minutes during recovery
- 4 Metabolism during walk, run and recovery are calculated as described in previous sections
- 5 Oxygen debt. This measurement is the difference in oxygen consumption for fifteen minutes following the run and the basal requirement for that period. To be complete the period should be 60 minutes

## Example of a Protocol of a Typical Experiment

| SUBJECT B<br>Grade 86% | DATE<br>8/3/41               | BAR. PRESSURE 61 mm Hg<br>GASOMETER TEMP 34 C |                 |         | TISSOT FACTOR 5.0<br>x CM SIP FACTOR 0.826<br>COMBINED FACTOR 4.68 |      |         |   |  |
|------------------------|------------------------------|---|-----------------|---------|--|------|---------|---|--|
|                        |                              | 1 F F   | 1 L L           | 1 A A V |  |      |         |   |  |
| WALK<br>3.5 m.p.h.     | SPIROMETER<br>READINGS<br>ml | cm  | VENTIL<br>L/MIN | CO<br>% | O <sub>2</sub>   | RQ   | TF<br>% | O<br>USED<br>L/MIN                      |  |
| 4 mins                 | 46                           |   |                 |         |  |      |         |   |  |
| 5                      | 131                          | 85  | 91              | 4.9     | 16.6   | 0.91 | 4.0     | 1.87                                    |  |
| 6                      | 210                          | 9   | 30.9            | 4.1     | 15.97  | 0.91 | 5.1     | 1.91                                    |  |
| Run 7.0 m.p.h.         |                              |   |                 |         |  |      |         |   |  |
| 0 mins.                | 45                           |   |                 |         |  |      |         |   |  |
| 1                      | 169                          |   |                 |         |  |      |         |   |  |
| 1½                     | 56                           |   |                 |         |  |      |         |   |  |
| 2                      | 340                          | 1.1   | 80.0            | 4.0     | 16.48  | 1.06 | 4.40    | 3.40                                    |  |
| 2½                     | 410                          |   |                 |         |  |      |         |   |  |
| 3                      | 34                           | 19.4  | 90.8            | 4.9     | 16.63  | 1.14 | 4.16    | 3.78                                    |  |
| 3½                     | 631                          |   |                 |         |  |      |         |   |  |
| 4                      | 731                          | 19  | 9               | 4.80    | 16.6   | 1.17 | 4.1     | 3.84                                    |  |
| 4½                     | 830                          |   |                 |         |  |      |         |   |  |
| 5                      | 85                           | 19.4  | 90.8            | 4.97    | 16.61  | 1.19 | 4.1     | 3.9                                     |  |
| Recovery               |                              |   |                 |         |  |      |         |   |  |
| Bag                    |                              |   |                 |         |  |      |         |   |  |
| 0 min                  | 47                           |   |                 |         |  |      |         |   |  |
| 3                      | 44.5                         | 19.8  | 90              | 4.04    | 18.09  | 1.0  | 2.5     | 1.9                                     |  |
| Gasometric             |                              |   |                 |         |  |      |         |   |  |
| 3 mins                 | 48                           |   |                 |         |  |      |         |   |  |
| 5                      | 165                          |   |                 |         |  |      |         |   |  |
| 7                      | 217                          |   |                 |         |  |      |         |   |  |
| 9                      | 25                           |   |                 |         |  |      |         |   |  |
| 11                     | 398                          |   |                 |         |  |      |         |   |  |
| 13                     | 408                          |   |                 |         |  |      |         |   |  |
| 15                     | 516                          | 46.8  | 19.1            | 3.53    | 1.59   | 1.06 | 3.0     | 0.60                                    |  |
|                        |                              |   |                 |         |  |      |         | Basal O <sub>2</sub> consumption = 0.30 |  |

Net O<sub>2</sub> debt in 15 minutes

$$3 \times (1.58 - 0.9) + 12 \times (0.60 - 0.30) = 7.48 \text{ liters}$$

## Precautions

- 1 The whole method depends on good sampling so be very careful. Run a test at least 3 times before drawing the final sample.
- 2 Accurate timing is almost as important as accurate sampling.
- 3 It will be noted that the two values for the walk should and do check quite closely and similarly for the last three minutes of the run. Generally a man must run at least two minutes at this grade of work before the maximum is obtained.
- 4 We have found that the oxygen debt for 15 minutes is only roughly proportional to the total oxygen debt. For precise measurement the period should be 60 minutes or until the O<sub>2</sub> consumption approaches the basal level.



- 5 Gasometer readings are taken every half minute during run
- 6 Expired air samples are taken over an entire minute for the second, third fourth and fifth minutes (the first minute is not of sufficient value) or for as long as the subject runs. Samples of less than 20 seconds are of doubtful value

#### **Procedure—Recovery After Exhausting Exercise**

- 1 When subject stops running by jumping to a straddling position (keeping the mouthpiece in and the nose clip on), the following must be done precisely at the time he stops running
  - a Valve at treadmill superstructure is turned to collect expired air in Douglas bags for three minutes
  - b Cardiotachometer tape is marked
  - c Time is noted
  - d Treadmill is stopped. Subject sits on tall stool in place on treadmill still breathing through mouthpiece
- 2 As soon as possible the following additional steps should be taken
  - a Final gasometer reading is taken
  - b Temperature of expired air in gasometer is noted
  - c Gasometer is emptied. Connection removed from mixing chamber and placed on direct connection to gasometer. Valve turned for collection, and reading noted
- 3 Blood pressure is taken as soon as possible preferably during the first minute after subject stops running
- 4 At the end of three minutes valve at treadmill superstructure is turned to gasometer
- 5 Gasometer readings are taken every minute during the remainder of the 15 minutes
- 6 Finger blood samples for blood sugar and lactate are taken at 5, 10 and 15 minutes during recovery
- 7 Cardiotachometer is left running at least for the first five minutes of recovery and may be turned off then and used only to take heart rate during 10th and 15th minutes of recovery
- 8 At 15 minutes gasometer valve is closed temperature is recorded and final reading is noted
- 9 Two expired air samples are taken from gasometer (3-15 minutes of recovery)
- 10 After gasometer is emptied valves are turned connecting Douglas bag to gasometer and initial reading is noted. Two samples of contents of bag are drawn from small side tube on bag as gas is forced to large gasometer for measurement of volume. Bag is finally rolled to expel all air connection clamped and gasometer volume and temperature are noted. (Three minutes of recovery)
- 11 Barometric pressure is recorded  
(The work experiment as carried out above required at least three and preferably four men to guarantee smooth operation)

#### **Calculations and Data Obtained from Work Experiments**

- 1 Blood pressure at rest (sitting) and immediately following exhausting run
- 2 Continuous heart rate record counted to obtain following
  - a Rest (sitting)
  - b Maximum during walk (sixth minute)
  - c Maximum during run (usually last minute)
  - d First five minutes of recovery. May be counted by six second intervals for first minute of recovery and longer intervals up to five minutes
  - e At 10 minutes during recovery
  - f At 15 minutes during recovery

- 3 Blood sugar and blood lactate in mg per 100 ml
  - a During walk
  - b At 5, 10 and 15 minutes during recovery
- 4 Metabolism during walk, run and recovery are calculated as described in previous sections
- 5 Oxygen debt. This measurement is the difference in oxygen consumption for fifteen minutes following the run and the basal requirement for that period. To be complete the period should be 60 minutes.

#### Example of a Protocol of a Typical Experiment

| SUBJECT B<br>Grasle 86%                 | DATE<br>8/3/49                | BAR PRESSURE 101 mm Hg<br>GAS METER TEMP 34 C |                           |                      | TISSOT FACTOR 5.2<br>× cm STP FACTOR 0.896<br>COMBINED FACTOR 4.68 |      |           |                                 |
|---|-------------------------------|---|---------------------------|----------------------|--|------|-----------|---------------------------------|
|   |                               | DIFF<br>cm                                    | PULMONARY<br>CAP<br>L MIN | CO <sub>2</sub><br>V | O <sub>2</sub>   | PQ   | TRCE<br>" | O <sub>2</sub><br>USED<br>l/min |
| WALK<br>3.5 mph                         | SPIRO METER<br>READING—<br>cm |   |                           |                      |  |      |           |                                 |
| 4 mins                                  | 46                            |   |                           |                      |  |      |           |                                 |
| 5                                       | 131                           | 8   | 39.7                      | 4.59                 | 16.6   | 0.97 | 4.0       | 1.97                            |
| 6                                       | 210                           | 19  | 36.9                      | 4.11                 | 15.97  | 0.91 | 5.17      | 1.91                            |
| Run 7.0 mph.                            |                               |   |                           |                      |  |      |           |                                 |
| 0 mins.                                 | 40                            |   |                           |                      |  |      |           |                                 |
| 1                                       | 169                           |   |                           |                      |  |      |           |                                 |
| 1½                                      | 36                            |   |                           |                      |  |      |           |                                 |
|   | 340                           | 11  | 40.0                      | 4.0                  | 16.48  | 1.06 | 4.40      | 3.5                             |
| ½                                       | 440                           |   |                           |                      |  |      |           |                                 |
| 3                                       | 534                           | 19.4  | 37.8                      | 4.3                  | 16.60  | 1.14 | 4.16      | 3.8                             |
| 3½                                      | 631                           |   |                           |                      |  |      |           |                                 |
| 4                                       | 31                            | 19  | 3                         | 4.89                 | 16.67  | 1.17 | 4.1       | 3.84                            |
| 4½                                      | 830                           |   |                           |                      |  |      |           |                                 |
| 5                                       | 805                           | 19.4  | 40.9                      | 4.97                 | 16.61  | 1.19 | 4.1       | 3.79                            |
| Recovery                                |                               |   |                           |                      |  |      |           |                                 |
| Basal                                   |                               |   |                           |                      |  |      |           |                                 |
| 0 min                                   | 47                            |   |                           |                      |  |      |           |                                 |
| 3                                       | 445                           | 9.8   | 6.0                       | 4.04                 | 16.09  | 1.56 | 0.50      | 1.59                            |
| Grasle                                  |                               |   |                           |                      |  |      |           |                                 |
| 3 mins                                  | 48                            |   |                           |                      |  |      |           |                                 |
| 5                                       | 165                           |   |                           |                      |  |      |           |                                 |
| 7                                       | 47                            |   |                           |                      |  |      |           |                                 |
| 9                                       | 375                           |   |                           |                      |  |      |           |                                 |
| 11                                      | 88                            |   |                           |                      |  |      |           |                                 |
| 13                                      | 458                           |   |                           |                      |  |      |           |                                 |
| 15                                      | 216                           | 46.8  | 19.3                      | 3.53                 | 1.9  | 1.06 | 0         | 0.60                            |
| Basal O <sub>2</sub> consumption = 0.30 |                               |   |                           |                      |  |      |           |                                 |

Net O<sub>2</sub> debt in 15 minutes

$$3 \times (1.58 - 0.30) + 1 \times (0.60 - 0.30) = 1.48 \text{ liters}$$

#### Precautions

- 1 The whole method depends on good sampling so be very careful. Rinse out at least 3 times before drawing the final sample.  
*Accurate timing is almost as important as accurate sampling.*
- 2 It will be noted that the two values for the walk should and do check quite closely and similarly for the last three minutes of the run. Generally a man must run at least at two minutes at this grade of work before the maximum is obtained.
- 3 We have found that the oxygen debt for 15 minutes is only roughly proportional to the total oxygen debt. For precise measurement the period should be 60 minutes or until the O<sub>2</sub> consumption approaches the basal level.

- 5 Gasometer readings are taken every half minute during run
- 6 Expired air samples are taken over an entire minute for the second, third fourth and fifth minutes (the first minute is not of sufficient value) or for as long as the subject runs. Samples of less than 20 seconds are of doubtful value

#### **Procedure—Recovery After Exhausting Exercise**

- 1 When subject stops running by jumping to a straddling position (keeping the mouthpiece in and the nose clip on), the following must be done precisely at the time he stops running
  - a Valve at treadmill superstructure is turned to collect expired air in Douglas bags for three minutes
  - b Cardiotachometer tape is marked.
  - c Time is noted
  - d Treadmill is stopped. Subject sits on tall stool in place on treadmill still breathing through mouthpiece
- 2 As soon as possible the following additional steps should be taken
  - a Final gasometer reading is taken
  - b Temperature of expired air in gasometer is noted.
  - c Gasometer is emptied. Connection removed from mixing chamber and placed on direct connection to gasometer. Valve turned for collection and reading noted
- 3 Blood pressure is taken as soon as possible, preferably during the first minute after subject stops running
- 4 At the end of three minutes valve at treadmill superstructure is turned to gasometer
- 5 Gasometer readings are taken every minute during the remainder of the 15 minutes
- 6 Finger blood samples for blood sugar and lactate are taken at 5, 10 and 15 minutes during recovery
- 7 Cardiotachometer is left running, at least for the first five minutes of recovery and may be turned off then and used only to take heart rate during 10th and 15th minutes of recovery
- 8 At 15 minutes gasometer valve is closed, temperature is recorded and final reading is noted
- 9 Two expired air samples are taken from gasometer (3, 15 minutes of recovery)
- 10 After gasometer is emptied valves are turned connecting Douglas bag to gasometer and initial reading is noted. Two samples of contents of bag are drawn from small side tube on bag as gas is forced to large gasometer for measurement of volume. Bag is finally rolled to expel all air, connection clamped and gasometer volume and temperature are noted. (Three minutes of recovery)
- 11 Barometric pressure is recorded  
(The work experiment as carried out above required at least three and preferably four men to guarantee smooth operation.)

#### **Calculations and Data Obtained from Work Experiments**

- 1 Blood pressure at rest (sitting) and immediately following exhausting run
- 2 Continuous heart rate record counted to obtain following
  - a Rest (sitting)
  - b Maximum during walk (sixth minute)
  - c Maximum during run (usually last minute)
  - d First five minutes of recovery. May be counted by six second intervals for first minute of recovery and longer intervals up to five minutes
  - e At 10 minutes during recovery
  - f At 15 minutes during recovery

- 3 Blood sugar and blood lactate in mg per 100 ml
  - a During walk
  - b At 5, 10 and 15 minutes during recovery
- 4 Metabolism during walk, run and recovery are calculated as described in previous sections
- 5 Oxygen debt. This measurement is the difference in oxygen consumption for fifteen minutes following the run and the basal requirement for that period. To be complete the period should be 60 minutes.

## Example of a Protocol of a Typical Experiment

| SUBJECT B<br>Grade 86 <sup>00</sup> | DATE<br>8/3/49               | BAR PRESSURE 761 mm Hg<br>GASOMETER TEMP 34 C |                                   |         | TISSUE FACTOR 0<br>× cm STP FACTOR 0.898<br>CORRECTED FACTOR 4.69 |      |           |                                 |
|-------------------------------------|------------------------------|---|-----------------------------------|---------|---|------|-----------|---------------------------------|
| WALK<br>3.5 mph                     | SPISOMETER<br>READINGS<br>cm | DIFF<br>cm                                    | PULMONARY<br>VENTILATION<br>L/MIN | CO<br>% | O <sub>2</sub>  | RQ   | TRUF<br>L | O <sub>2</sub><br>USED<br>L/MIN |
| 4 mins                              | 46                           |   |                                   |         |   |      |           |                                 |
| 5                                   | 131                          | 8.5   | 39                                | 4.59    | 16.6  | 0.97 | 4.70      | 1.8                             |
| 6                                   | 10                           | 9   | 74.9                              | 4.1     | 15.87   | 0.91 | 5.17      | 1.91                            |
| Run, 6.0 mph.                       |                              |   |                                   |         |   |      |           |                                 |
| 0 mins                              | 45                           |   |                                   |         |   |      |           |                                 |
| 1                                   | 169                          |   |                                   |         |   |      |           |                                 |
| 1½                                  | 56                           |   |                                   |         |   |      |           |                                 |
| 2                                   | 340                          | 1.1   | 80.0                              | 4.0     | 16.48   | 1.06 | 4.40      | 3.5                             |
| 2½                                  | 440                          |   |                                   |         |   |      |           |                                 |
| 3                                   | 534                          | 19.4  | 90.8                              | 4.9     | 16.65   | 1.14 | 4.16      | 3.78                            |
| 3½                                  | 631                          |   |                                   |         |   |      |           |                                 |
| 4                                   | 731                          | 19  | 9                                 | 4.89    | 16.6  | 1.17 | 4.17      | 3.84                            |
| 4½                                  | 830                          |   |                                   |         |   |      |           |                                 |
| 5                                   | 855                          | 18.4  | 90.8                              | 4.9     | 16.61   | 1.19 | 4.17      | 3.9                             |
| Recovery                            |                              |   |                                   |         |   |      |           |                                 |
| Bag                                 |                              |   |                                   |         |   |      |           |                                 |
| 0 min                               | 47                           |   |                                   |         |   |      |           |                                 |
| 3                                   | 445                          | 39.9  | 6.0                               | 4.04    | 18.08   | 1.50 | 2.5       | 1.58                            |
| Gasometer                           |                              |   |                                   |         |   |      |           |                                 |
| 3 mins                              | 48                           |   |                                   |         |   |      |           |                                 |
| 5                                   | 185                          |   |                                   |         |   |      |           |                                 |
| 7                                   | 147                          |   |                                   |         |   |      |           |                                 |
| 9                                   | 325                          |   |                                   |         |   |      |           |                                 |
| 11                                  | 388                          |   |                                   |         |   |      |           |                                 |
| 13                                  | 458                          |   |                                   |         |   |      |           |                                 |
| 15                                  | 516                          | 46.8  | 18.3                              | 3.3     | 1.59  | 1.06 | 3.30      | 0.60                            |
| Total O <sub>2</sub> consumption =  |                              |   |                                   |         |   |      |           | 0.30                            |

Net O<sub>2</sub> debt in 15 minutes

$$3 \times (1.58 - 0.30) + 1 \times (0.60 - 0.30) = 4.8 \text{ litres}$$

## Precautions

- 1 The whole method depends on good sampling so be very careful. Put out at least 3 times before drawing the final sample.
- 2 Accurate timing is almost as important as accurate sampling.
- 3 It will be noted that the two values for the walk should and do check quite closely and similarly for the last three minutes of the run. Generally a man must run at least two minutes at the grade of work before the maximum is obtained.
- 4 We have found that the oxygen debt for 15 minutes is only roughly proportional to the total oxygen debt. For precise measurement the period should be 60 minutes or until the O<sub>2</sub> consumption approaches the basal level.

## SECTION VI

# PHYSIOLOGICAL MEASUREMENTS (Continued)

---

### E PULMONARY MEASUREMENTS

#### 1 Alveolar Air in Rest and Work

##### A. IN REST

###### Reference

Dill D B, Hurxthal L. M., Van Caulaert, C, Folling A, and Bock, A. V. The Carbon Dioxide Equilibrium in Alveolar Air and Arterial Blood. J Biol. Chem 74 303-313 (Aug) 1927

###### Principle

Alveolar air is that portion of the air in the lungs that is in contact with the alveolar walls. According to the concept of Haldane, the last portion of expired air at the end of a forced expiration may be taken as representative of the alveolar air of the whole lung.

###### Apparatus


- 1 Alveolar air sampler (see line drawing Fig 50). This is a modification of the Haldane Priestley tube with a three way stopcock between the rubber tube and the mouthpiece. The smaller the dead space in the valve the better. The rubber tube must be at least 3 feet long.
- 2 A number of gas sampling tubes (Barcroft type). The tube completely filled with mercury including the bore above the stopcock is attached to a sampling outlet on the part of the cock carrying the tube.
- 3 A nose clip
- 4 A clean mouthpiece
- 5 A Haldane gas analysis apparatus

###### Reagents

- 1 Clean mercury
- 2 Reagents for the Haldane gas analyzer

###### Procedure

- 1 The mouthpiece and the nose clip are adjusted on the subject after he has rested for at least one half hour.
- 2 The cock has been turned so that the subject is breathing directly into the room.
- 3 Towards the end of a normal expiration the valve is turned into the rubber tube and the signal 'Blow' is given.
- 4 The subject expires forcibly and completely without inspiring after the signal. Blow is given. The tap is immediately turned back so that the subject may breathe normally. This maneuver leaves a sample of alveolar air in the tube.
- 5 A sample is withdrawn from the tube into a gas sampling tube.
- 6 The sample is then analyzed in the Haldane apparatus for carbon dioxide alone.




## ALVEOLAR AIR AT REST

C. E. SORDON

Fig 50—Arrangement of valves and sampling tubes for collecting alveolar air by the Haldane-Priestley technique

# PARTIAL PRESSURE OF A GAS FROM BAROMETRIC PRESSURE AND PERCENTAGE OF GAS FOUND



— G. E. GORDON —

Fig. 51.—Nomogram for estimating partial pressure of a gas from barometric pressure and percentage of gas found. Lay ruler to connect the two outside scales at the proper points for total barometric pressure (corrected for tension of aqueous vapor), and percentage of gas found. The partial pressure (mm Hg) of the gas will then be indicated at the point where the ruler crosses the center scale.


Calculation

Alveolar  $p_{CO}$  in mm Hg =

$$(\% CO \text{ from Haldane analysis}) \times \frac{(B - W)}{100}$$

Where B is barometric pressure in mm Hg and W is vapor pressure of water at 37 C  
i.e. 47 mm Hg

# ALVEOLAR AT WORK





## 350 *Physiological Measurements*

- 2 The line chart Fig 51 gives directly the partial pressure of a gas from the percentage and barometric pressure *corrected for vapor tension*. In the case of  $\text{CO}_2$  it is necessary to use some factor such as 10 to enable the right hand scale to be used. The reading in the middle scale must be corrected accordingly.

### Precautions

- 1 Draw the sample as soon as possible after the expiration
- 2 Be sure to give the signal *Blow* only after a normal expiration.
- 3 Subjects have to be trained in this procedure to give reliable results

## B IN WORK

### Reference

Henderson Y and Haggard, H W The Circulation and Its Measurement *Am J Physiol* 73 193 253 (June) 1905

### Principle

During work respiration is not under the same easy voluntary control as it is in rest. Hence the procedure for collecting alveolar air in rest cannot be used in work. A device of Y Henderson and associates collects samples automatically. It utilizes the slight negative pressure which is produced during inspiration to suck samples of the last of the preceding expiration into a collecting bottle.

### Apparatus

- 1 Y Henderson's apparatus for collecting alveolar air, during exercise (see Fig 50)
- 2 Gas sampling tubes (50 ml volume) filled with mercury
- 3 Other apparatus as for work metabolism

### Reagents

- 1 Clean mercury
- 2 Reagents for Haldane gas analysis
- 3 Water acidulated slightly with sulfuric acid for filling the apparatus

### Procedure

- 1 The mouthpiece and nose clip are adjusted on the subject when the metabolic collection period begins
- 2 All clamps on the alveolar air sampler are closed
- 3 When it is desired to collect a sample of alveolar air the steps are
  - a Attach a gas sampling tube to the rubber outlet
  - b Open clamps 1 2 and 4.
  - c The depth of the tube in the acidulated water is determined by the depth of the subject's respiration. It is usually about 2.5 mm
  - d The air is allowed to flush out the chamber for at least 2 minutes
  - e Clamps 1 and 4 are closed.
  - f Clamp 3 is opened
  - g A sample is obtained by opening the cock of the sampling tube
  - h Analysis of carbon dioxide is performed in the Haldane gas analysis apparatus.

### Calculation

Same as for alveolar air in rest

### Example

Same as for alveolar air in rest

### Precautions

Same as for alveolar air in rest and in addition

- 1 Acidulation of the water in the collecting apparatus is necessary to prevent absorption of  $\text{CO}_2$

## Functional Residual Air

### Reference

- Cournaud A Baldwin E DeF Darling R C and Richards D W Jr Studies on Intrapulmonary Mixtures of Gases IV The Significance of the Pulmonary Emptying Rate and a Simplified Open Circuit Measurement of Residual Air J Clin Investigation 20 681 689 (Oct) 1941

### Principle

Practically all the nitrogen in the lungs is washed out by seven minutes of breathing oxygen. The total expired air is collected, measured and analyzed for  $N_2$ . The volume of the gas space in the lungs is calculated from these data. An appropriate correction is made for nitrogen freed from solution in the blood.

### Apparatus

- 1 A 100 liter Tissot gasometer is used with connections as shown in Fig 53 in which the symbols mean the following

M—mouthpiece

Alv—Evacuated gas sampling tubes

V V V V—Three way respiratory valves

F F F F—One way rubber flutter valves

D<sub>r</sub>—Rubber anesthesia bag

T—One hundred liter (Tissot) gasometer

B—Valve and attachment for obtaining gasometer samples

- a The main circuit is essentially that used for open circuit measurement of basal metabolism except that the inlet arm comes from a reservoir bag B which is connected to an oxygen tank. The valve V is an essential feature but may be improvised in various ways. It must be near to the mouthpiece M so that there is negligible dead space but if the subject's breathing can be shifted from room air to oxygen in the room air or side circuit there should be any suitable alveolar sampling device (Pictured in the diagram as evacuated sampling tubes inserted through holes in the rubber tubing connection). The recording pen and drum D are useful but unessential parts of the apparatus.

- b The dead space under the bell of the gasometer T should be measured before using the apparatus. A satisfactory accuracy within 90 ml may be attained by geometrical calculation of the space equally easy and more accurate measurement can be made by adding a known volume of oxygen to the dead space volume of room air through the thermometer socket (not shown on diagram) and analyzing the mixture in the bell.

- c If V = vol O added DS = dead space under bell a = % N in oxygen supply b = % N in mixture under bell

$$DS = \frac{V(b - a)}{a - b}$$

A Van Slyke manometric apparatus should be available for analysis of the high oxygen mixtures. A Scholander gas analyzer is equally as good.

### Reagents

- 1 A tank of oxygen
- 2 Reagents for Van Slyke or Scholander apparatus

### Procedure

- 1 With V turned to the side circuit flush the main circuit and gasometer six consecutive times with oxygen. Then turn V to connect main circuit to outside air and maintain an oxygen flow approximately equal to the subject's ventilation.
- 2 Have the subject under basal conditions reclining with a one or two pillow head rest. Insert the mouthpiece and put on the nose clip.

- 3 Take the initial gasometer reading
- 4 Turn V to connect to gasometer
- 5 Observing the subject's respiratory rhythm and being sure that he is breathing easily turn V to the main circuit exactly at the end of a normal expiration
- 6 Have the subject continue breathing oxygen for 7 minutes adjust the oxygen flow as necessary to keep the bag B nearly full
- 7 During this time prepare the side circuit for alveolar sampling
- 8 At the beginning of the expiration closest to the seven minute point turn the valve V to the side circuit and instruct the subject to exhale forcibly Take the alveolar sample The subject may then be disconnected

### OPEN CIRCUIT APPARATUS FOR MEASUREMENT OF RESIDUAL AIR




Fig. 53—The apparatus used by Cournand, Baldwin, Darling and Richards for estimation of functional residual air.

- 9 Flush the main circuit with approximately 10 liters of oxygen and close valve 1.
- 10 Take second gasometer reading and temperature
- 11 Take two samples from the gasometer after first flushing the gasometer tubing with the first portion of the contents
- 12 Repeat steps 1 to 11 for a check determination after  $\frac{1}{2}$  hour interval
- 13 Estimate the nitrogen in the alveolar and spirometer gases using the Van Slyke apparatus or the Scholander gas analyzer

## Calculation

## Symbols

- FPA (dry) = functional residual air in ml (calculated as dry gas)  
 V = volume of spirometer gas corrected to 37°C (dry)  
 DS = volume of dead space under gasometer bell corrected to 37°C  
 a = % N (+ other inert gases) in oxygen supply (usually < 0.5%)  
 b = % N from gasometer analysis  
 Alv N = % N in alveolar sample  
 C = correction for nitrogen from blood (See formula below)  
 BS = body surface of subject in square meters  
 BP = barometric pressure mm Hg

## Equations

$$1 \text{ FPA (dry)} = \frac{(V + DS)(b - a) - C}{81 - \text{Alv N}}$$

2 Substitute value of C from following

$$C \text{ ml} = \frac{(BS \times 96.5) + 3}{6.66}$$

3 For calculation of FRA saturated with water vapor as exists (approximately) within the lungs

$$\text{FRA} = \text{FPA (dry)} \times \frac{\text{BP}}{\text{BP} - 4}$$

## Example

V was 0.21 l DS 0.21 a 0.3% b 8.30% Alv N 3.00% BS 1.75 sq m and BP 760 mm Hg

1 FPA liter dry =

$$\frac{(0.21 + 0.21)(8.30 - 0.32) - (96.5 \times 1.75 + 3)/6.66}{81 - 3}$$

$$= 0.6$$

FRA liters wet =

$$0.6 \times \frac{60}{60 - 47}$$

$$= 2.19$$

## Precautions

This method depends upon

- 1 Accurate sampling
- 2 Accurate timing
- 3 Extremely accurate analysis of samples of gas rich in oxygen

## 3 Cardiac Output (Grollman's Method)

## Reference

- Marshall E K Jr and Grollman A A Method for the Determination of the Circulatory Minute Volume in Man *Am J Physiol* 86 127-137 (Aug) 1953  
 Grollman A and Marshall E K Jr The Time Necessary for Rebreathing in a Lung Bag System to Attain Homogeneous Mixture *Am J Physiol* 86 110-116 (Aug) 1953

## Principle

The subject breathes acetylene which dissolves in plasma but does not combine with hemoglobin lipoids or other constituents of the blood. The quantity absorbed

is measured and from this measurement and the coefficient of solubility of acetylene in plasma the quantity of blood which has passed through the lungs is calculated. The blood flow through the lungs is affected by the rebreathing of acetylene. Hence in order to calculate the true cardiac output is that immediately preceding the breathing. Grollman computes arteriovenous difference from the oxygen absorbed during rebreathing and the blood flow during the rebreathing as determined by acetylene absorption. The total oxygen consumption immediately preceding rebreathing divided by the arteriovenous oxygen difference gives the cardiac output prior to rebreathing. In the analysis of the gas samples three gases are absorbed and a fourth is obtained by difference. The analysis is performed in a machine similar in principle to that of the Haldane previously described. There are two machines in common use the Grollman Haldane and the Krogh Haldane.

#### Apparatus

- 1 A Grollman Haldane apparatus. This apparatus has an extra absorption tube for the analysis of acetylene.
- 2 A rubber bag connected to a 2 way valve.
- 3 A rack of gas sampling tubes 30 ml volume.
- 4 Apparatus for determining metabolic rate.

#### Reagents

- 1 A tank of oxygen gas.
- 2 A tank of air.
- 3 A tank of acetylene gas.
- 4 Same reagents as for the Haldane plus an acetylene absorber. 20 gm of mercuric cyanide and 8 gm of sodium hydroxide are dissolved and made up to 100 ml of water.

#### Procedure—Estimation of Total Oxygen Consumption

- 1 Immediately preceding the rebreathing of acetylene estimate the oxygen consumption of the subject as described under metabolic measurements.

#### Procedure—Rebreathing of Acetylene

- 1 Depending on the subject's activity fill the bag as follows

|          | Oxygen     | Air        | Acetylene  |
|----------|------------|------------|------------|
| Rest     | 0          | 900 0 ml   | 500 700 ml |
| Activity | 500 700 ml | 100 500 ml | 500 700 ml |

- 2 It is essential that the bag mixture maintain near saturation of the arterial blood throughout the period of rebreathing.
- 3 Prepare a series of evacuated gas sampling tubes.
- 4 The mixture of O<sub>2</sub>, N<sub>2</sub> and acetylene is rebreathed from a bag equipped with a 2 way valve and also equipped with two small outlets between bag and valve for taking gas samples.
- 5 The subject with mouthpiece and nose clip in place is instructed to exhale completely, at which moment the valve is turned to connect him to the bag contents.
- 6 He is then instructed to rebreathe the bag contents rapidly and completely but without undue effort. The speed should be set so that 5-6 breaths are complete before sample I is taken. (See below).
- 7 It is essential to complete the rebreathing within one circulation of the blood. At rest total rebreathing time must be not more than 40 seconds and in moderate activity not more than 15 seconds.
- 8 With a resting subject collect Sample I after 10 seconds of rebreathing and Sample II after 18-20 seconds. In work take Sample I after 7 seconds and Sample II after 12 seconds.

## Procedure—Analysis of Samples

- 1 The procedure for analysis is essentially that of the Haldane machine with a few differences
  - a After the initial total volume reading the gases are passed into the KOH exactly 8 times no more and a single unchecked reading is obtained. The reason for this is that acetylene is slightly soluble in KOH
  - b The acetylene is then absorbed. Due to the formation of the white precipitate it is difficult to read the meniscus mark. This difficulty is avoided by having a series of marks on the tube and by using these successively beginning with the bottom and ascending as the mark becomes indistinct
  - c The oxygen is absorbed as in ordinary gas analysis.

## Calculations

- 1 Volume of lung bag system is reduced during rebreathing since the sum of the O and  $C_2H_2$  absorbed is greater than the CO eliminated. N takes no part in the reactions the increase in its percentage is proportional to the decrease in lung bag volume. The relative gas volumes at the times of Sample I and Sample II are then

$$V = V_1 \times \frac{\%N}{\%N_1} \quad (\text{Equation 1})$$

- 2 The amount of  $C_2H_2$  absorbed is proportional to the mean  $C_2H_2$  tension during rebreathing and to the solubility coefficient of  $C_2H_2$  in blood

$$\%C_2H_2 \text{ Av} = \frac{\%C_2H_2 + \%C_2H_2}{2} \quad (\text{Equation 2})$$

and ml  $C_2H_2$  absorbed/liter blood during interval between Sample I and Sample

$$II = \frac{740 \times \%C_2H_2 \times (B - 47)}{100 \times 100} \quad (\text{Equation 3})$$

where 740 is ml  $C_2H_2$  which will be dissolved by blood at 760 mm Hg  $C_2H_2$  pressure at 37 C. B is observed barometric pressure and 47 is vapor tension of water at 37 C.

- 3 The blood flow through the lungs is determined by the ratio between total  $C_2H_2$  absorbed during period of observation and the quantity of  $C_2H_2$  absorbed per liter of blood. Liters blood flow through lungs during rebreathing (i.e. liters cardiac output during rebreathing) =

$$\frac{(V_1 \times \%C_2H_2 - V \times \%C_2H_2) \times (60 \times 100)}{(740 \times \%C_2H_2) \times (B - 47)} \quad (\text{Equation 4})$$

- 4 To obtain cardiac output prior to rebreathing it is now necessary to calculate the arterio-venous oxygen difference during rebreathing the quantity of O absorbed by each liter of blood in its passage through the lungs

A - V oxygen difference in ml =

$$\frac{(V_1 \times \%O_1 - V \times \%O_2) \times (740 \times \%C_2H_2) \times (B - 47)}{(V_1 \times \%C_2H_2 - V \times \%C_2H_2) \times 760 \times 100} \quad (\text{Equation 5})$$

Combining the factors 740 760 and 100

A - V Difference =

$$\frac{(V_1 \times \%O_1 - V \times \%O_2) \times (B - 47) \times 0.00974 \times \%C_2H_2}{(V_1 \times \%C_2H_2 - V \times \%C_2H_2)}$$

Substituting values from Equation 1

A - V difference =

$$\frac{(V_1 \times \%O_1 - V_1 \times \%O_2 \times \frac{\%N}{\%N_1}) \times (B - 47) \times 0.00974 \times \%C_2H_2}{(V_1 \times \%C_2H_2 - V_1 \times \%C_2H_2 \times \frac{\%N}{\%N_1})}$$

And  $V_1$  cancels out leaving

A - V difference =

$$\frac{(\%O_1 - \%O_{II} \times \frac{\%N_I}{\%N_{II}}) \times \%C_2H_{2A} \times (B - 47) \times 0.0097}{(\%C_2H_{2I} - \%C_2H_{2II} \times \frac{\%N_I}{\%N_{II}})} \quad (\text{Equation 6})$$

5 Cardiac output, l/min prior to rebreathing =

$$\frac{\text{Total oxygen consumption ml/min}}{\text{A - V oxygen difference ml/l blood}} \quad (\text{Equation 7})$$

6 Stroke volume ml per heart beat =

$$\frac{\text{Cardiac output l/min}}{\text{Pulse rate beats/min}} \times 1000 \quad (\text{Equation 8})$$

#### Precautions

- 1 This is a difficult method. Much practice is required on the part of both the observer and the subject to obtain consistent and reproducible results
- 2 The calculations above derive step by step the final equations which are used in routine calculation Equations 6, 7 and 8

## 4 Pulmonary and Circulatory Measurements Miscellaneous References

### Acetylene

Coulson Smith C and Seyfang A P A Colorimetric Method for the Estimation of Small Quantities of Acetylene in Air, *Analyst* 67 39 41 (Feb) 1942

### Ballistocardiograph

Nickerson J L Symposium on Cardiac Output The Low Frequency Critically Damped Ballistocardiograph Federation Proc 4 (2) 201 '06 (June) 1945

Nickerson J L Some Observations on the Ballistocardiographic Pattern With Special Reference to the H and K Waves *J Clin. Investigation* 28 369 377 (March) 1949

### Cardiac Output

McMichael J and Sharpey Schafer E P Cardiac Output in Man by the Direct Fick Method Effects of Posture Venous Pressure Change Atropine and Adrenaline *Brit Heart J* 6 (1) 33 40 (July) 1944

Ring G C Greisheimer E M Baier H N Oppenheimer M J Sokalchuk A Ellis D and Friday S J Electrokymograph for Estimation of Heart Output Comparison With Direct Fick in Dogs *Am J Physiol* 161 231 35 (May) 1950

Ring G C, Sokalchuk A, Baier H N, Pudel H Oppenheimer M J, Friday S J, and Navis G Electrokymograph for Estimation of Heart Output Comparison With Stewart in Dogs *Am. J Physiol* 161 236 238 (May) 1950

Seely R D and Gregg D E Technique for Measuring Cardiac Output Directly by Cannulation of the Pulmonary Artery *Proc Soc Exper Biol & Med.* 73 960 2,0 (Feb) 1950

Warren J V Stead E A Jr and Brannon E S Cardiac Output in Man A Study of Some of the Errors in the Method of Right Heart Catheterization *Am. J Physiol.* 145 459 464 (Feb) 1948

Werko L Berseus S, and Lagerlof H A Comparison of the Direct Fick and the Grollman Methods for Determination of Cardiac Output in Man *J Clin Investigation* 28 516 520 (May) 1949

### Catheterization

Bing R J Vandam L D, Gregoire F, Handelsman J C Goodale W T and Eickenhoff J E Catheterization of the Coronary Sinus and Middle Cardiac Vein in Man *Proc Soc Exper Biol & Med* 66 (1) 239 240 (Oct.) 1947

Cournand A Measurement of the Cardiac Output in Man Using the Right Heart Catheterization Description of Technique Discussion of Validity and of Place in the Study of Circulation *Fed Proc* 4 207 212 (June) 1945

Cournand A. and Ranges H A Catheterization of the Right Auricle in Man *Proc Soc Exper Biol & Med.* 46 462-466 (March) 1941

Cournand A. Ranges H. A. and Riley R L Comparison of Results of the Normal Ballistocardiogram and a Direct Fick Method in Measuring the Cardiac Output in Man *J Clin. Investigation* 21 (3) 437 493 (May) 1942

- Corrao A. Riley R L Breed, E S Ballwin, E de F and Richards D W Jr  
Measurement of Cardiac Output in Man Using the Technique of Catheterization of the Right Auricle or Ventricle *J Clin Investigation* 24 (1) 106-116 (Jan) 1945
- Dexter L Haynes F W Barwell C S and Eppinger E C Sagereson R P and Evans J M Studies of Congenital Heart Disease II. The Pressure and O<sub>2</sub> Content in the Right Auricle Right Ventricle and Pulmonary Artery in Control Patients With Observations on the O<sub>2</sub> Saturation and Source of Pulmonary Capillary Blood *J Clin Investigation* 26 544-560 (May) 1947
- Dexter L Haynes F W Barwell C D Eppinger E C, Seibel R E and Evans J M Studies of Congenital Heart Disease I. Technique of Venous Catheterization as a Diagnostic Procedure *J Clin. Investigation* 26 547-553 (May) 1947

### Cerebral Circulation

- Kety S S and Schmidt C F The Nitrous Oxide Method for the Quantitative Determination of Cerebral Blood Flow in Man Theory Procedure and Normal Values *J Clin Investigation* 27 476-493 (July) 1948
- Kety S S The Quantitative Determination of Cerebral Blood Flow in Man Methods in Medical Research Chicago Vol I 1948 Yearbook Publishers
- Myerson A. Helloran R D and Hirsch H L. Technique for Obtaining Blood from the Internal Jugular Vein and Internal Carotid Artery *Arch Neurol and Psychiat* 17 (6) 878-893 (June) 1947

### Circulation Time

- Wexler J and Whittenberger J L Objective Method for Determining Circulation Time from Pulmonary to Systemic Capillaries by Use of the Oximeter *J Clin Investigation* 25 (3) 44-450 (May) 1946

### Coronary Blood Flow

- Eckenhoff J E Hafenschied J H Harmel M H Goodale W T Lubin M Bing R J, and Kety S S Measurement of Coronary Blood Flow by the Nitrous Oxide Method *Am J Physiol* 152 (1) 356-364 (Feb) 1946

### Electrocardiograph

- Graybiel A and White P D *Electrocardiography in Practice* ed 2 Philadelphia London 1946 W B Saunders Company
- Katz L N Exercises in Electrocardiographic Interpretation ed Philadelphia 1946 L A and Febiger

### Gas Tension in Blood

- Boos A., and Black H Direct Determination of Partial and Total Tensions of Respiratory Gases in Blood, *Am. J Physiol.* 160 1031-6 (Jan) 1950

### Lung Volume

- Herrald F J C and McMichael J Determination of Lung Volume A Simple Constant Volume Modification of Christie's Method *Proc Roy Soc London* 126 (B345) 431-501 (Feb) 1933
- McMichael J A Rapid Method of Determining Lung Capacity *Clin Sci* 4 (-) 167-173 (Dec) 1933
- Mencely G F and Faltreder N L Use of Helium for Determination of Pulmonary Capacity *Proc Soc Exper Biol & Med* 46 (1) 66-69 (Jan) 1941
- Penney G R and Haltreider N L The Volume of the Lung Determined by Helium Dilution Description of the Method and Comparison With Other Procedures *J Clin Investigation* 28 191-199 (Jan) 1949
- Bateman J B Boothby W M and Helmholz H F Jr Studies of Lung Volumes and Intrapulmonary Mixing Notes on Open Circuit Methods Including Use of a New Pivoted Type Gasometer for Lung Clearance Studies *J Clin Investigation* 28 69-89 (July) 1949

### Pressure Breathing

- Barach A L Fenn W O Ferris E L and Schmidt C F The Physiology of Pressure Breathing A Brief Review of Its Present Status *J Aviation Med* 18 (1) 13-18 (Feb) 1947

### Venous Pressure

- Lyons R H Kennedy J A and Barwell C S Measurement of Venous Pressure by the Direct Method *Am Heart J* 16 (6) 563-573 (Dec) 1938



## SECTION VI

### PHYSIOLOGICAL MEASUREMENTS (Continued)

---

#### F PLASMA VOLUME, BLOOD VOLUME, EXTRACELLULAR VOLUME

##### 1 Plasma and Blood Volumes (Dye Method)

###### Reference

Nitshe G A Jr and Cohen P I Simplified Determination of Blood Volume *Am J Clin Path* 17 239 243 (March) 1947

###### Principle

The procedure consists of the injection of a known amount of the dye T 1824 (Evans Blue) and the calculation of the plasma volume from the concentration of the dye in the plasma after complete mixing

###### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 10 x 75 mm
- 3 Specially calibrated syringes sterilized for dye injection
- 4 10 ml syringes
- 5 Two 20 gauge needles for venipuncture (sterile)
- 6 Heparin solution 110 unit/mg Dilute 100 mg dry heparin to 147 ml isotonic saline and use one drop for the preservation of 10 ml of whole blood
- 7 0.2 ml syringe pipette accurately calibrated.
- 8 4 ml hematocrit tubes with cork stoppers to fit
- 9 Electric centrifuge

###### Reagents

- 1 Sterile solution of T 1824 Dye The solution of the dye T 1824 is prepared by dissolving 2 gm in 500 ml distilled water It is passed through a 3 G 3 Jena sintered glass filter and divided between a suitable number of brown bottles and autoclaved Following this rubber diaphragm stoppers are inserted and the bottles stored in the dark Such solutions remain stable and sterile indefinitely when used with reasonable care Each new batch must be calibrated.

###### Procedure

- 1 After necessary equipment has been sterilized and thoroughly dried two hematocrit tubes are placed in a rack One drop (0.05 ml) of heparin is placed in each tube Care must be taken that the tubes are chemically clean and that the heparin is not deposited along the rim of the mouth If the tubes are not absolutely clean some hemolysis may result If heparin is deposited at the mouth it may be covered by the cork stopper that is to be used
- 2 The amount of dye to be injected should be approximately 0.3 mg/kg of body weight The calibrated syringe is filled with the correct quantity of dye
- 3 With the patient in a reclining position 5 ml of blood is removed from a cubital vein with a dry syringe It is important that there be no stasis of venous blood. Consequently the tourniquet must be removed at least one minute before the blood is drawn into the syringe

- 4 The measured amount of dye is then injected slowly through the same needle the syringe being filled and emptied at least five times with the patient's blood in order to wash out all traces of dye
- 5 A notation of the time must be made at the beginning of the injection of the dye. While the dye is being injected an assistant should take the freshly drawn blood and fill the first hematocrit tube to the 10 cm mark. This should be done quickly but with caution, in order that there be no frothing which may cause hemolysis. The cork stopper should be inserted and the tube gently inverted several times until good mixing with the heparin is effected.
- 6 Approximately nine minutes after beginning the injection of the dye a clean needle on a clean dry syringe is inserted into a cubital vein. The tourniquet is removed and time allowed for venous stasis to correct itself, and when exactly ten minutes have elapsed from the time of injection of dye 5 ml of blood is again withdrawn.
- 7 The needle and syringe are removed and the blood cautiously placed in the second hematocrit tube which is inverted and stoppered in the same manner as before. It will be noted that some blood will collect at the junction of the stoppers and the tubes. This should be carefully removed with a cotton tipped applicator, the corks likewise wiped off and then replaced.
- 8 Both tubes should then be centrifuged for thirty minutes at 3000 r.p.m.
- 9 The hematocrit is then read in the conventional manner.
- 10 The supernatant plasma is removed by means of capillary pipettes placed in 10 x 75 mm cuvettes and centrifuged to remove any suspended cells.
- 11 The first cuvette which contains the sample without dye is placed in the spectrophotometer with the wave length set at 6.0 millimicrons. The machine is adjusted to 0 density with this sample. The second sample which contains dye is then inserted and the density read and recorded.

#### Calibration Curves

Each batch of dye must be used to make new calibration curves. The K value is determined in the presence of plasma at a wave length of 6.0  $\mu$ . The procedure is

- 1 It is unwise to try to calculate exactly the dye strength in mg/ml because each batch of dye seems to have its own solubility and color characteristics. The stock solution contains approximately 4 mg/ml.
- 2 Into 50 ml volumetric flask pipette exactly 1 and 0.5 ml of stock solution and dilute with distilled water. The strengths are approximately 0.08 and 0.04 mg/ml.
- 3 Prepare a batch of 10-15 ml plasma for use in the standardization.
- 4 In separate cuvettes mix 1.8 ml of plasma with exactly 0.2 ml of diluted dye. This gives plasma samples of approximately 0.008 and 0.004 mg/ml.
- 5 Set a control cuvette containing 1.8 ml plasma and 0.2 ml isotonic saline at 0 density at wave length 6.0  $\mu$ .
- 6 Compute K from the formula  $K = D/C$  where D is the optical density and C is the dye equivalent in mg/ml. C is understood to be an arbitrary figure varying in reality with each batch of dye.

#### Calculation

$$1 \text{ Circulatory plasma volume ml} = \frac{(\text{mg dye injected})}{(\text{mg dye/ml plasma})}$$

Where mg dye injected is calculated from the equation (mg dye/ml solution) x (ml injected) and mg dye/ml plasma is calculated from the equation D/K.

- 2 The total circulating blood volume is calculated from the plasma volume and the hematocrit. One correction has to be made since the hematocrit of peripheral blood does not necessarily equal that of the total blood. Several authors have

# BLOOD VOLUME AND CIRCULATING RED CELL MASS FROM OBSERVED PLASMA VOLUME AND HEMATOCRIT




Fig. 54.—Line chart for calculating blood volume and circulating red cell mass from observed plasma volume and hematocrit. Note that the hematocrit is multiplied by an arbitrary correction factor of 0.92.

urged that for purposes of calculating blood volume the peripheral hematocrit must be corrected by a constant factor 0.92

$$\text{Total circulating blood volume ml} = \frac{(\text{ml plasma volume}) \times (100)}{(100) - (\text{Hematocrit} \times 0.92)}$$

- 3 From plasma volume and total circulating blood volume total circulating red cell volume ml = ml total blood volume - ml plasma volume

#### Example

Into a 6 kg man were injected 1 mg of dye. The plasma concentration was 0.0075 mg/ml. Hematocrit was 43 ml %

- 1 Circulating plasma volume =

$$\frac{1}{0.0075} = 2800 \text{ ml}$$

- 2 Corrected circulating total blood volume =

$$\frac{(2800) \times (100)}{100 - (43 \times 0.9)} = 4632 \text{ ml}$$

- 3 Corrected circulating red cell mass =

$$4632 - 2800 = 1832 \text{ ml}$$

- 4 These calculations are simplified by use of the line chart Fig. 54

#### Precautions

- 1 The calibration of each new batch of dye is essential for consistent results
- 2 Hemolysis vitiates the results. Although some authors have proposed means of correcting for hemolysis it is our practice to discard any samples in which hemolysis has occurred
- 3 While extrapolation to zero time by employing a series of successive samples would in all probability give more accurate values it was found that in order to obtain a reliable straight line too many samples were required to make the method practical. Consequently a single ten minute sample was used. The validity of using a single sample has been reported by Gregersen, Noble and Miller, Hopper, Tabor and Winkler and Elkinton, Gilder, Muller and Phillips have shown that there is complete mixing of the dye with the circulating plasma in this period of time

## 2 Total Circulating Red Cell, Blood and Plasma Volumes (Carbon Monoxide)

#### Reference

- Root W. S., Roughton F. J. W. and Gregersen M. I. Simultaneous Determinations of Blood Volume by CO and Dye (T 18 4) Under Various Conditions. *Am. J. Physiol.* 146: 739-755 (Aug.) 1946

#### Principle

A measured volume of CO is introduced into the blood by inhalation. The increase in CO content of the blood divided into the total amount of CO administered gives the blood volume. The individual estimations required are

- 1 Carbon monoxide content of the blood before and after inhalation of carbon monoxide
- 2 Loss of CO if any during and after inhalation
- 3 Carbon monoxide content and volume of inspired air
- 4 Residual unabsorbed CO in the expired air

#### Apparatus

- 1 Scholander Roughton apparatus for estimating carbon monoxide in blood.
- 2 Apparatus for generating carbon monoxide

- 3 Douglas bag apparatus for collecting expired air The method described under metabolic procedures is used with one exception inspiration of outside air is from a six foot length of hose attached to the mouthpiece inlet valve

#### Reagents

- 1 Those for the Scholander Roughton method of estimating carbon monoxide
- 2 Carbon monoxide (CO) Place formic acid ( $\text{HCOOH}$ ) in a flask From a dropping funnel add conc sulfuric acid ( $\text{H}_2\text{SO}_4$ ) Pass the escaping carbon monoxide through a trap filled with pellets of dry KOH and collect in a large bottle by displacing water The stored gas should be 96 to 98% CO

#### Procedure

- 1 Allow the subject to lie down for at least 30 minutes
- 2 Draw a sample of venous blood without stasis from an antecubital vein using heparin for the anticoagulant The blood is collected in a syringe protected from light and stored in the refrigerator
- 3 Place the nose clip in place and allow the subject to breathe outdoor air through a mouthpiece connected with a 6 foot length of rubber hose one inch in diameter The outlet valve is connected to a 100 liter Douglas bag, which has previously been rinsed well with outside air
- 4 A 240 ml tonometer containing CO over mercury is connected by means of a short tightly fitted rubber tube with a #20 hypodermic needle, which is thrust into the lumen of the intake rubber hose
- 5 Immediately after starting to collect the subject's expired air a measured volume of about 240 ml of CO of analyzed purity is introduced into the intake hose over a period of from 3 to 4 minutes
- 6 Collect expired air for 10 minutes During the last 7 minutes the subject must be inspiring CO free air and the residual CO in his lungs at the end of that time must be negligible compared to that collected in the Douglas bag
- 7 Allow the subject to rest for 60 minutes more and collect 4 blood samples at timed intervals of 15 minutes from the time the CO was administered.
- 8 Analyze the expired air and the blood samples for CO using the Scholander Roughton micro technique
- 9 Mix the gas thoroughly in the Douglas bag and push the volume through a calibrated gasometer The total volume of the expired air is obtained by the difference in gasometer readings
- 10 Estimate the hematocrit value of the four blood samples and take the average

#### Calculation

- 1 Increase in blood CO ml/100 ml =  
 $(\text{ml CO}/100 \text{ ml at zero time} - \text{ml CO}/100 \text{ ml before inhalation})$   
 where CO at zero time is obtained by extrapolation when the logarithms of the four timed blood specimens are plotted against time
- 2 Volume of CO inspired dry STP ml =  
 $(\text{ml injected into tube}) \times \text{STP factor} \times \text{purity factor}$   
 where ml injected is measured STP factor is obtained from the line chart in the section on respiratory metabolism and the purity factor is determined by analysis of the gas administered
- 3 Total volume of expired air dry STP is calculated according to the procedure given under the Douglas bag technique from volume wet at room temperature  $T$  and ambient barometric pressure
- 4 CO expired ml, dry STP =  

$$\frac{\% \text{ CO in expired air}}{100} \times \text{volume of expired air in ml}$$
 where % CO is determined by analysis and volume of expired air is calculated as in 3 above

5 CO absorbed ml dry STP = (ml CO inspired - ml CO exp red)

where these values are as found in 2 and 4 above

6 Total circulating red cell volume ml =

$$\frac{\text{ml CO absorbed} \times 100}{\text{increase in blood CO ml/100 ml}}$$

7 Total circulating blood volume ml =

$$\frac{\text{red cell volume} \times 100}{\text{average \% hematocrit}}$$

8 Total circulating plasma volume ml =

$$(\text{total blood volume} - \text{total red cell volume})$$

#### Example

.34 ml of CO was inhaled at 24.5°C and 753 mm Hg STP factor 0.881. The CO was 96% pure, total volume of wet expired air was 66.16 l with a CO content of 0.117%. The extrapolated zero CO content was 8.10 ml CO/100 ml, the content before inhalation was 3.00 ml CO/100 ml. The average hematocrit was 47.6 ml/100 ml.

1 Increase in blood CO ml/100 ml =

$$8.10 - 3.00 = 5.10$$

2 CO absorbed ml =

$$(.34 \times 0.881 \times 0.96) - (66.16 \times 0.881 \times 1000 \times \frac{0.117}{100}) = 128 - 68 = 130 \text{ ml}$$

3 Total circulating red cell volume ml =

$$\frac{130 \times 100}{5.10} = 50 \text{ ml}$$

4 Total circulating blood volume ml =

$$\frac{50 \times 100}{47.6} = 6000 \text{ ml}$$

5 Total plasma volume ml =

$$6000 - 50 = 5950 \text{ ml}$$

6 These calculations may be facilitated by use of the line chart (Fig. 54)

#### Precautions

1 With small subjects it may be desirable to avoid toxic effects by decreasing the inspired CO to between 150 and 200 ml.

### 3 True Blood Volume

#### Reference

Same as for dye method and carbon monoxide method for estimating blood volume.

#### Principle

The dye method gives a reliable estimate of circulating plasma volume; the carbon monoxide estimates reliably the circulating red cell mass. Calculations of blood volume by either method involves empirical corrections of the hematocrit of peripheral blood to make it represent total blood hematocrit. By addition of the plasma volume as determined by the dye method and the red cell volume as determined by the carbon monoxide method the true blood volume is calculated.

#### Apparatus

1. Apparatus for plasma volume by the dye method.

2. Apparatus for red cell volume by the carbon monoxide method.

#### Reagents

1. Reagents for the dye method.

2. Reagents for the carbon monoxide method.

**Procedure**

- 1 Estimate plasma volume by the dye method.
- 2 Estimate red cell volume by the carbon monoxide method

**Calculation**

' True blood volume ' = ml plasma volume + ml total red cell volume

**Example**

By the dye method, the plasma volume was 3780 ml By the carbon monoxide method the circulating red cell volume was 2550 ml.

' True blood volume ' =  $3780 + 2550 = 6330$

**Precautions**

Same as for individual methods

**4 Extracellular Fluid (Thiocyanate Space)****References**

- 1 Crandall L. A and Anderson M. X Estimation of the State of Hydration of the Body by the Amount of Water Available for the Solution of Sodium Thiocyanate *Am J Digest Dis & Nutrition* 1 126 131 1934
- 2 Barker M. H The Blood Cyanates in the Treatment of Hypertension *J. A. M. A* 106 762 767 (March) 1936

**Principle**

It is alleged that when thiocyanate is injected into the blood stream it becomes distributed uniformly in the interstitial spaces, but not in the cellular mass and that its concentration in plasma after equilibrium is reached is the same as that in the extracellular fluid.

**Apparatus**

- 1 Sterile needles and syringes for intravenous injection
- 2 Coleman Jr Spectrophotometer
- 3 Cuvettes 19 x 150 mm.
- 4 A quantity of 15 ml round bottom centrifuge tubes with No 1 solid rubber stoppers
- 5 Syringe pipettes 1 ml and 5 ml.
- 6 An electric centrifuge

**Reagents**

- 1 Trichloroacetic acid a 10% solution
- 2 Ferric nitrate solution. Dissolve 50 gm of ferric nitrate crystals in 600 ml of water Add 25 ml of concentrated nitric acid and make up to 1000 ml with distilled water
- 3 Thiocyanate standard solution (Use potassium or sodium) Dissolve 1 gm KCNS in 800 ml of distilled water Titrate a 20 ml portion of a silver nitrate solution (2.9195 gm/1000 ml water) acidified with 5 ml of conc HNO<sub>3</sub> with the KCNS solution using ferric ammonium sulfate as the indicator Calculate the amount of water necessary to add to the KCNS solution to make 20 ml equivalent to 20 ml of silver nitrate solution Add the amount of distilled water and mix thoroughly (1 mg/ml)
- 4 Thiocyanate for injection. Dissolve 30 gm KCNS in 1000 ml redistilled water Shake vigorously to dissolve and filter carefully through a sintered glass filter Store about 11 ml in a 15 ml ampoule and seal. 10 ml = 0.3 gm KCNS Each batch must be analyzed for KCNS
- 5 Redistilled water

**Procedure—Collection of Samples**

- 1 Draw 15 ml of blood for a control sample
- 2 Now slowly inject intravenously 10 ml of KCNS equivalent to 0.3 gm KCNS and record the time
- 3 Draw specimens of blood at 1 hour 3 hours and at 24 hours for analysis
- 4 Just before injection of KCNS have the subject empty his bladder. Discard this specimen. For the next 24 hours collect urine quantitatively in three separate bottles: the first for the period 0 to 1 hour, the second for the period 1 to 3 hours and the third for the period 3 to 24 hours
- 5 Allow each specimen of blood to clot, collect serum and analyze the serum and the corresponding urine sample for KCNS

**Procedure—Analysis of Serum**

- 1 Pipette 5 ml of trichloroacetic acid into a 15 ml round bottom centrifuge tube. Add exactly 5 ml of serum, shake vigorously and centrifuge
- 3 Pipette exactly 5 ml of filtrate into a cuvette
- 4 Add 1 ml ferric nitrate solution
- 5 For a blank, pipette 5 ml of trichloroacetic acid into a cuvette and add 1 ml ferric nitrate solution
- 6 Set blank at 100% T at 550 mμ and read the unknowns.

**Procedure—Analysis of Urine**

- 1 Use 5 ml of urine and add 5 ml of trichloroacetic acid
- 2 Mix and pipette 5 ml into a cuvette
- 3 Add 1 ml ferric nitrate and read at 550 mμ
- 4 If reading is too high use smaller volume of urine and add more trichloroacetic acid.

**Calibration Curve**

1. Make up a series of standards so that 5 ml contains from 0.10 to 0.50 mg of KCNS
- 2 Pipette exactly 2.5 ml into a cuvette
- 3 Add 2.5 ml of trichloroacetic acid and 1 ml of ferric nitrate. Mix and read at 550 mμ with the blank set at 100% T
- 4 Plot mg/cuvette vs % T and use this curve in the subsequent calculations.

**Calculation—Thiocyanate Concentration**

mg KCNS/100 ml serum or urine =

$$\frac{(\text{mg in cuvette}) \times (\text{total dilution}) \times 100}{(\text{mg aliquot}) \times (\text{ml serum or urine used})}$$

**Calculation—Extracellular Fluid**

kilograms water in extracellular fluid =

$$\frac{\text{ml KCNS injected} - \text{total mg KCNS excreted at time of blood}}{\text{mg KCNS/100 ml serum}}$$

**Example**

Ten ml of KCNS (0.3 gm) were injected. Blood levels at 1, 3 and 24 hours were 7, 4 and 3 mg/100 ml. Urine concentrations for the corresponding specimens of urine were 5, 15 and 5 mg/100 ml with volumes of 60, 90 and 900 ml respectively.

kg water in extracellular fluid =

$$\frac{300 - \left( 5 \times \frac{60}{100} \right)}{7 \times 10} = 4.00 \text{ kg}$$



Sample II

$$\frac{300 - \left(20 \times \frac{60}{100}\right) - \left(15 \times \frac{200}{100}\right)}{4 \times 10} = 6.40 \text{ kg}$$

Sample III

$$\frac{300 - \left(25 \times \frac{60}{100}\right) - \left(15 \times \frac{200}{100}\right) - \left(5 \times \frac{900}{100}\right)}{2 \times 10} = 10.5 \text{ kg}$$

**Precautions**

- 1 This method has never been proved to measure extracellular fluid.
- 2 Thiocyanate is a toxic drug. Do not inject more than 0.3 gm/day. Toxemia may show itself at blood levels over 15 mg/100 ml.

## 5 Plasma Volume, Blood Volume and Extracellular Fluid, Miscellaneous References

**Blood Volume**

- Kelly F J, Simonsen D H and Elman R. Blood Volume Determination in the Human With Red Cells Containing Radioactive Phosphorus ( $P^{32}$ ) and With Pure Human Albumin. *J Clin Investigation* 27: 795-804 (Nov) 1948
- Nachman H. M., James III G. W., Moore J. W. and Evans E. L. A Comparative Study of Red Cell Volumes in Human Subjects With Radioactive Phosphorus Tagged Red Cells and T 1824 Dye. *J Clin. Investigation* 29: 558-564 (Feb) 1950
- Nylin G. Blood Volume Determinations With Radioactive Phosphorus. *Brit Heart J* 7: 81-84 (April) 1945
- Reid Allen and Orr M. A Rapid Method for Determining Blood Volumes by the Use of  $P^{32}$  Labelled Red Cells. *J Clin. Investigation* 29: 313-316 (March) 1950

**Cell Volume**

- Grav S. The Use of Fibrinogen in a Rapid Method of Determining Cell Volume. *Am J Med. Sc* 207: 29-39 (Jan) 1944

**Extracellular Fluid**

- Schwartz I. L. Measurement of Extracellular Fluid by Means of Constant Infusion Technique Without Collection of Urine. *Am J Physiol* 160: 565-571 (March) 1950

**Sodium Test**

- Green D. M. and Farah A. Influence of Sodium Load on Sodium Excretion. *Am J Physiol* 158: 444-456 (Sept) 1949

**T 1824**

- Miller A. T. Jr. Excretion of the Blue Dye T 1824 in the Bile. *Am J Physiol* 151: 229-233 (Nov) 194

## SECTION VII

### STUDIES IN THE FIELD

---

#### A GENERAL CONSIDERATIONS

The conduct of scientific surveys is far easier now than it was during most of World War II. When one of the present authors wrote a paper in 1945 summarizing the experience of his colleagues during the war (Johnson R. E. *A Field Nutritional Laboratory War Medicine* 7: 222-226 April 1945) the scientist had to struggle against difficulties of transportation, lack of permanent construction in many areas to house a field laboratory, the necessity to be self-supporting in almost all respects, and the carrying out of chemical analyses under primitive conditions.

In addition, there was a justifiable suspicion in the minds of many military men in the field that research scientists were a hindrance, not a help, and their cooperation was often given only grudgingly. By comparison, the situation now is almost easy. Everywhere in the world are permanent installations which can be reached in a few days by air. Samples for analysis can be frozen and delivered to a well-equipped base laboratory. Supplies can be obtained by air from the United States very quickly. Finally, most military men realize that the scientists are a help, not a hindrance, and welcome their presence—provided they do not get in the way too much. The construction of well-equipped base laboratories in important areas has been an important step for the biological scientist. Examples are the Canadian laboratory at Fort Churchill, the U. S. Air Force Laboratory at Ladd Air Force Base, and the laboratory supported at Point Barrow by the U. S. Navy.

In spite of the great improvement in conditions, certain basic difficulties remain, of which separation from one's family and friends is one of the most important. The collection, storage, handling, and analysis of specimens is still not as easy in the field as it is in university laboratories. Recruitment of able workers presents another difficulty, because in many universities field studies are not regarded as legitimate academic pursuits. It is still the opinion of the present authors that some phenomena can be studied realistically only in the field, whatever the difficulties that must be overcome.

#### B LABORATORY AND CLINICAL METHODS

There is no handbook in existence covering field methods for the medical sciences. Each worker must develop his own methods. The present writers feel that field methods should be suitable for use in extremes of temperature.

and barometric pressure, should be as simple as is consistent with reliability should require the use only of rugged apparatus, and should be rapid enough to handle the large number of specimens which field surveys almost always accumulate. Listed below are selected references for chemical and clinical procedures which are adaptable for medical surveys in the field.

### 1 *Chemical and Clinical Laboratory Methods*

- a Many of the methods described in the present manual have been labelled specifically as field methods. Most of the other methods in the present manual have been used in the field with success.
- b Bessey and Lowry have devised a system of blood and urine analysis requiring only very small samples. Their excellent methods require highly trained technicians and rather expensive apparatus but are otherwise extremely useful for survey work because of their reliability and rapidity.

#### Use of the Beckman Spectrophotometer for Ultramicro Methods

Lowry O H and Bessey O A. The Adaptation of the Beckman Spectrophotometer to Measurements on Minute Quantities of Biological Materials. *J Biol Chem.* 163 633 639 (June) 1946

#### Serum Protein

Lowry O H and Hunter T H. Determination of Serum Protein Concentration With a Gradient Tube. *J Biol Chem.* 159 465 475 (July) 1945

#### Phosphatase

Bessey O A, Lowry O H and Brock M J. Method for the Determination of Alkaline Phosphatase With Five Cubic Millimeters of Serum. *J Biol Chem.* 164 3 13-9 (July) 1946

#### Inorganic Phosphate

Lowry O H and Lopez J. The Determination of Inorganic Phosphate in the Presence of Labile Phosphate Esters. *J Biol Chem.* 162 421 423 (March) 1946

#### Vitamin A

Bessey, O A, Lowry O H, Brock M J and Lopez J A. The Determination of Vitamin A and Carotene in Small Quantities of Blood Serum. *J Biol Chem.* 166 177 188 (Nov) 1946

#### Vitamin C

Lowry O H, Lopez J A., and Bessey O A. The Determination of Ascorbic Acid in Small Amounts of Blood Serum. *J Biol Chem.* 160 609 615 (Oct) 1945

Bessey O A. A Method for the Determination of Small Quantities of Ascorbic Acid and Dehydroascorbic Acid in Turbid and Colored Solutions in the Presence of Other Reducing Substances. *J Biol. Chem.* 126 771 784 (Dec) 1933

#### Riboflavin

Lowry O H, Bessey O A and Hull T Z. Method for the Determination of Five tenths to Two Millimicrograms of Riboflavin. *J Biol Chem.* 155 71 77 (Sept) 1944.

### 2 *Clinical Diagnosis*

The Food and Nutrition Board has produced a monograph entitled *Nutrition Surveys—Their Techniques and Value*. Bulletin of the National Research Council No 17 May 1949. The chapter on clinical diagnosis presents a critical discussion of the clinical study of populations.

# C TREADMILL TEST OF PHYSICAL FITNESS FOR HARD MUSCULAR WORK

## Reference

Johnson R E Brouha L and Darling R C A Test of Physical Fitness for Strenuous Exertion *Rev Canad de Biol* 1 491 53 (June) 1954

## Principle

This test is a simple rapid means of assessing the physical fitness of normal healthy men for hard muscular work Any such test must be judged by the following considerations

- 1 Athletes in training should on the average get much better scores than men not training  
For a large group of untrained men the scores should range widely from bad to good.
- 2 Normal men not in training should be able to improve their scores by following a program of hard physical work
- 3 In the case of any given man the score should in general agree with the judgment of trained observers with the results of other reliable tests and most important of all with the man's demonstrated physical capacities in actual performance of hard work

As judged by these criteria certain tests now in common use are of little value in assessing fitness for hard muscular work In our experience when hard muscular work is in question one must expose the subject to hard work in order to test him accurately This is chiefly because the easier the work the smaller and less regular are the differences between the fit and the unfit Satisfactory assessment of a man's fitness for hard and prolonged work such as marching long distances can be made by subjecting him to a standard exhausting exercise and taking into account only two factors the length of time he can sustain it and his pulse rate in recovery after the work is over

## Apparatus

- 1 A motor driven treadmill
- 2 An accurate stopwatch

## Procedure

- 1 The subject warms up by walking uphill for 5 minutes at 35 miles per hour at a grade of 8.6%
- 2 He sits on a chair for 5 minutes
- 3 At a signal he begins to run uphill at a speed of 7 miles per hour and a grade of 8.6% If not exhausted he is stopped at 5 minutes The duration of the run is noted to the nearest second
- 4 Commence men using the time of recovery from the time when he stops
- 5 Count the radial pulses from 1 to 15 from 2 to 25 and from 4 to 45 minutes of recovery

## Calculation

$$\begin{aligned} \text{Arbitrary score of physical fitness} = \\ \frac{(\text{seconds the subject ran}) \times (100)}{(\quad) \times (\text{sum of three half minute pulse rates})} \end{aligned}$$

## Example

The subject ran 3 minutes 45 seconds. His three half minute pulse rates were 85 90 and 60 respectively

- 1 Arbitrary score of physical fitness =

$$\frac{5 \times 100}{3 \times 15} = 5$$

— The line chart Fig 5 facilitates these computations.

# LINE CHART FOR CALCULATING INDEX OF PHYSICAL FITNESS! (HARVARD FATIGUE LABORATORY)




Fig 55—Calculation of index of physical fitness for hard muscular work in the treadmill, peck and step tests.

### Interpretation

From measurements on several hundreds of men of military age of all degrees of fitness the following interpretation of scores is made  
 A poor score is below 40  
 An average score is from 41 to 75  
 A good score is from 76 to 90  
 A superior score is above 90

### Precautions

1 Any test involving the cooperation of the subject requires the will to work anatomical adequacy to do the work and physiological fitness for the task. A poor score may mean that any one of the three factors or a combination of them has been weak

## D FIELD TESTS OF PHYSICAL FITNESS

### 1 The Pack Test

#### Reference

Same as for the treadmill test

#### Principle

For field purposes the treadmill is not usually available. Therefore a modified test is recommended based on the same general principle. A test that involves hard work. Some of the accuracy of the treadmill is necessarily sacrificed for the sake of field application.

#### Apparatus

- 1 Scales accurate to one pound
- 2 One watch preferably a stop watch
- 3 One haversack for each subject
- 4 A series of bags of metal rock or sand weighing 10 5 and 1 pound
- 5 Stepping platform similar to that of the McCordy Larson test. In practice standard gymnasium stools 16 high are very good when placed as close as possible to the standard gymnasium stall bars. Other platforms can be improvised from boxes or chairs 16 high. Ladders and broomsticks can be used for hand hold. One platform is needed for each subject.

#### Procedure for Testing One Subject

- 1 Weigh the subject stripped or if he is clothed subtract the approximate weight of clothes usually about pounds
- Load the pack according to the man's weight as follows

#### MAN'S WEIGHT STRIPPED

#### TOTAL PACK WEIGHT

|                 |           |
|-----------------|-----------|
| Below 10 pounds | 40 pounds |
| 10 135 pounds   | 45 pounds |
| 135 150 pounds  | 50 pounds |
| 150 165 pounds  | 55 pounds |
| 165 180 pounds  | 60 pounds |
| 180 195 pounds  | 65 pounds |
| 195 210 pounds  | 70 pounds |
| Over 210 pounds | 75 pounds |

Note that the pack weight is approximately 1/3 of the subject's body weight. Be sure that the subject wears a shirt or some other protect on to prevent chafing the shoulders and back. Help him put on the appropriate pack

# LINE CHART FOR CALCULATING INDEX OF PHYSICAL FITNESS: (HARVARD FATIGUE LABORATORY)




Fig 55—Calculation of index of physical fitness for hard muscular work in the treadmill, pack and step tests.

4. The times mentioned should be followed as strictly as possible
5. Observers unaccustomed to counting the pulse rate in exercise must practice counting rapid pulses before they can be expected to obtain accurate results.

## 2 The Step Test

### Reference

1. Brouha L. Graybiel A. and Heath C W. *Rev Canad de biol* 2: 86-91 1943
2. Johnson R E. Crowley L V. Dorsey E A. Jr. Swana H L. Price J I. Saffell W B. and Jones E O. *Test of Operational Ration Type E* pages 67-70 (March) 1947. Medical Nutrition Laboratory Report #12

### Principle

The present test is recommended for use only when it is impossible to use one of the more reliable tests. The step test was devised for speed and simplicity of use with large numbers of subjects but a great deal of accuracy has been sacrificed. The 30 inch step is so high that it tends to penalize short men and also to involve to some extent the factor of agility as well as general fitness. Consequently gross errors may sometimes be made in judging very short or clumsy men. However if a lower step is used the work is too easy which so increases the number of errors and impairs the value of the test.

### Apparatus

1. One watch preferably a stop watch for each observer
2. A stepping platform 30 inches high (plus or minus not more than 3/4 inch). This should be sturdily constructed so as not to give when stepped on. It should be at least 14 inches from front to back. When many men are to use the same step simultaneously it should be long enough to allow at least 30 inches for each man. These 30 inch spaces should be marked off. Standard gymnasium stools if adjusted to 30 inches in height are good.
3. A device for keeping rhythm. The subjects must go up and down every two seconds. If a metronome is used it must be a good one that will work without error for 5 minutes. We have used satisfactorily a simple pendulum consisting of a weight on a 39 inch string. It must be calibrated daily.

### Procedure

1. As many subjects as there are observe a line up in front of the stepping platform. They should be stripped to the underwear and should wear rubber soled shoes or no shoes at all. One observer calls the rhythm. At the signal Start! each subject places one foot on the platform steps up placing both feet on the platform straightens the legs and back and immediately steps down again bringing down first the same foot he placed up first. At exactly two second intervals the signal Up! is given and rhythm is maintained by giving the count UP! - 2 - 3 - 4 - UP! - 2 - 3 - 4 etc. The subject should lead off with the same foot each time and not try to alternate feet. However one or two changes of lead off during the test are immaterial. He must not touch anything with his hands but may move his arms freely.
2. Begin counting the time when the subject starts exercising and exercise him for five minutes continuously unless he stops before then from exhaustion. If he falls behind stop him after he has been unable to keep up the pace for 30 seconds. Note the duration of his effort to the nearest second. All men are stopped at 5 minutes if they can go that long.
3. When the subject stops start counting the time and have him sit quietly on a chair.



4. Exercise the subject as follows Have him grasp with both hands a horizontal bar at the height of his shoulders or lightly lower and stand with both feet on the floor The distance between hands should be about the width of his shoulders At the signal "UP" he places one foot on the stool, steps up placing both feet on the stool straightens his legs and back and immediately steps down again bringing down first the same foot he put up first, and retaining his hold on the horizontal bar Every two seconds the signal "UP" is given and the subject steps up and then down He should lead off" with the same foot each time and not try to alternate the feet However one or two changes of "lead off" in the course of the test are immaterial
5. Begin counting the time when the subject starts and exercise him for five minutes continuously unless he stops from exhaustion before then If he falls behind, stop him after he has been unable to keep the pace for twenty seconds Note the duration of the exercise to the nearest second.
6. When he stops start counting the time then take off the pack and have him sit on a stool or chair
7. Count and record the pulse during each of the following intervals 1 to 15 " to 25 and 4 to 45 minutes after he stops working The counts must be limited precisely to these intervals The whole test is completed in ten minutes No other observations are necessary

#### Procedure for Testing a Large Group

1. Procure as many stop watches as there are observers
2. Procure 14 packs and load them as follows
 

|  |
|--|
| 1 pack weighing 40 pounds for subjects weighing below 100 pounds |
| 2 packs weighing 45 pounds for subjects weighing 100-135 pounds  |
| 3 packs weighing 50 pounds for subjects weighing 135-150 pounds  |
| 3 packs weighing 55 pounds for subjects weighing 150-165 pounds  |
| 2 packs weighing 60 pounds for subjects weighing 165-180 pounds  |
| 1 pack weighing 65 pounds for subjects weighing 180-195 pounds   |
| 1 pack weighing 70 pounds for subjects weighing 195-210 pounds   |
| 1 pack weighing 75 pounds for subjects weighing above 210 pounds |
3. Assign two men the task of weighing the subjects giving them packs and acting as secretaries Start as many subjects exercising all at the same time as there are observers to count the pulse One observer stands behind each subject making sure that he has the proper pack seeing that he keeps the rhythm timing the duration of his effort taking off his pack having him sit down counting his pulse and telling them to the secretaries One of the observers gives the signal "Start!" and then calls "UP!" every two seconds With six observers it is possible to test about thirty six subjects every hour

#### Calculation, Example Interpretation

Same as for treadmill test

#### Precautions

1. The height of the platforms should be 16" with an outside tolerance of 0.5"
2. The handhold should be a horizontal bar at shoulder height If it is too high the subject will find the work too easy If it is too low he will find it awkward. Vertical bars can be used but the subject must keep his hands at the level of his shoulders throughout the exercise
3. The rhythm should be maintained as closely as possible It is easy to distinguish between those who are unable to keep the pace because of poor sense of rhythm and those who fall behind because of exhaustion. Only the latter should be stopped. Those with a poor sense of rhythm can be kept to the pace by a tap on the elbow every two seconds

## E SURVEYS OF FOOD CONSUMPTION

### Reference

- 1 Berryman G H and Chatfield C Short Method of Calculating Nutritive Value of Diets *J Nutrition* 25 43 27 (Jan.) 1943
- 2 Berryman G H and Howe P E Short Method of Calculating Nutritive Value of Diets *J Nutrition* 27 31 40 (Mar.) 1944

### Principle

In the field it is often necessary to estimate the average food and nutrient consumption of a group. For this purpose a short survey is done over a period of ten days. The total food consumption (by groups) is calculated from inventory data corrected for plate and kitchen waste. The total nutrient consumption is calculated by adding together the nutrients contributed by each food group.

### Apparatus and Supplies

- 1 Scales folding 300 lb capacity
- 2 Hand tally counter
- 3 Chip board
- 4 Ruler
- 5 At least ten 5 gallon cans or equivalent containers for segregation and weighing of waste both kitchen and plate
- 6 Slide rule
- 7 Adding or calculating machine
- 8 Essential books
  - a Handbook of Chemistry and Physics
  - b National Research Council Tables of Food Composition
  - c Table for converting weights to pounds per man per day
  - d Bowes and Church or other standard book on food composition
  - e Master Menu
  - f TMI 10 412 (Army Recipes)
- 9 Pads of paper 8 x 10½
- 10 Small notebooks one for each observer
- 11 Pencils assorted
- 1 Paper graph one pad
- 13 Clips paper
- 14 8 x 4¼ bibliography cards
- 15 Proper forms
- 16 Rubber bands
- 17 One dozen manila folders 8 x 10½

### Reagents

None unless it is proposed to collect specimens of food. If so see sections on collection of specimens and the individual methods contemplated.

### Procedure

- 1 It is recommended that the survey be planned to cover a complete 10 day period. Shorter periods may give misleading and erroneous results.
- 2 Inventory and Issue A complete inventory is made of food stores in the mess on the night before each survey and on the night after the last meal of each survey. All items issued during the survey are weighed and recorded.
- 3 Head Count At each meal one observer is stationed near the place where trays are returned. Each man is tallied on a hand counter as he hands in his tray. Meals at odd hours especially at night are recorded by night cook and baker.
- 4 Kitchen and Plate Waste Food groups are segregated and weighed in the following categories: meats, fish and poultry, eggs, milk products, butter, other.

- 5 Beginning exactly one minute after he stops count the number of heart beats for exactly 30 seconds The base of the neck is the easiest place to find the pulse after exercise
- 6 Record the duration of effort and the number of heart beats in the 30 second period beginning one minute after he stops No other observations are necessary

#### Calculation

A simplified table is used for scoring (Table 2.) From the duration of effort (left hand vertical column) and the heart beats (top horizontal column) read directly the arbitrary score

#### Example

A subject lasted 3 minutes and 35 seconds His pulse from 1 to 1½ minutes in recovery was 70 beats in the 30 seconds From the table his score was 53

#### Interpretation

This test is easier to perform than either the 'treadmill or pack' tests.

For normal healthy young men

Poor is below 50,

Average is 50 through 80,

Good is above 80

#### Precautions

- 1 The platform must be 20 inches high
- 2 The rhythm should be maintained as closely as possible It is easy to distinguish between those unable to keep pace from exhaustion and those who merely have a poor sense of rhythm Only exhausted men should be stopped before five minutes Those with a poor sense of rhythm can be kept to the pace by a tap on the elbow every 2 seconds
- 3 Times must be followed strictly
- 4 The subjects must straighten legs and back at each step
- 5 Observers unaccustomed to counting rapid pulses must practice before they can be expected to obtain accurate results

TABLE 2°

#### SCORING THE 'STEP TEST' OF PHYSICAL FITNESS

Instructions for using the table (1) Find the appropriate line for duration of effort (°) then find the appropriate column of the pulse count (3) read off the score where the line and the column intersect

| DURATION<br>OF<br>EFFORT<br>MINS | TOTAL HEART BEATS 1 MIN TO 1½ MIN IN RECOVERY |     |     |    |    |    |    |    |    |    |    |    |    |    |
|----------------------------------|---|-----|-----|----|----|----|----|----|----|----|----|----|----|----|
|                                  | 40  | 45  | 50  | 55 | 60 | 65 | 70 | 75 | 80 | 85 | 90 | 95 |    |    |
|                                  | TO  | TO  | TO  | TO | TO | TO | TO | TO | TO | TO | TO | TO | TO | TO |
|                                  | 44  | 49  | 54  | 59 | 64 | 69 | 74 | 79 | 84 | 89 | 94 | 99 |    |    |
| SCORE (ARBITRARY UNITS)          |   |     |     |    |    |    |    |    |    |    |    |    |    |    |
| 0 to ½                           | 6   | 6   | 5   | 5  | 4  | 4  | 4  | 4  | 3  | 3  | 3  | 3  |    |    |
| ½ to 1                           | 19  | 17  | 16  | 14 | 13 | 12 | 11 | 11 | 10 | 9  | 9  | 8  |    |    |
| 1 to 1½                          | 32  | 29  | 26  | 24 | 22 | 20 | 19 | 18 | 17 | 16 | 15 | 14 |    |    |
| 1½ to 2                          | 45  | 41  | 38  | 34 | 31 | 29 | 27 | 25 | 23 | 22 | 21 | 20 |    |    |
| 2 to 2½                          | 58  | 50  | 47  | 43 | 40 | 36 | 34 | 30 | 30 | 28 | 27 | 25 |    |    |
| 2½ to 3                          | 71  | 64  | 58  | 53 | 48 | 45 | 42 | 39 | 37 | 34 | 33 | 31 |    |    |
| 3 to 3½                          | 84  | 75  | 68  | 62 | 57 | 53 | 49 | 46 | 43 | 41 | 39 | 37 |    |    |
| 3½ to 4                          | 97  | 87  | 79  | 72 | 66 | 61 | 57 | 53 | 50 | 47 | 45 | 4  |    |    |
| 4 to 4½                          | 110   | 98  | 89  | 82 | 75 | 70 | 65 | 61 | 57 | 54 | 51 | 48 |    |    |
| 4½ to 5                          | 123   | 110 | 100 | 91 | 84 | 77 | 70 | 68 | 63 | 60 | 57 | 54 |    |    |
| 5                                | 129   | 116 | 105 | 96 | 88 | 80 | 76 | 71 | 67 | 63 | 60 | 56 |    |    |

## **E SURVEYS OF FOOD CONSUMPTION**

### **Reference**

- 1 Berryman G H and Chatfield C Short Method of Calculating Nutritive Value of Diets *J Nutrition* 25 3-7 (Jan) 1943
- 2 Berryman G H and Howe P E Short Method of Calculating Nutritive Value of Diets *J Nutrition* 27 31-40 (Mar) 1944.

### **Principle**

In the field it is often necessary to estimate the average food and nutrient consumption of a group. For this purpose a short survey is done over a period of ten days. The total food consumption (by groups) is calculated from inventory data corrected for plate and kitchen waste. The total nutrient consumption is calculated by adding together the nutrients contributed by each food group.

### **Apparatus and Supplies**

- 1 Scales folding 300 lb capacity
- 2 Hand tally counter
- 3 Clip board.
- 4 Ruler
- 5 At least ten 5 gallon cans or equivalent containers for segregation and weighing of waste both kitchen and plate
- 6 Slide rule
- 7 Adding or calculating machine
- 8 Essential books
  - a Handbook of Chemistry and Physics
  - b National Research Council Tables of Food Composition
  - c Table for converting weights to pounds per man per day
  - d Bowes and Church or other standard book on food composition
  - e Master Menu
  - f TM 10-412 (Army Recipes)
- 9 Pads of paper 8 x 10 1/2
- 10 Small notebooks one for each observer
- 11 Pencils assorted
- 1 Paper graph one pad
- 13 Clips paper
- 14 3 x 4 1/2 bibliography cards
- 15 Proper forms
- 16 Rubber bands
- 17 One dozen manila folders 8 x 10 1/2

### **Reagents**

None unless it is proposed to collect specimens of food. If so see sections on collection of specimens and the individual methods contemplated.

### **Procedure**

- 1 It is recommended that the survey be planned to cover a complete 10 day period. Shorter periods may give misleading and erroneous results.
- 2 Inventory and Issue. A complete inventory is made of food stores in the mess on the night before each survey and on the night after the last meal of each survey. All items issued during the survey are weighed and recorded.
- 3 Head Count. At each meal one observer is stationed near the place where trays are returned. Each man is tallied on a hand counter as he hands in his tray. Meals at odd hours especially at night are recorded by night cook and baker.
- 4 Kitchen and Plate Waste. Food groups are segregated and weighed in the following categories: meats, fish and poultry, eggs, milk products, butter, other.

- 5 Beginning exactly one minute after he stops count the number of heart beats for exactly 30 seconds The base of the neck is the easiest place to find the pulse after exercise
- 6 Record the duration of effort and the number of heart beats in the 30 second period beginning one minute after he stops No other observations are necessary

#### Calculation

A simplified table is used for scoring (Table 22) From the duration of effort (left hand vertical column) and the heart beats (top horizontal column) read directly the arbitrary score

#### Example

A subject lasted 3 minutes and 35 seconds His pulse from 1 to 1½ minutes in recovery was 79 beats in the 30 seconds From the table his score was 53

#### Interpretation

This test is easier to perform than either the "treadmill" or "pack" tests

For normal, healthy young men

Poor is below 50

Average is 50 through 80

Good is above 80

#### Precautions

- 1 The platform must be 20 inches high
- 2 The rhythm should be maintained as closely as possible It is easy to distinguish between those unable to keep pace from exhaustion and those who merely have a poor sense of rhythm. Only exhausted men should be stopped before five minutes Those with a poor sense of rhythm can be kept to the pace by a tap on the elbow every 2 seconds
- 3 Times must be followed strictly
- 4 The subjects must straighten legs and back at each step
- 5 Observers unaccustomed to counting rapid pulses must practice before they can be expected to obtain accurate results

TABLE 22  
SCORING THE STEP TEST OF PHYSICAL FITNESS

Instructions for using the table (1) Find the appropriate line for duration of effort (2) then find the appropriate column of the pulse count (3) read off the score where the line and the column intersect

| DURATION<br>OF<br>EFFORT<br>MINS | TOTAL HEART BEATS 1 MIN TO 1½ MIN IN RECOVERY |     |     |     |    |    |    |    |    |    |    |
|----------------------------------|---|-----|-----|-----|----|----|----|----|----|----|----|
|                                  | 40  | 45  | 50  | 55  | 60 | 65 | 70 | 75 | 80 | 85 | 90 |
|                                  | TO  | TO  | TO  | TO  | TO | TO | TO | TO | TO | TO | TO |
|                                  | 44  | 49  | 54  | 59  | 64 | 69 | 74 | 79 | 84 | 89 | 94 |
|                                  | SCORE (ARBITRARY UNITS)                       |     |     |     |    |    |    |    |    |    |    |
| 0 to ½                           | 6   | 6   | 5   | 5   | 4  | 4  | 4  | 4  | 3  | 3  | 3  |
| ½ to 1                           | 19  | 17  | 16  | 14  | 13 | 12 | 11 | 11 | 10 | 9  | 8  |
| 1 to 1½                          | 32  | 29  | 26  | 24  | 22 | 20 | 19 | 18 | 17 | 16 | 15 |
| 1½ to 2                          | 45  | 41  | 38  | 34  | 31 | 29 | 27 | 25 | 23 | 22 | 21 |
| 2 to 2½                          | 58  | 52  | 47  | 43  | 40 | 36 | 34 | 32 | 30 | 28 | 25 |
| 2½ to 3                          | 71  | 64  | 58  | 53  | 48 | 45 | 42 | 39 | 37 | 34 | 31 |
| 3 to 3½                          | 84  | 75  | 68  | 62  | 57 | 53 | 49 | 46 | 43 | 41 | 39 |
| 3½ to 4                          | 97  | 87  | 79  | 72  | 66 | 61 | 57 | 53 | 50 | 47 | 45 |
| 4 to 4½                          | 110   | 98  | 89  | 82  | 75 | 70 | 65 | 61 | 57 | 54 | 51 |
| 4½ to 5                          | 123   | 110 | 100 | 91  | 84 | 77 | 72 | 68 | 63 | 60 | 57 |
| 5                                | 136   | 123 | 116 | 105 | 96 | 88 | 81 | 76 | 71 | 67 | 63 |

## E SURVEYS OF FOOD CONSUMPTION

### Reference

- 1 Berryman G H and Chatfield C Short Method of Calculating Nutritive Value of Diets *J Nutrition* 25 23-27 (Jan) 1943
- 2 Berryman G H and Howe P E Short Method of Calculating Nutritive Value of Diets *J Nutrition* 27 31-40 (Mar) 1944

### Principle

In the field it is often necessary to estimate the average food and nutrient consumption of a group. For this purpose a short survey is done over a period of ten days. The total food consumption (by groups) is calculated from inventory data corrected for plate and kitchen waste. The total nutrient consumption is calculated by adding together the nutrients contributed by each food group.

### Apparatus and Supplies

- 1 Scales folding 300 lb capacity
- 2 Hand tally counter
- 3 Clip board.
- 4 Ruler
- 5 At least ten 5 gallon cans or equivalent containers for segregation and weighing of waste both kitchen and plate
- 6 Slide rule
- 7 Adding or calculating machine
- 8 Essential books
  - a Handbook of Chemistry and Physics
  - b National Research Council Tables of Food Composition
  - c Table for converting weights to pounds per man per day
  - d Bowes and Church or other standard book on food composition
  - e Master Menu
  - f TM 10-412 (Army Recipes)
- 9 Pads of paper 8 x 10½
- 10 Small notebooks one for each observer
- 11 Pencils assorted.
  - 1 Paper graph one pad
- 12 Clips paper
- 14 3 x 4½ bibliography cards
- 15 Proper forms
- 16 Rubber bands
- 17 One dozen manila folders 8 x 10½

### Reagents

None unless it is proposed to collect specimens of food. If so see sections on collection of specimens and the individual methods contemplated.

### Procedure

- 1 It is recommended that the survey be planned to cover a complete 10 day period. Shorter periods may give misleading and erroneous results.
- 2 Inventory and Issue. A complete inventory is made of food stores in the mess on the night before each survey and on the night after the last meal of each survey. All items issued during the survey are weighed and recorded.
- 3 Head Count. At each meal one observer is stationed near the place where trays are returned. Each man is tallied on a hand counter as he hands in his tray. Meals at odd hours especially at night are recorded by night cook and baker.
- 4 Kitchen and Plate Waste. Food groups are segregated and weighed in the following categories: meats, fish and poultry, eggs, milk products, butter, other

- Beginning exactly one minute after he stops count the number of heart beats for exactly 30 seconds. The base of the neck is the easiest place to find the pulse after exercise.
- Record the duration of effort and the number of heart beats in the 30 second period beginning one minute after he stops. No other observations are necessary.

#### Calculation

A simplified table is used for scoring (Table 22). From the duration of effort (left hand vertical column) and the heart beats (top horizontal column) read directly the arbitrary score.

#### Example

A subject lasted 3 minutes and 35 seconds. His pulse from 1 to 1½ minutes in recovery was 79 beats in the 30 seconds. From the table his score was 53.

#### Interpretation

This test is easier to perform than either the "treadmill" or "pack" tests.

For normal, healthy young men

Poor is below 50

Average is 50 through 80

Good is above 80

#### Precautions

- The platform must be 20 inches high.
- The rhythm should be maintained as closely as possible. It is easy to distinguish between those unable to keep pace from exhaustion and those who merely have a poor sense of rhythm. Only exhausted men should be stopped before five minutes. Those with a poor sense of rhythm can be kept to the pace by a tap on the elbow every 2 seconds.
- Times must be followed strictly.
- The subjects must straighten legs and back at each step.
- Observers unaccustomed to counting rapid pulses must practice before they can be expected to obtain accurate results.

TABLE 22

#### SCORING THE STEP TEST OF PHYSICAL FITNESS

Instructions for using the table: (1) Find the appropriate line for duration of effort ( ) then find the appropriate column of the pulse count (3) read off the score where the line and the column intersect.

| DURATION<br>OF<br>EFFORT<br>MINS. | TOTAL HEART BEATS 1 MIN. TO 1½ MIN. IN RECOVERY |     |     |    |    |    |    |    |    |    |    |    |
|-----------------------------------|---|-----|-----|----|----|----|----|----|----|----|----|----|
|                                   | 40  | 45  | 50  | 55 | 60 | 65 | 70 | 75 | 80 | 85 | 90 | 95 |
|                                   | TO  | TO  | TO  | TO | TO | TO | TO | TO | TO | TO | TO | TO |
|                                   | 44  | 49  | 54  | 59 | 64 | 69 | 74 | 79 | 84 | 89 | 94 | 99 |
|                                   | SCORE (ARBITRARY UNITS)                         |     |     |    |    |    |    |    |    |    |    |    |
| 0 to ½                            | 6   | 6   | 5   | 5  | 4  | 4  | 4  | 4  | 3  | 3  | 3  | 3  |
| ½ to 1                            | 19  | 17  | 16  | 14 | 13 | 12 | 11 | 11 | 10 | 9  | 9  | 8  |
| 1 to 1½                           | 30  | 29  | 26  | 24 | 20 | 20 | 19 | 18 | 17 | 16 | 15 | 14 |
| 1½ to 2                           | 45  | 41  | 38  | 34 | 31 | 29 | 27 | 25 | 23 | 22 | 21 | 0  |
| 2 to 2½                           | 58  | 52  | 47  | 43 | 40 | 36 | 34 | 30 | 30 | 28 | 7  | 25 |
| 2½ to 3                           | 71  | 64  | 58  | 53 | 48 | 45 | 42 | 39 | 37 | 34 | 33 | 31 |
| 3 to 3½                           | 84  | 75  | 68  | 60 | 57 | 53 | 49 | 46 | 43 | 41 | 39 | 37 |
| 3½ to 4                           | 97  | 87  | 79  | 70 | 66 | 61 | 57 | 53 | 50 | 47 | 45 | 4  |
| 4 to 4½                           | 110   | 98  | 89  | 80 | 75 | 70 | 65 | 61 | 57 | 54 | 51 | 48 |
| 4½ to 5                           | 13  | 110 | 100 | 91 | 84 | 77 | 70 | 68 | 63 | 60 | 57 | 54 |
| 5                                 | 13  | 110 | 100 | 91 | 84 | 77 | 70 | 68 | 63 | 60 | 57 | 54 |

## E SURVEYS OF FOOD CONSUMPTION

### Reference

- 1 Berryman G H. and Chatfield C Short Method of Calculating Nutritive Value of Diets *J Nutrition* 25 43-7 (Jan) 1943
- 2 Berryman G H and Howe P E Short Method of Calculating Nutritive Value of Diets *J Nutrition* 27 231-240 (Mar) 1944.

### Principle

In the field it is often necessary to estimate the average food and nutrient consumption of a group. For this purpose a short survey is done over a period of ten days. The total food consumption (by groups) is calculated from inventory data corrected for plate and kitchen waste. The total nutrient consumption is calculated by adding together the nutrients contributed by each food group.

### Apparatus and Supplies

- 1 Scales folding 300 lb capacity
- 2 Hand tally counter
- 3 Chip board.
- 4 Ruler
- 5 At least ten 5 gallon cans or equivalent containers for segregation and weighing of waste both kitchen and plate
- 6 Slide rule
- 7 Adding or calculating machine
- 8 Essential books
  - a Handbook of Chemistry and Physics
  - b National Research Council Tables of Food Composition
  - c Table for converting weights to pounds per man per day
  - d. Bowes and Church or other standard book on food composition
  - e Master Menu
  - f TM 10 412 (Army Recipes)
- 9 Pads of paper 8 x 10 1/4"
- 10 Small notebooks one for each observer
- 11 Pencils assorted.
- 1 Paper graph one pad
- 13 Clips paper
14. 3 x 4 1/2 bibliography cards
- 15 Proper forms
- 16 Rubber bands
- 17 One dozen manila folders 8 x 10 1/4"

### Reagents

None unless it is proposed to collect specimens of food. If so see sections on collection of specimens and the individual methods contemplated.

### Procedure

- 1 It is recommended that the survey be planned to cover a complete 10 day period. Shorter periods may give misleading and erroneous results.
- 2 **Inventory and Issue** A complete inventory is made of food stores in the mess on the night before each survey and on the night after the last meal of each survey. All items issued during the survey are weighed and recorded.
- 3 **Head Count** At each meal one observer is stationed near the place where trays are returned. Each man is tallied on a hand counter as he hands in his tray. Meals at odd hours especially at night are recorded by night cook and baker.
- 4 **Kitchen and Plate Waste** Food groups are segregated and weighed in the following categories: meats, fish and poultry, eggs, milk products, butter, other.



- 5 Beginning exactly one minute after he stops count the number of heart beats for exactly 30 seconds The base of the neck is the easiest place to find the pulse after exercise
- 6 Record the duration of effort and the number of heart beats in the 30 second period beginning one minute after he stops No other observations are necessary

#### Calculation

A simplified table is used for scoring (Table 22) From the duration of effort (left hand vertical column) and the heart beats (top horizontal column) read directly the arbitrary score

#### Example

A subject lasted 3 minutes and 35 seconds His pulse from 1 to 1½ minutes in recovery was 79 beats in the 30 seconds From the table his score was 53

#### Interpretation

This test is easier to perform than either the 'treadmill or pack' tests.

For normal healthy young men

Poor is below 50

Average is 50 through 80

Good is above 80

#### Precautions

- 1 The platform must be 20 inches high
- 2 The rhythm should be maintained as closely as possible It is easy to distinguish between the e unable to keep pace from exhaustion and those who merely have a poor sense of rhythm Only exhausted men should be stopped before five minutes Those with a poor sense of rhythm can be kept to the pace by a tap on the elbow every 2 seconds
- 3 Times must be followed strictly
- 4 The subjects must straighten legs and back at each step
- 5 Observers unaccustomed to counting rapid pulses must practice before they can be expected to obtain accurate results

TABLE 22  
SCORING THE 'STEP TEST' OF PHYSICAL FITNESS

Instructions for using the table (1) Find the appropriate line for duration of effort ( ) then find the appropriate column of the pulse count (3) read off the score where the line and the column intersect

| DURATION<br>OF<br>EFFORT<br>MINS | TOTAL HEART BEATS 1 MIN TO 1½ MIN IN RECOVERY |     |     |    |    |    |    |    |    |    |    |    |    |    |
|----------------------------------|---|-----|-----|----|----|----|----|----|----|----|----|----|----|----|
|                                  | 40  | 45  | 50  | 55 | 60 | 65 | 70 | 75 | 80 | 85 | 90 | 95 |    |    |
|                                  | TO  | TO  | TO  | TO | TO | TO | TO | TO | TO | TO | TO | TO | TO | TO |
|                                  | 44  | 49  | 54  | 59 | 64 | 69 | 74 | 79 | 84 | 89 | 94 | 99 |    |    |
|                                  | SCORE (ARBITRARY UNITS)                       |     |     |    |    |    |    |    |    |    |    |    |    |    |
| 0 to ½                           | 6   | 6   | 5   | 5  | 4  | 4  | 4  | 3  | 3  | 3  | 3  | 3  |    |    |
| ½ to 1                           | 19  | 17  | 16  | 14 | 13 | 12 | 11 | 11 | 10 | 9  | 9  | 8  |    |    |
| 1 to 1½                          | 3   | 29  | 26  | 24 | 20 | 20 | 19 | 18 | 17 | 16 | 15 | 14 |    |    |
| 1½ to 2                          | 45  | 41  | 38  | 34 | 31 | 29 | 27 | 25 | 23 | 22 | 21 | 20 |    |    |
| 2 to 2½                          | 58  | 52  | 47  | 43 | 40 | 36 | 34 | 31 | 30 | 28 | 27 | 25 |    |    |
| 2½ to 3                          | 71  | 64  | 58  | 53 | 48 | 45 | 41 | 39 | 37 | 34 | 33 | 31 |    |    |
| 3 to 3½                          | 84  | 75  | 68  | 62 | 57 | 53 | 49 | 46 | 43 | 41 | 39 | 37 |    |    |
| 3½ to 4                          | 97  | 87  | 79  | 72 | 66 | 61 | 57 | 53 | 50 | 47 | 45 | 43 |    |    |
| 4 to 4½                          | 110   | 98  | 89  | 82 | 75 | 70 | 65 | 61 | 57 | 54 | 51 | 48 |    |    |
| 4½ to 5                          | 123   | 110 | 100 | 91 | 84 | 77 | 72 | 68 | 63 | 60 | 57 | 54 |    |    |
| 5                                | 129   | 116 | 105 | 96 | 88 | 81 | 76 | 71 | 67 | 63 | 60 | 56 |    |    |

## E SURVEYS OF FOOD CONSUMPTION

### Reference

- 1 Berryman G H and Chatfield C Short Method of Calculating Nutritive Value of Diets *J Nutrition* 25 327 (Jan.) 1943
- 2 Berryman G H and Howe P E Short Method of Calculating Nutritive Value of Diets *J Nutrition* 27 31-40 (Mar.) 1944.

### Principle

In the field it is often necessary to estimate the average food and nutrient consumption of a group. For this purpose a short survey is done over a period of ten days. The total food consumption (by groups) is calculated from inventory data corrected for plate and kitchen waste. The total nutrient consumption is calculated by adding together the nutrients contributed by each food group.

### Apparatus and Supplies

- 1 Scales folding 300 lb capacity
- 2 Hand tally counter
- 3 Clip board
- 4 Ruler
- 5 At least ten 5 gallon cans or equivalent containers for segregation and weighing of waste both kitchen and plate
- 6 Slide rule
- 7 Adding or calculating machine
- 8 Essential books
  - a *Handbook of Chemistry and Physics*
  - b National Research Council Tables of Food Composition
  - c Table for converting weights to pounds per man per day
  - d Bowes and Church or other standard book on food composition
  - e Master Menu
  - f TM 10-412 (Army Recipes)
- 9 Pads of paper 8 x 10½
- 10 Small notebooks one for each observer
- 11 Pencils assorted.
- 1 Paper graph one pad
- 13 Clips paper
- 14 3 x 4½ bibliography cards
- 15 Proper forms
- 16 Rubber bands
- 17 One dozen manila folders 8 x 10½

### Reagents

None unless it is proposed to collect specimens of food. If so see sections on collection of specimens and the individual method contemplated.

### Procedure

- 1 It is recommended that the survey be planned to cover a complete 10 day period. Shorter periods may give misleading and erroneous results.  
Inventory and Issue. A complete inventory is made of food stores in the mess on the night before each survey and on the night after the last meal of each survey. All items used during the survey are weighed and recorded.
- 3 Head Count. At each meal one observer is stationed near the place where trays are returned. Each man is tallied on a hand counter as he hands in his tray. Meals at odd hours especially at night are recorded by night cook and baker.
4. Kitchen and Plate Waste. Food groups are segregated and weighed in the following categories: meats, fish and poultry, eggs, milk products, butter, other

fats sugars and syrups grain products dry legumes and nuts leafy green and yellow vegetables tomatoes citrus fruits white potatoes, other vegetables other fruits dried fruits condiments, and miscellaneous

For the measurement of plate waste a five gallon can is convenient to use for each food group represented in a particular meal. At the end of each meal, each man supervised by an observer segregates his own leftovers according to food groups scraping the food from his tray into the appropriate five gallon can. Kitchen waste is also segregated according to food groups and both plate and kitchen waste are weighed. Trap grease is weighed every day and is counted as fat other

### Calculation

Calculation although tedious is simple in principle

#### 1 Food intake by food groups =

(First inventory + Food issued) minus

(Last inventory + kitchen waste + Plate waste)

The breakdown of kitchen and plate waste into groups is made in accordance with the recipes of the various dishes. The food groups of canned and other items are determined from recipes or standard tables

#### 2 Nutrient intake by individual nutrients =

(Total food intake, by food groups)  $\times$  (Group factor for nutrient)

The food group factors may be calculated from tables of food composition such as Bowes A deP and Church, C F Food Values of Portions Commonly Used. A deP Bowes Philadelphia (1946), and US Department of Agriculture Miscellaneous Publication No 572 Tables of Food Composition in Terms of Eleven Nutrients prepared by the Bureau of Human Nutrition and Home Economics US Department of Agriculture in cooperation with National Research Council, 1945 Table 23 which follows is a table of nutrient conversion factors used by the Medical Nutrition Laboratory in one of its surveys. Cooking losses may be corrected for when the food group factors are computed

#### 3 Food or nutrient intake/man/day =

$$\frac{\text{Total food intake or nutrient intake}}{\text{Average number of men messing} \times \text{total days of survey}}$$

TABLE 23

### NUTRIENTS PER POUND OF FOOD GROUPS AS PURCHASED

(Based on food prescribed for the Army, May October 1944 and NRC Tables of Food, Composition Third Revised Edition.)

| NUTRIENT            | MEATS <sup>1</sup><br>FISH<br>POULTRY | EGGS<br>FPELH | MILK <sup>2</sup><br>MILK<br>PRODUCTS | BUTTER | FATS<br>OTHER | SUGARS<br>SYRUPS | CEREALS<br>GRAIN<br>PRODUCTS |
|---------------------|---------------------------------------|---------------|---------------------------------------|--------|---------------|------------------|------------------------------|
| Calories/lb         | 1158                                  | 636           | 46                                    | 33.7   | 3860          | 1730             | 146                          |
| Protein, gm/lb      | 59                                    | 52            | 21                                    | 3      | 2             | 2                | 49                           |
| Fat, gm/lb          | 102                                   | 46            | 24                                    | 368    | 4.6           | 0.1              | 14                           |
| Carbohydrate gm/lb  | 3                                     | 3             | 31                                    | 2      | 4             | 431              | 85                           |
| Calcium mg/lb       | 41                                    | 218           | 695                                   | 73     | 4             | 96               | 217                          |
| Phosphorus mg/lb    | 668                                   | 818           | 558                                   | 73     | 75            | 13               | 491                          |
| Iron mg/lb          | 9.0                                   | 10.9          | 0.5                                   | 1.0    | 1.2           | 2.2              | 11.9                         |
| Vit. A IU/lb        | 473                                   | 4590          | 1169                                  | 17500  | 263           | 4                | 31                           |
| Thiamin mg/lb       | 0.91                                  | 0.43          | 0.18                                  | 0.01   | 0.04          | 0.01             | 1.60                         |
| Riboflavin, mg/lb   | 0.84                                  | 1.35          | 1.00                                  | 0.05   | 0.04          | 0.01             | 1.70                         |
| Niacin, mg/lb       | 18.0                                  | 0.3           | 0.6                                   | 0.5    | 0             | 0.1              | 13.9                         |
| Ascorbic acid mg/lb | 6.0                                   | 0             | 6.0                                   | 0      | 0             | 2.0              | 0                            |

1. Meat values - a carcass basis lean pork cuts. Calculate them by the long method.
2. Excluding milk products - which is to a straight basis. Do not use milk equivalents for calculation. The values for milk and basis.
3. When liver 8.000 and 12.73 mg riboflavin/lb liver

TABLE 23 (CONT'D)

| NUTRIENT       | BEANS<br>OTHER DRY<br>LEGUMES<br>NUTS | VEGETABLES        |          |                              |       | FRUITS |                |       |
|----------------|---------------------------------------|-------------------|----------|------------------------------|-------|--------|----------------|-------|
|                |                                       | POTATOES<br>WHITE | TOMATOES | LEAFY<br>GREEN AND<br>YELLOW | OTHER | CITRUS | OTHER<br>FRESH | DRIED |
| Cal            | 18.3                                  | 2.5               | 143      | 1.5                          | 18.5  | 15     | 38             | 1300  |
| Prot           | 104                                   | 8                 | 5        | 6                            | 6     | 3      | 2              | 11    |
| Fat            | 55                                    | 0                 | 1        | 1                            | 1     | 1      | 0.8            | 2     |
| CHO            | 33                                    | 73                | 8        | 4                            | 38    | 35     | 56             | 308   |
| Ca             | 530                                   | 4                 | 49       | 1.4                          | 90    | 9      | 41             | 295   |
| P              | 1947                                  | 13                | 115      | 198                          | 16    | 61     | 57             | 516   |
| Fe             | 3.0                                   | 2.7               | 2.8      | 3.6                          | 2     | 11     | 1.5            | 18.2  |
| Vit. A         | 12                                    | 70                | 5160     | 11910                        | 309   | 3.6    | 1610           | 6815  |
| B              | 0.30                                  | 0.40              | 0.26     | 0.25                         | 0.15  | 0.19   | 0.11           | 0.13  |
| B <sub>2</sub> | 1.50                                  | 0.15              | 0.17     | 0.6                          | 0.18  | 0.06   | 0.30           | 0.69  |
| Niacin         | 23.9                                  | 4.4               | 3.9      | 1.7                          | 1.6   | 0.7    | 1.4            | 8.40  |
| C              | 60                                    | 570               | 770      | 610                          | 400   | 1460   | 80             | 80    |

Example

In Table 4 are listed the significant data for a recent survey in Alaska. The step by step calculations are not shown but only the final data for food consumption and nutrient intake.

Interpretation of Data

If it is desired to compare the results with recommended dietary allowances Table 5 shows the NRC recommended allowances taken from National Research Council Food and Nutrition Board Recommended Dietary Allowances Revised 1943 Reprint and Circular Series No. 19 (Oct.) 1948 and the Department of Army standards from U. S. Department of the Army Army Regulation No. 40-130 Medical Department Nutrition Washington, D. C. (23 Oct.) 1947.

Precautions

1. Administrative considerations are important in surveys of food consumption. The kitchen and mess personnel must be cooperative and honest enough to give accurate data of inventory and issue.

TABLE 4

EXAMPLE OF RESULTS OF A TYPICAL SURVEY OF FOOD CONSUMPTION AND NUTRIENT INTAKE IN A SUBARCTIC CAMP

| FOOD GROUP                        | LB. /MAN/DAY | NUTRIENT                      | AMOUNT |
|-----------------------------------|--------------|-------------------------------|--------|
| Meats fish poultry                | 1.8          | Calories                      | 1550   |
| Eggs fresh                        | 0.6          | Protein gm                    | 149    |
| Milk & Milk Products              | 1.13         | Fat gm                        | 2.8    |
| Butter                            | 10           | Carbohydrate gm               | 5.1    |
| Fats other                        | 10           | Calium gm                     | 1.2    |
| Sugars syrups etc                 | 44           | Phosphorus gm                 | 0.3    |
| Cereals                           | 49           | Iron mg                       | 6      |
| Legumes                           | 0.01         | Vitamin A IU                  | 14700  |
| Vegetables leafy green and yellow | 58           | Thiamine <sup>1</sup> mg      | 1.8    |
| Tomatoes                          | 37           | Riboflavin <sup>2</sup> mg    | 3.1    |
| Fruits citrus                     | 0            | Niacin <sup>3</sup> mg        | 0.9    |
| Potatoes white                    | 86           | Ascorbic acid <sup>3</sup> mg | 113    |
| Vegetables other                  | 0            |                               |        |
| Fruits other                      | 51           |                               |        |
| Fruits dried                      | 03           |                               |        |
| Total pounds                      | 660          |                               |        |
| Waste (lbs kitchen plus plate)    | 89           |                               |        |

Figures represent average/man/day not including intake outside the mess.

<sup>1</sup>Cooking losses calculated according to AR 40-150 page 637

fats sugars and syrups grain products, dry legumes and nuts leafy green and yellow vegetables tomatoes citrus fruits, white potatoes other vegetables other fruits, dried fruits condiments, and miscellaneous

For the measurement of plate waste, a five gallon can is convenient to use for each food group represented in a particular meal. At the end of each meal each man supervised by an observer segregates his own leftovers according to food groups scraping the food from his tray into the appropriate five gallon can. Kitchen waste is also segregated according to food groups and both plate and kitchen waste are weighed. Trap grease is weighed every day, and is counted as fat other.

### Calculation

Calculation although tedious is simple in principle

#### 1 Food intake, by food groups =

(First inventory + Food issued) minus  
(Last inventory + Kitchen waste + Plate waste)

The breakdown of kitchen and plate waste into groups is made in accordance with the recipes of the various dishes. The food groups of canned and other items are determined from recipes or standard tables.

#### 2 Nutrient intake by individual nutrients =

(Total food intake by food groups)  $\times$  (Group factor for nutrient)

The food group factors may be calculated from tables of food composition such as Bowes A. deP. and Church, C. F. Food Values of Portions Commonly Used A. deP. Bowes Philadelphia (1946), and US Department of Agriculture Miscellaneous Publication No. 572. Tables of Food Composition in Terms of Eleven Nutrients prepared by the Bureau of Human Nutrition and Home Economics, US Department of Agriculture in cooperation with National Research Council, 1945. Table 23 which follows is a table of nutrient conversion factors used by the Medical Nutrition Laboratory in one of its surveys. Cooking losses may be corrected for when the food group factors are computed.

#### 3 Food or nutrient intake/man/day =

$$\frac{\text{Total food intake or nutrient intake}}{\text{Average number of men messing} \times \text{total days of survey}}$$

TABLE 23

#### NUTRIENTS PER POUND OF FOOD GROUPS AS PURCHASED

(Based on food prescribed for the Army May October, 1944 and NRC Tables of Food, Composition Third Revised Edition.)

| NUTRIENT            | MEATS <sup>1</sup><br>FISH<br>POULTRY | EGGS<br>FRESH | MILK <sup>2</sup><br>MILK<br>PRODUCTS | BUTTER | FATS<br>OTHER | SUGARS<br>SYRUPS | CEREALS<br>GRAIN<br>PRODUCTS |
|---------------------|---------------------------------------|---------------|---------------------------------------|--------|---------------|------------------|------------------------------|
| Calories/lb         | 1158                                  | 636           | 46                                    | 3327   | 3860          | 1730             | 146                          |
| Protein gm/lb       | 59                                    | 52            | 21                                    | 3      | 2             | 2                | 49                           |
| Fat, gm/lb          | 102                                   | 46            | 24                                    | 363    | 46            | 0.1              | 14                           |
| Carbohydrate gm/lb  | 3                                     | 3             | 31                                    | 2      | 4             | 431              | 85                           |
| Calcium mg/lb       | 41                                    | 218           | 695                                   | 73     | 75            | 76               | 17                           |
| Phosphorus mg/lb    | 666                                   | 848           | 503                                   | 73     | 75            | 13               | 491                          |
| Iron mg/lb          | 90                                    | 10.9          | 0.5                                   | 10     | 12            | 2                | 11.9                         |
| Vit. A IU/lb *      | 473                                   | 4590          | 1169                                  | 17500  | 263           | 4                | 31                           |
| Thiamin, mg/lb      | 0.91                                  | 0.43          | 0.18                                  | 0.01   | 0.04          | 0.01             | 1.60                         |
| Riboflavin mg/lb *  | 0.84                                  | 1.35          | 1.00                                  | 0.05   | 0.04          | 0.01             | 1.0                          |
| Niacin mg/lb        | 180                                   | 0.3           | 0.6                                   | 0.5    | 0             | 0.1              | 13.9                         |
| Ascorbic acid mg/lb | 60                                    | 0             | 60                                    | 0      | 0             | 20               | 0                            |

- 1 Meat values exclude carcass beef and all lean pork cuts. Calculate them by the long method.
- 2 Excluding cheese which is to be calculated by the long method. The values for milk and milk products are on a straight poundage basis. Do not use milk equivalents for calculations.
- 3 When liver is used add 87000 IU Vitamin A and 1.73 mg riboflavin/lb liver.

TABLE 23 (CONT'D)

| NUTRIENT       | BEANS<br>OTHER DRY<br>LEGUMES<br>NUTS | VEGETABLES        |          |                              |       | FRUITS |                |       |
|----------------|---------------------------------------|-------------------|----------|------------------------------|-------|--------|----------------|-------|
|                |                                       | POTATOES<br>WHITE | TOMATOES | LEAFY<br>GREEN AND<br>YELLOW | OTHER | CITRUS | OTHER<br>FRESH | DRIED |
| Cal            | 183                                   | 35                | 143      | 18                           | 185   | 15     | 38             | 1310  |
| Prot           | 104                                   | 8                 | 5        | 6                            | 6     | 3      | ~              | 11    |
| Fat            | 56                                    | 0                 | 1        | 1                            | 1     | 1      | 0.8            | 2     |
| CHO            | 738                                   | 73                | ~8       | 24                           | 38    | 35     | 56             | 308   |
| Ca             | 530                                   | 4~                | 49       | 124                          | 90    | 79     | 41             | 95    |
| P              | 1967                                  | 213               | 115      | 124                          | 167   | 61     | 57             | 516   |
| Fe             | 350                                   | 27                | ~8       | 36                           |       | 11     | 15             | 182   |
| Vit. A         | 1                                     | 70                | 5160     | 11910                        | 308   | 356    | 1610           | 6915  |
| B              | ~30                                   | 0.40              | 0.6      | 0.25                         | 0.15  | 0.19   | 0.11           | 0.13  |
| B <sub>2</sub> | 150                                   | 0.15              | 0.17     | 0.6                          | 0.18  | 0.06   | 0.30           | 0.69  |
| Niacin         | 239                                   | 44                | 39       | 17                           | 16    | 07     | 14             | 840   |
| C              | 60                                    | 570               | 770      | 610                          | 400   | 1460   | ~80            | ~80   |

#### Example

In Table 24 are listed the significant data for a recent survey in Alaska. The step by step calculations are not shown but only the final data for food consumption and nutrient intake

#### Interpretation of Data

If it is desired to compare the results with recommended dietary allowances, Table 25 shows the NRC recommended allowances taken from National Research Council Food and Nutrition Board Recommended Dietary Allowances Revised 1948 Reprint and Circular Series No. 19 (Oct.) 1948 and the Department of Army standards from U. S. Department of the Army Army Regulation No. 40-20 Medical Department Nutrition Washington 23 D C (8 Oct.) 1947

#### Precautions

1. Administrative considerations are important in surveys of food consumption. The kitchen and mess personnel must be cooperative and honest enough to give accurate data of inventory and issue

TABLE 4

EXAMPLE OF RESULTS OF A TYPICAL SURVEY OF FOOD CONSUMPTION AND NUTRIENT INTAKE IN A SUBARCTIC CAMP

| FOOD GROUP                        | LBS/MAN/DAY | NUTRIENT         | AMOUNT |
|-----------------------------------|-------------|------------------|--------|
| Meats fresh poultry               | 1.8         | Calories         | 3150   |
| Eggs fresh                        | 6           | Protein gm       | 149    |
| Milk & Milk Products              | 1.13        | Fat gm           | 76     |
| Butter                            | 10          | Carbohydrate gm  | 51     |
| Fats other                        | 10          | Calcium gm       | 12     |
| Sugars syrups etc.                | 44          | Phosphorus gm    | 23     |
| Cereals                           | 49          | Iron mg          | ~6     |
| Legumes                           | 0.01        | Vitamin A IU     | 14.00  |
| Vegetables leafy green and yellow | .58         | Thiamine mg      | 18     |
| Tomatoes                          | 37          | Riboflavin mg    | 31     |
| Fruit citrus                      | ~           | Niacin mg        | 29     |
| Potatoes white                    | 86          | Ascorbic acid mg | 113    |
| Vegetables other                  | ~           |                  |        |
| Fruits other                      | .51         |                  |        |
| Fruits dried                      | .03         |                  |        |
| Total pounds                      | 660         |                  |        |
| Waste (lbs kitchen plus plate)    | 89          |                  |        |

Figures represent average/man/day not including intake outside the mess.

\*Cooking losses calculated according to ARJ 40-250 page 637

- 2 The trap grease is important being pure fat  
 3 The results at best must be regarded as rough estimates This is particularly true of nutrient intake

TABLE 25

RECOMMENDED DAILY DIETARY ALLOWANCES OF THE NATIONAL RESEARCH COUNCIL, FOOD AND NUTRITION BOARD AND ALLOWANCES PRESCRIBED BY AR 40 250

| NUTRIENT         | FOOD AND NUTRITION BOARD |           |        |               | AR 40 250 |        |
|------------------|--------------------------|-----------|--------|---------------|-----------|--------|
|                  | BOYS<br>16 20            | MEN       |        |               | SEDENTARY | ACTIVE |
|                  |                          | SEDENTARY | ACTIVE | HEAVY<br>WORK |           |        |
| Calories†        | 3800                     | 4000      | 3000   | 4500          | 3000      | 3600   |
| Protein gm       | 100                      | 70        | 70     | 70            | 100       |        |
| Calcium gm       | 14                       | 10        | 10     | 10            | 07        |        |
| Iron gm          | 15                       | 12        | 12     | 12            |           |        |
| Vit A IU         | 6000                     | 5000      | 5000   | 5000          | 5000      |        |
| Thiamine mg      | 17                       | 12        | 15     | 18            | 16        |        |
| Riboflavin mg    | 25                       | 18        | 18     | 18            | 2.2       |        |
| Niacin mg        | 17                       | 12        | 15     | 18            | 16        |        |
| Ascorbic acid mg | 100                      | 75        | 75     | 75            | 50        |        |
| Vit D IU         | 400                      | -         | -      | -             | -         |        |

Calorie allowances must be adjusted up or down to meet specific needs The proper calorie allowance is that which over an extended period of time will maintain body weight or rate of growth at the level most conducive to well being AR 40 250 prescribes 4100 Calories when prevailing temperatures are sub zero for more than 1 week.

## F HEAT BALANCE—CLOTHING

### Reference

- Belding, H S, Darling, R C, Robinson, S Turrell E and Griffin, D E Notes on Methods Used at the Fatigue Laboratory to Evaluate the Thermal Insulation Provided by Cold Weather Clothing, Report No 19 under OSRD Contract No OEM CMR 54 Sept 10 1943  
 Belding H S Chapter in Clothing Test Methods Edited by Newburgh L H, and Harris, M 1945  
 Lee D H K. and Lewsons, H Clothing for Global Man, Geographical Review 29 181-113 (No 2) 1943

### Importance of Clothing in Field Research

In any study of man in relation to his environment, heat balance is important, and clothing usually has to be taken into account in military studies Clothing affects heat balance through its effects on heat exchange by radiation conduction convection and evaporation convection In addition clothing may affect heat production by its weight or restrictive effects

In relation to radiation effects the clothing may be considered as a barrier between the body and the environment Clothing shuts off about 95% of the radiation from the body which emits in the far infrared Clothing reflects a variable proportion of light incident upon it

The thermal resistance of cloth to conduction is for the most part proportional to the amount of still air held in the substance of the cloth and between its layers For most purposes the nature of the fabric may be neglected and the thermal insulation may be considered as proportional to the total thickness of the clothing assembly Direct conduction to an external object may be increased by compression and wetting of the clothing





<sup>2</sup> The trip gressio is important being pure fat

<sup>3</sup> The results at best must be regarded as rough estimates This is particularly true of nutrient intake

TABLE 25

RECOMMENDED DAILY DIETARY ALLOWANCES OF THE NATIONAL RESEARCH COUNCIL FOOD AND NUTRITION BOARD AND ALLOWANCES PRESCRIBED BY AR 40 250

| NUTRIENT         | FOOD AND NUTRITION BOARD |           |        |               | AR 40 250 |        |
|------------------|--------------------------|-----------|--------|---------------|-----------|--------|
|                  | BOYS<br>16 00            | MEN       |        |               | SEDENTARY | ACTIVE |
|                  |                          | SEDENTARY | ACTIVE | HEAVY<br>WORK |           |        |
| Calories         | 3500                     | 4000      | 3000   | 4500          | 3000      | 3600   |
| Protein gm       | 100                      | 70        | 70     | 70            | 100       |        |
| Calcium gm       | 14                       | 10        | 10     | 10            |           | 07     |
| Iron, gm         | 15                       | 12        | 12     | 12            |           |        |
| Vit A IU         | 6000                     | 5000      | 5000   | 5000          | 5000      |        |
| Thiamine, mg     | 17                       | 17        | 15     | 18            |           | 18     |
| Riboflavin mg    | 5                        | 18        | 18     | 18            |           | 22     |
| Niacin mg        | 17                       | 13        | 15     | 18            |           | 16     |
| Ascorbic acid mg | 100                      | 75        | 75     | 75            | 50        |        |
| Vit D IU         | 400                      | -         | -      | -             | -         |        |

<sup>1</sup>Calorie allowances must be adjusted up or down to meet specific needs The proper calorie allowance is that which over an extended period of time will maintain body weight or rate of growth at the level most conducive to well being AR 40 250 prescribes 4400 Calories when prevailing temperatures are sub zero for more than 1 week

## F HEAT BALANCE—CLOTHING

### Reference

Belding H S, Darling, R C, Robinson S, Turrell E, and Griffin, D R Notes on Methods Used at the Fatigue Laboratory to Evaluate the Thermal Insulation Provided by Cold Weather Clothing, Report No 19 under OSRD Contract No OEM CMR 54, Sept 10, 1943

Belding, H S Chapter in *Clothing Test Methods*, Edited by Newburgh L H, and Harris M 1945

Lee, D H H and Lemons H. *Clothing for Global Man* Geographical Review 23 181 213 (No 2) 1949

### Importance of Clothing in Field Research

In any study of man in relation to his environment heat balance is important and clothing usually has to be taken into account in military studies Clothing affects heat balance through its effects on heat exchange by radiation conduction convection and evaporation convection In addition, clothing may affect heat production by its weight or restrictive effects

In relation to radiation effects the clothing may be considered as a barrier between the body and the environment Clothing shuts off about 95% of the radiation from the body which emits in the far infrared Clothing reflects a variable proportion of light incident upon it

The thermal resistance of cloth to conduction is for the most part proportional to the amount of still air held in the substance of the cloth and between its layers For most purposes the nature of the fabric may be neglected, and the thermal insulation may be considered as proportional to the total thickness of the clothing assembly Direct conduction to an external object may be increased by compression and wetting of the clothing

When the ambient temperature is above skin temperature evaporation is the only avenue of heat loss. At all ambient temperatures evaporation must be taken into account. When the skin is naked heat of vaporization is drawn from the skin itself and the effectiveness of such evaporation is nearly 100%. When clothing is added to the situation there are complications. Under hot conditions other avenues of heat are altered; the efficiency of convection over the skin is reduced; the sweat is evaporated away from the skin; part of the heat of vaporization being taken from the ambient air; the rate of effective evaporation (i.e. the evaporation that cools the skin) is reduced; and when the whole body has been covered with a film of moisture evaporative cooling has reached its maximum. In cold conditions water vapor may recondense in some layer of the clothing assembly, returning some of its heat to the clothing system and increasing the thermal conductivity of the fabric by wetting it. Under all conditions clothing impedes the diffusion of water vapor outward from the skin; textiles differing in this respect one from the other; and the movements of the body increase the transfer of heat in the layers of the clothing assembly by a bellows action.

#### The Thermal Insulation of Clothes—The Clo Unit

For historical reasons the Clo Unit is important since much of the data on clothing were expressed in terms of Clo units. The unit is defined as the thermal insulation necessary to maintain in comfort a sitting resting subject in a normally ventilated room (air movement 0 feet per minute) at a temperature of 70 F and a humidity of less than 50% (Gagge, A. I. and associates, *A Practical System of Units for the Description of the Heat Exchange of Man with His Environment*, Science 94: 48-430, June 1941). One Clo unit equals 0.18 C/kilocal/square meter/hour.

Experimentally the thermal insulation of a clothing assembly on a man is calculated from the equation

$$\text{Clo value} = \frac{3.09 (T_s - T_a)}{H_c} - I_a \quad (\text{Equation 1})$$

Where  $T_s$  is average skin temperature F

$T_a$  is ambient temperature F

$H_c$  is heat lost through the clothing and

$I_a$  is the insulation of the air in Clo units at the wind velocity of the test.

$T_s$  The mean skin temperature is determined by thermocouple measurements of the following areas: front of chest; lateral aspect of upper arm; medial aspect of first phalanx of great toe; medial aspect of thigh; medial aspect of lower arm; over shoulder blade; over gastrocnemius; edge of hand at base of first finger and forehead.

$I_a$  the insulation of air at the wind velocity of the test is obtained from tables compiled by A. C. Burton.

$$H_c = M - S - (F + A) \quad (\text{Equation 2})$$

where  $M$  is metabolic heat production;  $S$  is heat stored or lost in the body mass during the experiment;  $E$  is heat lost by vaporization from lungs and skin; and  $A$  is heat lost in warming inspired air.

$M$  is estimated by standard procedures as described in the section on respiratory methods.  $S$  is estimated by direct observation of the change in skin and internal temperatures and calculated according to Burton's equations.

(1) Mean temperature of body mass =

$$0.53 \times \text{mean skin temperature} + 0.67 \times \text{rectal temperature} \quad (\text{Equation 3})$$

$$(2) \text{ Stored heat} = \Delta \text{ mean temperature} \times 0.83 \quad (\text{Equation 4})$$

0.83 being the specific heat of the body in Cals/kg/°C and  $\Delta$  mean temperature is in °C

E is estimated from the change in the subject's weight (to the nearest gram) corrected for the excess weight of CO<sub>2</sub> expired compared with O consumed according to the equation

$$\frac{\text{Cals/square meter/hr vaporized from skin and lungs} = \frac{(\text{gm change in weight/hr}) - (0.30 \times O_2 \text{ l/hr})}{\text{surface area, square meters}} \times 0.586 \quad (\text{Equation 5})$$

where  $0.30 \times O_2$  l/hr represents the excess weight of CO<sub>2</sub> expired over O consumed at R Q 0.88 and 0.586 is the heat of vaporization of water

A is calculated from pulmonary ventilation as measured by techniques described in previous sections

$$A \text{ (cals/square meter/hour from warming inspired air)} = \frac{\text{pulmonary ventilation l/hr} \times 0.00031 \times (33^\circ\text{C} - \text{ambient T } ^\circ\text{C})}{\text{surface area square meters}}$$

where 0.00031 is the specific heat of air, kg Cal/cubic meter/°C

#### Thermal Significance of Clothing

A generalization was made by Burton to express the effect of clothing on heat loss from the surface of the body to a given environment (Burton A C An Analysis of the Physiological Effects of Clothing in Hot Environments Report to Aviation Medical Research Council of Canada C 2754 SPC 186, Nov 24 1944)

His equation is

$$H \text{ heat loss in kilocal/square meter/hour} = \frac{5 \times W \times (P - P_a)}{100 \times (r_1 + r)} + \frac{C(t - t_a - RI/^\circ\text{C})}{I_1 + I} \quad (\text{Equation 6})$$

where

W is wetted skin

P is vapor pressure of air in mm Hg at skin temperature obtained from tables,

P<sub>a</sub> is vapor pressure of ambient air in mm Hg from standard tables

r<sub>1</sub> and r are respectively the resistance of clothing and air expressed as cm of still air r<sub>1</sub> is determined experimentally or calculated from the measure thickness of the clothing r from a standard curve which gives the value corrected for the estimated radius of curvature and velocity of ambient air movement

C is a constant 5.55 for one °C

t and t<sub>a</sub> are respectively the mean skin temperature and ambient air temperature Measurement of mean skin temperature is described above under 'Clo' unit measurement

I<sub>1</sub> and I are respectively the insulation value of clothing and air, expressed in 'Clo' units I<sub>1</sub> is determined experimentally as described in the section on 'Clo' units or is calculated at the same time as r<sub>1</sub> I is determined from a standard curve R is absorbed radiant energy in kilocal/square meter/hour

Burton's equation has three limitations it is applicable only to equilibrium conditions it makes no allowances for subclothing convection, and some of the values are difficult to measure or calculate especially the radiant energy factor It is very useful for predicting the value of various clothing assemblies or predicting the thermal state of a man wearing a given assembly in a given environment

## G METEOROLOGICAL MEASUREMENTS

In all field studies certain meteorological data should be collected systematically and regularly. These include the altitude of the station its latitude and longitude barometric pressure maximum minimum and ambient temperatures precipitation hours of sunshine cloud conditions absolute and relative humidities and wind velocity. The charts in the present section are designed to assist in the calculation of such data.

### Temperature

The first chart (Fig. 6) gives direct conversion of Centigrade temperatures into Fahrenheit temperatures and vice versa according to the equation

$$^{\circ}\text{F} = (9/5 \times ^{\circ}\text{C}) + 32$$

$$^{\circ}\text{C} = (5/9 \times ^{\circ}\text{F}) - 32$$

### Wind Velocity

The second chart (Fig. 57) gives conversion for temperature barometric pressure and wind velocity. Wind velocity is converted from meters/second to miles/hour by the equations

$$^{\circ}\text{ miles/hour} = 360 \times (\text{meters per second})$$

$$^{\circ}\text{ meters/second} = 0.4470 \times (\text{miles per hour})$$

### Barometric Pressure Altitude and Related Variables

The third chart (Fig. 58) gives barometric conversion from millimeters of mercury to inches of mercury according to the equations

$$^{\circ}\text{ millimeters Hg} = 25.40 \times (\text{inches Hg})$$

$$^{\circ}\text{ inches Hg} = 0.03937 \times (\text{millimeters Hg})$$

The third chart gives various altitude equivalents. From Table 30 Altitude Pressure Tables Based on the United States Standard Atmosphere National Advisory Committee for Aeronautics Report No. 538 1935 are given the following altitude correlated equivalents altitude (feet and meters) barometric pressure (inches and millimeters of mercury) temperature (actual and mean) and oxygen (millimeters partial pressure and equivalent sea level  $\text{O}_2$ ).

### Humidity and Absolute Vapor Tension

The fourth chart (Fig. 59) gives the relative humidity from the dry bulb temperature and the wet bulb temperatures. The fifth chart (Fig. 60) for use in cold areas gives relative humidity from the dry bulb temperature and the difference between the dry bulb and the wet bulb temperatures.

It is sometimes necessary to compute absolute vapor tension instead of relative humidity. A line chart (Fig. 61) has been constructed for this purpose from the equation

$$E (\text{absolute vapor tension in mm Hg}) =$$

$$E_w + 0.0066 B (T_d - T_w) [1 + 0.00115 (T_d - T_w)]$$

where  $E_w$  is vapor tension mm Hg at wet bulb temperature

$B$  is barometric pressure mm Hg

$T_d$  is dry bulb temperature  $^{\circ}\text{C}$  and

$T_w$  is wet bulb temperature  $^{\circ}\text{C}$ .

The line chart gives vapor tension for various barometric pressure together with the vapor tension from the wet bulb temperature and the difference between the wet bulb and dry bulb temperatures. It may also be used to compute relative humidity as follows

# TEMPERATURE CONVERSION (CENTIGRADE - FAHRENHEIT)






Fig 56—Line chart for converting Centigrade and Fahrenheit scales

# CONVERSIONS, METEOROLOGY

## TEMPERATURE



## BAROMETER



## WIND VELOCITY




Fig. 57.—Conversion charts for temperature, barometric pressure and wind velocity

## ALTITUDE EQUIVALENTS

| ALTITUDE |        | BAROMETER |       | TEMPERATURE |         | OXYGEN  |       |
|----------|--------|-----------|-------|-------------|---------|---------|-------|
| FEET     | METERS | INCHES    | mm Hg | ACTUAL °C   | MEAN °C | EQUIV % | mm Hg |
| 45000    | 14000  | 4         | 100   |             | -30     | 3       |       |
|          |        |           | 110   |             |         |         | 25    |
|          | 13000  | 5         | 120   |             |         |         |       |
|          |        |           | 130   |             |         |         |       |
| 40000    | 12000  | 6         | 140   | -55         | -25     | 4       | 30    |
|          |        |           | 150   |             |         |         |       |
|          | 11000  | 7         | 160   |             |         |         | 35    |
|          |        |           | 170   |             |         |         |       |
| 35000    | 10000  | 8         | 180   | -50         | -20     | 5       | 40    |
|          |        |           | 190   |             |         |         |       |
|          |        |           | 200   |             |         |         |       |
| 30000    | 9000   | 9         |       | -45         | -15     | 6       | 45    |
|          |        |           |       |             |         |         | 50    |
|          | 8000   | 10        | 250   | -40         |         | 7       | 55    |
| 25000    | 7000   | 11        |       | -35         | -10     | 8       | 60    |
|          |        |           |       |             |         |         | 65    |
|          | 6000   | 12        | 300   | -30         |         | 9       | 70    |
| 20000    |        | 13        |       |             |         |         | 75    |
|          |        | 14        | 350   | -25         | -5      | 10      | 80    |
|          |        | 15        |       | -20         |         | 11      | 85    |
|          | 5000   | 16        | 400   | -15         | 0       | 12      | 90    |
| 15000    |        | 17        |       |             |         |         | 95    |
|          | 4000   | 18        | 450   | -10         |         | 13      | 100   |
|          |        | 19        |       |             |         |         | 105   |
|          |        | 20        | 500   | -5          | +5      | 14      | 110   |
| 10000    | 3000   | 21        |       |             |         |         | 115   |
|          |        | 22        | 550   | 0           |         | 15      | 120   |
|          |        | 23        |       |             |         |         | 125   |
|          | 2000   | 24        | 600   | +5          | +10     | 17      | 130   |
| 5000     |        | 25        |       |             |         |         | 135   |
|          |        | 26        | 650   | +10         |         | 18      | 140   |
|          | 1000   | 27        |       |             |         |         | 145   |
|          |        | 28        | 700   |             |         | 20      | 150   |
| 0        | 0      |           |       | +15         | +15     |         |       |

C. E. SIMON

Fig. 58—Chart constructed for correlated data at different altitudes.




Fig. 58.—Cartesian graph for relative humidity from psychrometric data.

—G. E. Horton—



## RELATIVE HUMIDITY BY SLING PSYCHROMETER





Fig. 60—Line chart for relative humidity

# ABSOLUTE VAPOR TENSION FROM PSYCHROMETRY



11 5000

Fig. 11—Line chart for computing absolute vapor tension. Measure barometric pressure and dry bulb temperature. (a) Draw line between wet and dry bulb temperatures on chart. (b) Stretch ruler between wet bulb line (L) and dry bulb minus wet bulb line (R). (c) Note the point where the ruler intersects the barometric pressure vertical line. (d) Turn ruler to run between (c) on the right line and the barometric pressure point. (e) Read absolute vapor tension on center line.

# "WINDCHILL"

$$K_0 = \left( \sqrt{WV \times 100 - WV + 10.5} \right) (33 - T_A)$$




Fig. 6.—Line chart showing "wind chill" and state of comfort.

(a) determine absolute vapor tension, (b) determine vapor tension corresponding to saturation at the dry bulb temperature using the wet bulb scale and a wet bulb-dry bulb difference

of  $\theta$  (c) calculate relative humidity R.H. in %  $= \frac{\text{absolute vapor tension}}{\text{vapor tension dry bulb equivalent}} \times 100$

#### Windchill

The concept of windchill has been very useful in expressing the relative discomfort of cold climates in relation to absolute temperature and wind velocity. The quantitative aspects of this concept have been discussed by Siple and Passel (Siple P. A. and Passel C. F. Measurements of Dry Atmospheric Cooling in Subfreezing Temperatures. *Proc Am Phil Soc* 29 1:710 1949) and Court [Court Arnold Windchill Bull. Am. Meteorol. Soc. 29 437-493 (Dec) 1949]. The mathematical validity of the equation of Siple and Passel is questionable but nevertheless the concept emphasizes the most important single avenue of body heat loss under windy conditions—convection to the air.

The accompanying line chart (Fig. 6) enables one to derive windchill from wind velocity and ambient temperature and to estimate at the same time the degree of discomfort of the environment. The equation used was that of Siple and Passel:

Windchill kologram Cal/square meter/hour =

$$(\sqrt{WV \times 100} - W + 10.45) \times (33 - T)$$

where WV is wind velocity meters/second

10.45 is an arbitrary constant

33 C is the average skin temperature and

T is the ambient dry bulb temperature C

## H PRESERVATION OF ANIMAL TISSUES FOR HISTOLOGICAL STUDIES

(This section was prepared by K. Sutherland and C. L. Pirani.)

#### References

- Lillie R. D. *Histopathologic Technique* ed. 6 Philadelphia 1948 The Blakiston Company
- Glick D. *Techniques of Histo and Cytochemistry* ed. 1 New York 1949 Interscience Publishers Inc.
- U. S. War Dept. *Methods for Laboratory Technicians TM 8-40* War Department Technical Manual 1946
- Cowdry P. V. *Laboratory Technique in Biology and Medicine* ed. 6 Baltimore Md. 1948 Williams and Wilkins Company
- Gersh I. *Application in Pathology of the Method of Fixation by Freezing and Drying of Tissues* Bull. Internat. Assoc. Med. Museums 23 1:9 195 1948

During field studies it may be necessary or desirable to study tissues from either animal or human sources. The tissues may have been obtained at autopsy or by biopsic or surgical methods. If proper precautions in the preservation and fixation of tissue specimens are observed histological, histochemical and chemical studies may be carried out and contribute substantially to the understanding of physiological and physiopathological processes. Application of such techniques undoubtedly could be very useful in metabolic and nutritional studies that might be conducted in the field.

It is beyond the purpose of this manual however to give detailed techniques which can be applied to investigations of this type. It is obvious that

in most cases equipment necessary to carry out these methods will not be available in field studies. The reader is referred to a number of excellent manuals and textbooks which are listed among the references for detailed histological and histochemical techniques. This section will therefore be limited to a brief review of the best available methods to preserve tissues which can be applied to field conditions and which will then permit either storage or shipment of tissue specimens.

*1 Preservation of Tissue Specimens for Histological Studies* Probably the best fixative for all around purposes is 10 per cent formalin. This fixative is prepared by diluting 1 volume of 40% Formaldehyde (USP) with 9 volumes of distilled water. If distilled water is not available tap water can be used satisfactorily.

Blocks of tissue should not be more than 0.4 to 0.5 cm in thickness and 2 to 3 cm in width to insure proper fixation and to prevent distortion. The proportion of the formalin solution to volume of tissue to be fixed should be no less than 10:1. In addition it is advisable to change the formalin solution after 24 hours. If these simple precautions are used tissues can be kept indefinitely in this fixative.

Formalin fixation can be used satisfactorily for most routine histological methods. It is not adequate however for many histochemical techniques except those for non water soluble substances such as lipids. According to Lallie proper treatment of formalin fixed tissue makes possible the use of additional special staining techniques.

In general the use of other fixatives is not advisable under field laboratory conditions unless tissues can be shipped to reach a processing laboratory within 24 hours. An important exception is 80% alcohol in which tissues can be kept indefinitely. In addition this fixative offers the advantage of freezing at a temperature lower than that of 10% formalin.

*2 Preservation of Tissue Specimens for Histochemical and Chemical Studies* Deep freezing of tissue specimens is best suited for this purpose. This can be accomplished either by storing the specimens in a deep freezing apparatus or by placing them in a thin walled container which is then surrounded by crushed dry ice. Tissues can be kept indefinitely in a deep freeze state without any significant chemical change taking place. In shipment of frozen specimens is desired the following precautions are necessary: (a) Place specimens in an insulated container such as a vacuum bottle and pack dry ice all around them. (b) The volume of dry ice should be at least 20 times greater than that of the specimen. This will insure its preservation for a period of no more than 24 hours. (c) Make sure that the container is not gas tight and will allow the expanding carbon dioxide to escape.

Defrosting of frozen tissue specimens at room temperature is satisfactory for most chemical analysis except when labile substances (for example ascorbic acid and some enzymes) are involved. In such case adequate precautions should be used as indicated by the different analytical methods.

For histochemical studies neither slow nor quick defrosting can prevent some degree of distortion and loss of histological detail. Paraffin embedding of frozen specimens under vacuum in the Altman Gersh apparatus is without doubt the method of choice should such apparatus be available. With this technique tissues undergo practically no chemical or histological changes and the blocks obtained can be permanently stored for future slides.

3 *Labeling and Packing of Tissue Specimens* Wide mouthed bottles sealed with paraffin should be generally used for tissues in fixing solution (wet tissue). The bottle should be adequately labeled on the outside. A second small label of heavy drawing paper written in soft pencil (avoid ink) should be put in the bottle with the specimen. Adequate packing into a mailing case marked *Fragile First Class Mail Rush Specimen for Diagnosis* is satisfactory for shipment. Reference is made to the Technical Manual 8 227 of the War Department for additional information on this important phase of the handling of specimens.

## SECTION VIII

# TECHNIQUES FOR METABOLIC WARDS

---

### A GENERAL CONSIDERATIONS

#### References

Reifenstein E C, Jr, Albright F and Wells S L The Accumulation, Interpretation and Presentation of Data Pertaining to Metabolic Balances Notably Those of Calcium Phosphorus and Nitrogen *J Clin Endocrinol* 5 367 395 (Nov) 1945

#### The Physical Arrangement of a Metabolic Ward

Numerous authors have discussed the ideal arrangement of a metabolic ward The paper of Reifenstein Albright and Wells gives a good bibliography on this subject Our own experience has been that the work is greatly facilitated by a well organized physical plant Ideally the patients should have private rooms and in addition a bathroom and common recreation room On the same floor and adjacent to the ward should be

- 1 A dining room
- 2 A diet kitchen
- 3 An office for the staff
- 4 A record office with a calculating machine
- 5 A room for preparing and storing specimens
- 6 A laboratory for routine clinical pathology
- 7 A separate laboratory for the many special procedures required in metabolic work

All space should be air conditioned because it is becoming increasingly clear that environmental control is extremely desirable for the most consistent results In addition, the comfort of the staff which is usually on a 12 month schedule, is very important for efficient work

#### Finances of a Metabolic Ward

Few persons not actively engaged in metabolic work realize the large expenses involved For patient care alone in 1949 the cost in Chicago was of the order of \$14 per patient per day Add to this the numerous and some times expensive laboratory procedures and the cost is such that metabolic work should not be started without an assured adequate budget

#### Collection Preparation and Storage of Specimens

These topics have been discussed in previous parts of this manual If at all feasible it is better to use analytical data for food consumption rather than calculated data Our experience with stool nitrogen is that for a given patient

within wide limits the daily stool nitrogen is quite constant but that different patients have greatly different characteristic excretions. Therefore it is unsound to use calculated values for stool nitrogen.

## B. DERIVED DATA— THEORETICAL BALANCES

### 1. Calcium, Phosphorus and Nitrogen

The concept of theoretical balances has been extensively used by Albright and colleagues (see reference previously cited). By a theoretical balance is meant the balance of an element such as phosphorus which would be expected from the observed balances of other elements such as nitrogen and calcium. The two chief virtues of such derived data are their assistance in detecting gross errors in the data and their emphasis on situations where the measured balances differ from the theoretical.

For computation of theoretical calcium, phosphorus and nitrogen balances one from the other the following factors are used:

$$\begin{aligned}\text{Calcium/phosphorus in bone} &= .23 \\ \text{Nitrogen/phosphorus in muscle} &= 14.7/15.5\end{aligned}$$

In computing the balances the following equations are used assuming that muscle is characteristic of total protoplasmic mass:

- a. Theoretical phosphorus from calcium and nitrogen =

$$\frac{\text{gm calcium retained}}{3} + \frac{\text{gm nitrogen retained}}{15}$$

- b. Theoretical calcium from phosphorus =  
(gm phosphorus retained  $\times$  .23)

- c. Theoretical nitrogen from phosphorus and calcium =

$$(\text{gm phosphorus retained} \times 14.7) - \left\{ \frac{\text{gm calcium retained}}{3} \times 14.7 \right\}$$

### 2. Potassium and Sulfur

From potassium, sulfur and nitrogen it is possible to calculate theoretical balances. The factors used are:

$$\begin{aligned}\text{mEq potassium} &= 2.1 \times \text{gm nitrogen} \\ \text{mg potassium} &= 10 \times \text{gm nitrogen} \\ \text{gm sulfur} &= \frac{\text{gm nitrogen}}{14.5}\end{aligned}$$

To calculate the theoretical balances the following equations are used (note that in these cases the results are not additive as they are for Ca, P and N):

- a. Theoretical nitrogen from potassium =  $\frac{\text{mEq K retained}}{7}$  or  $\frac{\text{mg K retained}}{10}$   
 b. Theoretical nitrogen from sulfur = gm S retained  $\times$  14.5  
 c. Theoretical mEq potassium from nitrogen = gm N retained  $\times$  2.7  
 d. Theoretical mg potassium from nitrogen = gm N retained  $\times$  105  
 e. Theoretical mEq potassium from sulfur = gm S retained  $\times$  14.5  $\times$  .7



## SECTION VIII

# TECHNIQUES FOR METABOLIC WARDS

---

### A GENERAL CONSIDERATIONS

#### References

- Reifenstein E C, Jr Albright F and Wells S L The Accumulation, Interpretation and Presentation of Data Pertaining to Metabolic Balance: Notably Those of Calcium Phosphorus and Nitrogen *J Clin Endocrinol* 5 367-395 (Nov) 1945

#### The Physical Arrangement of a Metabolic Ward

Numerous authors have discussed the ideal arrangement of a metabolic ward. The paper of Reifenstein, Albright and Wells gives a good bibliography on this subject. Our own experience has been that the work is greatly facilitated by a well organized physical plant. Ideally the patients should have private rooms and in addition a bathroom and common recreation room. On the same floor and adjacent to the ward should be

- 1 A dining room
- 2 A diet kitchen
- 3 An office for the staff
- 4 A record office with a calculating machine
- 5 A room for preparing and storing specimens
- 6 A laboratory for routine clinical pathology
- 7 A separate laboratory for the many special procedures required in metabolic work

All space should be air conditioned because it is becoming increasingly clear that environmental control is extremely desirable for the most consistent results. In addition the comfort of the staff which is usually on a 12 month schedule, is very important for efficient work.

#### Finances of a Metabolic Ward

Few persons not actively engaged in metabolic work realize the large expenses involved. For patient care alone in 1949 the cost in Chicago was of the order of \$14 per patient per day. Add to this the numerous and sometimes expensive laboratory procedures and the cost is such that metabolic work should not be started without an assured adequate budget.

#### Collection, Preparation and Storage of Specimens

These topics have been discussed in previous parts of this manual. If at all feasible it is better to use analytical data for food consumption rather than calculated data. Our experience with stool nitrogen is that for a given patient

within wide limits the daily stool nitrogen is quite constant but that different patients have greatly different characteristic excretions. Therefore it is unsound to use calculated values for stool nitrogen.

## B DERIVED DATA— THEORETICAL BALANCES

### 1 Calcium, Phosphorus and Nitrogen

The concept of theoretical balances has been extensively used by Albright and colleagues (*see* reference previously cited). By a *theoretical balance* is meant the balance of an element such as phosphorus which would be expected from the observed balances of other elements, such as nitrogen and calcium. The two chief virtues of such derived data are their assistance in detecting gross errors in the data and their emphasis on situations where the measured balances differ from the theoretical.

For computation of theoretical calcium, phosphorus and nitrogen balances one from the other the following factors are used:

$$\begin{aligned}\text{Calcium/phosphorus in bone} &= 2.3 \\ \text{Nitrogen/phosphorus in muscle} &= 14.155\end{aligned}$$

In computing the balances the following equations are used assuming that muscle is characteristic of total protoplasmic mass:

- a. Theoretical phosphorus from calcium and nitrogen —

$$\frac{\text{gm calcium retained}}{3} + \frac{\text{gm nitrogen retained}}{15}$$

- b. Theoretical calcium from phosphorus =

$$(\text{gm phosphorus retained} \times 2.23)$$

Theoretical nitrogen from phosphorus and calcium =

$$(\text{gm phosphorus retained} \times 14.7) - \left( \frac{\text{gm calcium retained}}{3} \times 14.7 \right)$$

### 2 Potassium and Sulfur

From potassium, sulfur and nitrogen it is possible to calculate theoretical balances. The factors used are:

$$\begin{aligned}\text{mEq potassium} &= 7 \times \text{gm nitrogen} \\ \text{mg potassium} &= 105 \times \text{gm nitrogen} \\ \text{gm sulfur} &= \frac{\text{gm nitrogen}}{14}\end{aligned}$$

To calculate the theoretical balances the following equations are used (note that in these cases the results are not additive as they are for Ca, P and N):

- a. Theoretical nitrogen from potassium =  $\frac{\text{mEq K retained}}{7}$  or  $\frac{\text{mg K retained}}{105}$

- b. Theoretical nitrogen from sulfur =  $\text{gm S retained} \times 14.5$

- c. Theoretical mEq potassium from nitrogen =  $\text{gm N retained} \times 2.7$

- d. Theoretical mg potassium from nitrogen =  $\text{gm N retained} \times 105$

- e. Theoretical mEq potassium from sulfur =  $\text{gm S retained} \times 14.5 \times 2.7$

$$f \text{ Theoretical mg potassium from sulfur} = \text{gm S retained} \times 14.5 \times 105$$

$$g \text{ Theoretical sulfur from nitrogen} = \frac{\text{gm nitrogen}}{14.5}$$

$$h \text{ Theoretical sulfur from potassium} = \frac{\text{mEq K retained}}{2.7 \times 14.5} \text{ or } \frac{\text{mg K retained}}{105 \times 14.5}$$

## C DERIVED DATA— THEORETICAL WEIGHT

Albright and colleagues have used the concept of "theoretical weight" as a means of discussing the causes of shifts of body weight. The conditions under which their calculations are valid are

- 1 The caloric output must be constant
- 2 The state of hydration and saturation with glycogen must be the same at the beginning and end of the study
- 3 In other words the patient must be in the same state of equilibrium at the beginning and end of the study. The calculations become more accurate the longer the study.

In principle one computes the following changes in weight from metabolic data

- 1 Change in protoplasmic mass associated with nitrogen, from nitrogen balance
- 2 Change in weight of extracellular fluid not accounted for with protoplasm from sodium or chloride balance
- 3 Change in weight of intracellular fluid glycogen or both not accounted for with protoplasm from potassium balance

The difference between the observed change in weight and the sum of the above three separate changes in weight is related to the "caloric discrepancy" or the weight of fat which would be necessary to account for the change in weight on the basis of energy alone.

To carry out the actual calculations Table 26 gives the necessary conversion factors

The steps in the calculation are

- 1 Net change in body weight gm =  
Body wt in gm at beginning - Body wt in gm at end
- 2 Change in weight of protoplasm fat free but not extracellular fluid related to nitrogen gm =  
Nitrogen balance gm  $\times$  32
- 3 Change in weight in grams of extracellular fluid not accounted for with protoplasm =

$$\frac{(\text{Sodium balance mEq})}{(151 \times 1000)} - (\text{Nitrogen balance, gm} \times 5)$$

$$\frac{(\text{Chloride balance mEq})}{(104 \times 1000)} - (\text{Nitrogen balance gm} \times 5)$$

- 4 Change in weight in grams of intracellular fluid and/or glycogen not accounted for with protoplasm =

$$\frac{(\text{Potassium balance mEq})}{(143 \times 1000)} - (\text{Nitrogen balance gm} \times 19)$$

- 5 Theoretical weight change assuming caloric balance grams =  
 Sum of changes in weight due to protoplasm extracellular fluid and intracellular fluid
- 6 Caloric discrepancy i.e. grams fat required to account for difference between observed and theoretical weight change =  
 (Observed net weight change) - (theoretical weight change)

TABLE 6

## FACTORS USED IN DERIVING CERTAIN COMPONENTS OF MUSCLE PROTOPLASM FROM NITROGEN BALANCE

(Taken from Reifenstein Albright and Wells J Clin Endocrinol 5 367 9, Nov 1945)

| COMPONENT TO BE DERIVED   | FINAL UNIT | FACTOR |
|---|------------|--------|
| 1 Protein in protoplasm   | gm         | 6.5    |
| 2 Protoplasm fat free but not extracellular fluid free                  | gm         | 3.0    |
| 3 Protoplasm fat free and extracellular fluid free (true muscle)        | gm         | .97    |
| 4 Intracellular fluid in protoplasm                                     | ml         | 19     |
| 5 Extracellular fluid in protoplasm                                     | ml         | 5      |
| 6 Potassium in extracellular fluid in protoplasm                        | mEq        | 2.7    |
| 7 Sodium in extracellular fluid in protoplasm                           | mEq        | 0.77   |
| 8 Phosphorus in protoplasm  | gm         | 0.066  |
| 9 Sulfur in protoplasm  | gm         | 0.069  |
| 10 Fat calorically equivalent to protoplasm                             | gm         | 9      |
| 11 Protoplasm (fat free) minus fat calorically equivalent to protoplasm | gm         | 9      |

To derive the component multiply nitrogen in gram by the factor

The paper of Reifenstein Albright and Wells gives the justification refinements and applications of this sort of theoretical calculation. Their method of computation is different from that presented above. To start their computation they calculate what the observed weight loss would have been had the subject been in nitrogen balance then express this in terms of fat and correct this fat figure for protoplasm (from Nitrogen balance) extracellular fluid (from Sodium balance) and intracellular fluid (from Potassium balance). Their method of calculation therefore involves one more assumption than the one presented in detail above i.e. the assumption that it is permissible to take the observed change in body weight as a function of fat. Their method of calculation has the advantage that it emphasizes discrepancies between observed and theoretical weight change whereas the method presented above derives fat by difference and therefore the theoretical always equals the observed change in weight.

## SECTION IX

# CLINICAL LABORATORY PROCEDURES

---

### A MISCELLANEOUS REFERENCES TO TEXTBOOKS

It is not the intention of this manual to reduplicate detailed instructions for clinical laboratory procedures which have been well described in many textbooks of clinical methods. The selected list of references below will give the reader a satisfactory working library in which may be found details of most of the methods commonly in use. In the sections which follow, the term "standard methods" refers to methods described in one or the other of the selected references. References to individual papers will be found in the text that follows this bibliography.

- Bray W E. *Synopsis of Clinical Laboratory Methods*, ed. 4. St. Louis 1950. The C V Mosby Company.
- Gradwohl, E B H. *Clinical Laboratory Methods and Diagnosis*, ed. 4. St. Louis 1948, The C V Mosby Company.
- Ham H. and others. *Syllabus for Second Year Course in Laboratory Diagnosis*. Preliminary Edition, 1948. Dept. of Medicine. Harvard Medical School, Longwood Publishing Company.
- Hawk P B, Osler B L, and Summerson W H. *Practical Physiological Chemistry*, ed. 12. Philadelphia-Toronto 1947. The Blakiston Company.
- Kolmer, J A. *Clinical Diagnosis by Laboratory Examinations*. 1st Edition revised. New York-London 1944. D Appleton Century Co. Inc.
- Levinson S. A., and McFate, R P. *Clinical Laboratory Diagnosis*, ed. 3. Philadelphia 1946. Lea and Febiger.
- Osgood E L. *A Textbook of Laboratory Diagnosis*, ed. 3. Philadelphia-Toronto, 1948. The Blakiston Company.
- Simmons J S. and Gentzlow C J. *Laboratory Methods of the United States Army*, ed. 5, Philadelphia 1944. Lea and Febiger.
- Stutt E R, Clough P W, and Clough M. C. *Practical Bacteriology, Hematology and Animal Parasitology*, ed. 9, Philadelphia 1945. The Blakiston Company.
- Todd J C. and Sanford A W. *Clinical Diagnosis With Laboratory Methods*, ed. 10. Philadelphia and London 1943. W B Saunders Company.
- U S War Department. *Methods for Laboratory Technicians*. Oct 1946, U S Government Printing Office.
- Wadsworth A B. *Standard Methods of the Division of Laboratories and Research of the New York State Department of Health*, ed. 3. Baltimore 1947, Williams and Wilkins Company.
- Whitby L E H. and Britton C J C. *Disorders of the Blood*, ed. 5, Philadelphia 1946. The Blakiston Company.
- Wintrobe M. M. *Clinical Hematology*, ed. 2 (thoroughly revised). Philadelphia 1947, Lea and Febiger.

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### B HEMATOLOGY

##### 1 Calculation of Mean Corpuscular Volume Mean Corpuscular Hemoglobin Content and Mean Corpuscular Hemoglobin Concentration

The calculation of mean corpuscular volume hemoglobin content and hemoglobin concentration has been popularized by Wintrobe. These derived values are useful in the classification of anemias. The data required for the computation are

- 1 Erythrocyte count millions/cu mm
- 2 Hematocrit ml packed red cells/100 ml blood
- 3 Hemoglobin concentration gm Hb/100 ml blood

The equations used are

- 1 Mean corpuscular volume cubic microns =

$$\frac{\text{Hematocrit} \times 10}{\text{Erythrocyte count}}$$

Mean corpuscular hemoglobin content micromicrograms =

$$\frac{\text{Hemoglobin concentration of blood} \times 10}{\text{Erythrocyte count}}$$

- 2 Mean cell hemoglobin concentration % Hb =

$$\frac{\text{Hemoglobin concentration} \times 100}{\text{Hematocrit}}$$

The line chart (Fig. 63) facilitates the computations

- 1 To obtain mean corpuscular volume stretch a straight edge between hematocrit (right hand) and red count (left hand). Read MCV on middle line.
- 2 To obtain mean corpuscular hemoglobin content stretch a straight edge between hemoglobin concentration (right hand) and red count (left hand). Read MCH on middle line.
- 3 To obtain mean corpuscular hemoglobin concentration stretch a straight edge between hemoglobin concentration (right hand) and hematocrit (left hand).

# MCV, MCHC AND MCH FROM HEMATOCRIT, HEMOGLOBIN AND RED BLOOD CELL COUNT




Fig. 63—Line chart for calculating cellular derived data using Wintrobe's equations.

## 2 Calculation of Sedimentation Rate

### References

- 1 Rourke M D and Ernestine A C A Method for Correcting the Erythrocyte Sedimentation Rate for Variations in the Cell Volume Percentage of Blood *J Clin Investigation* 8 545 559 (June) 1930
- 2 Wintrobe W M and Landsberg J W A Standardized Technique for the Blood Sedimentation Test *Am J M S* 169 107 115 (Jan) 1935

### Calculation of Sedimentation Rate

Most authors who have published data on the blood sedimentation agree that it is necessary to convert the readings obtained to a standard hematocrit of 45 to 47 ml/100 ml blood. Opinions vary on the exact form of the correction but it is agreed that the sedimentation rate becomes higher the lower the hematocrit. The chart of Rourke and Ernestine will suffice for most of the methods.

- 1 Obtain the uncorrected reading in mm per hour or mm per minute
- 2 Determine the hematocrit
- 3 On the chart (Fig 64) use the left hand line in mm per minute as it stands or multiply by 60 to convert to mm per hour
- 4 On the line for hematocrit as observed find the point for mm per minute as observed
- 5 Follow the empirical line down to hematocrit 45% and read the corrected sedimentation rate on the left hand line
- 6 Compare with the normal for hematocrit 45 (solid block)

## 3 Direct Eosinophil Count

### Reference

- 1 Randolph T G Blood Studies in Allergy I. The Direct Counting Chamber Determinations of Eosinophils by Propylene Glycol Aqueous Stains *J Allergy* 15 89 95 (March) 1944
- 2 Randolph T G Differentiation and Enumeration of Eosinophils in the Counting Chamber With a Glycol Stain a Valuable Technique in Appraising ACTH Dosage *J Lab & Clin Med* 34 1696 1698 (Dec) 1949

### Principle

Counting eosinophils from a differential smear is not accurate enough for many research purposes. A direct count can be made in a counting chamber after dilution of the blood with the appropriate solution.

### Apparatus

- 1 A quantity of white blood cell pipettes
- 2 Counting chamber (Levy)
- 3 A microscope
- 4 Petri dishes
- 5 Mechanical shaker for blood pipettes

### Reagents

- 1 Heparin solution Prepare heparin tubes by placing one drop of heparin solution in a glass vial of 10 ml capacity. Allow drop to evaporate to dryness before using tube.






Fig 64 —Cor ection of hematocrit values by the method of Rourke and Ernsten.

**2 Diluting fluid**

**Solution A** Prepare a 0.05 per cent solution of phloxin in propylene glycol and water (3:1 mixture) Store in a brown bottle

**Solution B** Prepare an 0.02 per cent solution of methylene blue in propylene glycol and water (1:1 mixture) Store in a brown bottle

The diluting fluid is prepared fresh daily by mixing three parts of solution A with one part of solution B

**Procedure**

- 1 Draw 0.5 ml of blood by venipuncture and transfer the blood carefully to the heparinized tube Mix gently
- 2 After thoroughly mixing draw blood to the 10 mark in a white blood count pipette
- 3 Bring the liquid level up to the 11 mark with the freshly prepared diluting fluid
- 4 Place the pipette in a mechanical mixing device and mix the cell suspension for 3 minutes
- 5 After discarding the first 4 to 5 drops fill two Levy counting chambers with fluid from the pipette
- 6 The filled counting chambers are placed under an inverted Petri dish containing a piece of moist filter paper to prevent evaporation Twenty to 30 minutes are allowed for the settling of cells
- 7 After sedimentation enumerate the eosinophils in the 9 square millimeters on both sides of the two counting chambers The eosinophils are identified by their deeply red stained granules Other granules take a faint blue stain

**Calculation**

Total eosinophil count = average eosinophils enumerated per side (9 sq mm)  $\times$  11 = cells per cubic millimeter

**Precautions**

- 1 Timing is important Try to arrange the work so that all counts are made at approximately the same time after dilution of the blood
- 2 The exact proportions of the diluting fluid are important

**4 Hematology Miscellaneous References****Anticoagulants**

Wintrobe M M J Lab & Clin Med 17 899-91 (June) 1933

Moore R M Blood Coagulation and Its Clinical Modification With Anticoagulants Texas Rep on Biol and Med 5 ( ) 1-31 1947

**Ashby Count**

Ebert R V and Pomeroy C P Jr A Clinical Study of Transfusion Reactions The Hemolytic Effect of Group O Blood and Pooled Plasma Containing Incompatible Isoagglutinins J Clin Investigation 25 67-83 (July) 1946

**Bone Marrow**

Berman L and Axelrod A R Aspiration of Sternal Bone Marrow Technique for Obtaining Volumetric Readings Smears Imprints and Histopathologic Sections Am J Clin Pathol 17 61-66 (Jan) 1947

**Clotting Time**

Pohlke F J and Taylor F H L Coagulation Defect in Hemophilia J Clin Investigation 18 741-747 (Sept) 1937

**Erythrocytes**

Joazez J H P Determination of Red Blood Corpuscles with the Aid of a Photoelectric Colorimeter Nederl tydschr v geneesk 92 20-4 0-6 1948

Piper A The Diffraction Method of Measuring Red Blood Cells J Lab & Clin Med 32 (7) 857-877 (July) 1947

## 402 *Clinical Laboratory Procedures*

### **Fibrinolysis**

Tagren, H J Levenson S M Davidson, C S and Taylor F H L Occurrence of Fibrinolysis in Shock, With Observations on Prothrombin Time and Plasma Fibrinogen During Hemorrhagic Shock *Am J M Sc* 211 88 96 (Jan ) 1946

### **Icteric Index**

Hosley R J Nomograph for Rapid Calculation of Icteric Index When Determined by Colorimeter *Am J Clin Path* 19 884 885 (Sept ) 1949

### **Leukocytes**

Minor A H, and Burnett L A Method for Obtaining Living Leukocytes From Human Peripheral Blood by Acceleration of Erythrocyte Sedimentation *Blood* 3 ,99 80° (July) 1948

Vallee B L, Hughes W L Jr and Gibson, J G II The Separation of Leucocytes From Whole Blood by Flotation on Serum Albumin, *Blood Special Issue* 1 8 87 (July) 1947

### **Platelets**

Pohle, F J The Blood Platelet Count in Relation to the Menstrual Cycle in Normal Women *Am J M Sc* 197 49 47 (Jan ) 1939

### **Prothrombin Time**

Munro F L and Munro M P The Preparation of Prothrombin by Adsorption on and Elution From Aluminum Hydroxide *Arch Biochem* 15 (2) 295 304 (Nov ) 1947

Honorato René A Simple Means to Determine Exact Moment of Clotting in Prothrombin or Thrombin Time Determinations, *Proc Soc Exper Biol & Med* 65 (1) 41 (May) 1947

Ware Arnold G Guest M Mason and Seegers Walter H A Factor in Plasma Which Accelerates the Activation of Prothrombin *J Biol Chem* 169 (1) 31 33 (June) 1947

Sternberger L A The Stabilisation of Thrombin in Plasma Development of a Simple Two stage Method for the Determination of Prothrombin *Brit J Exper Path* 28 (3) 165 177 (June) 1947

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### C TOLERANCE TESTS

##### 1 Glucose Tolerance Tests

###### A. ORAL GLUCOSE TOLERANCE

###### Reference

Janney N W and Isaacsen V I A Blood Sugar Tolerance Test Proc Soc  
Exper Biol & Med 15 1-16 1917

###### Principle

After an oral dose of 50 to 75 gm of glucose the blood sugar concentration normally comes back to the control level in two hour and little if any glucose appears in the urine. In diabetes mellitus the level stays high and glucose appears in the urine. In renal glycosuria the blood level behaves normally but glucose appears in the urine.

###### Apparatus

- 1 Micro pipettes 0.1 ml accurately calibrated
- 2 A quantity of watch glasses
- 3 A quantity of 15 ml round bottom centrifuge tubes
- 4 Apparatus as for blood sugar analysis
- 5 A lancet for getting blood

###### Reagents

- 1 Dilute tungstic acid solution for deproteinizing blood (see glucose method)
- 2 Reagents for quantitative estimation of glucose
- 3 Benedict's qualitative glucose reagent Weigh out separately 8.0 gm copper sulfate ( $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ ) 85.0 gm sodium citrate and 50.0 gm sodium carbonate ( $\text{Na}_2\text{CO}_3$  anhydrous). Dissolve the citrate and carbonate in 400 ml of water and filter into a large beaker. Dissolve copper sulfate in 50 ml of water and pour through a funnel into the bottom of the citrate carbonate solution. As much as possible avoid agitation and loss of  $\text{CO}_2$  to the air. Dilute to one liter. The solution keeps well at room temperature.
- 4 Glucose solution The dose is 1.0 gm glucose/kg body weight. Dissolve glucose in 4 parts of water chill and flavor with lemon.

###### Procedure

- 1 The subject comes to the laboratory without any breakfast. Measure the body weight.  
A blood sample is drawn from the ear or finger and 0.1 ml is diluted to 10 ml with acid tungstate solution.
- 2 Collect a urine specimen.
- 3 The subject then drinks the glucose mixture as soon as possible and time is noted (if the mixture is cold it is more palatable).

- 5 Blood samples are then collected at 30 minutes 60 minutes 120 minutes and 180 minutes and analyzed for glucose by one of the quantitative methods
- 6 Urine samples are collected at one, two and three hours A qualitative test is run using Benedict's reagent Heat 10 ml of Benedict's reagent in a water bath. Add 8 drops of urine Heat for 5 minutes Cool Read in terms of color of sediment and supernatant fluid The gradations are blue green, yellow and brown

**Calculation**

As for glucose estimation

**Precautions**

- 1 As for glucose estimation
- 2 For interpretation see standard textbooks According to Conn (Conn J W Am J M Sc 199 555 564 1940) a short period of carbohydrate restriction causes a marked delay in the utilization of carbohydrate as indicated by the glucose tolerance test Undernourished persons lose even more carbohydrate tolerance than do normals Erroneous diagnosis of diabetes mellitus may be made when the factor of previous carbohydrate restriction is neglected Hence the use of preparatory diet (300 gm of carbohydrate 80 gm protein at best maintenance caloric intake) for at least 3 days before the test is recommended.

**B INTRAVENOUS GLUCOSE TOLERANCE****Reference**

Lozner E L Winkler A W Taylor F H L and Peters J P Intravenous Glucose Tolerance Test J Clin Investigation 20 507 515 (Sept) 1941

**Principle**

A single sugar determination on a specimen of blood drawn two hours after intravenous injection of 5 gm of glucose is a reliable guide in distinguishing between benign glycosuria and diabetes mellitus

**Apparatus**

- 1 Sterile needles and syringes for injection of glucose and drawing of blood
- 2 Apparatus as for glucose analysis

**Reagents**

- 1 Glucose solution a 50% sterile solution in water
- 2 Reagents as for blood sugar

**Procedure**

- 1 Prepare the patient as for oral glucose tolerance
- 2 Collect a sample of blood before injection
- 3 Inject 50 ml of glucose solution over a period of 10 minutes
- 4 For screening purposes a sample of blood at 2 hours is used For more precise studies samples are taken at 15 minutes 30 minutes 60 minutes and 120 minutes together with specimens of urine for qualitative test

**Calculation, Precautions**

Same as for oral test

**Interpretation**

See reference of Lozner and associates

## 2 Vitamin Tolerance Tests

### A. ORAL VITAMIN TOLERANCE

#### References

- 1 Johnson P E Contreras L A Consolazio C F and Robinson P F A Comparison of Intravenous and Oral Vitamin Tolerance Tests *Am J Physiol* 144 54-63 (June) 1945
- 2 Goldsmith G A *Chemical Measurements in Relation to Physiological Evidence of Malnutrition Federation Proc* 8 553-567 (June) 1949

#### Principle

A dose of vitamins is administered and the urine is collected thereafter for four hours. If the body is well saturated with vitamins a considerable percentage of the test dose is excreted. If the body is unsaturated the vitamins are absorbed by the tissues and the urinary excretion is small.

#### Apparatus (assuming 50 subjects)

- 1 A two gallon coffee pot or other container from which pouring is easy
- 2 A large spoon for stirring
- 3 50 paper cups
- 4 Apparatus for collecting and storing urine as described in previous sections

#### Reagents

- 1 Water soluble tablets of thiamine hydrochloride   mg each
- 2 Water soluble tablets of riboflavin 5 mg each
- 3 Water soluble tablets of nicotinamide 50 mg each
- 4 Water soluble tablets of ascorbic acid 100 mg each
- 5 Chemicals required for collecting and storing urine

#### Procedure

- 1 Allow at least 2 hours for mixing the solutions before administration of the dose
- 2 For each subject expected add 150 ml of water to the pot and 1 tablet each of thiamine riboflavin and nicotinamide. Stir until dissolved.
- 3 Before administration add for each subject 5 tablets of vitamin C. Stir until dissolved. Add the pills just in time to insure complete solution. 1/2 hour is usually enough except in cold weather.
- 4 On a convenient table place one cup for each subject and pour 150 ml of solution. It is well to calibrate the cups with a heavy pencil line at the 150 ml mark.
- 5 Immediately after the voiding of the fasting specimen command the subjects to file past the table take a cup drink its contents in the presence of an observer, and deposit the empty cups in a garbage can or other disposal unit provided.
- 6 The doses of vitamins are so large that the subject may now be allowed to eat breakfast without interfering with the test.
- 7 All subjects should be briefed with instructions that the urine for the next 4 hours is to be collected. If they have to urinate an observer should provide a paper cup for the purpose.
- 8 Exactly 4 hours after the administration of the vitamins the specimens of urine should be collected measured and stored as described in section on storage of specimens.




#### Calculations

It is desired to determine how much of the dose of vitamins is excreted in 4 hours over and above the amount which would have been excreted if the subject had not taken any vitamins. The assumption is made that the hourly excretion measured

# EXCRETION OF VITAMINS FROM URINARY VOLUME AND CONCENTRATION OF VITAMINS

CONCENTRATION  
PER 100 ml URINE

TOTAL AMOUNT  
OF SUBSTANCE



in the fasting specimen would have changed very little if the subjects had not taken vitamins. The following general formula is used:

$$\text{Excess excretion in 4 hours} = \left( \frac{\text{urine volume in load period}}{100} \times \text{concentration of vitamin in loaded specimen} \right) - (\text{excretion per hour in fasting specimen} \times 4)$$

Calculations are facilitated by the line chart in Fig. 6.

#### Example

##### 1 Fasting urine

|                                    |        |      |          |     |
|------------------------------------|--------|------|----------|-----|
| Vitamin C                          | mg/hr  | 1    | mg/4 hr  | 5   |
| Vitamin B <sub>1</sub>             | mcg/hr | 12   | mcg/4 hr | 48  |
| N <sup>1</sup> Methyl nicotinamide | mg/hr  | 0.50 | mg/4 hr  | 3.4 |
| Vitamin B <sub>6</sub>             | mg/hr  | 80   | mg/4 hr  | 0   |

##### 2 Tolerance test

Volume of urine in 4 hrs 300 ml

|                |      |     |          |      |                  |      |
|----------------|------|-----|----------|------|------------------|------|
| C              | mg%  | 1.5 | mg/4 hr  | 45   | mg over fasting  | 40   |
| B <sub>1</sub> | mcg% | 80  | mcg/4 hr | 240  | mcg over fasting | 192  |
| NMN            | mg%  | 0.0 | mg/4 hr  | 0.75 | mg over fasting  | 3.8  |
| B <sub>6</sub> | mcg% | 500 | mcg/4 hr | 1800 | mcg over fasting | 1480 |

#### Precautions

- 1 The solution of vitamins must be protected from direct sunlight at all times because of the sensitivity of riboflavin. Vitamin C is so unstable that it must be added only just before administration.
- 2 Very strict supervision is required to insure that the subjects drink all the vitamins and also that they collect all of their urine.
- 3 The load test described above is a necessary compromise between practical conditions in the field and scientific desirability. In particular a longer collection period would be desirable owing to the relatively slow excretion of certain of the vitamins. Many workers feel that intravenous injection is preferable to oral administration.
- 4 It is advisable to buy vitamins in large batches and essential to assay each batch before use.
- 5 For every 50 subjects it is a good idea to make up enough solution for 5 extra doses in case of spillage or inaccurate measuring.

## B INTRAVENOUS VITAMIN TOLERANCE

#### Reference

Same as for oral vitamin tolerance.

#### Principle

Same as for oral vitamin tolerance.

#### Apparatus

- 1 Sterile needles and syringes for injecting vitamins.
- 2 Apparatus for collection and storage of urine specimens.

#### Reagents

- 1 Reagents for collecting and storing urine.
- 2 Ampoules of vitamin solutions for intravenous injection. Because of incompatibility of some of the vitamins it is necessary to have two sets of ampoules:
  - a Ascorbic acid 0 mg in 4 ml



- b In 5 ml 5 mg thiamine hydrochloride 3 mg riboflavin, 75 mg nicotinamide  
1 mg calcium pantothenate 1 mg of methylglucamine salt of folic acid 10  
mg of biotin and 10 mg of pyridoxine hydrochloride

**Procedure**

- 1 Obtain a fasting hour specimen as described under oral vitamin tolerance
- 2 In a 20 ml sterile syringe mix 4 ml of the ascorbic acid with 5 ml of B complex solution and inject the 9 ml quantitatively

**Calculation**

Same as for oral vitamin tolerance

**Precautions**

Vitamin tolerance tests are not well standardized. Their interpretation is well discussed in Bulletin No 117 National Research Council Nutrition Surveys Their Techniques and Value May 1949. A mathematical analysis of the time course of vitamin excretion is given by Johnson and colleagues.

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### D KIDNEY FUNCTION

##### 1 Titratable Acidity

###### Reference

Henderson L. J. and Palmer W. W. On Several Factors of Acid Excretion J Biol. Chem. 17: 305-315 (Apr.) 1914.

###### Principle

Titratable acidity is measured by the amount of alkali necessary to alter the pH of a biological fluid to 7.4

###### Apparatus

1. A quantity of 500 ml Erlenmeyer flasks and a burette calibrated in 0.05 ml divisions
2. Calibrated pipettes or syringe pipettes 2 and 10 ml
3. A 250 ml volumetric flask.

###### Reagents

1. Phosphate solution of pH 7.4 M/15 buffer solution prepared according to Sørensen's table on phosphate buffers. Use 80.8 ml of M/15  $\text{Na}_2\text{HPO}_4$  and 19.2 ml of M/15  $\text{KH}_2\text{PO}_4$ .
2. Sodium hydroxide (NaOH) a 0.1000 N solution accurately titrated
3. Phenol red 0.1% aqueous solution of sodium salt. Dissolve 0.5 gm phenol red in 500 ml of water and add 3.5 ml of 0.05 N NaOH
4. Standard comparison solution. In a 250 ml volumetric flask dilute 2 ml phenol red and 10 ml standard phosphate buffer to 100 ml with water

###### Procedure

1. To a 500 ml Erlenmeyer flask add exactly 10 ml of urine. Add approximately 150 ml of water and exactly 2 ml of phenol red solution
2. Dilute to the mark and shake well
3. Titrate the unknown with 0.1000 N NaOH until the color matches the standard. The end point may be checked with the spectrophotometer at a wave length of 511 m $\mu$  or the pH itself can be checked with a glass electrode

###### Calculation

Titratable acidity ml of 1 N titratable acid/liter of urine =  
$$\frac{(\text{ml alkali}) \times \text{normality of alkali} \times 1000}{\text{ml urine}}$$

###### Example

Ten ml of urine required 5.5 ml of 0.1000 N NaOH to change the pH to 7.4  
Titratable acidity ml 1 N titratable acid/liter urine =  
$$\frac{(5.5) \times 0.1000 \times 1000}{10} = 55.5$$

**Precautions**

Matching the end point colors may be difficult in the presence of a highly pigmented specimen of urine. In such cases it may be necessary to use a 5 ml sample.

**2 Phenolsulfonphthalein Test****Reference**

Rowntree L G and Geraghty J T. The Phthalein Test. An Experimental and Clinical Study of Phenolsulphonphthalein in Relation to Renal Function in Health and Disease. *Arch Int Med* 9: 284-338 (March) 1912.

**Principle**

The normal kidney clears injected phenolsulfonphthalein rapidly from the blood stream. The amount of dye excreted in a given time is measured spectrophotometrically.

**Apparatus**

- 1 Graduated cylinders 1000 ml
- 2 Coleman Jr Spectrophotometer
- 3 Cuvettes 19 x 150 mm.
- 4 Sterile 1 ml tuberculin syringes and #23 needles

**Reagents**

- 1 Sodium hydroxide approximately 1 N 40 gm/1000 ml water
- 2 Phenolsulfonphthalein (PSP) 1 ml equivalent to 6 mg of dye sterilized in ampoules
- 3 Standard solution of dye Dilute 1 ml of dye (= 6 mg dye) to exactly 2000 ml with water containing a few ml of NaOH.

**Procedure—Collection of Specimens**

- 1 The subject omits breakfast the morning of the test
- 2 Allow the subject to drink 300-400 ml of water
- 3 One half hour later instruct him to urinate and discard the specimen
- 4 Inject 1 ml of PSP (6 mg) intravenously and note the time
- 5 Collect urine samples quantitatively at 15, 30 and 60 minutes

**Procedure—Analysis of Specimens**

- 1 Pour the whole specimen of urine into a graduated liter cylinder
- 2 Add enough NaOH to make urine alkaline (10-15 ml)
- 3 Dilute to 1000 ml. Mix well
- 4 Read aliquot in cuvette setting water blank at 100% T at 520 mμ

**Calibration Curve**

- 1 Using the standard PSP solution as described in reagents prepare a calibration curve of the following dilutions

| STANDARD<br>ml | WATER<br>ml | DYE/1<br>mg | % OF DYE<br>INJECTED |
|----------------|-------------|-------------|----------------------|
| 100            | 0           | 3.0         | 50                   |
| 80             | 20          | 2.4         | 40                   |
| 60             | 40          | 1.8         | 30                   |
| 40             | 60          | 1.2         | 20                   |
| 20             | 80          | 0.6         | 10                   |
| 10             | 90          | 0.3         | 5                    |

\* Read these against a water blank set at 100% T at wave length 520 mμ.

- 3 Plot the calibration curve of % T against % dye excreted

### Calculation

- 1 Calculate "C" dye in each specimen by reference to the calibration curve
- 2 % dye excreted in 60 minutes =  
(Sum of "C" dye in specimens 1 and 2)

### Example

- 1 The "C" dye excreted in specimen 1 was 40 in specimen 25 and in specimen 10
- 2 % dye excreted in 60 minutes =  
 $45 + 25 + 10 = 80\%$
- 3 The time course of excretion was normal

### Precautions

- 1 New calibration curves should be run with every new batch of dye
- 2 In heavily pigmented urine, a urine blank instead of a water blank may be required.

## 3 Renal Clearance Tests Miscellaneous References

### General

- Smith H W The Physiology of the Kidney New York 1937, Oxford University Press  
 Smith H W Notes on the Interpretation of Clearance Methods J Clin Investigation  
 20 631-635 (Nov) 1941

### Diodrast

- White H L and Kolf D A Rapid Micro Method for Determining Diodrast and Inorganic Iodide Iodine in Blood and Urine Proc Soc Exper Biol & Med 43 1 (Jan) 1940

### Glomerular Filtration Rate

- Earle D P Jr and Berliner B W A Simplified Clinical Procedure for Measurement of Glomerular Filtration Rate and Renal Plasma Flow Proc Soc Exper Biol & Med 62 (2) 62-64 (June) 1946

### Inulin Clearance

- Shaunon J A and Smith H W The Excretion of Inulin Xylose and Urea by Normal Man J Clin Investigation 14 393-401 (July) 1935

### Phenolsulphonphthalein Test

- Scardino P L and Scott W W The Determination of Phenolsulphonphthalein With the Photoelectric Colorimeter and Its Application of the Phthalein Elimination Curve J Urology 58 (2) 143-150 (Aug) 1947

### Renal Blood Flow

- Poà P P and Poà V L A Simple Method for Determining Effective Renal Blood Flow and Tubular Excretory Mass in Man Proc Soc Exper Biol & Med 51 (3) 3-8 (Dec) 1944

### Urea Clearance

- Peters J P and Van Slyke D D Quantitative Clinical Chemistry Vol II Methods Baltimore 1935 Williams and Wilkins Company page 564-572

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### E LIVER FUNCTION

##### 1 Bromsulphalein Test

###### References

- 1 Rosenthal S M and White E C Clinical Application of the Bromsulphalein Test for Hepatic Function J A M A 84 1112 1114 (April) 1925
- 2 Mateer G Baltz J I Marion D F, and MacMillan, J M Liver Function Tests J A. M. A. 121 723 728 (March) 1943
- 3 Gaebler O H Determination of Bromsulphalein in Normal Turbid Hemolyzed or Icteric Serums Am J Clin Path 15 432 455 (Oct) 1945

###### Principle

The dye bromsulphalein is removed from the blood stream by the liver. After intravenous injection bromsulphalein in the serum is measured colorimetrically.

###### Reagents

- 1 Bromsulphalein dye 150 mg/3 ml ampoule
- 2 Hydrochloric acid a 5% solution
- 3 Sodium hydroxide, 0.1 N solution, 4 gm/1000 ml of water

###### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 x 150 mm
- 3 A quantity of 15 ml centrifuge tubes (round bottom)
- 4 Electric centrifuge
- 5 Syringes 5 ml and #20 gauge needles (both sterile)
- 6 Syringe pipette 10 ml

###### Procedure—Collection of Specimens

- 1 The patient is weighed and the dosage to be injected is calculated on a basis of 5 mg per kilogram of body weight
- 2 The calculated amount of dye to be injected is measured into a 5 ml syringe
- 3 A control specimen is drawn from a vein. The syringe is separated from the needle and the syringe containing the dye is connected in its place. The dye is slowly injected into the vein the time being approximately one minute. Care should be taken not to allow infiltration of the dye outside the vein, as it is very painful.
- 4 Record the time of injection and after 45 minutes draw a sample of blood (4 to 5 ml) preferably from the other arm. Allow the blood to run into a clean dry test tube
- 5 Allow the blood to coagulate and then centrifuge and separate the serum

###### Procedure—Analysis of Serum

- 1 Into each of two cuvettes pipette exactly 0.5 ml serum
- 2 To the control tube add 2.5 ml water and 3 ml of hydrochloric acid. Mix.

- 3 To the other tube add 2.0 ml of water and 3 ml of 0.1 N sodium hydroxide
- 4 Set the control at 100% T at 565 m $\mu$  and read the unknown

#### Calibration Curve

- 1 Make up dilution of the dye in dilute sodium hydroxide in the range 0 to 0.1 mg/5 ml
- 2 Using a water blank set at 100% T read the standards at a wave length of 565 m $\mu$ .

#### Calculation

- 1 A value of 10 mg bromsulphalein/100 ml of serum is arbitrarily taken as 100% retention when the dosage is 5 mg of bromsulphalein/kilogram body weight

$$\% \text{ retention} = \frac{(\text{mg in cuvette}) \times 100 \times 100}{(\text{ml serum in cuvette}) \times 10}$$

#### Example

- The % T was 91 % corresponding to 0.0 mg bromsulphalein in the cuvette
- % retention =

$$\frac{(0.0) \times 100 \times 100}{0.5 \times 10} = 40\%$$

#### Precautions

- 1 Painful sloughing will be caused by extravascular bromsulphalein

## 2 Serum Bilirubin

#### Reference

- Malloy H T and Evelyn K. A. The Determination of Bilirubin With the Photoelectric Colorimeter *J Biol. Chem* 119 431-490 (July) 1937

#### Principle

Bilirubin like aromatic amines combines with diazonium salts to form azo dyes. Addition of alcohol sets all the bilirubin free to react with the diazo reagent and is used to make the reaction quantitative for total bilirubin content of the serum.

#### Apparatus

- 1 Coleman Jr Spectrophotometer Model 6
- 2 Cuvettes 19 x 1.0 mm
- 3 A quantity of volumetric flasks 10 ml or 25 ml
- 4 A 10 ml syringe or pipette accurately calibrated

#### Reagents

- 1 Methyl alcohol absolute
- 2 Sulfanilic acid solution Dissolve 1 gm of sulfanilic acid in 15 ml of concentrated HCl and dilute to 1 liter with distilled water
- 3 Sodium nitrite solution Dissolve 0.5 gm of NaNO<sub>2</sub> (ACS) in a little water and dilute to 100 ml. The solution is to be prepared on the day it is used
- 4 Diazo reagent Add 0.3 ml of the nitrite solution to 10 ml of the sulfanilic acid solution. The diazo reagent must be prepared within a few minutes of the time when it is to be used
- 5 Dilute HCl for blanks Dilute 15 ml of concentrated HCl to 1 liter with water
- 6 Bilirubin stock standard. Transfer exactly 40 mg of purified bilirubin to a 100 ml glass stoppered volumetric flask. Dissolve bring to the mark with chloroform (USP) and mix by inverting 10 times (This solution is stable if stored in a brown bottle with properly fitted glass stopper)

- \* Bilirubin working standard Dilute 1 ml of the above stock solution to 100 ml with absolute methyl alcohol The working standard is 0.4 mg/100 ml

#### Procedure—Collection of Specimens

- 1 To avoid turbidity of serum by alimentary lipemia the blood should be taken when the patient is in the postabsorptive state Transfer the sample to a round bottom tube without anticoagulant close with a pledget of cotton and then keep at room temperature for approximately 15 minutes until the clot is completely formed Chill in ice water and place in the refrigerator until the clot has contracted after which it is gently freed from the top part of the tube if there is adherence and the tube then centrifuged After centrifuging transfer the clear serum to a clean dry tube The diazo reaction must be applied within 2 hours after the blood is clotted and if the separated serum is allowed to stand the reaction is altered After long standing only the indirect reactions occur

#### Procedure—Analysis of Serum

- 1 Pipette exactly 1 ml of serum into a 10 ml volumetric flask fill to the mark with distilled water stopper and mix by inverting 10 times
- 2 Into a cuvette pipette exactly 5 ml of absolute methyl alcohol and then 1 ml of fresh diazo reagent
- 3 Follow with exactly 4 ml of the (8 + 1) diluted serum mix stopper and hold 30 minutes at room temperature for the full color to develop (If bilirubin exceeds 5 mg/100 ml dilute the serum to 25 ml instead of 10 ml)
- 4 At the same time prepare a reagent blank with 4 ml of dilute hydrochloric acid 5 ml of absolute methyl alcohol and 1 ml of diazo reagent
- 5 Set the blank at 100% T at 550 mμ and read the unknown

#### Calibration Curve

- 1 Prepare a calibration curve with working bilirubin standard in the range 0 to 0.016 mg bilirubin per cuvette Mix in each case 5 ml of absolute methyl alcohol 1 ml of diazo reagent the requisite amount of working standard and dilute hydrochloric acid to make the total final volume 10 ml
- 2 Stopper all cuvettes tightly mix each by inverting 4 times Let stand together for 30 minutes at room temperature for the color to develop fully then measure the Transmittance (T) at 550 mμ

#### Calculation

$$\begin{aligned} \text{mg bilirubin/100 mg serum} = \\ \frac{(\text{mg in cuvette}) \times (\text{dilution}) \times 100}{(\text{ml aliquot}) \times (\text{ml serum})} \end{aligned}$$

#### Example

One ml serum was diluted to 10 ml of which 4 ml were taken for analysis A reading of 80% T was equivalent to 0.010 mg bilirubin in the cuvette  
mg total bilirubin/100 ml serum =

$$\frac{(0.010) \times 10 \times 100}{4 \times 1} = 2.5 \text{ mg bilirubin}$$

#### Precautions

- 1 Careful handling of the serum is essential to this method
- 2 If it is desired to determine direct and indirect bilirubin two stages are required (a) without methyl alcohol (b) with methyl alcohol See the original paper for this procedure

### 3 Urobilinogen in Urine and Feces

#### A. UROBILINOGEN IN URINE

##### Reference

Schwartz S Shorav V and Watson O J Urobilinogen IV The Determination of Urobilinogen by Means of the Evelyn Photoelectric Colorimeter *Am J Clin Path* 14 593 604 (Dec) 1944

##### Principle

A red color is formed by the reaction between urobilinogen and Ehrlich's aldehyde reagent

##### Apparatus

- 1 Separatory funnels 500 ml
- Graduated cylinders 100 ml
- 3 Cuvettes size 10 x 150 mm
- 4 Coleman Jr Sps photometer Model 6
- 5 A quantity of 175 ml Erlenmeyer flasks
- 6 An 0.5 ml syringe pipette accurately calibrated

##### Reagents

- 1 Sodium carbonate (reagent grade)
- 2 Petroleum ether (b.p. 30-60°C)
- 3 Ferrous sulfate a 0.6% solution (freshly prepared) 0 gm of ferrous sulfate are dissolved in 9 ml of water
- 4 Sodium hydroxide a 10% solution in water
- 5 Ehrlich's aldehyde reagent 0.7 gm para dimethyl amino benzaldehyde 150 ml of concentrated hydrochloric acid and 100 ml of water
- 6 Sodium acetate aqueous saturated solution
- 7 Glacial acetic acid
- 8 Ethyl alcohol 95%
- 9 Standard solution Dissolve 5 mg of Pontacyl Carmine "B" (Du Pont) and 95 mg of Pontacyl Violet 6B in 1000 ml of 0.5% acetic acid. A further dilution of one part in 6 gives a color intensity equivalent to 95 mg of urobilinogen per 100 ml. With a drop of CHCl<sub>3</sub> the color is stable indefinitely. (If available it is better to use pure crystalline urobilinogen)

##### Preservation of Sample

The procedure is applicable to an undisturbed 4 hour sample. If the latter the sample should be in a dark bottle containing approx. 5 gm of sodium carbonate under a layer of petroleum ether.

##### Procedure

- 1 To 50 ml of urine in a 125 ml Erlenmeyer flask add 5 ml of ferrous sulfate solution
- 2 Add 5 ml of 10% sodium hydroxide mix thoroughly and set away in the dark for one hour. Filter or centrifuge
- 3 At the end of the hour the filtrate should not exhibit any urobilin band at 68 mμ. To determine the amount of filtrate to use do a qualitative test first.
  - a To 1 ml of the filtrate add 1 ml of Ehrlich's reagent and 4 to 6 ml of sodium acetate solution
  - b If the resultant color is an intense red 0.5 ml of filtrate is enough for the quantitative procedure. If pale red use 5 to 10 ml. If faint pink use 15 to 30 ml and if there is no color use 0 ml.
- 4 The amount selected is placed in a pear shaped funnel and if small is diluted to bring the volume to 3 ml
- 5 Then add 50 ml of petroleum ether and acidify with 5 ml of glacial acetic acid



## 416 Clinical Laboratory Procedures

- 6 Shake the mixture vigorously and allow to separate. The emulsions are broken by the addition of more acetic acid and in some cases 12 ml of 95% ethyl alcohol is used.
- 7 Decant the petroleum ether from the top by running out as much of the aqueous fraction as possible through the lower stopcock. The aqueous layer is extracted twice more with 25 ml portions of petroleum ether.
- 8 The combined ether extractions are placed in a clean separatory funnel and washed once with a small amount of distilled water which is discarded. This serves to remove small amounts of mechanically entrained impurities.
- 9 Add 2 ml of Ehrlich's reagent and shake the solutions vigorously.
- 10 Now add 6 ml of the sodium acetate solution and shake the mixture again. (At times it is noted at this stage that the color of the aqueous fraction becomes somewhat less intense and at the same time a fine reddish brown precipitate appears at the interphase. This is believed to be indol or skatol aldehyde since it has been observed repeatedly that the same phenomenon occurs when indol is similarly handled with petroleum ether, Ehrlich's reagent and sodium acetate.)
- 11 The colored aqueous fraction after complete separation is removed into a 100 ml graduated cylinder.
- 12 Repeat the extraction of the petroleum ether with the Ehrlich's reagent and sodium acetate one or more times or until the final aqueous fraction is not more than very faint. Dilute to 100 ml.
- 13 Place 10 ml of the well mixed solution into a cuvette.
- 14 For a blank reading add 9 ml of sodium acetate to 3 ml of Ehrlich's reagent. Set at 100°C T at 508 mμ and read the unknown.

### Calibration Curve

The arbitrary standard is equivalent to 0.6 mg urobilinogen/100 ml. Set up a calibration curve in the range from 0 to 10 mcg in 10 ml and run at 508 mμ.

### Calculation

mg urobilinogen/liter of urine =

$$\frac{(\text{mcg in cuvette}) \times (\text{final dilution}) \times (\text{initial dilution}) \times (1000)}{(1000) \times (\text{ml final aliquot}) \times (\text{ml filtrate}) \times (\text{ml urine})}$$

### Example

A reading of 75% T was equivalent to 1.6 mcg in 10 ml in the cuvette. Ten ml aliquots were read of a final dilution of 100 ml. Ten ml of filtrate had been analyzed from an initial dilution of 50 ml of urine to 100.

mg urobilinogen/liter of urine =

$$\frac{(1.6) \times 100 \times 100 \times 1000}{1000 \times 10 \times 10 \times 50} = 32 \text{ mg/liter}$$

### Precautions

- 1 There are many places in this method at which losses will occur. Care is required at every step. The method itself is not as quantitative as one would desire.

## B UROBILINOGEN IN FECES

### Reference Principle

Same as for urobilinogen in urine.

### Apparatus

Same as for urobilinogen in urine and in addition

- 1 Quantitative filter paper
- 2 A number of funnels
- 3 A Waring Blender

**Reagents**

Same as for urobilinogen in urine

**Procedure—Processing of Sample**

- 1 All fecal matter for a period of 3 days is collected as described in the section of this manual on collection of feces
- 2 In Waring Blender the whole mass is homogenized and finally diluted to exactly 1000 ml

**Procedure—Preparation of the Alkaline Extract of Stool**

- 1 Dilute 25 to 100 ml of the emulsion up to exactly 300 ml with water
- 2 Add 100 ml of ferrous sulfate
- 3 Then add 100 ml of sodium hydroxide stopper the flask shake and place in the dark for one to two hours
- 4 Filter through quantitative filter paper

**Procedure—Analysis of Alkaline Extract**

- 1 Pipette 20 ml of the extract into a separatory funnel
- 2 Add 1 ml glacial acetic acid
- 3 Then add 30-40 ml of petroleum ether and shake for three minutes
- 4 Release the watery part into a second separatory funnel add 10-40 ml petroleum ether again and shake for three minutes Combine both ether extracts
- 5 Wash twice with distilled water and then add 1 ml of Ehrlich's aldehyde reagent
- 6 Shake for one minute and then add 4 ml of a cold solution of sodium acetate (saturated)
- 7 Shake for 3 minutes and release the colored phase into a 100 ml graduated cylinder
- 8 Repeat the process of adding the Ehrlich's reagent and sodium acetate until no colored compound is formed.
- 9 Combine all the colored washings dilute to 100 ml and read an aliquot in the spectrophotometer setting the blank at 100% T as for urine

**Calibration Curve**

Same as for urobilinogen in urine

**Calculation**

mg urobilinogen in feces/day =

$$\frac{(\text{mcg in cuvette}) \times (\text{dilution \# 1}) \times (\text{dilution \# 2}) \times (2000)}{(100) \times (\text{ml aliquot}) \times (\text{ml alkaline filtrate}) \times (\text{ml stool}) \times (\text{days collected})}$$

**Example**

A 3 day sample was diluted to 1000 ml. Of the emulsion 50 ml were diluted to 500 of which 20 ml aliquots were analyzed the final dilution being 100 ml. The concentration of urobilinogen in the cuvette was 15 mcg in 10 ml  
mg fecal urobilinogen/day =

$$\frac{15 \times 100 \times 500 \times 1000}{100 \times 10 \times 20 \times 50 \times 3} = 50 \text{ mg/day}$$

**Precautions**

Same as for urobilinogen in urine

**4 Liver Function, Miscellaneous References****General**

- Mateer, J G, Baltz, J L, Steele, H H, Brouwer, S W and Colvert, J P. Chronic Subclinical Impairment of the Liver. Early Diagnosis and Treatment. Further Improvements and Evaluation of Certain Liver Function Tests. J A M A 133: 909-916 (March) 1947
- Osgood, Edwin E. Interpretation of Liver Function Tests. J A M A 134 (7): 855-859 (June) 1944

## 418 *Clinical Laboratory Procedures*

### **Amyloidosis**

Unger P N Zuckerbrod M Beck G J and Steele J M Study of the Disappearance of Congo Red From the Blood of Non amyloid Subjects and Patients With Amyloidosis J Clin Investigation 27 111 117 (Jan) 1948

### **Bile in Urine**

Hawkinson V, Wat on C J and Turner R H Modification of Harrison's Test for Bilirubin in Urine Especially Suited for Mass and Serial Use J A M A 129 514 515 (Oct) 1945

### **Cephalin Flocculation**

Hanger F M Serological Differentiation of Obstructive From Hepatogenous Jaundice by Flocculation of Cephalin Cholesterol Emulsion J Clin Investigation 18 61 69 (May) 1939

### **Cholinesterase Activity**

Vorhaus L J II and Kark R M Measurement of Serum Cholinesterase Activity A Useful Laboratory Tool in the Study of Diseases of the Liver and Biliary System To be published Med cal Nutrition Laboratory Report No 56 June 8 1949

### **Galactose**

Colcher H Patek A J Jr and Kendall F E Galactose Disappearance From the Blood Stream Calculation of the Galactose Removal Constant and Its Application as a Test for Liver Function J Clin Investigation 25 768 775 (Sept) 1946

### **Hemosiderin in Urine**

Fous P Urinary Siderosis Hemosiderin Granules in the Urine as an Aid to Diagnosis of Pernicious Anemia Hemochromatosis and Other Diseases Causing Siderosis J Liver Med 28 645 658 (Nov) 1918

### **Hippuric Acid**

Quick A J The Clinical Application of the Hippuric Acid and Prothrombin Test Am J Clin Path 10 2-- 2-3 (March) 1940

### **Phosphatase**

Cutman A B and Hanger F M Jr Differential Diagnosis of Jaundice I Combined Serum Phosphatase Determination and Cephalin Flocculation Test M Clin North America 25 83 848 (May) 1941

### **Portal Pressure**

Bean W B Paul W D and Franklin M Preliminary Studies on an Indirect Method for Determining Portal Pressures J Clin Investigation 28 769 770 (July) 1949

### **Prothrombin Time**

Aggeler P M Howard J Lucia S P Clark W and Astaff A Standardization of the Quick Prothrombin Test Blood 1 290 L 3 (March) 1946

Quick A J The Clinical Significance of Prothrombin as a Factor in Hemorrhage Pennsylvania M J 43 1-5 130 (Oct) 1939

Unger P N and Shapiro S Prothrombin Response to the Parenteral Administration of Large Doses of Vitamin K in Subjects With Normal Liver Function and in Cases of Liver Disease A Standardized Test for the Estimation of Hepatic Function J Clin Investigation 27 39 41 (Jan) 1948

### **Thymol Turbidity**

Dekema H The Thymol Turbidity Test and the Cephalin Cholesterol Flocculation Test Nederl tijdschr v geneesk 93 (II) 19 4 82 (1949)

Huerga J dela and Popper H Standardized Reagent for Thymol Turbidity Test J Lab & Clin Med 34 (6) 847 850 (June) 1949

Shank R E and Hoagland C L A Modified Method for the Quantitative Determination of the Thymol Turbidity Reaction of Serum J Biol Chem 162 133 138 (Jan) 1946

### **Urobilinogen**

Wilson T M and Davidson L S P Ehrlich's Aldehyde Test for Urobilinogen British M J 1949 834 889 (May) 1949

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### FUNCTIONAL TESTS OF ENDOCRINE, GASTRIC AND OTHER SYSTEMS

##### 1 Pituitary Adrenocorticotrophic Hormone (ACTH) Test

###### Reference

Forsham P H, Thorn G W, Prunty F T G and Hills A G. Clinical Studies With Pituitary Adrenocorticotrophic Hormone. *J Clin Endocrinol* 8: 15-66 (Jan) 1948

###### Principle

Following the injection of pituitary adrenocorticotrophic hormone (ACTH) in normal individuals there is an immediate and striking fall in circulating eosinophils and a rise in the uric acid excretion. Patients with Addison's disease fail to show these changes.

###### Apparatus

- 1 Sterile syringe and needle for intramuscular injection of ACTH
- 2 Apparatus for direct eosinophil count
- 3 Apparatus for estimating uric acid and creatinine

###### Reagents

- 1 ACTH for Intramuscular Injection. The solubility varies with different batches of ACTH. The hormone requires alkalization to a pH between 8 and 9 (faint pink with phenolphthalein). To 5 ml of sterile saline add three drops of approximately N/10 NaOH. This is taken up in a syringe and added to the rubber capped vial containing 5 mg of ACTH. Shake gently until homogeneous. A somewhat cloudy solution is obtained. This should not be kept longer than twelve hours at 4°C or longer than two hours at room temperature. It should be frozen for long storage. The powder appears stable.
- 2 Reagents for direct eosinophil count
- Reagents for estimating uric acid and creatinine

###### Procedure

- 1 No food is allowed after 8:00 P.M. at night
- 2 On the following day give the patient 200 ml of water at 6:00, 8:00 and 10:00 A.M.
- 3 Collect urine from 6:00 to 8:00 A.M.
- 4 An eosinophil count is done at 8:00 A.M.
- 5 Immediately after the 8:00 A.M. specimen inject 5 mg ACTH intramuscularly deep into the gluteus maximus. A moderate degree of soreness often persists for a day or two.
- 6 Urine is collected from 9:00 A.M. to 1:00 noon and the eosinophil count is repeated at noon.
- 7 Analyze the urine specimen for uric acid and creatinine.

## 418 *Clinical Laboratory Procedures*

### **Amyloidosis**

Unger P N Zuckerbrod M Beck G J and Steele J M Study of the Disappearance of Congo Red From the Blood of Non amyloid Subjects and Patients With Amyloidosis J Clin Investigation 27 111 117 (Jan) 1948

### **Bile in Urine**

Hawkinson V Watson C J and Turner R H Modification of Harrison's Test for Bilirubin in Urine Especially Suited for Mass and Serial Use J A M A 129 514 515 (Oct) 1945

### **Cephalin Flocculation**

Hanger F M Serological Differentiation of Obstructive From Hepatogenous Jaundice by Flocculation of Cephalin Cholesterol Emulsion J Clin Investigation 18 661 69 (May) 1939

### **Cholinesterase Activity**

Vorhaus L J II and Kark R M Measurement of Serum Cholinesterase Activity A Useful Laboratory Tool in the Study of Diseases of the Liver and Biliary System To be published Medical Nutrition Laboratory Report No 56 June 8 1949

### **Galactose**

Colecher H Patek A J Jr and Kendall F E Galactose Disappearance From the Blood Stream Calculation of the Galactose Removal Constant and Its Application as a Test for Liver Function J Clin Investigation 25 63 75 (Sept) 1946

### **Hemosiderin in Urine**

Rous P Urinary Siderosis Hemosiderin Granules in the Urine as an Aid to Diagnosis of Pernicious Anemia Hemochromatosis and Other Diseases Causing Siderosis J Exper Med 28 645 658 (Nov) 1918

### **Hippuric Acid**

Quick A J The Clinical Application of the Hippuric Acid and Prothrombin Test Am J Clin Path 10 229 233 (March) 1940

### **Phosphatase**

Cutman A B and Hanger F M Jr Differential Diagnosis of Jaundice by Combined Serum Phosphatase Determination and Cephalin Flocculation Test M Clin North America 25 83 848 (May) 1941

### **Portal Pressure**

Bean W B Paul W D and Franklin M Preliminary Studies on an Indirect Method for Determining Portal Pressures J Clin Investigation 28 69 110 (July) 1949

### **Prothrombin Time**

Aggeler P M Howard J Lucia S P Clark W and Astaff A Standardization of the Quick Prothrombin Test Blood 1 20 223 (March) 1946

Quick A J The Clinical Significance of Prothrombin as a Factor in Hemorrhage Pennsylvania M J 43 125 130 (Oct) 1939

Unger P N and Shapiro S Prothrombin Response to the Parenteral Administration of Large Doses of Vitamin K in Subjects With Normal Liver Function and in Cases of Liver Disease A Standardized Test for the Estimation of Hepatic Function J Clin Investigation 27 23 41 (Jan) 1948

### **Thymol Turbidity**

Dekema H The Thymol Turbidity Test and the Cephalin Cholesterol Flocculation Test Nederl tijdschr Geneesk 93 (II) 19 480 (1949)

Huerga J dela and Popper H Standardized Reagent for Thymol Turbidity Test J Lab & Clin Med 34 (6) 87 880 (June) 1949

Shank R E and Hoagland C L A Modified Method for the Quantitative Determination of the Thymol Turbidity Reaction of Serum J Biol Chem 162 133 135 (Jan) 1946

### **Urobilinogen**

Wilson T M and Davidson L S P Ehrlich's Aldehyde Test for Urobilinogen British M J 1949 884 889 (May) 1949

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### F FUNCTIONAL TESTS OF ENDOCRINE GASTRIC AND OTHER SYSTEMS

##### 1 Pituitary Adrenocorticotrophic Hormone (ACTH) Test

###### Reference

Forsham P H, Thorn G W, Prunty F T G and Hills A G. Clinical Studies With Pituitary Adrenocorticotrophic Hormone. *J Clin Endocrinol* 8: 1, 66 (Jan) 1948

###### Principle

Following the injection of pituitary adrenocorticotrophic hormone (ACTH) in normal individuals there is an immediate and striking fall in circulating eosinophils and a rise in the uric acid excretion. Patients with Addison's disease fail to show these changes.

###### Apparatus

- 1 Sterile syringe and needle for intramuscular injection of ACTH
- 2 Apparatus for direct eosinophil count
- 3 Apparatus for estimating uric acid and creatinine

###### Reagents

- 1 ACTH for Intramuscular Injection. The solubility varies with different batches of ACTH. The hormone requires alkalization to a pH between 8 and 9 (faint pink with phenolphthalein). To 10 ml of sterile saline add three drops of approximately N/10 NaOH. This is taken up in a syringe and added to the rubber capped vial containing 5 mg of ACTH. Shake gently until homogeneous. A somewhat cloudy solution is obtained. This should not be kept longer than twelve hours at 4°C or longer than two hours at room temperature. It should be frozen for long storage. The powder appears stable.
- 2 Reagents for direct eosinophil count
- 3 Reagents for estimating uric acid and creatinine

###### Procedure

- 1 No food is allowed after 8:00 P.M. at night.
- 2 On the following day give the patient 100 ml of water at 6:00, 8:00 and 10:00 A.M.
- 3 Collect urine from 6:00 to 8:00 A.M.
- 4 An eosinophil count is done at 8:00 A.M.
- 5 Immediately after the 8:00 A.M. specimen inject 25 mg ACTH intramuscularly deep into the gluteus maximus. A moderate degree of soreness often persists for a day or two.
- 6 Urine is collected from 9:00 A.M. to 1:00 noon, and the eosinophil count is repeated at noon.
- 7 Analyze the urine specimens for uric acid and creatinine.

# **Amyloidosis**

Unger P N Zuckerbrod M Beck G J and Steele J M Study of the Disappearance of Congo Red From the Blood of Non amyloid Subjects and Patients With Amyloidosis J Clin Investigation 27 111 117 (Jan ) 1948

# **Bile in Urine**

Hawkinson V Watson C J and Turner R H Modification of Harrison's Test for Bilirubin in Urine Especially Suited for Mass and Serial Use J A M A 109 514 515 (Oct ) 1945

# **Cephalin Flocculation**

Hanger F M Serological Differentiation of Obstructive From Hepatogenous Jaundice by Flocculation of Cephalin Cholesterol Emulsion J Clin Investigation 18 661 669 (Mar) 1939

# **Cholinesterase Activity**

Vorhaus L J II and Kark R M Measurement of Serum Cholinesterase Activity A Useful Laboratory Tool in the Study of Diseases of the Liver and Biliary System To be published Medical Nutrition Laboratory Report No 56 June 8 1949

# **Galactose**

Colcher H Patek A J Jr and Wendall F E Galactose Disappearance From the Blood Stream Calculation of the Galactose Removal Constant and Its Application as a Test for Liver Function J Clin Investigation 25 768 775 (Sept ) 1946

# **Hemosiderin in Urine**

Rous P Urinary Siderosis Hemosiderin Granules in the Urine as an Aid to Diagnosis of Pernicious Anemia Hemochromatosis and Other Diseases Causing Siderosis J Exper Med 28 645 658 (Nov ) 1918

# **Hippuric Acid**

Quick A J The Clinical Application of the Hippuric Acid and Prothrombin Test Am J Clin Path 10 222 233 (March) 1940

# **Phosphatase**

Cutman A B and Hanger F M Jr Differential Diagnosis of Jaundice by Combined Serum Phosphatase Determination and Cephalin Flocculation Test M Clin North America 25 837 848 (Mar) 1941

# **Portal Pressure**

Bean W B Paul W D and Franklin M Preliminary Studies on an Indirect Method for Determining Portal Pressures J Clin Investigation 28 1069 1110 (July) 1949

# **Prothrombin Time**

Aggeler P M Howard J Lucia S P Clark W and Astaff A Standardization of the Quick Prothrombin Test Blood 1 220 223 (March) 1946

Quick A J The Clinical Significance of Prothrombin as a Factor in Hemorrhage Pennsylvania M J 43 105 130 (Oct ) 1939

Unger P N and Shapiro S Prothrombin Response to the Parenteral Administration of Large Doses of Vitamin K in Subjects With Normal Liver Function and in Cases of Liver Disease A Standardized Test for the Estimation of Hepatic Function J Clin Investigation 27 39 47 (Jan ) 1948

# **Thymol Turbidity**

Dekema H The Thymol Turbidity Test and the Cephalin Cholesterol Flocculation Test Nederl tijdschr v geneesk 93 (II) 19 48 (1949)

Huerga J dela and Popper H Standardized Reagent for Thymol Turbidity Test J Lab & Clin Med 34 (6) 877 880 (June) 1949

Shank R E and Hoagland C L A Modified Method for the Quantitative Determination of the Thymol Turbidity Reaction of Serum J Biol Chem 162 133 138 (Jan ) 1946

# **Urobilinogen**

Wilson T M and Davidson L S P Ehrlich's Aldehyde Test for Urobilinogen British M J 1949 884 889 (Mar) 1949

## SECTION IX

### CLINICAL LABORATORY PROCEDURES (Cont'd)

---

#### F FUNCTIONAL TESTS OF ENDOCRINE GASTRIC AND OTHER SYSTEMS

##### 1 Pituitary Adrenocorticotrophic Hormone (ACTH) Test

###### Reference

Forsham P H, Thorn G W, Prusky F T G and Hills A G. Clinical Studies With Pituitary Adrenocorticotrophic Hormone. J Clin Endocrinol 8: 15-66 (Jan.) 1948

###### Principle

Following the injection of pituitary adrenocorticotrophic hormone (ACTH) in normal individuals there is an immediate and striking fall in circulating eosinophils and a rise in the uric acid excretion. Patients with Addison's disease fail to show these changes.

###### Apparatus

- 1 Sterile syringe and needle for intramuscular injection of ACTH
- 2 Apparatus for direct eosinophil count
- 3 Apparatus for estimating uric acid and creatinine

###### Reagents

- 1 ACTH for Intramuscular Injection. The solubility varies with different batches of ACTH. The hormone requires alkalization to a pH between 8 and 9 (faint pink with phenolphthalein). To 5 ml of sterile saline add three drops of approximately N/10 NaOH. This is taken up in a syringe and added to the rubber capped vial containing 0.5 mg of ACTH. Shake gently until homogeneous. A somewhat cloudy solution is obtained. This should not be kept longer than twelve hours at 4°C or longer than two hours at room temperature. It should be frozen for long storage. The powder appears stable.
- Reagents for direct eosinophil count
- 3 Reagent for estimating uric acid and creatinine

###### Procedure

- 1 No food is allowed after 8:00 P.M. at night.
- On the following day give the patient 100 ml of water at 8:00, 8:00 and 10:00 A.M.
- 1 Collect urine from 6:00 to 8:00 A.M.
- 4 An eosinophil count is done at 8:00 A.M.
- 5 Immediately after the 8:00 A.M. specimen inject 0.5 mg ACTH intramuscularly deep into the gluteus maximus. A moderate degree of soreness often persists for a day or two.
- 6 Urine is collected from 9:00 A.M. to 1:00 noon and the eosinophil count is repeated at noon.
- 7 Analyze the urine specimens for uric acid and creatinine.



## 418 *Clinical Laboratory Procedures*

### **Amylodiosis**

Unger P N Zuckerbrod M Beck G J and Steele J M Study of the Disappearance of Congo Red From the Blood of Non amyloid Subjects and Patients With Amyloidosis *J Clin Investigation* 27 111 117 (Jan) 1948

### **Bile in Urine**

Hawkinson A Watson C J and Turner R H Modification of Harrison's Test for Bilirubin in Urine Especially Suited for Mass and Serial Use *J A M A* 129 514 515 (Oct) 1945

### **Cephalin Flocculation**

Hanger F M Serological Differentiation of Obstructive From Hepatogenous Jaundice by Flocculation of Cephalin Cholesterol Emulsion *J Clin Investigation* 18 61 69 (May) 1939

### **Cholinesterase Activity**

Vorhaus L J H and Kark R M Measurement of Serum Cholinesterase Activity A Useful Laboratory Tool in the Study of Diseases of the Liver and Biliary System To be published *Medical Nutrition Laboratory Report No 56 June 8 1949*

### **Galactose**

Coleher H Patek A J Jr and Kendall F L Galactose Disappearance From the Blood Stream Calculation of the Galactose Removal Constant and Its Application as a Test for Liver Function *J Clin Investigation* 25 768 775 (Sept) 1946

### **Hemosiderin in Urine**

Rous P Urinary Siderosis Hemosiderin Granules in the Urine as an Aid to Diagnosis of Pernicious Anemia Hemochromatosis and Other Diseases Causing Siderosis *J Exper Med* 28 645 658 (Nov) 1918

### **Hippuric Acid**

Quick A J The Clinical Application of the Hippuric Acid and Prothrombin Test *Am J Clin Path* 10 22-23 (March) 1940

### **Phosphatase**

Cusman A B and Hanger F M Jr Differential Diagnosis of Jaundice by Combined Serum Phosphatase Determination and Cephalin Flocculation Test *M Clin North America* 25 531 548 (Mar) 1941

### **Portal Pressure**

Bean W B Paul W D and Franklin M Preliminary Studies on an Indirect Method for Determining Portal Pressures *J Clin Investigation* 28 769 770 (July) 1949

### **Prothrombin Time**

Aggeler P M Howard J Lucia S P Clark W and Astaff A Standardization of the Quick Prothrombin Test *Blood* 1 2 0 203 (March) 1946

Quick A J The Clinical Significance of Prothrombin as a Factor in Hemorrhage *Pennsylvania M J* 43 125 130 (Oct) 1939

Unger P N and Shapiro S Prothrombin Response to the Parenteral Administration of Large Doses of Vitamin K in Subjects With Normal Liver Function and in Cases of Liver Disease A Standardized Test for the Estimation of Hepatic Function *J Clin Investigation* 27 29 41 (Jan) 1948

### **Thymol Turbidity**

Dekema H The Thymol Turbidity Test and the Cephalin Cholesterol Flocculation Test *Nederl tijdschr v geneesk* 93 (II) 194 82 (1949)

Ruerva J dela and Popper H Standardized Reagent for Thymol Turbidity Test *J Lab & Clin Med* 34 (6) 87 880 (June) 1949

Shank R E and Hoagland C L A Modified Method for the Quantitative Determination of the Thymol Turbidity Reaction of Serum *J Biol Chem* 162 133 138 (Jan) 1946

### **Urobilinogen**

Wilson T M and Davidson L S P Ehrlich's Aldehyde Test for Urobilinogen *British M J* 1949 884 889 (May) 1949

# Gold Pressor Tests

Markert N J and Lang J D Physicians Handbook ed 4 1410 Calif 1940  
University Medical Publishers page 28

## Porphyrins

Welker M L The Porphyrins New England J Med 32 11 19 (Jan) 1940  
Nesbitt S Acute Porphyrins J A M A 124 86 94 (Jan) 1944  
Katzon C J Schwartz S and Hawkenson V Studies of the Uroporphyrins II  
Further Studies of the Porphyrins of the Urine Feces Bile and Liver in Cases of  
Porphyrin With Particular Reference to a Waldenström Type Porphyrin Behaving  
as an Entity on the Tswett Column J Biol Chem 157 345 362 (Jan) 1940

## SUBSTANCES OF OCCASIONAL CLINICAL IMPORTANCE MISCELLANEOUS REFERENCES

### Alcohol

Harger R N Raney, B B Bridwell E G and Hatcher M F The Partition Ratio of  
Alcohol Between Air and Water Urine and Blood Estimation and Identification of  
Alcohol in These Liquids From Analysis of Air Equilibrated With Them J Biol  
Chem 123 104 14 (March) 1940

### Heparin

Astrup I On the Determination of Heparin in Blood Plasma and Urine Acta pharmacol  
et toxicol 3 (2) 165 178 1940

### Mercury

Van Susteren H Warmoltz N and Wiggerink G L Method for Determining the  
Mercury Content of Urine Phillips F H Re 31 119 120 1940

### Nicotine

Wolff William A, Hawkins Marina H and Giles W E The Spectrophotometric Estimation  
of Nicotine in Blood J Biol Chem 175 8 831 (9) 1948

### Phenols

Schmidt L O Urinary Phenol IV The Simultaneous Determination of Phenol and  
p Cresol in Urine J Biol Chem 179 11 16 (May) 1949  
Theis R C and Benedict S P The Determination of Phenols in the Blood J Biol  
Chem 61 64 74 (Aug) 1944  
Verra Mario Urinary Phenols I Method of Determination Urinary Phenols II The  
Significance in Normal and Pathological Conditions Am J Clin Path 12 5 533  
580 589 (Nov) 1940

### Procaine

Tung H C Coon J M and Conway A C A Spectrophotometric Method for Determination  
of Procaine and p Aminobenzoic Acid J Lab & Clin Med 34 (6) 8 7873  
(June) 1949

### Salicylates

Verra Mario and Jacobs Mildred D A Simple Method for the Determination of  
Salicylates in Blood J Lab & Clin Med 32 (10) 1 8 1 83 (Oct) 1941

### Streptomycin

Mashall E K, Jr Blanchard I C and Buhle Emmett I Colorimetric Methods for  
Determination of Streptomycin J Pharm & Exper Therap 90 (4) 367 374 (Aug)  
1944

### Sulfonamides

Bratton A C and Marshall F K A New Coupling Component for Sulfanilamide  
Determination J Biol Chem 178 375 0 (Mar) 1939

**Calculation**

- 1 Calculate eosinophils uric acid and creatinine as described under those methods
- 2 Calculate the uric acid creatinine ratio for both specimens of urine according to the formula

uric acid/creatinine ratio =

$$\frac{\text{mg uric acid/100 ml urine}}{\text{mg creatinine/100 ml urine}}$$

- 3 Calculate the percentage change of eosinophils and uric acid/creatinine ratio between the pre injection and post injection values according to the formula
- change =

$$\frac{(\text{Pre injection value} - \text{post injection value}) \times 100}{(\text{Pre injection value})}$$

**Interpretation**

**Eosinophil Count** Normal subjects have a decrease of 70% or more 50% is the lower limit of normal Patients with Addison's disease show little or no drop in eosinophil count

**Uric Acid Creatinine Ratio** Normal subjects will increase their uric acid creatinine ratio by approximately 100 percent average Patients with Addison's disease show approximately a 40 percent increase An increase of over 50 percent is strong evidence against adrenal insufficiency

**2 References to Endocrine Gastric and Other Functional Tests****Water and Salt Tests of Adrenal Sufficiency**

- Robinson, F J Power M H and Kepler E J Two New Procedures to Assist in the Recognition and Exclusion of Addison's Disease *Proc Staff Meet, Mayo Clin* 16 577 583 (Sept) 1947
- Cutler H H Power M H and Walder R M Concentration of Chloride Sodium and Potassium in Urine and Blood Their Diagnostic Significance in Adrenal Insufficiency *J A M A* 111 117 122 (Jan) 1938

**Adrenal Hormones**

- Fraser R W Forbes A P Albright F Sulkowitch S and Reifenshtein E C Jr Colorimetric Assay of 17 Ketosteroids in Urine *J Clin Endocrinol* 1 34 55 (March) 1941
- Talbot N B Saltzman F Wixon R L and Wolfe J K Colorimetric Assay of Urinary Corticosteroid like Substances *J Biol Chem* 160 535 546 (Oct) 1945

**Pancreas**

- Fraser R Albright F and Smith P H The Value of Glucose Tolerance Tests Insulin Tolerance Test and the Glucose Insulin Tolerance Test in the Diagnosis of Endocrinological Disorders of Glucose Metabolism *J Clin Endocrinol* 1 29, 306 (April) 1941
- Comfort M W and Osterberg A E Serum Amylase and Serum Lipase in the Diagnosis of Diseases of the Pancreas *M Clin North America* 24 1137 1149 (July) 1940
- Beazell J M Schmidt C R and Ivy A C The Diagnosis and Treatment of Achylia Pancreatica *J A M A* 116 736 739 (June) 1941

**Parathyroids**

- Albright F The Parathyroids—Physiology and Therapeutics *J A M A* 117 5 7 533 (Aug) 1941

**Gastric Function**

- Glenn P M Histamine Stimulated Fractional Gastric Analysis The Diagnostic Value of Total Secretion *Gastroenterology* 6 409 416 (May) 1946
- Hollander F The Insulin Test for the Presence of Intact Nerve Fibers After Vagotomy Operations for Peptic Ulcer *Gastroenterology* 7 607 614 (Dec) 1946
- Vanant F R and Alvarez W C Calculating the Diagnostic Value of Gastric Analysis A Study in the Methodology of Diagnosis *Am J Digest Dis & Nutrition* 2 466 472 (Oct) 1935

Yule C Udny An Introduction to the Theory of Statistics London England 1940  
Charles Griffin and Company Ltd  
*This text covers the theory of statistics has an excellent bibliography and gives answers to problems presented in the text*

## **2 General Handbooks**

Handbook of Chemistry and Physics Cleveland Ohio 1947 Chemical Rubber Publishing Company  
Handbook of Chemistry edited by Lange Adolph N Sandusky Ohio 1946 Handbook Publishers Inc  
Van Nostrand's Chemical Annual edited by Olsen John C New York 1934 D Van Nostrand and Company Inc  
Marks I one Mechanical Engineer's Handbook New York 1943 McGraw Hill Book Co Inc

## SECTION X

# STATISTICAL METHODS

---

### A REFERENCES TO STANDARD WORKS

#### 1 Textbooks of Statistics

- Buros Oscar Krisen Research and Statistical Methodology New Brunswick N J 1948  
Rutgers University Press  
An annotated bibliography
- Buros Oscar Krisen The Second Yearbook of Research and Statistical Methodology--  
Books and Reviews Highland Park N J 1941 The Gryphon Press  
An annotated bibliography
- Cramer H Mathematical Method of Statistics Princeton N J 1946 Princeton Uni-  
versity Press  
Very theoretical
- Croxton F F and Cowden D J Applied General Statistics New York 1946 Prentice  
Hall Inc  
Gives a good presentation of statistical methods and their application in various fields
- Ezekiel Mordecai Methods of Correlation Analysis New York 1941 John Wiley & Sons  
Inc  
Covers that part of the field of statistics which is concerned with studying the relations  
between variables
- Fisher Donald A Statistical Methods for Research Workers Edinburgh Scotland 1941  
Oliver and Boyd, Ltd. pp xiii 339  
A comprehensive text with emphasis on biology
- Fisher Ronald A and Yates Frank Statistical Tables for Biological Agricultural and  
Medical Research, Edinburgh Scotland 1938 Oliver & Boyd Ltd pp 90
- Goodwin H M Elements of the Precision of Measurements and Graphical Method  
McGraw Hill Book Co N Y 1939  
An excellent account of the subject of making measurements of errors involved, and of the  
statistical methods of handling instrumental data There is also a section dealing  
with graphical methods Illustrative problems are given together with examples
- Kenney John F Mathematics of Statistics New York 1939 D Van Nostrand Company  
Inc  
An excellent dictionary
- Kurtz Albert R and Edgeston Harold A Statistical Dictionary of Terms and Symbols  
New York 1939 John Wiley & Sons Inc pp xiii 191  
An excellent dictionary
- Mainland Donald The Treatment of Clinical and Laboratory Data An Introduction  
to Statistical Ideas and Method for Medical and Dental Workers 1948 Edinburgh  
Scotland Oliver & Boyd Ltd pp xi 40  
This book has been prepared for the express purpose of aiding the clinician untrained in  
statistics in the evaluation of the significance of laboratory observations It is an  
appeal to understanding rather than to mathematical derivation
- Rider, P R An Introduction to Modern Statistical Methods New York, 1939, John Wiley  
and Sons Inc  
A practical statistical book that has proved useful in handling biological data.
- Snedecor George W Statistical Methods Applied to Experiments in Agriculture and  
Biology ed 3 Ames Iowa, 1940, Iowa State College Press pp xv 400  
A standard textbook

## B CONVERSION BETWEEN DIFFERENT SYSTEMS OF WEIGHTS AND MEASURES

### 1 Length—Metric System

| Angstrom Units<br>Å | Millimicrons<br>mμ | Microns<br>μ     | Millimeters<br>mm | Centimeters<br>cm          |
|---------------------|--------------------|------------------|-------------------|----------------------------|
| 1                   | 0.1                | 10               | 10 <sup>-6</sup>  | 10 <sup>-5</sup>           |
| 10                  | 1                  | 10 <sup>-1</sup> | 10 <sup>-5</sup>  | 10 <sup>-4</sup>           |
| 10                  | 10 <sup>-1</sup>   | 1                | 10 <sup>-4</sup>  | 10 <sup>-3</sup>           |
| Millimeters<br>mm   | Centimeter<br>cm   | Meter<br>m       | Kilometers<br>km  | U S<br>Equivalents         |
| 1                   | 0.1                | 0.001            | 10 <sup>-6</sup>  | 0.03937 inch               |
| 10                  | 1                  | 0.01             | 10 <sup>-5</sup>  | 0.3937 inch                |
| 1000                | 100                | 1                | 10 <sup>-3</sup>  | 39.37 inch                 |
| 10                  | 10                 | 1000             | 1                 | 3.9093 feet<br>0.6137 mile |

### 2 Length—U S System

| Inches<br>in   | Feet<br>ft | Yards<br>yd | Miles<br>mi | Metric<br>Equivalent |
|----------------|------------|-------------|-------------|----------------------|
| 1              | 0.08333    | 0.000778    |             | 2.54001 cm           |
| 1 <sup>2</sup> | 1          | 0.333       |             | 0.91401 m            |
| 36             | 3          | 1           | 0.00069     | 0.9140 m             |
| 63 360         | 5 280      | 1 600       | 1           | 1.6093 km            |

### 3 Area—Metric System

| Square<br>Millimeters<br>mm <sup>2</sup> | Square<br>Centimeters<br>cm <sup>2</sup> | Square<br>Meter<br>m <sup>2</sup> | Square<br>Kilometers<br>km <sup>2</sup> | U S<br>Equivalents |
|--|--|-----------------------------------|---|--------------------|
| 1  | 0.01                                     | 10 <sup>-6</sup>                  | 10 <sup>-10</sup>                       | 0.00155 sq in.     |
| 100                                      | 1  | 10 <sup>-4</sup>                  | 10 <sup>-8</sup>                        | 0.15500 sq in.     |
| 10                                       | 10.000                                   | 1                                 | 10 <sup>-6</sup>                        | 10.764 sq ft       |
| 10                                       | 100                                      | 100                               | 1                                       |                    |

### 4 Area—U S System

| Square Inches<br>sq in | Square Feet<br>sq ft | Square Yards<br>sq yd | Acres<br>A | Square Miles<br>sq mi | Metric<br>Equivalent   |
|------------------------|----------------------|-----------------------|------------|-----------------------|------------------------|
| 1                      | 0.006944             |                       | —          | —                     | 6.4516 cm <sup>2</sup> |
| 144                    | 1                    | 0.1111                | —          | —                     | 0.093 m <sup>2</sup>   |
| 1 96                   | 9                    | 1                     | —          | —                     | 0.8361 m <sup>2</sup>  |
|                        | 47,900               | 4940                  | 1          | 0.001563              | 404.69 m <sup>2</sup>  |
|                        |                      | 309 60                | 640        | 1                     | 2.5999 km <sup>2</sup> |

### 5 Volume—Metric System

| Cubic Millimeters<br>mm <sup>3</sup> | Cubic Centimeters<br>cm <sup>3</sup> | Cubic Meters<br>m <sup>3</sup> | U S Equivalent               |
|--------------------------------------|--------------------------------------|--------------------------------|------------------------------|
| 1                                    | 0.001                                | 10 <sup>-6</sup>               | 0.000061 cu in.              |
| 1000                                 | 1                                    | 10 <sup>-3</sup>               | 0.0610 cu in.                |
| 10 <sup>6</sup>                      | 10                                   | 1                              | 3.7854 cu ft<br>1.3564 cu yd |

### 6 Volume—U S System

| Cubic Inches<br>cu in | Cubic Feet<br>cu ft | Cubic Yards<br>cu yd | Metric Equivalent       |
|-----------------------|---------------------|----------------------|-------------------------|
| 1                     | 0.00058             |                      | 16.39 cm <sup>3</sup>   |
| 17 28                 | 1                   | 0.0357               | 0.028317 m <sup>3</sup> |
| 46 086                | 27                  | 1                    | 0.646 m <sup>3</sup>    |

# SECTION XI

## MISCELLANEOUS DATA

### A ATOMIC WEIGHTS

#### References

Table of International Atomic Weights 1946

1947 Summary of Nuclear Data (Aug 1946 Nov 1947) Compiled by Nuclear Data Committee of Clinton National Laboratory Nucleonics 2 Part 2 (May) 1948

TABLE 2 — ATOMIC WEIGHTS

| ELEMENT    | SYMBOL | ATOMIC NO | ATOMIC WT | ELEMENT      | SYMBOL | ATOMIC NO | ATOMIC WT |
|------------|--------|-----------|-----------|--------------|--------|-----------|-----------|
| Actinium   | Ac     | 81        | 1 ca      | Neodymium    | Nd     | 60        | 144.7     |
| Alabamine  | Ab     | 85        | 1 ca      | Neon         | Ne     | 10        | 20.183    |
| Aluminum   | Al     | 13        | 26.98     | Neptunium    | Np     | 93        | 237 ca.   |
| Americium  | Am     | 95        | 241 ca    | Nickel       | Ni     | 28        | 58.69     |
| Antimony   | Sb     | 51        | 121.76    | Nitrogen     | N      | 7         | 14.008    |
| Arsenic    | As     | 33        | 74.92     | Omium        | Os     | 76        | 190       |
| Astatine   | At     | 85        | 210 ca    | Oxygen       | O      | 8         | 16.0000   |
| Barium     | Ba     | 56        | 137.36    | Palladium    | Pd     | 46        | 106       |
| Beryllium  | Be     | 4         | 9.01      | Phosphorus   | P      | 15        | 30.97     |
| Bismuth    | Bi     | 83        | 208.98    | Platinum     | Pt     | 78        | 195.08    |
| Boron      | B      | 5         | 10.81     | Plutonium    | Pu     | 94        | 239 ca    |
| Bromine    | Br     | 35        | 79.90     | Polonium     | Po     | 84        | 209 ca    |
| Cadmium    | Cd     | 48        | 112.41    | Potassium    | K      | 19        | 39.098    |
| Calcium    | Ca     | 20        | 40.08     | Praseodymium | Pr     | 59        | 140.91    |
| Carbon     | C      | 6         | 12.01     | Protactinium | Pa     | 91        | 231       |
| Cerium     | Ce     | 58        | 140.13    | Radium       | Ra     | 88        | 226       |
| Cesium     | Cs     | 55        | 132.91    | Radon        | Rn     | 86        | 222       |
| Chlorine   | Cl     | 17        | 35.457    | Rhenium      | Rh     | 45        | 186.21    |
| Chromium   | Cr     | 24        | 52.01     | Rhodium      | Rd     | 37        | 101.07    |
| Cobalt     | Co     | 27        | 58.94     | Rubidium     | Rb     | 37        | 85.47     |
| Columbium  | Cb     | 41        | 94.01     | Ruthenium    | Ru     | 44        | 101.07    |
| Copper     | Cu     | 29        | 63.55     | Samarium     | Sm     | 62        | 150.35    |
| Curium     | Cm     | 96        | 247 ca.   | Scandium     | Sc     | 21        | 44.96     |
| Deuterium  | D      | 2         | 2.014     | Selenium     | Se     | 34        | 78.96     |
| Erbium     | Er     | 68        | 167.27    | Silicon      | Si     | 14        | 28.086    |
| Europium   | Eu     | 63        | 151.96    | Silver       | Ag     | 47        | 107.87    |
| Fluorine   | F      | 9         | 18.998    | Sodium       | Na     | 11        | 22.989    |
| Francium   | Fr     | 87        | 223 ca    | Strontium    | Sr     | 38        | 87.62     |
| Gallium    | Gal    | 69        | 69.72     | Sulfur       | S      | 16        | 32.06     |
| Germanium  | Ge     | 32        | 72.64     | Tantalum     | Ta     | 73        | 180.95    |
| Gold       | Au     | 79        | 196.97    | Technetium   | Tc     | 43        | 98 ca     |
| Hafnium    | Hf     | 72        | 178.49    | (Manganese)  |        |           |           |
| Helium     | He     | 2         | 4.003     | Tellurium    | Te     | 52        | 127.6     |
| Holmium    | Ho     | 67        | 164.93    | Terbium      | Tb     | 65        | 158.93    |
| Hydrogen   | H      | 1         | 1.0080    | Thallium     | Tl     | 81        | 204.38    |
| Indium     | In     | 49        | 114.71    | Thorium      | Th     | 90        | 232.04    |
| Iodine     | I      | 53        | 126.90    | Thulium      | Tm     | 69        | 168.93    |
| Iridium    | Ir     | 77        | 192.22    | Tin          | Sn     | 50        | 118.71    |
| Iron       | Fe     | 26        | 55.85     | Titanium     | Ti     | 22        | 47.88     |
| Krypton    | Kr     | 36        | 83.80     | Tungsten     | W      | 74        | 183.85    |
| Lanthanum  | La     | 57        | 138.91    | Uranium      | U      | 92        | 238.03    |
| Lead       | Pb     | 82        | 207.2     | Vanadium     | V      | 23        | 50.94     |
| Lithium    | Li     | 3         | 6.941     | Virginium    | Vi     | 87        | 223 ca.   |
| Lutecium   | Lu     | 71        | 174.97    | Xenon        | Xe     | 54        | 131.30    |
| Magnesium  | Mg     | 12        | 24.31     | Ytterbium    | Yb     | 70        | 173.05    |
| Manganese  | Mn     | 25        | 54.94     | Yttrium      | Y      | 39        | 88.91     |
| Mercury    | Hg     | 80        | 200.59    | Zinc         | Zn     | 30        | 65.38     |
| Molybdenum | Mo     | 42        | 95.94     | Zirconium    | Zr     | 40        | 91.22     |
|            |        |           |           | Flement 61   | —      | 61        | 145 ca    |
|            |        |           |           | Flement 86   | —      | 86        | —         |

## B CONVERSION BETWEEN DIFFERENT SYSTEMS OF WEIGHTS AND MEASURES

### 1 Length—Metric System

| Angstrom Units<br>$\text{\AA}$ | Millimicrons<br>$\text{m}\mu$ | Microns<br>$\mu$    | Millimeters<br>$\text{mm}$ | Centimeters<br>$\text{cm}$ |
|--------------------------------|-------------------------------|---------------------|----------------------------|----------------------------|
| 1                              | 0.1                           | 10                  | 10 <sup>-3</sup>           | 10 <sup>-2</sup>           |
| 10                             | 1                             | 10 <sup>1</sup>     | 10 <sup>-2</sup>           | 10 <sup>-1</sup>           |
| 10 <sup>4</sup>                | 10 <sup>3</sup>               | 1                   | 10 <sup>-3</sup>           | 10 <sup>-2</sup>           |
| Millimeters<br>$\text{mm}$     | Centimeter<br>$\text{cm}$     | Meter<br>$\text{m}$ | Kilometers<br>$\text{km}$  | U. S. Equivalents          |
| 1                              | 0.1                           | 0.001               | 10 <sup>-4</sup>           | 0.03937 inch               |
| 10                             | 1                             | 0.01                | 10 <sup>-3</sup>           | 0.393 inch                 |
| 1000                           | 100                           | 1                   | 10 <sup>-3</sup>           | 39.3 inch<br>1.3123 feet   |
| 10                             | 10                            | 1000                | 1                          | 0.621372 mile              |

### 2 Length—U. S. System

| Inches<br>$\text{in}$ | Feet<br>$\text{ft}$ | Yards<br>$\text{yd}$ | Miles<br>$\text{mi}$ | Metric Equivalent |
|-----------------------|---------------------|----------------------|----------------------|-------------------|
| 1                     | 0.08333             | 0.078                |                      | 2.5401 cm         |
| 1                     | 1                   | 0.33333              |                      | 0.304801 m        |
| 36                    | 3                   | 1                    | 0.00068              | 0.91440 m         |
| 63.60                 | 5.30                | 1.60                 | 1                    | 1.6093 km         |

### 3 Area—Metric System

| Square Millimeters<br>$\text{mm}^2$ | Square Centimeters<br>$\text{cm}^2$ | Square Meters<br>$\text{m}^2$ | Square Kilometers<br>$\text{km}^2$ | U. S. Equivalents |
|-------------------------------------|-------------------------------------|-------------------------------|------------------------------------|-------------------|
| 1                                   | 0.01                                | 10 <sup>-6</sup>              | 10 <sup>-10</sup>                  | 0.00155 sq in     |
| 100                                 | 1                                   | 10 <sup>-4</sup>              | 10 <sup>-8</sup>                   | 0.1550 sq in      |
| 10 <sup>6</sup>                     | 10 000                              | 1                             | 10 <sup>-6</sup>                   | 10.7639 sq ft     |
| 10 <sup>10</sup>                    | 10 <sup>11</sup>                    | 10 <sup>4</sup>               | 1                                  |                   |

### 4 Area—U. S. System

| Square Inches<br>$\text{sq in}$ | Square Feet<br>$\text{sq ft}$ | Square Yards<br>$\text{sq yd}$ | Acres<br>$\text{A}$ | Square Miles<br>$\text{sq mi}$ | Metric Equivalent      |
|---------------------------------|-------------------------------|--------------------------------|---------------------|--------------------------------|------------------------|
| 1                               | 0.006944                      |                                | —                   | —                              | 6.4516 cm <sup>2</sup> |
| 144                             | 1                             | 0.1111                         | —                   | —                              | 0.0929 m <sup>2</sup>  |
| 1 440                           | 9                             | 1                              | —                   | —                              | 0.8361 m <sup>2</sup>  |
|                                 | 47.910                        | 4840                           | 1                   | 0.001563                       | 4046.86 m <sup>2</sup> |
|                                 |                               | 309 600                        | 640                 | 1                              | 2.59 km <sup>2</sup>   |

### 5 Volume—Metric System

| Cubic Millimeters<br>$\text{mm}^3$ | Cubic Centimeters<br>$\text{cm}^3$ | Cubic Meters<br>$\text{m}^3$ | U. S. Equivalent              |
|------------------------------------|------------------------------------|------------------------------|-------------------------------|
| 1                                  | 0.001                              | 10 <sup>-9</sup>             | 0.0000001 cu in               |
| 1000                               | 1                                  | 10 <sup>-6</sup>             | 0.0610 cu in                  |
| 10 <sup>9</sup>                    | 10                                 | 1                            | 35.3147 cu ft<br>1.3579 cu yd |

### 6 Volume—U. S. System

| Cubic Inches<br>$\text{cu in}$ | Cubic Feet<br>$\text{cu ft}$ | Cubic Yards<br>$\text{cu yd}$ | Metric Equivalent       |
|--------------------------------|------------------------------|-------------------------------|-------------------------|
| 1                              | 0.00058                      |                               | 16.39 cm <sup>3</sup>   |
| 17 28                          | 1                            | 0.0357                        | 0.028317 m <sup>3</sup> |
| 46 556                         | 27                           | 1                             | 0.6456 m <sup>3</sup>   |



## 7 Capacity—Metric System

| Milliliters<br>ml | Liters<br>l | U S<br>Equivalent   |
|-------------------|-------------|---|
| 1                 | 0.001       | $\left\{ \begin{array}{l} 16 \frac{31}{32} \text{ minims} \\ 0.0610 \text{ cu in} \\ 0.018 \text{ fl dram} \end{array} \right.$ |
| 1000              | 1           | $\left\{ \begin{array}{l} 33 \frac{8147}{100000} \text{ fl ounce} \end{array} \right.$  |

## 8 Capacity—U S System—Liquid

| Gills<br>gi | Pints<br>pt | Quarts<br>qt | Gallons | Metric<br>Equivalent |
|-------------|-------------|--------------|---------|----------------------|
| 1           | 0.5         | 0.125        | 0.03125 | 118.297 ml           |
| 4           | 1           | 0.5          | 0.125   | 0.473 l              |
| 8           | 2           | 1            | 0.5     | 0.946 l              |
| 32          | 8           | 4            | 1       | 3.785 l              |

## 9 Capacity—U S System—Dry

| Pints<br>pt | Quarts<br>qt | Peck<br>pk | Bushels<br>bu | Metric<br>Equivalent |
|-------------|--------------|------------|---------------|----------------------|
| 1           | 0.5          | 0.0625     | 0.015625      | 0.5506 l             |
| 2           | 1            | 0.125      | 0.03125       | 1.101 l              |
| 16          | 8            | 1          | 0.25          | 8.806 l              |
| 64          | 32           | 4          | 1             | 35.238 l             |

## 10 Mass—Metric System

| Milligrams<br>mg | Grams<br>gm | Kilograms<br>kg | U S<br>Equivalent  |
|------------------|-------------|-----------------|--|
| 1                | 0.001       | 0.001           | 0.01543 grain  |
| 1000             | 1           | 0.001           | 15.4336 grain  |
| 10 <sup>6</sup>  | 1000        | 1               | $\left\{ \begin{array}{l} 0.03224 \text{ oz av} \\ 2.2046 \text{ lb av} \end{array} \right.$ |

## 11 Mass—U S System—Avoirdupois Commercial

| Grain<br>gr | Drams<br>dr av | Ounces<br>oz av | Pounds<br>lb av | Ton<br>short tn | Metric<br>Equivalent |
|-------------|----------------|-----------------|-----------------|-----------------|----------------------|
| 1           | 0.015625       | —               | —               | —               | 0.0648 gm            |
| ~437.5      | 1              | 0.0625          | —               | —               | 1.13 gm              |
| 437.5       | 16             | 1               | 0.0625          | —               | 9.349 gm             |
| 7000        | 256            | 16              | 1               | 0.0005          | 0.4536 kg            |
| —           | —              | 16 000          | 1000            | 1               | 907.18 kg            |

## 12 U S Apothecaries Fluid Measure

| Minim<br>min | Fluid Drams<br>fl dr | Fluid Ounces<br>fl oz | Pints<br>pt | Metric<br>Equivalent |
|--------------|----------------------|-----------------------|-------------|----------------------|
| 1            | 0.01666              | 0.0009                | —           | 0.0616 ml            |
| 60           | 1                    | 0.125                 | —           | 3.69 ml              |
| 480          | 8                    | 1                     | 0.0625      | 9.573 ml             |
| 7680         | 128                  | 16                    | 1           | 37.85 l              |

## 13 U S Apothecaries Weight

| Grain<br>gr | Scruples<br>s ap | Drams<br>dr ap | Ounces<br>oz ap | Pounds<br>lb ap | Metric<br>Equivalent |
|-------------|------------------|----------------|-----------------|-----------------|----------------------|
| 1           | 0.03             | 0.0166         | —               | —               | 64.91 mc             |
| 60          | 1                | 0.125          | 0.04167         | —               | 1.9 gm               |
| 480         | 8                | 1              | 0.125           | 0.0104          | 3.89 gm              |
| 480         | 8                | 8              | 1               | 0.08333         | 31.103 gm            |
| ~60         | ~88              | 96             | 1               | 1               | 33.94 gm             |

## 14 Time

| Seconds<br>sec | Minutes<br>min | Hours<br>hr | Days<br>da. | Years<br>yr |
|----------------|----------------|-------------|-------------|-------------|
| 1              | 0.01667        | 0.000 778   |             | -           |
| 60             | 1              | 0.166       |             | -           |
| 3 600          | 60             | 1           | 0.011667    |             |
| 86 400         | 1440           | 24          | 1           | -           |
| -              |                |             | 365.4       | 1 common    |
| -              |                |             | 366.56      | 1 a leap    |

## 7 Capacity—Metric System

| M illi ters<br>ml | L iters<br>l | U S<br>Equ alent                     |
|-------------------|--------------|--------------------------------------|
| 1                 | 1001         | { 16 31 min ms<br>00610 cu           |
| 1000              | 1            | { 0.26 fl dram<br>33.8147 fl ounce s |

## 8 Capacity—U S System—Liquid

| G ill<br>gi | Flint<br>pt | Quarts<br>qt | Gallo s | Metre<br>Equ alent |
|-------------|-------------|--------------|---------|--------------------|
| 1           | 0.5         | 0.125        | 0.03125 | 118.9 ml           |
| 4           | 1           | 0.5          | 0.125   | 0.473 l            |
| 8           |             | 1            | 0.5     | 0.9463 l           |
| 32          | 8           | 4            | 1       | 8.45 l             |

## 9 Capacity—U S System—Dry

| Pints<br>pt | Quart<br>qt | Pecks<br>pk | Bushels<br>bu | Metre<br>Equ alent |
|-------------|-------------|-------------|---------------|--------------------|
| 1           | 0.5         | 0.0625      | 0.015625      | 0.5506 l           |
|             | 1           | 0.125       | 0.03125       | 1.1012 l           |
| 16          | 8           | 1           | 0.25          | 8.8096 l           |
| 64          | 32          | 4           | 1             | 35.233 l           |

## 10 Mass—Metric System

| Milli grams<br>mg | Grams<br>gm | Kilograms<br>kg  | U S<br>Equ alent              |
|-------------------|-------------|------------------|-------------------------------|
| 1                 | 1001        | 10 <sup>-6</sup> | 0.01543 gr n                  |
| 1000              | 1           | 0.001            | 1.543236 gr n                 |
|                   |             |                  | { 0.035 74 oz av<br>0.46 fl n |
| 10                | 1000        | 1                |                               |

## 11 Mass—U S System—Avoirdupois Commercial

| Grains<br>gr | Drams<br>dr av | Ounces<br>oz a | Pounds<br>lb a | Tons<br>ton | Metre<br>Equ alent |
|--------------|----------------|----------------|----------------|-------------|--------------------|
| 1            | 0.0625         |                |                |             | 0.0648 gm          |
| 7000         | 1              | 0.0625         |                |             | 17.18 gm           |
| 437.5        | 16             | 1              | 0.0625         |             | 93.47 gm           |
| 283.5        | 8              | 16             | 1              | 0.00045     | 0.476 kg           |
|              |                | 16000          | 1000           | 1           | 907.18 kg          |

## 12 U S Apothecaries Fluid Measure

| Minims<br>min | Fluid Drams<br>fl dr | Fluid Ounces<br>fl oz | Pints<br>pt | Metre<br>Equ alent |
|---------------|----------------------|-----------------------|-------------|--------------------|
| 1             | 1666                 | 00.08                 |             | 0.0616 ml          |
| 60            | 1                    | 0.125                 |             | 3.697 ml           |
| 480           | 8                    | 1                     | 0.0625      | 9.463 ml           |
| 640           | 16                   | 16                    | 1           | 3.785 l            |

## 13 U S Apothecaries Weight

| Grains<br>gr | Scruples<br>s ap | Drams<br>dr ap | Ounces<br>oz ap | Pounds<br>lb ap | Metre<br>Equ alent |
|--------------|------------------|----------------|-----------------|-----------------|--------------------|
| 1            | 0.03             | 0.0156         |                 |                 | 64.8 mg            |
| 60           | 1                | 0.333          | 0.0416          |                 | 1.90 gm            |
| 60           | 3                | 1              | 0.125           | 0.0104          | 3.494 g            |
| 480          | 24               | 8              | 1               | 0.03133         | 31.103 gm          |
| 5760         | 288              | 96             | 1               | 1               | 373.24 gm          |

Time

| Seconds<br>sec | Minutes<br>min | Hours<br>hr | Days<br>da | Years<br>yr |
|----------------|----------------|-------------|------------|-------------|
| 1              | 0.01667        | 0.000 78    | -          | -           |
| 60             | 1              | 0.0166      | -          | -           |
| 3 600          | 60             | 1           | 0.0166     | -           |
| 86 400         | 1 440          | 24          | 1          | -           |
| -              | -              | -           | 36 547     | 1 common    |
| -              | -              | -           | 36 526     | 1 sidereal  |

## 7 Capacity—Metric System

| Milliliters<br>ml | Liters<br>l | U S<br>Equivalent                      |
|-------------------|-------------|--|
| 1                 | 0.001       | { 16.31 minims<br>{ 0.0010 cu in       |
| 1000              | 1           | { 0.018 fl dram<br>{ 33.8147 fl ounces |

## 8 Capacity—U S System—Liquid

| Gill<br>gi | Pints<br>pt | Quart<br>qt | Gallons | Metric<br>Equivalent |
|------------|-------------|-------------|---------|----------------------|
| 1          | 0.5         | 0.125       | 0.03125 | 118.3 ml             |
| 4          | 1           | 0.5         | 0.125   | 0.43 l               |
| 8          |             | 1           | 0.25    | 0.9463 l             |
| 32         | 8           | 4           | 1       | 3.785 l              |

## 9 Capacity—U S System—Dry

| Pints<br>pt | Quarts<br>qt | Pecks<br>pk | Bushels<br>bu | Metric<br>Equivalent |
|-------------|--------------|-------------|---------------|----------------------|
| 1           | 0.5          | 0.0625      | 0.015625      | 0.5506 l             |
|             | 1            | 0.125       | 0.03125       | 1.101 l              |
| 16          | 8            | 1           | 0.25          | 8.8096 l             |
| 64          | 32           | 4           | 1             | 35.239 l             |

## 10 Mass—Metric System

| Milligrams<br>mg | Grams<br>gm | Kilograms<br>kg  | U S<br>Equivalent                 |
|------------------|-------------|------------------|-----------------------------------|
| 1                | 0.001       | 10 <sup>-6</sup> | 0.01543 grain                     |
| 1000             | 1           | 0.001            | 15.4336 grain                     |
|                  |             |                  | { 0.035274 oz av<br>{ 2.2046 l av |
| 10 <sup>6</sup>  | 1000        | 1                |                                   |

## 11 Mass—U S System—Avoirdupois Commercial

| Grain<br>gr | Dram<br>dr av | Ounce<br>oz av | Pound<br>lb av | Ton<br>short tn | Metric<br>Equivalent |
|-------------|---------------|----------------|----------------|-----------------|----------------------|
| 1           | 0.036         | 0.0625         | 0.00044328     | —               | 0.0648 gm            |
| 43.75       | 1             | 0.0625         | —              | —               | 1.774 gm             |
| 43          | 16            | 1              | 0.0625         | —               | 28.349 gm            |
| 7000        | 16            | 1              | 0.00044328     | —               | 0.4536 kg            |
| —           | —             | 32000          | 0.00044328     | 1               | 907.18 kg            |

## 12 U S Apothecaries Fluid Measure

| Minimum<br>min | Fluid Drams<br>fl dr | Fluid Ounces<br>fl oz | Pints<br>pt | Metric<br>Equivalent |
|----------------|----------------------|-----------------------|-------------|----------------------|
| 1              | 0.016666             | 0.002083              | —           | 0.0616 ml            |
| 60             | 1                    | 0.125                 | —           | 3.69 ml              |
| 480            | 8                    | 1                     | 0.0625      | 29.573 ml            |
| 7680           | 128                  | 16                    | 1           | 378.5 l              |

## 13 U S Apothecaries Weight

| Grain<br>gr | Scruple<br>s ap | Dram<br>dr ap | Ounces<br>oz ap | Pounds<br>lb ap | Metric<br>Equivalent |
|-------------|-----------------|---------------|-----------------|-----------------|----------------------|
| 1           | 0.03            | 0.0156        | —               | —               | 64.8 mg              |
| 60          | 1               | 0.3           | 0.04167         | —               | 1.996 gm             |
| 480         | 3               | 1             | 0.125           | 0.0104          | 3.889 gm             |
| 480         | 24              | 8             | 1               | 0.09333         | 31.103 gm            |
| 5760        | 288             | 96            | 16              | 1               | 373.4 gm             |

## 14 Time

| seconds<br>sec | Minutes<br>min | Hours<br>hr | Days<br>da | Years<br>yr |
|----------------|----------------|-------------|------------|-------------|
| 1              | 0.01667        | 0.000 778   | -          | -           |
| 60             | 1              | 0.0166      | -          | -           |
| 3 600          | 60             | 1           | 0.0166     | -           |
| 86 400         | 1440           | 4           | 1          | -           |
| -              | -              | -           | 36 240     | 1 common    |
| -              | -              | -           | 365.26     | 1 #dereal   |

## C ELEMENTS LISTED IN ORDER OF ATOMIC NUMBER

TABLE 28

## ELEMENTS LISTED IN ORDER OF ATOMIC NUMBER

(For mass numbers and  $\infty$  abundance of stable isotope see International Tables of Stable Isotopes. For mass number and half life of unstable isotopes see 1947 Summary of Nuclear Data Nucleonics 2 Part 2 May 1948)

| ATOMIC NO | NAME       | SYMBOL | ATOMIC NO | NAME          | SYMBOL |
|-----------|------------|--------|-----------|---------------|--------|
| 0         | Neutron    | n      | 48        | Cadmium       | Cd     |
| 1         | Deuterium  | D      | 49        | Indium        | In     |
| 1         | Hydrogen   | H      | 50        | Tin           | Sn     |
| 2         | Helium     | He     | 51        | Antimony      | Sb     |
| 3         | Lithium    | Li     | 52        | Tellurium     | Te     |
| 4         | Beryllium  | Be     | 53        | Iodine        | I      |
| 5         | Boron      | B      | 54        | Xenon         | Xe     |
| 6         | Carbon     | C      | 55        | Cesium        | Cs     |
| 7         | Nitrogen   | N      | 56        | Barium        | Ba     |
| 8         | Oxygen     | O      | 57        | Lanthanum     | La     |
| 9         | Fluorine   | F      | 58        | Cerium        | Ce     |
| 10        | Neon       | Ne     | 59        | Praseodymium  | Pr     |
| 11        | Sodium     | Na     | 60        | Neodymium     | Nd     |
| 12        | Magnesium  | Mg     | 61        | Element 61    | —      |
| 13        | Aluminum   | Al     | 62        | Samarium      | Sm     |
| 14        | Silicon    | Si     | 63        | Europium      | Eu     |
| 15        | Phosphorus | P      | 64        | Gadolinium    | Gd     |
| 16        | Sulfur     | S      | 65        | Terbium       | Tb     |
| 17        | Chlorine   | Cl     | 66        | Dysprosium    | Dy     |
| 18        | Argon      | A      | 67        | Holmium       | Ho     |
| 19        | Potassium  | K      | 68        | Erbium        | Er     |
| 20        | Calcium    | Ca     | 69        | Thulium       | Tm     |
| 21        | Scandium   | Sc     | 70        | Ytterbium     | Yb     |
| 22        | Titanium   | Ti     | 71        | Lutecium      | Lu     |
| 23        | Vanadium   | V      | 72        | Hafnium       | Hf     |
| 24        | Chromium   | Cr     | 73        | Tantalum      | Ta     |
| 25        | Manganese  | Mn     | 74        | Tungsten      | W      |
| 26        | Iron       | Fe     | 75        | Rhenium       | Re     |
| 27        | Cobalt     | Co     | 76        | Osmium        | Os     |
| 28        | Nickel     | Ni     | 77        | Iridium       | Ir     |
| 29        | Copper     | Cu     | 78        | Platinum      | Pt     |
| 30        | Zinc       | Zn     | 79        | Gold          | Au     |
| 31        | Gallium    | Ga     | 80        | Mercury       | Hg     |
| 32        | Germanium  | Ge     | 81        | Thallium      | Tl     |
| 33        | Arsenic    | As     | 82        | Lead          | Pb     |
| 34        | Selenium   | Se     | 83        | Bismuth       | Bi     |
| 35        | Bromine    | Br     | 84        | Polonium      | Po     |
| 36        | Krypton    | Kr     | 85        | Astatine      | At     |
| 37        | Rubidium   | Rb     | 86        | —             | —      |
| 38        | Strontium  | Sr     | 87        | Francium      | Fr     |
| 39        | Yttrium    | Y      | 88        | Radium        | Ra     |
| 40        | Zirconium  | Zr     | 89        | Actinium      | Ac     |
| 41        | Niobium    | Nb     | 90        | Thorium       | Th     |
| 42        | Molybdenum | Mo     | 91        | Protoactinium | Pa     |
| 43        | Technetium | Tc     | 92        | Uranium       | U      |
| 44        | Ruthenium  | Ru     | 93        | Neptunium     | Np     |
| 45        | Rhodium    | Rh     | 94        | Plutonium     | Pu     |
| 46        | Palladium  | Pd     | 95        | Americium     | Am     |
| 47        | Silver     | Ag     | 96        | Curium        | Cm     |

## **D COMPOSITION OF BODY FLUIDS**

In metabolic studies it is sometimes helpful to find in a convenient form and in one place average values for various components of the body. The four accompanying tables (Tables 29, 30, 31 and 32) list for whole blood, plasma, urine and feces various constituents in several general categories. These categories are (1) water, (2) hydrogen ion concentration, (3) minerals, (4) nitrogen and protein derivatives, (5) carbohydrates and derivatives, (6) fats and derivatives, (7) vitamins, (8) hormones, (9) enzymes, (10) pigments, (11) gases, and (12) measurements of clinical interest.

It is hardly necessary to warn the reader that the term "normal range" is notoriously difficult to define, that the concentration of any substance may vary from person to person, that the effects of many factors such as age, sex, disease, previous diet, exercise and season may be large, and that no table can hope to be complete.



TABLE 29 COMPOSITION OF WHOLE BLOOD

(Listed below are normal ranges for various substances in whole blood all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake and to exercise.)

| SUBSTANCE                                 | UNITS                       | USUAL RANGE         |
|---|-----------------------------|---------------------|
| <b>1 Water</b>                            |                             |                     |
| Specific gravity                          | vs water                    | 1.035-1.075         |
| Total solids                              | gm/100 ml                   | 18-22               |
| Water                                     | ml/100 ml                   | 78-82               |
| <b>2 Hydrogen Ion Concentration</b>       |                             |                     |
| Hydrogen ion concentration arterial       | pH                          | 7.35-7.45           |
| Hydrogen ion concentration venous         | pH                          | 7.35-7.38           |
| <b>3 Minerals</b>                         |                             |                     |
| Aluminum                                  | mg/100 ml                   | 0-10 <sup>o</sup>   |
| Bromine                                   | mg/100 ml                   | 0 <sup>o</sup> -0.3 |
| Chloride                                  | mEq/liter                   | 83-86               |
| Iodine                                    | mcg/100 ml                  | 0.6-8.4             |
| Magnesium                                 | mg/100 ml                   | 1-2.0               |
| Manganese                                 | mcg/100 ml                  | 0-1.03              |
| Rubidium                                  | mg/100 ml                   | 0 <sup>o</sup> -0.4 |
| Silicon                                   | mg SiO <sub>2</sub> /100 ml | 0-1                 |
| Zinc                                      | mg/100 ml                   | 0.6-1.4             |
| <b>4 Nitrogen and Protein Derivatives</b> |                             |                     |
| Alanine                                   | mg/100 ml                   | 0-8.50              |
| Ergothioneine                             | mg/100 ml                   | 3-9.17              |
| Glycine                                   | mg/100 ml                   | 1.8-2.3             |
| Nonprotein nitrogen                       | mg/100 ml                   | 30-40               |
| Urea nitrogen                             | mg/100 ml                   | 10-15               |
| Uric acid                                 | mg/100 ml                   | 2-4                 |
| <b>5 Carbohydrates and Derivatives</b>    |                             |                     |
| Acetaldehyde                              | mg/100 ml                   | 0 <sup>o</sup> -0.6 |
| Alpha ketoglutaric acid                   | mg/100 ml                   | 0-1.03              |
| Citric acid                               | mg/100 ml                   | 13-23               |
| Glucose                                   | mg/100 ml                   | 70-110              |
| Glycogen                                  | mg/100 ml                   | 1-16                |
| Lactic acid                               | mg/100 ml                   | 5-20                |
| Oxalic acid                               | mg/100 ml                   | 0.9-1.5             |
| Pyruvic acid                              | mg/100 ml                   | 0-1 <sup>o</sup>    |
| <b>6 Fats and Derivatives</b>             |                             |                     |
| Acetoacetic acid                          | mg/100 ml                   | 0-10 <sup>o</sup>   |
| Acetone                                   | mg/100 ml                   | 0-0.1               |
| Beta hydroxybutyric acid                  | mg/100 ml                   | 0.5-1.0             |
| Cholic acid                               | mg/100 ml                   | 0.2-3.0             |
| <b>7 Vitamins</b>                         |                             |                     |
| Ascorbic acid                             | mg/100 ml                   | 0.5-1.5             |
| Biotin                                    | mcg/100 ml                  | 1-2                 |
| Coccarboxylase                            | mcg/100 ml                  | 5-10                |
| Folic acid                                | mcg/100 ml                  | 0.5-1.3             |
| Niacin                                    | mg/100 ml                   | 0.6-0.8             |
| Pantothenic acid                          | mg/100 ml                   | 1-47                |
| Pyridoxine                                | mcg/100 ml                  | 20-60               |
| Thiamine                                  | mcg/100 ml                  | 3-9                 |
| Tocopherols                               | mg/100 ml                   | 0.5-1.8             |
| Tryptophane                               | mg/100 ml                   | 250-300             |
| Vitamin B                                 | mcg/ml                      | 0.6-1.4             |
| <b>8 Hormones</b>                         |                             |                     |
| Adrenaline                                | micro mcg/ml                | 20 <sup>o</sup> -0  |
| Histamine                                 | mcg/100 ml                  | 1-10                |

TABLE 9—CONT'D

|   |                   |  |           |
|---|-------------------|--|-----------|
| <b>9 Enzymes</b>                            |                   |  |           |
| Amylase                                     | units/ml          |  | 30 90     |
| Arginase                                    | units/ml          |  | 8 16      |
| Carbonic anhydrase                          | Roughton units/ml |  | 1 ° 6     |
| Phosphatase (acid)                          | Bodansky units    |  | 10 ° 0    |
| Zin anhydrase                               | mcg/ml            |  | , 10      |
| <b>10 Pigments</b>                          |                   |  |           |
| CO Hemoglobin                               | gm/100 ml         |  | 0 10 J    |
| Hemoglobin (men)                            | gm/100 ml         |  | 14 18     |
| Hemoglobin (women)                          | gm/100 ml         |  | 1 16      |
| Methemoglobin                               | gm/100 ml         |  | 0 03 0 13 |
| Sulfhemoglobin                              | gm/100 ml         |  | 0 1 0 3   |
| <b>11 Gases</b>                             |                   |  |           |
| Arterial saturation                         | per cent          |  | 93 98     |
| Carbon dioxide (arterial)                   | ml/100 ml STP     |  | 45 55     |
| CO capacity (T <sub>m</sub> )               | mEq/l ter         |  | 19 5 ° 3  |
| Carbon monoxide (arterial)                  | ml/100 ml STP     |  | 0 1 0 5   |
| Nitrogen (arterial)                         | ml/100 ml STP     |  | 1         |
| Oxygen (arterial)                           | ml/100 ml STP     |  | 18 °      |
| pCO   | mm Hg             |  | 40-4      |
| <b>12 Measurements of Clinical Interest</b> |                   |  |           |
| Hematocrit (men)                            | ml/100 ml         |  | 4 5°      |
| Hematocrit (women)                          | ml/100 ml         |  | 9 46      |
| Mean corpuscular hemoglobin concentration   | per cent          |  | 33 38     |
| Mean corpuscular volume                     | cu microns        |  | 90 94     |
| Phenols                                     | mg/100 ml         |  | 0 7 1 8   |
| Idenols (free)                              | mg/100 ml         |  | 0 05 0 09 |
| Sedimentation rate                          | ml/min            |  | 0 1 0 3   |
| Urethane                                    | mg/100 ml         |  | 0 1 1 4   |

TABLE 30 COMPOSITION OF PLASMA (OR SERUM)

(Listed below are normal ranges for various substances in human plasma all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake and to exercise.)

| SUBSTANCE                                 | UNITS              | USUAL RANGE |
|---|--------------------|-------------|
| <b>1 Water</b>                            |                    |             |
| Specific gravity                          | vs water           | 1.015-1.035 |
| Total solids                              | gm/100 ml          | 8-10        |
| Water                                     | ml/100 ml          | 90-9        |
| <b>2 Hydrogen Ion Concentration</b>       |                    |             |
| Hydrogen ion concentration arterial       | pH                 | 7.35-7.45   |
| Hydrogen ion concentration venous         | pH                 | 7.30-7.35   |
| <b>3 Minerals</b>                         |                    |             |
| Aluminum                                  | mg/100 ml          | 0.05-0.1    |
| Arsenic                                   | mcg/100 ml         | 0.50        |
| Bromine                                   | mg/100 ml          | 0.2-0.6     |
| Calcium                                   | mg/100 ml          | 90-115      |
| Chloride                                  | mEq/liter          | 95-110      |
| Cobalt                                    | mcg/100 ml         | 5-10        |
| Copper                                    | mcg/100 ml         | 50-250      |
| Iodine (protein bound)                    | mcg/100 ml         | 4-8         |
| Iron                                      | mcg/100 ml         | 100-150     |
| Lead                                      | mg/100 ml          | 0.01        |
| Magnesium                                 | mg/100 ml          | 1-3         |
| Manganese                                 | mcg/100 ml         | 0.1-0.3     |
| Phosphorus (as inorganic P)               | mg/100 ml          | 3-4         |
| Potassium                                 | mEq/liter          | 4.5-5.5     |
| Sodium                                    | mEq/liter          | 135-145     |
| Sulfate (inorganic)                       | mg/100 ml          | 0.9-1.1     |
| Thiocyanate                               | mg/100 ml          | 0.008       |
| Total base                                | mEq/liter          | 155-165     |
| Total sulphur                             | mg/100 ml          | 90-110      |
| Zinc                                      | mcg/100 ml         | 1.7-900     |
| <b>4 Nitrogen and Protein Derivatives</b> |                    |             |
| Albumin                                   | gm/100 ml          | 4.5-5.5     |
| Alpha amino nitrogen                      | mg/100 ml          | 32-55       |
| Alpha globulin                            | gm/100 ml          | 0.7-1.5     |
| Amino acid nitrogen                       | mg/100 ml          | 5-8         |
| Arginine                                  | mg/100 ml          | 5           |
| Beta globulin                             | gm/100 ml          | 0.6-1.3     |
| Creatine                                  | mg/100 ml          | 3-7         |
| Creatinine                                | mg/100 ml          | 1-2         |
| Cystine                                   | mg/100 ml          | 90-400      |
| Fibrinogen                                | gm/100 ml          | 0.2-0.4     |
| Free amino acid nitrogen                  | mg/100 ml          | 9           |
| Gamma globulin                            | gm/100 ml          | 1.0-1.4     |
| Globulin (total)                          | gm/100 ml          | 1.5-3.5     |
| Glutamic acid                             | mg/100 ml          | 0.4-1.7     |
| Glutamine                                 | mg/100 ml          | 4-11        |
| Glutathione                               | mg/100 ml          | 40-65       |
| Histidine                                 | mg/100 ml          | 1.1-2.1     |
| Isoleucine                                | mg/100 ml          | 1.0-2.3     |
| Leucine                                   | mg/100 ml          | 1.6-2.6     |
| Lysine                                    | mg/100 ml          | 1.3-3.0     |
| Methionine                                | mg/100 ml          | 100-120     |
| Nonprotein nitrogen                       | mg/100 ml          | 20-30       |
| Protein total                             | gm/100 ml          | 6.0-7.0     |
| Phenylalanine                             | mg/100 ml          | 0.01-0.03   |
| Sulphydryl groups                         | mg cysteine/100 ml | 6-7         |
| Threonine                                 | mg/100 ml          | 1.3-3       |
| Tyrosine                                  | mg/100 ml          | 1.0-2.0     |
| Tryptophane                               | mg/100 ml          | 1           |
| Urea nitrogen                             | mg/100 ml          | 10-25       |
| Uric acid                                 | mg/100 ml          | 3-6         |
| Valine                                    | mg/100 ml          | 3-4         |

TABLE 30—CONT'D

## 5 Carbohydrates and Derivatives

|                         |           |         |
|-------------------------|-----------|---------|
| Acetaldehyde            | mg/100 ml | 0.0004  |
| Alpha ketoglutaric acid | mg/100 ml | 0.1-0.3 |
| Citric acid             | mg/100 ml | 1.6 "   |
| Lactic acid             | mg/100 ml | 0.0     |
| Oxalic acid             | mg/100 ml | 1.0-0   |
| Pyruvic acid            | mg/100 ml | 0.0-0   |

## 6 Fats and Derivatives

|                          |             |         |
|--------------------------|-------------|---------|
| Acetoacetic acid         | mg/100 ml   | 0.103   |
| Acetone                  | mg/100 ml   | 0.01    |
| Beta hydroxybutyric acid | mg/100 ml   | 0.0-1.5 |
| Cholesterol esters       | mg/100 ml   | 100-155 |
| Cholesterol (free)       | mg/100 ml   | 3-100   |
| Cholesterol (total)      | mg/100 ml   | 130-330 |
| Esterified fatty acids   | mEq/liter   | 13      |
| Ethanolamine             | mg/100 ml   | 40-50   |
| Lecithin                 | mg/100 ml   | 8-10    |
| Linoleic acid            | mg/100 ml   | 0       |
| Lipid chlorine           | mEq/liter   | 7.57    |
| Phospholipids            | mg P/100 ml | 1-14    |
| Total acetone bodies     | mg/100 ml   | 0.8-5.0 |
| Total lipids             | mg/100 ml   | 500-600 |

## 7 Vitamins

|  |            |           |
|--|------------|-----------|
| Ascorbic acid                            | mg/100 ml  | 0.0-1.5   |
| Beta-carotene                            | mcg/100 ml | 100-50    |
| Biotin                                   | mcg/100 ml | 0.00-1.00 |
| Choline                                  | mg/100 ml  | 0.5-1.5   |
| Desoxyribonuclease                       | units      | 0.14-0.54 |
| Flavine adenine dinucleotide             | mg/100 ml  | 1.5-3.0   |
| Folic acid                               | mg/100 ml  | 1.0-1.9   |
| Free B <sub>12</sub> monophosphate (FBL) | mcg/100 ml | 0.3-1.3   |
| Inositol                                 | mg/100 ml  | 0.4-0.8   |
| Niacin                                   | mg/100 ml  | 0.6-0.78  |
| Pantothenic acid                         | mcg/100 ml | 6         |
| Riboflavin                               | mcg/100 ml | 1.9       |
| Thiamine                                 | mcg/100 ml | 1.0-1.9   |
| Tocopherols                              | mg/100 ml  | 0.0-1.8   |
| Tryptophane                              | mg/100 ml  | 106-1.0   |
| Vitamin A                                | mcg/100 ml | 0.0       |

## 8 Hormones

|            |           |      |
|------------|-----------|------|
| Adrenaline | mg/100 ml | 0.10 |
|------------|-----------|------|

## 9 Enzymes

|                      |                 |           |
|----------------------|-----------------|-----------|
| Amylase              | unit/100 ml     | 15-30     |
| Beta glucuronidase   | units/100 ml    | 0-150     |
| Catalase             | peroxide units  | 4-10      |
| Choline esterase     | 1/2 chel units  | 0.7-1.4   |
| Dialase              | unit/100 ml     | 9-64      |
| Lipase               | ml of N/10 NaOH | 1-1 ml    |
| Peptidase            | units/ml        | 6-8       |
| Phosphatase acid     | Bodansky unit   | 0 or 1.4  |
| Phosphatase alkaline | Bodansky units  | 0 or less |
| Tributyrylase        | Goldstein units | 90-150    |

## 10 Pigments

|                    |           |         |
|--------------------|-----------|---------|
| Bilirubin (direct) | mg/100 ml | 0.04    |
| Bilirubin (total)  | mg/100 ml | 0.4-0.9 |

## 11 Gases

|                |           |      |
|----------------|-----------|------|
| Carbon dioxide | mEq/liter | 0-10 |
|----------------|-----------|------|

## 12 Measurements of Clinical Interest

|                               |           |          |
|-------------------------------|-----------|----------|
| Ikenols (conjugated)          | mg/100 ml | 0.103    |
| Ikenols (free)                | mg/100 ml | 0.04-0.1 |
| Prothrombin (dilution method) | seconds   | 30-45    |
| Prothrombin (Quick's)         | seconds   | 1-1      |

TABLE 31 COMPOSITION OF URINE

(Listed below are normal ranges for various substances in human urine all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake and to exercise.)

| SUBSTANCE                                 | UNITS     | USUAL RANGE   |
|---|-----------|---------------|
| <b>1 Water</b>                            |           |               |
| Specific gravity                          | vs water  | 1.003-1.030   |
| Total solids                              | gm/100 ml | 17            |
| Water                                     | ml/100 ml | 95-99         |
| <b>2 Hydrogen Ion Concentration</b>       |           |               |
| Hydrogen ion concentration                | pH        | 4.5-5.5       |
| <b>3 Minerals</b>                         |           |               |
| Bromine                                   | mg/day    | 15            |
| Calcium                                   | gm/day    | 0.5-0.5       |
| Chloride (as NaCl)                        | gm/day    | 10-15         |
| Chromium                                  | mg/day    | 0.1-0.4       |
| Copper                                    | mg/day    | 0.1-0.5       |
| Fluorine                                  | mg/day    | 0.2-0.5       |
| Iodine                                    | mcg/day   | 35-75         |
| Iron                                      | mg/day    | below 10      |
| Lead                                      | mg/day    | 0.05-0.5      |
| Magnesium                                 | mg/day    | 10-20         |
| Molybdenum                                | mcg/day   | 10-30         |
| Nickel                                    | mcg/day   | 0.3-30        |
| Phosphorus                                | gm/day    | 0.5-2.0       |
| Potassium                                 | mEq/day   | 20-64         |
| Sodium                                    | mEq/day   | 150-197       |
| Solids (total)                            | gm/day    | 20-60         |
| Sulfates (total)                          | gm/day    | 0.6-1.0       |
| Thiocyanate                               | mg/day    | 5-8           |
| Water                                     | l/day     | 0.5-1.0       |
| Zinc                                      | mg/day    | 0.3-0.4       |
| <b>4 Nitrogen and Protein Derivatives</b> |           |               |
| Albumin                                   | mg/day    | less than 100 |
| Alpha amino nitrogen (total)              | mg N/day  | 300-700       |
| Alpha amino nitrogen (bound)              | mg N/day  | 100-200       |
| Alpha amino nitrogen (free)               | mg N/day  | 1.0-4.0       |
| Allantoin                                 | mg/day    | 5-30          |
| Arginine                                  | mg/day    | 50-150        |
| Creatine                                  | mg/day    | less than 100 |
| Creatinine                                | gm/day    | 1.0-1.5       |
| Glutamine                                 | mg N/day  | 10-15         |
| Glycine                                   | mg/day    | 1.0-30        |
| Histidine                                 | mg/day    | 180-300       |
| Isoleucine                                | mg/day    | 1-10          |
| Leucine                                   | mg/day    | 0-30          |
| Lysine                                    | mg/day    | 73-100        |
| Methionine                                | mg/day    | 8-12          |
| Nitrogen (total)                          | gm/day    | 10-18         |
| Threonine                                 | mg/day    | 50-65         |
| Tryptophane                               | mg/day    | 0.0-0.00      |
| Urea                                      | gm/day    | 10-40         |
| Uric acid                                 | gm/day    | 0.5-2.0       |
| Valine                                    | mg/day    | 18-25         |
| Xanthine                                  | mg/day    | 0-40          |
| <b>5 Carbohydrates and Derivatives</b>    |           |               |
| Alpha ketoglutaric acid                   | mg/day    | 1-44          |
| Citric acid                               | mg/day    | 0.1-0.4       |
| Fructose                                  | mg/day    | 0             |
| Pyruvic acid                              | mg/day    | 10-20         |
| Total reducing sugar                      | mg/day    | 0-100         |

TABLE 31—CONT'D

|   |                       |           |
|---|-----------------------|-----------|
| <b>6 Fats and Derivatives</b>               |                       |           |
| Acetone bodies (total)                      | mg/day                | 10 100    |
| Cholesterol                                 | mg/day                | 0 3 10    |
| <b>7 Vitamins</b>                           |                       |           |
| Ascorbic acid                               | mg/day                | 6 18      |
| B <sub>6</sub> otin                         | mcg/day               | 10 2      |
| Choline                                     | mg/day                | 2 4       |
| Niacin                                      | mg/day                | 0 2 10    |
| N-Methylnicotinamide                        | mg/day                | 4 1       |
| N-Methyl - pyridine 5 carboxylamide         | mg/day                | 5 10      |
| Pantothenic acid                            | mg/day                | 10 30     |
| Pyridoxine                                  | mg/day                | 0 0 0 4   |
| Riboflavin                                  | mg/day                | 0 4 1 50  |
| Thiamine                                    | mg/day                | 0 1 0 4   |
| Vitamin A                                   | IU/day                | 0         |
| <b>8 Hormones</b>                           |                       |           |
| Androsterone                                | mg/day                | 1 2       |
| Androgen                                    | mg/day                | 3 10      |
| Corticosterone                              | mg/day                | 0 1 0 3   |
| Cortin                                      | mg/day                | 0 2 0 5   |
| 11 Desoxycorticosteroid                     | mg/day                | 0 1 0 5   |
| Estrogens (female)                          | mcg/day               | 50 150    |
| Estrogens (male)                            | mg/day                | 0 6       |
| Glycogenic steroids                         | mouse units/day       | 40 80     |
| 17 Ketosteroids (female)                    | mg/day                | 5 14      |
| 17 Ketosteroids (male)                      | mg/day                | 8 20      |
| <b>9 Enzymes</b>                            |                       |           |
| Diastase                                    | unit/day              | 8 3       |
| Phosphatase (acid)                          | h. n. Armstrong units | 80 300    |
| Trypsin                                     | unit                  | 40 50     |
| Uropepsin                                   | unit                  | 1000 3000 |
| <b>10 Pigments</b>                          |                       |           |
| Coproporphyrin (Type I)                     | mg/day                | 15 90     |
| Coproporphyrin (Type III)                   | mg/day                | 1 24      |
| Indican                                     | mg/day                | 10 150    |
| Urobilinogen                                | mg/day                | 0 4       |
| Uroporphyrin                                | mg/day                | 0         |
| <b>11 Blood Gases</b>                       |                       |           |
| Carbon dioxide                              | mEq/day               | 0 50      |
| <b>12 Measurements of Clinical Interest</b> |                       |           |
| Phenols (conjugated)                        | mg/day                | 15 40     |
| Phenols (free)                              | mg/day                | 0 0 4     |
| Phenols (total)                             | mg/day                | 15 40     |

TABLE 3. COMPOSITION OF FECES

(Listed below are normal ranges for various substances in human feces all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake. In addition there are large differences from person to person.)

| SUBSTANCE                                 | UNITS                       | USUAL RANGE                  |
|---|-----------------------------|------------------------------|
| 1 <i>Water</i>                            |                             |                              |
| Specific gravity                          | vs water                    | 1.030-1.100                  |
| Total solids                              | gm/100 gm                   | 15-30                        |
| Water                                     | ml/100 gm                   | 63-85                        |
| 2 <i>Hydrogen Ion Concentration</i>       |                             |                              |
| Hydrogen ion concentration                | pH                          | 7.0-7.5                      |
| 3 <i>Minerals</i>                         |                             |                              |
| Aluminum                                  | mg/day                      | 1.0-9                        |
| Calcium                                   | gm/day                      | 0.1-1.0                      |
| Chloride                                  | mEq/day                     | traces except<br>in diarrhea |
| Copper                                    | mg/day                      | 15-211                       |
| Iron                                      | gm/day                      | 0.7-1.0                      |
| Lead                                      | mg/day                      | 0.3-0.4                      |
| Manganese                                 | mg/day                      | 19-24                        |
| Molybdenum                                | mg/day                      | 2-4                          |
| Nickel                                    | mg/day                      | 5-10                         |
| Phosphorus                                | gm P/day                    | 0.9-1.7                      |
| Potassium                                 | mEq/day                     | 19-226                       |
| Sodium                                    | mEq/day                     | traces except<br>in diarrhea |
| Tin                                       | mg/day                      | 0.5-1.7                      |
| Zinc                                      | mg/day                      | 5-10                         |
| 4 <i>Nitrogen and Protein Derivatives</i> |                             |                              |
| Arginine                                  | gm/day                      | 1-2.1                        |
| Histidine                                 | gm/day                      | 0.6-0.8                      |
| Isoleucine                                | gm/day                      | 14-23                        |
| Leucine                                   | gm/day                      | 18-29                        |
| Lysine                                    | gm/day                      | 19-24                        |
| Methionine                                | gm/day                      | 0.5-0.8                      |
| Nitrogen (total)                          | gm/day                      | 0.7-2.1                      |
| Threonine                                 | gm/day                      | 14-20                        |
| Valine                                    | gm/day                      | 1.5-2.6                      |
| 5 <i>Carbohydrates and Derivatives</i>    |                             |                              |
| Total reducing sugar                      | mg/day                      | 0                            |
| 6 <i>Fats and Derivatives</i>             |                             |                              |
| Total fat                                 | gm/day                      | 1-7                          |
| Total fat                                 | per cent by weight<br>(dry) | 10-25                        |
| Total fat (unsaponifiable)                | per cent by weight<br>(dry) | 0-5                          |
| 7 <i>Vitamins</i>                         |                             |                              |
| Beta carotene                             | mg/day                      | 1.7-3.3                      |
| Biotin                                    | mcg/day                     | 100-200                      |
| Niacin                                    | mg/day                      | 35-55                        |
| 1-antiothenic acid                        | mg/day                      | 18-38                        |
| Pyridoxine                                | mg/day                      | 0.1-0.5                      |
| 4-Pyridoxic acid                          | mg/day                      | 0.5-0.6                      |
| Proflavin                                 | mg/day                      | 0.4-1.20                     |
| Thiamine                                  | mg/day                      | 0-0.8                        |
| Vitamin A                                 | mg/day                      | 0.17-0.33                    |
| 8 <i>Hormones</i>                         |                             |                              |
| No data available                         |                             |                              |

TABLE 3 --CONT'D

|  |                             |  |         |
|--|-----------------------------|--|---------|
| 9 <i>Enzymes</i>                           |                             |  |         |
| No data available                          |                             |  |         |
| 10 <i>Pigments</i>                         |                             |  |         |
| Porphyrin total                            | mcg/day                     |  | 300 400 |
| Protoporphyrin                             | mcg/day                     |  | 0 300   |
| Urobilinogen                               | mg/day                      |  | 40 80   |
| 11 <i>Gases</i>                            |                             |  |         |
| No data available                          |                             |  |         |
| 1 <i>Measurements of Clinical Interest</i> |                             |  |         |
| Bacterial debris                           | per cent of weight<br>(dry) |  | 10 30   |



TABLE 3 COMPOSITION OF FECES

(Listed below are normal ranges for various substances in human feces all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake. In addition there are large differences from person to person.)

| SUBSTANCE                                 | UNITS                       | USUAL RANGE                  |
|---|-----------------------------|------------------------------|
| 1 <i>Water</i>                            |                             |                              |
| Specific gravity                          | vs water                    | 1.030-1.100                  |
| Total solids                              | gm/100 gm                   | 10-30                        |
| Water                                     | ml/100 gm                   | 60-80                        |
| 2 <i>Hydrogen Ion Concentration</i>       |                             |                              |
| Hydrogen ion concentration                | pH                          | 4.0-7.5                      |
| 3 <i>Minerals</i>                         |                             |                              |
| Aluminum                                  | mg/day                      | 10-9                         |
| Calcium                                   | gm/day                      | 0.1-1.0                      |
| Chloride                                  | mEq/day                     | traces except<br>in diarrhea |
| Copper                                    | mg/day                      | 15-11                        |
| Iron                                      | gm/day                      | 0.7-1.0                      |
| Lead                                      | mg/day                      | 0.3-0.4                      |
| Manganese                                 | mg/day                      | 10-4                         |
| Molybdenum                                | mg/day                      | 2-4                          |
| Nickel                                    | mg/day                      | 5-10                         |
| Phosphorus                                | gm P/day                    | 0.9-1.7                      |
| Potassium                                 | mEq/day                     | 19-226                       |
| Sodium                                    | mEq/day                     | traces except<br>in diarrhea |
| Tin                                       | mg/day                      | 0-1.7                        |
| Zinc                                      | mg/day                      | 5-10                         |
| 4 <i>Nitrogen and Protein Derivatives</i> |                             |                              |
| Arginine                                  | gm/day                      | 1-2.1                        |
| Histidine                                 | gm/day                      | 0.6-0.8                      |
| Isoleucine                                | gm/day                      | 14-23                        |
| Leucine                                   | gm/day                      | 18-29                        |
| Lysine                                    | gm/day                      | 19-29                        |
| Methionine                                | gm/day                      | 0.5-0.9                      |
| Nitrogen (total)                          | gm/day                      | 0.7-1                        |
| Threonine                                 | gm/day                      | 14-22                        |
| Valine                                    | gm/day                      | 1.5-2.6                      |
| 5 <i>Carbohydrates and Derivatives</i>    |                             |                              |
| Total reducing sugar                      | mg/day                      | 0                            |
| 6 <i>Fats and Derivatives</i>             |                             |                              |
| Total fat                                 | gm/day                      | 1-7                          |
| Total fat                                 | per cent by weight<br>(dry) | 10-5                         |
| Total fat (unsaponifiable)                | per cent by weight<br>(dry) | 0.5                          |
| 7 <i>Vitamins</i>                         |                             |                              |
| Beta carotene                             | mg/day                      | 1-33                         |
| Biotin                                    | mcg/day                     | 100-200                      |
| Niacin                                    | mg/day                      | 3.5-5.5                      |
| Pantothenic acid                          | mg/day                      | 18-38                        |
| Pyridoxine                                | mg/day                      | 0.1-0.5                      |
| 4-Pyridoxic acid                          | mg/day                      | 0.5-0.6                      |
| Riboflavin                                | mg/day                      | 0.4-1.20                     |
| Thiamine                                  | mg/day                      | 0-0.8                        |
| Vitamin A                                 | mg/day                      | 0.17-0.33                    |
| 8 <i>Hormones</i>                         |                             |                              |
| No data available                         |                             |                              |

TABLE 3 —CONT'D

|   |                             |  |         |
|---|-----------------------------|--|---------|
| 9 <i>Enzymes</i>                            |                             |  |         |
| No data available                           |                             |  |         |
| 10 <i>Pigments</i>                          |                             |  |         |
| Porphyrin total                             | mcg/day                     |  | 300-400 |
| Protoporphyrin                              | mcg/day                     |  | 0-300   |
| Urobilinogen                                | mg/day                      |  | 40-80   |
| 11 <i>Gases</i>                             |                             |  |         |
| No data available                           |                             |  |         |
| 12 <i>Measurements of Clinical Interest</i> |                             |  |         |
| Bacterial debris                            | per cent of weight<br>(dry) |  | 10-30   |

TABLE 32 COMPOSITION OF FECES

(Listed below are normal ranges for various substances in human feces all values are for healthy individuals in a resting condition at sea level in a temperate environment. The values for many substances are closely related to dietary intake. In addition there are large differences from person to person.)

| SUBSTANCE                                 | UNITS                       | USUAL RANGE                  |
|---|-----------------------------|------------------------------|
| 1 <i>Water</i>                            |                             |                              |
| Specific gravity                          | vs water                    | 1.030-1.100                  |
| Total solids                              | gm/100 gm                   | 10-30                        |
| Water                                     | ml/100 gm                   | 60-80                        |
| <i>Hydrogen Ion Concentration</i>         |                             |                              |
| Hydrogen ion concentration                | pH                          | 7.0-8.5                      |
| 3 <i>Minerals</i>                         |                             |                              |
| Aluminum                                  | mg/day                      | 10-29                        |
| Calcium                                   | gm/day                      | 0.1-1.0                      |
| Chloride                                  | mEq/day                     | traces except<br>in diarrhea |
| Copper                                    | mg/day                      | 15-211                       |
| Iron                                      | gm/day                      | 0.7-1.0                      |
| Lead                                      | mg/day                      | 0.3-0.4                      |
| Manganese                                 | mg/day                      | 19-24                        |
| Molybdenum                                | mg/day                      | 0.4                          |
| Nickel                                    | mg/day                      | 5-10                         |
| Phosphorus                                | gm P/day                    | 0.9-1.7                      |
| Potassium                                 | mEq/day                     | 19-200                       |
| Sodium                                    | mEq/day                     | traces except<br>in diarrhea |
| Tin                                       | mg/day                      | 0.5-1.7                      |
| Zinc                                      | mg/day                      | 5-10                         |
| 4 <i>Nitrogen and Protein Derivatives</i> |                             |                              |
| Arginine                                  | gm/day                      | 1-2.1                        |
| Histidine                                 | gm/day                      | 0.6-0.8                      |
| Isoleucine                                | gm/day                      | 14-23                        |
| Leucine                                   | gm/day                      | 18-29                        |
| Lysine                                    | gm/day                      | 19-29                        |
| Methionine                                | gm/day                      | 0.5-0.8                      |
| Nitrogen (total)                          | gm/day                      | 0.7-1.1                      |
| Threonine                                 | gm/day                      | 14-22                        |
| Valine                                    | gm/day                      | 15-26                        |
| 5 <i>Carbohydrates and Derivatives</i>    |                             |                              |
| Total reducing sugar                      | mg/day                      | 0                            |
| 6 <i>Fats and Derivatives</i>             |                             |                              |
| Total fat                                 | gm/day                      | 1-2                          |
| Total fat                                 | per cent by weight<br>(dry) | 10-25                        |
| Total fat (unsaponifiable)                | per cent by weight<br>(dry) | 0.5                          |
| 7 <i>Vitamins</i>                         |                             |                              |
| Beta carotene                             | mg/day                      | 1.7-3.3                      |
| Biotin                                    | mcg/day                     | 100-200                      |
| Niacin                                    | mg/day                      | 35-55                        |
| Pantothenic acid                          | mg/day                      | 18-38                        |
| Pyridoxine                                | mg/day                      | 0.1-0.5                      |
| 4-Pyridoxic acid                          | mcg/day                     | 0.5-0.6                      |
| Riboflavin                                | mg/day                      | 0.4-1.0                      |
| Thiamine                                  | mg/day                      | 0.2-0.8                      |
| Vitamin A                                 | mcg/day                     | 0.1-0.33                     |
| 8 <i>Hormones</i>                         |                             |                              |
| No data available                         |                             |                              |

TABLE 3<sup>a</sup>—CONT D

|       |  |                             |         |
|-------|--|-----------------------------|---------|
| <hr/> |  |                             |         |
| 9     | <i>Enzymes</i>                           |                             |         |
|       | No data available                        |                             |         |
| 10    | <i>Pigments</i>                          |                             |         |
|       | Porphyrin, total                         | meg/day                     | 300 400 |
|       | Protoporphyrin                           | meg/day                     | 0 300   |
|       | Urobilinogen                             | mg/day                      | 40 80   |
| 11    | <i>Gases</i>                             |                             |         |
|       | No data available                        |                             |         |
| 1     | <i>Measurements of Clinical Interest</i> |                             |         |
|       | Bacterial debris                         | per cent of weight<br>(dry) | 10 30   |
| <hr/> |  |                             |         |



## A

- Al orb at  
 Decalco 185  
 Supe orb 193  
 Zeolite 140  
 Al orpt on columns (Hennes ey type) 184  
 Al or plate a say for botin of 39  
 Angeler P M 418  
 Agr n G 148 149  
 A r lamped balance (Sartorius typ) 6 68  
 Alkanol solution 40  
 Alkanese A A 118 120 153  
 Allet A 19  
 Al riv P A 80  
 Al right F 19 0 39 400  
 Al u n n serum determination of 1 1  
 Al h i  
 apryl 6 81  
 l te mination n  
 blood 163 164  
 ur ne 163 164  
 amyl 161 17  
 butyl 185  
 normal butyl 08  
 Al ch l m Hancous references 4 1  
 Al ohol eth r mixture 1 4  
 Al ohol c  
 glycerol 111  
 pho photogati a it solution 137  
 p ta um hi brovide solat on 18  
 Aller lerg D 14  
 Al van le B 101  
 Alkaline earth 113  
 Alkaline phosphatase in serum determina  
 tion of ~ 3  
 Alkaline p ale solution preparation 139  
 Alkaline substrate for pho phata c prepara  
 tion 4  
 Allen F H J 1 3 9  
 Allerv V 1  
 Al la n n ulotide bound 19  
 Al ip A 150  
 Alpha alpha pyridine react 99  
 Al n F 116  
 Altitud equval nt 384  
 Altimet h apparatus f r embolling tis  
 sue 391  
 Alum num 115  
 Alva z W C 4 0  
 Alve la air ne r 316-350  
 apparatus 316-34  
 relation 319 3 n  
 principle 346  
 Alveolar air in work 3 0  
 apparatus 319 3 0  
 f incip 350  
 Al a ur d nicator 94  
 Al m nt p e ure 61  
 Al m nt l mperature 461  
 An 84 1 14 1 0

- Amino acetic acid standard solution, 14°  
 amino acids microbiological 259  
 Amino acids miscellaneous references 259  
 Amino acid nitrogen in urine determination of 14°  
 Amino naphtholsulfonic acid (ANSA) solution, 174  
 amino nitrogen manometric methods 299  
 Amino nitrogen miscellaneous reference 299  
 Ammonia free urine preparation of 141  
 Ammoniacal alcohol solution 97  
 Ammonium hydroxide solution 93  
 Ammonium molybdate solution preparation 111 14  
 Ammonium oxalate solution preparation 93  
 Ammonium potassium oxalate mixture preparation 131  
 Ammonium sulfate standard solution preparation 130 135 140  
 Amylase  
   in serum, determination of 226 27,  
   simplified calculation (Table) 227  
 Amylase miscellaneous references 27  
 Amylase activity of serum 227  
 Amylase units in serum 227  
 Amyloidosis miscellaneous references 418  
 Anaerobic 235  
   handling of blood 263  
   storage of blood (Scholander technique) 264  
 Andersch M A 223  
 Anderson C T 117  
 Anderson M X 364  
 Anderson W N 178  
 Androgens miscellaneous references 18  
 Androsterone standard solution preparation 217  
 Anemia 131  
 Angstrom units 29  
 Annersten S 218  
 ANSA (amino naphtholsulfonic acid)  
   recrystallization 109  
   solution preparation 174  
 Anson C R 259  
 Anson M L 147  
 Antecubital fossa 261  
 Anticoagulants for blood miscellaneous references 401  
 Anticoagulants for blood preparation  
   heparin 20 301  
   oxalate 131  
   sodium citrate solution 130  
 Antimony 115  
 Antimony trichloride solution preparation 179  
 Anti protease miscellaneous references 227  
 Anti tryptic activity miscellaneous reference 227  
 Apparatus  
   alveolar air at rest (Figure) 341  
   alveolar air at work (Figure) 349  
   balance 80  
   capillary syringe and sling for the combined CO and O method 296  
   electron microscope 80  
   electrophoresis 80  
   extraction unit for 17 ketosteroids 216  
   flame photometer 80  
   Apparatus—Cont d  
   gas analysis 80  
   gas mixture preparation 30  
   glass electrode (Daly Forbes) 239  
   Haldane gas analysis (Figure) 313  
   infrared spectrometer 80  
   Kjeldahl micro (Hayes) 121  
   Krogh syringe burette 13  
   lactic acid 158  
   manometric techniques 80  
   mass spectrometer 80  
   measurement of functional residues 350  
   microrespirometer 80  
   oximeter 80  
   pipettes 68 81  
   polarography 81  
   Pechberg microburette 12  
   respirometer 81  
   Scholander micrometer burette 14  
   Scholander micrometer gas analyzer  
   Tissot gasometer for collecting expired gas 328  
   total base electrolysis unit 84  
   treadmill 340  
   ultracentrifuge 81  
   ultrafiltration 81  
   ultra microanalysis 81  
   Van SlykeNeill 264  
 Apparatus volumetric 63 74  
   burettes 68 74  
   calibration 74 6  
   microburettes 70 4  
   syringe pipettes 63  
   syringe screw burettes 74  
 Aqua regia 115  
 Arcand G M 116  
 Archibald R M 127 130  
 Arginase miscellaneous reference 7  
 Arginine microbiological miscellaneous references 29  
 Armstrong A R 223  
 Arnold sterilizer 244  
 Arsenic pentoxide preparation 143  
 Arsenophosphotungstic acid solution preparation 143 145  
 Arterial blood technique for obtaining 65  
   Goldschmidt and Licht method 263  
   Lahenthal and Pilev method 263  
   Mook's technique 263  
   Scholander technique 63 6,  
   usual method 61 6  
 Arterio venous difference 354  
 Ascorbic acid field method determination serum 201 204  
   urine 201 204  
 Ascorbic acid macro method determination  
   serum 190 197  
   tissue 192 200  
   urine 198  
   whole blood 193 194  
 Ascorbic acid micro method determination  
   serum 198 199  
   whole blood 198 199  
 Ascorbic acid miscellaneous references

- Ascorbic acid standard solution for field use preparation °01  
 Ascorbic acid standard solution preparation 196  
 Ashby count miscellaneous reference 401  
 Ashworth J N 15  
 Asparagine solution preparation 245  
 Astaff A. 418  
 Atkin, L. 251  
 Atomic number of the elements (table) 4 8  
 Atomic weight definition 76  
 table 4 4  
 Austin J H 311  
 Autoclave °° 139  
 Automatic pipette  
 Brewer °57  
 Fisher Volumat °44  
 Autotransformer mercury vapor lamps 55  
 Average normal values of  
 hematocrit 133  
 hemoglobin 133  
 Avery B F 300  
 Awapara J. 148 149  
 Axelrod A E 401
- B
- Baier H N, 358  
 Balaubridge H W 116  
 Balances  
 air damped (Bartorius type) 47 68  
 calibration 67 68  
 Roller Smith torsion 5 63 67  
 calibration 65  
 operation 65  
 Baldwin E. deF 35 357  
 Bale W F, 117  
 Ballistocardiograph miscellaneous references 356  
 Bandelin F J 59  
 Baner B 178  
 Banerjee, S 213  
 Baptist J H 148  
 Barach A L 357  
 Barah J, 150  
 Barbour H G 15.  
 Barclay J A 15  
 Barcroft gas sampling tubes 3 7  
 Barcroft tonometer, 501  
 Barium chloride solution 113  
 Barium hydroxide 156  
 Barker M. H 364  
 Barker S D 117  
 Barnes H D 33  
 Barnett J 19  
 Barometric pressure conversion chart 334  
 Barrett, E. 15  
 Bartczak, E. 2 0  
 Basal medium for folio acid preparation 246  
 Basal metabolic rate (BMR) 337 341  
 calculation, 31 341  
 principle 337  
 procedure 337  
 Base bound by plasma proteins effects of temperature 311  
 Bave E. 115  
 Bascett, A. M. 219  
 Bastrup J T 178  
 Batchelor T M. 117  
 Bateman J B 357  
 Bates R. W 219  
 Bates Smith E C 150  
 Bauer F C Jr 178  
 Bauernfeind J C 213  
 Baumann C 4 31  
 Baxter D 130 13.  
 Beadle G W °9  
 Bean R. S 2 8  
 Bean W B 418  
 Beazell J M 4°9  
 Beck G J 418  
 Beckman flame photometer attachment 60  
 Beckman pH meter °42  
 Beckman Quartz spectrophotometer Model DU 38 50  
 absorption cells holders and compartments 38 42  
 accuracy 4.  
 calculation of solution concentrations 47  
 calibration 43 50  
 electrical system 38  
 focusing of hydrogen lamp 48 49  
 focusing of tungsten lamp 47 43  
 light sources 39  
 monochromator 39  
 nominal band width 46  
 operation of instrument 50  
 optical system 38 39  
 phototubes 42  
 sensitivity 42  
 slits 39  
 spectral range 4  
 Beckman W W 15  
 Behre J A. 178  
 Behrwaltes W H 117  
 Belding H S 3 8 378  
 Bell J C 219  
 Bell R. D 140  
 Benzoe B 213  
 Benedict B R 143 166 4°1  
 Benedict Hitchcock stock standard solution for uric acid preparation 143  
 Benedict's qualitative glucose reagent 403  
 Benedict's qualitative test for urine sugar 404  
 Benne E J 233  
 Bennett H B 91 118  
 Bentley H R. 150  
 Bentley O G 260  
 Benzene 145  
 Benzene acid standard solution preparation 78 85  
 Berg C P, 153  
 Bergdahl A M. °60  
 Bergman O °9  
 Berhner R. W., 411  
 Berman L. 401  
 Bernheim A. R., 168  
 Berry H K. 148  
 Berryman, G H 375  
 Bernsen S 3.6  
 Berthles A K C 300  
 Beatty O A., 198 °1 363  
 Bhattacharya, G 213



- Bicarbonate manometric method miscel-  
laneous references 299
- Bicknell F, 212
- Biehler R 213
- Bile in urine miscellaneous references 418
- Bilirubin in serum combining with diazon-  
ium salts 413  
determination of 413 414  
principle 413  
standard solution preparation 413
- Bills, C E 118
- Bing R J 356
- Binkley, F 116 152
- Biochemical procedures 82 60
- Biotin agar plate assay 256 69  
medium 258  
principle 257  
standard curve 258
- Biotin microbiological method 249 250  
determination in  
blood 249 250  
feces 249 250  
food 249 250  
plasma 249 250  
urine 249 250  
standard solution preparation 245
- Bisulfite combining capacity in blood 163
- Bird E J 117
- Bird O D 213
- Biuret reagent preparation 1 4
- Black A 255
- Black, H 357
- Blanchard K. C 421
- Blackfan K D 1 9
- Blatz J I 412 417
- Blanch M B 145
- Blendor Waring 19
- Blux G 80
- Block H 209
- Block P J 148 149
- Block W D 145 153
- Blodner H 163
- Blood anti-coagulant  
heparin solution 20 301  
potassium oxalate solution 164  
sodium citrate solution 135
- Blood arterial  
calculation  
hemoglobin 306  
percent saturation 306  
physically dissolved oxygen 305  
oxygen bound by hemoglobin 305  
oxygen capacity 306  
total oxygen 305  
dissociation curves 305 310  
equilibration 304  
normal values 311  
preparation for analysis 303
- Blood, arterial and venous  
analysis of gas phase 304 305  
calculation 305  
equilibration of in gas mixture of known  
composition 301 305  
gas mixture preparation 301 303  
principle, 301
- Blood determination of  
acetone bodies total 175 1 6  
alcohol 163 165  
ascorbic acid  
macro method 195 197  
micro method 198 199  
biotin 249 250  
carbon dioxide (S R) 291 295  
carotenoids (Carr Price) 179 180  
combined carbon dioxide and oxygen (Scho-  
lander) 295 299  
eosinophil count (direct) 399 401  
folic acid 244 247  
glucose  
colorimetric 154 155  
volumetric 156 157  
hemoglobin (Evelyn Malloy) 230 231  
hemoglobin (Sanford Sheard) 229 230  
lactic acid 157 161  
lipid phosphorus 173 174  
mean corpuscular  
concentration 397 398  
hemoglobin content 397 398  
volume 397 398  
methemoglobin 231 232  
nicotinic acid microbiological 48  
nitrogen dissolved 282 285  
nonprotein nitrogen 135 137  
oxyhemoglobin 231 232  
oxygen (Roughton Scholander) 277 288  
oxygen (Van Slyke) 267  
pantothenic acid 251  
pH measurement 235 242  
riboflavin 19 195  
sedimentation rate 399 400  
sulfhemoglobin 231 232  
thiamine 186 189  
total bisulfite binding capacity 165  
total circulating volume  
calculation 36 363  
carbon monoxide method 361 363  
tryptophan 255 256  
urea nitrogen 140 141  
vitamin A (Carr Price) 179 180
- Blood drawing in the field 21
- Blood gases physical chemistry 301 311
- Blood total carbon dioxide calculation  
carbon dioxide content of equilibrated  
blood 306  
glycolysis during equilibration 306  
total carbon dioxide in plasma 306
- Blood volume dye method 358 361  
calculation 359 361
- Blood volume miscellaneous references 368
- Blood volume True 363 364
- Bloom A, 219
- Bloor W R 168 169 178 218
- Blum A E 259 260
- Bock A V 327 346
- Bockus H L 221
- ten Bokkel Huinink H 178
- Bone marrow miscellaneous references 401
- Boothby W M 357
- Boothby Sandiford standards for BMR 340
- Borate buffer pH 9 245
- Borek E, 260
- Boric acid solution preparation 1 0
- Botty M. C 150

Bowen A deP 376  
 Bowes and Church book on food compos  
 tion 375  
 Bowman J M 25  
 Boyer F D 166  
 Boyle A J 117  
 Brach urinary blood drawing 81  
 Brand F C 171  
 Brannon E S 306  
 Bratton A C 41  
 Bray W F 396  
 Brecher A 250  
 Breed E S 30  
 Breh F 118  
 Brendler H 4  
 Bridwell E G 41  
 Briggs A P 97 11  
 Brimley R C 147  
 Brink N G 260  
 Britton, C J C 206  
 Brock M J 368  
 Bromeresol green and cator 10 185 191  
 09 235  
 Bromeresol purple indicator 235  
 Bromine 190  
 Bromine in seallaneous references 116  
 Bromine water preparat on 143  
 Bromphenol blue and cator 186 230  
 Bromphenol red and cator 235  
 Bromsulfale n test in serum determinat on  
 of 41 413  
 solut on preparat on 412  
 Bromthymol blue and cator 23 44  
 Brown L 369 373  
 Brower S W 417  
 Brown D M 80  
 Brown F 150  
 Brown F A 213  
 Brown S 167  
 Buchanan O H 140  
 Buffer  
 borate pH 2 145  
 borate (Clark's) pH 7 8 10 0 237  
 definit on 35  
 pH 4 58 89  
 pH 47 100  
 pH definit on 35  
 phosphate pH 65 30  
 phosphate pH 66 230 231  
 phosphate Strensen's pH 5 8 8 236  
 pH definit on 235  
 potassium citrate 2  
 sodum diethylbarbiturate stock solution  
 24  
 Strensen's glycine HCl mixture pH 11 35  
 38  
 Strensen's glycine NaOH mixture pH 8 5  
 1 5 38  
 Strensen's phosphate solut on at 20 C and  
 37 C pH 5 8 8 40  
 Walpole's acetate mixtures pH 3 6 5 8 237  
 Behle E L 41  
 Bulger H A 153  
 Bull H B 14 143  
 Bullen S 9 18  
 Bureau of standard filters, for Beckman in  
 strument 50

## Burettes

Haldane cal brat on 31 310  
 Krogh syringe type 3  
 micro Fehberg 70 4 91  
 micro Scholander 74 6  
 cal brat on 6  
 construct on 2  
 syringe a re Krogh type 74  
 Barnett I 407  
 Burrus R T 80  
 Burrus O K 4 2  
 Burton A C 3 2 390  
 Burton C J 150  
 Burton R B 151  
 Burton's equation 290  
 Burwell C S 3 7  
 Buskirk H 11 60  
 Butler A M 8 118 134 21 2 0

## C

Ca L 149  
 Calcum  
 calculat on of theoret al balance 393  
 determinat on 2  
 feces 96 3  
 food 96 9  
 serum 93 90  
 tissue 96 9  
 urine 93 96  
 standard solut on preparat on 94  
 Calcum in seallaneous references 116  
 Calcum hydroxide  
 powdered 87  
 solut on preparat on 175  
 Calcum pantothenate standard solut on 2  
 Calulat on of  
 carbonate n serum 306  
 RMI 337 341  
 carbon dioxide content of blood equilib  
 rated 306  
 Clo unit 379 380  
 free carbon dioxide of serum 306  
 gas volumes 261  
 general equation 61  
 glycolysis during equilibrat on 306  
 hemoglobin in blood 306  
 oxygen bound by hemoglobin 30  
 oxygen capacity 306  
 percent saturat on a blood 306  
 phys ally dissolved oxygen in blood 300  
 serum pH 306  
 temperature effects on base bound by  
 plasma protein 311  
 T<sub>m</sub> from PCO<sub>2</sub> (total CO<sub>2</sub>) and pH 306  
 total carbon dioxide in blood 306  
 total carbon dioxide in plasma 306  
 total oxygen in blood 303  
 true blood volume 363 364  
 urine volume per hour 3  
 Calibrat on of  
 balances  
 a dampel type d 63  
 tors on 60  
 Beckman spectrophotometer 43 09  
 burette (Haldane) 31 315

- Calibration of—Cont'd  
 Coleman, Jr, spectrophotometer 33 38  
 concentration transmittance graph 37 33  
 volumetric apparatus 74 76  
 wave length 29
- Callow N H 218  
 Callow, R K 218  
 Calomel half cell 237  
 Calorie definition 336  
 Caloric discrepancy definition 390  
 Calvery E 123  
 Camien M N 259  
 Capillary syringe and sling for the combined carbon dioxide and oxygen method (Scholander) 296
- Capps B F 60  
 Caprylic alcohol, 26 181 193 268 270  
 Carbohydrates determination of  
 alcohol 163 165  
 bisulfite binding capacity 165 166  
 blood glucose or blood sugar  
 colorimetric 154 155  
 volumetric 156 157  
 glycogen 161 163  
 lactic acid 157 161  
 Carbohydrates miscellaneous references 166 167  
 aldehydes 166  
 citric acid 166  
 glucose 166  
 glycogen 166  
 insulin 166  
 lactic acid 166  
 mannitol 167  
 methylpentose 167  
 pentose 167  
 pyruvic acid 167
- Carbon dioxide  
 absorption with sodium hydroxide 270 272  
 blood determination of 291 295  
 apparatus 292  
 calculation 294 295  
 principle 291  
 procedure 293 294  
 liberation of by lactic acid 270  
 manometric method, 299  
 production calculation of 336  
 serum determination of (Van Slyke) 270 272  
 solid 25
- Carbon dioxide and oxygen combined method (Scholander) 295 299  
 apparatus 296  
 calculation 299  
 principle 295 296  
 procedure 297 298
- Carbon monoxide  
 absorbent Winkler's solution 286  
 blood determination of 275 276  
 blood determination (Scholander Pough ton) 285 291  
 calculation 287  
 general method 285 288  
 procedure 286 287  
 combined method for estimation of in single samples of blood 288 289  
 calculation 288 289  
 blood containing small amounts 289 291
- Carbon monoxide—Cont'd  
 liberation with hydrosulfite borate solution, 275  
 manometric method 299  
 preparation of gas, 362  
 small amounts in blood determination of 289 291
- Carbon monoxide miscellaneous references 299
- Carbon tetrachloride 215
- Carboxyhemoglobin miscellaneous references 233
- Cardiac output (Grollman's method) 353 356  
 analysis of acetylene 355  
 apparatus 354  
 calculation, 355 356  
 principle, 353 354  
 procedure 354 355  
 rebreathing of acetylene 354
- Cardiac output miscellaneous references 356
- Cardini C E 212
- Cardiotachometer, 342 343
- Carmine for collection of feces 26
- Carotene miscellaneous references 212
- Carotene standard solution preparation 182
- $\beta$  carotene standard solution 179
- Carotenoids  
 carotene, 180  
 serum (Dichlorhydrin), 181 184  
 xanthophyll 180
- Carpenter A M, 153  
 Carpenter D C 153  
 Carpenter F H 256
- Carr Price reaction 179 180
- Carter, P 20
- Cartwright E 104
- Casein, acid hydrolyzed 245
- Casein hydrolysate solution 255
- Caswell M C 260
- Catalyst for micro Kjeldahl analysis 120  
 mercury 123  
 selenium 123
- Catheterization miscellaneous references 356 357
- Cavallini D 150
- Cell volume miscellaneous references 366
- Cellophane transparent, 83
- Cells determination of  
 cholinesterase 274 2 6  
 pH measurement 235 241  
 potassium 92  
 sodium 90 91  
 total base 83 87  
 total solids 82 83  
 total water 82 83
- Centrifuge Clay Adams 20
- Cephalin flocculation miscellaneous references 418
- Cerebral circulation miscellaneous references 357
- Cerecedo L R, 152 184, 213
- Cesium oxide phototube 42
- Chalkoff I L 151 212
- Chambers D C 153
- Chanutin A 151
- Charcoal activated Norit 249

- Chargaff E 148 151  
 Charney J " 3  
 Charts alignment, D Oague 18  
 Chatfield C 3 5  
 Cheldelin C 37  
 Chen G 2 9  
 Chen Y C 118  
 Chopard F P 99  
 Chloride determination in  
   serum 105 108  
   urine 105 109  
 Chloroform 1 0  
 Chlorophenol indicator 35  
 Chlorophyll miscellaneous references 33  
 Chloroplatinic acid solution 91  
 Cholesterol  
   acetone alcohol mixture 171  
   acetone ether mixture 171  
   calculation 173  
   determination of  
     combined 171 173  
     esterified 171 173  
     free, 171 1 3  
     total (Bloor's) 169 1 0  
   development of color 170  
   isolation 1  
   precipitation with digitonin 171  
   standard solution preparation 1 0 171  
 Cholesterol miscellaneous references 178  
 Choline microtological miscellaneous references, 39  
 Cholinesterase determination in  
   cells 4 4  
   serum 42 6  
 Cholinesterase miscellaneous references 418  
 Cholinesterase activity miscellaneous references 418  
 Chou, B F 151  
 Christ C L 1 0  
 Christman A A 145  
 Chromioproton miscellaneous reference 13  
 Church C F 376  
 Chymotrypsin miscellaneous references 8  
 Circulation true miscellaneous reference 3,7  
 Claassen A 116  
 Clarke  
   coating solution  
   hydrophobic coating of micro gas analyzer 3 0  
 Clark F P 93  
 Clark L G 39  
 Clark L C Jr " 0  
 Clark W 418  
 Clark W M 34 3  
 Clark's borate buffer mixture pH 7.6 10 0  
   37  
 Claassen S W 14  
 Clifton F F  
 Clinical laboratory procedures 206 4  
   hematology 39 398  
   miscellaneous references to textbooks 396  
   Clo unit  
     calculation 3 9 386  
     definition 3 9  
 Clotting time miscellaneous references 401  
 Clough M C 368  
 Clouh L W 206  
 Cobalt miscellaneous references 116  
 Cohen H 19  
 Cohen H R 117  
 Cohen J A 17  
 Cohen P P 358  
 Cohn C 103 151  
 Cohn E J 15  
 Colcher H 416  
 Cold pressor tests miscellaneous reference 4 1  
 Coleman Jr spectrophotometer  
   adjusting the mechanical galvanometer  
     zero 36 37  
   apparatus 34  
   calibration 36 38  
   general principles 33 36  
   monochromator 33 34  
   optimum transmittance range 36  
   photometer 34 36  
   routine use of the instrument 38  
 Coleman Photofluorometer Model 17B 1  
   57  
   adjusting the instrument 51  
   apparatus 5  
 Collection of expired air 331  
 Collection of feces using  
   carmine, 23  
   charcoal 6  
   malachite green 6  
 Collection and storage of feces with carmine 6  
 Collection and storage of specimens 19 6  
   blood drawing 21  
   feces 4  
   food (whole livers in the field) 19 0  
   general considerations 19  
   to use 25  
   urine in the field 24  
   whole blood and serum in the field 20 2  
 Coller H B 150  
 Collimating mirror 48  
 Collip J B 93  
 Convert J R 417  
 Combined carbon monoxide and oxygen method (1 B) 88 89  
 Combined reagent for liberation of oxygen  
   and carbon dioxide 73  
 Comfort M W 470  
 Composition of normal arterial human blood 311  
 Computation of theoretical balances 393 394  
   calcium 393  
   nitrogen 393  
   phosphorus 393  
   potassium 393 394  
   sulfur 393 394  
 Conroe J H Jr 81  
 Concentration transmittance calibration  
   graph 3 33  
 Condenser, reflector 1 5  
 Condit H 13  
 Conn J W 404  
 Connor A C 11  
 Connor F T 194 190  
 Conley H 144 149  
 Consalvo C F 01 401  
 Conslazio W V 83 8 90 91 9 15  
 Contreras L A 406

## Conversions

- centrifuge to Fahrenheit scales 38.  
 dehydroascorbic acid with Norit 195  
 hemoglobin to oxyhemoglobin 209  
 temperature, barometer wind velocity 383  
 weights and measures different systems 405 407  
   area—metric system, 425  
   area—U S system, 425  
   capacity—metric system 406  
   capacity—U S system dry 406  
   capacity—U S system liquid 406  
   length metric system, 425  
   length—U S system 425  
   mass—metric system, 426  
   mass—U S system avoirdupois commercial 406  
   time, 427  
   U S apothecaries fluid measure, 406  
   U S apothecaries' weight 406  
   volume—metric system 425  
   volume—U S system 405
- Conway A C 401  
 Cook, E F 193  
 Cooke W D 117  
 Coon J M 401  
 Coors A H 219  
 Copper in serum determination of, 100 100  
   standard solution preparation 101  
 Copper miscellaneous references 116  
 Copper sulfate solution of 175 198  
   standard solution preparation, 127 1 8  
 Corex cells Beckman Model DU 39  
 Corcoran A C 167 218  
 Corley R C 96 110  
 Corn starch suspension 206  
 Coronary blood flow miscellaneous references 357  
 Cortical steroids, miscellaneous references, 218 219  
 Coulson, R A, 189  
 Coulson Smith C 356  
 Counting chamber, Leav 399  
 Coupling reaction, with 2 4 dinitrophenyl hydrazine 197  
 Cournaud A 351, 356 35  
 Court A 389  
 Coward K H 210  
 Cowden D J 400  
 Cowdry E V, 389  
 Cramer H 422  
 Crammer, J L 150  
 Crandall L A, 364  
 Creatine in urine, determination of 139 140  
 Creatinine in urine determination of 138 139  
 Creatinine preformed 139  
 Creatinine zinc chloride standard solution preparation 133  
 Cresol phthalein indicator 235  
 Cresol red, indicator 235  
 Crowley L V 373  
 Croxton F E 402  
 Crucibles platinum 96  
 Cruger D 152  
 Cullen G E, 103, 299  
 Curtis G M 117

## Cultures standard

- Lactobacillus arabinosus* 48 250  
*Saccharomyces carlsbergensis* 250  
 Curves spectral transmittance 29 30  
 Cutler, H H, 400  
 Cyanmethemoglobin conversion 230 231  
 L-cystine standard solution, preparation, 43  
 Cytochrome C miscellaneous references 233  
 Cytochrome oxidase miscellaneous references 208

## D

- Dahr E 118  
 Dak H C 50  
 Daly C A 135 235  
 Damon C M 190  
 Danielson I S 150  
 Dann W J 212  
 Darling R C 061 26 273 351 369 378  
 Daughaday, W H 218  
 Davidson C S 402  
 Davidson L S P, 418  
 Davis V I 153  
 Decalzo, activation of, 180  
 Definition of  
   atomic weight 76  
   calorie 336  
   equivalent weight 76  
   fluorometry, 50  
   gram equivalent weight 76  
   milliequivalent 76  
   millimicron 29  
   molar solution 76  
   molecular weight 76  
   normal oxidizing solution 76  
   normal reducing solution 76  
   normal solution 76  
   pH 2 5  
   phosphatase unit 000 223  
   phosphorescence 50  
   pK 230  
   selective absorption 28  
   spectrophotometer 33  
   standard solutions 76  
   transmission, 27  
   transmittance 27 28  
   T C (volumetric apparatus) 69  
   T D (volumetric apparatus) 69  
   volumetric analysis 68  
 Degeydt E 150  
 Dehydroascorbic acid conversion with Norit 190  
 7 dehydrocholesterol 181  
 Dehydro l as orbic acid in tissue determination 000  
 Dekans, J, 219  
 Dekema H, 418  
 de Laat B M 100  
 Delor R A 060  
 De Luca H A, 116  
 De Ment J, 50  
 De Ment J law of absorption 50  
 De Michele G 007  
 Dennis W 96 110  
 Density of mercury at various temperatures 314  
 Denstedt O F 166

- Dent C L 148 149  
 Deproteinizing reagent trichloroacetic acid 25  
 De Rittter E 213  
 de Verdier C H 149  
 De Witt J B 214  
 Dexter L 357  
 Diazo reagent preparation 413  
 2,6-dichlorophenolindophenol buffered solution preparation 201  
 Dicken D M 260  
 Didymium calibrating filter for Coleman Jr spectrophotometer 37  
 Digestion mixture for NPV 135  
 Digestion tube Folin 135  
 Digestion tubes 186 187  
 Digonide total cholesterol 17  
 Digitonin solution preparation 171  
 Diketo glycolic acid in tissue determination 199 200  
 Dill D B 87 335 301 305 3 7 346  
 Dilute tungstic acid solution preparation 154  
 Digenanase E 150  
 4-Dinitrophenylhydrazine reagent preparation 196 199  
 Diadrasat miscellaneous reference 411  
 Dische Z 167  
 Dishes petri 37 38  
 Disodium hydrogen phosphate solution preparation 21  
 Disodium phosphate 143  
 Dispersion data for spectrophotometer model DU 45  
 Dissociation curves 305 310  
 blood 309  
 oxygen in human blood 309  
 Distilled water metal free preparation 101  
 Diuresis 3  
 Dixon M 80  
 Doboulor P 11  
 Dobriner K 0  
 Do Cagne al pigment hart 14  
 Dodge E 219  
 Doherty E A Jr 33  
 Dole V P 1 7 130  
 Domingo W P 80  
 Donnan ratio calculation 310  
 Dorfman A 8  
 Dornbush A C 60  
 Douglas bag for collection of expired air 330 333 36  
 Downie G I 7  
 Drabkin D L 33  
 Dreker I J 0  
 Drew R 148  
 Drury P E 178  
 Dryden E C 151  
 Du Bois F I 33  
 Du Bois & Meek formula for surface area 335  
 Dunabin J I 145  
 Dunford, R. A 19  
 Dunn M R 149 54 60  
 Dunstan, B 16  
 Durso J G 9  
 Dyer C G 118
- E**
- Earle D P Jr 411  
 Ebert R V 401  
 Eckenhoff J E 356 357  
 Eder H A 127  
 Edgeston H A 400  
 Edman P 150  
 Edsall J T 147  
 Edwards G A 79 8  
 Edwards H T 157  
 Egnas E 13  
 Egnard J 178  
 Ehrlich aldehyde reagent preparation 415  
 Eickhorn F 150  
 Egen E 151  
 Einstein al photochemical equivalence law 51  
 Ekman B 150  
 Elberh I I M 148  
 Electrocardiograph miscellaneous reference 357  
 Electron microscope 80  
 Element by atomic number 4 9  
 Electrophoresis 80  
 Ellinger P 189  
 Ellis D 356  
 Ellis G H 117  
 Ellis J W 118  
 Elman R 366  
 Elvehjem C A 289 244 247 49  
 Emerson C P Jr 401  
 Emerson K Jr 197 130  
 E. M. K. 40 241  
 Emulsion Law of Stokes 51  
 Emmons C W 18  
 Emmett A D 13  
 Emulifying agent ascara 2 1  
 Emulsion of olive oil 2 1  
 Energy expenditure calculation 336  
 Engel F L 8  
 Engel L L 2 0  
 Entenman C 1  
 Enzymes 2 1 9  
 acid phosphatase 203  
 alkaline phosphatase 203  
 amylase 6-2 7  
 anti protease 7  
 anti tryptic activity 7  
 arginase 2 7  
 cholinesterase 2 4 0  
 chymotrypsin 2 8  
 clarsae 187 190  
 cytochrome oxidase 2 8  
 hyaluronidase 2 9  
 lipase 2 1 0 9  
 lipoxidase 2 8  
 peptic activity 9  
 peroxidase 2 8  
 phosphatase 2 8  
 proteolytic activity 2 4  
 takadiastase for thiamine analysis 19  
 takadiastase solution 19  
 trypsin 2 8  
 urease 145  
 Enzymes miscellaneous references 2 9 9  
 Eosinophil count, direct method 399 401  
 calculation 401  
 diluting fluid preparation 401

- Eppinger E C 357  
 Epstein J H 186  
 Equilibration of blood against gas mixtures  
 of known composition 301 305  
 Equivalent weight definition 76  
 Erdos J 118  
 Ergosterol 181  
 Ernste A C 399  
 Erythrocytes miscellaneous references 401  
 Estrogens miscellaneous references 419  
 Ether alcohol mixture solution for extraction  
 of cholesterol preparation 170  
 Ethyl alcohol aldehyde free standard solu-  
 tion 161  
 Ethylene glycol for absorption of water  
 vapor 265  
 Etton J 213  
 Evans E I 366  
 Evans H J 80  
 Evans J M 357  
 Evans J S 419  
 Evans blue dye T 18<sup>4</sup> sterile solution prep-  
 aration 358  
 Evelyn K A 230 231 413  
 Expired air collection  
 Douglas bag 330 333  
 gasometer 377 379  
 Extracellular fluid  
 calculation 360 366  
 change in weight 394  
 thiocyanate space 364 366  
 Extracellular fluid miscellaneous references  
 366  
 Extraction of cholesterol by ether alcohol  
 mixture 170  
 Extraction unit for 17 ketosteroids 418  
 Ezekiel M 4 2

## F

- F fluorescent factor 40, 409  
 Farah A 366  
 Fats and derivatives 168 178  
 acetone bodies 175 177  
 free and combined cholesterol 171 173  
 lipid phosphorus 173 174  
 total cholesterol (Bloor's) 169 170  
 total lipids 168 169  
 Fats and derivatives miscellaneous refer-  
 ences 178  
 acetoacetic acid 178  
 acetone 178  
 acetyl acetone 178  
 cholesterol 1 8  
 fatty acids 178  
 formic acid 178  
 hydroxybutyric acid 178  
 phospholipid 178  
 total fat 178  
 Feces  
 collection and storage 26  
 collection and storage in the field 26  
 determination of  
 biotin 249 250  
 calcium 96 97  
 carotenoids (Carr Price) 181  
 folic acid, 244 247

- Feces determination of—Cont d  
 pantothenic acid 51  
 phosphorus inorganic 110 111  
 potassium 93  
 pyridoxine 51 53  
 riboflavin 19 195  
 sodium 90  
 thiamine 180 189  
 total nitrogen 123  
 total solids 82 83  
 total water 82 83  
 urobilinogen 416 417  
 vitamin A (Carr Price) 181  
 Feces normal values 436 437  
 Fee D A 152  
 Feldstein M 233  
 Fenn W O 357  
 Ferguson M H 166  
 Ferric alum solution indicator 106  
 Ferric iron solution (gum ghatti) 154 155  
 Ferric nitrate solution 364  
 Ferric sulfate 155  
 Ferricyanide solution for Roughton Scho-  
 lander method 277 286  
 Ferris E B 357  
 Ferrous sulfate solution 415  
 Fibrin 23  
 Fibrinogen in plasma 134 135  
 Fibrinolysis miscellaneous references 402  
 Field methods 201 211  
 ascorbic acid in serum and urine 401 404  
 N-methylmaleimide 207 211  
 riboflavin 304 307  
 thiamine 07 211  
 Fieser L F 40  
 Fieser's modified reagent for absorbing oxy-  
 gen 467  
 Filters  
 Bureau of Standards 50  
 carbon yellow 50  
 cobalt blue 50  
 copper green 50  
 selenium orange 50  
 calibrating didymium 37  
 Jena sintered glass 8  
 red purple Corax A 43  
 Wood type 57  
 Fine S D 117  
 Fineman A 218  
 Fink K 148  
 Fink R M 148  
 Finkelstein M 219  
 Fisher R A 422  
 Fisher R B 149 152  
 Fisher volstat 244  
 Fiske C H 109 119 15  
 Fitch H M 237  
 Fitz E 175  
 Flame photometer 80  
 Flame photometer Beckman 59 65  
 computations 6  
 discussion of sources of error 62 65  
 general construction 59  
 general principles 59  
 procedure for routine analysis 61 62  
 setting up the instrument, 60 61

- Flachka H 11  
 Fleischer B C 90  
 Flood A E 150  
 Florry L 151  
 Fluorocence laws 0 7  
 Fluorometric methods using the field visual  
   fluorometer 5 58  
   apparatus 5 58  
   calculation and examples 9  
   precautions 58  
   principle 7  
   procedure 58  
 Fluorometry 0 1  
   definition 0  
   general principles 51  
   laws of fluorescence 0 1  
   visual 5 58  
 Fok N L 411  
 Fok P P 411  
 Fodor P J 8  
 Folch J 118  
 Folic acid miscellaneous references 212  
 Folic acid pteroylglutamic acid 44 47  
   determination in  
     blood 44 47  
     feces 244 47  
     food 14 17  
     plasma 14 47  
     urine 14 17  
   standard solution preparation 1  
 Folin O 110, 114 115 136 139 140 140  
   153 154 166  
 Folin digestion tubes 135  
 Folkers K 00  
 Fung A 346  
 Fondarat J 11  
 Food  
   collection and storage 10 0  
   collection and storage in the field 10 0  
   determination of  
     botulin 049 0  
     calcium 96 97  
     folic acid 14 14  
     nicotinic acid microbiological 43  
     pantothenic acid 01  
     phosphorus inorganic 110 111  
     potassium 93  
     pyridoxine 0  
     riboflavin 10 100  
     sodium 90  
     thiamine 146 189  
     total nitrogen 1 3  
     total solids 8 81  
     total water 8 83  
     tryptophan 5 56  
 Food consumption surveys 3 3 3  
 Forbes A P 0 0 4 0  
 Forbes T P 00  
 Forbes W H 030  
 Formalin solution preparation 390  
 Formic acid miscellaneous references 1 8  
 Forsham I H 19 412  
 Foster J C 11  
 Fowweather F 8 178  
 Fox S H 13 60  
 Fränkel Conrad H 0 8  
 France W G 80  
 Franke E 143  
 Frankl W 260  
 Franklin A E 150  
 Franklin A L 060  
 Franklin M 418  
 Frankton J E 153  
 Frazer R W 0 0 4 0  
 Friedman H A 010  
 Frier C N 51  
 Friday S J 356  
 Friedemann and Graesser lithium lactate solu-  
   tion standard preparation 159  
 Friedemann T E 166 181 113  
 Frontish N 150  
 Functional tests of endocrine gastric and  
   other systems 419-4 1  
   pituitary ACTH test 419 4 0  
   miscellaneous references 400 401  
 Furchgott R F 19  
 Functional residual air 351 353  
   apparatus 351 35  
   calculation 353  
   principle 351  
  
 G  
 Gaebler O H 118 41  
 Gage A P 3 9  
 Galactose miscellaneous references 418  
 Gallagher T F 18  
 Galston A W  
 Garner R J 117  
 Gas analyzers 31 3 6  
   Haldane apparatus 310 319  
   cylinder and micrometer 3 0 3 6  
 Gas mixtures preparation 301 303  
 Gas tension in blood miscellaneous references  
   267  
 Gas volume  
   calculation 261  
   general equation 61  
 Gasometer Tissot 3 330  
 Gastric function miscellaneous references  
   4 0  
 Gates F M 117  
 Geffen C 1  
 Gebel C 153  
 Genghof D R 06  
 Gentzko C J 140 396  
 Gerabaty J T 410  
 Gerish I 289  
 Ghosh N C 213  
 Gibson J Q H 119 163 40  
 Gills J F 149  
 Gil W E 4 1  
 Gillam A L 16  
 Gillard J 1 0  
 Gissel L 101  
 Gittleman I 00  
 Given H 4 150  
 Glass distilled water preparation 80  
 Glass electrodes 020 04  
   apparatus (Daly Forbes) 039  
   principle 035





- Hawkins M H 4 1  
 Hawkins V 3 416 4 1  
 Haynes F W 357  
 Hazel F 117  
 Heard R D H 3  
 Heat balance clothing 376 381  
   thermal insulation of clothes the ' Clo  
   unit 379 380  
   thermal significance of clothing 380  
 Heath C W 373  
 Heinen A J 118  
 Helmert O M 13  
 Helmholtz H F Jr 357  
 Hematocrit  
   average normal values 133  
   calculation 131  
   values (copper sulfate method) 130  
 Hematology, 397 400  
   direct eosinophil count 399 401  
   mean corpuscular volume calculation of  
   hemoglobin 397 399  
   hemoglobin concentration, 397 398  
   sedimentation rate calculation 399 400  
 Hematology miscellaneous references 401  
   400  
   anticoagulant 401  
   Ashby count 401  
   bone marrow 401  
   clotting time 401  
   erythrocyte 401  
   fibrinolysis 40  
   icteric index 40  
   leukocytes 400  
   platelet 400  
   prothrombin time 40  
 Hematoma prevention of 41  
 Hemingway A 81  
 Hemoglobin  
   calculation 131  
   calculation from oxygen capacity 306  
   calibration chart 30  
   concentration average normal value 131  
   concentration in terms of grams hemo  
   globin 133  
   concentration in terms oxygen capacity  
   133  
   standard normal 30  
   whole blood determination of (copper sul  
   fate method) 1 7 133  
 Hemoglobin manometric method miscel  
   laneous references 399  
 Hemolysis prevention of 41 20  
 Hemolysis in urine miscellaneous refer  
   ences 418  
 Henderson I J 301 305 310 400  
 Henderson I 349 350  
 Henderson L apparatus for collecting alive  
   clear air during exercise 349 350  
 Hendlin D 40  
 Henly A A 419  
 Hennessy D J 141 413  
 Hennessy absorption lumens 184  
   in parietal  
   solution 301  
   solution preparation 40  
 Hepatic misdiagnosis references 4 1  
 Herrick F J 4 1
- Hess W C 100 18  
 Hessel F H 453  
 Hettig R A 11  
 Hier S W 39  
 Hill B R 2 0  
 Hill J M 118  
 Hill R M 79  
 Hiller A 197 1 99 300 311  
 Hills A G 219 419  
 Himmer P F 218  
 Hilly W W 8  
 Hippel J A 80  
 Hippuric acid miscellaneous references 418  
 Hirsch E F 178  
 Hirsch H L 357  
 Hirst E L 150  
 Histamine miscellaneous references 219  
 Histamine microbiological miscellaneous ref  
   erences 409  
 Hitchings G B 157  
 Hoag E H 30  
 Hoagland C L 418  
 Hobbs H L 460  
 Hoch F I 119  
 Hochberg H 413  
 Hochberg M 210 41  
 Hoffman W B 116  
 Hoffmann C E. 80  
 Hogues K R 151  
 Holden H F 433  
 Hollen M 189  
 Holiday E P 150  
 Hollander F 4 0  
 Holmes E 80  
 Holmes P 149  
 Holtorf A F 4 0  
 Homburger F 149  
 Homogenate 19  
 Homogenizer 149 45  
 Honnato R 400  
 Hooker C W 460  
 Hormones 1 8  
 Hormones in miscellaneous references 418 400  
   acetylcholine 419  
   adrenalin 18  
   androgen 18  
   cortical steroids 419  
   corticotropin 18 19  
   estrogens 419  
   gonadotropin 419  
   histamine 219  
   parathormone 419  
   pituitary adrenocorticotropin 19  
   protein bound iodine 219  
   steroid 419 4 0  
   thyroxine 0  
 Horn M J 450 460  
 Horowitz N H 39  
 Horrocks R H 1 0  
 Horvath J. 200  
 Horvath R M 4 0  
 Hoxley R J 400  
 Hough L 150  
 Howard J. 414  
 Howe P E. 1 0 3 3  
 Huang H T. 4 0

Huckabay W B 178  
 Huddleston B 153  
 Hudson P B, 2<sup>nd</sup> 8  
 Huelin F E 212  
 Huerga J, 418  
 Huff J W 213  
 Hughes W L Jr 152 402  
 Huis in t Velt, L G 150  
 Hull T Z 368  
 Humidity and absolute vapor tension charts, conversion 380 386  
 Hunter D 152  
 Hunter T H 152 368  
 Hurtado A 300  
 Hurxthal L M 346  
 Hutching B L 249 260  
 Hyaluronidase miscellaneous references 2-8  
 Hydrochloric acid standard solution preparation 77  
 Hydrogen ion concentration 234 243  
   buffer solutions 234 235  
   glass electrode 235 238  
   indicators 234  
   measurement of pH of solution 240 243  
 Hydrogen ion concentration miscellaneous references 243  
   plasma pH  
   urine pH 243  
 Hydrogen manometric method miscellaneous references 300  
 Hydrogen peroxide 3% solution preparation 190 193  
   superoxol 174  
 Hydrolysis of tissue 160  
 Hydrolyzer steam 215  
 Hydroxybutyric acid 178  
 $\beta$  Hydroxybutyric acid 176  
 Hypoproteinemia 131

## I

Iceland spar pure calcium carbonate 94  
 Ichiba F 123  
 Icteric index miscellaneous references 402  
 Indicators  
   alizarin red 96  
   bromocresol green 140 185 194 208 235  
   bromocresol purple 235  
   bromphenol blue 235  
   bromphenol red 205  
   bromthymol blue 188 235 244  
   chlorophenol 230  
   cresol phthalein 235  
   cresol red 235  
   ferric alum solution 106  
   metacresol purple 35  
   methyl red 140 168  
   mixed for micro Kjeldahl analysis 100  
   phenol red 93 235  
   phenolphthalein 142  
   phenolphthalein solution preparation 171  
     2 1  
   p nitrophenol 99  
   starch 156  
   thymol blue 145 160 160 235  
   universal 234  
   universal color change 234  
 Indol aldehyde 416

Infrared spectrometer 80  
 Inoculum for microbiological assay preparation for  
   biotin 249 250  
   folic acid 245  
   nicotinic acid 248  
   pantothenic acid, 251  
   pyridoxine 253  
   riboflavin 254  
 Inositol standard solution preparation 5  
 Interstitial spaces distribution of thiocyanate in 364  
 Instrumentation 27 81  
   balances 65 67  
   Beckman photoelectric quartz spectrophotometer Model DU 38 50  
   flame photometry Beckman 59 65  
   fluorometric methods (field visual) 57 58  
   fluorometry general principles 50 51  
   photofluorometry the Coleman Model 10B 51 57  
   spectrophotometer Coleman Jr Model 6 33 38  
   spectrophotometry general principles 27 33  
   volumetric analysis 68 79  
 Instrumentation miscellaneous references 79 81  
 Intracellular fluid change in weight 394  
 Inulin clearance miscellaneous references 411  
 Iodate manometric method 300  
 Iodate miscellaneous references 300  
 Iodine  
   preparation of solution 159  
   standard solution 79 159 165  
   reference 117  
 Iodoplatinate, 92  
 Irby V, 150  
 Iron inorganic in serum determination 95 100  
   reference 117  
   standard solution preparation 99  
 Irving L 79 80 90  
 Isaacson V I 403  
 Isbell E R 260  
 Isbister J 151  
 Iselin B M 228  
 Iseri T 117  
 Isocamyl alcohol 101 192  
 Isobutyl alcohol 185 204  
 Isolation of free cholesterol 170  
 Isoleucine microbiological miscellaneous references 259  
 Ivy A C 420

## J

Jacobs M D 401  
 Jacobson S D 117  
 Jaffe H 218  
 Jager B V 80  
 Jailer J W 219  
 James G W III 366  
 Janney N W 403  
 Jena sintered glass filter, 87 358  
 John, H E 153  
 Johnson, E A 150

Johnson R. E. 165 301 367 369 373 405  
 Johnson R. M. 112  
 Johnson T. A. 21  
 Jones D. D. 9 260  
 Jones F. O. 373  
 Jones F. 118  
 Jones J. K. N. 104 140  
 Jones L. M. 2  
 Jones M. E. 118  
 Jonxis J. H. P. 401  
 Jorgensen K. S. 218  
 Jukes T. H. 260

## K

Kagan B. M. 150  
 Kahn B. S. 116 118  
 Kallmann S. 115  
 Kalsbeek F. 7  
 Kaltreider N. L. 357  
 Kamen M. D. 117  
 Kann S. 14  
 Kaplan D. 214  
 Kark R. M. 418  
 Kass F. H. 9  
 Katz, L. B. 337  
 Kaufman S. 116  
 Kawasaki C. 214  
 Kazmun V. E. 19  
 Kelly F. J. 368  
 Kendall F. E. 418  
 Kennedy J. A. 307  
 Kenney J. F. 4  
 Kenney R. A. 15  
 Kepler E. J. 219 4 0  
 Kerwin L. 80  
 Keston A. S. 148 149  
 Keto acids total increased during  
 exercise 165  
 thiamine deficiency 163  
 ketosis 165  
 17 Ketosteroids in urine  
 collection and storage 17  
 determination of 15 18  
 standard curve 17  
 Kety S. S. 337  
 Keutmann E. H. 151  
 Keys A. B. 83 103 1 0 13  
 Kidney function 409 411  
 phenolsulphonphthalin test 410 411  
 renal clearance test. miscellaneous ref  
 erences 411  
 titratable acidity 409 410  
 King C. G. 1  
 King F. J. S  
 Kingsley G. R. 151 166 178  
 Kins R. L. W. 117  
 Kirby H. 144  
 Kirk P. L. 81 153  
 Kitchel M. P. 4 1  
 Kitzer G. 98  
 Kjeldahl micro analysis 1 0 1 3  
 catalyst preparation 1 1  
 microanalysis apparatus 1 1  
 Klein D. 167  
 Klement R. 116  
 Klendshoj N. C. 333

Klimenko V. G. 118  
 Klucner R. G. 150  
 Klyne W. 80  
 Kmiecik T. C. 113  
 Knippenberg E. 116  
 Kobro M. 218  
 Koch F. C. 118 145 219 220  
 Koch W. 214  
 Kodacek E. 14  
 Koditschek L. K. 60  
 Kohn E. 218  
 Kolmer J. A. 396  
 Kolthoff I. M. 81  
 Konusky F. P. 260  
 Korenman I. M. 116  
 Kramer B. 93 14 300  
 Kramer H. 161  
 Kramli A. 260  
 Krehl W. A. 47  
 Krogh's new syringe burette 4  
 Krogh syringe pipette 69  
 Kuetter C. A. 193  
 Kugel V. H. 300  
 Kun E. 117  
 Kunkel H. O. 117  
 Kurtz A. K. 4.2  
 Kuttner T. 11

## L

Labhart H. 80  
 LaCour L. F. 149  
 Lacourt A. 117 150  
 Lactic acid determination in  
 blood 157 161  
 muscle 161  
 serum 157 161  
 solution preparation 265 70  
 Lactic acid manometric method miscellane-  
 ous references 300  
 Lactobacillus arabinosus (#8014) 247 248  
 50  
 Lactobacillus casei factor 244  
 Lactobacillus stock culture (#7469) 245  
 Lagerlof H. 206  
 Lambert Beer law 31 30  
 Lampe I. 117  
 Landsberg J. W. 399  
 Landis A. J. 149  
 Landy, M. 60  
 Lang E. F. 80  
 Lange A. N. 4 3  
 Lange J. D. 401  
 Langley W. D. 15  
 Langerstroth G. O. 15  
 Larson H. W. 151  
 Laser H. 80  
 Lauren M. 118  
 Lavietes P. H. 51  
 Laws  
 absorption of De Mont 50  
 emulsion of Stokes 51  
 fluorescence 50 1  
 Lambert Beer 31 3  
 photochemical equivalence of Finstein 51  
 Lawton F. M. Jr., 116  
 Leclerc R., 118

- Leder I G, 213  
 Lederer M 150  
 Lee D H K, 348  
 Leiboff S L 118  
 Leim M 118 153  
 Lemons H, 378  
 Lerman J 200  
 Lester, D 178  
 Letonoff T V 118  
 Leukocytes miscellaneous reference 40-  
 Levenson S M 409  
 Levin N R 208  
 Levine C 148  
 Levine P 153  
 Levinson S A 396  
 Levy E D 213  
 Levy M 148 149  
 Levy counting chamber 390  
 Lewis J C 260  
 Liberation of  
   carbonate 270  
   carbon dioxide with lactic acid 270  
   carbon monoxide and carbon dioxide with  
   hydrosulfite borate solution prepa-  
   ration 275  
   oxygen with potassium ferricyanide solu-  
   tion preparation 267  
 Lieberman, S, 20  
 Lihenthal and Riley technique for obtaining  
   arterial blood 263  
 Lillie P D 389  
 Line charts  
   ascorbic acid 203  
   calculating blood volume and circulatory  
   mass 360  
   calculating BMR from heat production  
   age and sex 340  
   calculating (CO) - (CO) from T<sub>a</sub> and  
   oxygen capacity 307  
   calculating factors to reduce saturated gas  
   volumes to dry volumes at 0 C and  
   760 mm Hg 269 334  
   calculating heat production from surface  
   area and O consumption 339  
   calculating hematocrit 400  
   calculating index of physical fitness 370  
   calculating mean corpuscular volume mean  
   corpuscular Hb concentration and  
   mean corpuscular Hb content 398  
   calculating plasma protein hemoglobin and  
   hematocrit 130  
   calculating relative humidity from wet and  
   dry bulb temperatures 385  
   calculating R Q and true oxygen from an  
   analysis of expired air (insert in  
   back)  
   calculating surface area from height and  
   weight 338  
   calculating T<sub>a</sub> from pCO (total CO) and  
   Hb (in ert in back)  
   calculating (total CO) reduced - (total  
   CO) oxygenated from T<sub>a</sub> and O  
   capacity 308  
   calculating volume of urine per hour from  
   total volume and time of collection  
   23  
 Lineweaver H 28
- Lingane J J 81  
 Lipase in serum determination of 2012  
 Lipase miscellaneous references 23  
 Lipemia alimentary 414  
 Lipid phosphorus determination in  
   blood 173 174  
   serum, 173 174  
 Lipids, total in serum determination of 168  
   169  
 Lipoxida e miscellaneous references 8  
 Lithium carbonate solution preparation 145  
 Lithium chloride solution preparation 143  
 Lithium lactate standard solution prepara-  
   tion 159  
 Little J M 166  
 Littman B 118  
 Liu S H 299  
 Liver function 410 418  
   determination of  
     bilirubin 413 414  
     bromsulphalein test 41 413  
     urobilinogen 415 417  
 Liver function miscellaneous references 417  
   418  
   amyloidosis 418  
   bile in urine 418  
   cephalin flocculation 418  
   cholinesterase activity 418  
   galactose 418  
   general 417  
   hemosiderin in urine 418  
   hippuric acid 418  
   phosphatase 418  
   portal pressure 418  
   prothrombin time 418  
   thymol turbidity 418  
   urobilinogen 418  
 Locke W 219  
 Longnecker W H 148  
 Longworth L G 80  
 Longwell B 19  
 Looney J M 118  
 Lopez J A 198 368  
 Love W 300  
 Low B 100  
 Lowry O H 19 15 193 14 369  
 Loy H W Jr 213  
 Lozner E L 404  
 Lubin W T 357  
 Lubricating grease 67 301  
 Lucia S P 418  
 Leucke P W 259  
 Lugg J W H 148  
 Lugovoy J K 166  
 Lung volume miscellaneous references 307  
 Lutwack R 118  
 Lyons R H 357  
 Lysine microbiological miscellaneous ref-  
   erences 260
- M
- Ma T S 120  
 MacFadyen D A 299  
 MacMillan J M 412  
 McCall P J 18  
 McCance R A 118

- McCarthy E F 63  
 McCoord A B 214  
 McCready M H 8  
 McCreary J F 179  
 McCullagh D R 117 118  
 McCurdy-Larson test of physics 117 118  
 McDonald F G 118  
 McEathern D 18  
 McFate R F 396  
 McGarack T H 8  
 McGargue J S 117  
 McIntire F C 19  
 McKay F A 117  
 McLachlan E A 0  
 McLester J S 12  
 McMichael J 306 35  
 McNabb W M 117  
 McNaught K J 116  
 McOmie F F W 143  
 Magnesium 117  
   determination in serum 97 98  
   standard solution 97  
 Magnesium nitrate solution 101  
 Mainland D 4  
 Malloy H T 30 31 413  
 Malvros H 154  
 Malpre F H 19  
 Malt extract Bacto 5  
 Manganese miscellaneous references 11  
 Manganous sulfate solution 19  
 Manning G H 150  
 Manometric apparatus Van Slyke H 65  
   47  
 Manometric method methane use refer-  
   ence 8  
   amino acids 199  
   bicarbonate 199  
   carbon dioxide 199  
   carbon monoxide 199  
   hemoglobin 199  
   hydrogen 199  
   iodate 199  
   lactic acid 199  
   methemoglobin 199  
   nonprotein nitrogen 199  
   oxalic acid 199  
   potassium 199  
   sugar 199  
   total nitrogen 199  
   urea 199  
 Maren T H 113  
 Marensi A D 21  
 Margen R 117  
 Marinis T F 118  
 Marion D F 41  
 Markham R 150  
 Marks L 4  
 Marrian G F 19  
 Marriot W M 18  
 Marshall L H Jr 353 4  
 Martin A J 143 149  
 Marvin H N 8  
 Mason H 143  
 Mason J M 140  
 Mass spectrometry miscellaneous refer-  
   ence 80  
 Mateer G 41  
 Mateer J G 417  
 Maurer A P 14  
 Maw G A 150  
 Maxey W 166  
 Maximum stability range of  
   N-methylnicotinamide 3  
   riboflavin 3  
   thiamine 3  
 May C D 19  
 Mayer A M 14  
 Mean corpuscular hemoglobin concentration  
   calculation 397 398  
 Mean corpuscular hemoglobin content cal-  
   culation 397 398  
 Mean corpuscular volume calculation 397 398  
 Measurement of  
   boundary potentials by glass electrode 30  
   pH in solution 143  
   osmotic activity 143  
 Meier D J 116  
 Meiklejohn V P 13  
 Meinhart J E 148  
 Meitinen J E 148  
 Melin M 10  
 Melnick D 11 13  
 Membrane tubes 83  
 Menzies G I 30  
 Meniscus mercury 63  
 Menzies J T 10  
 Mercuric acetate sulfate 10  
 Mercuric sulfate solution 10  
 Mercury  
   base free 83  
   cleaning of 6  
 Mercury determination of miscellaneous ref-  
   erences 41  
 Mercury vapor lamp (H 4) 5  
 Merikanto B 116  
 Metabolic product of  
   niacin 180  
   N-methylnicotinamide 180  
 Metabolic wards technique for 392 393  
   general considerations 392 393  
   theoretical balances 392 394  
   theoretical weight 394 395  
 Metabolism during exhausting work 343 344  
 Metabolism during work 341 345  
   apparatus 341 341  
   calculation and data obtained 341 343  
   principles 341  
   procedure 343 344  
 Metabolism after exhausting exercise 344  
 Meta-cresol purple indicator 35  
 Meta-dinitrobenzene solution 215  
 Metal free distilled water 101  
 Metaphosphoric acid stabilizing reagent 5  
 Meteorological measurements 341 343  
   barometric pressure altitude and related  
     variable 381  
   humidity and absolute vapor tension 391  
     392  
   temperature 391  
   wind velocity 381  
   windchill 380  
 Methemoglobin in blood determination 31  
 Methemoglobin manometric method mis-  
   cellaneous references 300

- Methylene blue solution 401  
 Methyl red indicator, 1° 163  
 Mettler A. V 178  
 Meyer, L. M 227  
 Meyers G B 117  
 Michaels G D 117  
 Michel H O 234  
 Michel's method cholinesterase simplified calculation 225  
 Mickelsen O 213  
 Microbiological procedures 244 260  
   agar plate method for biotin 256 259  
   biotin 249 250  
   folic acid 244 247  
   nicotinic acid 247 248  
   pantothenic acid 250 251  
   pyridoxine 251 253  
   riboflavin 253 255  
   tryptophan 255 256  
 Microbiological procedures miscellaneous references 259 260  
   amino acids 259  
   p amino benzoic acid 260  
   arginine 259  
   choline 259  
   glutamic acid 259  
   histamine 259  
   isoleucine 259  
   lysine 260  
   manganese 260  
   pantothenic acid 260  
   phenylalanine 260  
   progesterone 260  
   pteronic acid 260  
   riboflavin 260  
   sterols 260  
   tryptophan 260  
   vitamin B<sub>12</sub> 260  
 Microburettes Rehberg type 70 74  
 Micro burners 135  
 Micro Kjeldahl apparatus 1, 1  
 Micrometer burette 75 76  
 Micro respirometer miscellaneous references 80  
 Miller A T Jr 366  
 Milnequivalent definition 76  
 Milkman G A 80  
 Millimicron definition 29  
 Mills M B 199 212  
 Milton R 153 166  
 Mindus R L 212  
 Minerals 82 119  
   determination of  
     calcium 93 97  
     chloride 105 109  
     copper 100 105  
     iron inorganic 98 100  
     magnesium 97 98  
     phosphorus inorganic 109 112  
     potassium 91 93  
     sodium 87 91  
     sulfur 112 115  
     total base, 83 87  
     total solids 82 83  
     total water 82 83  
 Minerals miscellaneous references 115 119  
 Minor A. H. 402  
 Minot A S 213  
 Mirkan A 300  
 Mirror collimating 48  
 Mitchell H K., 150 152 260  
 Molar solution definition 76  
 Mold prevention salicylic acid solution 156  
 Molecular weight definition, 76  
 Molybdate I solution 109  
 Molybdate II solution 109  
 Monge Hedde M F 117  
 Monochromator Coleman Jr 33 34  
 Montgomery, H 134  
 Moos's technique for obtaining arterial blood 263  
 Moore D H 80  
 Moore J W 366  
 Moore R M 401  
 Moore S 149  
 Morell D B 213  
 Morga P W 233  
 Morris C J O R 219  
 Morton, R A 153  
 Mosley V M 149  
 Moubasher P 151  
 Moyer L. S 80  
 Mueller A 213  
 Murrehead E E 118  
 Mulford D J 152  
 Muller P B 214  
 Mulligan G C 68  
 Munro F L 402  
 Munro M P, 40°  
 Munson P L 218  
 Murtaghan M F 116  
 Muscle lactic acid determination of 161  
 Myerson A, 357

## N

- Nachlas M & 2°  
 Nachman H M 366  
 Naftalin L 149  
 Nakamura K. 15°  
 Natelson B 166 178  
 Navis G 356  
 Neidermeyer W, 118  
 Neill J M, 265 267 270 271  
 Near A O 80  
 Nephritis 134  
 Nesbitt S 421  
 Nessler's reagent  
   dilute solution preparation 140  
   solution for NPN preparation 135 136  
 Newburgh L H 378  
 Newerly K 214  
 Newton C J 178  
 Nicholas R E H 233  
 Nicholls J R 178  
 Nicholson D E 148  
 Nickerson J L 356  
 Nickerson M 80  
 Nicotine miscellaneous references 421  
 Nicotinic acid microbiological method 447  
   248  
   determination in  
     blood 248  
     food 248  
     urine, 248  
   principle 247  
   vitamin solution preparation 247

Nicotinic acid miscellaneous references 213  
 Nicotinic acid standard solution preparation 47-48  
 Nielsen A T 213  
 Niemann C 28  
 Nilsson T 149  
 Nitrogen  
   calculation of theoretical balances 393  
   dissolved in blood 98-985  
     acid phosphate buffer 282  
     apparatus 4  
     bicarbonate solution 3  
     calculation 284  
     hydroxylate solution 283  
     principle 3  
     procedure 283-284  
   soluble 135  
   standard solution preparation 135  
 p-nitrophenol indicator 99  
 Nish G A Jr 38  
 N-methylthiopyridine  
   determination in  
     urine 149  
     urine field method 6-11  
   maximum stability range 23  
   standard solution preparation 159  
 Nomogram for  
   calculating excretion of vitamins from  
     urinary volume and concentration  
     of vitamins 406  
   calculating relative humidity from dry  
     bulb temperature and differences  
     between wet and dry bulb 386  
   calculating windchill from tempera-  
     ture and wind velocity 333  
   estimating partial pressure of a gas from  
     barometric pressure and percentage  
     of gas found 318  
 Noninterfering substances in ascorbic acid  
   analysis (hydrolysis method)  
     fructose 197  
     glucose 197  
     glucuronic acid 197  
     xylose 197  
 Nonprotein nitrogen  
   determination in  
     blood 133-137  
     serum 133-137  
 Nonprotein nitrogen manometric method mis-  
   cellaneous references 309  
 Norit  
   acid washed preparation 396  
   activated charcoal 15-17  
 Normal butyl alcohol 94  
 Normal oxidizing solution definition 6  
 Normal ranges  
   feces 436-442  
   plasma (or serum) 42-433  
   urine 431-435  
   whole blood 430-431  
 Normal reducing solution definition 6  
 Normal solution definition 76  
 Normal values  
   plasma specific gravity  
     men 131  
     women 131

Normal values—Cont'd  
   whole blood specific gravity  
     men 131  
     women 131  
 Normochromic blood, 131  
 Normocytic blood 131  
 Nutrient 19  
 Nutrient per pound of food groups as per  
   chapter 3:6377  
 Nylon G 366  
 O  
 Oastler E G 6  
 Obermer E 133  
 Oesterling M J 199  
 O'Leithorpe C C 116  
 Olbrich O 166  
 Oleic acid 169  
 Olive oil emulsion solution preparation 271  
 Olsen J C 473  
 Open circuit apparatus for measurement of  
   functional residual air 5  
 Opieska Blauth J 19  
 Oppenheimer M J 34  
 Optical density 33  
 Optimum transmission range Coleman Jr  
   spectrophotometer 36  
 Orr M 366  
 Osier B L 219-213 396  
 Osool E F 390 417  
 Osterberg A F 40  
 Ott M L 8  
 Ottolenghi L 162  
 Overell B T 148  
 Oxalate as an anti-oxidant 131  
 Oxalic acid  
   manometric method 306  
   preservative 19  
   solution preparation 95  
   standard solution preparation 7  
 Oxalic acid miscellaneous references 309  
 Oxidizing reagent for thiamine 185  
 Oximeter 89  
 Oxygen  
   determination in blood (Roughton Scho-  
     lander) 7-8  
   apparatus 7  
   calculation, 281  
   principle 7  
   procedure 29-31  
   disoxygenation curves of human blood 309  
   total in arterial blood calculation 305  
   total in blood  
     bound by hemoglobin 305  
     calculated 305-306  
     capacity calculation 306  
     percent saturation calculation 306  
     physically dissolved 305  
   whole blood determination (Van Slyke)  
     27  
   factors (Van Slyke) 279  
 Oxygen absorbents  
   Fieser's modified reagent 26  
   gas analysis 31  
   hydroxylate solution preparation 24  
   Scholander meter gas analyzer 3  
   sodium hydroxylate solution preparation  
     26



- Oxygen capacity in blood determination (Van Slyke) 272 275  
 Oxygen and carbon dioxide in single sample of blood 212 273  
 Oxygen consumption calculation 338  
 Oxygen debt  
   calculation 345  
   definition 345  
 Oxygen reagent for liberation of oxygen 267  
 Oxyhemoglobin in blood determination of 231
- P**
- Pack test of physical fitness procedure 31 373  
 Page I H 167 218  
 Palmer W W 409  
 Palmitic acid 169  
 Pancreas miscellaneous references 420  
 Pantothenic acid "50  
   determination in  
     blood 251  
     feces 251  
     food 251  
     plasma "51  
     urine "51  
   standard solution preparation 200  
 Pantothenic acid microbiological miscellaneous references 260  
 Papain 247  
 Paper chromatography 147 151  
   amino acids 148 149  
   carbohydrates 150  
   methodology 147 148  
 Paper chromatography miscellaneous references 150 151  
 Parahormone miscellaneous references 219  
 Parathyroids miscellaneous references 40  
 Parsons D S 148  
 Partial pressure of a gas calculation 304  
 Partridge C W H "56  
 Partridge S M 150  
 Passel C F 389  
 Patek A J Jr 418  
 Patterson H R "0  
 Patterson W I 151  
 Paul W D 418  
 Peacock G 213  
 Pearson F B 117 259  
 Pearson S "0  
 Peckford, L 150  
 Pecover V 212  
 Pedersen K O 80 81  
 Pedley F G 116  
 Peffer E L 68  
 Penicillin reader Fisher 257 258  
 Pennell R B 233  
 Peptic activity miscellaneous references 23  
 Peptone 249  
   alkali treated 250  
   supplement "49  
 Peripheral hematocrit correction factor 361  
 Perlzweig W A 213  
 Permutit 142  
   preparation 140  
 Peroxidase miscellaneous references 2  
 Petering H G 233  
 Petermann M L 151  
 Peters J P 118 404 411  
 Peters J T 233  
 Peterson W H, 249  
 pH  
   definition "35  
   measurement of blood plasma and RBC 235 242  
   measurement of solutions "4 243  
   procedure 242  
 Phenol red indicator 93 235  
   standard solution preparation 409  
 Phenolphthalein indicator 85 140 171  
   solution preparation 21  
 Phenols miscellaneous references 41  
 Phenolsulfonphthalein  
   standard solution preparation 410 411  
   test in urine 410 411  
 Phenolsulfonphthalein miscellaneous references 411  
 Phenylalanine microbiological miscellaneous references 260  
 Phillips D M P 148  
 Phillips P H 260  
 Phillips R A 127 130  
 Phloxin solution 402  
 Phosphatase alkaline and acid "23  
 Phosphatase miscellaneous references 418  
 Phosphate buffers  
   pH 65 134  
   pH 66 230 231  
 Phospholipid miscellaneous references 178  
 Phosphoric acid definition 50  
 Phosphoric acid syrupy 143 181  
 Phosphorus  
   calculation of thermal balances 393  
   determination (inorganic) in 109 110  
     feces 110 111  
     food 110 111  
     serum 109 112  
     urine 110  
   reference 117  
   standard solution preparation 110 111  
 Phosphotungstic acid 134 143  
 Photochemical equivalent law of Einstein 51  
 Photofluorometer Coleman Model 10B 5157  
   adjustment and replacement of parts 55 67  
   calculation of results 54  
   computation of results 54  
   general construction 51  
   routine measurement 53 54  
   setting up and adjusting 51 53  
 Photometer Coleman Jr spectrophotometer 34 36  
 Photorex solution preparation 21  
 Phototube, cesium oxide, 42  
 Physical chemistry of blood gases 301 311  
   calculation of arterial pCO 306 308  
   calculation involving carbon dioxide 306  
   calculation involving oxygen and hemoglobin 305 306  
   calculation (miscellaneous) of the physicochemical properties of blood 310 311

- Physical chemistry of blood gases—Cont d  
composition of normal human arterial blood 311  
dissociation curves of blood 306 310  
Donnan ratio calculation 310  
temporary effects on base found by plasma proteins 311
- Physical fitness self tests 371 374  
pack test 371 3 3  
step test 3 3 3 4
- Physical fitness test using treadmill 369 371  
calculation of score 369 371  
principle 369
- Physicochemical system 301
- Physiological measurements 261 300  
calculation of gas volumes 261  
techniques for obtaining arterial blood and for handling anaerobically 61 6  
Van Slyke-Neill manometric apparatus 60 67
- Pickels E G 81
- Pierce acid solution preparation 139
- Pigments  
  arotenoid 180  
  hemoglobin 29 31  
  methemoglobin 31 3  
  oxyhemoglobin 31 3  
  sulphemoglobin 31 3
- Pigments miscellaneous references 33  
  carboxyhemoglobin 33  
  chlorophyll 33  
  chromoproteins 33  
  cytochrome C 31  
  hemoglobin 33  
  porphyrins 33
- Piper A 401
- Pineu J B 166
- Pipette Brewer automatic 3  
Pipettes  
  Bureau of Standards 69  
  Fisk 69  
  Micro Polum 1  
  Normax 69  
  Ostwald 71  
  serological 71  
  syringe 71  
  Van Slyke-Neill for blood 1
- Pirani C L 399
- Purdie A 13
- Pittsers 1 0  
  statutory interference of horn test 419 420  
  calculation 4 0
- Putnam interference otropin hormone measurement 419
- pH indicator 33  
pH normal value of 306  
pH standardization of  
  luciferin 0 9  
  nitrogen 134 135  
  folic acid 44 4  
  pantothenic acid 31  
  pH measurement 35 4  
  total carbon dioxide calculation 306  
  total circulating volume of carbon monoxide method 361 367  
  tryptophan 3 36
- Plasma normal values 430 433  
Plasma pH miscellaneous references 43  
Plasma volume blood volume and extracellular fluid 355 36  
  extracellular fluid 361 365  
  plasma and blood volume (dye method) 359 361  
  calculation 3 9 361  
  red cell plasma volumes (carbon monoxide method) 361 363  
  true blood volume 363 364
- Platelets miscellaneous references 40
- Platinum crucibles 96
- Platinum electrode 83
- Platinum wire as a stirring device 160
- Plazin J 1 7 1 2
- Pohle F G 401 402
- Polarography miscellaneous references 81
- Pollard F H 148
- Pollay I F 117
- Polson A 149 151
- Pontacyl carmine B 415
- Pontacyl violet GR 415
- Pope C G 149
- Popjak G 80
- Popper H 418
- Porphyrin miscellaneous references 33
- Portal pressure miscellaneous references 418
- Potassium  
  calculation of theoretical value 383 394  
  determination in  
    diets 91 91  
    feces 91 93  
    serum 91 93  
    tissue 91 93  
    urine 91 93  
  manometric method miscellaneous references 300  
  reference 115
- Potassium chloride  
  saturated solution preparation 1  
  solution preparation 153 64
- Potassium citrate buffer 3
- Potassium dichromate solution 1 3  
  standard solution preparation 9 164 165
- Potassium dihydrogen phosphate solution 1
- Potassium ferricyanide  
  liberation of oxygen 26  
  solution 154 164 31
- Potassium hydrogen phthalate standard solution
- Potassium iodate standard solution 9
- Potassium iodide solution 91 164
- Potassium oxalate blood anti-coagulant 164
- Potassium permanganate standard solution 9 93
- Potassium phosphate monobasic solution 9
- Potassium thiovanate standard solution 9
- Potassium thiovanate standardization 106 10
- Potter V F
- Power M H 114 4 0
- Prater A 33

- Preformed creatinine 139
- Preparation of  
 basal medium for  
 biotin 249 250  
 folic acid 246  
 nicotinic acid 248  
 pyridoxine 252 253  
 riboflavin 254  
 tryptophan 256  
 carbon monoxide 62  
 gas mixtures for equilibration of blood  
 302  
 inoculum (microbiology) for  
 biotin 249  
 folic acid 245  
 nicotinic acid 248  
 pantothenic acid 251  
 pyridoxine 252  
 riboflavin 254  
 Norit acid washed 196  
 Permutil 140  
 samples for folic acid 246  
 serum in the field 21 22  
 standard solutions 77 79  
 standard solutions for field use  
 ascorbic acid 201  
 N<sup>1</sup> methylnicotinamide 210  
 riboflavin 206  
 thiamine 210  
 uricase 145 146  
 zeolite activated 207 208
- Prescott B A 260  
 Prescott F 212
- Preservation of animal tissues for histological  
 studies  
 dry ice 390  
 formalin solution 390
- Preservatives  
 benzoic acid 155  
 oxalic acid for urine 19  
 sodium benzoate solution 221  
 thymol crystals 126  
 toluene for  
 copper determination 23  
 urinary calcium 23
- Prevention of  
 hematomas 21  
 hemolysis, 21 22  
 molds using salicylic acid 156
- Price J I, 373  
 Price V E 208  
 Priestly J G, 312 333
- Primary standards 77 78 94  
 hydrochloric acid solution 78  
 potassium dichromate solution 79  
 potassium hydrogen phthalate 78  
 potassium iodate solution 78  
 sodium oxalate 78
- Procaine miscellaneous references, 421
- Progesterone microbiological miscellaneous  
 references 260
- Propylene glycol 401
- Protein bound iodine miscellaneous refer-  
 ences 219
- Protein metabolism during work 337
- Protein and nitrogen compounds determina-  
 tion in  
 albumin (Howe method), 124
- Protein and nitrogen compounds determina-  
 tion in—Cont d  
 albumin and globulin in serum (Biuret  
 method), 126 127  
 alpha globulin 125  
 beta globulin 125  
 gamma globulin 125  
 creatine 136 140  
 creatinine 136 139  
 fibrinogen 134 135  
 non protein nitrogen, 135 137  
 protein (Biuret method) 123 126  
 protein (copper sulfate method) 127  
 133  
 total nitrogen micro Kjeldahl 1 0 127  
 urea nitrogen 140 141  
 uric acid 143 144
- Protein and nitrogen compounds miscel-  
 laneous references 151 153  
 alanine 151  
 albumin and globulin 151  
 allantoin 151  
 amino acids 151  
 amino nitrogen 152  
 creatine 152  
 cystine 152  
 histidine 152  
 non protein nitrogen 152  
 nucleic acid 152  
 paper chromatography 147 151  
 proteins 152  
 purines 152  
 total nitrogen 153  
 tryptophan 153  
 tyrosine 153  
 urea 153  
 uric acid 153
- Proteolytic activity miscellaneous references  
 228
- Prothrombin time miscellaneous references  
 402 418
- Protoplasm mass change in 394
- Prunty F T G 219 419
- Prussian blue 154
- Pteric acid miscellaneous biological refer-  
 ences 460
- Pteroylglutamic acid (folic acid)  
 determination in  
 blood 244 247  
 feces 244 247  
 food 244 247  
 plasma 244 247  
 urine 244 247  
 standard solution, 245
- Pulmonary and circulatory measurements  
 miscellaneous references 356 357  
 acetylene 356  
 ballistocardiograph 356  
 cardiac output 356  
 catheterization 356 357  
 cerebral circulation 357  
 circulation time 357  
 coronary blood flow 357  
 electrocardiograph 357  
 gas tension in blood 357  
 lung volume 357  
 pressure breathing 357  
 venous pressure 357

Pulmonary measurements 340 357  
 alveolar air in rest 340 350  
 alveolar air in work 3 0  
 cardiac output 353 366  
 functional residual air 3 1 363  
 Pulmonary ventilation calculation 3 0 333  
 Purification of mercury 267  
 Pyridine 04  
 Pyridine of dimethylnicotinamide miscel-  
 laneous references 12  
 Pyridine solution 193  
 4 Pyridine acid miscellaneos r leren et  
 13  
 Pyrioxime (vitamin B) 51 503  
 d rmination in  
 feces 1  
 food 51  
 urine 51  
 sugar and salt solution 51  
 Pyrioxime miscellaneos referen s 13  
 Pyrioxime hydrochloride standard solution  
 52  
 Pyrogallol solution  
 as an oxygen absorber 4  
 Loughlin Scholander r method 7 56 0

## Q

Quante M L 218 224  
 Quartz tube 01  
 Quartz J H 150  
 Quack A J 418  
 Quinine 55  
 Quinine sulfate standard solution for thi-  
 amine preparation 280

## R

Railey J A 50  
 Ramsey L L 1 1  
 Randall F F 3  
 Randolph T G 30  
 Rancy B B 4 1  
 Ranges H 1 356  
 Rappaport 1 13 1 8  
 Raymond W D 113  
 Recommended daily dietary allowances of  
 National Research Council 278  
 Recovery of sodium uranyl nitrate acetate 01  
 Recrystallization of 4:5:1 103  
 Reile H  
 cholinesterase 1 1 6  
 pH measurement 225 217  
 potassium  
 sodium 90 91  
 total base 83 87  
 total circulant volume 361  
 total circulating volume (carbon monoxide  
 method) 361 363  
 total solid 8 81  
 total water 8 83  
 Red oxide of mercury 1 5  
 Red Purple Cores A filter 43  
 Rehberg P B 2  
 Rehberg microburette 074 01  
 Reichard I 100

Reid A 366  
 Reichenstein E C Jr 29 4 0  
 Reinhardt H L 200  
 Reinhold J G 118 166  
 Renal blood flow miscellaneous references  
 411  
 Renal clearance tests miscellaneous refer-  
 ences 411  
 diodrast 411  
 general 411  
 glomerular filtration rate 411  
 renal clearance 411  
 phenol sulfophthalein test 411  
 renal blood flow 411  
 urea clearance 411  
 Respiratory gases estimation 31-3 7  
 Haldane apparatus for measuring carbon  
 dioxide oxygen 31 319  
 Scholander micrometer gas analyzer 3 0  
 3 6  
 Respiratory metabolism 3 7 345  
 basal metabolic rate 337 341  
 calculation 343 337  
 carbon dioxide production 336  
 energy expenditure 336  
 example of all 330 337  
 oxygen consumption 336  
 protein metabolism during work 337  
 pulmonary ventilation 333  
 respiratory quotient (RQ) 333 336  
 true carbon dioxide 333 336  
 true oxygen 333 336  
 collection of expired air Douglas bag tech-  
 nique 330 333  
 collection of expired air Tissot gasometer  
 1 9  
 metabolism during work 341 345  
 Respirometer 81  
 Respirometry meter 80  
 Rhoads C P 0  
 Riboflavin  
 determination in  
 blood 10 103  
 feces 10 103  
 food 10 103  
 urine 10 103  
 urine cell method 04 0  
 maximum stability range 13  
 in subtilis methods determination  
 203 255  
 principle 3  
 principle of  
 clarification of sample 19  
 determination of fluorescence 19  
 liberation of bound albumin as an indicator  
 190  
 oxidation of impurities by adsorption  
 19  
 removal of impurities by adsorption 19  
 standard solution preparation 190 193  
 215 214  
 table for simplified calculation 190 191  
 Riboflavin miscellaneous references 12  
 Rice L W 167  
 Richard D W Jr 331 33  
 Riches L L 60  
 Rider, P E 422

- Preformed creatinine 139  
 Preparation of  
   basal medium for  
     biotin 249 250  
     folic acid 246  
     nicotinic acid 248  
     pyridoxine 252 253  
     riboflavin 254  
     tryptophan 256  
   carbon monoxide 62  
   gas mixtures for equilibration of blood 302  
   inoculum (microbiology) for  
     biotin 249  
     folic acid 245  
     nicotinic acid 248  
     pantothenic acid 251  
     pyridoxine 253  
     riboflavin 254  
   Nitric acid washed 196  
   Permutit, 140  
   samples for folic acid 246  
   serum in the field 21 22  
   standard solutions 77 79  
   standard solutions for field use  
     ascorbic acid 201  
     N<sup>1</sup> methylnicotinamide 210  
     riboflavin 206  
     thiamine 210  
     uricase 145 146  
   zeolite activated 207 208  
 Prescott B A 260  
 Prescott F 212  
 Preservation of animal tissues for histological studies  
   dry ice 390  
   formalin solution 390  
 Preservatives  
   benzoic acid 155  
   oxalic acid for urine 19  
   sodium benzoate solution 221  
   thymol crystals 126  
   toluene for  
     copper determination 23  
     urinary calcium 23  
 Prevention of  
   hematoma 21  
   hemolysis, 21 23  
   molds using salicylic acid 156  
 Price J I, 373  
 Price V E 228  
 Priestly J G 312 333  
 Primary standards 77 78 94  
   hydrochloric acid solution 78  
   potassium dichromate solution 79  
   potassium hydrogen phthalate 78  
   potassium iodate solution 78  
   sodium oxalate 78  
 Procaine miscellaneous references, 421  
 Progesterone microbiological miscellaneous references 60  
 Propylene glycol 401  
 Protein bound iodine miscellaneous references 219  
 Protein metabolism during work 337  
 Protein and nitrogen compounds determination in  
   albumin (Howe method), 1-4  
 Protein and nitrogen compounds determination in—Cont'd  
   albumin and globulin in serum (Buret method) 1-6 197  
     alpha globulin 125  
     beta globulin 1-5  
     gamma globulin 195  
   creatinine 136 140  
   creatinine 136 139  
   fibrinogen 134 135  
   non protein nitrogen, 135 137  
   protein (Buret method) 193 126  
   protein (copper sulfate method) 197 133  
   total nitrogen micro Kjeldahl 1-0 1 7  
   urea nitrogen 140 141  
   uric acid 143 144  
 Protein and nitrogen compounds miscellaneous references 151 153  
   alanine 151  
   albumin and globulin 151  
   allantoin 151  
   amino acids 151  
   amino nitrogen 153  
   creatinine 152  
   cystine 152  
   histidine 152  
   non protein nitrogen 152  
   nucleic acid 152  
   paper chromatography 147 151  
   proteins 152  
   purines 152  
   total nitrogen 153  
   tryptophan 153  
   tyrosine 153  
   urea 153  
   uric acid 153  
 Proteolytic activity miscellaneous references 298  
 Prothrombin time miscellaneous references 402 418  
 Protoplasm mass change in 394  
 Prunty F T G 219 419  
 Prussian blue 154  
 Pteric acid miscellaneous biological references 260  
 Pteroylglutamic acid (folic acid)  
   determination in  
     blood 244 247  
     feces 244 247  
     food 244 247  
     plasma 244 247  
     urine 244 247  
   standard solution 245  
 Pulmonary and circulatory measurements miscellaneous references 356 357  
   acetylene 356  
   ballistocardiograph 356  
   cardiac output 356  
   catheterization 356 357  
   cerebral circulation 357  
   circulation time 357  
   coronary blood flow 357  
   electrocardiograph 357  
   gas tension in blood 357  
   lung volume 357  
   pressure breathing 357  
   venous pressure 357

- Pulmonary measurements 345 357  
 alveolar air in rest 346 350  
 alveolar air in work 350  
 cardiac output 353 3 6  
 functional residual air 351 353  
 Pulmonary ventilation calculation 3 9 333  
 Purification of mercury 267  
 Pyridine 204  
 Pyridine of N-methylacetanamide miscel-  
 laneous references 213  
 Pyridine solution 193  
 4-Pyridoxic acid miscellaneous references  
 213  
 Pyridoxine (vitamin B<sub>6</sub>) 1 253  
 determination in  
 feces 251  
 food 21  
 urine 51  
 sugar and salt solution 12  
 Pyridoxine miscellaneous reference 13  
 Pyridoxine hydrochloride standard solution  
 2  
 Pyrocollol solution  
 as an oxygen absorbent 274  
 Roughton-Scholander method 277 286 27

## Q

- Quante M L 213 214  
 Quarts tubes 91  
 Quastel J H 150  
 Quick, A J 418  
 Quinine 59  
 Quinine sulfate standard solution for the  
 same preparation 193

## R

- Radley J A 50  
 Ramsey L L 151  
 Randall F E 23  
 Randolph T O 390  
 Roney R B 41  
 Ranges H A 2 6  
 Rappaport F 15 1 8  
 Raymond W D 213  
 Recommended daily dietary allowances of  
 National Research Council 378  
 Recovery of sodium uranyl oxalate 91  
 Rectification of ANSA 103  
 Red cell  
 cholinesterase 2 4 6  
 pH measurement 235 24  
 potassium 90  
 sodium 90 91  
 total base 83 87  
 total circulating volume 361  
 total circulating volume (carbon monoxide  
 method) 361 363  
 total solids 82 83  
 total water 8 83  
 Red oxide of mercury 175  
 Red Purple Corex A filter 43  
 Rehberg P B 9  
 Rehberg microburette 9 24 91  
 Reichard P 150

- Red A. 366  
 Reifenstein E C Jr 270 39 400  
 Reinhardt H L 20  
 Reinhold J O 118 166  
 Renal blood flow miscellaneous references  
 411  
 Renal clearance tests miscellaneous refer-  
 ences 411  
 diodras 411  
 general 411  
 glomerular filtration rate 411  
 insulin clearance 411  
 phenolsulfonphthalein test 411  
 renal blood flow 411  
 urea clearance 411  
 Respiratory gases estimation 31 327  
 Haldane apparatus for measuring carbon  
 dioxide oxygen 312 319  
 Scholander micrometer gas analyzer 3 0  
 246  
 Respiratory metabolism 3 7 345  
 basal metabolic rate 337 341  
 calculations 32 337  
 carbon dioxide production on 336  
 energy expenditure 336  
 example of all 336 337  
 oxygen consumption 336  
 protein metabolism during work 337  
 pulmonary ventilation 333  
 respiratory quotient (R.Q.) 333 336  
 true carbon dioxide 337 336  
 true oxygen 333 336  
 collection of expired air Douglas bag tech-  
 nique 330 333  
 collection of expired air Tissot gasometer  
 3 0  
 metabolism during work 341 345  
 Respirometer 81  
 Respirometer micro 80  
 Richards C F 0  
 Riboflavin  
 determination in  
 blood 192 192  
 feces 19 19  
 food 19 192  
 urine 190 192  
 uric acid method 04 207  
 maximum ability range 3  
 microchemical methods determination  
 23 255  
 principle 53  
 principle of  
 clarification of sample 100  
 determination of fluorescence 192  
 liberation of bound alloxazine nucleotide  
 192  
 oxidation of impurities by adsorption  
 192  
 removal of impurities by adsorption 19  
 standard solution preparation 190 193  
 245 234  
 table for simplified calculation 190 191  
 Riboflavin miscellaneous references 213  
 Rice F W 267  
 Richards D W Jr 331 33  
 Ruckes E L 60  
 Rider P R, 400

- Riley R L 356 357  
 Rimington C 233  
 Ring G C 356  
 Pingbom A 116  
 Rusley H A 148 150  
 Robert D 151  
 Roberts E 277 260  
 Roberts S 157 219  
 Robinson F J 420  
 Robinson H W 311  
 Robinson J C 100  
 Robinson P F 201 405  
 Robinson S 341 378  
 Pobschiet Robbins F S 299  
 Robson J S 166  
 Rochelle salts 156  
 Rockland L B 149 160  
 Roe J H 116 166 167 193 199 217 220  
 Rohse W G 153  
 Polf D 411  
 Rollier Smith torsion balance 25 66  
 Roos A 357  
 Root W S 273 361  
 Rose W C 152  
 Rosen F 213  
 Rosenberg H R 212  
 Rosenkrantz H 219  
 Rosenthal O 233  
 Rosenthal S M 178 411  
 Ross E C 119  
 Rotator mechanical 304  
 Roth L W 219  
 Rothera A C H 178  
 Rothschild L 80  
 Boughton F J W 273 275 276 287 285  
 291 301  
 Boughton Scholander syringe pipette 276  
 Hourke M D 399  
 Hourke Ernestine line chart for calculating  
 hematocrit 400  
 Rous P 418  
 Rowntree L G 410  
 Rubin S H 213  
 Rudel H 356  
 Rutter L 148
- S
- Saccharomyces carlbergensis etc k culture  
 352  
 Saffell W B 373  
 Saffran M 166  
 Sagerson, H P 351  
 Sahyun M 151 210  
 Salicylates miscellaneous references 471  
 Saline ammonium sulfate solution prepara-  
 tion 124  
 Saline solution isotonic preparation 20 22  
 Salt solution  
 A for microbiology preparation 245  
 B for microbiology preparation 245  
 isotonic preparation 245  
 Salter W T 117 219  
 Saltzman A H 219  
 Saltzman F 4 0
- Sanford A H 209  
 Sanford A W 396  
 Sanger F 149  
 Sapirstein L A 119  
 Saponification for vitamin A 18  
 Saponification with sodium ethylate 163  
 Saret H P 213 250  
 Sargent F H 201  
 Satterfield G H 112  
 Saur B 152  
 Sawitsky A 27  
 Sayers G 219  
 Sayers M A 219  
 Shorov V 415  
 Scardino P L 411  
 Schaffer N K 153  
 Schaffer P A 173 178  
 Schaffert R R 178  
 Schales O 168  
 Schales S S 166  
 Schunkel M 220  
 Schmidt C F 357  
 Schmidt C F 470  
 Schmidt E G 421  
 Schneider J J 218  
 Schoenheimer F 171  
 Schoff H F 260  
 Scholander P F 75 79 80 81 277 298  
 285 291 295 300  
 Scholander a micrometer burette 576  
 Scholander a micrometer gas analyzer 300  
 306  
 analysis of sample 323 324  
 apparatus 303 321  
 assembling apparatus 321  
 calculation 324 325  
 principle 320  
 Scholander a technique for handling blood  
 anaerobically 263 265  
 Schroeder W A 148  
 Schuette H A 98  
 Schultz A S 251  
 Schuman R L 213  
 Schwartz I L 366  
 Schwartz M C 118  
 Schwartz S 233 415 411  
 Schwartz T B 80 208  
 Schweigert A B 117 259  
 Schwenk E 278  
 Scott M L 233  
 Scott W W 218 411  
 Scoz G 207  
 Scrimshaw, N S 214  
 Senstone C V 228  
 Sebrall W H 255  
 Sedimentation rate calculation 399 400  
 Seegers W H 401  
 Seeley R D 356  
 Seibel, R E 357  
 Selective absorption definition 28  
 Seligman A L 208  
 Seligman A M 151  
 Sendroy J Jr 299 300 311  
 Sensitivity of flame photometer concerning  
 metals 63

Serum collection and storage in the field 20 3

determination of  
albumin  
  Bauret method 1<sup>st</sup>  
  Howe method 1 6 127  
amylase 2 6 7  
ascorbic acid  
  field method 201 204  
  macro method 185 197  
  micro method 198 199  
base total 83 8  
bicarbonate 306  
bilirubin 413 414  
bromsulfalein 41 413  
calcium 93 9  
carbon dioxide 0  
carbonic acid free 306  
carotenoids  
  Carr Price method 179 180  
  dichlorohydrin method 181 184  
chloride 106 109  
cholesterol Bloor's method 169 1 0  
cholesterol combined or esterified 171  
  173  
cholesterol free Sper y's method 1 1  
  173  
cholinesterase 24 8  
copper 100 105  
globulin  
  Bauret method alpha beta gamma  
  1 5  
  Howe method 1 6 1<sup>st</sup>  
iron inorganic 98 100  
lactic acid 15 161  
lipase 12  
lipid phosphorus 173 174  
lipids total 168 169  
magnesium 97 98  
nitrogen total 1 0 1 3  
nonprotein nitrogen 135 137  
pH 306  
phosphatase acid 3  
phosphatase alkaline 3  
phosphorus inorganic 109 11  
potassium 91 9  
protein  
  Bauret method 1 3  
  copper sulfate method 1<sup>st</sup> 7 133  
  Kjeldahl micro method 1 0 1 3  
sodium 87 89  
solids total 8 83  
sulfate  
  etheral, 115  
  inorganic 11 114  
  total 11 114  
urea nitrogen 140 141  
uric acid, 143  
vitamin A  
  Carr Price 179 180  
  dichlorohydrin 181 184  
  water total 8 83  
Weyg M G 133  
Weyfang A P 148 3 6  
Wfortunato T 15  
Waffer P A. 156  
Waffer Hartman Somogyi sugar reagent 16

Shales O 116  
Shales S S 116  
Shank R E 418  
Shankman S 59  
Shannon J A. 411  
Shapiro S 418  
Sharpey Schafer E P 356  
Shaw J L 219  
Sheard C 9  
Shen S S 118  
Shepard M 80  
Sherlock S P 9  
Sherman A. I 19  
Sherman H C 1  
Shettles L B 167  
Shevsky M C 133  
Shimer S R. 13  
Shinowara Q Y 2  
Shipp H L 118  
Shohl A T 81 91 116 118  
Shorb M S 60  
Shorr E 19  
Shull G M 49  
Shun Chuan Liang I W 118  
Siggia S 166  
Sinca pebble 173 1 4  
Silver dichromate solution preparation 168  
Silver nitrate  
  solution preparation 143  
  standard solution 28 106  
Simmons J S 396  
Simmons S M 19  
Simmons D H 366  
Simplified calculation tables  
  amylase activity of serum  
  blood glucose 157  
  cholinesterase 5  
  riboflavin in urine 190 191  
  serum sodium 88  
  step test of physical fitness scoring 3 4  
  total base 86  
  urine sodium 98  
  urine sulfur 114  
Sintered glass (Pyrex) crucibles 1 5  
Siple P A 380  
Siple and Lassell equation 380  
Skatol aldehyde 416  
Skiersky S 116  
Slater E C 113  
Smith B W 2 7  
Smith C G 116  
Smith E L 80  
Smith F A S  
Smith G N 8  
Smith H W 411  
Smith J D 150  
Smith P Jr 219 4 9  
Smith W E 233  
Snedecor G W 4  
Snell E. F 24 253 255 26 2 9 60  
Snow S D 181  
Snodgrass J. 1 8  
Sobel A L 116 181 14  
Sobel H 20  
Sobotka H 214



**Sodium**

- determination in
  - cells 90 91
  - diets 90
  - feces 90
  - serum 87 89
  - tissue 90
  - urine 89 90
- serum simplified calculation table 83
- urine simplified calculation table 89
- Sodium determination miscellaneous references 118
- Sodium acetate
  - saturated solution preparation 95
  - solution preparation 142 185
- Sodium aluminum silicate 140
- Sodium anthraquinone  $\beta$  sulfonate 267
- Sodium benzoate 142
- Sodium benzoate solution preservative 271
- Sodium bicarbonate solution preparation 19
- Sodium  $\beta$  glycerophosphate solution preparation 222
- Sodium  $\beta$  naphthoquinone 4 sulfonate solution preparation 142
- Sodium carbonate solution preparation 14 145
- Sodium chloride standard solution preparation 106
- Sodium citrate solution preparation, 135
- Sodium cyanide solution preparation 143 232
- Sodium diethylbarbiturate buffer preparation 224
- Sodium (mono) diethylbarbiturate preparation 2
- Sodium diethyldithiocarbamate solution preparation 101
- Sodium ethylate solution preparation 163
- Sodium fluorescein reference standard solution preparation 191
- Sodium hydrosulfate solution for absorption of oxygen 191 193 267
- Sodium hydrosulfite (Lycopon) 217
- Sodium hydrosulfite borate solution preparation 175 239
- Sodium hydroxide
  - 5 N solution preparation 175
  - 10 N solution preparation, 181
  - an absorber of carbon dioxide 270 2, 2
  - standard solution preparation 77 91
- Sodium meta bisulfite solution preparation 159 165
- Sodium nitrite solution preparation 413
- Sodium oxalate standard solution preparation 73 93
- Sodium phosphate monobasic 143
- Sodium pyrophosphate solution preparation 101
- Sodium sulfate solution preparation 91 1 6
- Sodium test miscellaneous references 366
- Sodium thiosulfate
  - solution preparation 142
  - standard solution preparation 79 156 163
- Sodium tungstate solution preparation 140 143
- Sodium uranyl zinc acetate (triple sodium salt) 87
- recovery of 91

- Sokalichuk A 356
- Solomon J D 259
- Soluble nitrogen 135
- Solution definition of
  - molar 76
  - normal 76
  - normal oxidizing 76
  - normal reducing 76
- Solvents anhydrous
  - alcohol absolute 171
  - alcohol methyl, 1:1
  - ethylene chloride 171
- Somogyi M 156 161 236
- Somogyi's phosphate tartrate copper iodide solution 156
- Sommerey's G 150
- Soodak M 152
- Sørensen's glycine HCl mixtures for pH 11 35 238
- Sørensen's glycine NaOH mixtures for pH 8.5 12.5 233
- Sørensen's phosphate buffer pH 7.8 8.2 236
- Sørensen's phosphate buffer solutions at 0 C and 37 C pH 5.8 8.2 38
- Sortwell M 166
- Soyenkov D 117
- Span ether reagent 1 4
- Spectral transmittance curves 29 30
- Spectral transmittance curves for copper 103
- Spectrometer miscellaneous references 60
- infrared 80
- mass 80
- Spectrophotometer Beckman Quartz Model DU 38 50
  - absorption cells holders 38 40
  - calibration 43 50
  - of cells 46 47
  - solution concentration 47
  - wavelength 49
  - electrical construction 38
  - filter slide 43
  - focusing
    - hydrogen lamp 48 49
    - tungsten lamp 47 48
  - general construction 38 42
  - light sources 38 43
  - monochromator, 38
  - normal band width 46
  - operation 50
  - optical systems 33
  - photometer calibration 49 50
  - power requirements 43
  - sensitivity control 44
  - slits 38
  - spectral range 42
- Spectrophotometer Coleman Jr, Model 6 33 38
  - general principles 33 36
  - monochromator 33
  - optimum transmittance range 36
  - photometer 34 36
  - spectrophotometer 33
  - instructions for routine use 38
  - setting up and calibration of instrument 36 38
  - adjusting the mechanical galvanometer zero 36 37

- Spectrophotometer setting up and adjusting the mechanical galvanometer—Cont'd  
 calibration of the instrument 37 38  
 mounting and connecting 36  
 verifying the wavelength calibration 37
- Spectrophotometer definition 33
- Spectrophotometry general principles 7 33  
 color and wave length 38 39  
 concentration transmittance calibration graph 3 33
- Lambert Beer law 31 32  
 optical density 33  
 selecting the proper reference solution 8  
 spectral transmittance curves 29 30  
 transmittance 7  
 transmittance 7 23
- Spencer W H 8
- Sperber I 118
- Sperry W M 171
- Spes J F 153
- Spes T D 212
- Spirometer for collection of expired air 352
- Sprague A L 33
- Stabilizing agent  
 acetic acid solution for urine calcium 92  
 metaphosphoric acid 2  
 oxalic acid 2  
 stannous chloride solution 5
- Stadler W C 119 311
- Stafford D D 133
- Standard solution definitions  
 molar 76  
 normal 76  
 normal oxidizing 78  
 normal reducing 76
- Standard solutions preparation of  
 adenine guanine and uracil (microbiology) 245  
 alcohol ethyl 161  
 amino acetic acid 14  
 ammonium sulfate 120 135  
 androsterone 217  
 ascorbic acid 196  
   for field use 91  
 benzoic acid 78 82  
 beta carotene 129  
 bilirubin 413  
 broth 245  
 calcium 94  
 calcium pantothenate 22  
 carotene 18  
 cholesterol 170 171  
 copper 161  
 copper sulfate 1 71 8  
 creatinine zinc chloride crystalline 138  
 L cystine 245  
 glycogen 162  
 glucose 155  
 hydrochloric acid 77  
 inositol 52  
 iodine 79 1 9  
 iron 99  
 lithium lactate 1 9  
 magnesium 97  
 N-methylglucosaminide 189  
 nicotinic acid 17 18  
 oxalic acid 77
- Standard solutions preparation of—Cont'd  
 pantothenic acid 250  
 phenol red 409  
 phenol ulfonphthalate 410 411  
 phosphate (KH<sub>2</sub>PO<sub>4</sub>) 110 111  
 potassium dichromate 9 164 168  
 potassium hydrogen phthalate 77  
 potassium iodate 9  
 potassium permanganate 9 93  
 potassium thiocyanate 78 106  
 pteroylglutamic acid 245  
 pyridoxine hydrochloride 22  
 quinacrine sulfate 180  
 riboflavin 190 193 145 4  
 silver nitrate 78 106  
 sodium chloride 106  
 sodium fluorescein 191  
 sodium hydroxide 7 1  
 sodium oxalate 78 93  
 sodium thiosulfate 79 156 168  
 sulfuric acid 78  
 thiamine 185  
 thiocyanate 364  
 total base 85  
 tryptophane 255  
 L tryptophane 45  
 uric acid 143 145  
 uricase 415  
 vitamin A alcohol 179 18
- Standardization example of 107
- Standardization of solutions  
 benzoic acid 78  
 hydrochloric acid 78  
 iodine 78  
 oxalic acid 8  
 potassium permanganate 79  
 potassium thiocyanate 78 106 107  
 silver nitrate 78  
 sodium thiosulfate 79  
 sulfuric acid 78
- Stannous chloride stabilizing reagent 5
- Stannous chloride solution preparation 199
- Starch indicator 156
- Starch suspension solution 168
- Statistical methods 1
- Statistical methods references to standard works 4 4 3
- Staub H 80
- Stauffer J F 80
- Stead E A Jr 354
- Steam hydrolyzer 115
- Stearic acid 169
- Steele H H 417
- Steele J M 418
- Steele R L
- Steers E 163
- Stehl F L 350
- Stein W H 149
- Step test of physical fitness 373 374
- Stephens I M 21
- Stepka W 149
- Sterilizer Arnold 714
- Stern, H. O 80
- Sternberger L A 40
- Steroids miscellaneous references 19 9
- Sterols microbiological miscellaneous references 60

## Tables—Cont d

- Haldane burette (sample calibration of) 315
- Indicators selected by Clark and Lubs (1917) supplemental by Cohn (1927) 35
- Mercury (density of) at various temperatures 314
- Michelson method for estimating cholinesterase activity corrected values for pH decrement per hour for pH 6.0 to pH 7.99 5
- Muscle protoplasm (factor used in deriving certain components of) from nitrogen balance 330
- Nutrients per pound of food groups as purchased 343 377
- Oxygen factors (Van Slyke) 69
- Partial pressure of water and temperature factor for equilibrated blood at different room temperatures 305
- Riboflavin in urine simplified calculation 190 191
- Serum sodium simplified calculation 88
- Sodium in urine simplified calculation 89
- Sørensen's glycolic HCl mixtures for pH 11.35 938
- Sørensen's glycine NaOH mixtures for pH 8.5 125 38
- Sørensen's phosphate buffer solutions at 0°C and 37°C for pH range 5.8 to 8.2 240
- Sørensen's phosphate mixture for pH 5.8 6 936
- Step test of physical fitness (scoring the) 374
- Sulfur in urine simplified calculation 114
- Tolerances established by the U.S. Bureau of standards for volumetric apparatus 67
- Total base simplified calculation 80
- Walpole's acetate buffer mixtures for pH 3.65 8 937
- Tabone J. 151
- Tagres H. J. 40
- Takadiastase solution preparation 45
- Talbot N. B. 18 19 0 4 0
- Talbot J. H. 83 91 15 3
- Tallen H. H. 7
- Taquini A. C. 303
- Taurer A. 11 91
- Tavinsky H. H. 166
- Taylor C. B. 81
- Taylor F. H. L. 401 40 404
- Taylor H. L. 1 7
- Taylor W. A. 116
- Technique for absorption of carbon monoxide in glass cup (Poughton & Scholander) 78
- Technique for absorption of excess carbon dioxide 78
- Technique for obtaining arterial blood 81 83
- Goldsmith and Light method 63
- Lilenthal and Riley method 67
- Mook's technique 63
- Scholander's technique 63 65
- usual method 61 6
- Techniques for metabolic wards 39 39a
- derived data theoretical balances 393 394
- derived data theoretical weights 394 39a
- general considerations 39 393
- Tee, A. E. 13
- Temperature conversion chart 381 38
- Temperature conversion table 38a
- Temperature equilibration of gas bubble in capillary 78
- Tepley L. J. 71 44 49
- Tesluk H. 151
- Thies R. 4 1
- Theoretical balance metabolic state 393 394
- calcium 393
- computation 393
- nitrogen 393
- phosphorus 393
- potassium 393 394
- sulfur 393 394
- Theoretical weight calculation 394 39a
- Thermal insulation of clothes Clo unit 3 9
- Thermal significance of clothing 380 calculation 380
- Thermobarometer Scholander micro gas analyzer 3 0
- Thermocouple measurements for skin temperature 3 0
- Thiamine
- collection of samples in
- blood 187
- feces 187
- food 187
- urine 18a
- determination in
- blood 186 189
- feces 186 189
- food 186 189
- urine 184 186
- urine field method 0 11
- enzymes for 187
- hydrolysis of bound thiamine 188
- maximum stability range 3
- precipitation of protein 188
- adsorption 188
- elution 189
- extraction 188
- oxidation 188
- preservative for (H<sub>2</sub>O) 18a
- standard curve preparation 186
- standard solution preparation 183
- Thiamine miscellaneous references 913
- Thiochrome
- conversion 184
- formation 54
- Thiocyanate
- preparation for reaction 364
- standard solution preparation 364
- Thiocyanate concentration calculation 36
- Thiocyanate space
- calculation 365 366
- extracellular fluid 364 366
- Thioglycolic acid solution preparation 99
- Thiourea 199
- solution preparation 196 198
- Thompson H. 0

Tungstate dilute solution preparation 154  
 Turner E H 418  
 Turrell E 378  
 Tuttle E 87 118  
 Tyner E H 116  
 Tzu Cheng C 118

## U

Udenfr end S 148 149  
 Ultracentrifuge 81  
 Ultrafiltration 81  
 Ultramicro methods of Dessey and Lowry  
 miscellaneous references 368  
 Ultramicroanalysis 81  
 Umbreit W W 80  
 Unger P N 418  
 Universal indicator and color change 34  
 Uranyl zinc acetate solution preparation  
 87  
 Urbach K F 151  
 Urea clearance miscellaneous references 411  
 Urea cyanide solution preparation 145  
 Urea manometer method miscellaneous  
 references 300  
 Urea nitrogen determination in  
 blood 140 141  
 serum 140 141  
 urine 141  
 Urea solution preparation 277 \*86  
 Urea solution preparation 140  
 Urea solution  
 determination in  
 serum 142  
 urine 144  
 standard solution preparation 143 145  
 Uricase  
 method for uric acid in urine 145 147  
 preparation 145 146  
 Uricolytic activity measurement 146  
 Urine collection and storage in the field 22  
 24  
 determination of  
 acetone bodies total 176 177  
 alcohol 163 164  
 amino acid nitrogen 14  
 ammonia free 141  
 ascorbic acid field method 81 \*04  
 ascorbic acid macro method 197  
 biotin 49 50  
 calcium 95 96  
 chloride, 108 109  
 creatine 139 140  
 creatinine 136 139  
 folic acid 44 47  
 17 ketosteroids \*15  
 succinic acid microbiological 45  
 N methylglutamine 169  
 N methylglutamine field method \*07  
 11  
 pantothenic acid 1  
 pH miscellaneous references \*43  
 phenolsulfonphthalein test 411  
 phosphorus inorganic 110  
 potassium 83  
 protein, 133 134

Urine determination of—Cont'd  
 pyridoxine 23  
 riboflavin 19 19  
 riboflavin field method \*04 \*07  
 sod am 83 90  
 sulfates  
 ethereal sulfates 11  
 inorganic sulfate 114 115  
 total sulfate, 112 114  
 thiamine 184 186  
 thiamine field method 07 211  
 titratable acidity 109 410  
 total base 83 87  
 total nitrogen 1 0 1 3  
 tryptophane \*3 236  
 urea nitrogen 141  
 urobilinogen 414 416  
 uric acid 144  
 Urine normal values 434 435

Urobilinogen  
 determination in  
 feces 416 417  
 urine 416 416  
 standard solution preparation 415  
 Urobilinogen miscellaneous references 418

## V

Vallee B L 119 40  
 Values normal  
 feces 436 437  
 plasma (or serum) 437 433  
 plasma specific gravity  
 men 131  
 women 131  
 urine 434 435  
 whole blood 430 431  
 whole blood specific gravity  
 men 131  
 women 131  
 Van Bruggen J T 19  
 Van Cappel C 316  
 Van de Kamer J H 178  
 Vandam L D 356  
 Vandebelt J M 213  
 Van Klooster H 8 118  
 Van Rysselberghe P J 118  
 Van Slyke D D 1 7 153 175, \*43 83  
 87 0 299 00 311 411  
 Van Slyke K K \*43  
 Van Slyke Neill manometric apparatus 83  
 \*67  
 air free solen on preparation \*6  
 illustration 86  
 machine corrections 26  
 moisture in the manometer 265  
 precautions 67  
 principle \*65  
 Van Soest H 4 1  
 Van Winkle Q 86  
 Vanzant F R 4 0  
 Vapor tension absolute conversion chart,  
 355 386  
 Varney R F \*19  
 Venning E H 216  
 Venous pressure miscellaneous references 357  
 Vials heparinized preparation 1 7

Virtanen A J 148  
 Vischer E 151  
 Visual fluorometry 57 58  
 Vitamin A alcohol standard solution preparation on 179 182  
 Vitamin A and carotene miscellaneous references 214  
 Vitamin A and carotenoids (Carr Price) 1 9 181  
     determination on  
     blood 179 180  
     feces 181  
     serum 1 9 180  
 Vitamin A in serum (dichloroethane) 181 184  
 Vitamin B 184 186 190 191  
 Vitamin B<sub>2</sub> 190 195 194 197  
 Vitamin B  
     pyridoxal 53  
     pyridoxine 193  
 Vitamin B<sub>6</sub> microbiological miscellaneous references 160  
 Vitamin C 190 199 191 194  
 Vitamin D miscellaneous references 14  
 Vitamin solution for intravenous injection 407 408  
 Vitamin tolerance tests  
     determination of  
     intravenous 407 408  
     oral 405 407  
 Vitamins and nutrition miscellaneous references 190 194  
     ascorbic acid 210  
     carotene 210  
     choline 212  
     folate 212  
     nicotinic acid 212  
     pyridine of N-methylpicolamide 213  
     pyridoxine 213  
     pyridoxine 193  
     riboflavin 213  
     thiamine 213  
     tocopherol 213  
     tryptophan 214  
     vitamin A and carotene 194  
     vitamin D 214  
 Volterra M 401  
 Volume of urine per hour from total volume and time of collection 13  
 Volumetric analysis 68 9  
     apparatus 68 4  
     burettes 68 4  
     pipettes 68  
     syringe pipettes 68  
     calibration of apparatus 74 76  
     definition 68  
     standard solutions preparation 6 9  
 Volustat Fisher 244  
 Vorhaus L J II 418  
 Wye tubes 91

# W

Wadman H 150  
 Wadsworth A B 396  
 Waelsch H 160  
 Wagner R 166  
 Watt R B 218

Wakefield E G 118  
 Waldvogel M 197  
 Walker P 81  
 Walpole's acetate buffer mixtures pH 3.65 3.87  
 Walshe V 18  
 Wang Y L 14  
 Ware A G 40  
 Waring blender 19  
 Warkentin J 4 1  
 Warmoltz N 421  
 Warren G H 208  
 Warren J V 306  
 Warringa M G P J 19  
 Warburton H W 80  
 Water  
     glass distilled  
     preparation on 8  
     reference to method, 119  
     salt tests of adrenal sufficiency miscellaneous references 420  
 Water vapor  
     absorption with ethylene glycol 165  
     tension 261  
 Watson C J 233 415 418 4 1  
 Wave length calibration 19  
 Weibelbaum T E 104  
 Weights and measures conversions between different systems 405 407  
 Weiland H 116  
 Weinbach A P 118  
 Weisger J E 43  
 Welcker M L 401  
 Wells B B 17  
 Wells S L 390  
 Werkheiser W C 133  
 Werko L 356  
 Weeselman H J 18  
 Westall R G 150  
 Westerveld W 116  
 Wexler J 357  
 Weyers H A 178  
 Whipple G H 117  
 Whyte L E H 396  
 White A 15  
 White E C 410  
 White F D 208  
 White H L 411  
 White J U 80  
 White J W Jr 151  
 White P D 357  
 Whitehead J K 150  
 Whitehead T 117  
 Whitehorn J C 116 118  
 Whittenberger J L 35  
 Whole blood  
     carbon dioxide determination of 2 0 27  
     carbon monoxide determination of 175 76  
     carbon monoxide (SR) 185 291  
     blood containing small amounts, 89 191  
     combined carbon monoxide and oxygen method 188 189  
     general method 185 188  
     collection and storage in the field 19  
     combined oxygen and carbon dioxide method 270 273

- Whole blood--Cont d  
 hemoglobin copper sulfate method 1 133  
 normal values 430 431  
 oxygen determination of 87 7 8  
 oxygen total arterial 303 306
- Wiggerink G L, 4 1
- Wilder R M 4<sup>00</sup>
- Wiley H F 80
- Wilhelmi A E 15<sup>0</sup>
- Willard H H 116
- Willebrands A F Jr 118
- Williams R H 218
- Williams R J 148 259
- Williams R P 12
- Williams W L 2<sup>01</sup>
- Wilson H 2 0
- Wilson T M 418
- Wundchill  
 calculation 389  
 definition 389
- Wind velocity conversion chart 381 383
- Winegard H M 149
- Winkler A W 401
- Winkler's solution for absorption of carbon  
 monoxide 286
- Winsten W A 148 151
- Winternitz, W 214
- Wintrobe M 104 33 394 399 401
- Wintrobe tubes 13
- Wiring diagram of Beckman Quartz spectro  
 photometer 41
- Wiss O 151
- Witter F F 178
- Wixon R L, 19 4 0
- Woisod A J 11 149
- Wolcott G H 166
- Wolf L M 11
- Wolfe J K, 219 4 0
- Wolff H 115
- Wolff W A 4<sup>01</sup>
- Wolfson W Q 1 3 151 1 3
- Wong S Y 151
- Wood filter (ultra violet) 5
- Wood L L 117
- Wood T R 2<sup>60</sup>
- Woodbury L A 219
- Woodman F J 148
- Woodruff H B 151
- Woodson, H W 259
- Woolley J G 255
- Wright L D 247 55 56
- Wu H 136 139 160
- Wuth A 116
- Wyatt P F 115
- Wyckoff R W 149
- Wyckoff, R W G 80
- Wynn V 148
- Y
- Yanthophyll 180
- Y
- Yeast dextrose agar 45
- Yeast dextrose broth 45
- Younts J B 1
- Young A R 34
- Young N F 149 151
- Youngsberg G E, 1 3
- Youngsberg M V 173
- Yule G U 4 3
- Z
- Zaffaroni A 151
- Zechmeister L 148
- Zeolite 140  
 activated 87 208
- Zimmerman W J 116
- Zinc reference to method 119
- Zinc sulfate solut on preparation 156
- Zustaga G 1 0
- Zuckerbrod M, 419
- Zuckerman J L 178
- Zymaris M C 178



